

Alfred Butler

Allyson

2000-2005

DURHAM UNIVERSITY ARCHIVES	
LOCATOR	5169
BOX	11
PAGE	3

SERVING THE SCIENTIFIC COMMUNITY FOR OVER 30 YEARS

Aldrichimica ACTA

VOL. 33, NO. 1 • 2000



The Chemical
Vapor Deposition
of Metal Nitride Films
Using Modern
Metalorganic
Precursors

Synthetic
Applications of
Indium Trichloride
Catalyzed Reactions


Purpald[®]: A Reagent
That Turns Aldehydes
Purple!



 **ALDRICH**[®]
CHEMISTS HELPING CHEMISTS




New Products

The *tert*-butyldimethylsilyl analog of Danishefsky's diene is now available. It has been used to prepare a variety of bicyclic enones and 2,3-dihydro-4*H*-pyran-4-ones.¹⁻³



(1) Uchida, H. et al. *Tetrahedron Lett.* **1999**, 40, 113. (2) Pudukulathan, Z. et al. *J. Am. Chem. Soc.* **1998**, 120, 11953. (3) Annunziata, R. et al. *J. Org. Chem.* **1992**, 57, 3605.

51,536-1 *trans*-3-(*tert*-Butyldimethylsilyloxy)-1-methoxy-1,3-butadiene, 95%
5g \$53.00



A variety of nucleophiles have been used to open the oxirane ring of these compounds. Examples include lithium acetylides^{1,2} and lithium dithianes.³


(1) Arista, L. et al. *Heterocycles* **1998**, 48, 1325. (2) Maguire, R. J. et al. *J. Chem. Soc., Perkin Trans. 1* **1998**, 2853. (3) Smith, A. B., III; Boldi, A. M. *J. Am. Chem. Soc.* **1997**, 119, 6925.

52,599-5 *tert*-Butyldimethylsilyl (*R*)-(+)-glycidyl ether, 98%
5g \$65.00

52,490-5 *tert*-Butyldimethylsilyl (*S*)-(-)-glycidyl ether, 98%
5g \$65.00

52,613-4 *tert*-Butyldimethylsilyl glycidyl ether, 98%
10g \$42.00

Building blocks for *trans*-2-methyl-6-substituted piperidines via deprotonation with *sec*-butyllithium followed by reaction with an electrophile.¹⁻³



52,288-0



52,290-2

(1) Chackalamannil, S. et al. *J. Am. Chem. Soc.* **1996**, 118, 9812. (2) Beak, P.; Lee, W. K. *J. Org. Chem.* **1993**, 58, 1109. (3) *Idem ibid.* **1990**, 55, 2578.

52,288-0 (*S*)-(+)-*N*-(*tert*-Butoxycarbonyl)-2-methylpiperidine, 98%
1g \$31.50

52,290-2 (*S*)-(+)-2-Methylpiperidine, 97%
1g \$68.00


Precursors for a wide variety of substituted indoles.¹⁻³



(1) Battaglia, S. et al. *Eur. J. Med. Chem.* **1999**, 34, 93. (2) Joseph, B. et al. *J. Heterocycl. Chem.* **1997**, 34, 525. (3) Le Borgne, M. et al. *Bioorg. Med. Chem. Lett.* **1999**, 9, 333.


51,874-3 5-Bromoindole-3-carboxaldehyde, 98%
5g \$26.40

51,112-9 1-Methylindole-2-carboxaldehyde, 97%
5g \$45.70



51,520-5 3-Indoleglyoxylyl chloride, 98%
5g \$40.60


α,β -Unsaturated *N*-methoxy-*N*-methylamides are prepared from this reagent through deprotonation with sodium hydride, reaction with an alkyl halide, and in situ heating.



Benev, C. et al. *Tetrahedron Lett.* **1998**, 39, 5779.

51,139-0 *N*-Methoxy-*N*-methyl-2-(phenylsulfinyl)acetamide, 96%
1g \$30.00; 5g \$100.00


2-Substituted picolines have been prepared from this compound. Examples include (2-pyridyl)indoles and endothelin receptors.^{1,2}



(1) Amat, M. et al. *J. Org. Chem.* **1997**, 62, 3158. (2) Kourounakis, A. et al. *Bioorg. Med. Chem. Lett.* **1997**, 7, 2223.

51,894-8 2-Chloro-3-methylpyridine, 97%
5mL \$35.40; 25mL \$117.90

This carbonate resin is used to bind amines or amino acids as urethanes.



Dipeptides and hydantoins have been prepared from these polymer-bound urethanes.¹⁻³


(1) Dixit, D. M.; Leznoff, C. C. *J. Chem. Soc., Chem. Commun.* **1977**, 798. (2) Dressman, B. A. et al. *Tetrahedron Lett.* **1996**, 37, 937. (3) Gouilleux, L. et al. *ibid.* **1996**, 37, 7031.

51,529-9 4-Nitrophenyl carbonate, polymer-bound on Wang Resin
1g \$35.60; 10g \$199.00


Starting material for the preparation of 4-substituted imidazoles.^{1,2}

(1) Lange, J. H. M. et al. *Tetrahedron* **1995**, 51, 13447. (2) Singh, B. et al. *J. Med. Chem.* **1992**, 35, 4858.

47,869-5 4-Bromo-1*H*-imidazole, 97%
1g \$31.50



Naphthoquinones are prepared from this compound via palladium-catalyzed coupling reactions with tributylstannyl-heteroaromatics or by nucleophilic displacement of one or both bromides.^{1,2}




(1) Yoshida, S. et al. *Chem. Lett.* **1996**, 139. (2) Falling, S. N.; Rapoport, H. *J. Org. Chem.* **1980**, 45, 1260.

52,342-9 2,3-Dibromo-1,4-naphthoquinone, 97%
25g \$58.00

Precursor for 3-substituted-2-methyl-2-cyclopenten-1-ones.^{1,2}

(1) Cossy, J. et al. *Tetrahedron Lett.* **1997**, 38, 4069. (2) Junga, H.; Bleichert, S. *ibid.* **1993**, 34, 3731.

51,440-3 3-Ethoxy-2-methyl-2-cyclopenten-1-one, 97%
5g \$36.00



Aldrichimica ACTA

VOL. 33, NO. 1 • 2000

Aldrich Chemical Co., Inc.

1001 West Saint Paul Ave.
Milwaukee, WI 53233 USA

To Place Orders

Telephone 800-558-9160 (USA)
or 414-273-3850
FAX 800-962-9591 (USA)
or 414-273-4979
Mail P.O. Box 2060
Milwaukee, WI 53201 USA

General Correspondence

William J. Seitz, Sharbil J. Firsan, or Jennifer L. Botic
P.O. Box 355, Milwaukee, WI 53201 USA

Customer & Technical Services

Customer Inquiries 800-558-9160
Technical Service 800-231-8327
Sigma-Aldrich Fine Chemicals 800-336-9719
Custom Synthesis 800-336-9719
Flavors & Fragrances 800-227-4563
International 414-273-3850
24-Hour Emergency 414-273-3850
Web Site <http://www.sigma-aldrich.com>
E-Mail aldrich@sial.com

To request your **FREE** subscription to the
Aldrichimica Acta,

please call: **800-558-9160** (USA)

or write: **Attn: Mailroom**
Aldrich Chemical Co., Inc.
P.O. Box 355
Milwaukee, WI 53201-9358

International customers, please contact your
local Sigma-Aldrich office.

The *Aldrichimica Acta* is also available on the
Internet at <http://www.sigma-aldrich.com>.

All prices listed in this publication are subject to change
without notice.

Aldrich brand products are sold through Sigma-Aldrich, Inc.
Sigma-Aldrich, Inc. warrants that its products conform to the
information contained in this and other Sigma-Aldrich
publications. Purchaser must determine the suitability of the
product for its particular use. See reverse side of invoice or
packing slip for additional terms and conditions of sale.

For worldwide mailing/contact information, please see
the inside back cover.

Aldrichimica Acta is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich family. © 2000 by Sigma-Aldrich Co.
Printed in the United States.

About Our Cover

Our first cover of 2000 is a montage of various art works that have graced some of our past publications. The top three and number 6 are part of the extensive personal art collection of Dr. Alfred Bader, founder of Aldrich Chemical Co. and formerly the *Acta's* "Chemist Collector". These are *The Farm Scale* by John Whalley [no. 1], *King David* by Govaert Flinck [no. 2], *Trompe L'oeil* [no. 3], and *Lady in Black* [no. 6].

The other two are *The Dreamer* by Pierre Auguste Renoir [no. 4] and *The Column* by Hubert Robert [no. 5]. Both of these paintings are in the permanent collection of the Saint Louis Art Museum, Saint Louis, MO.

A reproduction (11x14 in.) of each of these six paintings is available for only a nominal \$4.00 shipping-and-handling charge. Request your own copy by calling **800-558-9160** (USA) or your local Sigma-Aldrich office, e-mailing us at aldrich@sial.com, or visiting our Web site at www.sigma-aldrich.com.



"Please Bother Us."

Clint Lane, President

Dear Fellow Chemists:

It is a great pleasure for me to again provide a contribution to our *Aldrichimica Acta*, this time as President of Aldrich. During the years 1973–1977, I contributed nine reviews to this publication. Now, 23 years later, I am back on these pages. Much has changed over the years, but we never stopped publishing what we here at Aldrich simply call "the ACTA".

Year 2000 was viewed as a perfect time to roll out a new look for our classic. However, we have not deleted any of the features that you have come to expect. Art on the cover (reprints available), excellent and timely reviews, lab notes, featured new products, etc. are all here for you to enjoy. Our goal is to continue to be *Chemists Helping Chemists* as we move into the new millennium. We know you will continue to find new ideas and information to help you in your research efforts.

Thank you for staying with us over all these years. This is the first issue in the year 2000. We are up to Volume 33 (33 years) and still counting. As always, "Please Bother Us" with comments or contributions.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page or on the inside back cover.




Lab Notes

Aldrich Rotary Evaporator Antisplash Adapters as Solvent Traps in Recrystallizations

For practicing synthetic chemists, recrystallization of reaction intermediates is routine. For maximal product recovery, it is important to prepare a saturated solution of the compound. In many instances, this is difficult to judge, since excess solvent is sometimes needed to completely dissolve the solute, or since filtration subsequent to charcoal treatment is followed by washing of the residue with hot solvent. These steps result in dilution of the crystallizing solution, and, therefore, concentration of such a solution has to be performed. This is usually accomplished by simply heating the solution to evaporate some of the solvent. However, this leads to two problems: (a) how much solvent is actually left in the flask cannot be accurately estimated, and (b) solvent vapors are discharged into the hood and are not collected for recycling or disposal. We have devised a simple solution to these two problems that uses the Aldrich rotary evaporator antisplash adapter (without return holes, Cat. No. Z17,604-4 to Z20,329-7).

The rotary evaporator antisplash adapter doubles as a solvent trap for the recovery of the crystallization solvent. The compound is placed in

a round-bottom flask, and heated while the solvent is gradually added until the compound



dissolves completely. The antisplash adapter is then attached to the flask and heating is continued. The solvent that evaporates at this

stage condenses and collects in the antisplash adapter. A couple of modifications of this setup are possible. For larger-scale operations, the solutions can be magnetically stirred and heated. Also, for better recoveries of volatile solvents, a reflux condenser can be attached to the top of the antisplash adapter. After the solvent has collected in the adapter, the adapter is disconnected, and the solvent is removed and its volume measured. Besides effecting the desired concentration of the crystallization solution, this procedure allows one to calculate the amount of solvent left in the crystallization flask, and to recycle or properly dispose of the condensed solvent.

We routinely use this setup and find it very convenient. We hope that other researchers will find it equally helpful.

Mahesh K. Lakshman, Assistant Professor
Department of Chemistry
University of North Dakota
Grand Forks, ND 58202-9024


E-mail: mlakshman@mail.chem.und.nodak.edu

Please turn to page 12 for more Lab Notes.

Aldrich splash-guard adapters


Antisplash, without return holes

Top \$ joint	Bottom \$ joint	Cap. (mL)	Cat. No.	Each
24/40	24/40	100	Z17,604-4	\$49.00
24/40	24/40	250	Z14,779-6	55.00
24/40	24/40	500	Z14,781-8	59.95
24/40	14/20	250	Z14,782-6	54.50
29/32	29/32	100	Z20,344-0	49.00
29/32	29/32	250	Z20,327-0	55.00
29/32	29/32	500	Z20,328-9	59.95
29/32	14.5/23	250	Z20,329-7	54.50



*For custom configurations, e-mail the
Aldrich Glass Shop at aldglass@sial.com.*

Do you have an innovative shortcut or unique laboratory hint you'd like to share with your fellow chemists? If so, please send it to Aldrich (attn: Lab Notes, Aldrichimica Acta). For submitting your idea, you will receive a complimentary, laminated periodic table poster (Cat. No. Z15,000-2). If we publish your lab note, you will also receive an Aldrich periodic table turbo mouse pad (Cat. No. Z24,409-0). It is Teflon®-coated, 8½ x 11 in., with a full-color periodic table on the front. We reserve the right to retain all entries for future consideration.



Teflon is a registered trademark of E.I. du Pont de Nemours & Co., Inc.

Sigma-Aldrich Worldwide Locations

Argentina

Sigma-Aldrich de Argentina, S.A.
Av. Pueyrredon 2446 Piso 5 "B"
C.P. C1119ACU C.A. de Buenos Aires
Argentina
Tel: (54-11) 4807-0321
Fax: (54-11) 4807-0346
E-mail: cservice@sigma-aldrich.com.ar

Australia

Sigma-Aldrich Pty., Limited
P.O. Box 970
Castle Hill, N.S.W. 1765
Free Tel: 1-800-800-097
Free Fax: 1-800-800-096
Tel: 61-2-9841-0555
Fax: 61-2-9841-0500
E-mail: ausmail@sial.com

Austria

Sigma-Aldrich Handels GmbH
Hebbelplatz 7
1100 Wien
Tel: 43 1 605 8110
Fax: 43 1 605 8120
E-mail: sigma@sigma.co.at

Belgium

Sigma-Aldrich nv/sa
K.Cardiijnplein 8
B-2880 Bornem
Free Tel: 0800-147.47
Free Fax: 0800-147.45
Tel: 03 899.13.01
Fax: 03 899.13.11
E-mail: becustsv@eurnotes.sial.com

Brazil

Sigma-Aldrich Quimica Brasil Ltda.
Rua Sabará, 566- Cj. 53
01239-010 São Paulo, SP
Tel: 55-11-231-1866
Fax: 55-11-257-9079
E-mail: sigmabr@ibm.net

Canada

Sigma-Aldrich Canada Ltd.
2149 Winston Park Drive
Oakville, Ontario L6H 6J8
Free Tel: 800-565-1400
Free Fax: 800-265-3858
Tel: 905-829-9500
Fax: 905-829-9292
E-mail: canada@sial.com

Czech Republic

Sigma-Aldrich s.r.o.
Pobřežní 46
186 00 Praha 8
Tel: 00-420-2-2176 1300
Fax: 00-420-2-2176 3300
E-mail: CZECustSV@eurnotes.sial.com

Denmark

Sigma-Aldrich Denmark A/S
Vejlegaardsvej 65 B
2665 Vallensbaek Strand
Tel: +45 43 56 59 10
Fax: +45 43 56 59 05
E-mail: denorder@eurnotes.sial.com

Eire

Sigma-Aldrich Ireland Ltd.
Airton Road
Tallaght
Dublin 24
Free Tel: 1 800 200 888
Free Fax: 1 800 600 222
Tel: 01 404 1900
Fax: 01 404 1910
E-mail: EICustsv@eurnotes.sial.com

Finland

Sigma-Aldrich Finland
YA-Kemia Oy
Teerisuonkuja 4
00700 Helsinki
Tel: (09) 350-9250
Fax: (09) 350-92555
E-mail: finorder@eurnotes.sial.com

France

Sigma-Aldrich Chimie S.a.r.l
L'Isle D'Abeau Chesnes
B.P. 701
38297 St. Quentin Fallavier
Cedex
Tel Numero Vert: 0800 21 14 08
Fax Numero Vert: 0800 03 10 52
E-mail: fradsv@eurnotes.sial.com

Germany

(also SE Europe, the Baltics, Africa and the Middle East)

Sigma-Aldrich Chemie GmbH
Eschenstr. 5
82024 Taufkirchen
Freecall Tel: 0800-5155 000
Freecall Fax: 0800-6490 000
Tel: +49(0)89-6513-0
Fax: +49(0)89-6513-1169
Export Department:
Tel: +49(0)89-6513-1807
Fax: +49(0)89-6513-1877
E-mail: DeOrders@eurnotes.sial.com

Greece

Sigma-Aldrich (o.m.) Ltd.
72 Argonafon Str.
16346 Ilioupoli, Athens
Tel: 30-1-9948010
Fax: 30-1-9943831
E-mail: GRCustSV@eurnotes.sial.com

Hungary

Sigma-Aldrich Kft
1399 Budapest
Pf. 701/400
Magyarország
Tel ingyenes: 06-80-355-355
Fax ingyenes: 06-80-344-344
Tel: 06-1-235-9055
Fax: 06-1-235-9050
E-mail: info@sigma.sial.hu

India

Sigma-Aldrich Corporation
Bangalore location:
Survey No. 31/1, Sitharamapalaya
Mahadevapura P.O.,
Bangalore 560048
Phone: 91 80 852 4222
Fax: 91 80 852 4214
New Delhi location:
Flat No. 4082, Sector B 5/6
Vasant Kunj, New Delhi 110070
Phone: (011) 689 9826
Fax: (011) 689 9827
E-mail: sigma@del2.vsnl.net.in

Israel

Sigma-Aldrich Israel Ltd.
Park Rabin, Rehovot 76100
Toll Free Tel: 1-800-70-2222
Tel: 972-8-948-4222
Fax: 972-8-948-4200
E-mail: sigisr@sigma.co.il

Italy

Sigma-Aldrich S.r.l.
Via Gallarate 154
20151 Milano
Numero Verde: 800-827018
Tel: 02-33417.310
Fax: 02-38010.737
E-mail: itorder@eurnotes.sial.com

Japan

Sigma-Aldrich Japan K.K.
Higashi Nihonbashi Sky Bldg.
1-1-7 Higashi Nihonbashi, Chuo-ku
Tokyo 103-0012
Tel: 81-3-5821-3111
Fax: 81-3-5821-3170

Korea

Sigma-Aldrich Korea
Samhan Camus Annex, 10th Floor
17-26 Yoido-dong, Yungdeungpo-ku
Seoul, South Korea
Toll Free Tel: 080-023-7111
Toll Free Fax: 080-023-8111
Tel: 82-2-783-5211
Fax: 82-2-783-5011
E-mail: order@sial.co.kr

Malaysia

Sigma-Aldrich (M) Sdn Bhd
60, 60-1 Jalan Awan Jawa
Taman Yarl
Off Jalan Kelang Lama
58200 Kuala Lumpur
Malaysia
Tel: 03-7824181
Fax: 03-7824067
E-mail: sigalm@pojaring.my

Mexico

Sigma-Aldrich Química, S.A. de C.V.
Calle 6 Norte No. 107
Parque Industrial Toluca 2000
50200 Toluca, México
Teléfono sin costo: 01-800-007-5300
Teléfono en Toluca: (72)-76-1600
Fax sin costo: 01-800-712-9920
Fax en Toluca: (72)-76-1601
E-mail: mexico@sial.com

Netherlands

Sigma-Aldrich Chemie BV
Stationsplein 4 E
Postbus 27
NL-3330 AA Zwijndrecht
Tel Gratis: 0800 022 9088
Fax Gratis: 0800 022 9089
Tel: 078-620 54 11
Fax: 078-620 54 21
E-mail: nlcustsv@eurnotes.sial.com

Norway

Sigma-Aldrich Norway AS
Tevlingveien 23
1081 Oslo, Norway
Tel: 47 2317 6000
Fax: 47 2317 6010
E-mail: nororder@sial.com

Poland

Sigma-Aldrich Sp. z o.o.
Szelągowska 30
61-626 Poznań
Tel: (48-61) 823-24-81
Fax: (48-61) 823-27-81
E-mail: plcustsv@eurnotes.sial.com

Portugal

Sigma-Aldrich Quimica S.A.
(Sucursal em Portugal)
Aptdo. 131
2710 Sintra
Free Tel: 0800-202180
Free Fax: 0800-202178
E-mail, encomendas:
poorders@eurnotes.sial.com

Russia

Sigma-Aldrich Russia
TechCare Systems, Inc.
Makarenko Str. 2/21 Bldg. 1 Flat 22
Moscow 103062
Tel: 7-095-975-3321
Fax: 7-095-975-4792
E-mail: techcare@online.ru

Singapore

Sigma-Aldrich Pte., Ltd.
102E Pasir Panjang Road
#08-01, Citilink Warehouse
Singapore, 118529
Tel: (65) 271 1089
Fax: (65) 271 1571
E-mail: sapl@sialsea.com

South Africa

Sigma-Aldrich South Africa (Pty) Ltd.
Southern Life Industrial Park
Unit 16 & 17
CNR Kelly & Ackerman Streets
Jet Park, Boksburg 1459
Tel: (011) 397 8886
Fax: (011) 397 8859
E-mail: sigald@icon.co.za

Spain

Sigma-Aldrich Química S.A.
Aptdo. 161
28100 Alcobendas (Madrid)
Free Tel: 900-101376
Free Fax: 900-102028
Tel: 91-6619977
Fax: 91-6619642
E-mail, pedidos:
esorders@eurnotes.sial.com

Sweden

Sigma-Aldrich Sweden AB
Solkraftsvagen 14C
135 70 Stockholm
Tel (in Sweden only): 020-35 05 10
Fax (in Sweden only): 020-35 25 22
Outside Sweden Tel: 46-8-742 4200
Outside Sweden Fax: 46-8-742 4243
E-mail: sweorder@eurnotes.sial.com

Switzerland

Fluka Chemie GmbH
Postfach 260
CH-9471 Buchs
Tel: 081-755 27 23
Fax: 081-755 28 40
Swiss Free Call: 0800 80 00 80
E-mail: fluka@sial.com

United Kingdom

Sigma-Aldrich Company Ltd.
The Old Brickyard
New Road
Gillingham
Dorset
SP8 4XT
Free Tel: 0800 717181
Free Fax: 0800 378785
Tel: (01747) 822211
Fax: (01747) 823779
E-mail: ukcustsv@eurnotes.sial.com

United States

Aldrich Chemical Company, Inc.
1001 West Saint Paul Avenue
Milwaukee, Wisconsin 53233
Toll Free Tel: 1-800-558-9160
Toll Free Fax: 1-800-962-9591
Tel: 414-273-3850
Fax: 414-273-4979
E-mail: aldrich@sial.com
Sigma Chemical Company
3050 Spruce Street
St. Louis, Missouri 63103
Toll Free Tel: 1-800-521-8956
Toll Free Fax: 1-800-325-5052
Tel: 314-771-5765
Fax: 314-771-5757
E-mail: sigma@sial.com

High-Hazard Chemicals



Such as
Chlorosulfonic
Acid,
Epichlorohydrin,
Allyl Chloride,
and more...

Handled **Safely** and **Conveniently**
18 to 850L Returnable Containers

MINI-BULKTM

MINI-BULK is a trademark of Sigma-Aldrich Co.

For details, call **1-800-213-1206** (USA) or
visit our Web site at **www.sigma-aldrich.com**

ALDRICH CHEMICAL COMPANY, INC.
P.O. BOX 355
MILWAUKEE, WISCONSIN 53201 USA



CHANGE SERVICE REQUESTED


PRESORTED STANDARD
U.S. POSTAGE
PAID
MILWAUKEE, WISCONSIN
PERMIT NO. 552

014

AN EXCLUSIVE REVIEW BY DR. REUBEN D. RIEKE

Aldrichimica ACTA

VOL. 33, NO. 2 • 2000




Catalytic
Synthesis
of Thiocrowns

New
Organometallic
Reagents

Tellurium
Chemistry

 **ALDRICH**[®]
CHEMISTS HELPING CHEMISTS


NEW PRODUCTS



Quinoline building block.¹ Used as an amino component in a novel series of bladder-selective diaminocyclobutenedione potassium channel openers.²

(1) Armer, R.E. et al. *Bioorg. Med. Chem. Lett.* **1999**, 9, 2430. (2) Bulera, J.A. *J. Med. Chem.* **2000**, 43, 1187.


52,809-9 8-Amino-6-methoxyquinoline hydrobromide, 97%
1g \$57.10; 10g \$317.00



Employed in the synthesis of biologically active naphthalenes.¹ The aldehyde moiety has been elaborated into an oxetane in a photocyclization² and to a butanone in the preparation of the NSAID nabumetone.³

(1) Bosca, F. et al. *Photochem. Photobiol.* **2000**, 71, 173. (2) Abe, M. et al. *J. Chem. Soc., Perkin Trans. 1* **1998**, 3261. (3) Prabhakar, C. et al. *Org. Process Res. Dev.* **1999**, 3, 121.


53,308-4 6-Methoxy-2-naphthaldehyde, 98%
5g \$138.00



Used recently in the synthesis of oligothiophenes^{1,2} and oligothiophenylferrocene complexes.³

(1) Jiang, B. et al. *J. Am. Chem. Soc.* **1999**, 121, 9744. (2) Ibrahim, M. A. et al. *Synth. Met.* **1999**, 105, 35. (3) Zhu, Y.; Wolf, M. O. *Chem. Mater.* **1999**, 11, 2995.


52,285-6 5-Bromo-2,2'-bithiophene, 96%
1g \$40.90; 5g \$171.20



Direct copper catalyzed substitution of the bromines with nitrogen led to chiral bipyridine ligands based on camphor sulfam or a prolinol ether.¹ The bis(3-pyrazolyl) derivative was studied as a tetradentate chelating ligand.²

(1) Kandzia, C. et al. *Tetrahedron: Asymmetry* **1993**, 4, 39. (2) Couchman, S. M. et al. *Polyhedron* **1999**, 18, 2633.


51,388-1 6,6'-Dibromo-2,2'-dipyridyl, 90%
500mg \$45.70



The BPO initiated graft polymerizations of this N-substituted acrylamide and N-isopropylacrylamide (41,532-4) onto ethylene-propylene-diene terpolymer were conducted. Thermal and mechanical properties were studied along with changes in morphology upon irradiation.

Park, K.G. et al. *J. Appl. Polym. Sci.* **1999**, 74, 3259.

53,004-2 N-Phenylacrylamide, 99%
10g \$84.60



Used for preparing the enantiopode of Garner's aldehyde,¹ peptide nucleic acid (PNA) analogs,² and chiral gadolinium complexes that are employed as MRI contrasting agents.³

(1) Avenoza, A. et al. *Synthesis* **1997**, 1146. (2) Ramasamy, K. S.; Seifert, W. *Bioorg. Med. Chem. Lett.* **1996**, 6, 1799. (3) Sajiki, H. et al. *Synth. Commun.* **1996**, 26, 2511.

47,375-8 (R)-(+)-2-Amino-3-benzyloxy-1-propanol, 97% (99% ee/GLC) (1)
1g \$38.60

47,376-6 (R)-(+)-3-Benzyloxy-2-(tert-butoxycarbonyl-amino)-1-propanol, 98% (2)
5g \$169.60

Used to prepare fluoro-substituted pyrrolidinyl-2-pyridinones as antibacterial agents¹ and in the study of Protein Kinase C (PKC) inhibitors.^{2,3}




1, R=CBZ
 2, R=BOC

(1) Li, Q. et al. *Bioorg. Med. Chem. Lett.* **1998**, 8, 1953. (2) Jagdmann, G.E. et al. *ibid.* **1995**, 5, 2015. (3) Lai, Y. et al. *ibid.* **1995**, 5, 2151.

49,412-7 Benzyl 3-pyrroline-1-carboxylate, 90% (1)
1g \$16.50; 5g \$54.90

47,751-6 tert-Butyl 2,5-dihydro-1H-pyrrole-1-carboxylate, 97% (2)
1g \$67.60


Bis-protected guanylating agent used in the preparation of guanidinophthyl amino acid esters, "inverse" substrates for trypsin and trypsin-like enzymes.^{1,2}



(1) Tanizawa, K. et al. *Chem. Pharm. Bull.* **1999**, 47, 104. (2) Sekizaki, H. et al. *Tetrahedron Lett.* **1997**, 38, 1777.

43,942-8 N,N'-Bis(benzyloxycarbonyl)-1H-pyrazole-1-carboxamide, 97%
1g \$24.60; 5g \$81.90


The highly reactive anhydride moiety makes this a commonly used scavenger resin to remove primary and secondary aliphatic amines from solution-phase reactions.



Coppola, G.M. *Tetrahedron Lett.* **1998**, 39, 8233.

51,437-3 Isatoic anhydride, polymer-bound, 2.0-2.5 mmol N/g
5g \$88.40


A pheromone intermediate was prepared by displacement of the bromide with an acetylene followed by selective deprotection of the silyl group with neutral alumina in the presence of a THP ether.



Capdevila, A. et al. *Org. Lett.* **1999**, 1, 845.

51,314-8 (6-Bromohexyloxy)-tert-butyl dimethylsilane, 99%
5mL \$38.40; 25mL \$127.90


Pharma building block used in the preparation of topomerase inhibitors,¹ neuroprotective agents,² and novel antagonists for vascular 5HT_{1B}-like receptors.³



(1) Dudouit, F. et al. *Bioorg. Med. Chem. Lett.* **2000**, 10, 553. (2) Di Fabio, R. et al. *J. Med. Chem.* **1999**, 42, 3486. (3) Moloney, G.P. et al. *ibid.* **1999**, 42, 2504.

52,562-6 4-[(N-BOC)aminomethyl]aniline, 97%
1g \$22.00

BOC protection makes this aniline a good substrate for Suzuki coupling. Two examples are the biaryls from reaction with pyridine¹ or 3-benzaldehyde² boronic acids.



(1) Lamothe, M. et al. *J. Med. Chem.* **1997**, 40, 3542. (2) Coutis, I.G.C. et al. *Tetrahedron Lett.* **1997**, 38, 5563.

52,724-6 N-(tert-Butoxycarbonyl)-4-bromoaniline, 97%
10g \$40.00

Aldrichimica ACTA

VOL. 33, NO. 2 • 2000

Aldrich Chemical Co., Inc.

1001 West Saint Paul Ave.
Milwaukee, WI 53233 USA

To Place Orders

Telephone 800-558-9160 (USA)

or 414-273-3850

FAX 800-962-9591 (USA)

or 414-273-4979

Mail P.O. Box 2060

Milwaukee, WI 53201 USA

General Correspondence

William J. Seitz, Sharbil J. Firsan, or Jennifer L. Botic

P.O. Box 355, Milwaukee, WI 53201 USA

Customer & Technical Services

Customer Inquiries 800-558-9160

Technical Service 800-231-8327

Sigma-Aldrich Fine Chemicals 800-336-9719

Custom Synthesis 800-336-9719

Flavors & Fragrances 800-227-4563

International 414-273-3850

24-Hour Emergency 414-273-3850

Web Site <http://www.sigma-aldrich.com>

E-Mail aldrich@sial.com

Aldrichimica Acta is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich family. © 2000 by Sigma-Aldrich Co. Printed in the United States.

To request your **FREE** subscription to the *Aldrichimica Acta*,

please call: 800-558-9160 (USA)

or write: **Attn: Mailroom**
Aldrich Chemical Co., Inc.
P.O. Box 355
Milwaukee, WI 53201-9358

International customers, please contact your local Sigma-Aldrich office.

The *Aldrichimica Acta* is also available on the Internet at <http://www.sigma-aldrich.com>.

All prices listed in this publication are subject to change without notice


Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale

For worldwide mailing/contact information, please see the inside back cover.

"Please
Bother
Us!"

Clint Lane

Clint Lane, President



Professor Jonathan S. Lindsey (NCSU, Raleigh) kindly suggested that we make 4-(trimethylsilyl)ethynylbenzaldehyde. This protected arylacetylene is a useful precursor to porphyrin building blocks, and, subsequently, porphyrin-based molecular devices.

Ravikanth, M.; Strachan, J.-P.; Li, F.; Lindsey, J.S. *Tetrahedron* 1998, 54, 7721.

52,338-0

**4-(Trimethylsilyl)ethynyl-
benzaldehyde, 97%**

5g \$57.00

Naturally, we made this useful building block. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page or on the inside back cover.

Table of Contents

The Catalytic Synthesis of Thiocrowns from Thietanes and Thiranes by Metal Carbonyl Complexes39
-R. D. Adams

The Preparation of Highly Reactive Metals and the Development of Novel Organometallic Reagents.....52
-R. D. Rieke

Add a Little Tellurium to Your Synthetic Plans!.....66
-J. V. Comasseto &
R. E. Barrientos-Astigarraga


About Our Cover

The City from Greenwich Village (oil on canvas, 26 in. x 33³/₄ in.) was painted in 1922 by the American artist John Sloan, one of a number of artists who came to be called the "Ash Can School" because of their use of busy city streets, tenements, and back alleys as subjects. Sloan, who was born in 1871 in Lock Haven, Pennsylvania, worked during the 1890s as a newspaper illustrator in Philadelphia, where he came in contact with the artist Robert Henri. Henri had returned from several years' study in Paris determined to lead a revolt against academic art and with a strong belief that works of art should reflect the reality of contemporary life.

After 1904, Sloan settled permanently in New York, and though he continued to work as an illustrator, he devoted himself increasingly to the pursuit of subjects of human interest and the depiction of scenes of everyday life among ordinary people in the city. He was one of the original members of a group of artists called The Eight, who were not united by a common style but by a rejection of academic aestheticism. This group came into being in 1907, when the National Academy of Design rejected the work of Sloan, George Luks, and William Glackens, and Henri withdrew his own paintings in protest from the Academy's annual exhibition.

The City from Greenwich Village shows the view from the artist's Washington Place studio out over the street on a rainy winter evening. Automobiles and figures bearing umbrellas make their way under the elevated trains, which ran then along Sixth Avenue. Electric lights reflect from the rainy streets and cast a glow between the buildings, and the lights of the skyscrapers at the far left are like a celestial vision in the distance. A much more romantic image than is commonly found in Sloan's earlier works, it is a celebration of the energy and excitement of a great modern city.

This painting is a gift of Helen Farr Sloan to the National Gallery of Art, Washington, D.C.





Lab Notes

Analysis of Volatile Compounds in Resins by GC/MS

A convenient and simple way to analyze volatile compounds in resins used in polymerization processes is by GC/MS. The power of GC/MS as an analytical tool is well developed and most laboratories are now equipped with bench-top GC/MS units. Therefore, it is quite convenient to use GC/MS to analyze the amount of monomers present in a resin, the volatile products from polymerization (condensation), or the impurities in a resin. The difficulty in analyzing a resin arises from the fact that it is usually a thick and viscous liquid. The normal procedure of using a solvent to dissolve the resin is tedious and difficult and requires the judicious choice of solvent. In addition, the solvent peak in GC/MS often obscures the small amount of volatile compounds to be analyzed. The use of fixed-volume micro syringes allows simple, direct analysis of the sample without contamination. The micro syringe (DynaTech Mini-injector, from 0.01 to 0.10 microliters, Aldrich Cat. No. Z17,096-8 and Z17,099-2) is dipped directly into the resin sample. Normal injection is made into a GC/MS. The injection port is kept at the appropriate temperature, and the needle is left in the injection port for sufficient time to allow the vaporization of volatile compounds from the needle (about one minute). The needle can be cleaned with appropriate solvents and any residue burned off using a burner. We have used this method successfully to analyze the amount of allyl alcohol (a potent irritant) in a Vibrin resin.

Ambrose Leong, Ph.D.
Billy Sue Hurst Professor of Chemistry
Department of Chemistry
Emory & Henry College
Emory, VA 24327-0947
E-mail: aleong@ehc.edu



Preventing Cross-Contamination in Solvent Stills Connected in Series to the Same Inert Gas Manifold

We have always had the problem of solvents getting contaminated with other solvents, when several stills were connected to a common inert gas source. The flow of the inert gas was designed to entrain the solvent vapors by entering the still with the low-vapor-pressure solvent and ending in a bubbler next to the high-vapor-pressure solvent. But, often, solvents were contaminating each other during the distillation due to changes in the rate of inert gas flow, water flow cutoff, or power failure.

To minimize this contamination, we have designed a new piece of glassware called a Protecting Head, as shown in **Figure 1**. Since the solvent vapors are carried deep inside the side arm, they condense, allowing only the inert gas to escape. In a typical setup, each still is fitted with a protecting head and connected to the common inert gas manifold through T-connectors. Any contaminated solvent accumulated in the side arm is periodically removed through the bottom stopcock, making sure that the accumulated solvent does not reach the bottom of the inner tube.

We have successfully used this setup in the distillation of such solvents as diethyl ether, dichloromethane, hexane, toluene, acetonitrile, and THF. We hope that chemists the world over will find it equally useful.

Maravanji S. Balakrishna, Assistant Professor
Department of Chemistry
Indian Institute of Technology
Powai, Mumbai 400 076, INDIA
E-mail: krishna@chem.iitb.ernet.in






Figure 1
Solvent Still with a Still Protecting Head.

Do you have an innovative shortcut or unique laboratory hint you'd like to share with your fellow chemists? If so, please send it to Aldrich (attn: Lab Notes, Aldrichimica Acta). For submitting your idea, you will receive a complimentary, laminated periodic table poster (Cat. No. **Z15,000-2**). If we publish your lab note, you will also receive an Aldrich periodic table turbo mouse pad (Cat. No. **Z24,409-0**). It is Teflon®-coated, 8½ x 11 in., with a full-color periodic table on the front. We reserve the right to retain all entries for future consideration.



Teflon is a registered trademark of E.I. du Pont de Nemours & Co., Inc.

Update your bookshelf
with the
new 81st edition of the



CRC
Handbook
of Chemistry
& Physics

Z42,457-9 \$129.95

Sigma-Aldrich Worldwide Locations

Argentina

Sigma-Aldrich de Argentina, S.A.
Tel: 54 11 4807 0321
Fax: 54 11 4807 0346
E-mail: info@sigma-aldrich.com.ar

Australia

Sigma-Aldrich Pty., Limited
Free Tel: 1-800-800-097
Free Fax: 1-800-800-096
Tel: 61-2-9841-0555
Fax: 61-2-9841-0500
E-mail: ausmail@sial.com

Austria

Sigma-Aldrich Handels GmbH
Tel: 43 1 605 8110
Fax: 43 1 605 8120
E-mail: sigma@sial.com

Belgium

Sigma-Aldrich nv/sa
Free Tel: 0800-147.47
Free Fax: 0800-147.45
Tel: 03 899.13.01
Fax: 03 899.13.11
E-mail: sigmaldr@sial.be

Brazil

Sigma-Aldrich Quimica Brasil Ltda.
Tel: 55-11-231-1866
Fax: 55-11-257-9079
E-mail: sigmabr@ibm.net

Canada

Sigma-Aldrich Canada Ltd
Free Tel: 800-565-1400
Free Fax: 800-265-3858
Tel: 905-829-9500
Fax: 905-829-9292
E-mail: canada@sial.com

Czech Republic

Sigma-Aldrich s.r.o.
Tel: 00-420-2-2176 1300
Fax: 00-420-2-2176 3300
E-mail: CZECustSV@eurnotes.sial.com

Denmark

Sigma-Aldrich Denmark A/S
Tel: +45 43 56 59 10
Fax: +45 43 56 59 05
E-mail: denorder@eurnotes.sial.com

Eire

Sigma-Aldrich Ireland Ltd
Free Tel: 1 800 200 888
Free Fax: 1 800 600 222
Tel: 01 404 1900
Fax: 01 404 1910
E-mail: EICustsv@eurnotes.sial.com

Finland

Sigma-Aldrich Finland
Tel: (09) 350-9250
Fax: (09) 350-92555
E-mail: finorder@eurnotes.sial.com

France

Sigma-Aldrich Chimie S.a.r.l.
Tel Numero Vert: 0800 21 14 08
Fax Numero Vert: 0800 03 10 52
E-mail: fradvsv@eurnotes.sial.com

Germany

(also SE Europe, the Baltics, Africa and the Middle East)

Sigma-Aldrich Chemie GmbH
Free-call Tel: 0800-5155 000
Free-call Fax: 0800-6490 000
Tel: +49(0)89-6513-0
Fax: +49(0)89-6513-1169
E-mail: deorders@eurnotes.sial.com
For Slovakia, Africa, Middle East and other European countries:
Tel: +49(0)89-6513-1807
Fax: +49(0)89-6513-1877

Greece

Sigma-Aldrich (o.m) Ltd.
Tel: 30-1-9948010
Fax: 30-1-9943831
E-mail: sigald@acci.gr

Hungary

Sigma-Aldrich Kft
Tel: 06-1-269-6474
Fax: 06-1-235-9050
E-mail: info@sigma.sial.hu

India

Sigma-Aldrich Corporation
Bangalore location:
Phone: 91 80 851 8797
Fax: 91 80 851 8358

New Delhi location

Phone: (011) 689 9826
Fax: (011) 689 9827
E-mail: sigma@del2.vsnl.net.in

Israel

Sigma Israel Chemicals Ltd
Toll Free Tel: 1-800-70-2222
Tel: 972-8-948-4222
Fax: 972-8-948-4200
E-mail: sigisr@sigma.co.il

Italy

Sigma-Aldrich S.r.l.
Numero Verde: 800-827018
Tel: 02-33417.310
Fax: 02-38010.737
E-mail: itorder@eurnotes.sial.com

Japan

Sigma-Aldrich Japan K.K.
Tel: 81-3-5821-3111
Fax: 81-3-5821-3170

Korea

Sigma-Aldrich Korea
Toll Free Tel: 080-023-7111
Tel: 03-7824181, 03-7808140
Toll Free Fax: 080-023-8111
Tel: 82-2-783-5211
Fax: 82-2-783-5011
E-mail: order@salkor.co.kr

Malaysia

Sigma-Aldrich (M) Sdn Bhd
Tel: 03-7824181, 03-7808140
Fax: 03-7824067
E-mail: sigalm@po.jaring.my

Mexico

Sigma-Aldrich Quimica, S.A. de C.V.
Tel Gratis: 91-800-0075300
Tel: (72) 76 1600
Fax: (72) 76 1601
E-mail: mexico@sial.com

Netherlands

Sigma-Aldrich Chemie BV
Tel Gratis: 0800 022 9088
Fax Gratis: 0800 022 9089
Tel: 078-620 54 11
Fax: 078-620 54 21
E-mail: nlcustsv@eurnotes.sial.com

Norway

Sigma-Aldrich Norway AS
Tel: 47 23 176 000
Fax: 47 23 176 010
E-mail: nororder@eurnotes.sial.com

Poland

Sigma-Aldrich Sp.z o.
Tel: (48-61) 823-24-81
Fax: (48-61) 823-27-81
E-mail: plcustsv@eurnotes.sial.com

Portugal

Sigma-Aldrich Quimica S.A.
Free Tel: 0800-202180
Free Fax: 0800-202178
E-mail: encomendas:
poorders@eurnotes.sial.com

Russia

Sigma-Aldrich Russia
Tel: 7-095-975-3321
Fax: 7-095-975-4792
E-mail: techcare@online.ru

Singapore

Sigma-Aldrich Pte., Ltd
Tel: (65) 271 1089
Fax: (65) 271 1571
E-mail: sup@sial.com

South Africa

Sigma-Aldrich South Africa (Pty) Ltd
Tel: (011) 397 8886
Fax: (011) 397 8859
E-mail: sigald@icon.co.za

Spain

Sigma-Aldrich Quimica S.A
Free Tel: 900-101376
Free Fax: 900-102028
Tel: 91-6619977
Fax: 91-6619642
E-mail, pedidos:
esorders@eurnotes.sial.com

Sweden

Sigma-Aldrich Sweden AB
Tel: 020-35 05 10
Fax: 020-35 25 22
Outside Sweden Tel: 46-8-742 4200
Outside Sweden Fax: 46-8-742 4243
E-mail: sweorder@eurnotes.sial.com

Switzerland

Aldrich Chemie
A division of Fluka Chemie AG
Tel: 081-755 27 23
Fax: 081-755 28 40
Swiss Free Call: 0800 80 00 80
E-mail: Fluka@sial.com

United Kingdom

Sigma-Aldrich Company Ltd.
Free Tel: 0800 717181
Free Fax: 0800 378538
Tel: 01202 733114
Fax: 01202 715460
E-mail: ukcustsv@eurnotes.sial.com

United States


Aldrich Chemical Company, Inc.
Toll Free Tel: 1-800-558-9160
Toll Free Fax: 1-800-962-9591
Tel: 414-273-3850
Fax: 414-273-4979
E-mail: aldrich@sial.com

Sigma Chemical Company

Toll Free Tel: 1-800-521-8956
Toll Free Fax: 1-800-325-5052
Tel: 314-771-5765
Fax: 314-771-5757
E-mail: sigma@sial.com

Each Dean-Stark trap comes standard with a Teflon® stopcock.

The 12mL and 25mL traps use a 2mm bore; all the rest use a 4mm.



Aldrich Dean-Stark Traps

"Providing Ample Clearance with Thick-Walled Aluminum Heating Mantles"

Designed specifically for azeotropic distillation using solvents lighter than water.

Cap. (mL)	Overall H (mm)	Side-Arm Clearance (mm)	Suggested Flask Size	Joint	Cat. No.	Each
12	205	80	50mL – 5L	14/20	Z20,254-1	\$97.00
25	255	80	50mL – 5L	14/20	Z42,307-6	137.00
25	280	80	50mL – 5L	24/40	Z42,302-5	137.00
100	350	100	5L – 12L	24/40	Z42,303-3	219.00
300	440	170	22L – 50L	24/40	Z42,304-1	249.00
500	455	170	22L – 50L	24/40	Z42,306-8	295.00

For custom Dean-Stark configurations, call the Aldrich Glass Shop or e-mail us at aldglass@sial.com.

FINE CHEMICALS

purity, PLUS.

with **sigma-Aldrich** as your product development partner, high-purity materials are only the **beginning**.

You get a partner with more than 50 years of experience in manufacturing, who can provide raw materials, scale-up, and complex multistep custom synthesis—all with rigorous attention to purity.

You get a partner who can develop your novel chemistry and scale it to commercial volumes while preserving the confidentiality of your ideas and research.

You get a partner who will deliver these quality, high-purity materials to your required specifications—even for the most demanding applications.

- Air-sensitive manufacturing
- Low-temperature synthesis: down to -110°C
- Ultrahigh purification: 99.999+%
- Worldwide manufacturing at flexible production scales: from kilograms to tons

No matter what your high technology application may be, Sigma-Aldrich Fine Chemicals stands ready to be your quality high-purity materials partner. For more information contact us at 800 336-9719 (USA) or visit us at www.sigmaaldrich.com, safe.



www.sigmaaldrich.com

SIGMA-ALDRICH FINE CHEMICALS • 280 SPRING ST. • ST. LOUIS • MISSOURI • 63103 • USA



SIGMA-ALDRICH
FINE CHEMICALS

ALDRICH CHEMICAL COMPANY, INC.
P.O. BOX 355
MILWAUKEE, WISCONSIN 53201 USA



CHANGE SERVICE REQUESTED

PRESORTED STANDARD
U.S. POSTAGE
PAID
MILWAUKEE, WISCONSIN
PERMIT NO. 552

014

SPECIAL SECTION: RESINS FOR COMBINATORIAL CHEMISTRY

Aldrichimica ACTA

VOL. 33, NO. 3 • 2000
(Last issue in 2000)



Flexible
Cross-Linkers
in SPS Resins


Highlights of
Synthetic Radical
Chemistry



ALDRICH[®]

CHEMISTS HELPING CHEMISTS


New Products



Used in the synthesis of nitrogen-containing heterocycles via vinylogous Pummerer reactions of amido sulfoxides.

Padwa, A.; Kuethe, J.T. *J. Org. Chem.* **1998**, *63*, 4256.


53,147-2 2-Iodophenylacetic acid, 97%



Nitrogen-containing heterocycles have recently been prepared from this chiral amine in a tandem aminocarbonylation/Pauson-Khand reaction.

Baisells, J. et al. *Org. Lett.* **1999**, *1*, 1981.


52,554-5 (S)-N-Allyl-α-methylbenzylamine, 97%



This silane is a key intermediate in the synthesis of a silicon linker, useful for the traceless synthesis of pyridine-based tricyclics.

Woolard, F.X. et al. *J. Org. Chem.* **1997**, *62*, 6102.


51,104-8 Allyl(4-methoxyphenyl)dimethylsilane, 98%



Useful for the synthesis of *m*-terphenyls, and, subsequently, novel cyclophanes.^{1,2}

(1) Rajakumar, P.; Srisailas, M. *Tetrahedron Lett.* **1997**, *38*, 5323. (2) Kannan, A. et al. *J. Org. Chem.* **1996**, *61*, 5090.


51,045-9 1,3-Dichloro-2-iodobenzene, 98%



A simple, one-step reaction employing this reagent leads to the direct regioselective dehydrogenation of aldehydes and cyclic ketones.

Shvo, Y.; Arisha, A.H.I. *J. Org. Chem.* **1998**, *63*, 5640.


52,585-5 Diethyl allyl phosphate, 98%



Employed in the synthesis of novel tetrazoles via tetrazolylbenzyl radicals.

Biadatti, T. et al. *Tetrahedron Lett.* **1998**, *39*, 19.


53,199-5 2-Chloro-N-benzylacetamide, 97%



This chiral auxiliary provides practical access to a wide range of α -fluoro-carboxylic acids in high enantiomeric excess.

Myers, A.G. et al. *Tetrahedron Lett.* **1998**, *39*, 1335.


52,561-8 (1S,2S)-Pseudoephedrine α -fluoroacetamide, 98%



Utilized in the catalytic synthesis of optically active allyl α -amino acids.

Fang, X. et al. *J. Org. Chem.* **1999**, *64*, 4844.


51,026-2 Methallyltrimethylsilane, 97%



Synthons for the synthesis of the immunosuppressive agent FK506¹ and the thrombin inhibitor inogatran,² and for preparing new cyclosporin derivatives.³

(1) Baker, R.K. et al. *Tetrahedron Lett.* **1998**, *39*, 229. (2) Preville, P. et al. *Bioorg. Med. Chem. Lett.* **1997**, *7*, 1563. (3) Eberle, M.K. et al. *J. Med. Chem.* **1995**, *38*, 1853.

51,636-8 (S)-(-)-1-(tert-Butoxycarbonyl)-2-piperidine-carboxylic acid, 98% (98% ee/GLC)



51,634-1 (R)-(+)-1-(tert-Butoxycarbonyl)-2-piperidine-carboxylic acid, 98% (98% ee/GLC)

Aldrichimica ACTA

VOL. 33, NO. 3 • 2000

Aldrich Chemical Co., Inc.

1001 West Saint Paul Ave.
Milwaukee, WI 53233 USA

To Place Orders

Telephone 800-558-9160 (USA)
or 414-273-3850
FAX 800-962-9591 (USA)
or 414-273-4979
Mail P.O. Box 2060
Milwaukee, WI 53201 USA

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201 USA

Customer & Technical Services

Customer Inquiries	800-558-9160
Technical Service	800-231-8327
Sigma-Aldrich Fine Chemicals	800-336-9719
Custom Synthesis	800-336-9719
Flavors & Fragrances	800-227-4563
International	414-273-3850
24-Hour Emergency	414-273-3850
Web Site	http://www.sigma-aldrich.com
E-Mail	aldrich@sial.com

To request your **FREE** subscription to the *Aldrichimica Acta*.

please call: 800-558-9160 (USA)

or write: **Attn: Mailroom**
Aldrich Chemical Co., Inc.
P.O. Box 355
Milwaukee, WI 53201-9358

International customers, please contact your local Sigma-Aldrich office.


The *Aldrichimica Acta* is also available on the Internet at <http://www.sigma-aldrich.com>.

Aldrichimica Acta is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich family. © 2000 by Sigma-Aldrich Co. Printed in the United States

“Please
Bother
Us.”

Clint Lane

Clint Lane, President



Professor Janusz Pawliszyn of the University of Waterloo, Waterloo, Canada, kindly suggested that we make sodium tetraethylborate- d_{20} , which is used as a derivatization reagent for the determination of alkyllead and inorganic lead in aqueous samples. Dr. Pawliszyn's method can be adapted for field measurements, and has the potential of becoming a general method for the speciation of inorganic and organometallic compounds.

Yu, X.; Pawliszyn, J. *Anal. Chem.* **2000**, *72*, 1788.

49,504-2

Sodium tetraethylborate- d_{20} ,
98 atom % D

Naturally, we made this useful reagent. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page or on the inside back cover.

About Our Cover

Street in Venice (oil on wood panel, 17 $\frac{3}{4}$ in. x 21 $\frac{1}{4}$ in.) was painted by the American artist John Singer Sargent. Sargent was born in Florence in 1856 to affluent expatriate American parents and spent most of his life in Europe. He lived mainly in Paris until 1884, and then in London, where he died in 1925. His formal instruction in art was in Rome, Florence, and Paris, and, from 1877, he exhibited at the official salons in the French capital. Sargent traveled extensively, copying Velázquez at the Prado in Madrid and Frans Hals in Belgium and Holland, and developing a style that combines close observation with a fluent and painterly technique derived at least partly from these earlier masters. In 1881, he met James McNeill Whistler in Venice, where he painted this picture, probably in the fall of 1882. By the end of the nineteenth century, he had become the most celebrated international society portrait painter of his time, with a wealthy, aristocratic, and fashionable clientele; but he sought every chance he could to escape the constraints of formal portraiture for genre subjects.




Sargent does not paint the colorful, idealized views of scenes like the Piazza San Marco or the Grand Canal, commonly portrayed by many of his contemporaries, but captures a casual moment in one of the back streets of a working class district of Venice. Observed by two men standing near a doorway, a young woman hurries along, her expression preoccupied and aloof, her shawl drawn about her against the chill. The limited range of color, the deep perspective view of the narrow street, and the dark silhouettes of the men, one of whom glances at the woman, create a mood that is almost sinister, and certainly ambiguous. Who is this young woman? Who are the two men? Do they know her? We cannot be sure, but this mystery is one reason that Sargent's painting is so compelling and so memorable.

This painting is a gift of the Avalon Foundation to the National Gallery of Art, Washington, D.C.

Table of Contents

Tailoring Polystyrene Solid-Phase Synthesis Resins: Incorporation of Flexible Cross-Linkers	87
-P. H. Toy, T. S. Reger, & K. D. Janda	
Highlights from Two Decades of Synthetic Radical Chemistry	104
-D. P. Curran	



Lab Notes

Effective Cooling System in Bulb-to-Bulb Distillation

Bulb-to-bulb distillation is a useful laboratory technique for the small-scale purification of liquid substances. In this microscale, short-path technique, a small quantity of material is redistilled using several glass bulbs connected in line. In some Kugelrohr models, heating is carried out by a built-in oven, and control of both temperature and stirring speed is performed by an integrated electronic controller. The flask containing the product to be distilled and the collecting bulbs are connected to a rotary device (**Figure 1**).

Traditionally, cooling of the collecting bulbs has been carried out with dry ice in acetone (-78°C), especially when working with low-boiling substances or when distilling under vacuum. This cooling technique has several drawbacks:

- Acetone is usually applied by a piece of cotton held with tweezers and may accidentally drip on the Kugelrohr apparatus and attack the plastic cover. It may also drip inside the housing, through the control buttons, and damage the electronic circuit boards. Acetone can also dissolve the plastic parts of the front panel and lead to the fusion of the buttons to the housing.
- Acetone, because of its low flash point (-18°C), could ignite if it contacts the electrical system since the rotor can generate sparks.
- From an operational point of view, the ground joints, which connect the glass bulbs, usually are not lubricated with vacuum grease to avoid contaminating the final distillate during its recovery. As a result, and even though the bulb-to-bulb seal is usually good, a reabsorption of acetone through the joints could happen principally when the distillation is carried out under vacuum. In this case, contamination of the distillate with acetone could become a problem, especially when dealing with low-boiling substances.

To avoid these problems, we have designed and used extensively a new, safe, and simple cooling system. In this system, the distillation bulbs are cooled with dry ice without the use of acetone. As containers of solid dry ice (coolers), we use cheap and readily available plastic or cardboard cylindrical tubs of 4 to 5 cm in diameter and 10 cm in length (**Figure 2a**).

To adapt the cylindrical container to the shape of the cooling bulb and to increase the cooling surface, it is recommended that two circular cuts be made on opposite sides of the bottom of the cylinder (**Figure 2b**).

Ángel M. Montaña, Ph.D., and Pedro M. Grima

Departamento de Química Orgánica, Universidad de Barcelona
Martí i Franquès 1-11, 08028-Barcelona, Spain
E-mail: ammontana@qo.ub.es

Do you have an innovative shortcut or unique laboratory hint you'd like to share with your fellow chemists? See the inside back cover for details.




Figure 1: Bulb-to-bulb distillation apparatus with a cooling system of open cylinders filled with dry ice.




Figure 2(a): Detail of the cylindrical dry-ice container of the cooling system.




Figure 2(b): Cross section of the cylindrical dry-ice container showing the circular cut made to adapt the cylinder to the rotating bulb.

Sigma-Aldrich Worldwide Locations

Argentina

Sigma-Aldrich de Argentina, S.A.
Tel: 54 11 4807 0321
Fax: 54 11 4807 0346
E-mail: info@sigma-aldrich.com.ar

Australia

Sigma-Aldrich Pty., Limited
Free Tel: 1-800-800-097
Free Fax: 1-800-800-096
Tel: 61-2-9841-0555
Fax: 61-2-9841-0500
E-mail: ausmail@sial.com

Austria

Sigma-Aldrich Handels GmbH
Tel: 43 1 605 8110
Fax: 43 1 605 8120
E-mail: sigma@sial.co.at

Belgium

Sigma-Aldrich nv/sa
Free Tel: 0800-147 47
Free Fax: 0800-147 45
Tel: 03 899.13.01
Fax: 03 899.13.11
E-mail: sigmaldr@sial.be

Brazil

Sigma-Aldrich Quimica Brasil Ltda.
Tel: 55-11-231-1866
Fax: 55-11-257-9079
E-mail: sigmabr@ibm.net

Canada

Sigma-Aldrich Canada Ltd
Free Tel: 800-565-1400
Free Fax: 800-265-3858
Tel: 905-829-9500
Fax: 905-829-9292
E-mail: canada@sial.com

Czech Republic

Sigma-Aldrich s.r.o.
Tel: 00-420-2-2176 1300
Fax: 00-420-2-2176 3300
E-mail: CZECustSV@eurnotes.sial.com

Denmark

Sigma-Aldrich Denmark A/S
Tel: +45 43 56 59 10
Fax: +45 43 56 59 05
E-mail: denorder@eurnotes.sial.com

Eire

Sigma-Aldrich Ireland Ltd
Free Tel: 1 800 200 888
Free Fax: 1 800 600 222
Tel: 01 404 1900
Fax: 01 404 1910
E-mail: EICustsv@eurnotes.sial.com

Finland

Sigma-Aldrich Finland
Tel: (09) 350-9250
Fax: (09) 350-9255
E-mail: finorder@eurnotes.sial.com

France

Sigma-Aldrich Chimie S.a.r.l
Tel Numero Vert: 0800 21 14 08
Fax Numero Vert: 0800 03 10 52
E-mail: fradvsv@eurnotes.sial.com

Germany

(also SE Europe, the Baltics, Africa and the Middle East)

Sigma-Aldrich Chemie GmbH
FreeCall Tel: 0800-5155 000
FreeCall Fax: 0800-6490 000
Tel: +49(0)89-6513-0
Fax: +49(0)89-6513-1169
E-mail: deorders@eurnotes.sial.com
For Slovakia, Africa, Middle East and other European countries:
Tel: +49(0)89-6513-1807
Fax: +49(0)89-6513-1877

Greece

Sigma-Aldrich (o.m) Ltd.
Tel: 30-1-9948010
Fax: 30-1-9943831
E-mail: sigald@acci.gr

Hungary

Sigma-Aldrich Kft
Tel: 06-1-269-6474
Fax: 06-1-235-9050
E-mail: info@sigma.sial.hu

India

Sigma-Aldrich Corporation
Bangalore location:
Phone: 91 80 851 8797
Fax: 91 80 851 8358

New Delhi location:

Phone: (011) 689 9826
Fax: (011) 689 9827
E-mail: sigma@del2.vsnl.net.in

Israel

Sigma Israel Chemicals Ltd.
Toll Free Tel: 1-800-70-2222
Tel: 972-8-948-4222
Fax: 972-8-948-4200
E-mail: sigisr@sigma.co.il

Italy

Sigma-Aldrich S.r.l.
Numero Verde: 800-827018
Tel: 02-33417.310
Fax: 02-38010.737
E-mail: itorder@eurnotes.sial.com

Japan

Sigma-Aldrich Japan K.K.
Tel: 81-3-5821-3111
Fax: 81-3-5821-3170

Korea

Sigma-Aldrich Korea
Toll Free Tel: 080-023-7111
Toll Free Fax: 080-023-8111
Tel: 82-2-783-5211
Fax: 82-2-783-5011
E-mail: order@safkor.co.kr

Malaysia

Sigma-Aldrich (M) Sdn Bhd
Tel: 03-7824181, 03-7808140
Fax: 03-7824067
E-mail: sigalm@po.jaring.my

Mexico

Sigma-Aldrich Quimica, S.A. de C.V.
Tel Gratis: 91-800-0075300
Tel: (72) 76 1600
Fax: (72) 76 1601
E-mail: mexico@sial.com

Netherlands

Sigma-Aldrich Chemie BV
Tel Gratis: 0800 022 9088
Fax Gratis: 0800 022 9089
Tel: 078-620 54 11
Fax: 078-620 54 21
E-mail: nlcustsv@eurnotes.sial.com

Norway

Sigma-Aldrich Norway AS
Tel: 47 23 176 000
Fax: 47 23 176 010
E-mail: nororder@eurnotes.sial.com

Poland

Sigma-Aldrich Sp.zo.o
Tel: (48-61) 823-24-81
Fax: (48-61) 823-27-81
E-mail: plcustsv@eurnotes.sial.com

Portugal

Sigma-Aldrich Quimica S.A.
Free Tel: 0800-202180
Free Fax: 0800-202178
E-mail, encomendas:
poorders@eurnotes.sial.com

Russia

Sigma-Aldrich Russia
Tel: 7-095-975-3321
Fax: 7-095-975-4792
E-mail: techcare@online.ru

Singapore

Sigma-Aldrich Pte., Ltd
Tel: (65) 271 1089
Fax: (65) 271 1571
E-mail: sup@sial.com

South Africa

Sigma-Aldrich South Africa (Pty) Ltd
Tel: (011) 397 8886
Fax: (011) 397 8859
E-mail: sigald@icon.co.za

Spain

Sigma-Aldrich Quimica S.A.
Free Tel: 900-101376
Free Fax: 900-102028
Tel: 91-6619977
Fax: 91-6619642
E-mail, pedidos:
esorders@eurnotes.sial.com

Sweden

Sigma-Aldrich Sweden AB
Tel: 020-35 05 10
Fax: 020-35 25 22
Outside Sweden Tel: 46-8-742 4200
Outside Sweden Fax: 46-8-742 4243
E-mail: sweorder@eurnotes.sial.com

Switzerland

Aldrich Chemie
A division of Fluka Chemie AG
Tel: 081-755 27 23
Fax: 081-755 28 40
Swiss Free Call: 0800 80 00 80
E-mail: Fluka@sial.com


United Kingdom

Sigma-Aldrich Company Ltd.
Free Tel: 0800 717181
Free Fax: 0800 378538
Tel: 01202 733114
Fax: 01202 715460
E-mail: ukcustsv@eurnotes.sial.com

United States

Aldrich Chemical Company, Inc.
Toll Free Tel: 1-800-558-9160
Toll Free Fax: 1-800-962-9591
Tel: 414-273-3850
Fax: 414-273-4979
E-mail: aldrich@sial.com
Sigma Chemical Company
Toll Free Tel: 1-800-521-8956
Toll Free Fax: 1-800-325-5052
Tel: 314-771-5765
Fax: 314-771-5757
E-mail: sigma@sial.com

New edition available in December 2000




March's Advanced Organic Chemistry: Reactions, Mechanisms, and Structure

5th ed., M. B. Smith & J. March, John Wiley & Sons, New York, NY, 2001, 1,824pp. Hardcover.

Contains over 600 reactions, including 50 new sections with detailed descriptions of the basic principles underlying those reactions. Describes the structure of organic compounds, including chemical bonding and stereochemistry. Features a reaction index classified by type of compound synthesized. Contains more than 20,000 references to original papers—5,000 new to this edition.

Z53,719-5

Do you have an innovative shortcut or unique laboratory hint you'd like to share with your fellow chemists? If so, please send it to Aldrich (attn: Lab Notes, Aldrichimica Acta). For submitting your idea, you will receive a complimentary, Aldrich periodic table poster (Cat. No. **Z54,320-9**). If we publish your lab note, you will also receive an Aldrich periodic table turbo mouse pad (Cat. No. **Z24,409-0**). It is Teflon®-coated, 8½ x 11 in., with a full-color periodic table on the front. We reserve the right to retain all entries for future consideration.



Teflon is a registered trademark of E. I. du Pont de Nemours & Co., Inc.

At
Sigma-Aldrich,
our
chemists
understand
there's one
precious
ingredient
you can't
afford to
waste...



...Time

Reduce product development time with Sigma-Aldrich. In all stages of custom manufacturing, time is an extremely valuable commodity. That's why you need a proven partner like Sigma-Aldrich.

Over the past 50 years, we've built an extensive technology database. We have over 33,000 preparations on file for both organic and biochemical processes. So, more than any other company, we have the resources, expertise and know-how you need to help you bring your product to market on time and on budget.

Our organic manufacturing facilities in five countries specialize in complex, multistep organic syntheses of:

- High-quality intermediates
- Air-sensitive organometallic compounds
- cGMP active pharmaceutical ingredients

To ensure consistent high quality and product purity, we also offer complete analytical services, including method development and validation, as well as long-term stability testing.

Looking to control project time—and costs? Call us at Sigma-Aldrich Fine Chemicals. *We'll help you put time on your side.*



SIGMA-ALDRICH
FINE CHEMICALS

ALDRICH CHEMICAL COMPANY, INC.
P.O. BOX 355
MILWAUKEE, WISCONSIN 53201 USA



CHANGE SERVICE REQUESTED

1951–2001: FIFTY YEARS OF CHEMISTS HELPING CHEMISTS

Aldrichimica ACTA

VOL. 34, NO. 1 • 2001



Preparation of Optically
Active α -Amino Acids

Alkoxy methylenemalonates
in Organic Synthesis



ALDRICH[®]

CHEMISTS HELPING CHEMISTS

NEW PRODUCTS

53,182-0 *N,N'*-Di-(*tert*-butoxycarbonyl)thiourea, 97%

This diprotected thiourea is widely used in the synthesis of heterocycles, including a pentacyclic guanidine system as an intermediate to ptilomycin A,¹ and, recently, in the preparation of *p-N,N'*-bis-Boc-guanidophenol, a key intermediate for the preparation of a series of aryl *o*-aroylbenzoates as serine protease inhibitors.²

(1) Nagasawa, K. et al. *Tetrahedron* **2000**, *56*, 187. (2) Jones, P.B.; Porter, N.A. *J. Am. Chem. Soc.* **1999**, *121*, 2753.

54,036-6 1,1-Diethoxy-3-methyl-2-butene, 97%

This acetal has been utilized in the synthesis of the related 1-alkoxy-3-phenylselenoalkenes and 3-phenylselenoalkanes,¹ and in the preparation of 2,2-dimethylchromenes from electron-deficient phenols.²

(1) Nishiyama, Y. et al. *Tetrahedron Lett.* **1998**, *39*, 8685. (2) North, J.T. et al. *J. Org. Chem.* **1995**, *60*, 3397.

54,345-7 Dimethyl L-tartrate, 99% (99+% ee/GLC)

Recent citations for this versatile chiral building block include the enantioselective synthesis of α -bromo carbonyl and carboxylic acid derivatives,¹ and the synthesis of enantiomerically pure spirane porphyrazines.²

(1) Boyes, S.A.; Hewson, A.T. *J. Chem. Soc., Perkin Trans. 1* **2000**, 2759. (2) Hachiya, S. et al. *Tetrahedron* **2000**, *56*, 6565.

54,944-0 (*R*)-(+)-1-(*tert*-Butoxycarbonyl)-2-pyrrolidinemethanol, 97%

This protected pyrrolidinemethanol was used in a recent study of 3-substituted indoles as potential antimigraine drugs.¹ It was also utilized in the preparation of a novel, potent, and selective 5-HT₁ antagonist via elongation of the side chain followed by transformation into the protected piperidineethylpyrrolidine, subsequent deprotection and arylsulfonylation.²

(1) Sternfeld, F. et al. *J. Med. Chem.* **1999**, *42*, 677. (2) Lovell, P.J. et al. *ibid.* **2000**, *43*, 342.

55,673-3 2-Fluoro-6-methoxybenzaldehyde, 98%

On conversion to a dinitro ester, the aryl fluorine acts as a useful leaving group early in the synthesis of the macrolactam portion of the ansamycin antibiotic (+)-thiazinotrienomycin E.

Smith, A.B., III; Wan, Z. *J. Org. Chem.* **2000**, *65*, 3738.

54,331-4 2-Fluoro-5-iodobenzaldehyde, 97%

The related iodobenzo[b]thiophenecarboxylate ester has been synthesized in moderate yield from this dihalo aldehyde.

Bridges, A.J. et al. *Tetrahedron Lett.* **1992**, *33*, 7499.

54,540-1 (1*R*,2*S*)-1-Phenyl-2-(1-pyrrolidinyl)-1-propanol, 98%

54,897-9 (1*S*,2*R*)-1-Phenyl-2-(1-pyrrolidinyl)-1-propanol, 98%

These norephedrine derivatives are useful chiral mediators and have applications in the enantioselective addition of acetylides to carbonyl compounds. Examples include the synthesis of the HIV-1 reverse transcriptase inhibitor DMP 266¹ and the enantioselective addition of diethylzinc to aldehydes.²

(1) Pierce, M.E. et al. *J. Org. Chem.* **1998**, *63*, 8536. (2) Soai, K. et al. *ibid.* **1991**, *56*, 4264.

54,174-5 Tri-*O*-acetyl- β -D-arabinosylbromide, 95%

Among the many applications for this protected bromosugar, the solid-state reaction to make glycopyranosyl pyrimidine nucleosides,¹ the novel synthesis of thioglycosides,² and the high-yield preparation of pyranoid glycals³ are some of the ones that have been cited.

(1) Im, J. et al. *Tetrahedron Lett.* **1997**, *38*, 451. (2) Pákulski, Z. et al. *Tetrahedron* **1994**, *50*, 2975. (3) Kovács, G. et al. *ibid.* **1999**, *55*, 5253.

52,869-2 2-(6-Bromohexyloxy)tetrahydro-2*H*-pyran, 97%

This THP-protected bromoalcohol is an important building block used in the synthesis of alkaloids,¹ symmetrical olefins,² and 14-membered macrocyclic ethers.³

(1) Kaiser, A. et al. *J. Org. Chem.* **1999**, *64*, 3778. (2) Poulain, S. et al. *Tetrahedron Lett.* **1996**, *37*, 7703. (3) Clyne, D.S.; Weiler, L. *Tetrahedron* **1999**, *55*, 13659.

51,587-6 Imidazole, trifluoromethanesulfonate salt, 97%

This triflate has been used as a coupling agent for oligonucleotide synthesis¹ and as a reagent for the synthesis of aryl triflates.²

(1) Hayakawa, Y.; Kataoka, M. *J. Am. Chem. Soc.* **1998**, *120*, 12395. (2) Effenberger, F.; Mack, K.E. *Tetrahedron Lett.* **1970**, 3947.

54,039-0 (3,4-Dihydro-1-naphthyl)oxy)trimethylsilane, 97%

The titanium tetrachloride promoted reaction of this trimethylsilyl enol ether with ethylene oxide affords the homoaldol-type product in moderate yield.¹ In other studies, *C*₂-symmetric copper complexes have been shown to catalyze the enantioselective amination of,² or ketomalonate addition to,³ this enol silane.

(1) Lalic, G. et al. *Tetrahedron Lett.* **2000**, *41*, 763. (2) Evans, D.A.; Johnson, D.S. *Org. Lett.* **1999**, *1*, 595. (3) Reichel, F. et al. *J. Chem. Soc., Chem. Commun.* **1999**, 1505.

Aldrichimica ACTA

VOL. 34, NO. 1 • 2001

Aldrich Chemical Co., Inc.

1001 West Saint Paul Ave.
Milwaukee, WI 53233 USA

To Place Orders

Telephone 800-558-9160 (USA)
or 414-273-3850
FAX 800-962-9591 (USA)
or 414-273-4979
Mail P.O. Box 2060
Milwaukee, WI 53201 USA

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201 USA

Customer & Technical Services

Customer Inquiries 800-558-9160
Technical Service 800-231-8327
Sigma-Aldrich Fine Chemicals 800-336-9719
Custom Synthesis 800-336-9719
Flavors & Fragrances 800-227-4563
International 414-273-3850
24-Hour Emergency 414-273-3850
Web Site <http://www.sigma-aldrich.com>
E-Mail aldrich@sial.com

To request your FREE subscription to the *Aldrichimica Acta*,

please call: 800-558-9160 (USA),
write: Attn: Mailroom
Aldrich Chemical Co., Inc.
P.O. Box 355
Milwaukee, WI 53201-9358

or email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office.

The *Aldrichimica Acta* is also available on the Internet at <http://www.sigma-aldrich.com>.

Aldrichimica Acta is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich family. © 2001 by Sigma-Aldrich Co. Printed in the United States

About Our Cover

Edgar Degas's *Ballet Scene* (pastel on cardboard, 30¼ in. x 43¾ in.) was executed ca. 1907 when the artist was in his seventies, and not long before his career was cut short by his growing blindness. The subject of the picture is one indelibly identified with Degas who, unlike his impressionist friends, more often chose to represent urban scenes than the country landscapes they preferred. Life in the cafés, at the racetrack, in the theatre, at the opera, and in particular at the ballet, was what most attracted him.



However, it was not the glamour, the drama, and the sumptuous spectacle of a ballet performance that interested him, but the behind-the-scenes work that goes on before the dancers go on stage. His goal was to distill from the constantly changing movements and postures of the dancers a single moment which would capture the essence of motion, and it was the exercises, the rehearsals, the minutes before a performance that gave him the opportunity to do this.


This goes some way towards explaining the medium he used for this work, which is not an oil painting on canvas but a picture executed in pastel chalks, a dry medium similar to the charcoal he often used to sketch quickly the varied life of the city of Paris all around him. Through his depiction of varied postures and gestures, Degas leads the eye along a curving line from the lower right corner of the picture to the left, contrasting the more static attitudes of the three figures in the right foreground with the animated movements of the dancers further back on the left, to capture in an essentially motionless work of art a sense of the dynamic actions of a group of dancers at work.

This painting is in the Chester Dale Collection at the National Gallery of Art, Washington, D.C.

“Please
Bother
Us.”

Clint Lane

Clint Lane, President



Professor AHCÈNE Boumendjel of the Université Joseph Fourier (La Tronche, France) kindly suggested that we make *N*-methoxy-*N*-methyl-2-(phenylsulfonyl)-acetamide.¹ This amide is useful for the homologation of alkyl halides to α,β -unsaturated *N*-methoxy-*N*-methylamides, which are valuable intermediates in the synthesis of natural products and therapeutics.²

(1) Beney, C.; Boumendjel, A.; Mariotte, A.-M. *Tetrahedron Lett.* **1998**, *39*, 5779. (2) Sibi, M. *P. Org. Prep. Proced. Int.* **1993**, *25*, 15.

51,139-0

N-Methoxy-*N*-methyl-2-(phenylsulfonyl)acetamide, 96%

Naturally, we made this useful reagent. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page or on the inside back cover.

Table of Contents

The Preparation of Optically Active α -Amino Acids from the Benzophenone Imines of Glycine Derivatives.....	3
—M. J. O'Donnell	
Alkoxy methylenemalonates in Organic Synthesis.....	20
—V. Milata	



Lab Notes


Inexpensive and Easily Constructed Space-Saving "Spaghetti Tubing Apparatus" for Creating Inert Atmospheres within a Large Number of Reaction Vessels

I read a recent Lab Notes article,¹ which prompted me to share another innovative method for creating and maintaining inert atmospheres within a large number of reaction vessels or reagent bottles. The setup is easily constructed from inexpensive equipment generally found in a laboratory, and occupies less space than a conventional manifold line.

The system has been used in our laboratories for several years with great success, supplying inert atmospheres to multiple reactions and easily coping with the higher pressures required for cannulation. Based on the "spaghetti tubing" manifold described by Casey et al.,² this apparatus consists of a single source of inert gas connected to the side-arm of a B29/32 cone adapter, which has been fitted at each end with rubber septa. The single source of inert gas is split into multiple supplies (comfortably up to six) by needle-tipped, 3-mm (o.d.) spaghetti lines piercing the septa. A positive pressure of inert gas in the apparatus is maintained by means of an outflow passing through a bubbler assembled from a boiling tube, an inverted Pasteur pipette, a needle and a rubber septum. The various outflow supply lines are easily moved to individual reaction vessels allowing each reaction to be individually purged with inert gas.

The present trend toward the efficient synthesis of compounds via parallel synthesis should make this apparatus invaluable for the pharmaceutical chemist. In addition, its cost effectiveness and small size make it ideal for university chemistry departments, where restricted space within fume cupboards and tight budgetary constraints may exist.

References: (1) Flemer, S., Jr. *Aldrichimica Acta* **1998**, *31*, 34.
(2) Casey, M.; Leonard, J.; Lygo, B.; Procter, G. *Advanced Practical Organic Chemistry*; Blackie & Son: London, UK, 1990; pp 108-109.



Materials for Constructing the Spaghetti Line and Bubbler:

- B29/32 socket/cone adapter with "T" connection (equivalent to Aldrich cat. no. **Z41,580-4**) (fitted with detachable plastic connector for easy coupling of tubing)
- Suba-Seal® No. 45 (Aldrich cat. no. **Z16,732-0**)
- Suba-Seal® No. 57 (Aldrich cat. no. **Z16,735-5**)
- Silicone tubing (Versilic) (o.d. 3 mm; bore x wall, 1x1 mm)
- Pasteur pipette (9.0 in. long) (Aldrich cat. no. **Z19,061-6**)
- Disposable needles (16 ga. x 1 in.). Using grips to hold the metal needle, pliers are used to remove the purple plastic sheath. The blunt end of the needle is then pushed into the spaghetti tubing) (Aldrich cat. no. **Z19,256-2**).

Suba-Seal is a registered trademark of William Freeman & Co., Ltd.

David Bebbington, D.Phil.

Previous Address:
Senior Scientist
Department of Chemistry
Cerebrus plc
Oakdene Court
613 Reading Road, Winnersh
Wokingham RG41 5UA
United Kingdom

Current Address:
Staff Investigator
Vertex Pharmaceuticals Ltd.
88 Milton Park
Abingdon
Oxfordshire OX14 4RY
United Kingdom
E-mail: david_bebbington@vpharm.com

Do you have an innovative shortcut or unique laboratory hint you'd like to share with your fellow chemists? See the inside back cover for details.

Sigma-Aldrich Worldwide Locations

Argentina

Sigma-Aldrich de Argentina, S.A
Tel: 54 11 4807 0321
Fax: 54 11 4807 0346
E-mail: info@sigma-aldrich.com.ar

Australia

Sigma-Aldrich Pty., Limited
Free Tel: 1-800-800-097
Free Fax: 1-800-800-096
Tel: 61-2-9841-0555
Fax: 61-2-9841-0500
E-mail: ausmail@sial.com

Austria

Sigma-Aldrich Handels GmbH
Tel: 43 1 605 8110
Fax: 43 1 605 8120
E-mail: sigma@sial.com.at

Belgium

Sigma-Aldrich nv/sa
Free Tel: 0800-147.47
Free Fax: 0800-147.45
Tel: 03 899.13.01
Fax: 03 899.13.11
E-mail: sigmaldr@sial.be

Brazil

Sigma-Aldrich Quimica Brasil Ltda.
Tel: 55-11-231-1866
Fax: 55-11-257-9079
E-mail: sigmabr@ibm.net

Canada

Sigma-Aldrich Canada Ltd.
Free Tel: 800-565-1400
Free Fax: 800-265-3858
Tel: 905-829-9500
Fax: 905-829-9292
E-mail: canada@sial.com

Czech Republic

Sigma-Aldrich s.r.o.
Tel: 00-420-2-2176 1300
Fax: 00-420-2-2176 3300
E-mail: CZEcustSV@eurnotes.sial.com

Denmark

Sigma-Aldrich Denmark A/S
Tel: +45 43 56 59 10
Fax: +45 43 56 59 05
E-mail: denorder@eurnotes.sial.com

Eire

Sigma-Aldrich Ireland Ltd
Free Tel: 1 800 200 888
Free Fax: 1 800 600 222
Tel: 01 404 1900
Fax: 01 404 1910
E-mail: ElCustsv@eurnotes.sial.com

Finland

Sigma-Aldrich Finland
Tel: (09) 350-9250
Fax: (09) 350-9255
E-mail: finorder@eurnotes.sial.com

France

Sigma-Aldrich Chimie S.a.r.l
Tel Numero Vert: 0800 21 14 08
Fax Numero Vert: 0800 03 10 52
E-mail: fradsvs@eurnotes.sial.com

Germany

(also SE Europe, the Baltics, Africa and the Middle East)

Sigma-Aldrich Chemie GmbH
Free-call Tel: 0800-5155 000
Free-call Fax: 0800-6490 000
Tel: +49(0)89-6513-0
Fax: +49(0)89-6513-1169
E-mail: deorders@eurnotes.sial.com
For Slovakia, Africa, Middle East and other European countries:
Tel: +49(0)89-6513-1807
Fax: +49(0)89-6513-1877

Greece

Sigma-Aldrich (o.m) Ltd.
Tel: 30-1-9948010
Fax: 30-1-9943831
E-mail: sigald@accr.gr

Hungary

Sigma-Aldrich Kft
Tel: 06-1-269-6474
Fax: 06-1-235-9050
E-mail: info@sigma.sial.hu

India

Sigma-Aldrich Corporation
Bangalore location:
Phone: 91 80 851 8797
Fax: 91 80 851 8358

New Delhi location:

Phone: (011) 689 9826
Fax: (011) 689 9827
E-mail: sigma@del2.vsnl.net.in

Israel

Sigma Israel Chemicals Ltd.
Toll Free Tel: 1-800-70-2222
Tel: 972-8-948-4222
Fax: 972-8-948-4200
E-mail: sigisr@sial.com

Italy

Sigma-Aldrich S.r.l.
Numero Verde: 800-827018
Tel: 02-33417.310
Fax: 02-38010.737
E-mail: itorder@eurnotes.sial.com

Japan

Sigma-Aldrich Japan K.K.
Tel: 81-3-5821-3111
Fax: 81-3-5821-3170
E-mail: sialjp@sial.com

Korea

Sigma-Aldrich Korea
Toll Free Tel: 080-023-7111
Toll Free Fax: 080-023-8111
Tel: 82-2-783-5211
Fax: 82-2-783-5011
E-mail: order@safkor.co.kr

Malaysia

Sigma-Aldrich (M) Sdn Bhd
Tel: 03-7824181, 03-7808140
Fax: 03-7824067
E-mail: sigalm@po.jaring.my

Mexico

Sigma-Aldrich Quimica, S.A. de C.V.
Tel Gratis: 91-800-0075300
Tel: (72) 76 1600
Fax: (72) 76 1601
E-mail: mexico@sial.com

Netherlands

Sigma-Aldrich Chemie BV
Tel Gratis: 0800 022 9088
Fax Gratis: 0800 022 9089
Tel: 078-620 54 11
Fax: 078-620 54 21
E-mail: nlcustsv@eurnotes.sial.com

Norway

Sigma-Aldrich Norway AS
Tel: 47 23 176 000
Fax: 47 23 176 010
E-mail: nororder@eurnotes.sial.com

Poland

Sigma-Aldrich Sp.z.o.
Tel: (48-61) 823-24-81
Fax: (48-61) 823-27-81
E-mail: plcustsv@eurnotes.sial.com

Portugal

Sigma-Aldrich Quimica S.A.
Free Tel: 0800-202180
Free Fax: 0800-202178
E-mail, encomendas:
poorders@eurnotes.sial.com

Russia

Sigma-Aldrich Russia
Tel: 7-095-975-3321
Fax: 7-095-975-4792
E-mail: techcare@online.ru

Singapore

Sigma-Aldrich Pte., Ltd.
Tel: (65) 271 1089
Fax: (65) 271 1571
E-mail: sup@sial.com

South Africa

Sigma-Aldrich South Africa (Pty) Ltd
Tel: (011) 397 8886
Fax: (011) 397 8859
E-mail: sigald@icon.co.za

Spain

Sigma-Aldrich Quimica S.A.
Free Tel: 900-101376
Free Fax: 900-102028
Tel: 91-6619977
Fax: 91-6619642
E-mail, pedidos:
esorders@eurnotes.sial.com

Sweden

Sigma-Aldrich Sweden AB
Tel: 020-35 05 10
Fax: 020-35 25 22
Outside Sweden Tel: 46-8-742 4200
Outside Sweden Fax: 46-8-742 4243
E-mail: sweorder@eurnotes.sial.com

Switzerland

Aldrich Chemie
A division of Fluka Chemie AG
Tel: 081-755 27 23
Fax: 081-755 28 40
Swiss Free Call: 0800 80 00 80
E-mail: Fluka@sial.com

United Kingdom

Sigma-Aldrich Company Ltd.
Free Tel: 0800 717181
Free Fax: 0800 378538
Tel: 01202 733114
Fax: 01202 715460
E-mail: ukcustsv@eurnotes.sial.com

United States

Aldrich Chemical Company, Inc.
Toll Free Tel: 1-800-558-9160
Toll Free Fax: 1-800-962-9591
Tel: 414-273-3850
Fax: 414-273-4979
E-mail: aldrich@sial.com

Sigma Chemical Company


Toll Free Tel: 1-800-521-8956
Toll Free Fax: 1-800-325-5052
Tel: 314-771-5765
Fax: 314-771-5757
E-mail: sigma@sial.com

Aldrich High-Vacuum Valve Adapters

Features PTFE J. Young valves for high-vacuum applications. For use with 1/4 in. i.d. tubing.

Joint	Cat. No.
14/20	Z54,967-3
19/22	Z54,968-1
24/40	Z54,970-3
29/32	Z54,971-1
34/45	Z54,973-8

Do you have an innovative shortcut or unique laboratory hint you'd like to share with your fellow chemists? If so, please send it to Aldrich (attn: Lab Notes, Aldrichimica Acta). For submitting your idea, you will receive a complimentary, Aldrich periodic table poster (Cat. No. **Z54,320-9**). If we publish your lab note, you will also receive an Aldrich periodic table turbo mouse pad (Cat. No. **Z24,409-0**). It is Teflon®-coated, 8½ x 11 in., with a full-color periodic table on the front. We reserve the right to retain all entries for future consideration.



Teflon is a registered trademark of E.I. du Pont de Nemours & Co. Inc.

ALDRICH CELEBRATES ITS

1951 + 2001

*What a Remarkable Half Century
This Has Been!*

From simple beginnings in a rented garage on Milwaukee's East Side, Aldrich has grown into a leading chemical company, offering over 40,000 products for research and industrial applications. Aldrich has also become one of the most recognizable brand names in chemistry and related fields worldwide.

As we celebrate reaching this milestone, Aldrich would like to thank our **customers**, whose support and loyalty over the years has made this phenomenal success possible.

Aldrich also owes a debt of gratitude to our **employees**, past and present, whose loyalty and dedication to the company have made it what it is today, and to the many **scientists** whose collaboration has permitted us to remain on the leading edge of chemistry research and development.

As Aldrich looks to the future, we would like to pledge to our customers that our commitment to providing top-notch service and the highest-quality products will continue to be our focus in the next half century. So will our unwavering support to the advancement of the chemical sciences.

Aldrich is proud of the last 50 great years of achievements & service to the chemistry community, and looks forward to the next 50 years with confidence and enthusiasm!



FIFTIETH ANNIVERSARY

ALDRICH CHEMICAL COMPANY, INC.
P.O. BOX 14508
ST. LOUIS, MISSOURI 63178 USA



1951-2001: FIFTY YEARS OF CHEMISTS HELPING CHEMISTS

Aldrichimica ACTA


VOL. 34, NO. 2 • 2001



 **ALDRICH**[®]
CHEMISTS HELPING CHEMISTS

New Products

21,502-3 3-Iodothiophene, 97%




Has been used recently in the stereoselective synthesis of (*E*)-(3-thienyl)vinylgermanes via a Stille coupling reaction,¹ and in the preparation and study of 3-(diacetoxyiodo)thiophene.²

(1) David-Quillot, F. et al. *Tetrahedron Lett.* **2000**, *41*, 9981. (2) Togo, H. et al. *J. Org. Chem.* **2000**, *65*, 8391.

1g \$33.90

54,332-2 2-Phenylindole-3-carboxaldehyde, 97%




Has been employed in a study of the inhibition of tubulin polymerization and the growth of breast cancer cells.

Gaspard, R. et al. *J. Med. Chem.* **1998**, *41*, 4965.

1g \$14.60; 5g \$48.50

55,075-2 *N* α ,*N*-(*im*)-Di-Boc-L-histidine methyl ester, 97%




Was employed in the preparation of the corresponding imidazole-2-thione during the development of the first synthesis of L-(+)-ergothioneine.

Xu, J.; Yadan, J.C. *J. Org. Chem.* **1995**, *60*, 6296.

1g \$14.00; 5g \$49.60

53,944-9 1-(*N*-Boc-aminomethyl)-4-(aminomethyl)benzene, 97%



1g \$31.50; 5g \$105.10

Has been utilized for the solid- and solution-phase syntheses of oligomeric thioureas,¹ the preparation of a bipyrrrole-based [2]catenane,² and, very recently, for the synthesis of model receptors for dicarboxylates and monosaccharides.³

(1) Smith, J. et al. *J. Org. Chem.* **1996**, *61*, 8811. (2) Andrievsky, A. et al. *J. Am. Chem. Soc.* **1998**, *120*, 9712. (3) Benito, J.M. et al. *J. Org. Chem.* **2001**, *66*, 1366.

55,065-5 Heptylmagnesium bromide, 1.0M solution in diethyl ether

100mL \$115.20; 1L \$640.00



In situ generation of the corresponding organozinc compound was a key step in the recent use of this Grignard reagent in the stereoselective synthesis of protected alkylhydroxypyrrolidinones.¹ In another report, this reagent was used in the enantioselective opening of an epoxide en route to (S)-(+)-8-hydroxyhexadecanoic acid.²

(1) Huang, P.Q. et al. *Synth. Commun.* **2000**, *30*, 2259. (2) Shimajo, M. et al. *Tetrahedron* **2000**, *56*, 9281.

54,330-6 2,7-Di-*tert*-butylfluorene, 98%

25g \$131.00




Two recent examples of the use of this building block include the preparation of new group 4 metal complexes containing aminofluorenyl ligands,¹ and a novel, soluble analog of the Fmoc group.²

(1) Miller, S.A.; Bercaw, J.E. *Organometallics* **2000**, *19*, 5608. (2) Stigers, K.D. et al. *J. Org. Chem.* **2000**, *65*, 3858.

54,151-6 2-Chloro-5,5-dimethyl-1,3,2-dioxaphosphorinane, 97%

25g \$74.00



54,152-4 5,5-Dimethyl-1,3,2-dioxaphosphorinan-2-one, 96%


1g \$9.10; 5g \$30.30

Reaction of the chlorophosphorinane with Cp_2NbH_3 proceeds via a carbene-like insertion mechanism to give a derived cationic diphosphite complex in high yield.¹ Other workers have prepared tungsten (Z)-1,2-diphosphite-alkene tetracarbonyl and (*E*)-vinyl-phosphite pentacarbonyl complexes from this reagent via phosphite Fischer carbene complexes.²

(1) Nikonov, G.I. et al. *J. Organomet. Chem.* **1997**, *547*, 183. (2) Barluenga, J. et al. *Organometallics* **1997**, *16*, 3732.

49,740-1 (*R*)-(+)-2-Methyl-2-propanesulfonamide, 98%

1g \$123.35




Developed by Ellman, this chiral sulfonamide has found widespread use in diastereoselective enolate alkylations,¹ the synthesis of protected chiral amines,² the Strecker synthesis of α -alkyl α -amino acids,³ and the preparation of novel ligands for asymmetric Lewis acid catalysis.⁴

(1) Backes, B.J. et al. *J. Org. Chem.* **1999**, *64*, 5472. (2) Borg, G. et al. *Tetrahedron Lett.* **1999**, *40*, 6709. (3) Davis, F.A. et al. *J. Org. Chem.* **2000**, *65*, 8704. (4) Owens, T.D. et al. *J. Am. Chem. Soc.* **2001**, *123*, 1539.

53,645-8 2-Chloro-4-methoxypyrimidine, 98%

1g \$24.85; 5g \$82.80




The introduction of the 2-pyrimidinyl functionality, in the synthesis of methylene phosphonate analogs of thymidine 3'-phosphate, was achieved via a nucleophilic substitution reaction by an intermediate thioglycoside on this chloropyrimidine.¹ A number of related pyrazolopyrimidines were prepared for study as potential cytoprotective antiulcer agents.²

(1) Yokomatsu, T. et al. *Heterocycles* **1999**, *50*, 21. (2) Ikeda, M. et al. *Chem. Pharm. Bull.* **1996**, *44*, 1700.

52,007-1 (*R*)-(+)-*N*-(Boc)-O-(TBDMS)serinol, 97%

1g \$32.40; 5g \$108.00



The TEMPO oxidation of this protected serinol yields the α -amino aldehyde in good yield and without racemization.¹ This amino alcohol has also been employed in the diastereoselective synthesis of 1,2-diamines.²

(1) Jurczak, J. et al. *Tetrahedron* **1998**, *54*, 6051. (2) Gonda, J. et al. *Synthesis* **1993**, 729.

Aldrichimica ACTA

VOL. 34, NO. 2 • 2001

Aldrich Chemical Co., Inc.

1001 West Saint Paul Ave.
Milwaukee, WI 53233 USA

To Place Orders

Telephone 800-558-9160 (USA)
or 414-273-3850
FAX 800-962-9591 (USA)
or 414-273-4979
Mail P.O. Box 2060
Milwaukee, WI 53201 USA

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201 USA

Customer & Technical Services

Customer Inquiries	800-558-9160
Technical Service	800-231-8327
Sigma-Aldrich Fine Chemicals	800-336-9719
Custom Synthesis	800-336-9719
Flavors & Fragrances	800-227-4563
International	414-273-3850
24-Hour Emergency	414-273-3850
Web Site	http://www.sigma-aldrich.com
E-Mail	aldrich@sial.com

To request your **FREE** subscription to the *Aldrichimica Acta*,

please call: **800-558-9160** (USA),

write: **Attn: Mailroom**
Aldrich Chemical Co., Inc.
P.O. Box 355
Milwaukee, WI 53201-9358

or email: **sams-usa@sial.com**

International customers, please contact your local Sigma-Aldrich office.

The *Aldrichimica Acta* is also available on the internet at <http://www.sigma-aldrich.com>.

All prices listed in this publication are subject to change without notice.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

For worldwide contact information, please see the inside back cover.

Aldrichimica Acta is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich family. © 2001 by Sigma-Aldrich Co. Printed in the United States

About Our Cover


Paintings of all kinds have graced the covers of our *Aldrichimica Acta* over the years. Most were old masters, a few (like the one on the cover of Vol. 1, No. 1) were alchemical, and a few others were modern. Perhaps, none has a more direct connection to chemistry than the one reproduced on the cover of this issue. It (oil on canvas, 40 in. x 50 in.) was painted by Thomas Phillips in London in 1816, and depicts a 24-year-old Michael Faraday watching his teacher, Professor W. T. Brande, make **Prussian Blue**. At the time, Faraday was not yet the great chemist and physicist he later became—arguably the ablest British scientist of the 19th century.

What a fitting cover for the *Acta* issue celebrating Aldrich's 50th birthday!

The interested reader can learn more about this painting in: (a) Bader, A. End of the Mystery. *Chem. Br.*, July 2001, in press. (b) Bader, A. Out of the Blue. *Chem. Br.*, November 1997, p 24.

This painting is in the collection of Alfred Bader Fine Arts, Milwaukee, WI.


Full-color reproductions (11 in. x 14 in.) of this painting are available for a nominal fee to cover postage and handling. Please call 800-558-9160 (USA) and specify the *Prussian Blue* painting (Z52,866-8) from the cover of *Aldrichimica Acta*, Volume 34, Number 2.



“Please
Bother
Us.”

Clint Lane

Clint Lane, President



Professor David Crich of the Department of Chemistry, University of Illinois at Chicago, kindly suggested that we offer the highly sterically hindered base, 2,4,6-tri-*tert*-butylpyrimidine (TTBP), as an alternative to 2,6-di-*tert*-butylpyrimidine, and its 4-methyl and 4-*tert*-butyl derivatives, in glycosylation reactions and in the preparation of vinyl triflates.

Crich, D.; Smith, M.; Yao, Q.; Picione, J. *Synthesis* 2001, 323.

54,996-7

2,4,6-Tri-*tert*-butylpyrimidine,
97% **1g \$21.15; 5g \$70.50**

Naturally, we made this useful reagent. It was no bother at all, just a pleasure to be able to help.


Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page or on the inside back cover.

Table of Contents

A Half Century of Chemists Helping Chemists: Aldrich from 1951 to 2001...35
—S. J. Firsan

Activated 2*H*-Azirines as Dienophiles and Electrophiles.....51
—T. L. Gilchrist

Reactive Chemical Hazard Evaluation in the Scale-Up of Chemical Processes....61
—G. Amery



Lab Notes

1-Chlorobutane—A Useful Solvent for Recrystallizations

When carrying out a recrystallization, there are sometimes problems with the choice of solvent and recourse is made to mixed solvents. Of these, the combination ethyl acetate–petroleum ether is perhaps the most common pairing. Usually, one obtains a hot solution by addition of the more polar solvent to a hot suspension of the would-be solute in the less polar solvent. However, on cooling there is frequently development of two phases followed by precipitation of the solute as a solid mass rather than as crystals of good texture.

After some experiences as outlined above, we sought an alternative approach. What we were seeking was a solvent of low but definite polarity, with a convenient boiling point. We chose 1-chlorobutane, bp 77–78 °C, and we have found this to be a convenient solvent for the recrystallization of oxime ethers,¹ nitrogen ylides,

such as $\text{Me}_3\text{N}^+-\text{N}^--\text{CN}$ (with no butylation of the ylide),² and 3-pyrrolylpentanoic acid.³ In unpublished work, we have used 1-chlorobutane to advantage in the recrystallization of the nitrogen ylide $\text{Me}_3\text{N}^+-\text{N}^--\text{C}(=\text{O})\text{Me}$ (again without butylation), and the bicyclic compounds 1-nitrocamphene and camphor-4-carboxylic acid.

Although the bulk dielectric constant does not tell the whole story about a solvent, 1-chlorobutane, interestingly, has a higher dielectric constant (7.28) than ethyl acetate (6.08); that of *n*-hexane is 1.89.⁴ Specific solute–solvent interactions also need to be taken into account.

We recommend that 1-chlorobutane be considered as a solvent for recrystallization in appropriate cases subject to a couple of caveats: (a) manipulations with the solvent be carried out in a fume cupboard, and (b) one is mindful that

1-chlorobutane may be attacked by nucleophiles, though its reactivity is, at best, sluggish.

For recrystallizations that require a higher temperature, 1-chloropentane is worth consideration.

References: (1) Bradley, G. W.; Morris, D. G. *J. Chem. Res. (S)* **1993**, 220. (2) Chardin, A.; Berthelot, M.; Laurence, C.; Morris, D. G. *J. Phys. Org. Chem.* **1995**, 8, 626. (3) Ryder, K. S.; Morris, D. G.; Cooper, J. M. *Langmuir* **1996**, 12, 5681. (4) *CRC Handbook of Chemistry and Physics*, 79th ed.; Lide, D. R., Ed.; CRC Press LLC: Boca Raton, FL, 1998–1999; Section 6, pp 139–161.

David G. Morris, Ph.D.
Department of Chemistry
Joseph Black Building
University of Glasgow
Glasgow G12 8QQ, Scotland
United Kingdom
E-mail: d.morris@chem.gla.ac.uk

Making Optimal Use of Peristaltic Pumps

Laboratories involved in process development (e.g., fermentation research) frequently need to pump liquids at constant rates for specified periods of time. These liquids may be culture media, acids, alkalis, or special feed supplements. Peristaltic pumping is the chief means by which this is achieved. However, it may occasionally happen that the rate at which a liquid needs to be delivered is much lower than what can be achieved with the available pump head and tubing. This is especially the case if the pump head cannot accommodate tubings of different diameters. One simple and inexpensive solution is to pass a tubing of a smaller diameter through a short length of the pump head tubing and fix this tube-in-a-tube arrangement into the pump head. The pressure that is applied to the outer tubing by

the pump rollers is transferred to the outer walls of the inner tubing, facilitating peristaltic action. Sometimes, the inner tubing may slip out of the pump head tubing when the pump is operating. To prevent this, the inner tubing may be fastened to physical supports using rubber bands or pieces of string near the two points where it emerges from the pump head.

If it is found difficult to pass the thinner tubing through the pump head tubing, smearing its outer walls with a lubricant (an oil, for instance) may help. Furthermore, the thinner tubing may be easily passed through the thicker one if twisting motions are used. If suction of liquid does not occur, a slight adjustment of the gap between the rollers and the wall of the pump head may be required.

We have successfully used this arrangement in our laboratory. It does not appear to have any drawbacks, either with regard to the pump or the tubing. In fact, it has a distinct advantage: the wear on the inner tubing is reduced, since it is shielded by the thicker pump head tubing.

Nelliah Hariharan, Ph.D.
Senior Research Scientist (Biotechnology)
Wockhardt Research Centre
D-4 MIDC, Chikalthana
Aurangabad - 431 210 (M.S.)
India
E-mail: nelmail@yahoo.com

Do you have an innovative shortcut or unique laboratory hint you'd like to share with your fellow chemists? See the inside back cover for details.

A Half Century of Chemists Helping Chemists

Aldrich from 1951 to 2001

Sharbil J. Firsan
Aldrich Chemical Company, Inc.
1001 W. Saint Paul Avenue
Milwaukee, WI 53233, USA
E-mail: sfirsan@sial.com

Outline

1. Introduction
2. What's in a Name?
 - 2.1. Alfred R. Bader
 - 2.2. Aldrich Chemical Company, Inc.
 - 2.3. MNNG and Other Products
3. Critical Years: 1954 and 1955
 - 3.1. A Crucial Decision
 - 3.2. Suberic Acid and D-Penicillamine
4. The Remainder of the 1950s
 - 4.1. Rapid Growth
 - 4.2. The Rare Chemical Library
 - 4.3. Aldrich Chemical Co Ltd
 - 4.4. Aldrich Chemie KG
 - 4.5. Custom Synthesis
5. The 1960s: A Decade of Transformation
 - 5.1. Early to Mid-1960s
 - 5.2. Business, Art, and Chemistry
 - 5.3. Joint Ventures
 - 5.3.1. ALFA Inorganics, Inc.
 - 5.3.2. Aldrich-Europe
 - 5.3.3. Riedel-deHaën® Laboratory Chemicals
 - 5.3.4. Aldrich-APL, L.L.C. (AAPL™)
6. Great Opportunities and Profound Changes (1970s)
 - 6.1. Aldrich-Boranes, Inc.
 - 6.2. Stable Isotopes
 - 6.3. Sigma-Aldrich Corporation
 - 6.4. Floyd Green's Stains and Dyes
 - 6.5. Nonchemical Products
 - 6.5.1. Laboratory Equipment (Techware)
 - 6.5.2. The Aldrich Glass Shop
 - 6.6. New Lines of Business
 - 6.7. Craftsmen in Chemistry?
7. Post-Merger Era
 - 7.1. The Succession
 - 7.2. The 80s and 90s
 - 7.3. The Present and the Future
 - 7.3.1. "A Member of the Sigma-Aldrich Family"
 - 7.3.2. www.sigma-aldrich.com/aldrich
8. The Role of Science and Scientists
 - 8.1. Scientists' Contributions
 - 8.2. Herbert C. Brown
 - 8.3. Rewarding Excellence
9. Valued Customers, Dedicated Employees
10. Acknowledgments
11. References and Notes

1. Introduction

This August, Aldrich Chemical Company turns fifty. What a remarkable fifty-year period this has been! While it is today a "household" name among chemical scientists and practitioners worldwide, it is hard to imagine that Aldrich is only a half century old. On the following pages, I will take a brief look back at the past fifty years. While the story of Aldrich has been told and retold in one form or another many times,^{1,2,3} it is my sincere hope that the fresh approach I am using will prove to be of most interest to our *Aldrichimica Acta* readers. This approach traces the development of Aldrich through the key chemical products and business ventures that played a crucial role in the uninterrupted success that the company has enjoyed over the past fifty years. Since the success of any great enterprise is dependent on the contributions of a legion of dedicated individuals who believe in its mission, the roles that key people played in the development of Aldrich will also be highlighted.¹

2. What's in a Name?

2.1. Alfred R. Bader


Alfred Robert Bader, a young Austrian immigrant and a chemistry graduate student at Harvard University, first entertained the idea of starting a company to sell research chemicals in 1949—on the suggestion of Warren Stockwood, the storeroom supervisor at Harvard's Department of Chemistry³—after being disappointed with the service he received from the leading supplier of research chemicals at the time. Acting on the premise that chemists needed a wider array of research chemicals and better service, he and Jack Nathan Eisendrath, a Milwaukee attorney, founded Aldrich Chemical Company, Inc. on August 17, 1951. Jack became the company's first president.⁵ In a curious twist of fate, the company was named following a coin toss,



not after either of the founders, but after Eisendrath's fiancée, Bettie Mae Aldrich.⁷ Had Bader won the coin toss, the company would have been named Daniels Chemical Company, after Helen Ann Daniels, Bader's fiancée and future wife.²

2.2. Aldrich Chemical Company, Inc.

In 1951, the company operated first from Eisendrath's office on 161 W. Wisconsin Avenue and, later that year, from a rented garage located on N. Farwell Avenue in Milwaukee's East Side. It had three part-time employees: Alfred, Jack, and Lorraine Worby (née Neau). Lorraine first worked for Aldrich part-time, 4–5 hours on Wednesday nights, then became the first full-time Aldrich employee in August of 1954.⁸ From 1951 to 1954, Alfred sowed the seeds of what later became very important collaborations and acquisitions through visits to small chemical producers in Continental Europe and the United Kingdom. Noteworthy are two visits in 1952 to Fluka AG Chemische Fabrik in St. Gallen, Switzerland, and Heidenheimer Chemisches Laboratorium (HCL) in Heidenheim, Germany. These two companies as well as many others served as



Scheme 2. Large-Scale Preparation of Suberic Acid.

Later on, Marvin was instrumental in guiding the growth of Aldrich, both as legal counsel and director of Aldrich and of Sigma-Aldrich. As Alfred put it, Marvin "was the moving spirit in the growth of Aldrich".¹⁹ Marvin also played a key role in the negotiations leading up to the merger of Aldrich Chemical Company, Inc. and Sigma International, Ltd. in 1975; but more on that later! Beginning in January 1955, Helen Bader, Alfred's wife since July 1952, began working full-time for Aldrich.^{20,21} On September 01, 1955, George Skeff, also a former PPG laboratory technician, was hired as a full-time laboratory technician, following the departure of Anthony D. Kontowicz. George retired from Aldrich on April 14, 1989, after close to 34 years of employment!

3.2. Suberic Acid and D-Penicillamine


More importantly, Aldrich delivered its first bulk order, 500 lb of suberic acid, to Du Pont toward the end of 1955 (Scheme 2).³

Furthermore, one of the important products that Aldrich began to offer that same year was D-penicillamine or 3-mercapto-D-valine (**5**), an antirheumatic α -amino acid and, at the time, a promising orphan drug, which chelates and removes the copper accumulated in the liver of patients suffering from the rare disorder known as Wilson's Disease.²²

4. The Remainder of the 1950s

4.1. Rapid Growth

The remainder of the 1950s was characterized by rapid growth in sales and the number of products offered. This growth in sales led to the expansion of the physical facilities and the hiring of additional staff. In February 1956, Beverly Horick, the fourth non-owner employee, was hired on a part-time basis.²³ The next employee hired was Stella Ward. By 1958, Aldrich had about a dozen employees. In the same year, Aldrich purchased and moved to a three-story, 27,000-sq ft building at 2369 N. 29th Street in Milwaukee's inner city. Two years later (1960), a larger, 70,000-sq ft building was purchased nearby at 2371 N. 30th Street. The 29th Street building housed the R&D and Production departments, while the 30th Street



building had administrative offices and the QC, Packaging, and Warehousing departments.²⁴ The 29th and 30th Street buildings remained Aldrich's headquarters until 1969, when they were condemned by the Milwaukee County Expressway and Transportation Commission in preparation for building the east-west Park Freeway (which was never built). In 1969, Aldrich relocated to 940 W. St. Paul Avenue into an eight-story building that previously belonged to General Electric Company.²⁵

The 7th edition of the Aldrich catalog (December 1955) and its supplements featured over 1,600 products.²⁶ Some of the products offered in the late fifties became success stories. Among the most interesting ones are dicyclohexylcarbodiimide (DCC, **6**), a useful reagent in peptide coupling reactions; *p*-tolylsulfonylethylmethyl nitrosamide (Diazald[®], **7**), a diazomethane precursor that is safer than MNNG; and *p*-phenylazomaleinanil (**8**), a reagent employed for the characterization of conjugated dienes (Figure 2). Other interesting products from this period include lithium borohydride (LiBH₄), sodium tetraphenyl boron (NaBPh₄), diketene (C₄H₆O₂), and triallylamine (C₆H₉N).

The births of several important company entities, as well as other significant developments, also took place in the late fifties.

4.2. The Rare Chemical Library

The Rare Chemical Library (RCL) grew out of the collecting and salvaging of valuable research samples of retiring or deceased academic researchers and from other sources. While the RCL was initially part of a separate company, the Alfred Bader Chemical Corporation, it was sold to Aldrich Chemical Co. on December 20, 1965, just prior to a public offering of 100,000 Aldrich shares.² Over the years, large-scale contributions of samples to the library came from, among others, the personal collection of Joe Karabinos (who founded Carbolabs.



Photo courtesy of the Helen Bader Foundation

Helen A. Bader (née Daniels), 1982.



Photo courtesy of Bettie Aldrich Eisendrath

Bettie Aldrich Eisendrath (née Aldrich), 1951.



Photo courtesy of A.R. Bader

Left to right: Lorraine, George, Stella, and Beverly—1st, 3rd, 5th, and 4th Aldrich employee, respectively (1978).





Figure 2. Best-Selling Products from the Period 1955–1959.



Sigma-Aldrich facility in Gillingham, Dorset, United Kingdom.



Sigma-Aldrich facility in Steinheim am Albuch, Germany.

Inc.) and from the laboratories of such chemical luminaries as Henry Gilman, George Wittig, Robert Woodward, and Louis and Mary Fieser. The RCL currently boasts over 90,000 listings of hard-to-find chemicals. Over the years, the RCL has been invaluable to researchers in the chemical sciences and has led to the discovery and commercialization, by others, of some very valuable chemical commodities, e.g., Roundup® (Monsanto Co.), based on lead compounds obtained from the RCL.

4.3. Aldrich Chemical Co Ltd

In 1959, Aldrich's British subsidiary formally began as Ralph N. Emanuel, Ltd., and was equally owned by the Emanuel and Bader families. It had less than \$1,000 in sales in its first year,²⁶ but sales grew briskly leading to a rapid expansion of the physical facilities and the number of employees. Between 1969 and 1973, Aldrich began gradually acquiring Ralph N. Emanuel, Ltd., and turning it into a wholly owned subsidiary. It was then renamed Aldrich Chemical Co. Ltd. Ralph Emanuel became an Aldrich and, afterwards, a Sigma-Aldrich director. Business growth led to the company moving to the current site in Gillingham, Dorset, U.K., where a distribution center and a new cGMP manufacturing facility are presently located. In 1986, Sigma-Aldrich purchased Bristol Organics, a small manufacturer and long-time Aldrich supplier of fluoroaromatics, and shortly thereafter integrated its operations with those of the larger British subsidiary. Three years after the acquisition of Fluka Chemie by Sigma-Aldrich Corp., the operation of Fluka's subsidiary in Glossop, Derbyshire, was transferred to the Gillingham facility in September 1992. The British subsidiary has a special significance to the corporation, not only because it has become Britain's largest supplier of research chemicals, but also because it has been a source of a number of individuals in corporate leadership positions. Moreover, success of the British subsidiary encouraged the company to open other branches in Continental Europe, where strong competition from well-established fine chemicals companies had been a concern.

4.4. Aldrich Chemie KG

The story of Aldrich's subsidiary in Germany is somewhat different. It starts out with Heidenheimer Chemisches Laboratorium (HCL) in Heidenheim, Germany, acting as Aldrich's best supplier for most of the 50s. HCL was then operated by Dr. Ernst Reif, a chemist, and Gerhard Keppler, a businessman. Following, an industrial accident at HCL and a legal setback for the company, Aldrich became involved in the restructuring and refinancing of its operations. It was renamed

EGA-Chemie KG²⁹ and moved to Steinheim am Albuch—with Aldrich owning about 80% of it. Later (1971),³⁰ Aldrich bought the remaining 20% of EGA and renamed it Aldrich-Chemie GmbH & Co. KG. Like its British counterpart, the German subsidiary has been a tremendous success story, and has grown steadily in capabilities, personnel, and facilities. It now manufactures a range of products [e.g., 2,4-dimethylbenzaldehyde (9), 6-hydroxydihydrotheaspirane (10), and 6-acetoxydihydrotheaspirane (11)—all three are important flavoring raw materials for the food industry (Figure 3)] and serves a corporate warehousing function for all markets in Continental Europe from warehouses in three German towns: Schnelldorf, Steinheim, and Seelze. In 1975, it started producing a full German language edition of the Aldrich catalog concurrently with the English language editions. Two decades later, the activities of all Sigma-Aldrich brands in Germany were combined in one legal entity, Sigma-Aldrich Chemie GmbH, in order to streamline their operations. Perhaps one name stands out more than any other and is credited for most of the early success of Aldrich Chemie—that of Dr. Alfred Griesinger. Dr. Griesinger joined the company in March 1963, owned an interest in EGA-Chemie (1965–1970), and served in various important capacities. He later became a director of Aldrich-Chemie KG, and remained with the company until his untimely death in August 1997.³¹

4.5. Custom Synthesis

The late fifties and early sixties also witnessed the growing importance of custom synthesis and bulk sales. In these early days, custom synthesis was formally one of the business activities of the separate corporation, Alfred Bader Chemical Corp., which was sold to Aldrich on December 20, 1965.²⁷ Over the years, custom synthesis became an important function within the Aldrich Production department, and, together with bulk sales, evolved into Sigma-Aldrich Fine Chemicals (SAFC), currently one of four strategic business units within Sigma-Aldrich Corporation. SAFC concentrates on worldwide large-scale manufacturing and sales. Past and present custom synthesis customers are some of the best-known chemical and pharmaceutical companies in the world. Perhaps some of the more interesting custom synthesis projects that Aldrich worked on in the sixties involved the preparation of *tert*-butoxycarbonyl azide (12), BSA (13), and acryloyl chloride (14).³² Some of the hard-to-find products offered in bulk (100–1,000 lb) in the late fifties included several dimethylphenols (15),

trimethylphenols (**16**), 2-methylresorcinol (**17**), and dihydroxybenzoic acids (**18**) (Figure 4). Along with the R&D group, the custom synthesis team routinely carries out significant process improvement and scale-up projects as well as the manufacture of newly introduced products.

5. The 1960s: A Decade of Transformation

Firstly, Aldrich went from being a privately owned company to being a publicly traded one. Secondly, the growth in sales began to occur, not only as a result of the addition of new products, but also as a result of joint ventures and acquisitions. Finally, important developments took place within the company, such as the birth of this magazine (1967) and the transformation of the Aldrich catalog from a simple list of products and prices into a "handbook", as a consequence of the inclusion of useful factual information about the compounds being offered.

5.1. Early to Mid-1960s

1962 was, in many respects, a watershed year for the young company. Annual sales reached the \$1 million mark for the first time, and the Aldrich catalog grew to 303 pages, as the number of products offered swelled to 10,000. Aldrich became Janssen Pharmaceutica's sales agent in the US, and started ALFA Inorganics, a joint venture with Metal Hydrides, Inc. Dr. John Biel joined the company as Director of Research, replacing Edmund (Pete) J. Eisenbraun, Harvey B. Hopps was hired as an R&D group leader, and Bernard (Bernie) E. Edelstein joined as a chemist. Bernie went on to become one of the company's directors, its secretary, and then its first Executive Vice President (1974). In 1962, William Buth was Aldrich's General Manager; he later became the first Aldrich Vice President. In the same year (1962), J. T. Baker Chemical Co. attempted to buy Aldrich for \$1.5 million, but was rebuffed.³¹

In the mid-sixties, the prior practice of listing only products that were not offered by Eastman Kodak's Fine Chemicals division was abandoned in favor of listing products based on their usefulness and marketability. An interesting offering from this period is 9-amino-1,2,3,4-tetrahydroacridine hydrochloride hydrate (Tacrine hydrochloride; **19**), which was first introduced as an Aldrich product in 1963, and is now sold by Warner-Lambert, a division of Pfizer Inc., under the trade name COGNEX® for the treatment of mild to moderate dementia of the Alzheimer's type.^{34,35}




Figure 3. Important Flavoring Raw Materials Manufactured by the German Subsidiary.





Figure 4. Custom Synthesis (12–14) and Hard-to-Find Products Offered in Bulk (15–18) in the Late Fifties.

5.2. Business, Art, and Chemistry

In 1965 and then in 1966, the Bader family sold some of their shares in the company first to a select group of chemists and friends and then to the public at large. By 1965, sales had almost doubled to about \$1.8 million, from \$1 million in 1962, and the number of employees had grown to over 100, of which 15 were chemists (7 with a Ph.D.).³³ The Aldrich Catalog/Handbook took on its now familiar name and look, as the 1967/1968 edition was the first to have a painting, *The Quill Cutter* by Paulus de Lesire, on its cover. Contrary to popular belief, the idea for placing a painting on the cover of the catalog came, not from Bader—the art lover and collector—but from Bernie Edelstein, an Aldrich employee.³⁶

Following the public offerings, sales and the share price rose steadily as Aldrich grew by expanding into new lines of business and entering into (exclusive) distributorship agreements with a number of commodity chemical producers and other companies. Aldrich also acquired stakes in a number of U.S. based, small chemical producers, such as Hexagon and Kaplop Laboratories, but later divested itself of these stocks after the companies ran into business difficulties.

The *Aldrichimica Acta* evolved from the *Kardindex Sheets* that Aldrich used to mail to its best customers to keep them informed of its newest product offerings.³⁷ A preview issue was printed in the fall of 1967, and then



publication on a regular basis started in 1968. Richard K. Vitek, Aldrich's Director of Marketing at the time, became its first editor—albeit for a very short period of time. A classical alchemical Dutch painting by




Figure 5. Ultrahigh-Purity Products Manufactured by AAPL™.



Aldrich's manufacturing site near Sheboygan, Wisconsin.

Thomas Wyck (17th century) from the Alfred Bader art collection was reproduced on the cover of the first issue of 1968. Of this very same issue, 10,000 copies were printed, whereas, today, over 130,000 copies of each issue are distributed worldwide free of charge. Alfred Bader, who coined the name *Aldrichimica Acta*, came up with the idea of placing a painting on its cover,¹⁷ following the precedent-setting reproduction of a painting on the cover of the 1967/1968 edition of the Aldrich Catalog/Handbook (vide supra).

5.3. Joint Ventures

5.3.1. ALFA Inorganics, Inc.

This 50/50 joint venture with Metal Hydrides, Inc. (later became Ventron, Inc. and then a part of Thiokol Corp.) began in 1962 and ended in 1967. It was created to market inorganics, organometallics, and others such as organoboron and organo-arsenic reagents.³⁸ Even though the joint venture lasted only five years, it helped Aldrich learn a great deal about inorganics and set the stage for Aldrich to expand into this market sector in the 1970s.

5.3.2. Aldrich-Europe

In the late fifties and early sixties, Janssen Pharmaceutica (JP) of Beerse, Belgium, had become one of Aldrich's better suppliers. It was no surprise then that, in 1962, Aldrich became JP's sales agent in the United States. Thus started the good relationship between the two companies. It led in May of 1970 to the creation of Aldrich-Europe, a wholly owned division of JP, charged with distributing Aldrich products in Continental Europe. The joint venture lasted until June 1982. Its dissolution opened the way for Sigma-Aldrich to start subsidiaries and sales offices in many countries in Continental Europe. Today, these number 18.

5.3.3. Riedel-deHaën® Laboratory Chemicals

Fast-forwarding to the present, Sigma-Aldrich first formed a joint venture, Riedel-deHaën® Laborchemikalien GmbH & Co. KG, with Allied Signal in 1997 to market and sell mainly analytical reagents and solvents carrying the Riedel-deHaën® brand name. Three years later, the joint venture became a wholly owned subsidiary of Sigma-Aldrich Corporation and is currently offering close to 4,000 products belonging to three general types: Karl Fischer reagents for water determination, standards for environmental analysis, and high-purity solvents.

5.3.4. Aldrich-APL, LLC. (AAPL™)

Also in the present, Aldrich has a majority stake in a joint venture with APL Engineered Materials to produce select inorganics of extremely high purity and low moisture content (ultradry) intended for the high-tech market. This collaboration began in September 1995 and presently operates from a manufacturing facility in Urbana, Illinois. Dr. John Long, who had been hired by Aldrich in the 1970s to head the inorganics production laboratories in Milwaukee, was initially charged with managing the Urbana facility. Some of the more popular,

ultrahigh-purity products that are presently manufactured by AAPL™ include anhydrous gallium(III) chloride (20), phosphoric acid (21), anhydrous cesium iodide (22), and anhydrous sodium iodide (23) (Figure 5).

6. Great Opportunities and Profound Changes (1970s)

6.1. Aldrich-Boranes, Inc.

The greatest opportunity came when Professor H. C. Brown of Purdue University convinced Aldrich to further develop and commercialize the hydroboration technology and organoborane chemistry that he had developed and patented. This led, in May 1972, to the establishment of Aldrich-Boranes, Inc., a wholly owned Aldrich subsidiary, created to manufacture, among others, hydroboration reagents and products. Aldrich-Boranes, Inc. began operation in September 1972, with Dr. Harvey B. Hopps as its manager, Professor Brown as one of its directors, and a small development group of chemists headed by Clinton F. Lane, an enthusiastic young Ph.D., who had trained with Professor Brown at Purdue. Some of the first compounds manufactured by Aldrich-Boranes were borane-THF (24), 9-BBN (25), borane-methyl sulfide (26), and compounds 27-30 (Figure 6).³⁹ A multitude of others followed in rapid succession. The early development of this chemistry has been described in several reviews by Lane.⁴⁰ In the early 80s, Aldrich-Boranes, Inc. was integrated into Aldrich.

In the first five years, Aldrich-Boranes, Inc. operated from the production laboratories on West St. Paul Avenue in Milwaukee's city center. Following the purchase, in December 1977, of a laboratory building and a property in the town of Wilson, Sheboygan County, Wisconsin, its operations were moved to this site in March 1978.⁴¹ This became the nucleus of Sigma-Aldrich's current 513-acre manufacturing site at 5485 County Road V. The range of products manufactured at the site has long expanded to include pharmaceutical intermediates, air-sensitive reagents, various organometallics, cGMP products, high-purity solvents, and many other compound types.

6.2. Stable Isotopes

Also in 1972, Aldrich acquired Diaprep, an Atlanta based manufacturer of deuterated solvents, and with it two experienced chemists—Tom Wickersham and Bob Askins. Tom remained with Aldrich until his retirement in 1999; he spent most of his career in the Stains & Dyes division, which he helped grow into an important part of the business (vide infra). With the acquisition of Diaprep, Aldrich started the production of

deuterated solvents; this production has expanded rapidly and considerably to the point that Aldrich is presently perhaps the world's largest producer of deuterated solvents. Moreover, through special agreements with companies such as Isotec Inc., Aldrich has also been able to significantly expand its offerings of products labeled with other stable isotopes, such as ^6Li , ^7Li , ^{10}B , ^{11}B , ^{12}C , ^{13}C , ^{15}N , and ^{18}O . In February 2001, Sigma-Aldrich purchased Isotec, Inc., thus becoming the leader in the stable isotopes market.

6.3. Sigma-Aldrich Corporation

The other profound event occurred three years later, in August 1975, when Aldrich Chemical Co. merged with Sigma International, Ltd. of St. Louis, Missouri, to form Sigma-Aldrich Corporation. At the time, both companies were publicly owned, with Aldrich the leading supplier of organic research chemicals and Sigma the leading supplier of research biochemicals. Having dodged several takeover attempts, Aldrich had, in 1967, approached Sigma—then a privately owned company—with an offer to merge, but was rebuffed. By 1975, however, changing trends in chemical research and the synergy to be realized from their complementary business practices and product offerings convinced the two companies to finally merge.

Dan Broida, Sigma's president at the time of the merger, became chairman of the board, while Alfred Bader, Aldrich's president, became president of Sigma-Aldrich Corp. At the time, the new company did not have a CEO. As anticipated, the merged companies drew on each other's strengths—Sigma's emphasis on quality and service and Aldrich's emphasis on introducing new products and maintaining good relationships with suppliers.

A little over a year prior to the merger (April 23, 1974), four Aldrich departmental managers—Robert Gorzek, Irwin (Ike) Klundt, Charles (Chuck) Pouchert, and Edward Segrin—were promoted to vice presidents. On the same date, Vice President Bernard Edelstein was promoted to Executive Vice President.⁴² The first and second Aldrich vice presidents, William Buth and John Biel, had left the company in 1973 and 1968, respectively.⁴³ Also in 1974, David R. Harvey (see Section 7.1) became Aldrich's Vice President of European Operations.⁴⁴

6.4. Floyd Green's Stains and Dyes

In 1973, Dr. Floyd J. Green, a widely respected authority on biological stains and dyes, founded Aristo Custom Chemicals, Inc. in Cincinnati, Ohio, to manufacture




Figure 6. Some of the Early, and Still Popular, Organoboranes.

biological stains. He sold Aristo to Sigma-Aldrich in 1977 and moved to Milwaukee to become an Aldrich vice president in charge of the company's newly created Stains and Dyes division. This division currently boasts ca. 1,300 products, and offers an attractive handbook on stains, dyes, and indicators that was originally prepared by Dr. Green.⁴⁵ Two of the early dyes sold by the division were pararosaniline chloride and acetate. Presently, the Stains and Dyes division operates from a 110,000-sq ft facility, at 230 S. Emmer Lane in Milwaukee, that was purchased in the early seventies.


6.5. Nonchemical Products

6.5.1. Laboratory Equipment (Techware)


Charles J. Pouchert, a long-time Aldrich employee and Manager of its Quality Control department in the early 70s, promoted the idea that Aldrich should sell books that are useful to chemists. Not surprisingly, the first nonchemical item that Aldrich offered (cat. no. Z10,000-5)⁴⁶ was the first edition (1970) of *The Aldrich Library of Infrared Spectra*, that Pouchert edited. The idea of offering a laboratory equipment item was advanced by Dr. Harvey B. Hopps, at the time Aldrich's Manager of Technical Services. This first laboratory equipment item was a Diazald[®] kit (cat. no. Z10,025-2) and was listed for sale for the first time in the very early 70s.²⁴ Thus started the Techware division of Aldrich—with much input and encouragement from Edward J. Segrin, Aldrich's Sales Manager at the time. This division expanded rapidly over the years into such product areas as books, electronic media products, and glassware. It currently offers over 13,000 nonchemical items⁴⁷—everything the practicing chemist needs in addition to chemicals.

6.5.2. The Aldrich Glass Shop

The Aldrich glass shop grew simply enough from the need to repair laboratory glassware in-house. In the early 70s, and by arrangement with Daytime Vacuum Products, a glass shop was set up on the 4th floor of the



2001 edition of the Aldrich Stable Isotopes catalog.



940 W. St. Paul Avenue building. When the lone glass blower retired, Aldrich purchased the shop and looked to hire a glass blower to operate it. An early hire in 1974 left after a few months.⁴⁸ However, in 1975, an able



Techware CD-ROM product.



Sophisticated piece of glassware manufactured by the Aldrich Glass Shop.



Company logo from the early days.

glass blower, Dieter Damrow, was hired to operate the glass shop. Dieter stayed with the shop until his retirement in 1999, and has had a greater impact on the growth and development of the shop than anyone else. During Dieter's tenure, the glass shop grew steadily and its mission expanded to include the manufacture of new glass apparatus and

custom glassware, and the taking on of outside repair jobs. Back in 1985, the glass shop was relocated to newer quarters in the 2905 W. Hope Avenue facility.

6.6. New Lines of Business

As is clearly evident from the preceding discussion, the 70s were an exciting period for Aldrich. In addition to those mentioned in the previous sections, Aldrich availed itself of other business opportunities. Thus, in 1971, Aldrich launched its line of biochemical products with Irwin Klundt as its technical manager.⁴² However, this line of products had a brief, independent existence up until Aldrich's merger with Sigma in 1975. In 1978—free from the restrictions of the separation agreement that led to the dissolution of the Aldrich-Ventron joint venture (see Section 5.3.1)—Aldrich hired Dr. John Long, a promising young inorganic chemist, to spearhead its development and production of inorganic products. It wasn't long before Aldrich was competing successfully with ALFA Inorganics, Inc. The inorganic line of products is thriving today, enhanced by the establishment of the joint venture with APL Engineered Materials in 1995 (see Section 5.3.4). Together with the line of organometallics, it offers close to 10,000 products for, among others, the high-technology markets.

6.7. Craftsmen in Chemistry?

Just as many other companies have had to, Aldrich also had to change with the times. The "Craftsmen in Chemistry" slogan, which was in use in the late 60s, was abandoned a few years later (1975), since it was deemed politically incorrect. It wasn't until early 1978 that a new slogan, "chemists helping chemists in research and industry", began appearing in company literature. Presently, this motto is slowly being phased out in favor of its shorter version, "chemists helping chemists". Aldrich's logo also underwent a face-lift; an earlier version is depicted here. Fortunately, and unlike many other prominent chemical companies, Aldrich did not succumb to the pressure of dropping the word "chemical" from its name, during a period of time when this word had become a public relations liability.

7. Post-Merger Era

7.1. The Succession

With Aldrich now a company within Sigma-Aldrich Corporation, Alfred Bader became President of the merged company and remained President of Aldrich until early 1981. He was succeeded as President of

Aldrich by Dr. David R. Harvey (1981–1986), Dr. Jai P. Nagarkatti (1987–1999), and Dr. Clinton F. Lane (1999–Present).⁴⁹

David R. Harvey, an Oxford University graduate, started out in Aldrich-UK in August 1974 as Vice President of European Operations. As president of Aldrich, he oversaw, among others, the relocation of the distribution center to the Hope Avenue facility and the construction of the laboratory building at the Sheboygan County site. He is also credited with starting the company's Flavors & Fragrances division. David went on to become President of Sigma-Aldrich Corp. and then, in 2000, its Chairman, President, and CEO.

A graduate of East Texas State University, Jai P. Nagarkatti started his career at Aldrich in 1976 as a Process Development Chemist in the Production Laboratory. He progressed through the ranks to become Vice President of Production in 1985 and President of Aldrich in 1987. Late in 1999, he was promoted to President of Sigma-Aldrich Fine Chemicals, one of four strategic business units of Sigma-Aldrich Corp. During his tenure as president, Aldrich experienced a substantial growth in business. He oversaw a major expansion of the manufacturing plant (Pro II) at the Sheboygan County site, and worked tirelessly to integrate Aldrich more closely into Sigma-Aldrich Corp.

Clinton F. Lane, the current president, is a Purdue University graduate who was the first chemist hired for the Aldrich-Boranes, Inc. venture back in September 1972. After working as a bench chemist in Milwaukee, he moved, along with Aldrich-Boranes, Inc. to the Sheboygan County site soon after its purchase in 1977. After 13 years as Plant Manager of the Sheboygan County site, he was promoted to Vice President, Executive Vice President, and then President of Aldrich in 1999. Clint is credited with the substantial growth that both the Sheboygan County site and the line of boron-containing products have experienced. The interested reader should review Section 6.1 for more details.

7.2. The 80s and 90s

Following the merger of Aldrich and Sigma, business expectations for the merged company were soon realized. Annual double-digit growth has since been the norm, and other companies^{50,51}—Floyd Green's Aristo Custom Chemicals (USA, 1977), Makor Chemicals (Israel, 1978), Pathfinder (USA, 1984), Bio Yeda (Israel, 1986), Bristol Organics (UK, 1986), Fluka Chemie AG (Switzerland, 1989), Supelco (USA, 1993), LabKemi AB (Sweden, 1994), Research



David R. Harvey, Aldrich President (1981-1986).



Jai P. Nagarkatti, Aldrich President (1987-1999).



Clinton F. Lane, Aldrich President (1999-Present).

Biochemicals International (USA, 1997), Carbolabs (USA, 1997), Genosys Biotechnologies, Inc. (USA, 1998), Riedel-deHaën[®] Laborchemikalien GmbH & Co. KG (Germany, 1999), ARK Scientific GmbH

(Germany, 2000), First Medical, Inc. (USA, 2000), Amelung GmbH (Germany, 2000), Isotec, Inc. (USA, 2001)—have also become part of Sigma-Aldrich Corp. In 2000, Sigma-Aldrich corporate sales were over one billion dollars!⁵² The brands that make up Sigma-Aldrich Corp. are now well-known and trusted worldwide. The total number of products they offer is about 85,000 of which about 40,000 are produced.

Aldrich continued its spectacular growth in these two decades as evidenced by the purchase of several large buildings in the mid-1980s: 2905 W. Hope Avenue (currently holds the RCL collection), 1101 W. St. Paul Avenue (contains the Flavors & Fragrances products), and 1001 W. St. Paul Avenue (houses administrative offices and support departments). In 1986, Aldrich added 284 acres to the Sheboygan County site, and, in the early 90s, purchased and then added to a site at 6000 N. Teutonia Avenue in Milwaukee. In the 1990s, groupings of products, that Aldrich had offered since its early days, evolved into distinct product lines with their own technical managers, e.g., chiral, nonracemic products (~2,400 listings) and monomers & polymers (~3,100 products).

7.3. The Present and the Future 7.3.1. A Vision of the Future (2000)

Today, Aldrich continues to thrive within Sigma-Aldrich Corp. and has expanded into new market sectors—such as combinatorial chemistry, active high-purity metals and inorganics, and high-purity gases—and new overseas markets. In the past decade, it has also upgraded and enlarged its Milwaukee and East Coast distribution centers to enable the company to become even more responsive to its customers.

What does the future hold for Aldrich? Aldrich's future is intimately tied to that of Sigma-Aldrich. In December 2000, Sigma-Aldrich launched a new strategic plan clearly focusing the company on "leadership in Life Science and High Technology". In fact, 75% of the company's current sales are for Life Science applications, while the remaining 25% are in a variety of High Technology areas. The other key initiatives undertaken focus on service and process improvements.⁵³ The strategic plan made Aldrich a part of the Scientific Research business unit of Sigma-Aldrich Corp. The startup of a state-of-the-art, \$25 million production plant (Pro I) near Sheboygan Falls, WI, a combinatorial chemistry product line, a wide array of active and high-purity metals, and a strong Internet presence are but



Aldrich's Teutonia Avenue site.



Pro I building at Aldrich's Sheboygan County site.

a few of the ways in which Aldrich is implementing this strategy and is continuing to evolve.

Sigma-Aldrich's corporate vision is embodied in its motto: "We Are Committed to the Success of Our Customers, Employees and Shareholders through Leadership in Life Science, High Technology and Service."

Over the past five years, Sigma-Aldrich Corp. has been implementing its strategic plan for a strong, independent presence on the Internet with the goal of eventually transacting at least 50% of sales via the Internet. In 2000, corporate Internet sales represented over 10% of sales in the USA and about 5% of sales worldwide.⁵² Ongoing updates and a vast array of useful information, coupled with a focus on ease of use and visual appeal, have characterized the corporation's Web site.

Aldrich has vigorously participated in this effort by making available free of charge not only its catalog/handbook and various specialty catalogs and promotional materials,



Photo courtesy of H.C. Brown

Professor Herbert C. Brown (1999).



Photo courtesy of B.J. Horick

Beverly J. Horick (2001), recently retired Aldrich employee with the most years of service (February 17, 1956 to February 28, 2001).

but also its vast store of MSDSs, CofAs, and product technical data. One of its goals for the near future is to make available, free of charge, its IR, UV, and NMR spectral libraries on its Web site.

8. The Role of Science and Scientists

This is a topic that is so dear to the "heart" of the company that it merits a separate treatment.

8.1. Scientists' Contributions

Aldrich is a science-based company. It was cofounded by a Ph.D. chemist, currently employs several dozen Ph.D. chemists and many hundreds of college-educated chemistry professionals. Most of its customers are scientists of all walks of life. Early on, Aldrich recognized the importance of scientists to the growth and health of the company: Scientists were the originators of many of the successful Aldrich products, as well as the main consumers of its products. As discussed in Section 4.2, a large portion of the company's collection of hard-to-find research samples, known as the Rare Chemical Library, comes from the laboratories of these scientists.

Aldrich has been fortunate to have had long-standing professional collaborations with many of the leading chemists of the second half of the twentieth century. These relationships proved advantageous not only to the company, but also to the community of chemists by making available reagents that are now indispensable for chemistry research: $\text{Me}_3\text{S}\cdot\text{BH}_3$, NaBH_4CN , and the family of Selectride[®] reducing agents, to name a few. In the early days, it was natural for Alfred Bader to turn to Louis Fieser, his Ph.D. advisor at Harvard, or Martin Ettlinger, his graduate school contemporary, for ideas on what compounds would be of interest to researchers. Subsequently, Aldrich has had significant collaborations with many other leading chemists. What follows is only a partial list (in alphabetical order):⁵⁴ E. J. Corey, Henry Gilman (deceased), Eric N. Jacobsen, Kim D. Janda, Richard Lerner, Andrew G. Myers, K. C. Nicolaou, Martin J. O'Donnell, David O'Hagan, Siegfried Pickholz (deceased), Reuben Rieke, Ian P. Rothwell, I. Herbert Scheinberg, Barry K. Sharpless, John C. Sheehan (deceased), Gilbert Stork, and Robert B. Woodward (deceased).

Perhaps more than any other factor, it was the development of ideas and technologies, invented by researchers and developed or commercialized by Aldrich, that propelled Aldrich (and later Sigma-Aldrich) to the prominent position that it is presently in. It is unquestionably the vigorous pursuit of these contacts and collaborations that will keep Sigma-Aldrich a leading technology company.

8.2. Herbert C. Brown

Aside from Alfred Bader, perhaps no other single chemist has had a greater impact on the success of Aldrich than Professor Herbert C. Brown of Purdue University. A Nobel laureate and a towering figure in chemistry, Brown not only was the catalyst and a driving force for Aldrich-Boranes, Inc., as explained in Section 6.1, but he also served on the Aldrich Board of Directors (1972-1975) and the Sigma-Aldrich Board of Directors (1975-1979). In recognition of his lasting contributions, Aldrich not only pays royalties to the Purdue Research Foundation, but also co-sponsors the *Herbert C. Brown Award for Creative Research in Synthetic Methods* that is administered by the American Chemical Society.

8.3. Rewarding Excellence

Aldrich has also had a tradition of rewarding excellence in chemistry research by sponsoring or co-sponsoring prestigious professional awards, symposia, and student fellowships. A few examples come to mind: *ACS Award for Creative Work in Synthetic Organic Chemistry*, *ACS Award in Inorganic Chemistry*, *Herbert C. Brown Award for Creative Research in Synthetic Methods*, Project SEED, 32nd Organosilicon Symposium (Milwaukee, 1999), Asymmetric Synthesis Symposium (Milwaukee, 1998), Boron-USA meetings, and various Gordon Conferences. Moreover, Alfred Bader, Aldrich's cofounder, personally sponsors the *Alfred Bader Award in Bioinorganic or Bioorganic Chemistry*, and has helped over the years many deserving academic chemists by underwriting some of their research.²

9. Valued Customers, Dedicated Employees

Aldrich recognized very early on the importance of establishing strong relationships with its customers, and the necessity to provide them with valuable information related to the products that they were purchasing. Thus, it was no surprise that the Aldrich catalog was transformed from a listing of available products and prices to a "handbook" containing a wealth of information, which has made it an indispensable desk reference in many academic and industrial laboratories and libraries. Chemists active in research were not only customers, but also partners in the chemical enterprise, who were also invited to share their insights with others through such widely circulated and free Aldrich publications as the *Aldrichimica Acta* (which has been in existence for 34 years). It is also no secret that ideas for new products often

Table 2. Present and Former Aldrich Employees with 25 or More Years of Continuous Service^{ab}

Adler, Wayne J.	Gorzek, Robert J.	Leitner, Lorraine	Rochwerger, Leonard L.
Ahmed, Waheeduddin	Griesinger, Alfred	Lenga, Robert E.	Roper, Mattie D.
Bader, Alfred R.	Griffiths, David W.	Lent, Mary A.	Saladin, Barbara L.
Benson, Christine F.	Gunther, Patricia A.	Lewis, Robert J.	Schreiber, Peter L.
Borenstein, Mark	Harvey, David R.	Lisztwan, Emilia M.	Settingsgaard, Jacqueline L.
Bourgeois, Shirley R.	Helmin, William T.	Malone, Rosie L.	Shortridge, Nelgene
Branski, Robert A.	Holm, Phillip L.	Mehta, Milan N.	Shuder, Diane L.
Brien, Diana M.	Horick, Beverly J.	Metz, Marian E.	Siegel, Brian S.
Brien, James J.	Jenkins, Dolores H.	Mititch, Jacqueline	Skeff, George
Bruesewitz, Richard J.	Kasprzak, Russell J.	Nagarkatti, Jai P.	Smith, Andrew P.
Creighton, Anthony J.	Kett, Jeffrey A.	Napiorkowski, Anna M.	Smith, Robert W.
Daniels, John J.	Koppel, Henry C.	Podd, Rodney L.	Stanton, Genevieve L.
Edelstein, Sara	Kopperud, Cynthia A.	Poth, Donna L.	Wallace, Kenneth J.
Farrell, Richard T.	Korthoff, Kristine L.	Pouchert, Charles J.	Ward, Stella L.
Feustel, Barbara L.	Kratzer, Phyllis C.	Pruss, Judith R.	Weber, Roger O.
Fox, Lyle G.	Kreinus, Timothy M.	Pykett, Jonathan R.	Wells, Sheila E.
Freeman, Roland P.	Kurzynski, Alice J.	Rebarchik, Joseph A.	Wickersham, Thomas W.
Gallaspy, Barbara A.	Lane, Clinton F.	Riedmaier, John E.	Wondra, Carl T.

^a As of June 2001. ^b Admittedly, length of service is an imperfect measure of an employee's contribution to the company; however, this author was at a loss to come up with a fair, objective way of recognizing those employees who have given so much to the company, but may not be mentioned in the text. I offer my apologies to those employees, whose names belong in the table, but were inadvertently left out.

came from customers. To appreciate this fact more fully, it is sufficient to consult the "Please Bother Us." section of any recent *Aldrichimica Acta* issue.

Aldrich's phenomenal success is a tribute to the vision and determination of its cofounder, Dr. Alfred R. Bader, and the dedication and hard work of thousands of former and present employees. **Table 2** is only a modest attempt at acknowledging their contributions.

10. Acknowledgments

In addition to the sources cited, I wish to acknowledge the specific assistance of (in alphabetical order): (i) Bettie Aldrich Eisendrath, Gerd Backes, Alfred Bader, Jim Brien, Tom Gandia, David Harvey, Don Hobbs, Harvey Hopps, Beverly Horick, Peter Hyland, Linda Kehren, Clinton Lane, Lorraine Leitner, Edward Niemiec, Judith Pruss, Robert Smith, Joan Suda, Robert Wandler, Tom Wickersham, and LaShannon Wilson, who provided me with valuable information and recollections either via personal communications or via source materials that they supplied to me; (ii) Brian Case, Robert Gorzek, David Harvey, Chris Hewitt, Harvey Hopps, Peter Hyland, Anthony J. La Loggia, Clinton Lane, Jai Nagarkatti, and Craig Recatto, who proofread the manuscript or sections thereof, and offered helpful comments; (iii) Jennifer L. Botic, who laid out the manuscript and this

issue of the *Acta*, and who assisted me in locating some of the photographs used; and (iv) Rebecca Zelenka, who helped with contacting former Aldrich employees. Finally, I would like to thank my wife for her patience and understanding while I was preparing the manuscript.

11. References and Notes

- (1) This brief tour of the past fifty years cannot possibly do justice to the topic. For more details and anecdotes, the interested reader is directed to the very readable book by one of the founders of Aldrich: Bader, A. *Adventures of a Chemist Collector*; Weidenfeld and Nicolson: London, U.K., 1995.
- (2) Cori, T.; Emanuel, R. N.; Harvey, D.; Klitsner, M. E. *Aldrichimica Acta* **1984**, *17*, 3.
- (3) (a) Bader, Alfred. The Building of Aldrich. My Advice to Entrepreneurs. *The Chemist*, November/December 1997, pp 1-5. (b) Buchan, P. Bruce. Three Boards and "A Bet Against the Company". *The Chemical Intelligencer*, October 1996, pp 24-29 and 41. (c) Edward, J. T. *Can. Chem. News* **1992**, *44*(6), 23. (d) Bohning, James J. Crystallizing Hamburger: Alfred Bader and the Aldrich Chemical Company, Part I. *Beckman Center News*, Spring 1991, pp 1 and 8-9. (e) Bohning, James J. Crystallizing Hamburger: Alfred Bader and the Aldrich Chemical Company, Part II. *Beckman Center News*, Fall 1991, pp 3-4. (f) Bader, A. R. *CHEMTECH* **1990** (March), 138. (g) A Chemical Company in Your Garage. An Interview with Dr. Alfred Bader. *The DEL-CHEM BULLETIN*, May 1974, pp 5-10.
- (4) While it is unquestionably the contributions of a great many dedicated employees that have made Aldrich what it is today, it is not possible in such a short overview to mention them all. The author regrets any inadvertent or necessary omissions.
- (5) Reference 1, pp 70, 101, and 185.
- (6) (a) The 500 shares of stock issued were owned 50% by Alfred R. Bader and 50% by Jack N., Frank N., and Bettie Mae Eisendrath. Professor A. F. McKay of the University of Toronto, Canada, owned only one share of stock: Bader, A. R. Alfred Bader Fine Arts, Milwaukee, WI. Personal communication, April 09, 2001. (b) Prior to cofounding Aldrich, Jack Eisendrath had attempted to start and run a number of mail order/catalog businesses (e.g., selling moccasins): Eisendrath, B. A. Washington, DC. Personal communication, April 03, 2001.
- (7) (a) Reference 1, p 70. (b) Up until the Eisendraths sold their 50% stake in the company, Bettie acted as the company's (unpaid) secretary. As of the writing of this review, Bettie was still a remarkably energetic and socially active octogenarian living in Washington, DC.
- (8) (a) Leitner, L. (née Neau; retired) Aldrich Chemical Co., Inc., Milwaukee, WI. Personal communication, March 16, 2001. (b) Horick, B. J. (retired) Aldrich Chemical Co., Inc., Milwaukee, WI. Personal communication, February 01, 2001.
- (9) (a) McKay, A. F.; Wright, G. F. *J. Am. Chem. Soc.* **1947**, *69*, 3028. (b) McKay, A. F. *J. Am. Chem. Soc.* **1948**, *70*, 1974. (c) McKay, A. F.; Ott, W. L.; Taylor, G. W.; Buchanan, M. N.; Crooker, J. F. *Can. J. Res., Sec. B* **1950**, *28*, 683.

- (10) Aldrich Chemical Co., Inc. *Diazald*[®], *MNNG*, and *Diazomethane Generators*. Aldrich Technical Information Bulletin No. AL-180; Milwaukee, WI, 1993.
- (11) Reference 1, pp 71–72.
- (12) Bader, A. R. Alfred Bader Fine Arts, Milwaukee, WI. Personal communication, March 19, 2001.
- (13) While at PPG, Anthony co-authored a paper with Alfred Bader on the easy preparation of phenyl esters of carboxylic acids by heating the acid and phenol in the presence of PPA: Bader, A. R.; Kontowicz, A. D. *J. Am. Chem. Soc.* **1953**, *75*, 5416.
- (14) Reference 1, p 74.
- (15) Jack Eisendrath continued his practice of general, family, and consumer law in Milwaukee for about 40 years. He died on November 06, 1997 at the age of 85: Knoche, Eldon. Consumers Had Advocate in Attorney Eisendrath. *The Milwaukee Journal Sentinel* [Online], November 9, 1997, main page (www.jsonline.com).
- (16) For a glimpse of Eisendrath's perspective of these events, see: (a) Lank, Avrum D. Chemist Mixes Knowledge, Savvy. *The Milwaukee Sentinel*, October 01, 1985, Part 4, pp 1–2. (b) Gillespie, Scott. Alfred Bader's Diverse Talents Led Firm to International Role. *The Business Journal Special Report* (Milwaukee), Week of September 22, 1986, p 9.
- (17) Bader, A. R. Alfred Bader Fine Arts, Milwaukee, WI. Personal communication, March 09 and April 09, 2001.
- (18) Kenney, Ray. Aldrich Firing Up All Burners. *The Milwaukee Sentinel*, September 02, 1974, Part 2, p 9.
- (19) Reference 1, p 77.
- (20) Bader, A. R. Alfred Bader Fine Arts, Milwaukee, WI. Personal communication, March 23, 2001.
- (21) Helen Bader (née Daniels) worked for Aldrich in various capacities on and off for over twenty years, whenever her family obligations permitted her to. In addition to being one of the owners of Aldrich, she was also a company director and treasurer for several years and the third company president for a very short period of time (1964). Following her death in 1989, her family established the Helen Bader Foundation, Inc. to honor her memory and continue the charitable work that she had started. The interested reader can find out more by accessing the Foundation's Web site at www.hbf.org.
- (22) Reference 1, pp 80–81.
- (23) Horick, Beverly J. Then & Now. *The Aldrich Reporter*, January 2001, p 10. When Beverly J. Horick retired from Aldrich on February 28, 2001 (after a little over 45 years of continuous service!) she set a record as the employee with the longest service to the company.
- (24) Hopps, H. B. Amarillo College, Amarillo, TX. Personal communication, February 21, 2001.
- (25) (a) Aldrich Chemical Co., Inc. The Aldrich Annual Report; Milwaukee, WI, September 03, 1969. (b) Aldrich Chemical Co., Inc. Annual Report; Milwaukee, WI, September 30, 1970.
- (26) Reference 1, p 98.
- (27) Aldrich Chemical Co., Inc. Prospectus; Milwaukee, WI, December 30, 1965; p 2.
- (28) Reference 1, p 79.
- (29) EGA from the names Ernst, Gerhard, and Alfred.
- (30) Aldrich Chemical Co., Inc. Annual Report 1971; Milwaukee, WI, March 24, 1972.
- (31) (a) Mack, H. Sigma-Aldrich Chemie GmbH, Steinheim, Germany. Personal communication, March 06, 2001. (b) Backes, G. Sigma-Aldrich Chemie GmbH, Steinheim, Germany. Personal communication, March 09, 2001.
- (32) Koppel, H. *Aldrichimica Acta* **1968**, *1*, 3.
- (33) Reference 1, p 88.
- (34) Medical Economics Company, Inc. *Physicians' Desk Reference*, 52nd ed.; Montvale, NJ, 1998; pp 2082–2086.
- (35) The statements in this paragraph should not be construed to imply that COGNEX[®] is formulated with a material obtained from Aldrich. Aldrich brand products are sold mainly for research or industrial applications and are not intended for drug or household use, unless specifically designated for that purpose.
- (36) Reference 1, pp 82–83.
- (37) Aldrich Chemical Co., Inc. *Aldrichimica Acta* (Preview Issue); Milwaukee, WI, Fall 1967.
- (38) Reference 1, p 133.
- (39) Lane, Clinton F. Greetings from the Future Leaders in Organoborane Chemicals! *Aldrich-Boranes, Inc. Newsletter*, February 1973, p 2.
- (40) For some examples, see: (a) Lane, C. F. *Aldrichimica Acta* **1973**, *6*, 21. (b) Lane, C. F. *Aldrichimica Acta* **1973**, *6*, 51. (c) Lane, C. F. *Aldrichimica Acta* **1974**, *7*, 7. (d) Lane, C. F. *Aldrichimica Acta* **1974**, *7*, 32. (e) Lane, C. F. *Synthesis* **1975**, *135*. (f) Lane, C. F. *Aldrichimica Acta* **1975**, *8*, 3. (g) Lane, C. F. *Aldrichimica Acta* **1975**, *8*, 20. (h) Lane, C. F.; Kabalka, G. W. *Tetrahedron* **1976**, *32*, 981. (i) Lane, C. F. *Chem. Rev.* **1976**, *76*, 773. (j) Lane, C. F. *Aldrichimica Acta* **1976**, *9*, 31. (k) Lane, C. F.; Kramer, G. W. *Aldrichimica Acta* **1977**, *10*, 11. (l) Brown, H. C.; Lane, C. F. *Heterocycles* **1977**, *7*, 453. (m) Lane, C. F. *Aldrichimica Acta* **1977**, *10*, 41.
- (41) Aldrich Chemical Co., Inc. *Welcome to Aldrich Sheboygan Ste.* Aldrich Informational Bulletin; Sheboygan Falls, WI, 1999.
- (42) Aldrich Chemical Co., Inc. The Aldrich Semimonthly Newsletter (Extra); Milwaukee, WI, April 25, 1974.
- (43) (a) Aldrich Chemical Co., Inc. Annual Report 1973; Milwaukee, WI, February 15, 1974. (b) Reference 1, p 167.
- (44) Aldrich Chemical Co., Inc. Annual Report 1974; Milwaukee, WI, February 21, 1975.
- (45) Green, F. J. *The Sigma-Aldrich Handbook of Stains, Dyes and Indicators*; Aldrich Chemical Co., Inc.: Milwaukee, WI, 1990.
- (46) Has since been replaced by the third edition, cat. no. Z10.750-6: Pouchert, C. J. *The Aldrich Library of Infrared Spectra*, 3rd ed.; Aldrich Chemical Co., Inc.: Milwaukee, WI, 1981.
- (47) Brien, J. J. Aldrich Chemical Co., Inc., Milwaukee, WI. Personal communication, March 08, 2001.
- (48) Damrow, D. (retired) Aldrich Chemical Co., Inc., Milwaukee, WI. Personal communication, February 22, 2001.
- (49) For the record, Aldrich has had the following presidents in the following order: Jack Eisendrath, Alfred Bader, Helen Bader, Alfred Bader, David Harvey, Jai Nagarkatti, and Clinton Lane.
- (50) Suda, Joan. *The Sigma-Aldrich Story. Overview of a Multinational Family of Companies*. Sigma-Aldrich Bulletin; Sigma-Aldrich Corp.: St. Louis, MO, 2001.
- (51) Sigma-Aldrich Corp. *Sigma-Aldrich Acquires Isotec to Enhance Life Science Capabilities*. Press Release; St. Louis, MO, February 16, 2001.
- (52) Suda, Joan. *An Interview with David R. Harvey*. Sigma-Aldrich Bulletin; Sigma-Aldrich Corp.: St. Louis, MO, January 2001; p 1.
- (53) Harvey, D. R. Sigma-Aldrich Corp., St. Louis, MO. Personal communication, May 05, 2001.
- (54) Many more scientists have had fruitful collaborations or a significant impact on the growth of Aldrich. The author regrets not being able to acknowledge all of them.

Monsanto Company: Roundup[®]; **Riedel-deHaën GmbH:** Riedel-deHaën[®]; **Rieke Metals:** Rieke[®]; **Sigma-Aldrich Co.:** AAPL[™], ALDRICH[®], Diazald[®], FLUKA[®], Selectride[®], SIGMA[®]; **Warner-Lambert Co.:** COGNEX[®].

Happy 50th Anniversary Aldrich!

About the Author

Sharbil J. Firsan was born and raised in Lebanon. He completed his undergraduate studies at the American University of Beirut and his graduate work on acyclic imidate and thioimidate *N*-oxides with Professor Robert M. Coates at the University of Illinois in Urbana-Champaign (Ph.D., 1986). He did postdoctoral work at the University of Oregon in Eugene, OR, and then moved to Oklahoma State University in Stillwater, OK, to become a Research Associate and then a Visiting Assistant Professor. In 1996, he joined Aldrich Chemical Co., Milwaukee, WI, as a Promotions and Publications Specialist. He is currently a Senior Promotions and Publications Specialist and Editor of the *Aldrichimica Acta*. With his wife, Leah (Leila), Sharbil enjoys outdoor activities, gardening, and travel.



Sigma-Aldrich Worldwide Locations

Argentina

Sigma-Aldrich de Argentina, S.A.
Tel: 54 11 4807 0321
Fax: 54 11 4807 0346
E-mail: info@sigma-aldrich.com.ar

Australia

Sigma-Aldrich Pty, Limited
Free Tel: 1-800-800-097
Free Fax: 1-800-800-096
Tel: 61-2-9841-0555
Fax: 61-2-9841-0500
E-mail: ausmail@sial.com

Austria

Sigma-Aldrich Handels GmbH
Tel: 43 1 605 8110
Fax: 43 1 605 8120
E-mail: sigma@sial.co.at

Belgium

Sigma-Aldrich nv/sa
Free Tel: 0800-147.47
Free Fax: 0800-147.45
Tel: 03 899.13.01
Fax: 03 899.13.11
E-mail: sigmaldr@sial.be

Brazil

Sigma-Aldrich Quimica Brasil Ltda.
Tel: 55-11-231-1866
Fax: 55-11-257-9079
E-mail: sigmabr@ibm.net

Canada

Sigma-Aldrich Canada Ltd.
Free Tel: 800-565-1400
Free Fax: 800-265-3858
Tel: 905-829-9500
Fax: 905-829-9292
E-mail: canada@sial.com

Czech Republic

Sigma-Aldrich s.r.o.
Tel: 00-420-2-2176 1300
Fax: 00-420-2-2176 3300
E-mail: CZEcustSV@eurnotes.sial.com

Denmark

Sigma-Aldrich Denmark A/S
Tel: +45 43 56 59 10
Fax: +45 43 56 59 05
E-mail: denorder@eurnotes.sial.com

Eire

Sigma-Aldrich Ireland Ltd.
Free Tel: 1 800 200 888
Free Fax: 1 800 600 222
Tel: 01 404 1900
Fax: 01 404 1910
E-mail: EICustsv@eurnotes.sial.com

Finland

Sigma-Aldrich Finland
Tel: (09) 350-9250
Fax: (09) 350-92555
E-mail: finorder@eurnotes.sial.com

France

Sigma-Aldrich Chimie S.a.r.l
Tel Numero Vert: 0800 21 14 08
Fax Numero Vert: 0800 03 10 52
E-mail: fradvsv@eurnotes.sial.com

Germany

(also SE Europe, the Baltics, Africa and the Middle East)

Sigma-Aldrich Chemie GmbH
FreeCall Tel: 0800-5155 000
FreeCall Fax: 0800-6490 000
Tel: +49(0)89-6513-0
Fax: +49(0)89-6513-1169
E-mail: deorders@eurnotes.sial.com
For Slovakia, Africa, Middle East and other European countries
Tel: +49(0)89-6513-1807
Fax: +49(0)89-6513-1877

Greece

Sigma-Aldrich (o.m.) Ltd
Tel: 30-1-9948010
Fax: 30-1-9948331
E-mail: sigald@acci.gr

Hungary

Sigma-Aldrich Kft
Tel: 06-1-269-6474
Fax: 06-1-235-9050
E-mail: info@sigma.sial.hu

India

Sigma-Aldrich Corporation
Bangalore location
Phone: 91 80 851 8797
Fax: 91 80 851 8358

New Delhi location:

Phone: (011) 689 9826
Fax: (011) 689 9827
E-mail: sigma@del2.vsnl.net.in

Israel

Sigma Israel Chemicals Ltd
Toll Free Tel: 1-800-70-2222
Tel: 972-8-948-4222
Fax: 972-8-948-4200
E-mail: sigisr@sigma.co.il

Italy

Sigma-Aldrich S.r.l
Numero Verde: 800-827018
Tel: 02-33417.310
Fax: 02-38010.737
E-mail: itorder@eurnotes.sial.com

Japan

Sigma-Aldrich Japan K.K.
Tel: 81-3-5821-3111
Fax: 81-3-5821-3170
E-mail: sialjp@sial.com

Korea

Sigma-Aldrich Korea
Toll Free Tel: 080-023-7111
Toll Free Fax: 080-023-8111
Tel: 82-2-783-5211
Fax: 82-2-783-5011
E-mail: order@safkor.co.kr

Malaysia

Sigma-Aldrich (M) Sdn Bhd
Tel: 603-56353321
Fax: 603-56354116
E-mail: sam@sial.com

Mexico

Sigma-Aldrich Química, S.A. de C.V.
Tel Gratis: 91-800-0075300
Tel: (72) 76 1600
Fax: (72) 76 1601
E-mail: mexico@sial.com

Netherlands

Sigma-Aldrich Chemie BV
Tel Gratis: 0800 022 9088
Fax Gratis: 0800 022 9089
Tel: 078-620 54 11
Fax: 078-620 54 21
E-mail: nlcustsv@eurnotes.sial.com

Norway

Sigma-Aldrich Norway AS
Tel: 47 23 176 000
Fax: 47 23 176 010
E-mail: nororder@eurnotes.sial.com

Poland

Sigma-Aldrich Sp.z o.o.
Tel: (48-61) 823-24-81
Fax: (48-61) 823-27-81
E-mail: plcustsv@eurnotes.sial.com

Portugal

Sigma-Aldrich Quimica S.A.
Free Tel: 0800-202180
Free Fax: 0800-202178
E-mail, encomendas:
poorders@eurnotes.sial.com

Russia

Sigma-Aldrich Russia
Tel: 7-095-975-3321
Fax: 7-095-975-4792
E-mail: techcare@online.ru

Singapore

Sigma-Aldrich Pte., Ltd
Tel: (65) 271 1089
Fax: (65) 271 1571
E-mail: sup@sial.com

South Africa

Sigma-Aldrich South Africa (Pty) Ltd
Tel: (011) 397 8886
Fax: (011) 397 8859
E-mail: sigald@icon.co.za

Spain

Sigma-Aldrich Química S.A.
Free Tel: 900-101376
Free Fax: 900-102028
Tel: 91-6619977
Fax: 91-6619642
E-mail, pedidos:
esorders@eurnotes.sial.com

Sweden

Sigma-Aldrich Sweden AB
Tel: 020-35 05 10
Fax: 020-35 25 22
Outside Sweden Tel: 46-8-742 4200
Outside Sweden Fax: 46-8-742 4243
E-mail: sweorder@eurnotes.sial.com

Switzerland

Aldrich Chemie
A division of Fluka Chemie AG
Tel: 081-755 27 23
Fax: 081-755 28 40
Swiss Free Call: 0800 80 00 80
E-mail: Fluka@sial.com

United Kingdom

Sigma-Aldrich Company Ltd.
Free Tel: 0800 717181
Free Fax: 0800 378538
Tel: 01202 733114
Fax: 01202 715460
E-mail: ukcustsv@eurnotes.sial.com

United States

Aldrich Chemical Company, Inc
Toll Free Tel: 1-800-558-9160
Toll Free Fax: 1-800-962-9591
Tel: 414-273-3850
Fax: 414-273-4979
E-mail: aldrich@sial.com

Sigma Chemical Company
Toll Free Tel: 1-800-521-8956
Toll Free Fax: 1-800-325-5052
Tel: 314-771-5765
Fax: 314-771-5757
E-mail: sigma@sial.com


Aldrich Vacuum Manifolds with Joints

This practical design provides clearance to accommodate 250mL flasks with the snap of a KECK® clip! Manifolds are available with either 14/20 or 24/40 joints. Glass stopcocks have a 4mm bore.

KECK is a registered trademark of Cole-Parmer Instrument Co.

Positions	Overall L (mm)	14/20 Cat. No.	24/40 Cat. No.	Each
3	300	Z53,066-2	Z53,069-7	\$435.00
4	400	Z53,067-0	Z53,070-0	525.00
5	500	Z53,068-9	Z53,071-9	595.00

Do you have an innovative shortcut or unique laboratory hint you'd like to share with your fellow chemists? If so, please send it to Aldrich (attn: Lab Notes, Aldrichimica Acta). For submitting your idea, you will receive a complimentary, Aldrich periodic table poster (Cat. No. **Z54,320-9**). If we publish your lab note, you will also receive an Aldrich periodic table mouse pad (Cat. No. **Z54,323-3**). It is Teflon®-coated, 8½ in. x 11 in., with a full-color periodic table on the front. We reserve the right to retain all entries for future consideration.



Teflon is a registered trademark of E.I. du Pont de Nemours & Co., Inc.

"A Member of the Sigma-Aldrich Family"

In 2001, Aldrich Chemical Company, Inc. is celebrating the 50th anniversary of its founding on August 17, 1951. As most of our customers and *Aldrichimica Acta* readers are aware of, Aldrich has been a part of Sigma-Aldrich Corporation since August 1975.

Sigma-Aldrich Corporation has built its reputation as a leader in life science and high technology on the strengths of its many outstanding member companies. Companies that are leaders in their own markets, and that have widely known and trusted names such as Aldrich, Fluka, Isotec, Riedel-deHaën, Sigma, and SUPELCO.

Outstanding customer and technical services, the largest selection of top-notch products, a wide range of manufacturing and sourcing capabilities, and a strong Internet presence are some of what Sigma-Aldrich Corporation offers its customers.

Whether you need a few milligrams of material for research or a few tons for manufacturing, Sigma-Aldrich Corporation has the infrastructure in place and the skilled people to meet your needs and exceed your expectations. Do you have a custom synthesis project or would like to scale up a process? Just contact our fine chemicals division. Sigma-Aldrich Fine Chemicals (SAFC) provides large-scale organic and biochemical manufacturing, custom manufacturing, and sourcing services. It can help you from development through commercialization with over 50,000 gallons of total reactor capacity worldwide, an extensive technology toolbox, and several cGMP and ISO compliant plants.

We, at Aldrich, would like to say that it has been our privilege to serve you, our customers, for the past fifty years. We also look forward to providing you with the same top-quality service in the next fifty years.

Riedel-deHaën® is a registered trademark under license from Riedel-Gesellschaft GmbH

50 YEARS

1951
to
2001

SIGMA-ALDRICH
CORPORATION

ALDRICH CHEMICAL COMPANY, INC.
P.O. BOX 14508
ST. LOUIS, MISSOURI 63178 USA



CHANGE SERVICE REQUESTED

PRESORTED STANDARD
U.S. POSTAGE
PAID
MILWAUKEE, WISCONSIN
PERMIT NO. 552

014

1951–2001: FIFTY YEARS OF CHEMISTS HELPING CHEMISTS

Aldrichimica ACTA

VOL. 34, NO. 3 • 2001

(Last issue in 2001)




Solid-Phase
Dendrimer Chemistry

[Hydroxy(tosyloxy)-
iodo]benzene



 **ALDRICH**[®]
CHEMISTS HELPING CHEMISTS


54,167-2 4-(*tert*-Butyldimethylsilyloxy)-1-butyne, 97%
5mL \$22.50; 25mL \$75.00



Has been utilized for a transition-metal-catalyzed regioselective cyclotrimerization with enones¹ and as a starting material in the synthesis of (±)-asteriscanolide.²

(1) Ikeda, S.-I. et al. *Chem. Commun.* **2000**, 815. (2) Krafft, M. E. et al. *Synthesis* **2000**, 1020.


55,628-9 Methyl picolinate, 99%
25g \$63.00; 100g \$175.00



Useful as a ligand for the formation of metal complexes¹ and as a starting material for the synthesis of biologically active aza-anthraquinones.²

(1) Britovsek, G. et al. *J. Mol. Catal. A: Chem.* **1996**, *110*, 77. (2) Epszajn, J. et al. *Tetrahedron* **1996**, *52*, 11025.


55,406-5 Ethyl α-bromophenylacetate, 97%
25g \$78.00; 100g \$235.00



This versatile alkylating reagent has been utilized in the preparation of multifunctional phosphonates,¹ oxazinones,² and α-alkoxycarbonyl compounds.³

(1) Prager, R. H. et al. *Aust. J. Chem.* **1997**, *50*, 813. (2) Remuzon, P. et al. *Tetrahedron* **1997**, *53*, 17711. (3) Berglund, P. et al. *Tetrahedron: Asymmetry* **1999**, *10*, 4191.


55,540-1 Monomethyl isophthalate, 97%
1g \$21.90; 5g \$72.90



Was recently used to introduce structural variations in benzazepines,¹ and to synthesize a monomer for the preparation of an authentic H-T poly(amide-ester).²

(1) Murakami, Y. et al. *J. Med. Chem.* **1999**, *42*, 2621. (2) Li, L. et al. *Macromolecules* **1999**, *32*, 3851.


51,810-7 *tert*-Butyl 1-indolecarboxylate, 97%
25mL \$59.00



Allows the easy preparation of indole-2-boronic acid¹ and the corresponding 3-nitroindole.²

(1) de Koning, C. B. et al. *J. Chem. Soc., Perkin Trans. 1* **2000**, 1705. (2) Pelkey, E. T.; Gribble, G. W. *Synthesis* **1999**, 1117.


54,735-2 Mono-*tert*-butyl succinate, 97%
1g \$12.75; 5g \$42.45



Key building block for the construction of a chiral lactone that is a versatile intermediate for possible HIV-1 protease inhibitors,¹ and for the construction of the carbacephem β-lactam framework.²

(1) Solladié-Cavallo, A. et al. *Tetrahedron: Asymmetry* **1996**, *7*, 1797. (2) Guzzo, P.R.; Miller, M. J. *J. Org. Chem.* **1994**, *59*, 4862.


52,117-5 4-Ethynylbiphenyl, 97% **5g \$68.00**



This terminal acetylene has been converted to enamines with resin-bound 2° amines.¹ It is also suitable for Pd(0)-mediated coupling on solid supports.²

(1) Aznar, F. et al. *Tetrahedron Lett.* **2000**, *41*, 5683. (2) Berteina, S. et al. *Synlett* **1998**, 676.


51,296-6 Z-L-prolinol, 97% **1g \$48.80**



Useful for synthesizing peptidomimetics. A relatively recent application is in the synthesis of an ester mimic for the generation of catalytic antibodies.

Anderson, G. T. et al. *J. Org. Chem.* **1996**, *61*, 125.

55,415-4 4-*tert*-Butyl-2,6-diformylphenol, 96%
1g \$15.00; 5g \$48.90



Reported applications for this diformylphenol include the synthesis of binucleating macrocyclic ligands¹ and homooxocalix[n]arenes.²

(1) Lindoy, L. F. et al. *Synthesis* **1998**, 1029. (2) Komatsu, N. *Tetrahedron Lett.* **2001**, *42*, 1733.


55,587-8 Chloromethyl(dimethyl)silane, 97%
5mL \$47.05; 25mL \$156.75




Hydrosilylation of olefins using this reagent provides the corresponding chloromethyl(dimethyl)silyl (CMDMS) derivatives.^{1,2} This method was employed to introduce CMDMS groups on the terminal branches of carbosilane dendrimers.²

(1) Thibon, J. et al. *J. Org. Chem.* **1997**, *62*, 4635. (2) Krška, S. W.; Seyferth, D. *J. Am. Chem. Soc.* **1998**, *120*, 3604.

55,360-3 2,3-Dibromomaleimide, 97%
1g \$12.00; 5g \$40.00



55,778-1 *N*-Benzyl-2,3-dibromomaleimide, 97%
1g \$13.50; 5g \$45.00



Dihalogenated maleimides have applications in a broad spectrum of synthetic transformations, such as mono- or disubstitutions, Suzuki-type couplings, Grignard reactions, and Diels-Alder cyclizations. They are key intermediates in the synthesis of arcyriarubin,¹ staurosporine,² and pyrroloquinoxaline systems.³

(1) Mahboobi, S. et al. *Pharmazie* **1999**, *54*, 820. (2) Joyce, R. P. et al. *J. Org. Chem.* **1987**, *52*, 1177. (3) Hanaineh-Abdelnour, L. et al. *Tetrahedron* **1999**, *55*, 11859.

Aldrichimica ACTA

VOL. 34, NO. 3 • 2001

Aldrich Chemical Co., Inc.

1001 West Saint Paul Ave.
Milwaukee, WI 53233 USA

To Place Orders

Telephone 800-558-9160 (USA)

or 414-273-3850

FAX 800-962-9591 (USA)

or 414-273-4979

Mail P.O. Box 2060

Milwaukee, WI 53201 USA

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.

P.O. Box 355, Milwaukee, WI 53201 USA

Customer & Technical Services

Customer Inquiries 800-558-9160

Technical Service 800-231-8327

Sigma-Aldrich Fine Chemicals 800-336-9719

Custom Synthesis 800-336-9719

Flavors & Fragrances 800-227-4563

International 414-273-3850

24-Hour Emergency 414-273-3850

Web Site <http://www.sigma-aldrich.com>

E-Mail aldrich@sial.com

To request your **FREE** subscription to the *Aldrichimica Acta*,

please call: **800-558-9160** (USA),

write: **Attn: Mailroom**
Aldrich Chemical Co., Inc.
P.O. Box 355
Milwaukee, WI 53201-9358

or email: **sams-usa@sial.com**

International customers, please contact your local Sigma-Aldrich office.

The *Aldrichimica Acta* is also available on the Internet at <http://www.sigma-aldrich.com>.

All prices listed in this publication are subject to change without notice.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

For worldwide contact information, please see the inside back cover.

Aldrichimica Acta is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich family. © 2001 by Sigma-Aldrich Co. Printed in the United States.

About Our Cover

Edouard Manet's *Masked Ball at the Opera* (oil on canvas, 23¼ in. x 28½ in.), signed and dated 1873, is a reminder of one of the many sophisticated entertainment forms available to Parisians in the third quarter of the nineteenth century: the balls held at the opera on Saturday nights in March during Lent. The foyer of the old Opera House in the Rue Le Peletier is shown crowded with figures. The flowing line of the




silk hats of a large group of elegantly dressed gentlemen creates a graceful rhythm across the center of the scene, helping to unify an otherwise somewhat fragmented composition. The informality and great variety of poses of the figures and the seemingly random cropping on either side and across the top of the image give a sense of immediacy and realism to the scene. The figures of Polichinelle on the left, of the gentleman on the right edge of the canvas, and of those on the balcony are all cut off by the frame, suggesting that only a small part of what is going on is visible in the picture. Manet included portraits of himself and several of his friends in the painting, which is remarkable for the absence of color. Manet saw the men's silk hats and evening clothes as a unique opportunity to depict every variation of ebony.

This painting is a gift of Mrs. Horace Havemeyer to the National Gallery of Art in Washington, D.C.

“Please
Bother
Us.”

Clint Lane

Clint Lane, President



55,692-0

55,693-9

Professor Steven V. Ley of the University of Cambridge, U.K., kindly suggested that we offer these enantiomeric 2,3-butanediol-derived dimethyl tartrates. These compounds have been used in the synthesis of protected chiral tetrol derivatives¹ and for the transfer of chirality in enantioselective reactions.²

(1) Dixon, D. J.; Foster, A. C.; Ley, S. V.; Reynolds, D. J. *J. Chem. Soc., Perkin Trans. 1* **1999**, 1635. (2) Barlow, J. S.; Dixon, D. J.; Foster, A. C.; Ley, S. V.; Reynolds, D. J. *J. Chem. Soc., Perkin Trans. 1* **1999**, 1627.

55,692-0

(2*R*,3*R*,5*R*,6*R*)-5,6-Dimethoxy-5,6-dimethyl-[1,4]dioxane-2,3-dicarboxylic acid dimethyl ester, 98%
25g \$56.85

55,693-9

(2*S*,3*S*,5*S*,6*S*)-5,6-Dimethoxy-5,6-dimethyl-[1,4]dioxane-2,3-dicarboxylic acid dimethyl ester
1g \$10.45; 5g \$34.80

Naturally, we made these useful reagents. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page or on the inside back cover.

Table of Contents

Solid-Phase Dendrimer Chemistry:
Synthesis and Applications75
—S. Lebreton, S. Monaghan,
and M. Bradley*

[Hydroxy(tosyloxy)iodo]benzene and
Closely Related Iodanes: The Second
Stage of Development89
—G. F. Koser



Lab Notes

Rapid-Dry TLC Spotter

Flash chromatography can very quickly generate a lot of samples that may require spotting on a TLC plate for analysis. If you want to get small spots at the origin you have to add the solvent slowly while blowing on the spot to evaporate the solvent. This can lead to hyperventilation or solvent inhalation and, depending on the number of fractions, can take longer than running the column.





Figure 1

To simplify this operation and speed up the process, I devised a simple apparatus (Figure 1) to facilitate spotting of the TLC plate. The device is made from a cylinder of polymer in which the top is bored out to receive a flat-nosed, 10-microliter syringe, which is held in place by an O-ring (a trace of silicone grease on the O-ring facilitates the insertion and removal of the syringe barrel). A hose connector is attached to the side of the cylinder and used as an inlet for nitrogen. A hole in the bottom of the cylinder is fitted with a 25-gauge needle through which the tip of the 10-microliter syringe just emerges when inserted through the top. Nitrogen flows into the body of the device, down through the 25-gauge annulus and around the syringe needle, where it evenly and quickly evaporates the solvent as it is applied by depression of the syringe plunger.

This device has saved me countless hours while spotting TLC plates with flash chromatography fractions, especially those in higher boiling solvents. You also get narrower spots and thus achieve better resolution upon development of the TLC plate.

Dennis K. Klipa, Ph.D.
The Dow Chemical Company
Custom & Fine Chemicals
1710 Building
Midland, MI 48674
E-mail: klipadk@dow.com

*The new
13th edition
is now available!*




The Merck Index

13th ed., Merck & Co.,
Rahway, NJ, 2001,
1,741pp. Hardcover.

Update your bookshelf with one of the best-known and respected sources of chemical information. The thirteenth edition contains over 10,000 monographs on chemicals, biologicals, and drugs, and provides a variety of information such as CAS registry numbers, molecular weights, physical properties, toxicity data, descriptions of commercial uses, and references to the scientific and patent literature. Includes an organic name reaction section and a name index with over 50,000 names.

Z52,876-5 \$59.95

Do you have an innovative shortcut or unique laboratory hint you'd like to share with your fellow chemists? If so, please send it to Aldrich (attn: Lab Notes, Aldrichimica Acta). For submitting your idea, you will receive a complimentary, Aldrich periodic table poster (Cat. No. Z54,320-9). If we publish your lab note, you will also receive an Aldrich periodic table mouse pad (Cat. No. Z54,323-3). It is Teflon®-coated, 8½ in. x 11 in., with a full-color periodic table on the front. We reserve the right to retain all entries for future consideration.



Teflon is a registered trademark of E. I. du Pont de Nemours & Co., Inc.

Sigma-Aldrich Worldwide Locations

Argentina
Sigma-Aldrich de Argentina, S.A
Tel: 54 11 4807 0321
Fax: 54 11 4807 0346
E-mail: info@sigma-aldrich.com.ar

Australia
Sigma-Aldrich Pty. Limited
Free Tel: 1-800-800-097
Free Fax: 1-800-800-096
Tel: 61-2-9841-0555
Fax: 61-2-9841-0500
E-mail: ausmail@sial.com

Austria
Sigma-Aldrich Handels GmbH
Tel: 43 1 605 8110
Fax: 43 1 605 8120
E-mail: sigma@sial.com.at

Belgium
Sigma-Aldrich nv/sa
Free Tel: 0800-147.47
Free Fax: 0800-147.45
Tel: 03 899.13.01
Fax: 03 899.13.11
E-mail: sigmaldr@sial.be

Brazil
Sigma-Aldrich Quimica Brasil Ltda.
Tel: 55-11-231-1866
Fax: 55-11-257-9079
E-mail: sigmabr@bm.net

Canada
Sigma-Aldrich Canada Ltd.
Free Tel: 800-565-1400
Free Fax: 800-265-3858
Tel: 905-829-9500
Fax: 905-829-9292
E-mail: canada@sial.com

Czech Republic
Sigma-Aldrich s.r.o.
Tel: 00-420-2-2176 1300
Fax: 00-420-2-2176 3300
E-mail: CZECUSTSV@eurnotes.sial.com

Denmark
Sigma-Aldrich Denmark A/S
Tel: +45 43 56 59 10
Fax: +45 43 56 59 05
E-mail: denorder@eurnotes.sial.com

Eire
Sigma-Aldrich Ireland Ltd.
Free Tel: 1 800 200 888
Free Fax: 1 800 600 222
Tel: 01 404 1900
Fax: 01 404 1910
E-mail: EICUSTSV@eurnotes.sial.com

Finland
Sigma-Aldrich Finland
Tel: (09) 350-9250
Fax: (09) 350-92555
E-mail: finorder@eurnotes.sial.com

France
Sigma-Aldrich Chimie S.a.r.l
Tel Numero Vert: 0800 21 14 08
Fax Numero Vert: 0800 03 10 52
E-mail: fradvsv@eurnotes.sial.com

Germany
(also SE Europe, the Baltics, Africa and the Middle East)

Sigma-Aldrich Chemie GmbH
FreeCall Tel: 0800-5155 000
FreeCall Fax: 0800-6490 000
Tel: +49(0)89-6513-0
Fax: +49(0)89-6513-1169
E-mail: deorders@eurnotes.sial.com

For Slovakia, Africa, Middle East and other European countries:
Tel: +49(0)89-6513-1807
Fax: +49(0)89-6513-1877

Greece
Sigma-Aldrich (o.m) Ltd.
Tel: 30-1-9948010
Fax: 30-1-9943831
E-mail: sigald@accigr

Hungary
Sigma-Aldrich Kft
Tel: 06-1-269-6474
Fax: 06-1-235-9050
E-mail: info@sigma.sial.hu

India
Sigma-Aldrich Corporation
Bangalore location:
Phone: 91 80 851 8797
Fax: 91 80 851 8358

New Delhi location:
Phone: (011) 689 9826
Fax: (011) 689 9827
E-mail: sigma@del2.vsnl.net.in

Israel
Sigma Israel Chemicals Ltd.
Toll Free Tel: 1-800-70-2222
Tel: 972-8-948-4222
Fax: 972-8-948-4200
E-mail: sigisr@sigma.co.il

Italy
Sigma-Aldrich S.r.l.
Tel Numero Verde: 800-827018
Tel: 02-33417.310
Fax: 02-38010.737
E-mail: itorder@eurnotes.sial.com

Japan
Sigma-Aldrich Japan K.K
Tel: 81-3-5821-3111
Fax: 81-3-5821-3170
E-mail: sialjp@sial.com

Korea
Sigma-Aldrich Korea
Toll Free Tel: 080-023-7111
Toll Free Fax: 080-023-8111
Tel: 82-2-783-5211
Fax: 82-2-783-5011
E-mail: order@salkor.co.kr

Malaysia
Sigma-Aldrich (M) Sdn Bhd
Tel: 603-56353321
Fax: 603-56354116
E-mail: sam@sial.com

Mexico
Sigma-Aldrich Quimica, S.A. de C.V.
Tel Gratis: 91-800-0075300
Tel: (72) 76 1600
Fax: (72) 76 1601
E-mail: mexico@sial.com

Netherlands
Sigma-Aldrich Chemie BV
Tel Gratis: 0800 022 9088
Fax Gratis: 0800 022 9089
Tel: 078-620 54 11
Fax: 078-620 54 21
E-mail: nlcustsv@eurnotes.sial.com

Norway
Sigma-Aldrich Norway AS
Tel: 47 23 176 000
Fax: 47 23 176 010
E-mail: nororder@eurnotes.sial.com

Poland
Sigma-Aldrich Sp.z.o
Tel: (48-61) 823-24-81
Fax: (48-61) 823-27-81
E-mail: plcustsv@eurnotes.sial.com

Portugal
Sigma-Aldrich Quimica S.A.
Free Tel: 0800-202180
Free Fax: 0800-202178
E-mail, encomendas:
poorders@eurnotes.sial.com

Russia
Sigma-Aldrich Russia
Tel: 7-095-975-3321
Fax: 7-095-975-4792
E-mail: techcare@online.ru

Singapore
Sigma-Aldrich Pte., Ltd
Tel: (65) 271 1089
Fax: (65) 271 1571
E-mail: sup@sial.com

South Africa
Sigma-Aldrich South Africa (Pty) Ltd
Tel: (011) 397 8886
Fax: (011) 397 8859
E-mail: rsa@sial.com

Spain
Sigma-Aldrich Quimica S.A.
Free Tel: 900-101376
Free Fax: 900-102028
Tel: 91-6619977
Fax: 91-6619642
E-mail, pedidos:
esorders@eurnotes.sial.com



Sweden
Sigma-Aldrich Sweden AB
Tel: 020-35 05 10
Fax: 020-35 25 22
Outside Sweden Tel: 46-8-742 4200
Outside Sweden Fax: 46-8-742 4243
E-mail: sweorder@eurnotes.sial.com

Switzerland
Aldrich Chemie
A division of Fluka Chemie AG
Tel: 081-755 27 23
Fax: 081-755 28 40
Swiss Free Call: 0800 80 00 80
E-mail: Fluka@sial.com

United Kingdom
Sigma-Aldrich Company Ltd.
Free Tel: 0800 717181
Free Fax: 0800 378538
Tel: 01202 733114
Fax: 01202 715460
E-mail: ukcustsv@eurnotes.sial.com

United States
Aldrich Chemical Company, Inc.
Toll Free Tel: 1-800-558-9160
Toll Free Fax: 1-800-962-9591
Tel: 414-273-3850
Fax: 414-273-4979
E-mail: aldrich@sial.com

Sigma Chemical Company
Toll Free Tel: 1-800-521-8956
Toll Free Fax: 1-800-325-5052
Tel: 314-771-5765
Fax: 314-771-5757
E-mail: sigma@sial.com



Modified Dean-Stark Apparatus

Designed specifically for azeotropic distillation using solvents heavier than water (e.g., chloroform). The 8-ml reservoir is specifically designed for small-scale azeotropic distillation of valuable solvents using heavy water. 2mm bore PTFE stopcock, 300mm overall height.

The modified Dean-Stark technique

- Solvent vapors condense and collect in the receiver.
- The modified Dean-Stark trap separates the two components.
- The heavier-than-water layer sinks to the bottom of the trap leaving the water on top.
- The solvent is forced back through the inside tube and returned to the boiling flask.
- The upper water layer is drained out through the 2mm PTFE stopcock.

	Reservoir Cap.	₹ 24/40 Cat. No.	₹ 29/32 Cat. No.	Each
A.	8mL	Z51,617-1	Z51,629-5	\$256.55
B.	35mL	Z24,959-9	Z51,630-9	263.25

Ready to Scale Up?

Reduce Your Reduction Times

Versatile borane reducing agents are now available from Sigma-Aldrich Fine Chemicals in large quantities to meet your growing needs. These products are manufactured and packaged in multiton quantities at our new state-of-the-art manufacturing facility for air-sensitive compounds.


- In stock—Ready to ship
- Packaged for safe and easy handling
- Competitive pricing
- Multiton quantities

L-Selectride[®], 1.0M in THF Steel cylinders
17,849-7 [CAS# 38721-52-7]
Super-Hydride[®], 1.0M in THF Steel cylinders
17,972-8 [CAS# 22560-16-3]
Sodium triacetoxyborohydride, 95% Polylined steel drums
31,639-3 [CAS# 56553-60-7]
Sodium cyanoborohydride, 95% Polylined steel drums
15,615-9 [CAS# 25895-60-7]

Selectride and Super-Hydride are registered trademarks of Sigma Aldrich Co.

www.sigma-aldrich.com/safc

Call us to place an order at 314-534-4900 (USA), or contact your local Sigma-Aldrich office.


SIGMA-ALDRICH
FINE CHEMICALS

ALDRICH CHEMICAL COMPANY, INC.
P.O. BOX 14508
ST. LOUIS, MISSOURI 63178 USA



CHANGE SERVICE REQUESTED

PRESORTED STANDARD
U.S. POSTAGE
PAID
MILWAUKEE, WISCONSIN
PERMIT NO. 552

014

DEDICATED TO PROFESSOR H. C. BROWN ON HIS 90TH BIRTHDAY

Aldrichimica ACTA

VOL. 35, NO. 1 • 2002



Herbert C. Brown
and Aldrich

Asymmetric
Reduction of Ketones
and Imines


Pinane-Based
“Allyl”boranes



ALDRICH[®]

CHEMISTS HELPING CHEMISTS

NEW Products



These Boc-protected piperidines have been utilized extensively in the drug discovery arena. Recent applications include their use in the synthesis of several types of GPIIb/IIIa antagonists^{1,4} and of potent and selective inhibitors of acetylcholinesterase.⁵


(1) Hoekstra, W. J. et al. *J. Med. Chem.* 1999, 42, 5254. (2) Fisher, M. J. et al. *Bioorg. Med. Chem. Lett.* 2000, 10, 835. (3) Egbertson, M. S. et al. *J. Med. Chem.* 1999, 42, 2409. (4) Xue, C.-B. et al. *Bioorg. Med. Chem. Lett.* 1998, 8, 3499. (5) Villalobos, A. et al. *J. Med. Chem.* 1994, 37, 2721.

55,601-7

N-Boc-4-piperidinemethanol, 97% 1g \$20.20; 5g \$67.25

54,724-7

N-Boc-4-piperidineethanol, 97% 1g \$22.80; 5g \$76.00




This versatile hydroquinone has been employed as a building block for the synthesis of biologically active marine metabolites such as *ent*-chromazonarol,¹ 9,11-drimen-8 α -ol,² and fulvanin.² It has also found applications in nonlinear optics as a monomer for phenylene-ethynylene oligomers.³

(1) Barrero, A. F. et al. *Bioorg. Med. Chem. Lett.* 1999, 9, 2325. (2) Barrero, A. F. et al. *Synlett* 2000, 11, 1561. (3) Wautelet, P. et al. *Macromolecules* 1996, 29, 446.

16,713-4

2-Bromohydroquinone, 97% 1g \$17.10; 10g \$95.00




This convenient building block allows the facile preparation of 4-(dimethylamino)benzoyl-substituted compounds.

Purchase, T. S. et al. *Bioorg. Med. Chem.* 1997, 5, 739.

52,611-8

4-(Dimethylamino)benzoyl chloride, 97%
1g \$6.00; 5g \$19.00




Utilized as a milder alternative to dichlorophenylphosphine (Aldrich catalog number D7,198-4).

(1) Nettekoven, U. *J. Org. Chem.* 1999, 64, 3996. (2) Tollefson, M. B. et al. *J. Am. Chem. Soc.* 1996, 118, 9052.

55,470-7

Bis(diethylamino)phenylphosphine, 97%
1g \$8.00; 5g \$25.50




Employed in the synthesis of amino-substituted vinylsilanes, which were smoothly cyclized to the corresponding pyrrolidines.

Miura, K. et al. *Org. Lett.* 2000, 2, 385.

56,237-8

Benzylchlorodimethylsilane, 97% 1mL \$15.00; 5mL \$50.00




Has been used in the synthesis of lycorine-type alkaloids¹ and substituted phthalazines.²

(1) Lida, H. et al. *J. Chem. Soc., Perkin Trans. I* 1975, 2502. (2) Watanabe, N. et al. *J. Med. Chem.* 1998, 41, 3367.

56,504-0

3-Chloro-4-methoxybenzaldehyde, 97%
1g \$27.55; 5g \$91.75




This bromo aldehyde was employed in total syntheses of (-)-lycoridine,¹ benzo[c]phenanthridines,² and the diaza analog of ellagic acid.³

(1) Keck, G. E. *J. Am. Chem. Soc.* 1999, 121, 5176.

(2) Bernabe, P. *Tetrahedron Lett.* 1998, 39, 9785. (3) Kanojia, R. M. *ibid.* 1995, 36, 8553.

56,301-3

6-Bromopiperonal, 97% 50g \$17.80



Important intermediate for the preparation of macrocyclic ketones, such as (*R*)-(-)-muscone,¹

and biologically active macrocyclic lactones.²

(1) Kamat, V. P. et al. *Tetrahedron* 2000, 56, 4397. (2) Kobayashi, Y.; Okui, H. *J. Org. Chem.* 2000, 65, 612.

56,085-5

10-Bromo-1-decene, 97% 1g \$19.25; 5g \$64.05



This versatile reagent has been used in Sonogashira¹ and Suzuki coupling reactions in the presence of carbon monoxide to produce α -pyridyl ketones,² as well as in coupling reactions with aryl thiols to produce heteroaromatic thioethers.³

(1) Koseki, Y. et al. *Tetrahedron Lett.* 2000, 41, 2377. (2) Couve-Bonnaire, S. et al. *ibid.* 2001, 42, 3689. (3) Schopfer, U.; Schlapbach, A. *Tetrahedron* 2001, 57, 3069.

55,876-1

2-Iodopyridine, 98% 1g \$16.00; 5g \$59.55

For competitive quotes on larger quantities contact www.sigma-aldrich.com/safc. SIGMA-ALDRICH FINE CHEMICALS

Aldrichimica ACTA

VOL. 35, NO. 1 • 2002

Aldrich Chemical Co., Inc.

1001 West Saint Paul Ave.
Milwaukee, WI 53233 USA

To Place Orders

Telephone 800-558-9160 (USA)
or 414-273-3850
FAX 800-962-9591 (USA)
or 414-273-4979
Mail P.O. Box 2060
Milwaukee, WI 53201 USA

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201 USA

Customer & Technical Services

Customer Inquiries	800-558-9160
Technical Service	800-231-8327
Sigma-Aldrich Fine Chemicals	800-336-9719
Custom Synthesis	800-336-9719
Flavors & Fragrances	800-227-4563
International	414-273-3850
24-Hour Emergency	414-273-3850
Web Site	http://www.sigma-aldrich.com
E-Mail	aldrich@sial.com

Aldrichimica Acta is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich family. © 2002 by Sigma-Aldrich Co. Printed in the United States.

To request your **FREE** subscription to the *Aldrichimica Acta*,

please call: **800-558-9160** (USA),

write: **Attn: Mailroom**
Aldrich Chemical Co., Inc.
P.O. Box 355
Milwaukee, WI 53201-9358

or email: **sams-usa@sial.com**

International customers, please contact your local Sigma-Aldrich office.

The *Aldrichimica Acta* is also available on the Internet at <http://www.sigma-aldrich.com>.


All prices listed in this publication are subject to change without notice.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

For worldwide contact information, please see the inside back cover.

“Please
Bother
Us.”

Clint Lane, President



Dr. Indrapal Singh Aiden and his research group at the Indian Institute of Technology Madras (Chennai, India) have developed *N*-methoxy-*N*-methyl-2-phenylsulfonylacetamide as a versatile reagent for the two-carbon homologation of alkyl halides under mild conditions (K_2CO_3 in DMF). Dr. Aiden kindly suggested that we offer this valuable amide to our customers.

Satyamurthi, N.; Singh, J.; Aiden, I. S. *Synthesis* 2000, 375.
Satyamurthi, N.; Aiden, I. S. *Carbohydr. Lett.* 1999, 3, 355.

56,108-8

N-Methoxy-*N*-methyl-2-phenylsulfonylacetamide, 97% **1g \$13.00; 5g \$42.90**

Naturally, we made this useful reagent. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page or on the inside back cover.

Table of Contents

Herbert C. Brown and Aldrich2
—S. J. Firsan

Boron-Based Reducing Agents for the Asymmetric Reduction of Functionalized Ketones and Ketimines3
—B. T. Cho

Pinane-Based Versatile “Allyl”boranes23
—P. V. Ramachandran


About Our Cover

Le *Château Noir* (oil on canvas, 73.7 cm x 96.6 cm) was painted by Paul Cézanne between 1900 and 1904. At the center of the painting stands the chateau, with its Gothic windows, red door, and terra cotta and ochre walls contrasted with the deep green of the surrounding vegetation and the intense blue of the Provençal sky. The title, *The Black Chateau*, comes from the fact that the original owner, a manufacturer of lampblack (a pigment made from soot), had the house's interior walls painted black.

Cézanne was born in nearby Aix-en-Provence, and as a youth took long walks to the chateau and the surrounding estate. His father wanted him to study law, but in 1861 he abandoned his legal studies and went to Paris to devote himself to painting. There, he met the artist Camille Pissarro and became acquainted with the impressionists. He exhibited with these artists in 1874 and 1877, but his artistic aims were never entirely compatible with impressionism, which sought to capture through the analysis of color and light the effect of transitory visual sensation. The impressionist technique tends to dissolve any sense of solid form, but Cézanne declared that he wished “to make of impressionism something solid and durable, like the art of the museums”. He wanted to “do Poussin over again, after nature”, which is to say that he wanted to reconcile two seeming opposites, the immediacy of observed nature and the timeless order of the compositions of his seventeenth-century predecessor. He increasingly thought in terms of simple planes and volumes modeled in color, which he understood as the embodiment of the underlying forms of visible objects.

In *Le Château Noir*, the walls of the house, the path, the tree trunks and leafy branches, even the sky, are painted in planes of rich color. Cézanne had no interest in specific or picturesque details, but believed that a painting must succeed both as an illusion of depth and as a textural surface of heavy and vibrant pigments. The result of this approach was an ever-greater abstraction, and his influence on artists like Picasso and Braque, who in time almost completely abandoned any connection to the illusionism of traditional representation, was incalculable.

This painting is a gift of Eugene and Agnes Meyer to the National Gallery of Art, Washington, D.C.



Herbert C. Brown and Aldrich

There is hardly any chemist today, who is not aware of Professor Brown's accomplishments in, and impact on, the field of chemistry. Perhaps no one is more keenly aware of his lasting contributions than we are at Aldrich Chemical Company.

When Professor Brown finally convinced a small chemical company (Aldrich) of the merits of commercializing his hydroboration technology,¹ a 30-year, highly successful collaboration began between Dr. Brown and Aldrich Chemical Company. Thus, the Aldrich-Boranes "acorn" was planted just over 30 years ago and has now grown into a tall, sturdy "oak". Some of the fruits of this collaboration have been the wide array of boron compounds that Aldrich has made available to chemists worldwide, and the financial rewards that the Purdue Research Foundation² and the recipients of the annual Herbert C. Brown Award³ have reaped.

Brown served on the Aldrich Board of Directors (1972–1975), the Sigma-Aldrich Board of Directors (1975–1979), and has been retained as a consultant by Aldrich (1972–Present). A few of Brown's students came to work for Aldrich and made important contributions to the growth of the business, perhaps none more so than Clint Lane, Aldrich's current president.

No write-up about Professor Brown, however brief, would be complete without mentioning his lifelong passion for *chemistry*.⁴ As attendees of professional meetings such as the ACS or the Boron Americas meetings can attest to, Nobel Laureate Brown still attends as many talks as he possibly can and unassumingly walks around the exhibition halls checking out exhibits, visiting vendors, and picking up chemical literature—just as an eager, young graduate student would do, except that Brown is nearly 90 years old! Professor Brown's other lifelong love is Sarah Baylen Brown, his wife of 65 years. As Brown put it "I take care of the chemistry, and Sarah takes care of the money and everything else". Brown relates a pertinent story that is dear to his heart: not long ago both of them were attending a professional function and were observed to treat each other affectionately. This prompted a couple, who didn't know them well, to ask them if they were newlyweds! After family and chemistry in general, Brown is perhaps most fond of boron and its compounds; he was particularly thrilled when he bought a laser pointer that emitted a green light—the color associated with the combustion of boron compounds. In fact, those very close to him know that just about the only way to get him off the subject of boron chemistry is to ask him about his two granddaughters.

In writing this piece about Professor Brown, we were hard-pressed to learn about any interests he might have outside of chemistry. The only one we could identify is photography. Of the few things he dislikes, doing chores that require physical exertion appears to be foremost.⁵

Professor Brown.

Your friends and admirers at Aldrich thank you for your lasting contributions to chemistry and to the success of Aldrich, and wish you a happy ninetieth birthday.⁶ Many of us look forward to celebrating your birthday—with you, Sarah, your son Charles, your two granddaughters, Tamar and Ronni, and many of the more than 500 former Brown students—at a special symposium to be held on May 23–25, 2002 at Purdue University.



Sarah and Herbert Brown in front of the new building devoted to the production of air-sensitive compounds at the Aldrich plant near Sheboygan Falls, WI (October 14, 1999).

Photo courtesy of H.C. Brown

References and Notes: (1) Firsan, S. J. *Aldrichimica Acta* 2001, 34, 35. (2) Reference 1, p 44. (3) *The Herbert C. Brown Award for Creative Research in Synthetic Methods*, administered by the American Chemical Society (ACS) and co-sponsored by Aldrich Chemical Company. More information on this award is found on the ACS Web site at www.chemistry.org. (4) (a) Many anecdotes and fascinating information, particularly about Brown's early years, have been compiled in a book celebrating Brown's 66th birthday and retirement: *Remembering HCB: Memoirs of Colleagues and Students of Herbert C. Brown*; Bank, S., Ed.; Purdue University: West Lafayette, IN, 1978. (b) An essay honoring Professor Brown on the occasion of his 75th birthday has appeared in an earlier issue of this magazine: Brewster, J. H. *Aldrichimica Acta* 1987, 20, 3. (5) Reference 4a, p 5. (6) Herbert C. Brown's ninetieth birthday falls on May 22, 2002.

Sigma-Aldrich Worldwide Locations

Argentina

Sigma-Aldrich de Argentina SA
Phone: 54 11 4556 1472
Fax: 54 11 4552 1698
E-mail: info-argentina@sial.com

Australia

Sigma-Aldrich Pty. Ltd.
Phone: 61 2 9841 0555
Fax: 61 2 9841 0500
E-mail: ausmail@sial.com

Austria

Sigma-Aldrich Handels GmbH
Phone: 43 1 605 8110
Fax: 43 1 605 8120
E-mail: sigma@sial.com.at

Belgium

Sigma-Aldrich N.V./S.A.
Phone: 32 3 899 1301
Fax: 32 3 899 1311
E-mail: becustsv@eurnotes.sial.com

Brazil

Sigma-Aldrich Brasil Ltda.
Phone: 55 11 3732 3100
Fax: 55 11 3733 5151
E-mail: sigmabr@sial.com

Canada

Sigma-Aldrich Canada Ltd.
Phone: (905) 829-9500
Fax: (905) 829-9292
E-mail: canada@sial.com

Czech Republic

Sigma-Aldrich s.r.o.
Phone: 420 2 2176 1111
Fax: 420 2 2176 3300
E-mail: czecustsv@eurnotes.sial.com

Denmark

Sigma-Aldrich Denmark A/S
Phone: 45 43 565900
Fax: 45 43 565905
E-mail: denorder@eurnotes.sial.com

Finland

Sigma-Aldrich Finland YA-Kemia Oy
Phone: 358 9 350 9250
Fax: 358 9 350 92555
E-mail: finorder@eurnotes.sial.com

France

Sigma-Aldrich Chimie S.a.r.l.
Phone: 33 4 74 822888
Fax: 33 4 74 956808
E-mail: fradsvs@eurnotes.sial.com

Germany

Sigma-Aldrich Chemie GmbH
Phone: 49 89 6513 0
Fax: 49 89 6513 1169
E-mail: deorders@eurnotes.sial.com

Greece

Sigma-Aldrich (O.M.) Ltd.
Phone: 30 10 994 8010
Fax: 30 10 994 3831
E-mail: grcustsv@eurnotes.sial.com

Hungary

Sigma-Aldrich Kft
Phone: 36 1 269 6474
Fax: 36 1 235 9050
E-mail: info@sigma.sial.hu

India

Sigma-Aldrich Foreign Holding Company
Phone: 91 80 852 4222
Fax: 91 80 852 4214
E-mail: india@sial.com

Ireland

Sigma-Aldrich Ireland Ltd.
Phone: 353 1 404 1900
Fax: 353 1 404 1910
E-mail: eicustsv@eurnotes.sial.com

Israel

Sigma-Aldrich Israel Ltd.
Phone: 972 8 948 4222
Fax: 972 8 948 4200
E-mail: hexter@sigma.co.il

Italy

Sigma-Aldrich S.r.l.
Phone: 39 02 33417 310
Fax: 39 02 38010 737
E-mail: itorder@eurnotes.sial.com

Japan

Sigma-Aldrich Japan K.K.
Phone: 81 3 5821 3111
Fax: 81 3 5821 3170

Korea

Sigma-Aldrich Korea
Phone: 82 31 329 9000
Fax: 82 31 329 9090
E-mail: hcjoo@sial.com

Malaysia

Sigma-Aldrich (M) Sdn. Bhd.
Phone: 60 3 5635 3321
Fax: 60 3 5635 4116
E-mail: sam@sial.com

Mexico

Sigma-Aldrich Quimica, S.A. de C.V.
Phone: 52 722 276 1600
Fax: 52 722 276 1601
E-mail: mexico@sial.com

Netherlands

Sigma-Aldrich Chemie B.V.
Phone: 31 78 620 5411
Fax: 31 78 620 5421
E-mail: nlcustsv@eurnotes.sial.com

Norway

Sigma-Aldrich Norway AS
Phone: 47 23 176 000
Fax: 47 23 176 010
E-mail: nortech@sial.com

Poland

Sigma-Aldrich Sp. z o.o.
Phone: 48 61 823 2481
Fax: 48 61 823 2781
E-mail: plcustsv@eurnotes.sial.com

Portugal

Sigma-Aldrich Quimica, S.A.
Phone: 351 1 924 2555
Fax: 351 1 924 2610
E-mail: poorders@eurnotes.sial.com

Rep. of Singapore

Sigma-Aldrich Pte., Ltd.
Phone: 65 271 1089
Fax: 65 271 1571
E-mail: sapl@sial.com

Russia

OOO SAF-LAB
Phone: 7 095 975 3321
Fax: 7 095 975 4792
E-mail: techcare@sial.com

South Africa

Sigma-Aldrich (Pty.) Ltd
Phone: 27 11 397 8886
Fax: 27 11 397 8859
E-mail: rsa@sial.com

Spain

Sigma-Aldrich Quimica S.A.
Phone: 34 91 661 9977
Fax: 34 91 661 9642
E-mail: esorders@eurnotes.sial.com

Sweden

Sigma-Aldrich Sweden AB
Phone: 46 8 742 4200
Fax: 46 8 742 4243
E-mail: sweorder@eurnotes.sial.com

Switzerland

Fluka Holding AG
Phone: 41 81 755 2511
Fax: 41 81 756 5449
E-mail: fluka@sial.com

United Kingdom

Sigma-Aldrich Company Ltd.
Phone: 44 1202 733114
Fax: 44 1202 715460
E-mail: ukcustsv@eurnotes.sial.com

United States

Aldrich
Phone: (414) 273-3850
Fax: (414) 273-4979
E-mail: aldrich@sial.com

Sigma

Phone: (314) 771-5765
Fax: (314) 771-5757
E-mail: sigma@sial.com

Phosphorus Pentasulfide

Available in 55-gallon drums
Shipped from Milwaukee, WI



SIGMA-ALDRICH
FINE CHEMICALS

Product No. 54,462-0 • CAS No. [1314-80-3]

Call Dick Farrell or Ted Eslinger at 800-213-1206 (USA)

Air-Sensitive and Low-Temperature Reactions *in* Large Scale

To meet the growing need for large-scale air-sensitive compounds in pharmaceutical and industrial applications, Sigma-Aldrich's Pro-T plant in Sheboygan, WI, is ready to serve.

Examples of products we offer:

- Arylboronic acids
- Boranes and boron hydrides
- Grignards
- Lithium compounds
- Metallocenes
- Organoaluminum
- Organometallics
- Organotin
- Silicon reagents
- Custom synthesis

We perform reactions under inert atmospheres in glass-lined, stainless steel, and Hastelloy® vessels (up to 4,000 gal), and low temperature reactions down to -110°C (Hastelloy® reactors up to 500 gal). We manufacture an extensive list of air-sensitive reagents, and can use them to synthesize custom organic and inorganic materials for you.

More than 25 years of experience with air-sensitive reactions, combined with large-scale manufacturing and packaging capabilities, make Sigma-Aldrich a partner you can count on.

Contact us today to get your next air-sensitive or low-temperature project on the way!

Visit our Web site
for a Virtual Tour of
our Facilities

DOT-E 11167

www.sigma-aldrich.com/safc

SIGMA-ALDRICH
FINE CHEMICALS

France	Germany	Italy	Japan	Switzerland	U.K.	USA
tel. 33 4 74 82 28 82	tel. 49 89 6513 1900	tel. 3902 33 41 7350	tel. 81 3 5821 3171	tel. 41 81 755 2732	tel. 44 1202 712222	tel. 314 534 4900
fax 33 4 74 82 28 90	fax 49 89 6513 1919	fax 3902 33 41 7201	fax 81 3 5614 6279	fax 41 81 755 2770	fax 44 1202 740171	fax 314 652 0000

ALDRICH CHEMICAL COMPANY, INC.
P.O. BOX 14508
ST. LOUIS, MISSOURI 63178 USA



CHANGE SERVICE REQUESTED

PRESORTED STANDARD
U.S. POSTAGE
PAID
SIGMA-ALDRICH

014

WE PROVIDE BUILDING BLOCKS FOR SUCCESSFUL R&D

Aldrichimica ACTA

VOL. 35, NO. 2 • 2002



Asymmetric
Epoxidation of
 α,β -Unsaturated
Ketones
Catalyzed by
Poly(amino acids)


Borane–Amine
Complexes for
Hydroboration



ALDRICH[®]

CHEMISTS HELPING CHEMISTS

NEW products




Was employed in the synthesis of nonlinear optical polymers,¹ and served as a starting material for the construction of complex organometallics.²

(1) Buchmeiser, M.; Schrock, R.R. *Macromolecules* **1995**, *28*, 6642. (2) Buchmeiser, M.; Schottenberger, H. *J. Organomet. Chem.* **1992**, *441*, 457.

44,265-8

Ethynylferrocene, 97% **1g \$59.00**



These chiral amines have been utilized in the asymmetric synthesis of β-lactams¹ and β-amino acids.²


(1) Bull, S. D. et al. *J. Chem. Soc., Perkin Trans. I* **2001**, 3106. (2) Davies, S. G.; Ichihara, O. *Tetrahedron Lett.* **1999**, *40*, 9313.

52,554-5

(S)-(-)-N-Allyl-α-methylbenzylamine, 97% **1mL \$9.90; 10mL \$55.00**

55,903-2

(R)-(+)-N-Allyl-α-methylbenzylamine, 97% **1g \$58.15; 5g \$193.80**




Utilized in the platinum-catalyzed hydrosilylation of olefins for the generation of dendrons with a chlorosilyl functional group.

Cuadrado, I. et al. *J. Am. Chem. Soc.* **1997**, *119*, 7613.

56,923-2

Chlorophenylsilane, 97% **1g \$16.95; 5g \$56.40**




Employed in the synthesis of cyclopentano-1,2,3,4-tetrahydroisoquinolines¹ and a novel 5-HT_{2c} receptor agonist.²

(1) Mathison, I. W. et al. *J. Org. Chem.* **1974**, *39*, 2852. (2) Adams, D. R.; Dunton, M. A. *J. Synth. Commun.* **2001**, *31*, 2029.

56,825-2

4,5-Dimethoxy-1-indanone, 97% **1g \$28.80; 5g \$96.00**




Has been employed to produce isoxazolidines through 1,3-dipolar cycloaddition reactions.

Chiacchio, U. et al. *Tetrahedron* **1994**, *50*, 6671.

56,306-4

2-Vinylanisole, 98% **1g \$22.60; 5g \$75.30**




Basic building block for the synthesis of 1-aryl-2,3,4,5-tetrahydro-1H-3-benzazepines.^{1,2}

(1) Weinstock, J. et al. *J. Med. Chem.* **1986**, *29*, 2315. (2) Gallagher, G., Jr. et al. Ger. Patent 2,802,087, Jul 1978; *Chem. Abstr.* **1978**, *89*, 163433v.

56,939-9

3-Chloro-4,5-dimethoxybenzaldehyde, 97% **1g \$25.60; 5g \$100.00**




Used in the preparation of 1,3-dienes with high *E*-isomer selectivity through the Horner-Wadsworth-Emmons reaction.^{1,2}

(1) West, F.G.; Wang, Y. *Synthesis* **2002**, 99. (2) Maryanoff, B.E.; Reitz, A. B. *Chem. Rev.* **1989**, *89*, 863.

56,541-5

Diethyl allyl phosphonate, 98% **1g \$13.35; 5g \$44.40**




Tetralones have proven their synthetic versatility in the development of numerous biologically active compounds. 5,8-Dimethoxy-1-tetralone, for example, is a key intermediate in the synthesis of the antifungal and anticancer natural product (±)-diepoxin σ and related isomers.

Wipf, P.; Jung, J.-K. *J. Org. Chem.* **2000**, *65*, 6319.

56,965-8

5,8-Dimethoxy-1-tetralone, 99% **1g \$42.55; 5g \$141.80**




Utilized in palladium-catalyzed coupling reactions that lead to vinylsilanes stereoselectively.

Chatani, N. et al. *J. Org. Chem.* **1995**, *60*, 1834.

51,941-3

3-Ethynylanisole, 96% **5g \$74.00**



An excellent glycoside for the stereoselective synthesis of complex oligosaccharides and glycoconjugates.^{1,2}

(1) Merayala, H.B.; Reddy, G.V. *Tetrahedron* **1991**, *47*, 6435. (2) Garcia, B.A.; Gin, D.Y. *J. Am. Chem. Soc.* **2000**, *122*, 4269.

57,004-4

Isopropyl 2,3,4,6-tetra-O-benzyl-β-D-glucopyranoside, 98% **1g \$32.35; 5g \$107.70**

Aldrichimica ACTA

VOL. 35, NO. 2 • 2002

Aldrich Chemical Co., Inc.

1001 West Saint Paul Ave.
Milwaukee, WI 53233 USA

To Place Orders

Telephone 800-558-9160 (USA)

or 414-273-3850

FAX 800-962-9591 (USA)

or 414-273-4979

Mail P.O. Box 2060

Milwaukee, WI 53201 USA

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.

P.O. Box 355, Milwaukee, WI 53201 USA

Customer & Technical Services

Customer Inquiries 800-558-9160

Technical Service 800-231-8327

Sigma-Aldrich Fine Chemicals 800-336-9719

Custom Synthesis 800-336-9719

Flavors & Fragrances 800-227-4563

International 414-273-3850

24-Hour Emergency 414-273-3850

Web Site <http://www.sigma-aldrich.com>

E-Mail aldrich@sial.com

To request your **FREE** subscription to the *Aldrichimica Acta*,

please call: **800-558-9160** (USA),

write: **Attn: Mailroom**
Aldrich Chemical Co., Inc.
P.O. Box 355
Milwaukee, WI 53201-9358

or email: **sams-usa@sial.com**

International customers, please contact your local Sigma-Aldrich office.

The *Aldrichimica Acta* is also available on the Internet at <http://www.sigma-aldrich.com>.

All prices listed in this publication are subject to change without notice.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

For worldwide contact information, please see the inside back cover.

Aldrichimica Acta is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich family. © 2002 by Sigma-Aldrich Co.

About Our Cover

Baby at Play (oil on canvas, 81.9 x 122.9 cm) was painted by Thomas Eakins in 1876. Eakins was born and grew up in Philadelphia. After taking courses at the Pennsylvania Academy of the Fine Arts, and studying anatomy there and at the Jefferson Medical College, he went in 1866 to Paris, where he attended the École des Beaux-Arts and studied independently with several artists, including Jean-Léon Gérôme. During the five or six years after his return to Philadelphia three years later, Eakins painted a number of both individual and group portraits of members of his family and that of his sister Frances,

who married Eakins's boyhood friend, William J. Crowell. The subject of *Baby at Play* is Eakins's two-year-old niece, Ella Crowell. The low vantage point from which we see the figure, which is life-size, helps to draw attention to the child's almost solemn concentration on what she is doing. The painting is a good example of what the artist called "genre portraits", i.e., representations of people shown engaged in typical activities in their usual settings. It is probable that it is more than that, however.

The picture was painted just at the time that educational reformers in Europe, especially the English social philosopher Herbert Spencer, whose *Education: Intellectual, Moral, and Physical* was published in New York in 1866, were having a great effect on educational thought and practice in America. Eakins, who was well-known for his serious dedication to the disciplined teaching of art, mentioned Spencer in a letter to his father just two years later. It is quite likely that this painting reflects Spencer's statement that "men are at last seeing that the spontaneous activity of the observing faculties in children has a meaning and a use. What was once thought mere purposeless action, or play, or mischief, as the case might be, is now recognized as the process of acquiring a knowledge upon which all after-knowledge is based." Eakins's acquaintance with Spencer's ideas means that this painting probably should not be considered as simply a representation of light-hearted childish play, but as a depiction of a more serious endeavor, the learning a child gains through experience and observation.

This painting is in the John Hay Whitney Collection at the National Gallery of Art, Washington, D.C.




Photograph © Board of Trustees, National Gallery of Art, Washington

"Please
Bother
Us"

Clint Lane

Clint Lane, President



Professor Matthew H. Todd of the Department of Chemistry at Queen Mary, University of London, kindly suggested that we offer diisopropylzinc. This organometallic reagent undergoes an asymmetric autocatalytic reaction with substituted pyrimidines and a chiral initiator.

Shibata, T. et al. *J. Am. Chem. Soc.* **1998**, *120*, 12157

56,811-2


Diisopropylzinc, 1.0M solution in toluene
5mL \$66.00; 25mL \$222.00

Naturally, we made this useful reagent. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page or on the inside back cover.

Table of Contents

Asymmetric Epoxidation of α,β -Unsaturated Ketones Catalyzed by Poly(amino acids)	47
-C. Lauret and S. M. Roberts	
Borane-Amine Complexes for Hydroboration	57
-J. V. B. Kanth	



Lab Notes

A Simple Procedure for Polishing Circular and Rectangular Infrared Crystal Windows¹




Figure 1. IR Crystal Window Before and After Polishing with the Current Procedure.

towel in such a way as to create a smooth, even layer of polishing compound. Grip the sample-free window firmly between the thumb and index finger, and rub its bottom side on the polishing compound using light pressure and tracing a figure 8. The number of times a figure 8 is traced will depend on how bad the window surface is. Lift the window off of the paper towel and trace figure 8's with it on the laboratory towel in order to dry it and rid it of excess ferric oxide.

Inspect the polished surface. If it looks transparent enough, repeat the above procedure with the opposite face of the window. When both surfaces are done, and only if needed, they can be rinsed with an anhydrous solvent such as anhydrous acetone. When not in use, the polished windows should be stored in a desiccator over a drying agent such as Drierite[®].

This polishing procedure is routinely used in our laboratories, and is recommended mainly for NaCl or KBr windows—which are the most widely used types.

Notes: (1) If the IR crystal window is badly scratched or pitted, or if its surfaces are no longer flat from uneven polishing, the damaged window can be reclaimed as follows: It is first lapped on an ultrafine 600 grit sandpaper (e.g., 3M's TRI-M-ITE[™] silicon carbide 600 TN4) to remove surface damage or flatten the entire window face. The window is then polished as described above. (2) For example, Aldrich cat. no. Z11,186-4.


KIMTOWELS and Drierite are registered trademarks of Kimberly-Clark Corp. and W. A. Hammond Drierite Co., respectively. TRI-M-ITE is a trademark of 3M Co.

Editor's Notes: (1) Other, less detailed or related procedures for polishing cloudy infrared crystal windows have been reported: (a) Winston, A. *J. Chem. Educ.* **1991**, *68*, A124 (NaCl). (b) Gallego G., J. M. *J. Chem. Educ.* **1978**, *55*, 681 (NaCl, KCl). (c) Nemo, T. E. *Aldrichimica Acta* **1977**, *10*, 22 (AgCl). (d) Fearheller, W. R., Jr.; DuFour, H. *Appl. Spectrosc.* **1967**, *21*, 45 (CsI, CsBr). (e) Levine, M. A. *Appl. Opt.* **1966**, *5*, 1957 (CsI). (2) Customary lab safety precautions, such as the donning of suitable gloves, should be taken when using this polishing procedure.

James J. Brien
Laboratory Equipment Group
Aldrich Chemical Company, Inc.
1001 West Saint Paul Avenue
Milwaukee, WI 53233, USA
E-mail: jbrien@sial.com

Do you have an innovative shortcut or unique laboratory hint you'd like to share with your fellow chemists? If so, please send it to Aldrich (attn: Lab Notes, Aldrichimica Acta). For submitting your idea, you will receive a complimentary, Aldrich periodic table poster (cat. no. **Z54,320-9**). If we publish your lab note, you will also receive an Aldrich periodic table mouse pad (cat. no. **Z54,323-3**). It is Teflon[®]-coated, 8½ in. x 11 in., with a full-color periodic table on the front. We reserve the right to retain all entries for future consideration.

Teflon is a registered trademark of E.I. du Pont de Nemours & Co., Inc.



Sigma-Aldrich Worldwide Locations

Argentina

Sigma-Aldrich de Argentina SA
Phone: 54 11 4556 1472
Fax: 54 11 4552 1698
E-mail: info-argentina@sial.com

Australia

Sigma-Aldrich Pty. Ltd.
Phone: 61 2 9841 0555
Fax: 61 2 9841 0500
E-mail: ausmail@sial.com

Austria

Sigma-Aldrich Handels GmbH
Phone: 43 1 605 8110
Fax: 43 1 605 8120
E-mail: sigma@sial.com.at

Belgium

Sigma-Aldrich N.V./S.A.
Phone: 32 3 899 1301
Fax: 32 3 899 1311
E-mail: becustsv@eurnotes.sial.com

Brazil

Sigma-Aldrich Brasil Ltda.
Phone: 55 11 3732 3100
Fax: 55 11 3733 5151
E-mail: sigmabr@sial.com

Canada

Sigma-Aldrich Canada Ltd
Phone: (905) 829-9500
Fax: (905) 829-9292
E-mail: canada@sial.com

Czech Republic

Sigma-Aldrich s.r.o.
Phone: 420 2 2176 1310
Fax: 420 2 2176 3300
E-mail: czechustsv@eurnotes.sial.com

Denmark

Sigma-Aldrich Denmark A/S
Phone: 45 43 565910
Fax: 45 43 565905
E-mail: denorder@eurnotes.sial.com

Finland

Sigma-Aldrich Finland YA-Kemia Oy
Phone: 358 9 350 9250
Fax: 358 9 350 92555
E-mail: linorder@eurnotes.sial.com

France

Sigma-Aldrich Chimie S.a.r.l
Phone: 0800 21 14 08
Fax: 0800 03 10 52
E-mail: fradsvs@eurnotes.sial.com

Germany

Sigma-Aldrich Chemie GmbH
Phone: 49 89 6513 0
Fax: 49 89 6513 1169
E-mail: deorders@eurnotes.sial.com

Greece

Sigma-Aldrich (O.M.) Ltd.
Phone: 30 10 994 8010
Fax: 30 10 994 3831
E-mail: grcustsv@eurnotes.sial.com

Hungary

Sigma-Aldrich Kft
Phone: 36 1 269 6474
Fax: 36 1 235 9050
E-mail: info@sigma.sial.hu

India

Sigma-Aldrich Foreign Holding Company
Phone: 91 80 852 4222
Fax: 91 80 852 4214
E-mail: india@sial.com

Ireland

Sigma-Aldrich Ireland Ltd.
Phone: 353 1 404 1900
Fax: 353 1 404 1910
E-mail: eicustsv@eurnotes.sial.com

Israel

Sigma-Aldrich Israel Ltd.
Phone: 972 8 948 4222
Fax: 972 8 948 4200
E-mail: hextex@sigma.co.il

Italy

Sigma-Aldrich S.r.l.
Phone: 39 02 33417 310
Fax: 39 02 38010 737
E-mail: itorder@eurnotes.sial.com

Japan

Sigma-Aldrich Japan K.K.
Phone: 81 3 5821 3111
Fax: 81 3 5821 3170

Korea

Sigma-Aldrich Korea
Phone: 82 31 329 9000
Fax: 82 31 329 9090
E-mail: hcjoo@sial.com

Malaysia

Sigma-Aldrich (M) Sdn. Bhd.
Phone: 60 3 5635 3321
Fax: 60 3 5635 4116
E-mail: sam@sial.com

Mexico

Sigma-Aldrich Quimica, S.A. de C.V.
Phone: 52 722 276 1600
Fax: 52 722 276 1601
E-mail: mexico@sial.com

Netherlands

Sigma-Aldrich Chemie B.V.
Phone: 31 78 620 5411
Fax: 31 78 620 5421
E-mail: nlcustsv@eurnotes.sial.com

Norway

Sigma-Aldrich Norway AS
Phone: 47 23 176 000
Fax: 47 23 176 010
E-mail: nortech@sial.com

Poland

Sigma-Aldrich Sp. z o.o.
Phone: 48 61 823 2481
Fax: 48 61 823 2781
E-mail: plcustsv@eurnotes.sial.com

Portugal

Sigma-Aldrich Quimica, S.A.
Phone: 351 1 924 2555
Fax: 351 1 924 2610
E-mail: poorders@eurnotes.sial.com

Russia

OOO SAF-LAB
Phone: 7 095 975 3321
Fax: 7 095 975 4792
E-mail: techcare@sial.com

Singapore

Sigma-Aldrich Pte., Ltd.
Phone: 65 6271 1089
Fax: 65 6271 1571
E-mail: sapl@sial.com

South Africa

Sigma-Aldrich (Pty.) Ltd
Phone: 27 11 979 1188
Fax: 27 11 979 1119
E-mail: rsa@sial.com

Spain

Sigma-Aldrich Quimica S.A
Phone: 34 91 661 9977
Fax: 34 91 661 9642
E-mail: esorders@eurnotes.sial.com

Sweden

Sigma-Aldrich Sweden AB
Phone: 46 8 742 4200
Fax: 46 8 742 4243
E-mail: sweorder@eurnotes.sial.com

Switzerland

Fluka Chemie GmbH
Phone: 41 81 755 2721
Fax: 41 81 756 2840
E-mail: fluka@sial.com

United Kingdom

Sigma-Aldrich Company Ltd.
Phone: 44 1747 833000
Fax: 44 1747 833313
E-mail: ukcustsv@eurnotes.sial.com

United States

Aldrich
Phone: (414) 273-3850
Fax: (414) 273-4979
E-mail: aldrich@sial.com

Sigma
Phone: (314) 771-5765
Fax: (314) 771-5757
E-mail: sigma@sial.com

Phosphorus Pentasulfide

Available in 55-gallon drums
Shipped from Milwaukee, WI



SIGMA-ALDRICH
FINE CHEMICALS

Product No. 54,462-0 • CAS No. [1314-80-3]

Call Dan Zagrodnik or Ted Eslinger at 800-213-1206 (USA)

A Modern Classic



The 2003–2004 Aldrich Handbook of Fine Chemicals and Laboratory Equipment is now available. In it, you will find more than 4,700 new products, including solvents, reagents, tools, books, and computer software products.

More than just a catalog, our Handbook contains extensive technical data, application information, and helpful reference materials, making it an indispensable resource. No chemistry lab should be without it.

If you have not yet received our Handbook, request your FREE copy today by emailing sams-usa@sial.com, or call our Technical Service Department at 800-231-8327 (USA).

ALDRICH CHEMICAL COMPANY, INC.
P.O. BOX 14508
ST. LOUIS, MISSOURI 63178 USA



CHANGE SERVICE REQUESTED

PRESORTED STANDARD
U.S. POSTAGE
PAID
MILWAUKEE, WISCONSIN
PERMIT NO. 552

014

ALDRICH CONGRATULATES THE 2003 ACS AWARD WINNERS

Aldrichimica ACTA

VOL. 35, NO. 3 • 2002

(Last issue in 2002)



Applications of Ionic
Liquids in Organic Synthesis

Novel Methods of Resolving
Racemic Diols and Amino Alcohols




ALDRICH[®]

CHEMISTS HELPING CHEMISTS

New PRODUCTS

46,068-0 Tricyclohexyltin hydride, 97%

250mg \$44.85; 1g \$124.55




Utilized in the hydrostannation of a variety of functionalized alkenes in good yields.

Jousseau, B. et al. *Organometallics* 1995, 14, 685.

51,158-7 Perruthenate, polymer-bound

5g \$82.00; 25g \$340.00




Has been used in the oxidation of pyridylcarbinols,¹ hydroxylamines,² benzylic alcohols,³ and primary alcohols to their corresponding aldehydes.

(1) Habermann, J. et al. *J. Chem. Soc., Perkin Trans. I* 1999, 1253. (2) Hinzen, B.; Ley, S.V. *J. Chem. Soc., Perkin Trans. I* 1998, 1. (3) Haunert, F. et al. *ibid.* 1998, 2235.

57,553-4 1-(3,5-Dimethoxyphenyl)heptan-1-one, 96%

1g \$25.75; 5g \$85.75




An intermediate in the synthesis of functionalized cannabinoids.^{1,2}

(1) Papahatjis, D. P. et al. *J. Med. Chem.* 1998, 41, 1195. (2) Harrington, P. E. et al. *J. Org. Chem.* 2000, 65, 6576.

57,669-7 Triethyl 1,3,5-triazine-2,4,6-tricarboxylate, 97%

1g \$24.30; 5g \$81.00




A reactive azadiene that is very useful in inverse-electron-demand Diels-Alder reactions for the efficient synthesis of highly functionalized pyrimidines.¹ Also utilized in the synthesis of purines and purine analogs via Diels-Alder reactions with aromatic dienophiles.²

(1) Dang, Q. et al. *J. Org. Chem.* 1996, 61, 5204. (2) Dang, Q. et al. *J. Am. Chem. Soc.* 1999, 121, 5833.

57,893-2 2-Allyl-1,1,1,3,3,3-hexamethyldisilazane, 97%

5g \$18.80; 25g \$62.50




Was utilized to generate aminoalkylboronic acids.

Goeller, B. et al. *Main Group Metal Chemistry* 1997, 20, 795.

57,878-9 5,6-Epoxy-5,6-dihydro[1,10]phenanthroline, 98%

1g \$28.00; 5g \$100.00




Versatile reagent that has been employed in syntheses of 5-substituted 1,10-phenanthrolines¹ and in the preparation of the marine alkaloid acididemin.²

(1) Riklin, M. et al. *J. Chem. Soc., Dalton Trans.* 2001, 1813. (2) Moody, C. J. et al. *Tetrahedron* 1992, 48, 3589.

57,944-0 Bis(trifluoroethyl) methylphosphonate, 99%

5g \$25.00; 25g \$85.00



An excellent starting material for the stereoselective synthesis of Z α,β -unsaturated esters¹ and ketones.²


(1) Still, W.C.; Gennari, C. *Tetrahedron Lett.* 1983, 24, 4405. (2) Yu, Y. et al. *ibid.* 1999, 40, 6725.

57,217-9 *cis*-Propenylboronic acid 1g \$31.30; 5g \$104.20


57,663-8 *trans*-Propenylboronic acid 1g \$59.00; 5g \$188.00

57,135-0 α -Phenylvinylboronic acid 1g \$25.75; 5g \$85.75


57,887-8 5-(4,4,5,5-Tetramethyl-1,3,2-dioxaborolan-2-yl)-2,2'-bithiophene 1g \$30.00; 5g \$88.00




57,217-9



57,663-8



57,135-0



57,887-8


Vinyl boronic acids, like arylboronic acids, undergo facile Suzuki-Miyaura coupling with aromatic halides in the presence of palladium catalysts.¹

Several years ago, Miyaura and co-workers demonstrated the utility of cyclic pinacol esters of arylboronic acids in Suzuki-Miyaura coupling reactions.^{2,3} Very recently, α -phenylvinylboronic acid was used for the preparation of the corresponding nonracemic alcohols via asymmetric hydrogenation, in the presence of a chiral rhodium catalyst, followed by oxidative cleavage.⁴

(1) Miyaura, N.; Suzuki, A. *Chem. Rev.* 1995, 95, 2457. (2) Ishiyama, T. et al. *J. Org. Chem.* 1995, 60, 7508. (3) Ishiyama, T. et al. *Tetrahedron Lett.* 1997, 38, 3447. (4) Ueda, U. et al. *J. Organomet. Chem.* 2002, 642, 145.

57,670-0 2,2'-Bithiophene-5-carbaldehyde, 98%

1g \$22.00; 5g \$80.00



Useful in medicinal chemistry^{1,2} and materials science.^{3,4}

(1) Rodriguez, M. J. et al. *J. Antibiot.* 1998, 51, 560. (2) Xu, W.-C. et al. *Bioorg. Med. Chem. Lett.* 1999, 9, 2279. (3) Kamal, M. R. et al. *Phosphorus, Sulfur Silicon Relat. Elem.* 1997, 126, 65. (4) Soudan, P. et al. *J. Mater. Chem.* 2001, 11, 773.

Please see pages 102-103 for additional new products.

For competitive quotes on larger quantities contact www.sigma-aldrich.com/safc.

SIGMA-ALDRICH
FINE CHEMICALS

Aldrichimica ACTA

VOL. 35, NO. 3 • 2002

Aldrich Chemical Co., Inc.
1001 West Saint Paul Ave.
Milwaukee, WI 53233 USA

To Place Orders

Telephone: 800-558-9160 (USA)
or 414-273-3850
FAX: 800-962-9591 (USA)
or 414-273-4979
Mail: P.O. Box 2060
Milwaukee, WI 53201 USA

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201 USA

Customer & Technical Services

Customer Inquiries 800-558-9160
Technical Service 800-231-8327
Sigma-Aldrich Fine Chemicals 800-336-9719
Custom Synthesis 800-336-9719
Flavors & Fragrances 800-227-4563
International 414-273-3850
24-Hour Emergency 414-273-3850
Web Site <http://www.sigma-aldrich.com>
E-Mail aldrich@sial.com

To request your **FREE** subscription to the *Aldrichimica Acta*,

please call: 800-558-9160 (USA),

write: **Attn: Mailroom**
Aldrich Chemical Co., Inc.
P.O. Box 355
Milwaukee, WI 53201-9358

or email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office.

The *Aldrichimica Acta* is also available on the Internet at <http://www.sigma-aldrich.com>.

All prices listed in this publication are subject to change without notice.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

For worldwide contact information, please see the inside back cover.

Aldrichimica Acta is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich family. © 2002 by Sigma-Aldrich Co. Printed in the USA

About Our Cover

Cattleya Orchid and Three Brazilian Hummingbirds (oil on wood panel, 34.8 x 45.6 cm), signed and dated 1871, was painted by the American artist Martin Johnson Heade. Heade's earliest works were portraits, but by the early 1860s he had turned to landscape, a subject more attuned to his artistic personality. Using a limited number of pictorial elements: sky, clouds, water, perhaps some trees or rocks in an essentially flat landscape, Heade created an art of varied and shifting moods. Often his portrayal of natural phenomena such as shifting sunlight, approaching rain, lightning, dark clouds, and fog gives his paintings a dramatic and even disquieting character.



Photograph © Board of Trustees, National Gallery of Art, Washington


During the 1860s, Heade turned to painting objects at close range, and produced a series of remarkably sensuous still life paintings of flowers. In 1863, he sailed to Rio de Janeiro to study and paint the major species of tropical hummingbirds for a book. The book was never published, but he made two other trips to Central and South America in 1866 and 1870, fascinated by the wildlife and the landscape. His approach was different from that of his friend Frederick Edwin Church, who also traveled to Latin America, but who sought to capture the grandeur of vast tropical landscapes, and it was different from that of naturalists like John James Audubon, whose purpose was to create an objective record of the birds and plants he saw. In *Cattleya Orchid and Three Brazilian Hummingbirds*, Heade carefully represents a specific kind of pink orchid and two particular species of hummingbird, one Sappho Comet, green with a yellow throat and red tail feathers, and two green-and-pink Brazilian Amethysts, but he sets these subjects in an evocative and mysterious tropical setting full of mist and diffuse light, combining the two kinds of painting at which he excelled, still life and landscape.

This painting is a gift of the Morris and Gwendolyn Cafritz Foundation to the National Gallery of Art, Washington, DC.

“Please
Bother
Us.”

Clint Lane

Clint Lane, President



Professor James S. Nowick of the University of California, Irvine, kindly suggested that we offer the amine-protecting group 2,7-di-(*tert*-butyl)-9-fluorenylmethyl chloroformate (Fmoc-Cl) as an alternative to Fmoc-Cl. Dr. Nowick and coworkers have demonstrated that Fmoc-protected amines have a greater solubility in common volatile organic solvents and fewer problems associated with byproduct removal than their Fmoc-protected counterparts.

Stigers, K. D.; Koutroulis, M. R.; Chung, D. M.; Nowick, J. S. *J. Org. Chem.* 2000, 65, 3858.


55,100-7
2,7-Di-(*tert*-butyl)-9-fluorenylmethyl
chloroformate, 97% (Fmoc-Cl)
1g \$31.10

Naturally, we made this valuable protecting group. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page or on the inside back cover.

Table of Contents

Applications of Ionic Liquids in Organic Synthesis.....	75
—H. Zhao and S. V. Malhotra*	
Novel Methods of Resolving Racemic Diols and Amino Alcohols.....	89
—M. Periasamy	




Lab Notes

(1*R*,2*S*,5*R*)-(-)-Menthol: Proposed Calibration Standard for Polarimetry

Since its discovery in the early part of the nineteenth century,¹ polarimetry has been employed extensively in the sugar² and pharmaceutical industries,³ and for the analysis of chiral compounds⁴ and asymmetric synthesis products⁵ and catalysts.⁶ Polarimetry is also used in the teaching laboratory.^{7,8}

While glucose and sucrose solutions are the most commonly utilized international standards for calibrating polarimeters,⁹ their use suffers from the following disadvantages:

- Their solutions can become easily contaminated with bacteria or fungi.
- Sucrose can hydrolyze to a mixture of glucose and fructose, with a resulting reversal of the optical rotation.¹⁰
- Mutarotation is observed for solutions of glucose, which give constant rotations only after some hours.¹¹



Thus, we thought it desirable to search for a new polarimetry standard, which ought to be stable, inexpensive, and easy to prepare and purify. We conducted preliminary tests on several candidates [(+)-tartaric acid, quinine sulfate, ephedrine, and (-)-menthol] and compared the results with those obtained with glucose and sucrose:

- Aqueous solutions of (+)-tartaric acid¹² and quinine sulfate racemized slowly and did not give stable readings.
- Aqueous solutions of ephedrine hydrochloride¹³ (generated in situ from the free base and hydrochloric acid) gave stable readings, but, after three months, these readings changed due to racemization.

In contrast, we obtained good results with (1*R*,2*S*,5*R*)-(-)-menthol,¹⁴ which may be obtained by resolving racemic synthetic menthol,⁸ by asymmetric synthesis,¹⁵ or, more commonly, from the essential oils of several species of mint, e.g., *Mentha piperita* or *Mentha herbensis*.¹⁶ It is also commercially available,¹⁷ inexpensive (much less so than ephedrine or quinine sulfate), and easy to purify by recrystallization or sublimation¹⁴—even easier than sucrose or glucose. We found ethanolic solutions of (-)-menthol to be stable for years.

We used samples of enantiomerically pure natural (-)-menthol¹⁹ of constant melting point (41–42 °C), and checked their chemical purity by the method of Ligor and Buszewski.^{19,20} The purities obtained fell in the range of >99.9 to 100%. We then measured the absolute optical rotation of a 10% solution of (-)-menthol in absolute ethanol on a JASCO DIP 370 digital polarimeter. 27 measurements were carried out at 25 °C over a period of 72 days. The average of these measurements, $[\alpha]_D^{25} = -50.40 \pm 0.01$ ($c = 10$, C₂H₅OH), compares very favorably with the reported value of $[\alpha]_D^{25} = -50$ ($c = 10$, C₂H₅OH).^{8,14,17}

We have been routinely using such an ethanolic solution of natural menthol (stored in a tightly closed vial) as a polarimeter calibration standard for over three years. We have obtained very stable readings during this period, and found no evidence of decomposition or loss by evaporation.

References and Notes

- (1) (a) Eliel, E. L.; Wilen, S. H.; Mander, L. N. *Stereochemistry of Organic Compounds*; John Wiley & Sons, Inc.: New York, 1994; pp 2–3. (b) The first practical polarimeter was constructed in 1840. Biot, J. B. *Ann. Chim. Phys.* **1840**, *74*, 401; <http://museum.nist.gov/panels/bates/intro.htm> (accessed May 14, 2002). (2) Ulber, R.; Faurie, R.; Sosnitzka, P.; Fischer, L.; Stärk, E.; Arbeck, C.; Scheper, T. *J. Chromatogr., A* **2000**, *882*, 329. (3) Maier, N. M.; Pilar, F.; Lindner, W. *J. Chromatogr., A* **2001**, *906*, 3. (4) Goodall, D. M. *Trac-Trend Anal. Chem.* **1993**, *12*, 177; *Chem. Abstr.* **1993**, *119*, 155032m. (5) Huerta, F. F.; Minidis, A. B. E.; Backvall, J. E. *Chem. Soc. Rev.* **2001**, *30*, 221. (6) Richards, C. J.; Locke, A. J. *Tetrahedron: Asymmetry* **1998**, *9*, 2377. (7) Mosher, M. D.; Kelly, C. O.; Mosher, M. W. *J. Chem. Educ.* **1996**, *73*, 567. (8) Poiré, C.; Rabiller, C.; Chon, C.; Hudhomme, P. *J. Chem. Educ.* **1996**, *73*, 93. (9) National Institute of Standards and Technology (NIST). *Standard Reference Materials*; SRM No. 17e and 917b. <http://www.srmcatalog.nist.gov/srmcatalog/tables/204-6.htm> (accessed May 8, 2002). (10) Eggleston, G.; Vercelotti, J. R. *J. Carbohydr. Chem.* **2000**, *19*, 1305. (11) Yamabe, S.; Ishikawa, T. *J. Org. Chem.* **1999**, *64*, 4519. (12) Dippy, J. F. J.; Hughes, S. R. C.; Rozanski, A. *J. Chem. Soc.* **1959**, 2492. (13) Goss, C. A.; Wilson, D. C.; Welsler, W. E. *Anal. Chem.* **1994**, *66*, 3093. (14) Perrin, D. D.; Armarego, W. L. F. *Purification of Laboratory Chemicals*, 3rd ed.; Pergamon Press: Oxford, U.K., 1988; p 215. (15) Milone, C.; Gangemi, C.; Neri, G.; Pistone, A.; Galvagno, S. *Appl. Catal. A, Gen.* **2000**, *199*, 239. (16) Veronese, P.; Li, X.; Niu, X. M.; Weller, S. C.; Bressan, R. A.; Hasegawa, P. M. *Plant Cell. Tiss. Organ Cult.* **2001**, *64*, 133. (17) For example, Aldrich Chemical Company. *Handbook of Fine Chemicals and Laboratory Equipment*, U.S. ed.; Milwaukee, WI, 2003–2004; p 1161, cat. no. M278-0. Also available at www.sigma-aldrich.com. (18) Brazilian standard pharmaceutical grade, which is equivalent to USP grade. The samples were obtained from Drogasil or Henrifarma. (19) Ligor, M.; Buszewski, B. *J. Chromatogr., A* **1999**, *847*, 161. (20) HP 5890 II gas chromatograph; HP5, cross-linked 5% methyl silicone capillary glass column (30 m x 0.32 mm x 0.25 μm); nitrogen carrier gas; the injector temperature (200 °C), the FID temperature (250 °C), and the heating program were the same as those reported by Ligor and Buszewski.


Alberto Federman Neto (Ph.D.),*^a Áurea Donizete Lanchote Borges,^b Fernando Costa Archanjo,^{a,b} and Alessandra Caramori Pellegrino^b

^a Universidade de São Paulo
Faculdade de Ciências Farmacêuticas de Ribeirão Preto
Av. Zeferino Vaz, sem num., Campus Universitário
14040-903, Ribeirão Preto, SP
Brasil

^b Universidade de São Paulo
Faculdade de Filosofia, Ciências e Letras de Ribeirão Preto
Depto. De Química
E-mail: albnfneto@icfrp.usp.br

Do you have an innovative shortcut or unique laboratory hint you'd like to share with your fellow chemists? If so, please send it to Aldrich (attn: Lab Notes, Aldrichimica Acta). For submitting your idea, you will receive a complimentary Aldrich periodic table poster (cat. no. **Z54,320-9**). If you publish your lab note, you will also receive an Aldrich periodic table mouse pad (cat. no. **Z54,323-3**). It is Teflon[®]-coated, 8½ in. x 11 in., with a full-color periodic table on the front. We reserve the right to retain all entries for future consideration.

Teflon is a registered trademark of E. I. du Pont de Nemours & Co., Inc.



Sigma-Aldrich Worldwide Locations

Argentina

Sigma-Aldrich de Argentina SA
Phone: 54 11 4556 1472
Fax: 54 11 4552 1698
E-mail: info-argentina@sial.com

Australia

Sigma-Aldrich Pty. Ltd.
Phone: 61 2 9841 0555
Fax: 61 2 9841 0500
E-mail: ausmail@sial.com

Austria

Sigma-Aldrich Handels GmbH
Phone: 43 1 605 8110
Fax: 43 1 605 8120
E-mail: sigma@sigma.co.at

Belgium

Sigma-Aldrich N.V./S.A.
Phone: 32 3 899 1301
Fax: 32 3 899 1311
E-mail: becustsv@eurnotes.sial.com

Brazil

Sigma-Aldrich Brasil Ltda.
Phone: 55 11 3732 3100
Fax: 55 11 3733 5151
E-mail: sigmabr@sial.com

Canada

Sigma-Aldrich Canada Ltd.
Phone: (905) 829-9500
Fax: (905) 829-9292
E-mail: canada@sial.com

China

Sigma-Aldrich China Inc.
Phone: 86 21 6386 2766
Fax: 86 21 6386 3966
E-mail: china@sial.com

Czech Republic

Sigma-Aldrich s.r.o.
Phone: 420 2 2176 1310
Fax: 420 2 2176 3300
E-mail: czecustsv@eurnotes.sial.com

Denmark

Sigma-Aldrich Denmark A/S
Phone: 45 43 565910
Fax: 45 43 565905
E-mail: denorder@eurnotes.sial.com

Finland

Sigma-Aldrich Finland YA-Kemia Oy
Phone: 358 9 350 9250
Fax: 358 9 350 92555
E-mail: finorder@eurnotes.sial.com

France

Sigma-Aldrich Chimie S.a.r.l.
Phone: 0600 21 14 08
Fax: 0600 03 10 52
E-mail: fradsv@eurnotes.sial.com

Germany

Sigma-Aldrich Chemie GmbH
Phone: 49 89 6513 0
Fax: 49 89 6513 1169
E-mail: deorders@eurnotes.sial.com

Greece

Sigma-Aldrich (O.M.) Ltd.
Phone: 30 10 994 8010
Fax: 30 10 994 3831
E-mail: grcustsv@eurnotes.sial.com

Hungary

Sigma-Aldrich Kft
Phone: 36 1 269 6474
Fax: 36 1 235 9050
E-mail: info@sigma.sial.hu

India

Sigma-Aldrich Foreign Holding Company
Phone: 91 80 852 4222
Fax: 91 80 852 4214
E-mail: india@sial.com

Ireland

Sigma-Aldrich Ireland Ltd.
Phone: 353 1 404 1900
Fax: 353 1 404 1910
E-mail: eicustsv@eurnotes.sial.com

Israel

Sigma-Aldrich Israel Ltd.
Phone: 972 8 948 4222
Fax: 972 8 948 4200
E-mail: hexter@sigma.co.il

Italy

Sigma-Aldrich S.r.l.
Phone: 39 02 33417 310
Fax: 39 02 38010 737
E-mail: itorder@eurnotes.sial.com

Japan

Sigma-Aldrich Japan K.K.
Phone: 81 3 5821 3111
Fax: 81 3 5821 3170

Korea

Sigma-Aldrich Korea
Phone: 82 31 329 9000
Fax: 82 31 329 9090
E-mail: hcjoo@sial.com

Malaysia

Sigma-Aldrich (M) Sdn. Bhd.
Phone: 60 3 5635 3321
Fax: 60 3 5635 4116
E-mail: sam@sial.com

Mexico

Sigma-Aldrich Quimica, S.A. de C.V.
Phone: 52 722 276 1600
Fax: 52 722 276 1601
E-mail: mexico@sial.com

Netherlands

Sigma-Aldrich Chemie B.V.
Phone: 31 78 620 5411
Fax: 31 78 620 5421
E-mail: nlcustsv@eurnotes.sial.com

Norway

Sigma-Aldrich Norway AS
Phone: 47 23 176 000
Fax: 47 23 176 010
E-mail: nortech@sial.com

Poland

Sigma-Aldrich Sp. z o.o.
Phone: 48 61 823 2481
Fax: 48 61 823 2781
E-mail: plcustsv@eurnotes.sial.com

Portugal

Sigma-Aldrich Quimica, S.A.
Phone: 351 1 924 2555
Fax: 351 1 924 2610
E-mail: poorders@eurnotes.sial.com

Russia

OOO SAF-LAB
Phone: 7 095 975 3321
Fax: 7 095 975 4792
E-mail: techcare@sial.com

Singapore

Sigma-Aldrich Pte., Ltd.
Phone: 65 6271 1089
Fax: 65 6271 1571
E-mail: sapl@sial.com

South Africa

Sigma-Aldrich (Pty) Ltd
Phone: 27 11 979 1188
Fax: 27 11 979 1119
E-mail: rsa@sial.com

Spain

Sigma-Aldrich Quimica S.A.
Phone: 34 91 661 9977
Fax: 34 91 661 9642
E-mail: esorders@eurnotes.sial.com

Sweden

Sigma-Aldrich Sweden AB
Phone: 46 8 742 4200
Fax: 46 8 742 4243
E-mail: sweorder@eurnotes.sial.com

Switzerland

Fluka Chemie GmbH
Phone: 41 81 755 2721
Fax: 41 81 756 2840
E-mail: fluka@sial.com

United Kingdom

Sigma-Aldrich Company Ltd
Phone: 44 1747 833000
Fax: 44 1747 833313
E-mail: ukcustsv@eurnotes.sial.com

United States

Aldrich
Phone: (414) 273-3850
Fax: (414) 273-4979
E-mail: aldrich@sial.com

Sigma

Phone: (314) 771-5765
Fax: (314) 771-5757
E-mail: sigma@sial.com

Acetaldehyde

Manufactured by Eastman Chemical Company

Custom Packaged by Sigma-Aldrich in 18-Liter
to 1000-Liter Returnable MINI-BULK™ Containers

Ask for product number **W20036-0** (FEMA # 2003)

Call Dan Zagrodnik or Ted Eslinger at 800-213-1206 (USA)

New Content

- ▶ Monthly new products updates
- ▶ Online glassware catalog
- ▶ Useful resources: pressure-temperature nomograph, wire gauge chart, common solvent properties, ring numbering charts
- ▶ Catalog of fine art pigments
- ▶ Library of Rare Chemicals updates
- ▶ Spectral Viewer software tool

Commerce Features

- ▶ Current pricing
- ▶ Stock status
- ▶ Order tracking
- ▶ Product shelf for frequently ordered products

Search Functions


- ▶ Search for products by name, molecular formula, CAS number, and substructure
- ▶ Coming soon: keyword site search

Technical Information

- ▶ Aldrichimica Acta magazine
- ▶ Technical information bulletins
- ▶ Technical brochures, including *Analytix* and *ChemFiles*

Product Information

- ▶ Physical properties
- ▶ Extended specifications
- ▶ Application descriptions
- ▶ Hazard listings
- ▶ MSDS's
- ▶ Certificates of analysis
- ▶ Structure images
- ▶ Literature references



Scientists all over the world regard the *Aldrich Handbook of Fine Chemicals and Laboratory Equipment* as an indispensable resource. What many people don't yet realize is that our website, www.sigma-aldrich.com/aldrich, contains all the information found in the Handbook, and more!

www.sigma-aldrich.com/aldrich

ALDRICH CHEMICAL COMPANY, INC.
P.O. BOX 14508
ST. LOUIS, MISSOURI 63178 USA



CHANGE SERVICE REQUESTED

PRESORTED STANDARD
U.S. POSTAGE
PAID
MILWAUKEE, WISCONSIN
PERMIT NO. 552

DEDICATED TO PROFESSOR D. SEEBACH ON HIS 65TH BIRTHDAY

Aldrichimica ACTA

VOL. 36, NO. 1 • 2003



Iron Lewis Acid
[$(\eta^5\text{-C}_5\text{H}_5)_2\text{Fe}(\text{CO})_2(\text{THF})$]
Catalyzed Organic Reactions


Applications of
Dialkylaminopyridine (DMAP)
Catalysts in Organic Synthesis



ALDRICH[®]

CHEMISTS HELPING CHEMISTS

new products




- 55,883-4** (R)-(+)-N-p-Tolylsulfinyltryptamine, 97%
1g \$32.00; 5g \$106.65
- 58,851-2** (S)-(-)-N-p-Tolylsulfinyltryptamine, 97%
1g \$32.00; 5g \$106.65

These chiral tryptamine-derived sulfinamines can be condensed with aldehydes to form N-sulfinyliminium ions with high diastereomeric ratios. These ions undergo Picket-Spengler cyclization to generate biologically active tetrahydro- β -carbolines.

Gremmen, C. et al. *Org. Lett.* **2000**, *2*, 1955.


- 57,926-2** 3,5-Dimethoxyphenol trifluoromethanesulfonate, 97%
5g \$44.40; 25g \$148.00



This triflate-activated dimethoxyphenol readily undergoes Suzuki cross-couplings under mild conditions. It has been utilized in the synthesis of naturally occurring 5-alkylresorcinols¹ and Δ^8 -tetrahydrocannabinol analogs.²

(1) Fürstner, A.; Seidel, G. *J. Org. Chem.* **1997**, *62*, 2332. (2) Crocker, P. J. et al. *Tetrahedron* **1999**, *55*, 13907.


- 59,276-5** N-Boc-thiourea 1g \$45.90; 5g \$153.00



This monoprotected thiourea has been used in the synthesis of 2-aminothiazole intermediates for the study of marine metabolites.

Schiavi, B. et al. *Synth. Commun.* **2002**, *32*, 1671.


- 57,108-3** 1-Benzoylpiperazine, 97% 1g \$22.15; 5g \$73.80



Employed in the synthesis of various biologically active compounds such as nitric oxide donors,¹ pyruvate dehydrogenase kinase inhibitors,² and serine protease factor Xa inhibitors.³

(1) Mu, L. et al. *Chem. Pharm. Bull.* **2000**, *48*, 808. (2) Aicher, T. D. et al. *J. Med. Chem.* **2000**, *43*, 236. (3) Jones, S. D. et al. *Bioorg. Med. Chem. Lett.* **2001**, *11*, 773.


- 52,285-6** 5-Bromo-2,2'-bithiophene, 96%
1g \$40.90; 5g \$171.20



Building block for the synthesis of soluble thiophene oligomers possessing electronic and electro-optical properties.

Sotgiu, G. et al. *Tetrahedron* **2002**, *58*, 2245.


- 15,293-5** 1-Piperidinepropanol, 97% 1g \$11.00; 5g \$38.00



Employed in the preparation of new nonimidazole histamine H₃ receptor ligands designed to increase histaminergic neurotransmission.

Apelt, J. et al. *J. Med. Chem.* **2002**, *45*, 1128.


- 58,928-4** 1-[(Trimethylsilyl)methyl]benzotriazole, 99%
1g \$16.65; 5g \$55.50



Reagent for the one-carbon homologation of acyl chlorides to the corresponding acids or esters.

Katritzky, A. R. et al. *J. Org. Chem.* **2001**, *66*, 5606.
Katritzky, A. R. et al. *ibid.* **2002**, *67*, 7526.


- 58,885-7** 3-(Ethoxydimethylsilyl)propylamine, 97%
1g \$39.75; 5g \$132.40



Was utilized to develop surface tension arrays for the synthesis of oligonucleotide arrays on glass surfaces.

Butler, J. H. et al. *J. Am. Chem. Soc.* **2001**, *123*, 8887.


- 57,896-7** 1-Hydroxybenzotriazole hydrate, 15 wt. % solution in dimethylformamide
25g \$28.75; 100g \$79.80



A key reagent for peptide¹ and polyamide² synthesis.

(1) *Encyclopedia of Reagents for Organic Synthesis*; Paquette, L. A., Ed.; Wiley: West Sussex, U.K., 1995; pp 2752-2755. (2) Xiao, J. et al. *J. Org. Chem.* **2000**, *65*, 5506.


- 55,414-6** (6-Bromopyridin-2-yl)methanol, 96%
1g \$38.25; 5g \$127.40



Has been used in the synthesis of oligopyridine-functionalized aza-crown ethers¹ and 4-pyridyl-1,4-dihydropyridines.²

(1) Tsukube, H. et al. *J. Org. Chem.* **1993**, *58*, 4389.
(2) Ashimori, A. et al. *Chem. Pharm. Bull.* **1990**, *38*, 2446.


- 54,944-5** 1,4-Phenylenebis(chlorodimethylsilane), 95%
1g \$13.90; 5g \$46.25



This organosilane was utilized in the synthesis of zirconocene polymers and macrocycles.

Mao, S. S. H. et al. *J. Am. Chem. Soc.* **1998**, *120*, 1193.

- 15,942-5** 9-Chloroacridine 1g \$35.00; 5g \$120.00



A main building block for many biologically active compounds such as anti-inflammatory,¹ antibacterial,¹ and antimalarial² substances.

(1) Yeh-Long, C. et al. *J. Med. Chem.* **2002**, *45*, 4689.
(2) Kimura, M.; Okabayashi, I. *J. Heterocycl. Chem.* **1986**, *23*, 965.

Please see pages 30-31 for additional new products.

For competitive quotes on larger quantities contact www.sigma-aldrich.com/safc.

SIGMA-ALDRICH
FINE CHEMICALS

Aldrichimica ACTA

VOL. 36, NO. 1 • 2003

Aldrich Chemical Co., Inc.

1001 West Saint Paul Ave.
Milwaukee, WI 53233 USA

To Place Orders

Telephone 800-558-9160 (USA)
or 414-273-3850
FAX 800-962-9591 (USA)
or 414-273-4979
Mail P.O. Box 2060
Milwaukee, WI 53201 USA

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201 USA

Customer & Technical Services

Customer Inquiries 800-558-9160
Technical Service 800-231-8327
Sigma-Aldrich Fine Chemicals 800-336-9719
Custom Synthesis 800-336-9719
Flavors & Fragrances 800-227-4563
International 414-273-3850
24-Hour Emergency 414-273-3850
Website <http://www.sigma-aldrich.com>
Email aldrich@sial.com

To request your **FREE** subscription to the *Aldrichimica Acta*, please

call: 800-558-9160 (USA),
write: Attn: Mailroom
Aldrich Chemical Co., Inc.
P.O. Box 355
Milwaukee, WI 53201-9358
or email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office.

The *Aldrichimica Acta* is also available on the Internet at <http://www.sigma-aldrich.com>.

All prices listed in this publication are subject to change without notice.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

For worldwide contact information, please see the inside back cover.

Aldrichimica Acta is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich family. © 2003 by Sigma-Aldrich Co. Printed in the USA.

About Our Cover

Landscape at Vétheuil (oil on canvas, 115 x 165 cm) was painted by the French impressionist painter Auguste Renoir around 1890. Vétheuil, on the river Seine northwest of Paris, was chosen by Renoir's fellow impressionist Claude Monet as a location sufficiently distant from Paris to avoid urban distractions and allow concentration on painting. The landscape around Vétheuil, like that near Argenteuil and other popular boating and picnicking spots along the Seine, was an important source of inspiration for a number of these artists.



Photograph © Board of Trustees, National Gallery of Art, Washington.

The impressionists were not actually a homogeneous group united by clearly defined principles. Some can be called impressionists only at certain times during their careers, and among the group there is considerable ideological and stylistic variety. Generally speaking, however, they were unhappy with academic training, rejected the notion central to romanticism that a work of art should convey emotional content, and were uninterested in the message of social realism. They felt that the proper purpose of a work of art is not to be a vehicle for the artist's imagination or to represent historical or literary subjects, but to record dispassionately the actuality of nature and contemporary life.


Renoir's *Landscape at Vétheuil* is characteristic of impressionist paintings in the artist's rejection of a balanced and contrived composition in favor of what might be called a snapshot view, and in his use of a loose technique of small strokes of paint to convey the immediacy of visual experience. The irony is that this quasiscientific approach is often misinterpreted by modern viewers, who see in impressionism a romantic and idealized vision of the world, rather than the direct record of experience it was meant to be.

This painting is in the Ailsa Mellon Bruce Collection at the National Gallery of Art, Washington, DC.

“Please
Bother
Us.”

Joe Porwoll

Joe Porwoll, President



Professor Philip Eaton of the Department of Chemistry at The University of Chicago kindly suggested that we market a solution of butylmagnesium diisopropylamide. This thermally stable reagent deprotonates cyclopropane carboxamides and weak CH acids stoichiometrically and regioselectively.

Zhang, M.-X.; Eaton, P. E. *Angew. Chem., Int. Ed.* **2002**, *41*, 2169

59,045-2

Butylmagnesium diisopropylamide.

1.0M solution in tetrahydrofuran

25mL \$21.10; 100mL \$58.50

Naturally, we made this valuable reagent. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page or on the inside back cover.

Table of Contents

Iron Lewis Acid [(η⁵-C₅H₅)Fe(CO)₂(THF)]
Catalyzed Organic Reactions 3
—M. D. Redlich, M. F. Mayer,
and M. M. Hossain³

Applications of
Dialkylaminopyridine (DMAP)
Catalysts in Organic Synthesis 21
—R. Murugan and E. F. V. Scriven³




Lab Notes

Purification of 1 α -Bromotriacetylglucuronate Methyl Ester


Glycoside linkages are most commonly synthesized by the Koenigs–Knorr reaction,¹ through coupling of an acetylated 1-bromosugar with an alcohol or phenol group in the presence of a salt such as silver carbonate. The key to the success of this reaction is the purity of the bromosugar, which is prepared by reacting the fully acetylated and esterified sugar with an acetic acid solution of hydrobromic acid. The original method for preparing 1 α -bromotriacetylglucuronate methyl ester recommended recrystallizing the product from absolute alcohol.² Over the years, observations in our labs have led us to the realization that residual alcohol in the purified bromosugar generates hydrobromic acid and causes the rapid degradation of the crystals.³ In fact, after just a few hours at room temperature, the material rapidly turns brown. After 1–3 days, the material is definitely degraded, turning an oily black, and must be recrystallized before it can be used. This degradation is so rapid that the material is not widely available commercially.⁴ Even when stored under anhydrous conditions in the freezer, the material has a relatively short life span.

Difficulties with maintaining large quantities of the pure bromosugar have led us to search for a better recrystallization solvent. The requisite solvent(s) must be anhydrous and not contain reactive groups such as hydroxyls. We found that benzene–hexane or, preferably, toluene–hexane meets these requirements. The crude bromosugar is dissolved in warm (not boiling) toluene and briefly treated with decolorizing charcoal. The solution is filtered while warm through a sintered glass funnel (do not use filter paper) containing a Celite[®] pad, and hexane is added with stirring to the colorless filtrate until it begins to turn cloudy. The white crystals, which form after refrigeration overnight, are collected on a sintered glass funnel and thoroughly dried under high vacuum at room temperature. Material prepared in this fashion may be stored under anhydrous conditions in the freezer for an indefinite period, and is even stable at room temperature under anhydrous conditions or under vacuum for several days.

Purity can be readily determined visually (white to off-white) or by odor, since HBr can be easily detected in the decomposed material. The melting point must be carefully determined. When rapidly heated, the material melts at 82–83 °C, but, if allowed to cool and resolidify and is then heated slowly, the sample melts with decomposition at 102–104 °C. A slowly heated sample also melts with decomposition at 102–104 °C. Infrared and NMR analyses are also useful for purity determination. The IR (KBr) spectrum has two sharp carbonyl absorptions at 1767 and 1748 cm⁻¹ and a set of C–O absorptions centered at 1229 cm⁻¹. The ¹H NMR (CDCl₃, 60 MHz) spectrum consists of overlapping acetate methyls centered at 2.05 ppm (9 H), the ester methoxy at 3.77 ppm (3 H), the C-2, C-3, C-4, and C-5 ring hydrogens between 4.5 and 5.8 ppm (4 H), and the C-1 hydrogen centered at 6.65 ppm (1 H, doublet).⁵



IR scan



NMR scan

References and Notes: (1) Conrow, R. B.; Bernstein, S. *J. Org. Chem.* **1971**, *36*, 863. (2) Bollenback, G. N.; Long, J. W.; Benjamin, D. G.; Lindquist, J. A. *J. Am. Chem. Soc.* **1955**, *77*, 3310. (3) Hadd, H. E.; Slikker, W., Jr.; Helton, E. D. *J. Steroid Biochem.* **1980**, *13*, 1107. (4) 1 α -Bromotriacetylglucuronate methyl ester (CAS Registry NumberSM 21085-72-3) is currently available from a very small number of companies, including the Sigma (cat. no. A8292) and Fluka (cat. no. 00533F) brands of Sigma-Aldrich Corp., and is shipped in dry ice. (5) The resonances at 7.27 and 7.37 ppm arise from a trace of CHCl₃ (from CDCl₃) and benzene (recrystallization solvent), respectively. While it is difficult to pump off all of the aromatic solvent, it should not affect further use of the material in the Koenigs–Knorr reaction, since this reaction is most often carried out in benzene or toluene.


Celite is a registered trademark of Celite Corp. CAS Registry Number is a registered service mark of the American Chemical Society.

John M. Kokosa^{*}, Professor
Department of Science and Mathematics
Kettering University
1700 W. Third Avenue
Flint, MI 48504
Email: jkokosa@kettering.edu

Harry E. Hadd
Associate Professor Emeritus (Deceased)
Department of Biochemistry and Molecular Biology
Indiana University School of Medicine, Northwest Center
3400 Broadway
Gary, IN 46408

Do you have an innovative shortcut or unique laboratory hint you'd like to share with your fellow chemists? If so, please send it to Aldrich (attn: Lab Notes, Aldrichimica Acta). For submitting your idea, you will receive a complimentary, Aldrich periodic table poster (cat. no. **Z54,320-9**). If we publish your lab note, you will also receive an Aldrich periodic table mouse pad (cat. no. **Z54,323-3**). It is Teflon[®]-coated, 8½ in. x 11 in., with a full-color periodic table on the front. We reserve the right to retain all entries for future consideration.

Teflon is a registered trademark of E.I. du Pont de Nemours & Co., Inc



Sigma-Aldrich Worldwide Locations

Argentina

Sigma-Aldrich de Argentina SA
Phone: 54 11 4556 1472
Fax: 54 11 4552 1698
Email: info-argentina@sial.com

Australia

Sigma-Aldrich Pty. Ltd.
Phone: 61 2 9841 0555
Fax: 61 2 9841 0500
Email: ausmail@sial.com

Austria

Sigma-Aldrich Handels GmbH
Phone: 43 1 605 8110
Fax: 43 1 605 8120
Email: sigma@sigma.co.at

Belgium

Sigma-Aldrich N.V./S.A.
Phone: 32 3 899 1301
Fax: 32 3 899 1311
Email: becustsv@eurnotes.sial.com

Brazil

Sigma-Aldrich Brasil Ltda.
Phone: 55 11 3732 3100
Fax: 55 11 3733 5151
Email: sigmabr@sial.com

Canada

Sigma-Aldrich Canada Ltd.
Phone: (905) 829-9500
Fax: (905) 829-9292
Email: canada@sial.com

China

Sigma-Aldrich China Inc.
Phone: 86 21 6386 2766
Fax: 86 21 6386 3966
Email: china@sial.com

Czech Republic

Sigma-Aldrich s.r.o.
Phone: 420 2 2176 1310
Fax: 420 2 2176 3300
Email: czecustsv@eurnotes.sial.com

Denmark

Sigma-Aldrich Denmark A/S
Phone: 45 43 565910
Fax: 45 43 565905
Email: denorder@eurnotes.sial.com

Finland

Sigma-Aldrich Finland YA-Kemia Oy
Phone: 358 9 350 9250
Fax: 358 9 350 92555
Email: finorder@eurnotes.sial.com

France

Sigma-Aldrich Chimie S.a.r.l.
Phone: 0800 21 14 08
Fax: 0800 03 10 52
Email: fradvsv@eurnotes.sial.com

Germany

Sigma-Aldrich Chemie GmbH
Phone: 49 89 6513 0
Fax: 49 89 6513 1169
Email: deorders@eurnotes.sial.com

Greece

Sigma-Aldrich (O.M.) Ltd.
Phone: 30 10 994 8010
Fax: 30 10 994 3831
Email: grcustsv@eurnotes.sial.com

Hungary

Sigma-Aldrich Kft
Phone: 36 1 269 6474
Fax: 36 1 235 9050
Email: info@sigma.sial.hu

India

Sigma-Aldrich Foreign Holding Company
Phone: 91 80 852 4222
Fax: 91 80 852 4214
Email: india@sial.com

Ireland

Sigma-Aldrich Ireland Ltd.
Phone: 353 1 404 1900
Fax: 353 1 404 1910
Email: eicustsv@eurnotes.sial.com

Israel

Sigma-Aldrich Israel Ltd.
Phone: 972 8 948 4222
Fax: 972 8 948 4200
Email: sigisr@sigma.co.il

Italy

Sigma-Aldrich S.r.l.
Phone: 39 02 33417 310
Fax: 39 02 38010 737
Email: itorder@eurnotes.sial.com

Japan

Sigma-Aldrich Japan K.K.
Phone: 81 3 5821 3111
Fax: 81 3 5821 3170
Email: sialjp@sial.com

Korea

Sigma-Aldrich Korea
Phone: 82 31 329 9000
Fax: 82 31 329 9090
Email: hcjoo@sial.com

Malaysia

Sigma-Aldrich (M) Sdn. Bhd.
Phone: 60 3 5635 3321
Fax: 60 3 5635 4116
Email: sam@sial.com

Mexico

Sigma-Aldrich Quimica, S.A. de C.V.
Phone: 52 722 276 1600
Fax: 52 722 276 1601
Email: mexico@sial.com

Netherlands

Sigma-Aldrich Chemie B.V.
Phone: 31 78 620 5411
Fax: 31 78 620 5421
Email: nlcustsv@eurnotes.sial.com

Norway

Sigma-Aldrich Norway AS
Phone: 47 23 176 000
Fax: 47 23 176 010
Email: nortech@sial.com

Poland

Sigma-Aldrich Sp. z o.o.
Phone: 48 61 823 2481
Fax: 48 61 823 2781
Email: plcustsv@eurnotes.sial.com

Portugal

Sigma-Aldrich Quimica, S.A.
Phone: 351 1 924 2555
Fax: 351 1 924 2610
Email: poorders@eurnotes.sial.com

Russia

OOO SAF-LAB
Phone: 7 095 975 3321
Fax: 7 095 975 4792
Email: techcare@sial.com

Singapore

Sigma-Aldrich Pte., Ltd.
Phone: 65 6271 1089
Fax: 65 6271 1571
Email: sapl@sial.com

South Africa

Sigma-Aldrich (Pty.) Ltd
Phone: 27 11 979 1188
Fax: 27 11 979 1119
Email: rsa@sial.com

Spain

Sigma-Aldrich Quimica S.A.
Phone: 34 91 661 9977
Fax: 34 91 661 9642
Email: esorders@europe.sial.com

Sweden

Sigma-Aldrich Sweden AB
Phone: 46 8 742 4200
Fax: 46 8 742 4243
Email: sweorder@eurnotes.sial.com

Switzerland

Fluka Chemie GmbH
Phone: 41 81 755 2721
Fax: 41 81 756 2840
Email: fluka@sial.com

United Kingdom

Sigma-Aldrich Company Ltd.
Phone: 44 1747 833000
Fax: 44 1747 833313
Email: ukcustsv@eurnotes.sial.com

United States

Aldrich
Phone: (414) 273-3850
Fax: (414) 273-4979
Email: aldrich@sial.com

Sigma


Phone: (314) 771-5765
Fax: (314) 771-5757
Email: sigma@sial.com

Chromatography Columns

- Low-pressure (max. 7 psig), medium-wall columns with reservoir
- Greaseless PTFE insert with screw-thread attached flow control valve to column
- Precision all-PTFE valve and drip-tip

I.D. x L (mm)	Vol. (mL)	Cat. No.	Each	With fritted glass disc	
				Cat. No.	Each
11 x 250	22	Z40,058-0	172.60	Z40,063-7	174.40
15 x 250	40	Z40,059-9	178.00	Z40,064-5	180.20
19 x 300	80	Z40,060-2	179.10	Z40,065-3	193.00
25 x 300	150	Z40,061-0	182.60	Z40,066-1	195.00
50 x 300	560	Z40,062-9	198.70	Z40,068-8	201.70





"It is gratifying to have so much useful and informative material... concentrated in one highly readable, easily accessible source."

Journal of Medicinal Chemistry
Book Reviews
 (2002, 45, 992-993)

The Sigma-RBI Handbook of Receptor Classification and Signal Transduction, 4th edition, is a highly respected and widely used 231-page reference book that contains over 90 separate entries on a range of receptors, ion channels and enzymes involved in cell signaling mechanisms.

Over 150 respected research scientists worked with Sigma-RBI to compile the handbook. A peer-reviewed overview is provided for each entry,

along with a chart summarizing a variety of information including over 1,200 key research tools that are available from Sigma-RBI.

Order your copy today and see for yourself why scientists depend on **The Sigma-RBI Handbook of Receptor Classification and Signal Transduction** (Prod. No. S 8188). Call 800-325-3010 (USA), Fax 800-325-5052 (USA), www.sigma-aldrich.com, or contact your local Sigma-Aldrich sales office.

availability, ordering & tracking
sigma-aldrich.com
 scientific tools, products & information

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
 SIGMA-ALDRICH CORPORATION • BOX 14508 • ST. LOUIS, MO 63178 USA

SIGMA
RBI

ALDRICH CHEMICAL COMPANY, INC.
 P.O. BOX 14508
 ST. LOUIS, MISSOURI 63178 USA



CHANGE SERVICE REQUESTED

PRESORTED STANDARD
 U.S. POSTAGE
 PAID
 MILWAUKEE, WISCONSIN
 PERMIT NO. 552

SPECIAL FOCUS ON NONRACEMIC AZIRIDINES AND OXAZOLINES

Aldrichimica ACTA

VOL. 36, NO. 2 • 2003



*Aziridines and Oxazolines:
Valuable Intermediates in the
Synthesis of Unusual Amino Acids*


**Highlights of the Chemistry
of Enantiomerically Pure
Aziridine-2-carboxylates**

sigma-aldrich.com


 **ALDRICH**[®]
CHEMISTS HELPING CHEMISTS

New Products from Aldrich R&D

Diethyl *trans*-cinnamylphosphonate, 98%

59,438-5  5g
25g

Diethyl (2-methylallyl)phosphonate, 97%

59,309-5  1g

These phosphonates were utilized in the Horner–Wadsworth–Emmons reaction to form conjugated carbon–carbon double bonds.¹ They were also employed as starting materials in an efficient and regioselective synthesis of 4-oxo-2-alkenylphosphonates, which can serve as building blocks for the construction of polyethylenic chains.²

(1) Oestreich, M.; Hoppe, D. *Tetrahedron Lett.* **1999**, *40*, 1881. (2) Lee, B. S. et al. *J. Org. Chem.* **2000**, *65*, 4175.


Ethyl [Bis(2,2,2-trifluoroethoxy)phosphinyl]acetate

59,557-8  5g
10g

This compound was exploited in the Horner–Wadsworth–Emmons reaction to synthesize α,β -unsaturated esters derived from 6-methoxytetrahydropyran-3-one.

López Tudanca, P. L. et al. *J. Chem. Soc., Perkin Trans. 1* **1992**, 533.

2,6-Dichloropyridine-1-oxide, 99%

59,405-9  1g
5g

It oxidizes alkenes to epoxides¹ and alkanes to alcohols² in the presence of ruthenium catalysts.

(1) Zhang, J.-L.; Che, C.-M. *Org. Lett.* **2002**, *4*, 1911. (2) Yamaguchi, M. et al. *Chem. Lett.* **2002**, 434.

1,5-Naphthyridine hydrochloride

59,416-4  1g
5g

Serves as a precursor of diaza-*cis*-decalins, a structurally novel class of diamine ligands.¹ Has also been used in the synthesis of one member of a series of antimicrobial parenteral 3'-quaternary ammonium cephalosporins.²

(1) Li, X. et al. *Org. Lett.* **2000**, *2*, 875. (2) Brown, R. F. et al. *J. Med. Chem.* **1990**, *33*, 2114.


1-(1,1-Dimethylheptyl)-3,5-dimethoxybenzene, 97%

59,522-5  1g
5g


Employed as a starting material in the synthesis of a number of THC analogs that were evaluated for their binding affinity towards cannabinoid receptors.

Gareau, Y. et al. *Bioorg. Med. Chem. Lett.* **1996**, *6*, 189.


2,3-Dibromo-*N*-methylmaleimide

59,593-4  1g
5g

2,3-Dibromomaleimide, 97%

55,360-3  1g
5g


N-Benzyl-2,3-dibromomaleimide, 97%

55,778-1  1g
5g

Dihalogenated maleimides can be used either as dienophiles or as electrophiles. 2,3-Dibromo-*N*-methylmaleimide is a key starting material in the synthesis of rebeccamycin¹ and 7-azarebeccamycin analogs.² These analogs were then evaluated for their antitumor activities.

Marminon, C. et al. *Bioorg. Med. Chem.* **2003**, *11*, 679. (2) Marminon, C. et al. *J. Med. Chem.* **2003**, *46*, 609.


Tris[(methylamino)ethyl]amine, 97%

46,353-1  5g
10g
25g

A tripodal metal chelating agent that has been employed in the preparation of *N*-methyl superbase (Aldrich Cat. No. 46,355-8),¹ and its stilbene and bismuth azaatran analogs:² *N,N',N''*-trimethylazastibatrane and *N,N',N''*-trimethylazabismatrane.

(1) Tang, J.-s.; Verkade, J. G. *Tetrahedron Lett.* **1993**, *34*, 2903. (2) Shutov, P. L. et al. *Inorg. Chem.* **2002**, *41*, 6147.

N,N-Diethyl-1,1-dimethylsilylamine, 97%

58,624-2  1g
5g

Complements NaBH_3CN , and has been used in the Lewis acid catalyzed reductive amination of carbonyl compounds.

Miura, K. et al. *Synlett* **2001**, 1617.

Please see pages 54–55 for additional new products.

For competitive quotes on larger quantities contact www.sigma-aldrich.com/safc.

SIGMA-ALDRICH
FINE CHEMICALS

Aldrichimica ACTA

VOL. 36, NO. 2 • 2003

Aldrich Chemical Co., Inc.

1001 West Saint Paul Ave.
Milwaukee, WI 53233 USA

To Place Orders

Telephone 800-558-9160 (USA)
or 414-273-3850
FAX 800-962-9591 (USA)
or 414-273-4979
Mail P.O. Box 2060
Milwaukee, WI 53201 USA

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201 USA

Customer & Technical Services

Customer Inquiries	800-558-9160
Technical Service	800-231-8327
Sigma-Aldrich Fine Chemicals	800-336-9719
Custom Synthesis	800-336-9719
Flavors & Fragrances	800-227-4563
International	414-273-3850
24-Hour Emergency	414-273-3850
Website	http://www.sigma-aldrich.com
Email	aldrich@sial.com

To request your FREE subscription to the *Aldrichimica Acta*, please

call: 800-558-9160 (USA),
write: Attn: Mailroom
Aldrich Chemical Co., Inc.
P.O. Box 355
Milwaukee, WI 53201-9358
or email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office.

The *Aldrichimica Acta* is also available on the Internet at <http://www.sigma-aldrich.com>.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

For worldwide contact information, please see the inside back cover.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich family. © 2003 by Sigma-Aldrich Co. Printed in the USA.

About Our Cover

The Railway (oil on canvas, 93.3 x 111.5 cm) was painted in 1873 by the French painter Edouard Manet. When it was exhibited in the following year, it was severely criticized by both the critics and the public, who were greatly puzzled by the subject of the picture, or rather by the fact that it seemed to have no real subject. In late 19th-century France, the most highly valued subjects in art were religious, mythological, historical, or literary. At the same time, a contrary naturalistic movement, paralleled in literature by the writings of reformist authors such as Émile Zola, favored subjects that portrayed the lower classes, like scenes of peasants working in the fields.

Manet's picture, however, does not represent an imaginary literary subject or a glorious historical event, nor does it portray the travails of the poor or idealize the dignity of manual labor. It simply shows a young woman, who is neither rich nor poor, pausing to rest on a bench with a puppy in her lap, accompanied by a little girl with her back to us who grasps the bars of an iron fence. It does not even seem to show what is indicated by the title of the painting, and the only clue to this is the steam rising in the background.

The clear outdoor light and bright color and the broad brushstrokes of his technique seem to link Manet with the impressionists, and indeed, a year after he painted *The Railway*, he was at Argenteuil painting in the open air alongside Renoir and Monet. Manet is not truly an impressionist painter, however. This is not a quickly executed representation of a chance moment, captured by the painter as a photographer makes a snapshot, but a carefully planned work. Manet sketched in the basic composition before carrying the canvas outdoors to work directly from the models. Such details as the placement of one figure facing out and the other into the picture, and the color scheme of the dresses, one white on blue and the other blue on white, show the calculation that underlies his representation of *The Railway*, a phenomenon common to modern life.

This painting is a gift of Horace Havemeyer to the National Gallery of Art, Washington, DC, in memory of his mother, Louise W. Havemeyer.




Photograph © Board of Trustees, National Gallery of Art, Washington

“Please
Bother
Us.”

Joe Porwoll

Joe Porwoll, President



R = Me, *i*-Pr, *t*-Bu

Professor John G. Verkade of the Department of Chemistry at Iowa State University kindly suggested that we provide the following three proaza-phosphatrane nonionic bases. This family of superbases has broad applications¹ including recently as ligands in the Pd-catalyzed amination of aryl bromides and iodides.²

(1) Wroblewski, A. E.; Bansal, V.; Kisanga, P.; Verkade, J. G. *Tetrahedron* **2003**, 59, 561. (2) Urganekar, S.; Nagarajan, M.; Verkade, J. G. *J. Org. Chem.* **2003**, 68, 452

46,355-8
2,8,9-Trimethyl-2,5,8,9-tetraaza-1-phosphabicyclo[3.3.3]undecane

55,695-5
2,8,9-Triisopropyl-2,5,8,9-tetraaza-1-phosphabicyclo[3.3.3]undecane


56,588-1
2,8,9-Triisobutyl-2,5,8,9-tetraaza-1-phosphabicyclo[3.3.3]undecane

Naturally, we made these valuable superbases. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page or on the inside back cover.

Table of Contents

Aziridines and Oxazolines: Valuable Intermediates in the Synthesis of Unusual Amino Acids.....	39
—G. Cardillo* L. Gentilucci, and A. Tolomelli	
Highlights of the Chemistry of Enantiomerically Pure Aziridine-2-carboxylates.....	57
—W. K. Lee and H.-J. Ha	



Lab Notes



Use of the Microwave Oven for Sublimation: Flash Sublimation

Sublimation is a useful technique for the purification¹⁻³ or isolation^{2,4} of some organic, inorganic, or organometallic compounds. Generally, if a compound can be sublimed, sublimation can be a good alternative to recrystallization or distillation. Sublimation has been known since alchemical times and, in the past, was carried out by simply heating the compound in a porcelain dish covered with a common filter paper.⁵ Nowadays, a sublimation apparatus or, sometimes, a Kugelrohr oven^{3,4} is used under ambient or reduced pressure.

We have recently developed an improved method for the synthesis of arenocenium salts using a simple assembly for reactions under microwave conditions.^{6a} It consists of a crystallizing dish and a 250-mL, tall beaker (Berzelius flask) that is covered with dry ice. Dry ice does not absorb microwaves and, therefore, does not vaporize under microwave irradiation conditions.^{6b} We found that this simple device may also be used for sublimations. Microwave sublimation has been utilized to manufacture and isolate carbon nanotubes⁷ and essential powders from fresh animal, plant, or microbial matter.⁸

We have carried out the sublimation, under microwave heating, of some representative inorganic, organometallic, and organic compounds in the apparatus shown here. The sublimations were fast and easy to carry out. Collection of the sublimate with a spatula was also straightforward. The compounds tested and the "yields" of the corresponding sublimes are presented in **Table 1**. Even certain slightly air-sensitive compounds (Table 1, entries 7, 8, and 10), that are generally purified by sublimation under reduced pressure, may be purified by this method.

Acetyl ferrocene,⁹ decadeuterioferrocene,⁴ and (cyclopentadienyl)manganese tricarbonyl¹ were prepared by published procedures. Bromopentacarbonylmanganese was prepared by reaction of dimanganese decacarbonyl (Aldrich Cat. No. 24,526-7) and bromine. We tested all the recommended solvents for this reaction: CS₂,¹⁰ dichloromethane,¹⁰ carbon tetrachloride,¹¹ and hexane (used for the rhenium analog¹²), but found that benzene¹³ was the best solvent. Mn(CO)₅Br was obtained in 96% yield, in practically pure form, without formation of manganese(II) bromide as side product.¹¹ All other compounds in Table 1 were obtained from commercial sources.



References: (1) *Purification of Laboratory Chemicals*, 4th ed.; Perrin, D. D., Armarego, W. L., Eds.; Butterworth Publishers: New York, 1996. (2) Verberne, M. C.; Brouwer, N.; Delbianco, F.; Linthorst, H. J. M.; Bol, J. F.; Verpoorte, R. *Phytochem. Anal.* **2002**, *13*, 45. (3) Federman Neto, A.; Borges, A. D. L.; Miller, J.; Darin, V. A. *Synth. React. Inorg. Met.-Org. Chem.* **1997**, *27*, 1299. (4) Federman Neto, A.; Borges, A. D. L.; de Arruda Campos, I. P.; Miller, J. *Synth. React. Inorg. Met.-Org. Chem.* **1997**, *27*, 1543. (5) (a) Engel, R. *Traité Élémentaire de Chimie*; Librairie J.-B. Baillière et Fils (Publisher): Paris, France, 1896; p 13. (b) Furniss, B. S.; Hannaford, A. J.; Smith, P. W. G.; Tatchell, A. R. *Vogel's Textbook of Practical Organic Chemistry*, 5th ed.; Longman Scientific & Technical: Harlow, U.K., 1989; pp 153-155. (6) (a) Federman Neto, A. *An. Acad. Bras. Cienc.* **2003**, in press. (b) Dabirmanesh, Q.; Fernando, S. I. S.; Roberts, R. M. G. *J. Chem. Soc., Perkin Trans. 1* **1995**, 743. (7) Mochizuki, T.; Yoshizawa, H. *Jpn. Patent 2000 272,913*, October 3, 2000; *Chem. Abstr.* **2000**, *133*, 254543z. (8) Maghami, P. *Fr. Patent 2,618,450*, January 27, 1989; *Chem. Abstr.* **1989**, *111*, 150119v. (9) Darin, V. A.; Federman Neto, A.; Miller, J.; de Freitas Afonso, M. M.; Fonsatti, H. C.; Borges, A. D. L. *J. Prakt. Chem.* **1999**, *341*, 588. (10) Quick, M. H.; Angelici, R. J. *Inorg. Synth.* **1979**, *19*, 160. (11) Abel, E. W.; Wilkinson, G. *J. Chem. Soc.* **1959**, 1501. (12) Schmidt, S. P.; Troglor, W. C.; Basolo, F.; Urbancic, M. A.; Shapley, J. R. *Inorg. Synth.* **1985**, *23*, 41. (13) Federman Neto, A.; Oliveira, W. J. S.; Borges, A. D. L. University of São Paulo, Ribeirão Preto, Brazil. Unpublished work, 2003.

Alberto Federman Neto (Ph.D.),* Paloma Los Angeles Gazola Cordo, Wagner José Dos Santos Oliveira, and Aurea Donizete Lanchote Borges

Universidade de São Paulo
Faculdade de Ciências Farmacêuticas de Ribeirão Preto
Av. Zeferino Vaz, sem num., Campus Universitário
14040-903, Ribeirão Preto, SP
Brasil
Email: albfneto@fcfrp.usp.br

Editor's Note: *Caution.* Sigma-Aldrich scientists have not tested this procedure in-house and do not have any experience with it. Its publication in this magazine should not be construed as being endorsed or recommended by Sigma-Aldrich. The user should base his/her decision to use this technique solely on the claims made by the authors.

Microwave Synthesis: Chemistry at the Speed of Light by Dr. B. L. Hayes (Aldrich Cat. No. Z55,386-7). See pages 50 and 72 for more details.

Table 1. Sublimation Using a Microwave Oven^a

Entry	Compound Sublimed		Microwave Setting (%) ^b	Heating Time (s)	"Yield" of Sublimate (%)
	Formula or Name	Amount (g)			
1	I ₂	1.0	60	360	>99
2	AlCl ₃	0.5	60	180	64
3	Hg ₂ Cl ₂	0.5	60 and 80	360	— ^c
4	(C ₁₀ H ₁₀)Fe	0.5	30	60	92
5	(C ₁₀ D ₁₀)Fe	0.1	30	45	96
6	Acetyl ferrocene	0.5	60	60	15 ^d
7	Mn ₂ (CO) ₁₀	0.1	40	40	58
8	Mn(CO) ₅ Br	0.1	30	40	33 ^e
9	Mo(CO) ₆	0.3	30	60	81
10	(C ₅ H ₅)Mn(CO) ₃	0.1	30	30	72 ^e
11	(-)-Menthol	0.5	30	180	98
12	(±)-Camphor	0.5	20	60	84
13	Vanillin	0.5	10	30	— ^f
14	Piperonal	0.5	80	20	66
15	Biphenyl	0.5	40	150	78
16	Naphthalene	0.5	40	120	96
17	Anthracene	0.5	40	360	89
18	Salicylic acid	0.5	40	120	97
19	Benzophenone	0.5	80	90	71
20	Benzoic acid	0.5	60	60	89

^aSharp microwave oven, model Carousel III; manufactured by SANYO®: Manaus, Amazonas State, AM, Brazil. Of all the conditions tested, the best ones are shown in this table. ^bSetting as a percent of maximum power of 800 W. ^cHg₂Cl₂ appears to sublime only at higher temperatures. ^dWith extensive decomposition. ^eWith some decomposition. ^fThe material sublimed easily, but the vapors were lost without good condensation, even when low power was used.

Sigma-Aldrich Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA, S.A.
Tel: 54 11 4556 1472
Fax: 54 11 4552 1698

Australia

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: 612 9841 0555
Fax: 612 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: 43 1 605 81 10
Fax: 43 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA.
Free Tel: 0800-14747
Free Fax: 0800-14745
Tel: 03 899 13 01
Fax: 03 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Tel: 55 11 3732-3105
Fax: 55 11 3733-5151

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 800-565-1400
Free Fax: 800-265-3858
Tel: 905-829-9500
Fax: 905-829-9292

China

SIGMA-ALDRICH CHINA INC.
Tel: 86-21-6386 2766
Fax: 86-21-6386 3966

Czech Republic

SIGMA-ALDRICH s.r.o.
Tel: 246 003 200
Fax: 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: 43 56 59 10
Fax: 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND
Tel: 358-9-350-92 50
Fax: 358-9-350-92 555

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Tel appel gratuit: 0800 211 408
Fax appel gratuit: 0800 031 052

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800-51 55 000
Free Fax: 0800-649 00 00

Greece

SIGMA-ALDRICH (O.M.) LTD.
Tel: 30 210 9948010
Fax: 30 210 9943831

Hungary

SIGMA-ALDRICH Kft
Tel: 06-1-235-9054
Fax: 06-1-269-6470
Ingyenes zöld telefon: 06-80-355-355
Ingyenes zöld fax: 06-80-344-344

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: 91-80-852 4222 / 4150
Hyderabad:
91-40-5531 5548 / 2784 2378
Mumbai: 91-22-2579 7588 / 2570 2364
New Delhi:
91-11-2616 5477 / 2619 5360
Fax
Bangalore: 91-80-852 4214
Hyderabad: 91-40-5531 5466
Mumbai: 91-22-2579 7589
New Delhi: 91-11-2616 5611

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222

Israel

SIGMA-ALDRICH ISRAEL LTD.
Tel: 08-948-4100
Fax: 08-948-4200

Italy

SIGMA-ALDRICH S.r.l.
Tel: 02 33417310
Fax: 02 38010737
Numero Verde: 800-827018

Japan

SIGMA-ALDRICH JAPAN K.K.
Tokyo Tel: 03 5821 3111
Tokyo Fax: 03 5821 3170

Korea

SIGMA-ALDRICH KOREA
Tel: 031-329-9000
Fax: 031-329-9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD.
Tel: 603-56353321
Fax: 603-56354116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01-800-007-5300
Free Fax: 01-800-712-9920

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Tel Gratis: 0800-0229088
Fax Gratis: 0800-0229089
Tel: 078-6205411
Fax: 078-6205421

New Zealand

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 0800 936 666
Free Fax: 0800 937 777

Norway

SIGMA-ALDRICH NORWAY AS
Tel: 23 17 60 60
Fax: 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: (+61) 829 01 00
Fax: (+61) 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 20 21 80
Free Fax: 800 20 21 78

Russia

SIGMA-ALDRICH RUSSIA
TechCare Systems, Inc.
(SAF-LAB)
Tel: 095-975-1917/3321
Fax: 095-975-4792

Singapore

SIGMA-ALDRICH PTE. LTD.
Tel: 65-6271 1089
Fax: 65-6271 1571

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD.
Tel: 27 11 979 1188
Fax: 27 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA S.A.
Free Tel: 900101376
Free Fax: 900102028

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: 020-350510
Fax: 020-352522
Outside Sweden Tel: +46 8 7424200
Outside Sweden Fax: +46 8 7424243

Switzerland

FLUKA CHEMIE GmbH
Swiss Free Call: 0800 80 00 80
Tel: +41 81 755 2828
Fax: +41 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717181
Free Fax: 0800 378785
Tel: 01747 833000
Fax: 01747 833313

United States

SIGMA-ALDRICH CO.
P.O. Box 14508
St. Louis, Missouri 63178
Toll-Free: 800-325-3010
Call Collect: 314-771-5750
Toll-Free Fax: 800-325-5052
Tel: 314-771-5765
Fax: 314-771-5757

Internet:

sigma-aldrich.com

Solutions for Oligonucleotide Synthesis

Precisely Blended Capping, Deblocking, and Activator solutions for oligonucleotide synthesis.

Catalog No.	Product Name	Unit	Catalog No.	Product Name	Unit
55,534-7	Cap Mix A (Contains 90% tetrahydrofuran: 10% acetic anhydride)	1L 2L	55,535-5	Cap Mix B, with pyridine (Contains 80% tetrahydrofuran: 10% pyridine: 10% 1-methylimidazole)	1L 2L
55,531-2	Cap Mix A, with 2,6-lutidine (Contains 80% tetrahydrofuran: 10% acetic anhydride: 10% 2,6-lutidine)	1L 2L	56,193-2	Deblock (Contains 3% trichloroacetic acid in dichloromethane)	1L 2L
55,533-9	Cap Mix A, with pyridine (Contains 80% tetrahydrofuran: 10% acetic anhydride: 10% pyridine)	1L 2L	55,404-9	Activator (1-H-Tetrazole, 3 wt. % solution in acetonitrile)	1L 2L
55,532-0	Cap Mix B (Contains 84% tetrahydrofuran: 16% 1-methylimidazole)	1L 2L			

For additional pre-blended solvents or larger container options, email us at labchem@sial.com.

 SIGMA-ALDRICH

Fine Chemicals for Full-Scale Reliability

Ready to Scale Up?

Process development chemists at Sigma-Aldrich Fine Chemicals take pride in helping our customers succeed. Combining decades of process development experience with the newest process technologies, we can meet your development, scale-up, and manufacturing needs with our broad range of capabilities.

Organic capabilities include:

- Air-sensitive reactions
- Friedel-Crafts
- Grignards
- High-pressure reactions
- Hydrogenation
- Low-temperature reactions to $-110\text{ }^{\circ}\text{C}$
- Phosgenation

Biochemical capabilities include:

- Serum-free/animal-protein-free cell culture media
- Isolation and purification of proteins from transgenic and natural sources
- Bio-organic synthesis of lipids, poly(amino acids), carbohydrates, and buffers

When your project moves to large-scale manufacturing, Sigma-Aldrich's worldwide facilities can support many specialized techniques and reactions in scales ranging from liters to thousands of gallons.

To obtain a quote or place an order, please call 800-336-9719, visit us at www.sigma-aldrich.com/safc, or email us at safinechem@sial.com.

sigma-aldrich.com/safc

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
SIGMA ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI • USA


SIGMA-ALDRICH
FINE CHEMICALS

ALDRICH CHEMICAL COMPANY, INC.
P.O. BOX 14508
ST. LOUIS, MISSOURI 63178 USA



CHANGE SERVICE REQUESTED

F

ALDRICH CONGRATULATES THE 2004 ACS AWARD WINNERS

Aldrichimica **ACTA**

VOL. 36, NO. 3 • 2003

(Last issue in 2003)



**Organosilicon
Reagents: Synthesis and
Application to
Palladium-Catalyzed
Cross-Coupling
Reactions**

**Cross Metathesis of
Nitrogen-Containing
Systems**

sigma-aldrich.com



ALDRICH®

CHEMISTS HELPING CHEMISTS

New Products from Aldrich R&D


2-(4-Biphenyl)-5-phenyloxazole, 98%

21,698-4  1g \$13.80
5g 45.90

Widely used as an activator in C_6F_6 -based scintillators for the study of α -particles,¹ γ -rays,² and solar neutrinos from β -decay.³

(1) Britvich, G.I. et al. *Nucl. Instrum. Methods Phys. Res., Sect. A* **1993**, A326, 483; *Chem. Abstr.* **1993**, 118, 200826r. (2) Barabanov, I.R. et al. *Prib. Tekh. Eksp.* **1995**, 75; *Chem. Abstr.* **1995**, 123, 125307v. (3) Barabanov, I.R. et al. *ibid.* **1996**, 41; *Chem. Abstr.* **1996**, 125, 98054t.


α -Sexithiophene

59,468-7  1g \$80.00

This thiophene oligomer is known for its semiconducting, electrochemical, and photoelectric properties.

Bungs, M.; Tributsch, H. *J. Appl. Electrochem.* **2002**, 32, 91.


1,10-Phenanthroline-5-ylamine, 97%

63,150-7  5g \$195.00

Employed in the formation of luminophores and metal-ligand complexes for the detection of chemical and biochemical materials.^{1,2}

(1) Lecomte, J.-P. et al. *J. Chem. Soc., Faraday Trans.* **1993**, 89, 3261. (2) Meggers, E. et al. *Helv. Chim. Acta* **1997**, 80, 640.


tert-Butyl N-(2-oxiranylmethyl)carbamate

63,066-7  1g \$58.00


This functionalized epoxide can undergo hydrolytic kinetic resolution.¹ It has also been used as a building block for the construction of HIV proteinase inhibitors.^{2,3}

(1) Schaus, S.E. et al. *J. Am. Chem. Soc.* **2002**, 124, 1307. (2) Rocheblave, L. et al. *J. Med. Chem.* **2002**, 45, 3321. (3) Kitchin, J. et al. *J. Med. Chem.* **1994**, 37, 3707.


tert-Butyl phenyl carbonate, 98%

12,430-3  25g \$37.40
100g 106.40

Allyl phenyl carbonate, 97%

63,065-9  1g \$11.00
5g 36.40


Benzyl phenyl carbonate, 97%

63,064-0  5g \$29.90

These carbonates provide a practical and versatile method for selective Boc, Alloc, and Cbz protection of primary amines in simple symmetrical aliphatic diamines, and can selectively protect primary amines in the presence of secondary amines.

Pittelkow, M. et al. *Synthesis* **2002**, 2195.


2,4,6-Trichloropyridine, 97%

63,353-4  1g \$35.00
5g 120.00

Building block for a variety of biologically active compounds such as some cephalosporins,¹ anticancer agents,² and herbicides.³

(1) D'Andrea, S.V. et al. *Tetrahedron* **2000**, 56, 5687. (2) Kyoji, T. et al. Int. Patent Appl. WO 9534,559, Dec 21, 1995; *Chem. Abstr.* **1996**, 124, 289559a. (3) Hisashi, K. et al. Eur. Patent Appl. EP 693,490, Jan 24, 1996; *Chem. Abstr.* **1996**, 124, 289242s.


β -D-Galactose pentapivalate, 98%

63,247-3  1g \$18.00
5g 55.00

This fully protected carbohydrate is a convenient precursor for 2,3,4,6-tetra-O-pivaloyl- β -D-galactopyranosylamine, a useful chiral auxiliary for chiral α -amino nitrile preparation via the Strecker synthesis.

Kunz, H. et al. *Tetrahedron Lett.* **1988**, 29, 4397.


6-Chloropyran-2-one, 97%

63,299-6  1g \$24.50
5g 84.50

Has been used as a starting material for alkynylpyranones,¹ and for a wide spectrum of biologically active substrates including pretetramides.²

(1) Biagetti, M. et al. *Tetrahedron Lett.* **2003**, 44, 607. (2) Gilbreath, G. S. et al. *J. Am. Chem. Soc.* **1988**, 110, 6172.


1-Benzenesulfinylpiperidine, 97%

63,023-3  1g \$13.80
5g 45.90

Novel reagent for the synthesis of glycosides from thioglycoside donors in high yields and excellent stereoselectivities.

Crich, D.; Li, H. *J. Org. Chem.* **2002**, 67, 4640.


4-(Trimethylsilylethynyl)morpholine, 97%

63,277-5  1g \$10.00
10g 29.90

Employed in the Lewis acid mediated ring opening of terminal epoxides leading to the corresponding γ -butanolides.

Movassaghi, M.; Jacobsen, E. N. *J. Am. Chem. Soc.* **2002**, 124, 2456.

Trichloro(4-chlorophenyl)silane, 97%

63,045-4  1g \$18.50
5g 61.65

Employed as a starting material in the synthesis of tripod-shaped oligophenylenes designed for thin-film applications.

Deng, X.; Cai, C. *Tetrahedron Lett.* **2003**, 44, 815.

Please see pages 90–91 for additional new products.

For competitive quotes on larger quantities contact www.sigma-aldrich.com/safc.

SIGMA-ALDRICH
FINE CHEMICALS

Aldrichimica ACTA

VOL. 36, NO. 3 • 2003

Aldrich Chemical Co., Inc.
1001 West Saint Paul Ave.
Milwaukee, WI 53233 USA

To Place Orders

Telephone 800-558-9160 (USA)
800-231-8327
FAX 800-962-9591 (USA)
or 414-273-4979
Mail P.O. Box 2060
Milwaukee, WI 53201 USA

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201 USA

Customer & Technical Services

Customer Inquiries 800-558-9160
Technical Service 800-231-8327
Sigma-Aldrich Fine Chemicals 800-336-9719
Custom Synthesis 800-336-9719
Flavors & Fragrances 800-227-4563
International 414-273-3850
24-Hour Emergency 414-273-3850
Website <http://www.sigma-aldrich.com>
Email aldrich@sial.com

To request your **FREE** subscription to the *Aldrichimica Acta*, please

call: 800-558-9160 (USA),
write: Attn: Mailroom
Aldrich Chemical Co., Inc.
P.O. Box 355
Milwaukee, WI 53201-9358
or email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office.

The *Aldrichimica Acta* is also available on the Internet at <http://www.sigma-aldrich.com>.

All prices listed in this publication are subject to change without notice.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

For worldwide contact information, please see the inside back cover.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich family. ©2003 by Sigma-Aldrich Co


About Our Cover

Autumn—On the Hudson River (oil on canvas, 151.8 x 274.9 cm) was signed and dated by the American painter Jasper Francis Cropsey in 1860. This enormous painting represents a panoramic view of the Hudson River valley about 60 miles north of New York City near West Point and Storm King Mountain.

At first glance, this landscape, painted in brilliant autumnal colors under a magnificent sunlit sky, may appear to show nature in a completely wild state. Mankind is not entirely absent, however. In the foreground and left of center, three hunters and their dogs have stopped to rest; a log cabin sits among the trees in the middle distance on the right; a town can be seen along the bank of the river; and a number of boats, including a steamer, are on the river itself. The setting is neither completely untouched by man nor overly domesticated. Man may at first seem dwarfed by nature, but is shown here to take his place harmoniously in the natural world.

More than simply a visual record of a certain time and place, however, this painting embodies certain ideas that were current in nineteenth-century America. The natural world was thought to be the most profound manifestation of the Divine order. Moreover, the magnificence of the American landscape came to signify the expansionist ideal and the opportunities and potential greatness of this new country. The critics praised the picture extravagantly when it was exhibited in London, where it was painted from memory and from sketches brought from America during the second of two study trips the artist made to Europe. Viewers, however, questioned the brilliant colors of the foliage represented in the painting, which were more intense than anything they had ever seen. Cropsey, however, had thought to bring from America samples of brightly colored autumn leaves pasted on cardboard to demonstrate that his painting was not an exaggeration, but was quite true to nature, at least in America.

This painting is a gift of the Avalon Foundation to the National Gallery of Art, Washington, DC.




Photograph © Board of Trustees, National Gallery of Art, Washington

“Please
Bother
Us.”

Joe Porwoll

Joe Porwoll, President



Professor Donal F. O'Shea of the Department of Chemistry at University College Dublin kindly suggested that we offer 2,4,6-trivinylcyclotriboroxane-pyridine complex. This stable vinylboronic acid equivalent undergoes facile Suzuki cross-coupling with aryl halides to provide valuable functionalized styrene derivatives.

Kenns, F.; O'Shea, D. F. *J. Org. Chem.* 2002, 67, 4968

63,799-8
Vinyl boronic anhydride-pyridine complex
1g \$21.00; 5g \$73.50
56,814-7
Vinylboronic acid dibutyl ester, 97%
1g \$25.70; 5g \$85.40
63,334-8
Vinylboronic acid pinacol ester, 95%
1g \$20.40; 10g \$112.50

Naturally, we made not only this reagent but also two other stable vinylboronic acid esters—vinylboronic acid dibutyl ester and vinylboronic acid pinacol ester—which are useful for the Suzuki coupling, Heck coupling, and Grubb's olefin cross-metathesis reactions. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page or on the inside back cover.

Table of Contents

Organosilicon Reagents: Synthesis and Application to Palladium-Catalyzed Cross-Coupling Reactions	75
—S. E. Denmark* and M. H. Ober	
Cross Metathesis of Nitrogen-Containing Systems	93
—A. J. Vernall and A. D. Abell	



Lab Notes

A Funnel Assembly for the Safe Disposal of Nitrogen Discharged from an NMR Magnet during Cryogen Refilling

In NMR laboratories, where large volumes of cryogenics are routinely dispensed, a condition of displacement or deprivation of atmospheric oxygen may occur if the released gases are not efficiently removed. There is at least one reported case of displacement of oxygen by nitrogen during the installation of a magnetic resonance imaging system causing the death of a worker by asphyxiation.¹ In addition, similar incidents of asphyxiation resulting from unsafe handling of liquid nitrogen have been reported in other laboratory settings.^{2,3} In setting up a Bruker 400 MHz NMR instrument in our laboratory, we were constrained to install the magnet as well as the ancillary LC-NMR components in a 6 1/2'-deep well and the computers at the main level of the room. The well and the main level of the room were equipped with ZoneGuard sensors (Biosystems Inc., Middletown, CT) for a constant monitoring of the well and room oxygen levels. The monitoring systems were set to activate an audible alarm if the oxygen level dropped to 19.5% from the normal reading of 20.9%. During refilling of the magnet with liquid nitrogen, the oxygen level in the well dropped below the danger level of 18.0% and the alarm rang continuously.⁴ To minimize the risk of oxygen displacement, we explored the possibility of diverting the flow of cold nitrogen gas exiting the magnet into an area outside the NMR room. To accomplish this goal, we have fabricated a funnel assembly (Figure 1) for efficiently capturing and disposing of nitrogen gas and for maintaining safe oxygen levels in the room.

The assembly consists of a funnel mounted on a tripod, a clamp supported by the funnel, and a flexible hose attached to the stem of the funnel—as shown in the figure.⁵ The mouth of the funnel is sealed with a lid that has a circular opening (3" diam) at its center, and through which the nitrogen stream enters the funnel's chamber. The internal diameter of the funnel at the lid is 9" and that of the stem 2 1/4" and, at 2" into its depth, the funnel's body gradually begins to narrow along the 8"-long curvature. The length of the funnel, including the two-inch-long stem, is 12" and its capacity is adequate to capture the nitrogen gas emitted from the NMR magnet during refilling. The clamp (4" long x 1/4" thick) is supported by a stand (7" high x 1/4" thick) welded to the rim of the funnel. A copper tube (7" long x 3/8" i.d.) is clamped in an upright position to deliver the nitrogen gas at the center of the funnel. The tip of the copper tube is positioned approximately 3" above the opening in the funnel lid, and its top end is attached to a latex tube (5' long x 3/8" i.d.) that receives nitrogen gas from the magnet. The supporting tripod consists of a partial O-ring (7" i.d.) mounted on top of 3 legs (16" high), which are fastened to a heavy circular base (14" diam). The O-ring has a 3/4"-wide slit through which the funnel along with its attached stand can be readily dismounted from the tripod whenever needed. The material used for the fabrication of the funnel is nonmagnetic stainless steel, whereas that of the tripod is aluminum. The stem of the funnel is connected to a cloth hose (14' long x 3" i.d.) for carrying nitrogen gas into an exhaust line. The figure shows a segment of the latex tubing carrying nitrogen into the funnel and of the cloth hose discharging the gas.

Before liquid-nitrogen refilling, the funnel assembly is placed near the magnet and the funnel's latex tubing is connected to the left-hand-side turret. During refilling of the magnet, the funnel collects the gas discharge with no significant diffusion into the air. As it is denser than air, cold nitrogen gas settles well in the funnel and flows through the hose without applying suction. It also appears that the force of nitrogen streaming through the circular opening prevents back diffusion and exerts pressure sufficient to cause expulsion of the gas from the funnel through the hose. The gas is discharged




Figure 1. Schematic Diagram of the Funnel Assembly.

near an exhaust line located in the well. When the magnet is full, liquid nitrogen spraying into the funnel is clearly visible from several feet away. There is practically no drop in the oxygen level in the well or in the rest of the room during the refilling procedure, suggesting efficient trapping and disposal of nitrogen by this device. Occasionally, water condenses in the first two-foot segment (from the stem of the funnel) of the exterior of the cloth hose carrying the cold gas; placing an absorbent pad underneath the hose takes care of this problem.

The funnel assembly helps in maintaining the room oxygen level during refilling of the NMR magnet with liquid nitrogen, especially if the room is not spacious or not adequately ventilated. In addition, it minimizes the chances of spillage of the cryogen and thus accidental freezing of the O-rings in the base plate and top flange of the magnet cryostat. (If freezing occurs, the sealing O-rings will be hardened leading to loss of vacuum between the casings of the magnet cryostat.) Certain modifications of the funnel assembly may be needed depending on the room's configuration or on whether added efficiency in the removal of nitrogen gas is warranted. For instance, if an exhaust line is not accessible, the delivery end of the hose may be placed outside the room for safe disposal of the gas. If improved gas flow is needed, a suction fan may be attached at the delivery end of the hose.

References and Notes: (1) Gill, J.R.; Ely, S.F.; Hua, Z. Environmental gas displacement: three accidental deaths in the workplace. *Am. J. Forensic Med. Pathol.* **2002**, *23*, 26; abstract available at the National Library of Medicine website at <http://www.ncbi.nlm.nih.gov/entrez/query.fcgi> (accessed March 2003). (2) Kernbach-Wighton, G.; Kijewski, H.; Schwanke, P.; Saur, P.; Sprung, R. Clinical and morphological aspects of death due to liquid nitrogen. *Int. J. Legal Med.* **1998**, *111*, 191. (3) Tabata, N.; Funayama, M.; Ikeda, T.; Azumi, J.-i.; Morita, M. On an accident by liquid nitrogen—histological changes of skin in cold. *Forensic Sci. Int.* **1995**, *76*, 61. (4) Aside from recommending the installation of oxygen sensors, the manufacturer does not provide any device or recommend any procedure for handling the nitrogen discharged from the magnet during cryogen refilling. (5) The funnel, stand, and clamp were fabricated by Atlantic Sheet Metal Manufacturing (Essex, MD). The tripod was fabricated at the Division of Engineering Services, National Institutes of Health (Bethesda, MD). The entire assembly can be readily fabricated in any workshop using nonmagnetic metal sheets and rods (e.g., aluminum, copper, or nonmagnetic stainless steel) and cloth hose.

H. Umeshia Shetty (Ph.D.),* **Jinsoo Hong**, and **Victor W. Pike**
PET Radiopharmaceutical Sciences Section
Molecular Imaging Branch
National Institute of Mental Health
National Institutes of Health
10 Center Drive, Room B3C351, MSC 1003
Bethesda, MD 20892-1003, USA
Email: shettyu@intra.nimh.nih.gov

Editor's Note: *Caution.* The potential user of this note should carefully evaluate its suitability for a particular application, as well as familiarize himself/herself with any potential hazards associated with the construction and use of this assembly. The procedure described in this note has not been tested by Sigma-Aldrich scientists. Its publication in this magazine should not be construed as being endorsed or recommended by Sigma-Aldrich.

Sigma-Aldrich Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA, S.A.
Tel: 54 11 4556 1472
Fax: 54 11 4552 1698

Australia

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 1800 800 097
Free Fax: 1800 800 196
Tel: 612 9841 0555
Fax: 612 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: 43 1 605 81 10
Fax: 43 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA.
Free Tel: 0800-14747
Free Fax: 0800-14745
Tel: 03 899 13 01
Fax: 03 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Tel: 55 11 3732-3100
Fax: 55 11 3733-5151

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 800-565-1400
Free Fax: 800-265-3858
Tel: 905-829-9500
Fax: 905-829-9292

China

SIGMA-ALDRICH CHINA INC.
Tel: 86-21-6386 2766
Fax: 86-21-6386 3966

Czech Republic

SIGMA-ALDRICH S.r.o.
Tel: 246 003 200
Fax: 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: 43 56 59 10
Fax: 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND
Tel: 358-9-350-92 50
Fax: 358-9-350-92 555

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Tel appel gratuit: 0800 211 408
Fax appel gratuit: 0800 031 052

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800-51 55 000
Free Fax: 0800-649 00 00

Greece

SIGMA-ALDRICH (O.M.) LTD
Tel: 30 210 9948010
Fax: 30 210 9943831

Hungary

SIGMA-ALDRICH Kft
Tel: 06-1-235-9054
Fax: 06-1-269-6470
Ingyenes zold telefon: 06-80-355-355
Ingyenes zold fax: 06-80-344-344

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: 91-80-5112-7272
Hyderabad: 91-40-5531 5548 / 2784 2378
Mumbai: 91-22-2579 7588 / 2570 2364
New Delhi: 91-11-2616 5477 / 2619 5360
Fax
Bangalore: 91-80-5112-7473
Hyderabad: 91-40-5531 5466
Mumbai: 91-22-2579 7589
New Delhi: 91-11-2616 5611

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222

Israel

SIGMA-ALDRICH ISRAEL LTD.
Tel: 08-948-4100
Fax: 08-948-4200

Italy

SIGMA-ALDRICH S.r.l.
Telefono: 02 33417310
Fax: 02 38010737
Numero Verde: 800-827018

Japan

SIGMA-ALDRICH JAPAN K.K.
Tokyo Tel: 03-5796-7300
Tokyo Fax: 03-5796-7315

Korea

SIGMA-ALDRICH KOREA
Tel: 031-329-9000
Fax: 031-329-9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: 603-56353321
Fax: 603-56354116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01-800-007-5300
Free Fax: 01-800-712-9920

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Tel Gratis: 0800-0229088
Fax Gratis: 0800-0229089
Tel: 078-6205411
Fax: 078-6205421

New Zealand

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 0800 936 666
Free Fax: 0800 937 777

Norway

SIGMA-ALDRICH NORWAY AS
Tel: 23 17 60 60
Fax: 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: (+61) 829 01 00
Fax: (+61) 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 20 21 80
Free Fax: 800 20 21 78

Russia

SIGMA-ALDRICH RUSSIA
TechCare Systems, Inc.
(SAF-LAB)
Tel: 095-975-1917/3321
Fax: 095-975-4792

Singapore

SIGMA-ALDRICH PTE. LTD.
Tel: 65-6271 1089
Fax: 65-6271 1571

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD.
Tel: 27 11 979 1188
Fax: 27 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA S.A.
Free Tel: 900101376
Free Fax: 900102028

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: 020-350510
Fax: 020-352522
Outside Sweden Tel: +46 8 7424200
Outside Sweden Fax: +46 8 7424243

Switzerland

FLUKA CHEMIE GmbH
Swiss Free Call: 0800 80 00 80
Tel: +41 81 755 2828
Fax: +41 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717181
Free Fax: 0800 378785
Tel: 01747 833000
Fax: 01747 833313

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-free: 800-325-3010
Call Collect: 314-771-5750
Toll-Free Fax: 800-325-5052
Tel: 314-771-5765
Fax: 314-771-5757
Internet:
sigma-aldrich.com


Solutions for Oligonucleotide Synthesis

Precisely Blended Capping, Deblocking, and Activator solutions for oligonucleotide synthesis.

Catalog No.	Product Name	Unit	Price
55,534-7	Cap Mix A (Contains 90% tetrahydrofuran: 10% acetic anhydride)	1L	\$62.60
		2L	96.30
55,531-2	Cap Mix A, with 2,6-lutidine (Contains 80% tetrahydrofuran: 10% acetic anhydride: 10% 2,6-lutidine)	1L	98.60
		2L	157.60
55,533-9	Cap Mix A, with pyridine (Contains 80% tetrahydrofuran: 10% acetic anhydride: 10% pyridine)	1L	78.90
		2L	124.10
55,532-0	Cap Mix B (Contains 84% tetrahydrofuran: 16% 1-methylimidazole)	1L	108.40
		2L	174.30

Catalog No.	Product Name	Unit	Price
55,535-5	Cap Mix B, with pyridine (Contains 80% tetrahydrofuran: 10% pyridine: 10% 1-methylimidazole)	1L	\$105.00
		2L	168.40
56,193-2	Deblock (Contains 3% trichloroacetic acid in dichloromethane)	1L	82.60
55,404-9	Activator (1H-Tetrazole, 3 wt. % solution in acetonitrile)	1L	221.30
		2L	344.20

For additional pre-blended solvents or larger container options, email us at labchem@sial.com.

 SIGMA-ALDRICH

DISCOVERY

Drug Discovery

PRECLINICAL

Toxicology

CLINICAL

Phase I

Phase II

Phase III

Phase IV

Drug Approval

COMMERCIALIZATION

Sigma-Aldrich Fine Chemicals provides products and contract manufacturing services from Discovery through Commercialization to the Pharmaceutical, Biopharmaceutical, and High Technology industries.

We manufacture over 40,000 products.

What can we manufacture for you?



Take a Virtual Tour of our MANUFACTURING facilities at
<http://www.sigma-aldrich.com/tours>

To discuss how our expertise can benefit your next scale-up project or to obtain a quote, call us at 800-336-9719 (USA) or visit us at www.sigma-aldrich.com/safc

France: tel. 04 74 82 28 82 - fax 04 74 82 28 90 • Germany: tel. 089 6513 1900 - fax 089 6513 1919
• Italy: tel. 02 33 41 7350 - fax 02 33 41 7201 • Japan: tel. 03 5821 3171 - fax 03 5821 3170 • Switzerland: tel. 081 755 2732 - fax 081 755 2770
U.K.: tel. 01202 712222 - fax 01202 740171 • USA: tel. 800 336 9719 or 314 534 4900 - fax 314 652 0000 • Email: safinechem@sial.com

sigma-aldrich.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
SIGMA-ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI • USA


SIGMA-ALDRICH
FINE CHEMICALS

ALDRICH CHEMICAL COMPANY, INC.
P.O. BOX 14508
ST. LOUIS, MISSOURI 63178 USA



CHANGE SERVICE REQUESTED

PRESORTED STANDARD
U.S. POSTAGE
PAID
SIGMA-ALDRICH
CORPORATION

GEP

PROAZAPHOSPHATRANES: SUPERIOR BASES AND CATALYSTS

Aldrichimica ACTA

VOL. 37, NO. 1 • 2004



**Recent Applications of
Proazaphosphatranes
in Organic Synthesis**


**Chemistry Without
Reagents: Synthetic
Applications of Flash
Vacuum Pyrolysis**

sigma-aldrich.com


 **ALDRICH®**

New Products from Aldrich R&D

Diethyl 2-butenylphosphonate, 95%, predominantly trans

59,759-7  1g \$14.65
5g 48.75

Diethyl 3-butenylphosphonate

64,052-2  1g \$23.50
10g 130.00

Diethyl cyclopropylmethylphosphonate, 97%

63,582-0  1g \$47.00
5g 163.00

Phosphonates have been utilized extensively in the Horner–Wadsworth–Emmons reaction for the construction of carbon–carbon double bonds. These phosphonates have been employed in the preparation of 3-acetoxy-1-alkenyl phosphonates.¹ They were also exploited in an efficient and regioselective synthesis of 4-oxo-2-alkenylphosphonates.²

(1) Principato, B. et al. *Tetrahedron* **1996**, 52, 2087. (2) Lee, B.S. et al. *J. Org. Chem.* **2000**, 65, 4175

Methyl 2-((succinimidoxy)carbonyl)benzoate, 97%

63,392-5  5g \$32.50

This reagent efficiently undergoes N-phthaloylation with α -amino acids, α -amino alcohols, α -amino carboxamides, α -amino esters, and dipeptides in high yields without racemization.

Casimir, J. et al. *J. Org. Chem.* **2002**, 67, 3764


Benzyl thioacetimidate hydrochloride, 95%

63,711-4  1g \$14.50
5g 51.00

This functionalized thioimide is a useful reagent for the conversion of anilines to N-substituted acetamidines.^{1,2}

(1) Doise, M. et al. *Tetrahedron Lett.* **1990**, 31, 1155. (2) Collins, J. L. et al. *J. Med. Chem.* **1998**, 41, 2858


3,5-Dimethoxyphenylmagnesium chloride

63,762-9  50mL \$69.50

This Grignard reagent was used recently to make materials for studies in crystal engineering.

Tanaka, T. et al. *J. Am. Chem. Soc.* **2002**, 124, 12453


2-Bromo-3-thiophenecarboxylic acid, 97%

63,812-9  1g \$39.00
5g 135.00

This bifunctional building block was recently utilized in the preparation of oligothiophenes for materials science applications.^{1,2}

(1) Zhang, O. T.; Tour, J. M. *J. Am. Chem. Soc.* **1997**, 119, 9624. (2) Pomerantz, M. et al. *J. Org. Chem.* **2002**, 67, 6931


4-(tert-Butyldimethylsilyloxy)cyclohexanone

63,815-3  1g \$11.50
5g 39.00

As a protected 4-hydroxycyclohexanone, this compound has been exploited in the synthesis of butenolides,¹ (\pm)-cocculolide,² and the enyne A-ring synthon of 1 α -hydroxy vitamin D₃.³

(1) Majewski, M. et al. *Tetrahedron: Asymmetry* **1995**, 6, 1837. (2) Kawasaki, T. et al. *Tetrahedron Lett.* **2001**, 42, 8003. (3) Parker, K. A.; Dermatakis, A. *J. Org. Chem.* **1997**, 62, 6692


N-Boc-2-naphthalenesulfonamide, 97%

63,921-4  1g \$21.00
5g 72.00

This novel reagent allows the stepwise synthesis of secondary aliphatic amines, such as selectively protected derivatives of spermidine. Cleavage of the sulfonylcarbamate group can be achieved with catalytic amounts of Mg(ClO₄)₂ without affecting other Boc-protected amines.


Grehn, L.; Ragnarsson, U. *J. Org. Chem.* **2002**, 67, 6557

Diazald[®], 20 wt. % mixture with sodium sulfate


64,053-0  25g \$22.00
20 wt. % mixture with sodium sulfate

This well-known reagent has been used for the safe production of diazomethane. Aldrich now proudly presents this mixture to address shipping and handling concerns.

6-Quinolyl trifluoromethanesulfonate, 97%

63,344-5  1g \$12.00
5g 39.50

9H-Carbazol-2-yl trifluoromethanesulfonate, 97%

63,920-6  1g \$17.00
5g 58.50

These heterocyclic triflates are convenient building blocks for the synthesis of the corresponding halides,¹ thiols,² and alkenes.¹ They have also been exploited in cross-coupling reactions with Grignard reagents.

(1) Ritter, K. *Synthesis* **1993**, 735. (2) Arnould, J. C. et al. *Tetrahedron Lett.* **1996**, 37, 4523
(3) Furstner, A. *J. Am. Chem. Soc.* **2002**, 124, 13856

Please see pages 32–33 for additional new products.

For competitive quotes on larger quantities contact www.sigma-aldrich.com/safc.

SIGMA-ALDRICH
FINE CHEMICALS

Diazald is a registered trademark of Sigma-Aldrich Biotechnology, L.P.

Aldrichimica ACTA

VOL. 37, NO. 1 • 2004

Aldrich Chemical Co., Inc.
1001 West Saint Paul Ave.
Milwaukee, WI 53233, USA

To Place Orders

Telephone 800-558-9160 (USA)
FAX 800-962-9591 (USA)
or 414-273-4979
Mail P.O. Box 2060
Milwaukee, WI 53201 USA

Customer & Technical Services

Customer Inquiries 800-558-9160
Technical Service 800-231-8327
Sigma-Aldrich Fine Chemicals 800-336-9719
Custom Synthesis 800-336-9719
Flavors & Fragrances 800-227-4563
International 414-273-3850
24-Hour Emergency 414-273-3850
Website <http://www.sigma-aldrich.com>
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201, USA

Subscriptions

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-558-9160 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
P.O. Box 355
Milwaukee, WI 53201-9358
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the inside back cover.

The *Aldrichimica Acta* is also available on the Internet at <http://www.sigma-aldrich.com>.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.



Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich family. ©2004 by Sigma-Aldrich Co.

"PLEASE BOTHER US."

Joe Porwoll, President

Professor Masaru Kobayashi of Hokkaido University, Japan, kindly suggested we offer 1,2,3,5-tetra-O-acetyl- α -D-arabinofuranose and 1,2,3,5-tetra-O-acetyl- α -L-arabinofuranose. These fully protected carbohydrates are convenient reagents for determining the absolute configuration of secondary alcohols by ^{13}C NMR spectroscopy.

Kobayashi, M. *Tetrahedron* 2002, 58, 9365



63,480-8 1,2,3,5-Tetra-O-acetyl- α -D-arabinofuranose 1g \$34.50; 5g \$125.00
63,744-0 1,2,3,5-Tetra-O-acetyl- α -L-arabinofuranose 1g \$34.50; 5g \$125.00

Naturally, we made these valuable reagents. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page or on the inside back cover.

TABLE OF CONTENTS

An Inexpensive, Manually Operated, Solid-Phase, Parallel Synthesizer 2
Balvinder S. Vig and Jane V. Aldrich, University of Maryland*

Recent Applications of Proazaphosphatranes in Organic Synthesis 3
John G. Verkade and Philip B. Kisanga, Iowa State University and Albany Molecular Research, Inc.*

Chemistry Without Reagents: Synthetic Applications of Flash Vacuum Pyrolysis 19
Hamish McNab, The University of Edinburgh

ABOUT OUR COVER

Flower Beds in Holland (oil on canvas, 48.9 x 66.0 cm) was painted by Vincent van Gogh, probably in 1883. Born in 1853, he left school at the age of sixteen to work for the Goupil art firm in The Hague. This was followed by a period of wandering and spiritual anguish that took him to London (1873), Paris (1875), and Belgium (1878), after which he returned to Holland determined to become an artist.

Van Gogh studied at the Antwerp Academy, but was impatient with formal training and, in 1886, left for Paris. His earliest works portray sympathetically the lives of peasants and workmen, but after he met Pissarro, Degas, Gauguin, Seurat, and Toulouse-Lautrec, he developed an obsessive interest in the symbolic and expressive possibilities of colors. In a frenzy of creation during the last years of his life, he produced over 800 paintings and 850 drawings.

Although its date makes it one of Van Gogh's earliest known works, one can already see in *Flower Beds in Holland* the interest in surface texture and expressive color so prominent in his later works. The bright stripes of color of the flower beds in the foreground contrast strongly with the relatively flat outlines of the tree branches against the gray sky. The low horizon follows the tradition of Dutch landscape painting, which reflects the natural geography of the country. Van Gogh's handling of paint and use of color, however, show the influence of his contacts with contemporary French artists, and point the way to the ever more impassioned and brilliantly colored works of his later career.

This painting is in the Collection of Mr. and Mrs. Paul Mellon at the National Gallery of Art, Washington, DC.



Photograph © Board of Trustees, National Gallery of Art, Washington

Lab Notes


An Inexpensive, Manually Operated, Solid-Phase, Parallel Synthesizer

Solid-phase synthesis is routinely used for the preparation of peptides and small molecules.¹⁻⁶ Because solid-phase synthesis involves sequential mixing and draining steps, it is frequently performed in automated solid-phase synthesizers. Although solid-phase synthesis offers many advantages over its solution-phase counterpart, the cost of instrumentation is a limiting factor that prevents many laboratories from venturing into this methodology.⁷ Reaction vessels for manual synthesis are commercially available, but each usually requires its own setup for mixing and draining. This makes it more complicated to run multiple reactions at a time. To address this shortcoming, we have designed an inexpensive, manually operated, solid-phase, multiple synthesizer ("CHOIR")⁸ that is easy to use and maintain, and that should be affordable (~\$250 for glassware) to a wide variety of researchers and educators, especially in cost-conscious academic laboratories.

CHOIR, in its current design, permits up to six syntheses to be performed simultaneously and independently. It consists of a modified vacuum manifold (Chemglass custom Cat. No. UM-2008-301D) with 8 outlets (see figure). The left end of the manifold serves as an inlet for the inert gas, which is bubbled through the reaction vessels in order to mix the reactants and resins. A regulator controls the gas flow and maintains enough gas pressure through the reaction vessels to prevent premature drainage of the solvents. The right end has a three-way stopcock that is opened only when the solvents are being drained. Solvents are drained into a solvent trap with the aid of a vacuum pump or aspirator. The remaining six outlets have LUER® connectors, which are used to attach the reaction vessels via two-way stopcocks. Solid-phase extraction tubes with polyethylene frits, available in a variety of sizes (e.g., Supelco Cat. No. 57176), are used as reaction vessels. These tubes are inexpensive and can be discarded after each synthesis. The inert gas, the solvent trap, and the vacuum pump are connected to CHOIR by means of chemically resistant tubing (e.g., Aldrich Cat. No. Z27,986-2). CHOIR is very easy to operate, since each reaction vessel can be controlled separately by closing and opening the individual stopcock attached to it. Once the parallel syntheses are completed, the reaction vessels are removed from CHOIR, and the peptides individually cleaved from the resins in the same reaction vessels.

We have used CHOIR for the synthesis, partial or complete, of a large number (100–150) of linear and cyclic peptides of various sizes (5–11 residues) in high yields and purities, or for the optimization of the reaction conditions leading to these peptides. Examples of peptides synthesized include cyclic and linear analogs of the opioid peptide dynorphin A:^{5,6,9-11}

- *cyclo*[D-Asp²,Xxx³,Dap⁵]Dyn A-(1–11)NH₂ and linear [D-Asp²,Xxx³,Dap⁵]Dyn A-(1–11)NH₂, where Xxx = Gly, Ala, D-Ala, Trp, D-Trp, or Pro.
- *cyclo*[D-Asp²,Xxx⁴,Dap⁵]Dyn A-(1–11)NH₂ and linear [D-Asp²,Xxx⁴,Dap⁵]Dyn A-(1–11)NH₂, where Xxx = Phe, D-Phe, HomoPhe, or D-HomoPhe.
- *cyclo*^x[Trp³,Trp⁴,Glu⁵]Dyn A-(1–11)NH₂ and its cyclic and linear analogs.
- *cyclo*^x[COCH²,Tyr¹,Lys⁵]Dyn A-(1–11)NH₂ and analogs.
- *cyclo*^x[COCH²,Tyr¹,Lys⁵]Dyn A-(1–11)NH₂ and analogs.



References and Notes: (1) For a review, see *Solid-Phase Synthesis, A Practical Guide*; Kates, S. A., Albericio, F., Eds.; Marcel Dekker, Inc.: New York, 2000. (2) Gorman, J. J. *Anal. Biochem.* **1984**, *136*, 397. (3) Stewart, J. M.; Young, J. D. *Solid Phase Peptide Synthesis*, 2nd ed.; Pierce Chemical Co.: Rockford, IL, 1984. (4) Knapp, D. R.; Oatis, J. E., Jr.; Papac, D. I. *Int. J. Pept. Protein Res.* **1993**, *42*, 259. This report describes a similar albeit less practical synthesizer, and gives a brief history of systems devised to carry out parallel syntheses of peptides. (5) Bennett, M. A.; Murray, T. F.; Aldrich, J. V. *J. Med. Chem.* **2002**, *45*, 5617. (6) Vig, B. S.; Murray, T. F.; Aldrich, J. V. *J. Med. Chem.* **2003**, *46*, 1279. (7) For example, the Bohdan MiniBlock® system, which is convenient to use, is roughly 14 times more expensive than the system described here (\$3,500 vs \$250). (8) CHOIR = Cheap Homemade Organic Inline Reactor. (9) Vig, B. S.; Aldrich, J. V. Strategies for the Synthesis of Novel Head-to-Side Chain Cyclic Peptides: Application to Dynorphin A Analogs. In *Peptides: The Wave of the Future*; Lebl, M., Houghten, R. A., Eds.; American Peptide Society: San Diego, CA, 2001; pp 144–145. (10) Vig, B. S.; Murray, T. F.; Aldrich, J. V. *J. Med. Chem.* **2004**, *47*, 446. (11) Vig, B. S.; Murray, T. F.; Aldrich, J. V. *Biopolymers* **2004**, in press.

LUER-LOK is a registered trademark of Becton, Dickinson and Company. MiniBlock is a registered trademark of METTLER-TOLEDO Bohdan, Inc.

Balvinder S. Vig*† and Jane V. Aldrich†

Department of Pharmaceutical Sciences
School of Pharmacy
University of Maryland
Baltimore, MD 21201, USA

Current Addresses: (‡) Dr. Balvinder S. Vig, Pharmaceutical Research Institute, Bristol-Myers Squibb Co., One Squibb Drive, P. O. Box 191, New Brunswick, NJ 08903-0191, USA; Email: balvinder.vig@bms.com. (†) Dr. Jane V. Aldrich, Department of Medicinal Chemistry, School of Pharmacy, University of Kansas, 4050 Malott Hall, 1251 Wescoe Hall Drive, Lawrence, KS 66045-7582, USA; Email: jaldrich@ku.edu.

Sigma-Aldrich Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA, S.A.
Tel: 54 11 4556 1472
Fax: 54 11 4552 1698

Australia

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: 612 9841 0555
Fax: 612 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: 43 1 605 81 10
Fax: 43 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA.
Free Tel: 0800-14747
Free Fax: 0800-14745
Tel: 03 899 13 01
Fax: 03 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Tel: 55 11 3732-3100
Fax: 55 11 3733-5151

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 800-565-1400
Free Fax: 800-265-3858
Tel: 905-829-9500
Fax: 905-829-9292

China

SIGMA-ALDRICH CHINA INC.
Tel: 86-21-6386 2766
Fax: 86-21-6386 3966

Czech Republic

SIGMA-ALDRICH s.r.o.
Tel: 246 003 200
Fax: 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: 43 56 59 10
Fax: 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND
Tel: 358-9-350-92 50
Fax: 358-9-350-92 555

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Tel appel gratuit: 0800 211 408
Fax appel gratuit: 0800 031 052

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800-51 55 000
Free Fax: 0800-649 00 00

Greece

SIGMA-ALDRICH (O.M.) LTD
Tel: 30 210 9948010
Fax: 30 210 9943831

Hungary

SIGMA-ALDRICH Kft
Tel: 06-1-235-9054
Fax: 06-1-269-6470
Ingyenes zöld telefon: 06-80-355-355
Ingyenes zöld fax: 06-80-344-344

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: 91-80-5112-7272
Hyderabad:
91-40-5531 5548 / 2784 2378
Mumbai: 91-22-2579 7588 / 2570 2364
New Delhi:
91-11-2616 5477 / 2619 5360
Fax
Bangalore: 91-80-5112-7473
Hyderabad: 91-40-5531 5466
Mumbai: 91-22-2579 7589
New Delhi: 91-11-2616 5611

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222

Israel

SIGMA-ALDRICH ISRAEL LTD.
Tel: 08-948-4100
Fax: 08-948-4200

Italy

SIGMA-ALDRICH S.r.l.
Telefono: 02 33417310
Fax: 02 38010737
Numero Verde: 800-827018

Japan

SIGMA-ALDRICH JAPAN K.K.
Tokyo Tel: 03 5796 7300
Tokyo Fax: 03 5796 7315

Korea

SIGMA-ALDRICH KOREA
Tel: 031-329-9000
Fax: 031-329-9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: 603-56353321
Fax: 603-56354116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01-800-007-5300
Free Fax: 01-800-712-9920

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Tel Gratis: 0800-0229088
Fax Gratis: 0800-0229089
Tel: 078-6205411
Fax: 078-6205421

New Zealand

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 0800 936 666
Free Fax: 0800 937 777

Norway

SIGMA-ALDRICH NORWAY AS
Tel: 23 17 60 60
Fax: 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: (+61) 829 01 00
Fax: (+61) 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 20 21 80
Free Fax: 800 20 21 78

Russia

SIGMA-ALDRICH RUSSIA
OOO SAF-LAB
Tel: 095 975-1917/3321
Fax: 095 975-4792

Singapore

SIGMA-ALDRICH PTE. LTD.
Tel: 65-6271 1089
Fax: 65-6271 1571

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD.
Tel: 27 11 979 1188
Fax: 27 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA S.A.
Free Tel: 900101376
Free Fax: 900102028

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: 020-350510
Fax: 020-352522
Outside Sweden Tel: +46 8 7424200
Outside Sweden Fax: +46 8 7424243

Switzerland

FLUKA CHEMIE GmbH
Swiss Free Call: 0800 80 00 80
Tel: +41 81 755 2828
Fax: +41 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717181
Free Fax: 0800 378785
Tel: 01747 833000
Fax: 01747 833313

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-Free: 800-325-3010
Call Collect: 314-771-5750
Toll-Free Fax: 800-325-5052
Tel: 314-771-5765
Fax: 314-771-5757

Internet:

sigma-aldrich.com

Solutions for Oligonucleotide Synthesis

Precisely Blended Capping, Deblocking, and Activator solutions for oligonucleotide synthesis.

Catalog No.	Product Name	Unit	Price
55,534-7	Cap Mix A (Contains 90% tetrahydrofuran; 10% acetic anhydride)	1L	\$65.10
		2L	100.00
55,531-2	Cap Mix A, with 2,6-lutidine (Contains 80% tetrahydrofuran; 10% acetic anhydride; 10% 2,6-lutidine)	1L	102.50
		2L	164.00
55,533-9	Cap Mix A, with pyridine (Contains 80% tetrahydrofuran; 10% acetic anhydride; 10% pyridine)	1L	82.10
		2L	129.00
55,532-0	Cap Mix B (Contains 84% tetrahydrofuran; 16% 1-methylimidazole)	1L	112.50
		2L	181.50

Catalog No.	Product Name	Unit	Price
55,535-5	Cap Mix B, with pyridine (Contains 80% tetrahydrofuran; 10% pyridine; 10% 1-methylimidazole)	1L	\$109.00
		2L	175.00
56,193-2	Deblock (Contains 3% trichloroacetic acid in dichloromethane)	1L	85.90
		2L	135.50
55,404-9	Activator (1 <i>H</i> -Tetrazole, 3 wt. % solution in acetonitrile)	1L	237.00
		2L	368.50
63,457-3	Activator (5-Ethylthio-1 <i>H</i> -tetrazole, 0.25 <i>M</i> in acetonitrile)	1L	230.00
		2L	347.50

For additional preblended solvents or larger quantities in 18-, 200-, or 400-L returnable stainless steel Pure-Pac™ tanks, please email dzagrodnik@sial.com.

Pure-Pac is a trademark of Sigma-Aldrich Biotechnology, L.P.


SIGMA-ALDRICH
A MILLIPORE COMPANY

DISCOVERY

Drug Discovery

PRECLINICAL

Toxicology

CLINICAL

Phase I

Phase II

Phase III

Phase IV

Drug Approval

COMMERCIALIZATION

Sigma-Aldrich Fine Chemicals provides products and contract manufacturing services, from Discovery through Commercialization, to the Pharmaceutical, Biopharmaceutical, and High-Technology industries.

We manufacture over 40,000 products.
What can we manufacture for you?



Take a Virtual Tour of our MANUFACTURING facilities at

<http://www.sigma-aldrich.com/tours>

To discuss how our expertise can benefit your next scale-up project or to obtain a quote, call us at 800-336-9719 (USA) or visit us at www.sigma-aldrich.com/safc

France: tel. 04 74 82 28 82 - fax 04 74 82 28 90 • Germany: tel. 089 6513 1900 - fax 089 6513 1919
Italy: tel. 02 33 41 7350 - fax 02 33 41 7201 • Japan: tel. 03 5821 3171 - fax 03 5821 3170 • Switzerland: tel. 081 755 2732 - fax 081 755 2770
U.K.: tel. 01202 712222 - fax 01202 740171 • USA: tel. 800 336 9719 or 314 534 4900 - fax 314 652 0000 • Email: safinechem@sial.com

sigma-aldrich.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
SIGMA-ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI • USA


SIGMA-ALDRICH
FINE CHEMICALS

ALDRICH CHEMICAL COMPANY, INC.
P.O. BOX 14508
ST. LOUIS, MISSOURI 63178 USA



CHANGE SERVICE REQUESTED

PRSR STD
U.S. POSTAGE PAID
MADISON, WI
PERMIT NO. 2223

GVC

DENDRIMERS: BUILDING BLOCKS FOR NANOSCALE SYNTHESIS

Aldrichimica ACTA

VOL. 37, NO. 2 • 2004



Birth of a New Macro-
molecular Architecture:
Dendrimers as Quantized
Building Blocks for Nanoscale
Synthetic Organic Chemistry

Recent Advances in
Microwave-Assisted
Synthesis


sigma-aldrich.com



ALDRICH®

New Products from Aldrich R&D

2-Pyridinylboronate (2,2'-bipyridylimino)diethanol ester

64,728-4		1 g	\$18.00
		5 g	71.30
		10 g	128.60


A stable 2-pyridinylboronate that is suitable for Suzuki-Miyaura coupling reactions.¹ This air-stable solid is a practical alternative to unstable 2-pyridinylboronic acids or esters.

(1) Hodgson, P. B.; Salingue, F. H. *Tetrahedron Lett.* **2004**, *45*, 685. (2) Fischer, F. C.; Havinga, E. *Recl. Trav. Chim. Pays-Bas* **1974**, *93*, 21. (3) Thompson, W. J. et al. *J. Org. Chem.* **1988**, *53*, 2052.

N,N-Diethylsacilamide, 97%


64,423-4		5g	\$27.50
		25g	95.50

N-(2-Hydroxybenzoyl)pyrrolidine, 97%

64,293-2		5g	\$19.00
		25g	65.00

These ligands have recently been utilized in copper-catalyzed aminations of aryl bromides with primary alkylamines under mild, solvent-free conditions and are compatible with a variety of functional groups. High yields have been obtained regardless of the ortho-, meta-, or para-substitution pattern on the aryl bromide. Kwong, F. Y.; Buchwald, S. L. *Org. Lett.* **2003**, *5*, 793.


Di-tert-butylmethylphosphine tetrafluoroborate

64,377-7		1g	\$28.00
		5g	98.50


This reagent is a convenient substitute for the air-sensitive di-tert-butylmethylphosphine (Aldrich Cat. No. **64,262-9**), and displays a remarkable reactivity in palladium-catalyzed Suzuki cross-coupling reactions.

Kirchoff, J. H. et al. *J. Am. Chem. Soc.* **2002**, *124*, 13662.

(1*R*,2*S*)-1-[(2-Hydroxy-3,5-di-tert-butylbenzylidene)]-aminoindan-2-ol

64,295-9		1g	\$16.50
		5g	58.00


(1*R*,2*S*)-1-[(2-Hydroxy-3,5-di-tert-butylbenzylidene)]-aminoindan-2-ol

64,414-5		1g	\$16.50
		5g	58.00

These chiral ligands have recently been used in the vanadium-catalyzed oxidation of sulfides in conjunction with hydrogen peroxide. Optically active sulfoxides are obtained in high yields and selectivities.

Pelotier, B. et al. *Synlett* **2002**, 1055.


4-[2-(4-Bromophenylsulfanyl)ethyl]pyridine, 97%

64,223-1		1g	\$16.50
		5g	58.00

The pyridinylethyl group protects the phenylthio group for subsequent modification. The later-freed thio group in phenylethynyl oligomers was used to anchor self-assembled monolayers (SAMs)¹ and molecular wires.²


(1) Collman, J. P. et al. *Synthesis* **2001**, 367. (2) YU, C. J. et al. *J. Org. Chem.* **1999**, *64*, 2070.

N,N',N''-Trihydroxyisocyanuric acid dimethylformamide complex, 97%

64,341-6		1g	\$35.00
		5g	119.50


An efficient catalyst for the aerobic oxidation of alkylbenzenes. Hirai, N. et al. *J. Org. Chem.* **2003**, *68*, 6587.

Bromobis(triphenylphosphine)(N-succinimide)palladium(II)


64,374-2		250mg	\$39.00
		1g	97.50

This bromosuccinimide-based palladium catalyst is used for Stille cross-coupling of allylic and benzylic halides. Crawforth, C. M. et al. *Chem. Commun.* **2003**, 2194.


4-Methoxy-N,N-diphenylaniline

64,612-1		1g	\$14.00
		10g	77.50


3-Methoxy-N,N-diphenylaniline, 97%

64,054-9		1g	\$12.00
		5g	37.50


4-Bromo-N,N-diphenylaniline, 97%

64,383-1		5g	\$34.50
		25g	115.50

4-(N,N-Diphenylamino)phenylboronic acid, 95%

64,729-2		1g	\$52.00
		5g	205.00

4-(N,N-Diphenylamino)benzaldehyde

64,720-9		1g	\$22.00
		5g	76.50

A substantial amount of research has been undertaken in the development of organic-based materials exhibiting metallic properties. Triarylamines have been extensively employed in electroluminescent materials as hole-transport materials and hole-transport emitters.¹ Triarylamines have also been used in the development of molecule-based magnets.^{2,3} The compounds above are some of the new triarylamines that Aldrich now offers.

(1) Mitschke, U.; Bäuerle, P. *J. Mater. Chem.* **2000**, *10*, 1471. (2) Miller, J. S.; Epstein, A. J. *MRS Bull.* **2000**, November, 21. (3) Veciana, J.; Iwamura, H. *MRS Bull.* **2000**, November, 41.

Please see pages 62-63 for additional new products.

For competitive quotes on larger quantities contact www.sigma-aldrich.com/safc.

SIGMA-ALDRICH
FINE CHEMICALS

Aldrichimica ACTA

VOL. 37, NO. 2 • 2004

Aldrich Chemical Co., Inc.
1001 West Saint Paul Ave.
Milwaukee, WI 53233, USA

To Place Orders

Telephone 800-558-9160 (USA)

FAX 800-962-9591 (USA)

or 414-273-4979

Mail P.O. Box 2060

Milwaukee, WI 53201 USA

Customer & Technical Services

Customer Inquiries 800-558-9160

Technical Service 800-231-8327

Sigma-Aldrich Fine Chemicals 800-336-9719

Custom Synthesis 800-336-9719

Flavors & Fragrances 800-227-4563

International 414-273-3850

24-Hour Emergency 414-273-3850

Website <http://www.sigma-aldrich.com>

Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.

P.O. Box 355, Milwaukee, WI 53201, USA

Subscriptions

To request your FREE subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-558-9160 (USA)

Mail: Attn: Mailroom

Aldrich Chemical Co., Inc.

P.O. Box 355

Milwaukee, WI 53201-9358

Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the inside back cover.

The *Aldrichimica Acta* is also available on the Internet at <http://www.sigma-aldrich.com>.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich family. ©2004 by Sigma-Aldrich Co.


"PLEASE BOTHER US."

Donald A. Tomalia

Joe Porwoll, President

Dr. Karl Pichler of Nanosolar, Inc. kindly suggested that we offer α,ω -dihexylsexithiophene (DH-6T). Oligothiophenes have attracted a lot of attention because of their high charge mobilities and on/off ratios as p-type semiconductors.¹⁻⁵ Garnier et al. fabricated the first all-organic transistor based on α -sexithiophene (6T).¹ Substituting the terminal α -hydrogen atoms on sexithiophene with alkyl chains increases the charge-carrier mobility further, as well as the solubility and thus the ease of processing

¹ Garnier, F. et al. *Adv. Mater.* **1990**, 2, 592. (2) Dodabalapur, A. et al. *Science* **1995**, 268, 270. (3) Garnier, F. *Acc. Chem. Res.* **1999**, 32, 239. (4) Sakamoto, Y. et al. *J. Am. Chem. Soc.* **2001**, 123, 4643. (5) Locklin, J. et al. *Lanhamur* **2002**, 18, 877. (6) Garnier, F. et al. *Adv. Mater.* **1993**, 115, 8716.



59,468-7

63,321-6

59,468-7 α -Sexithiophene (6T)

1g \$80.00

63,321-6 α,ω -Dihexylsexithiophene (DH-6T)

500mg \$93.50

Naturally, we made both α -sexithiophene and α,ω -dihexylsexithiophene for applications in various fields. It was no bother at all, just a pleasure to be able to help

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page or on the inside back cover.

TABLE OF CONTENTS

Birth of a New Macromolecular Architecture: Dendrimers as Quantized Building Blocks for Nanoscale Synthetic Organic Chemistry 39
Donald A. Tomalia, Dendritic NanoTechnologies, Inc.

Recent Advances in Microwave-Assisted Synthesis 66
Brittany L. Hayes, CEM Corporation

ABOUT OUR COVER

The Harbor at Lorient (oil on canvas, 43.5 × 73.0 cm) was painted by the French artist Berthe Morisot in 1869. In the right foreground of the picture, a young woman holding a parasol is shown seated on a low stone wall bordering the waters of the harbor, in which ships and boats of various sizes are anchored. The strong diagonal of the wall leads our eye straight to the figure of the woman, making her the obvious focal point of the painting. In the background can be seen the buildings of Lorient, a busy port on the south coast of Brittany on the estuary of the river Blavet



Photograph © Board of Trustees, National Gallery of Art, Washington

Morisot came to visit the town in 1869 to stay with her sister Edma, whose husband was a naval officer stationed there, and it is Edma who posed for the painting.

Morisot had studied with the painter Camille Corot, whose influence can be seen in the delicate palette and soft, atmospheric effects of this painting. In 1868, she began to study with the painter Edouard Manet. She came to know the impressionist painters, with whom she shared a preference for subjects from everyday life and for painting directly from nature to capture as much as possible the immediacy of visual experience. Beginning with the first group exhibition of these artists in 1874, she exhibited in all but one of the eight impressionist exhibitions held between that date and the last, in 1886. She married Manet's younger brother Eugène in 1874, but despite this close family tie, Morisot never convinced Manet to exhibit with the impressionist group. She had a uniquely important effect on the art of her brother-in-law, however, as it was under her influence that Manet abandoned the use of black, and adopted a lighter and more colorful impressionist palette.

This painting is in the Ailsa Mellon Bruce Collection at the National Gallery of Art, Washington, DC.

Queen's University Honors **Alfred Bader on His 80th Birthday**

Alfred Bader—chemist, art collector, philanthropist, and cofounder of Aldrich Chemical Company—has had a long association with Queen's University (Kingston, Ontario). He was a student there from 1941 to 1947, and has since been one of its most loyal supporters and generous benefactors.

It was no surprise, therefore, that Queen's University invited a group of distinguished guests and longtime friends of Alfred to participate in a celebration of his 80th birthday. The two-day gala, on May 12 and 13 of this year, combined festive social events with scholarly lectures organized by the chemistry and art departments at the university, as well as a public lecture by Dr.

Bader (*The Aldrich Chemical Company Story*) and a student awards luncheon. Nobel Laureate Barry Sharpless and Columbia Professor Emeritus Gilbert Stork were speakers at the chemistry symposium. In addition, the

university renamed Queen's Crescent, a street running through the center of campus, Bader Lane in gratitude to Alfred and his family's many contributions to the university over the years.

Alfred, his wife Isabel, and his two sons, Daniel and David, attended the celebrations. Sigma-Aldrich was represented by Henry van Oudenaren (Country Manager, Canada). Some highlights of these events are captured in the accompanying photographs.



Alfred Bader with mementos of the event



Alfred Bader at the student awards luncheon



Left to right: Barry Sharpless, Alfred Bader, Gilbert Stork, and Victor Snieckus at the chemistry symposium



Left to right: Henry van Oudenaren (Sigma-Aldrich), Gilbert Stork, and Alfred Bader

**Alfred, Your Sigma-Aldrich Family Wishes You a Happy 80th Birthday.
We Look Forward to Celebrating Your 90th!**

Sigma-Aldrich Worldwide Locations

Argentina
SIGMA-ALDRICH DE ARGENTINA, S.A.
Tel: 54 11 4556 1472
Fax: 54 11 4552 1698

Australia
SIGMA-ALDRICH PTY., LIMITED
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: 612 9841 0555
Fax: 612 9841 0500

Austria
SIGMA-ALDRICH HANDELS GmbH
Tel: 43 1 605 81 10
Fax: 43 1 605 81 20

Belgium
SIGMA-ALDRICH NV/SA.
Free Tel: 0800-14747
Free Fax: 0800-14745
Tel: 03 899 13 01
Fax: 03 899 13 11

Brazil
SIGMA-ALDRICH BRASIL LTDA
Tel: 55 11 3732-3100
Fax: 55 11 3733-5151

Canada
SIGMA-ALDRICH CANADA LTD.
Free Tel: 800-565-1400
Free Fax: 800-265-3858
Tel: 905-829-9500
Fax: 905-829-9292

China
SIGMA-ALDRICH CHINA INC.
Tel: 86-21-6386 2766
Fax: 86-21-6386 3966

Czech Republic
SIGMA-ALDRICH s.r.o.
Tel: 246 003 200
Fax: 246 003 291

Denmark
SIGMA-ALDRICH DENMARK A/S
Tel: 43 56 59 10
Fax: 43 56 59 05

Finland
SIGMA-ALDRICH FINLAND
Tel: 09-350-9250
Fax: 09-350-92555

France
SIGMA-ALDRICH CHIMIE S.à.r.l.
Tel appel gratuit: 0800 211 408
Fax appel gratuit: 0800 031 052

Germany
SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800-51 55 000
Free Fax: 0800-649 00 00

Greece
SIGMA-ALDRICH (O.M.) LTD
Tel: 30 210 9948010
Fax: 30 210 9943831

Hungary
SIGMA-ALDRICH Kft
Tel: 06-1-235-9054
Fax: 06-1-269-6470
Ingyenes zold telefon: 06-80-355-355
Ingyenes zold fax: 06-80-344-344

India
SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone:
Bangalore: 91-80-5112-7272
New Delhi: 91-11-2616 5477
Mumbai: 91-22-2570 2364
Hyderabad: 91-40-5531 5548
Fax:
Bangalore: 91-80-5112-7473
New Delhi: 91-11-2616 5611
Mumbai: 91-22-2579 7589
Hyderabad: 91-40-5531 5466

Ireland
SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: 353 1 4041900
Fax: 353 1 4041910

Israel
SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1-800-700-2222
Tel: 08-948-4100
Fax: 08-948-4200

Italy
SIGMA-ALDRICH S.r.l.
Telefono: 02 33417310
Fax: 02 38010737
Numero Verde: 800-827018

Japan
SIGMA-ALDRICH JAPAN K.K.
Tokyo Tel: 03 5796 7300
Tokyo Fax: 03 5796 7315

Korea
SIGMA-ALDRICH KOREA
Tel: 031-329-9000
Fax: 031-329-9090

Malaysia
SIGMA-ALDRICH (M) SDN. BHD
Tel: 603-56353321
Fax: 603-56354116

Mexico
SIGMA-ALDRICH QUÍMICA, S.A. de C.V
Free Tel: 01-800-007-5300
Free Fax: 01-800-712-9920

The Netherlands
SIGMA-ALDRICH CHEMIE BV
Tel Gratis: 0800-0229088
Fax Gratis: 0800-0229089
Tel: 078-6205411
Fax: 078-6205411

New Zealand
SIGMA-ALDRICH PTY., LIMITED
Free Tel: 800 936 666
Free Fax: 800 937 777

Norway
SIGMA-ALDRICH NORWAY AS
Tel: 23 17 60 60
Fax: 23 17 60 50

Poland
SIGMA-ALDRICH Sp. z o.o.
Tel: (+61) 829 01 00
Fax: (+61) 829 01 20

Portugal
SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 20 21 80
Free Fax: 800 20 21 78

Russia
SIGMA-ALDRICH RUSSIA
OOO SAF-LAB
Tel: 095 975-1917/3321
Fax: 095 975-4792

Singapore
SIGMA-ALDRICH PTE. LTD.
Tel: 65-6271 1089
Fax: 65-6271 1571

South Africa
SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: 27 11 979 1188
Fax: 27 11 979 1119

Spain
SIGMA-ALDRICH QUÍMICA S.A.
Free Tel: 900101376
Free Fax: 900102028

Sweden
SIGMA-ALDRICH SWEDEN AB
Tel: 020-350510
Fax: 020-352522
Outside Sweden Tel: +46 8 7424200
Outside Sweden Fax: +46 8 7424243

Switzerland
FLUKA CHEMIE GmbH
Swiss Free Call: 0800 80 00 80
Tel: +41 81 755 2828
Fax: +41 81 755 2815

United Kingdom
SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717181
Free Fax: 0800 378785
Tel: 01747 833000
Fax: 01747 833133

United States
SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-free: 800-325-3010
Call Collect: 314-771-5750
Toll-Free Fax: 800-325-5052
Tel: 314-771-5765
Fax: 314-771-5757
Internet:
sigma-aldrich.com

Solutions for Oligonucleotide Synthesis

Precisely Blended Capping, Deblocking, and Activator solutions for oligonucleotide synthesis.

Catalog No.	Product Name	Unit	Price	Catalog No.	Product Name	Unit	Price
55,534-7	Cap Mix A (Contains 90% tetrahydrofuran: 10% acetic anhydride)	1L 2L	\$65.10 100.00	55,535-5	Cap Mix B, with pyridine (Contains 80% tetrahydrofuran: 10% pyridine; 10% 1-methylimidazole)	1L 2L	\$109.00 175.00
55,531-2	Cap Mix A, with 2,6-lutidine (Contains 80% tetrahydrofuran: 10% acetic anhydride; 10% 2,6-lutidine)	1L 2L	102.50 164.00	56,193-2	Deblock (Contains 3% trichloroacetic acid in dichloromethane)	1L 2L	85.90 135.50
55,533-9	Cap Mix A, with pyridine (Contains 80% tetrahydrofuran: 10% acetic anhydride; 10% pyridine)	1L 2L	82.10 129.00	55,404-9	Activator (1 <i>H</i> -Tetrazole, 3 wt. % solution in acetonitrile)	1L 2L	237.00 368.50
55,532-0	Cap Mix B (Contains 84% tetrahydrofuran: 16% 1-methylimidazole)	1L 2L	112.50 181.50	63,457-3	Activator (5-Ethylthio-1 <i>H</i> -tetrazole, 0.25 <i>M</i> in acetonitrile)	1L 2L	230.00 347.50

For additional preblended solvents or larger quantities in 18-, 200-, or 400-L returnable stainless steel Pure-Pac™ tanks, please call 800-213-1206 (USA) or email dzagrodnik@sial.com.

Pure-Pac is a trademark of Sigma-Aldrich Biotechnology, L.P.


SIGMA-ALDRICH
FINE CHEMICALS

DISCOVERY

Drug Discovery

PRECLINICAL

Toxicology

CLINICAL

Phase I

Phase II

Phase III

Phase IV

Drug Approval

COMMERCIALIZATION

Sigma-Aldrich
acquires
Ultrafine
to enhance
Fine Chemicals
capabilities

Sigma-Aldrich Fine Chemicals provides products and contract manufacturing services, from Discovery through Commercialization, for Pharmaceutical, Biopharmaceutical, and Materials Science applications.

**We manufacture over 40,000 products
and can provide you with:**

NEW!

Rapid, Pre-Clinical Custom Synthesis
Clinical cGMP Custom Manufacturing
Commercial cGMP Custom Manufacturing

We look forward to discussing your manufacturing needs!

To discuss how our expertise can benefit your next scale-up project or to obtain a quote,
call us at 800-336-9719 (USA) or visit us at www.sigma-aldrich.com/safc

France: tel. 04 74 82 28 82 - fax 04 74 82 28 90 • Germany: tel. 089 6513 1900 - fax 089 6513 1919
Italy: tel. 02 33 41 7350 - fax 02 33 41 7201 • Japan: tel. 03 5821 3171 - fax 03 5821 3170 • Switzerland: tel. 081 755 2732 - fax 081 755 2770
U.K.: tel. 01202 712222 - fax 01202 740171 • USA: tel. 800 336 9719 or 314 534 4900 - fax 314 652 0000 • Email: safinechem@sial.com

sigma-aldrich.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
SIGMA-ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI • USA



SIGMA-ALDRICH
FINE CHEMICALS

ALDRICH CHEMICAL COMPANY, INC.
P.O. BOX 14508
ST. LOUIS, MISSOURI 63178 USA



CHANGE SERVICE REQUESTED

PSRST STD
U.S. POSTAGE PAID
MADISON, WI
PERMIT NO. 2223

HEX

ALDRICH CONGRATULATES THE 2005 ACS AWARD WINNERS

Aldrichimica **ACTA**

VOL. 37, NO. 3 • 2004

(Last issue in 2004)



**Palladium-Mediated
Synthesis of Aldehydes and
Ketones from Thiol Esters**

**Hydroalumination Reactions
in Organic Chemistry**

sigma-aldrich.com

 **ALDRICH®**

New Products from Aldrich R&D


(1S)-2-Methyl-CBS-oxazaborolidine

64,930-9		1 g	\$35.00
		10 g	175.00

(1S)-2-Methyl-CBS-oxazaborolidine, 1M solution in toluene

45,770-1		5 mL	\$58.30
		25 mL	187.00

(1R)-2-Methyl-CBS-oxazaborolidine

64,931-7		1 g	\$38.00
		10 g	184.00


(1R)-2-Methyl-CBS-oxazaborolidine, 1M solution in toluene

45,769-8		5 mL	\$58.30
		25 mL	187.00

Excellent catalysts for the asymmetric reduction of prochiral ketones.¹ Reduction of aralkyl ketones with borane sources such as BH₃-THF, BMS, amine-boranes, catecholborane, in the presence of these catalysts (0.02 to 0.1 mol %), provides the corresponding secondary alcohols in high (>95% ee) enantiomeric excesses.²⁻⁵ These reagents are also effective for kinetic resolutions.⁶

(1) Corey, E. J., Helal, C. J. *Angew. Chem., Int. Ed.* **1998**, 37, 1986. (2) Corey, E. J. et al. *J. Am. Chem. Soc.* **1987**, 109, 5551. (3) Corey, E. J. et al. *J. Am. Chem. Soc.* **1987**, 109, 7925. (4) Cho, B. T.; Kim, D. *J. Tetrahedron Asymmetry* **2001**, 12, 2043. (5) Cho, B. T. et al. *J. Chem. Soc., Perkin Trans. 1* **2001**, 1204. (6) Lebsack, A. D. et al. *J. Am. Chem. Soc.* **2001**, 123, 4851.


2-Pyridylboronic acid N-phenyldiethanolamine ester

64,728-4		1 g	\$18.00
		5 g	71.30
		10 g	128.60


3-Pyridylboronic acid N-phenyldiethanolamine ester

64,884-1		5 g	\$63.00
-----------------	--	-----	---------


4-Pyridylboronic acid N-phenyldiethanolamine ester

64,882-5		500 mg	\$27.50
		1 g	46.00
		5 g	181.00

5-Pyridylboronic acid N-phenyldiethanolamine ester

64,881-7		1 g	\$16.50
		5 g	66.00


6-Methoxy-4-pyridylboronic acid N-phenyldiethanolamine ester

64,915-5		1 g	\$43.50
		5 g	172.00

Novel 2-pyridylboronic acid derivatives that are stable for storage and suitable for Suzuki-Miyaura cross-coupling.

(1) Fischer, F. C.; Havinga, E. *Recl. Trav. Chim. Pays-Bas* **1974**, 93, 21. (2) Bouillon, A. et al. *Tetrahedron* **2003**, 59, 10043. (3) Hodgson, P. B.; Salinque, F. H. *Tetrahedron Lett.* **2004**, 45, 685.


2-Chloro-6-(methylamino)purine, 97%

64,459-5		5 g	\$59.50
-----------------	---	-----	---------

Has been employed in the preparation of selective P2Y₁ receptor agonists and antagonists as pharmacological probes.¹⁻³

(1) Kim, Y.-C. et al. *J. Med. Chem.* **2000**, 43, 746. (2) Nandanani, E. et al. *J. Med. Chem.* **2000**, 43, 829. (3) Kim, H. S. et al. *J. Med. Chem.* **2001**, 44, 3692.


Benzofurazan, 97%

65,013-7		1 g	\$17.50
		10 g	95.00

Modifications around this key heterocyclic nucleus have led to the preparation of fluorescence probes¹ and inhibitors of nucleic acid and protein biosyntheses.^{2,3}

(1) Uchiyama, S. et al. *J. Chem. Soc., Perkin Trans. 2* **1999**, 2525. (2) Crampton, M. R. et al. *J. Chem. Soc., Perkin Trans. 2* **2002**, 257. (3) Ghosh, P. B. et al. *J. Med. Chem.* **1972**, 15, 255.


2,3-Dibromobenzo[b]thiophene, 97%

64,737-3		1 g	\$29.50
		5 g	97.50

Served as a basic starting material for the synthesis of novel, chiral benzothio-phenone ligands via Pd-catalyzed allylation.¹ Selective displacement of each bromine makes this compound a very useful building block in Suzuki and Stille coupling reactions.^{2,3}

(1) Tietze, L. F.; Lohmann, J. K. *Synlett* **2002**, 2083. (2) Heynderickx, A. et al. *Synthesis* **2002**, 213. (3) Yamamura, K. et al. *Org. Biomol. Chem.* **2004**, 2, 1413.


2-Chloro-1,1,1-triethoxyethane

64,975-9		1 g	\$13.50
		10 g	75.00

An effective reagent for the large-scale preparation of 2-chloromethylbenzothiazoles and benzoxazoles.¹ Used also in the synthesis of 21-desoxy-21-chlorocorticosteroids with functionalized 17 α -ester groups.²

(1) Mylari, B. L. et al. *Synth. Commun.* **1989**, 19, 2921. (2) Ueno, H. et al. *J. Med. Chem.* **1991**, 34, 2468.


(Chloromethyl)dimethylvinylsilane, 97%

64,911-2		1 g	\$17.30
		10 g	95.00

Utilized in the synthesis of silylated β -lactams possessing antibacterial properties.

Altamura, M. et al. *J. Org. Chem.* **1995**, 60, 8403.


2-(Trimethylsilyl)pyridine, 98%

64,985-6		1 g	\$25.50
		10 g	140.00

Utilized in the preparation of a novel hydroxymethylating reagent, which was then used in the stereoselective synthesis of vinylsilanes.

Itami, K. et al. *Tetrahedron* **2001**, 57, 5045.

4,4'-Bis(triethoxysilyl)biphenyl, 95%

63,810-2		1 g	\$23.50
		10 g	128.50

Employed in the preparation of arylene-bridged polysilsesquioxanes.^{1,2}

(1) Kapoor, M. P. et al. *Chem. Lett.* **2003**, 32, 914. (2) Schaefer, D. W. et al. *Chem. Mater.* **2004**, 16, 1402.

Please see pages 100-101 for additional new products.

For competitive quotes on larger quantities, please contact safcglobal.com.

SAFC

Aldrichimica ACTA

VOL. 37, NO. 3 • 2004

Aldrich Chemical Co., Inc.
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-558-9160 (USA)
FAX 800-962-9591 (USA)
or 414-273-4979
Mail P.O. Box 2060
Milwaukee, WI 53201 USA

Customer & Technical Services

Customer Inquiries 800-558-9160
Technical Service 800-231-8327
SAFC 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 414-273-3850
24-Hour Emergency 414-273-3850
Web site sigma-aldrich.com
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201, USA

Subscriptions

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-558-9160 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
P.O. Box 355
Milwaukee, WI 53201-9358
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the inside back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.


Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich family. ©2004 by Sigma-Aldrich Co.

"PLEASE BOTHER US."

Editorial President

Professor Yi-Ping Wang of Shanghai Jiao Tong University kindly suggested that we offer the CLD-1 chromophore in a more convenient-to-use form. CLD-1 is a high $\mu\beta$ chromophore that possesses excellent thermal stability, low optical loss, and low modulation voltage.^{1,2} In addition to the powder form of CLD-1, Sigma-Aldrich now offers CLD-1 as a ready-to-cast solution with an amorphous polycarbonate (APC) support

¹ Zhang, C., Todorova, G., Wang, C., Londergan, T., Dalton, L. R. *Proc SPIE* **2000**, 4114, 77. ² Oh, M.-C., Zhang, H., Zhang, Erlig, H., Chang, Y., Tsap, Y., Chang, D., Szep, A., Steier, W. H., Fetterman, H. R., Dalton, L. R. *IEEE J Sel Top Quantum Electron* **2001**, 7, 100.



65,158-3 CLD-1-APC polymer, 10 wt. % solution in dichloroethane **5 mL \$385.00**
64,413-7 CLD-1 **100 mg \$310.00**

Naturally, we made this useful solution. It was no bother at all, just a pleasure to be able to help

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page or on the inside back cover.

TABLE OF CONTENTS

Palladium-Mediated Synthesis of Aldehydes and Ketones from Thiol Esters..... 87

Tohru Fukuyama* and Hidetoshi Tokuyama, The University of Tokyo

Hydroalumination Reactions in Organic Chemistry.....103

Herbert W. Roesky, Institut für Anorganische Chemie der Universität Göttingen

ABOUT OUR COVER

The Race Track (oil on canvas, 80.5 × 114.9 cm) was painted about 1891 by the French artist Jean-Louis Forain. Just to the right of center in the foreground of the painting, two men and a woman converse intently. A little further back and near the left edge of the picture, several figures appear to be rushing towards the horses and riders. The relative position and differing size—some larger, some smaller—of these groups of figures help to create a sense of space in the picture, which otherwise is dominated by large relatively flat areas of brilliant color, the light green of the grass, the white of the clouds, and the intense blue of the sky



Photograph © Board of Trustees, National Gallery of Art, Washington

At first glance, one might identify the artist of this painting as Degas, an artist well known for horseracing subjects. Not only the subject, however, but also the style of the painting may remind us of Degas, as it exhibits the same casual, seemingly unplanned composition that is usual in the work of the more famous artist. It is as if the figures in the painting were caught unawares in a snapshot shortly before the beginning of the race. Indeed, Forain had a chameleon-like tendency to imitate other artists he admired. Not only did he paint racehorses and ballet dancers in a manner that reminds one of Degas, but his wily lawyers and scheming politicians can also be easily mistaken for those of Daumier. And like Daumier, he was a graphic artist as well as a painter. In fact, it was his etchings that brought Forain his greatest fame. He was what we now might call a political cartoonist. Fine as his graphic technique was, however, his paintings almost always remind us of other, better-known artists. All his life, he remained something of a follower, rather than an inventive and original creator of his own unique style

This painting is in the Collection of Mr. and Mrs. Paul Mellon at the National Gallery of Art, Washington, DC.

PEPscreen™: Custom Peptide Libraries

High-Throughput Peptide Solution for Generating Fast Leads in Screening Applications

Sigma-Genosys has developed a proprietary peptide synthesis platform using state-of-the-art technology. This technology allows the fastest and most efficient high-throughput parallel synthesis of milligram quantities of peptide libraries. The PEPscreen™ service provides researchers with the ability to order large numbers of extremely affordable custom peptides delivered in less than 7 working days.

Features and Benefits

- Compatible with robotic automation for resuspension and dispensing in screening applications
- Ideal for performing solution-phase assays
- Suitable for robotic production of peptide microarrays
- N-Terminal modifications and nonstandard residue availability allows for maximum flexibility in project design
- Delivery of product in less than 7 working days allows for faster target identification

Modifications Available

- N-Terminal acetylation
- N-Terminal biotin
- N-Terminal fluorescein
- D-Amino acids
- Nonstandard residues
- Phosphorylated residues
- C-Terminal acid or amide

Product Specifications

- Quantity: 0.5–2 mg
- Length: 6–20 residues
- Purity: ~70% (based on 15 residues)
- Analysis: MALDI-TOF mass spectrometry on 100% of the peptides
- Format: lyophilized in 96-well format (2-D bar-coded tubes)
- Shipped with electronic PDF files of QC data and peptide sequence map
- Minimum order size is 48 peptides

Ordering Information

For more information, please refer to www.sigma-genosys.com/peptide_libraries.asp.
For technical inquiries, please contact peptides@sial.com.

Place your order by emailing your sequences in an 8-channel orientation (peptide #1 = A1, peptide #2 = B1, peptide #8 = H1, peptide #9 = A2, etc.).



PEPscreen™ is a trademark of Sigma-Aldrich Biotechnology, L.P.

sigma-genosys.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
SIGMA ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI 63178 • USA

SIGMA
GENOSYS

Sigma-Aldrich Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA, S.A.
Tel: 54 11 4556 1472
Fax: 54 11 4552 1698

Australia

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: 612 9841 0555
Fax: 612 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: 43 1 605 81 10
Fax: 43 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA
Free Tel: 0800-14747
Free Fax: 0800-14745
Tel: 03 899 13 01
Fax: 03 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Tel: 55 11 3732-3100
Fax: 55 11 3733-5151

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 800-565-1400
Free Fax: 800-265-3858
Tel: 905-829-9500
Fax: 905-829-9292

China

SIGMA-ALDRICH CHINA INC
Tel: 86-21-6386 2766
Fax: 86-21-6386 3966

Czech Republic

SIGMA-ALDRICH s.r.o.
Tel: 246 003 200
Fax: 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: 43 56 59 10
Fax: 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND
Tel: 09-3509250
Fax: 09-350-92555

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Tel appel gratuit: 0800 211 408
Fax appel gratuit: 0800 031 052

Germany

SIGMA-ALDRICH CHEMIE GmbH
Tel: 0800-51 55 000
Free Fax: 0800-649 00 00

Greece

SIGMA-ALDRICH (O.M.) LTD
Tel: 30 210 9948010
Fax: 30 210 9943831

Hungary

SIGMA-ALDRICH Kft
Tel: 06-1-235-9054
Fax: 06-1-269-6470
Ingyenes zold telefon: 06-80-355-355
Ingyenes zold fax: 06-80-344-344

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: 91-80-5112-7272
New Delhi: 91-11-2616 5477
Mumbai: 91-22-2570 2364
Hyderabad: 91-40-5531 5488
Fax
Bangalore: 91-80-5112-7473
New Delhi: 91-11-2616 5611
Mumbai: 91-22-2579 7589
Hyderabad: 91-40-5531 5466

Ireland

SIGMA-ALDRICH IRELAND LTD
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: 353 1 4041900
Fax: 353 1 4041910

Israel

SIGMA-ALDRICH ISRAEL LTD
Free Tel: 1-800-70-2222
Tel: 08-948-4100
Fax: 08-948-4200

Italy

SIGMA-ALDRICH S.r.l.
Telefono: 02 33417310
Fax: 02 38010737
Numero Verde: 800-827018

Japan

SIGMA-ALDRICH JAPAN K.K.
Tokyo Tel: 03 5796 7300
Tokyo Fax: 03 5796 7315

Korea

SIGMA-ALDRICH KOREA
Tel: 031-329-9000
Fax: 031-329-9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: 603-56353321
Fax: 603-56354116

Mexico

SIGMA-ALDRICH QUIMICA, S.A. de C.V.
Free Tel: 01-800-007-5300
Free Fax: 01-800-712-9920

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Tel Gratis: 0800-0229088
Fax Gratis: 0800-0229089
Tel: 078-6205411
Fax: 078-6205421

New Zealand

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 0800 936 666
Free Fax: 0800 937 777

Norway

SIGMA-ALDRICH NORWAY AS
Tel: 23 17 60 60
Fax: 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: (+61) 829 01 00
Fax: (+61) 829 01 20

Portugal

SIGMA-ALDRICH QUIMICA, S.A.
Free Tel: 800 202180
Free Fax: 800 202178
Tel: 21 9242555
Fax: 21 9242610

Russia

SIGMA-ALDRICH RUSSIA
OOO SAF-LAB
Tel: 095 975-1917/3321
Fax: 095 975-4792

Singapore

SIGMA-ALDRICH PTE. LTD
Tel: 65-67791200
Fax: 65-67791822

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: 27 11 979 1188
Fax: 27 11 979 1119

Spain

SIGMA-ALDRICH QUIMICA S.A.
Free Tel: 900 101376
Free Fax: 900 102028
Tel: 91 661 99 77
Fax: 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: 020-350510
Fax: 020-352522
Outside Sweden Tel: +46 8 7424200
Outside Sweden Fax: +46 8 7424243

Switzerland

FLUKA CHEMIE GmbH
Swiss Free Call: 0800 80 00 80
Tel: +41 81 755 2828
Fax: +41 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD
Free Tel: 0800 717181
Free Fax: 0800 378785
Tel: 01747 833000
Fax: 01747 833313

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-free: 800-325-3010
Call Collect: 314-771-5750
Toll-Free Fax: 800-325-5052
Tel: 314-771-5765
Fax: 314-771-5757

Internet

sigma-aldrich.com



Solutions for Oligonucleotide Synthesis

Precisely Blended Capping, Deblocking, and Activator solutions for oligonucleotide synthesis.

Catalog No.	Product Name	Unit	Price	Catalog No.	Product Name	Unit	Price
55,534-7	Cap Mix A (Contains 90% tetrahydrofuran: 10% acetic anhydride)	1 L	\$68.30	55,535-5	Cap Mix B, with pyridine (Contains 80% tetrahydrofuran: 10% pyridine; 10% 1-methylimidazole)	1 L	\$113.40
		2 L	105.00			2 L	184.00
55,531-2	Cap Mix A, with 2,6-lutidine (Contains 80% tetrahydrofuran: 10% acetic anhydride; 10% 2,6-lutidine)	1 L	106.60	56,193-2	Deblock (Contains 3% trichloroacetic acid in dichloromethane)	1 L	90.20
		2 L	172.00			2 L	142.50
55,533-9	Cap Mix A, with pyridine (Contains 80% tetrahydrofuran: 10% acetic anhydride; 10% pyridine)	1 L	86.20	55,404-9	Activator (1 <i>H</i> -Tetrazole, 3 wt. % solution in acetonitrile)	1 L	249.00
		2 L	135.50			2 L	387.00
55,532-0	Cap Mix B (Contains 84% tetrahydrofuran: 16% 1-methylimidazole)	1 L	118.00	63,457-3	Activator (5-Ethylthio-1 <i>H</i> -tetrazole, 0.25 <i>M</i> in acetonitrile)	1 L	241.50
		2 L	190.50			2 L	365.00

For additional preblended solvents or larger quantities in 18-, 200-, or 400-L returnable stainless steel Pure-Pac tanks, please call 800-325-3010 or visit us online at www.sigma-aldrich.com.

Pure-Pac is a trademark of Sigma-Aldrich Biotechnology, L.P.

SAFC

Introducing science beyond the ordinary



SAFC

Science. Advanced. Chemical.

SAFC requires only a few minutes of your time to help you understand the value of our services. We offer a range of services, including the full range of chemical synthesis, process development, and process optimization through our global network.

Over 75 years of chemistry know-how, process development and pilot-scale experience to bring your products to market.

Whether you need a small batch or large quantities.

Be inspired—by SAFC.

www.safcglobal.com

Member of the Sigma-Aldrich Group

ALDRICH CHEMICAL COMPANY, INC.
P.O. BOX 14508
ST. LOUIS, MISSOURI 63178 USA



CHANGE SERVICE REQUESTED

PRSR STD
U.S. POSTAGE PAID
MADISON, WI
PERMIT NO. 2223

HKY

NOVEL REAGENTS AND CATALYSTS FOR FACILITATING SYNTHESIS

Aldrichimica ACTA

VOL. 38, NO. 1 • 2005



ROM Polymerization
in Facilitated Synthesis

Polyurea-Encapsulated
Palladium Catalysts

sigma-aldrich.com




ALDRICH[®]

Advancing Science

New Products from Aldrich R&D

4-Diphenylphosphanybenzoic acid, 2-(trimethylsilyl)ethyl ester, 0.5 M in THF

63,801-3




5 mL \$28.50
25 mL 99.50

This reagent was developed to replace triphenylphosphine in the Mitsunobu reaction. The resulting phosphine oxide is easily removed from the reaction during workup.

Yoakim, C. et al. *Synlett* **2003**, 473

4-Iodotoluene difluoride

65,111-7




1 g \$31.50
5 g 108.00

A new and efficient fluorinating reagent that is recommended for the fluorination of the α carbon of sulfanyl amides under mild conditions.

Motherwell, W. B. et al. *J. Chem. Soc., Perkin Trans. 1* **2002**, 2816

mono-Methyl fumarate, 97%

65,141-9




1 g \$10.50
10 g 55.00

Employed in the design of affinity labels for opioid receptors¹ and in a highly stereoselective synthesis of optically active furfuryl fumarates.

(1) Chang, A.-C. et al. *J. Med. Chem.* **1994**, *37*, 4490. (2) Butz, T.; Sauer, J. *Tetrahedron Asymmetry* **1997**, *8*, 703

5-Bromo-6-bromomethylbenzo[1,3]dioxole

65,374-8




5 g \$29.50
25 g 105.00

Has been utilized in the preparation of toddaquinoline, an unusual medicinal alkaloid,¹ and lignans, which are known for their biological activity and effectiveness as antineoplastic agents.

(1) Harrowen, D. C. et al. *Tetrahedron* **2001**, *57*, 4447. (2) Cochran, J. E.; Padwa, A. *J. Org. Chem.* **1995**, *60*, 3938

Chlorodicyclopentylphosphine

64,906-6




1 g \$24.50
5 g 95.00
25 g 375.00

A phosphine precursor for ligand preparation in Negishi and Suzuki coupling reactions.

2,8,9-Triisopropyl-2,5,8,9-tetraaza-1-phosphabicyclo[3,3,3]undecane, 1 M in toluene

65,435-3




10 mL \$159.00

This highly useful and air-sensitive reagent is now available in an easier-to-handle toluene solution in Sure/Seal™ bottles.

Verkade, J. G.; Kisanga, P. B. *Aldrichimica Acta* **2004**, *37*, 3

N-(3-Bromophenyl)aniline, 97%

65,424-8




1 g \$13.50
10 g 74.50

A useful building block in the preparation of triaryl amines, which have been extensively employed in electroluminescent materials as hole-transport materials and hole-transport emitters,¹ and in the development of organic-based magnets.^{2,3}

(1) Mitschke, U.; Bäuerle, P. *J. Mater. Chem.* **2000**, *10*, 1471. (2) Miller, J. S.; Epstein, A. *J. MRS Bull.* **2000**, *25*, 21. (3) Veciana, J.; Iwamura, H. *MRS Bull.* **2000**, *25*, 41

Anilinium hypophosphite

65,411-6




25 g \$19.50

An easy-to-handle and relatively nonhygroscopic reagent that has been used in the preparation of monosubstituted phosphinic acids via Pd-catalyzed reaction with aryl halides or triflates.¹ A practical triethylborane-initiated radical addition of anilinium hypophosphite to olefins has been reported.²

(1) Montchamp, J.-L.; Dumond, Y. R. *J. Am. Chem. Soc.* **2001**, *123*, 510. (2) Deprele, S.; Montchamp, J.-L. *J. Org. Chem.* **2001**, *66*, 6745

2-(2'-Di-tert-butylphosphine)biphenylpalladium(II) acetate

65,541-4




250 mg \$28.50
1 g 75.00

A novel, air- and moisture-stable pre-catalyst for the amination of aryl chlorides.

Zim, D.; Buchwald, S. L. *Org. Lett.* **2003**, *5*, 2413

3,2':5',3''-Terthiophene

65,138-9




1 g \$105.00

An efficient singlet-oxygen sensitizer. This compound is also associated with photo-antibiotic and phototoxic properties.¹

(1) Beny, J.-P. et al. *J. Org. Chem.* **1982**, *47*, 2201. (2) Moriarty, R. M. et al. *Synth Commun.* **1985**, *15*, 789

N-Boc-2-aminomethylpyridine, 97%


65,157-5



5 g \$25.00
25 g 87.50

N-Boc-3-aminomethylpyridine, 97%


63,444-1



1 g \$21.00
5 g 73.50

N-Boc-4-aminomethylpyridine, 96%

64,976-7



1 g \$20.50
10 g 111.50

These Boc-protected pyridines can be used as pharmaceutical building blocks in thrombin¹ and PKC inhibitor² research, and can be readily hydrogenated to the protected piperidines.

(1) Hilpert, K. et al. *J. Med. Chem.* **1994**, *37*, 3889. (2) Shearer, B. G. et al. *J. Med. Chem.* **1991**, *34*, 2928

Please see pages 20–21 for additional new products.

For competitive quotes on larger quantities, contact safcglobal.com.

SAFC

Aldrichimica ACTA

VOL. 38, NO. 1 • 2005

Aldrich Chemical Co., Inc.
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-4979
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-231-8327
SAFC™ 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 414-438-3850
24-Hour Emergency 414-438-3850
Web site sigma-aldrich.com
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201, USA

Subscription

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
P.O. Box 355
Milwaukee, WI 53201-9358
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the inside back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich family. © 2005 Sigma-Aldrich Co.

"PLEASE BOTHER US."



Professor Ei-ichi Negishi, President

Professor Ei-ichi Negishi of Purdue University kindly suggested that we offer *trans*-1-bromo-2-iodoethene, a useful synthon in Pd-catalyzed cross-coupling reactions. Its utility has recently been demonstrated in the synthesis of xerulin, an inhibitor of cholesterol biosynthesis,¹ *cis*- and *trans*-bupleurynol,² in the preparation of 1-halo-1,3-dienes,³ and in coupling reactions with alkynylzinc reagents.⁴

(1) Negishi, E., Almardanov, A., Xu, C. *Org. Lett.* **2000**, *2*, 65 (2) Ghasemi, H., Antunes, L. M., Organ, M. G. *Org. Lett.* **2004**, *6*, 2913 (3) Qian, M., Huang, Z., Negishi, E. *Org. Lett.* **2004**, *6*, 1531. (4) Negishi, E., Qian, M.; Zeng, F.; Anastasia, L.; Babinski, D. *Org. Lett.* **2003**, *5*, 1597



57,780-4 trans-1-Bromo-2-iodoethene, 97% 5 g \$36.70
25 g 122.50

Naturally, we made this versatile building block. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the inside back cover.

TABLE OF CONTENTS

ROM Polymerization in Facilitated Synthesis 3
Andrew M. Harned, Mianji Zhang, Punitha Vedantham, Shubhashish Mukherjee, Russell H. Herpel, Daniel L. Flynn, and Paul R. Hanson,** The University of Kansas and Deciphera Pharmaceuticals, LLC

Polyurea-Encapsulated Palladium Catalysts: The Development and Application of a New and Versatile Immobilized-Homogeneous-Catalyst Technology 23
David A. Pears, Avecia Pharmaceuticals; and Stephen C. Smith, Syngenta

ABOUT OUR COVER

The Fence (oil on canvas, 37.8 × 45.7 cm) was signed and dated by the French painter Camille Jacob Pissarro in 1872. During the 1860s, Pissarro moved his family from Paris to various small villages in the French countryside. He was committed to the principles of socialism, and felt a strong affinity for the peasant farmers who worked the land. Like the other impressionists, Pissarro chose to represent subjects from modern life; but, while they often painted the pleasures of the urban bourgeoisie, scenes from the theatre, the opera, the ballet, or the racetrack, Pissarro was more likely to portray the rustic life of farmers working in the fields with whom he identified



Photograph © Board of Trustees, National Gallery of Art, Washington

The warm colors of the trees and shrubbery in the picture show that it is the fall of the year. At the left in the foreground is a large bent and broken tree, whose almost leafless branches are silhouetted against the light sky. Near the lower right corner of the painting, a peasant couple chat together on either side of a rustic fence. In the distance, one can make out the buildings of the local village, towards which a woman moves along the path near the right edge of the picture. As in a snapshot, nothing in nature has been rearranged to create a more pleasing, harmonious, or picturesque effect. The rapidly executed brushwork is variegated to suggest the different textures of diverse objects in the painting, which was almost certainly painted in a single session on the site. One might say, in fact, that a central purpose of all the impressionist artists was not to create an invented world on canvas, but to capture the immediacy of the unique conditions of a specific moment in time in a particular place

This painting is a part of the Collection of Mr. and Mrs. Paul Mellon at the National Gallery of Art, Washington, DC.

Introducing the New Sigma-Aldrich Interactive Virtual Stockroom...

Where products and information are literally at your fingertips!



Looking for information on solvents, LabBasics, EquipBasics, or general chemistry products? Interested in learning more about our E-commerce and on-site stocking options? The Virtual Stockroom has easily navigated screens that provide you with all of your product and information needs.

Navigate the Virtual Stockroom just as you would a real stockroom. Along with finding solvents in the flammables cabinet and reference materials on the bookshelf, you'll also find the following:

Products:


- Solvents arranged by grade and container type
- LabBasics products (commonly used chemical reagents)
- Aldrich custom glassware
- EquipBasics (commonly used laboratory equipment)

Technical Information:

- *Aldrichimica Acta* current and back issues
- Aldrich technical bulletins

Programs and Services:

- Chemical structure search
- Creating and maintaining a list of favorite products
- Current promotions
- Supplier product cross-matching
- On-site stocking and inventory management options
- Quick Quote
- E-commerce solutions
- New Lab Start-Up program



View the Virtual Stockroom at:
sigma-aldrich.com/virtualstockroom

Bookmark this page and allow us to become your personal stockroom!

Sigma-Aldrich Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA, S.A.
Tel: 54 11 4556 1472
Fax: 54 11 4552 1698

Australia

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: 612 9841 0555
Fax: 612 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: 43 1 605 81 10
Fax: 43 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA
Free Tel: 0800-14747
Free Fax: 0800-14745
Tel: 03 899 13 01
Fax: 03 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA
Tel: 55 11 3732-3100
Fax: 55 11 3733-5151

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 800-565-1400
Free Fax: 800-265-3858
Tel: 905-829-9500
Fax: 905-829-9292

China

SIGMA-ALDRICH CHINA INC
Tel: 86-21-6386 2766
Fax: 86-21-6386 3966

Czech Republic

SIGMA-ALDRICH s.r.o.
Tel: 246 003 200
Fax: 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: 43 56 59 10
Fax: 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND
Tel: 09-3509250
Fax: 09-350-92555

France

SIGMA-ALDRICH CHIMIE S.à.r.l
Tel appel gratuit: 0800 211 408
Fax appel gratuit: 0800 031 052

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800-51 55 000
Free Fax: 0800-649 00 00

Greece

SIGMA-ALDRICH (O.M.) LTD
Tel: 30 210 9948010
Fax: 30 210 9943831

Hungary

SIGMA-ALDRICH Kft
Tel: 06-1-235-9054
Fax: 06-1-269-6470
Ingyenes zöld telefon: 06-80-355-355
Ingyenes zöld fax: 06-80-344-344

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: 91-80-5112-7272
New Delhi: 91-11-2616 5477
Mumbai: 91-22-2570 2364
Hyderabad: 91-40-5531 5488
Fax
Bangalore: 91-80-5112-7473
New Delhi: 91-11-2616 5611
Mumbai: 91-22-2579 7589
Hyderabad: 91-40-5531 5466

Ireland

SIGMA-ALDRICH IRELAND LTD
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: 353 1 4041900
Fax: 353 1 4041910

Israel

SIGMA-ALDRICH ISRAEL LTD
Free Tel: 1-800-70-2222
Tel: 08-948-4100
Fax: 08-948-4200

Italy

SIGMA-ALDRICH S.r.l
Telefono: 02 33417310
Fax: 02 38010737
Numero Verde: 800-827018

Japan

SIGMA-ALDRICH JAPAN K.K.
Tokyo Tel: 03 5796 7300
Tokyo Fax: 03 5796 7315

Korea

SIGMA-ALDRICH KOREA
Tel: 031-329-9000
Fax: 031-329-9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: 603-56353321
Fax: 603-56354116

Mexico

SIGMA-ALDRICH QUIMICA, S.A. de C.V.
Free Tel: 01-800-007-5300
Free Fax: 01-800-712-9920

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Tel Gratis: 0800-0229088
Fax Gratis: 0800-0229089
Tel: 078-6205411
Fax: 078-6205421

New Zealand

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 0800 936 666
Free Fax: 0800 937 777

Norway

SIGMA-ALDRICH NORWAY AS
Tel: 23 17 60 60
Fax: 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: (+61) 829 01 00
Fax: (+61) 829 01 20

Portugal

SIGMA-ALDRICH QUIMICA, S.A.
Free Tel: 800 202180
Free Fax: 800 202178
Tel: 21 9242555
Fax: 21 9242610

Russia

SIGMA-ALDRICH RUSSIA
OOO SAF-LAB
Tel: 095 975-1917/3321
Fax: 095 975-4792

Singapore

SIGMA-ALDRICH PTE. LTD
Tel: 65-6779-1200
Fax: 65-6779-1822

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: 27 11 979 1188
Fax: 27 11 979 1119

Spain

SIGMA-ALDRICH QUIMICA S.A.
Free Tel: 900 101376
Free Fax: 900 102028
Tel: 91 661 99 77
Fax: 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: 020-350510
Fax: 020-352522
Outside Sweden Tel: +46 8 7424200
Outside Sweden Fax: +46 8 7424243

Switzerland

FLUKA CHEMIE GmbH
Swiss Free Call: 0800 80 00 80
Tel: +41 81 755 2828
Fax: +41 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD
Free Tel: 0800 717181
Free Fax: 0800 378785
Tel: 01747 833000
Fax: 01747 833313

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-free: 800-325-3010
Call Collect: 314-771-5750
Toll-Free Fax: 800-325-5052
Tel: 314-771-5765
Fax: 314-771-5757

Internet

sigma-aldrich.com

NEW Aldrich® Zipper-Lock AtmosBag™

Fast, easy access to contents with a reliable, gas tight seal

AtmosBag™ is a flexible, inflatable polyethylene chamber with built-in gloves that lets you work in a totally isolated and controlled environment. It is an inexpensive, portable alternative to the glove box, and is well suited to a variety of tasks and materials. It comes with the new **Technical Bulletin AL-211**, which contains instructional information and a listing of AtmosBag™ accessories

- Manipulation of air- and moisture-sensitive materials
- Sampling and weighing operations
- Desiccating or constant-humidity chamber
- Controlled habitat for botany and life sciences
- Field sampling and testing

One-Hand, Zipper-Lock AtmosBag™

Size	Opening (in.)	Inflation Vol. (L)	Cat. No.	Each
S	12	50	Z53,020-4	\$36 60
M	24	280	Z53,021-2	45 60
L	36	520	Z53,022-0	61 10



sigma-aldrich.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
ALDRICH • BOX 355 • MILWAUKEE • WISCONSIN • USA

ALDRICH
Advancing Science

AtmosBag and Aldrich are a trademark and a registered trademark, respectively, of Sigma-Aldrich Biotechnology, L.P.

Introducing science beyond the ordinary



SAFC

Inspiring Science™

SAFC is a leading provider of laboratory reagents and chemicals. We are committed to providing the highest quality products and services to our customers. Our products are used in a wide range of scientific and industrial applications. We are proud to be a part of the scientific community and to support the advancement of science.

Be inspired—by SAFC™.

SAFC, the SAFC logo, and Inspiring Science are trademarks of Sigma-Aldrich Biotechnology, L.P.

www.safcglobal.com

Member of the Sigma-Aldrich Group

ALDRICH CHEMICAL COMPANY, INC.
P.O. BOX 14508
ST. LOUIS, MISSOURI 63178, USA



RETURN SERVICE REQUESTED

PRSR STD
U.S. POSTAGE PAID
MADISON, WI
PERMIT NO. 2223

HVX

HONORING THE LIFE AND CONTRIBUTIONS OF H. C. BROWN

Aldrichimica ACTA

VOL. 38, NO. 2 • 2005



Bazille.

Herbert C. Brown and Aldrich:
Advancing Borane Chemistry for 32 Years

Organotrifluoroborates:
Expanding Organoboron Chemistry

Recent Advances in the Chemistry
of Lithium Aminoborohydrides

sigma-aldrich.com




ALDRICH[®]
Advancing Science

New Products from Aldrich R&D


REAGENTS FOR ASYMMETRIC SYNTHESIS

65,429-9 **1-(1S,2S)-2-(2-methylphenyl)pyrrolidin-1-yl-2-phenyl-1-oxo-1-phenylethane-1-thiolate** triflate, 0.5 M in toluene
5 mL \$31.50
25 mL 109.00



Excellent reagents for catalytic asymmetric transformations such as the Diels–Alder reaction¹ and the cyanosilylation of aldehydes.^{1,2} Chiral cationic oxazaborolidines, prepared from *o*-tolylloxazaborolidine and triflic acid, catalyze Diels–Alder reactions to provide high ee's.³ Recently, Corey and co-workers used these catalysts to synthesize several chiral natural products, otherwise synthesized in racemic form.⁴


65,430-2 **1-(1S,2S)-2-(2-methylphenyl)pyrrolidin-1-yl-2-phenyl-1-oxo-1-phenylethane-1-thiolate** triflate, 0.5 M in toluene
5 mL \$31.50
25 mL 109.00



(1) Corey, E. J. et al. *J. Am. Chem. Soc.* **2002**, *124*, 3808. (2) Ryu, D. H.; Corey, E. J. *J. Am. Chem. Soc.* **2004**, *126*, 8106. (3) Ryu, D. H. et al. *J. Am. Chem. Soc.* **2002**, *124*, 9992. (4) Hu, Q.-Y. et al. *J. Am. Chem. Soc.* **2004**, *126*, 13708.


REAGENTS FOR CROSS-COUPLING REACTIONS

64,848-5 **5-Chloroquinolin-8-yl trifluoromethanesulfonate, 5**
1 g \$25.00
5 g 78.00



Activated triflates of this type have served as alternatives to aryl halides in a variety of palladium-catalyzed coupling reactions

64,850-7 **5-Chloroquinolin-8-yl trifluoromethanesulfonate, 5**
1 g \$42.00
5 g 124.00



(1) Matthews, D. P. et al. *Tetrahedron Lett.* **1994**, *35*, 5177. (2) Tilley, J. W.; Zawoiski, S. *J. Org. Chem.* **1988**, *53*, 386. (3) Ellingboe, J. W. et al. *J. Med. Chem.* **1994**, *37*, 542.


LIGANDS FOR CATALYSIS

65,532-5 **tert-Butyldicyclohexylphosphine, 97%**
10 mL \$66.00
50 mL 230.00



This low-melting, air-sensitive, and pyrophoric solid is now available as an easier-to-use toluene solution in a Sure/Seal™ bottle

65,160-5 **tert-Butyldicyclohexylphosphine, 97%**
1 g \$49.00
5 g 172.00



This novel phosphine was used to prepare catalysts for the ring-opening metathesis polymerization of cyclooctene.

Jan, D. et al. *J. Organomet. Chem.* **2000**, *606*, 55.

REAGENT FOR C–X BOND FORMATION

65,531-7 **Dimethylphosphonothionic acid chloride**
1 g \$62.00
5 g 216.00




Reacts with imines under mild conditions and without demethylation to form α -aminophosphonothionates.

Tongcharoensirikul, P. et al. *J. Org. Chem.* **2004**, *69*, 2322.

BUILDING BLOCK

65,632-1 **3-Chloro-4-methoxyphenethylamine hydrochloride**
5 g \$80.00



A useful starting material for the synthesis of many natural products, such as 1,2,3,4-tetrahydroisoquinolines¹ and 2,3,4,5-tetrahydro-1H-3-benzazepines.²

(1) Charifson, P. S. et al. *J. Med. Chem.* **1988**, *31*, 1941. (2) Chumpradit, S. et al. *J. Med. Chem.* **1991**, *34*, 877.

For more information, please see pages 56, 58, and 59 for additional new products.

For competitive quotes on larger quantities, contact safcglobal.com.

SAFC

Aldrichimica ACTA

VOL. 38, NO. 2 • 2005

Aldrich Chemical Co., Inc.
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services:

Customer Inquiries	800-325-3010
Technical Service	800-231-8327
SAFC™	800-244-1173
Custom Synthesis	800-244-1173
Flavors & Fragrances	800-227-4563
International	414-438-3850
24-Hour Emergency	414-438-3850
Web site	sigma-aldrich.com
Email	aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201, USA

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
P.O. Box 355
Milwaukee, WI 53201-9358
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the inside back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.


All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich family. © 2005 Sigma-Aldrich Co.

"PLEASE BOTHER US."

Dr. Ian J. S. Fairlamb of the University of York, U.K., kindly suggested that we offer this Suzuki-Miyaura cross-coupling palladium catalyst in our catalog. Palladium-catalyzed cross-coupling reactions constitute some of the most important and applied transformations for the synthesis of natural products, pharmaceutical targets, and conjugated materials.¹ This palladium catalyst has shown a higher activity over Pd₂(dba)₃ in the Suzuki-Miyaura cross-coupling of organic halides with arylboronic acids.²

(1) *Handbook of Organopalladium Chemistry for Organic Synthesis*, Negishi, E., Ed., Wiley: Hoboken, NJ, 2002
(2) Fairlamb, I. J. S., Kapdi, A. R.; Lee, A. F. *Org Lett* **2004**, *6*, 4435



Ar = 3,5-dimethoxyphenyl

65,693-3 Bis(3,5,3',5'-dimethoxydibenzylideneacetone)-palladium(0) [Pd(dm-dba)₂] 1 g \$19.65
5 g 65.50

Naturally, we made this useful catalyst. It was no bother at all, just a pleasure to be able to help

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the inside back cover.

TABLE OF CONTENTS

Herbert C. Brown and Aldrich: Advancing Borane Chemistry for 32 Years	43
<i>Clinton F. Lane, Northern Arizona University</i>	
Organotrifluoroborates: Expanding Organoboron Chemistry	49
<i>Gary A. Molander* and Ruth Figueroa, University of Pennsylvania</i>	
Recent Advances in the Chemistry of Lithium Aminoborohydrides	61
<i>Lubov Pasumansky and Bakthan Singaram,* University of California, Santa Cruz, and Christian T. Goralski, CTG Consulting</i>	

ABOUT OUR COVER

The Ramparts at Aigues-Mortes (oil on canvas, 60 × 100 cm) was painted and signed by the artist Frédéric Bazille in 1867. Probably few tourists who visit the famous classical ruins at Nîmes and Arles in southern France are aware that nearby is the best-preserved medieval fortress in Europe, the subject of this painting. In 1246, King Louis IX began the tower seen just to the left of the center of this picture and had a five-mile channel dug through the estuary of the Rhone river to the Mediterranean for ships embarking on the Seventh and Eighth Crusades. His son, Philippe III, started the enormous adjoining battlements in 1272, but, within less than a century, silting closed the passage to the sea and the structure was abandoned.



Photograph © Board of Trustees, National Gallery of Art, Washington

Bazille was born not far away at Montpellier to a well-to-do family and was expected to become a doctor. While still in medical school, however, he convinced his father to subsidize his artistic education in Paris. He soon became fast friends with Alfred Sisley, Auguste Renoir, and Claude Monet, who suggested to Bazille and the others that they paint out-of-doors directly in front of their subjects, thus initiating the essential aspect of impressionist painting. When in the spring of 1867 Bazille went to paint *The Ramparts at Aigues-Mortes*, he wrote to his mother that he intended to paint "a very simple painting" of "the walls of the city reflected in a pond at sunset". The composition of the picture is indeed very simple, defined by the diagonal of the sandbar in the foreground and the almost equal areas of water and sky on either side of the long horizontal of the ramparts. The view is to the east, with the late afternoon sun striking the walls of the fortress. A tower near the right edge of the painting is reflected in the water, and the pink tinge of the clouds suggests the coming sunset.

This painting is a part of the Collection of Mr. and Mrs. Paul Mellon at the National Gallery of Art, Washington, DC.

Organoboranes for Asymmetric Synthesis

From Sigma-Aldrich

DIP-Chloride™ and DIP-Bromide™


Since its introduction 20 years ago, DIP-Chloride™ (*B*-chlorodiisopinocampheylborane, *lpc*₂BCl, DIP-Cl) has become one of the most effective reagents for the asymmetric reduction of prochiral ketones, especially aralkyl and α -hindered ketones.¹ The reagent has found popular use in the asymmetric synthesis of ligands,² β -amino alcohols,³ and alkaloids,⁴ and is also capable of effecting asymmetric aldol reactions.⁵ Furthermore, asymmetric allylations can be achieved when DIP-Cl is reacted with allylmagnesium bromide, allowing for further elaboration of the reduced product. We are now offering this valuable reagent also as an easier-to-handle solution in three different solvents: heptane, hexane, and α -pinene.

DIP-Bromide™ has been used for the asymmetric reduction of ketones, in enolboration reactions, and for the ring opening of meso epoxides.⁶ Sigma-Aldrich is pleased to offer these excellent reagents for asymmetric synthesis to our customers.

- (1) (a) Chandrasekharan, J. et al. *J. Org. Chem.* **1985**, *50*, 5446. (b) Srebnik, M.; Ramachandran, P. V. *Aldrichimica Acta* **1987**, *20*, 9
 (2) Drury, W. J., III et al. *Angew. Chem., Int. Ed.* **2004**, *43*, 70
 (3) Beardsley, D. A. et al. *Tetrahedron Lett.* **1994**, *35*, 1511
 (4) Felpin, F.-X.; Lebreton, J. *J. Org. Chem.* **2002**, *67*, 9192
 (5) Stocker, B. L. et al. *Eur. J. Org. Chem.* **2004**, 330
 (6) Dhar, R. K. *Aldrichimica Acta* **1994**, *27*, 43

(+)-DIP-Chloride™


31,701-2



5 g	\$33.50
25 g	108.00
100 g	300.50
1 kg	2145.00

(-)-DIP-Chloride™

31,702-0



5 g	\$35.00
25 g	106.00
100 g	289.50

(-)-DIP-Chloride™ solution, 65–75 wt. % in α -pinene

65,529-5

25 mL	\$90.00
100 mL	225.00

(-)-DIP-Chloride™ solution, 50–65 wt. % in heptane

64,841-8

25 mL	\$90.00
100 mL	225.00


(-)-DIP-Chloride™ solution, 50–65 wt. % in hexanes

64,842-6

25 mL	\$90.00
100 mL	225.00

(+)-DIP-Bromide™, 95%


41,427-1



5 g	\$74.10
25 g	270.00

(-)-DIP-Bromide™, 95%


41,099-3



5 g	\$177.00
25 g	590.00

L-Selectride® solution, 1.0 M in tetrahydrofuran


17,849-7



100 mL	\$37.40
800 mL	277.00
8 L	931.00
18 L	1865.00

LS-Selectride® solution, 1.0 M in tetrahydrofuran


22,592-4



100 mL	\$105.50
--------	----------

K-Selectride® solution, 1.0 M in tetrahydrofuran


22,076-0



100 mL	\$83.30
800 mL	562.00

KS-Selectride® solution, 1.0 M in tetrahydrofuran


22,077-9



100 mL	\$77.40
--------	---------

N-Selectride® solution, 1.0 M in tetrahydrofuran

21,340-3



100 mL	\$48.00
800 mL	257.50

The Selectride® Family

The Selectride® family of boranes is a group of powerful reagents for the asymmetric reduction of ketones on complex scaffolds. These reagents have been used in the synthesis of many natural products,⁷ in the stereoselective reduction of hemiacetals in various nucleoside syntheses,⁸ and are capable of effecting reductive cyclization⁹ and rearrangement reactions.¹⁰

- (7) Some recent examples: (a) Bahia, P. S.; Sneath, J. S. *J. Org. Chem.* **2004**, *69*, 3226. (b) Thede, K. et al. *Org. Lett.* **2004**, *6*, 4595. (c) Banwell, M. G. et al. *Org. Lett.* **2004**, *6*, 2737. (d) Lichtenthaler, F. W. et al. *Tetrahedron: Asymmetry* **2003**, *14*, 727. (e) Nelson, T. D. et al. *Tetrahedron Lett.* **2004**, *45*, 8917
 (8) (a) Hanessian, S.; Machaalani, R. *Tetrahedron Lett.* **2003**, *44*, 8321. (b) Navarre, J.-M. et al. *Tetrahedron Lett.* **2003**, *44*, 2199
 (9) Nagaoka, Y. et al. *Tetrahedron Lett.* **2002**, *43*, 4355
 (10) (a) Chen, W. et al. *J. Org. Chem.* **2003**, *68*, 1929. (b) Appendino, G. et al. *Eur. J. Org. Chem.* **2003**, 4422

For competitive quotes on larger quantities of all of these organoboranes, please contact us at 800-244-1173 (USA) or at safglobal.com.

sigma-aldrich.com

LEADERSHIP IN LIFE SCIENCE. HIGH TECHNOLOGY AND SERVICE
 ALDRICH • BOX 355 • MILWAUKEE • WISCONSIN • USA

ALDRICH
 Advancing Science

DIP-Chloride and DIP-Bromide are trademarks and Selectride a registered trademark of Sigma-Aldrich Biotechnology, L.P.

Sigma-Aldrich Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA, S.A.
Tel: 54 11 4556 1472
Fax: 54 11 4552 1698

Australia

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: 612 9841 0555
Fax: 612 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: 43 1 605 81 10
Fax: 43 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA.
Free Tel: 0800-14747
Free Fax: 0800-14745
Tel: 03 899 13 01
Fax: 03 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Tel: 55 11 3732-3100
Fax: 55 11 3733-5151

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 800-565-1400
Free Fax: 800-265-3858
Tel: 905-829-9500
Fax: 905-829-9292

China

SIGMA-ALDRICH CHINA INC.
Tel: 86-21-6386 2766
Fax: 86-21-6386 3966

Czech Republic

SIGMA-ALDRICH s.r.o.
Tel: +420 246 003 200
Fax: +420 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: 43 56 59 10
Fax: 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND
Tel: (09) 350 9250
Fax: (09) 350 9255

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Tel appel gratuit: 0800 211 408
Fax appel gratuit: 0800 031 052

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800-51 55 000
Free Fax: 0800-649 00 00

Greece

SIGMA-ALDRICH (O.M.) LTD
Tel: 30 210 9948010
Fax: 30 210 9943831

Hungary

SIGMA-ALDRICH Kft
Tel: 06-1-235-9054
Fax: 06-1-269-6470
Ingyenes zöld telefon: 06-80-355-355
Ingyenes zöld fax: 06-80-344-344

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: 91-80-5112-7272
New Delhi: 91-11-5165 4255
Mumbai: 91-22-2570 2364
Hyderabad: 91-40-5584 5488
Fax
Bangalore: 91-80-5112-7473
New Delhi: 91-11-5165 4266
Mumbai: 91-22-2579 7589
Hyderabad: 91-40-5584 5466

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: 353 1 4041900
Fax: 353 1 4041910

Israel

SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1-800-70-2222
Tel: 08-948-4100
Fax: 08-948-4200

Italy

SIGMA-ALDRICH S.r.l.
Telefono: 02 33417310
Fax: 02 38010737
Numero Verde: 800-827018

Japan

SIGMA-ALDRICH JAPAN K.K.
Tokyo Tel: 03 5796 7300
Tokyo Fax: 03 5796 7315

Korea

SIGMA-ALDRICH KOREA
Tel: 031-329-9000
Fax: 031-329-9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: 603-56353321
Fax: 603-56354116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01-800-007-5300
Free Fax: 01-800-712-9920

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Tel Gratis: 0800-0229088
Fax Gratis: 0800-0229089
Tel: 078-6205411
Fax: 078-6205421

New Zealand

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 0800 936 666
Free Fax: 0800 937 777

Norway

SIGMA-ALDRICH NORWAY AS
Tel: 23 17 60 60
Fax: 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: +4861 829 01 00
Fax: +4861 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202180
Free Fax: 800 202178
Tel: 21 9242555
Fax: 21 9242610

Russia

SIGMA-ALDRICH RUSSIA
OOO SAF-LAB
Tel: +7 (095) 975-1917/3321
Fax: +7 (095) 975-4792

Singapore

SIGMA-ALDRICH PTE. LTD.
Tel: 65-6779-1200
Fax: 65-6779-1822

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD.
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: 27 11 979 1188
Fax: 27 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA S.A.
Free Tel: 900 101376
Free Fax: 900 102028
Tel: 91 661 99 77
Fax: 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: 020-350510
Fax: 020-352522
Outside Sweden Tel: +46 8 7424200
Outside Sweden Fax: +46 8 7424243

Switzerland

FLUKA CHEMIE GmbH
Swiss Free Call: 0800 80 00 80
Tel: +41 81 755 2828
Fax: +41 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717181
Free Fax: 0800 378785
Tel: 01747 833000
Fax: 01747 833313
SAFC (UK): 01202 712305

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-free: 800-325-3010
Call Collect: 314-771-5750
Toll-Free Fax: 800-325-5052
Tel: 314-771-5765
Fax: 314-771-5757

Internet:

sigma-aldrich.com

Introducing science beyond the ordinary

Process Development and Production-Scale Raw Materials in

Returnable Containers

A few of the chemicals that SAFC™ offers in bulk quantities are:

- Allyl Alcohol
- Allyl Chloride
- Chlorosilanes
- Chlorosulfonic Acid
- Dimethyl Sulfate
- Epichlorohydrin
- Solvents, Anhydrous
- Solvents, HPLC-Grade

Packaged under an inert nitrogen atmosphere in a variety of cylinders and portable tanks in 1B- to 1,000-L sizes.

A safe, easy, and cost-effective way to purchase and transfer high-hazard, air-sensitive, and high-purity chemicals to your process vessels.

SAFC

For more information, please call Dan Zagrodnik
at 800-292-1206 (TJGA) or 314-438-3868.

Introducing science beyond the ordinary



SAFC

Inspiring Science™

SAFC™—where innovative scientists come to find the science they need to develop innovative products. Organic or biochemical, we offer integrated solutions throughout the product development plan, supporting you right through to commercialization. Applying years of chemistry know-how, process development and scale-up experience to bring your products to market faster.

No ordinary science here, but science beyond the ordinary.
Be inspired—by SAFC™.

www.safglobal.com

Member of The Sigma-Aldrich Group



ALDRICH®

P.O. Box 14508
St. Louis, MO 63178
USA

PRSR STD
U.S. POSTAGE PAID
MADISON, WI
PERMIT NO. 2223

Return Service Requested

THE NEGISHI COUPLING: AN UPDATE

Aldrichimica ACTA

VOL. 38, NO. 3 • 2005



Palladium-Catalyzed Alkenylation
by the Negishi Coupling

Enantiopure Sulfoxides and Sulfenamides

sigma-aldrich.com




ALDRICH[®]

Advancing Science

New Products from Aldrich R&D

TRIFLUOROBORATES


Potassium 4-(hydroxymethyl)phenyltrifluoroborate

659762		1 g	\$24.00
		5 g	80.00


Potassium 3-carboxyphenyltrifluoroborate

659789		1 g	\$24.00
		5 g	80.00


Potassium 3,4-(methylenedioxy)phenyltrifluoroborate

659754		1 g	\$24.00
		5 g	80.00


Potassium 3-fluorophenyltrifluoroborate

659770		1 g	\$24.00
		5 g	80.00


Potassium 3-hydroxyphenyltrifluoroborate

659746		1 g	\$24.00
		5 g	80.00


Potassium 4-*tert*-butylphenyltrifluoroborate, 95%

654728		1 g	\$14.50
		10 g	81.00

Potassium 2-methoxyphenyltrifluoroborate

654930		1 g	\$24.00
		5 g	80.00

Potassium vinyltrifluoroborate, 95%


655228		1 g	\$24.00
		5 g	80.00

Trifluoroborates are air-stable alternatives to boronic acids in the palladium-catalyzed Suzuki-Miyaura cross-coupling reaction.¹⁻³ They are more robust, easier to handle, and less prone to protodeboronation.¹ They display a remarkably uniform behavior.²

(1) Molander, G. A.; Biolatto, B. *J. Org. Chem.* **2003**, *68*, 4302. (2) Molander, G. A. et al. *J. Org. Chem.* **2003**, *68*, 5534. (3) Molander, G. A.; Figueroa, R. *Aldrichimica Acta* **2005**, *38*, 49

BUILDING BLOCKS


2,6-Dichloronicotinic acid, 97%

658219		1 g	\$33.00
		5 g	115.00

This pharmaceutical research^{1,2} key building block is now available in 97% purity. See **479586** for a technical grade (90%) of this material.

(1) Kelly, T. P. et al. *J. Labelled Compd. Radiopharm.* **2001**, *44*, 451. (2) Laeckmann, D. et al. *Bioorg. Med. Chem.* **2002**, *10*, 1793


2-Bromo-1-(*p*-toluenesulfonyl)pyrrole, 90%

650900		1 g	\$33.00
		5 g	115.00

This product is a stable and easily handled alternative for the labile bromopyrrole, and is an excellent substrate for the Suzuki coupling reaction. The tosyl group is readily cleaved under basic conditions to generate the 2-arylpiperole coupling product.

Knight, L. W. et al. *Synlett* **2003**, 1993

N-Boc-piperidine, 98%


655872		5 g	\$30.00
		25 g	105.00

A wide variety of substituted piperidines can be synthesized from this Boc-protected building block. Lithiation with BuLi leads to substitution at the 2 position with aryl, vinyl, and alkynyl halides.¹ Its use in the stereoselective synthesis of *D*-*threo*-methylphenidate (Ritalin[®]) and similar substrates has been reported.^{2,3}


(1) Dieter, R. K.; Li, S. J. *J. Org. Chem.* **1997**, *62*, 7726. (2) Axten, J. M. et al. *J. Am. Chem. Soc.* **1999**, *121*, 6511. (3) Beak, P.; Lee, W. K. *J. Org. Chem.* **1993**, *58*, 1109

BENZOFURAZANS


Benzofurazan, 97%

650137		1 g	\$17.50
		10 g	95.00

5-Methoxybenzofurazan, 97%

656925		1 g	\$26.50
		5 g	92.50

5-Chlorobenzofurazan

659118		1 g	\$24.00
		5 g	80.00

Benzofurazans have been utilized as effective fluorescence probe motifs¹ and as medicinal building blocks in the synthesis of antileukemic agents.^{2,3}

(1) Uchiyama, S. et al. *J. Chem. Soc., Perkin Trans. 2* **1999**, 2525. (2) Ghosh, P. B.; Whitehouse, M. W. *J. Med. Chem.* **1968**, *11*, 305. (3) Ghosh, P. B.; Whitehouse, M. W. *J. Med. Chem.* **1969**, *12*, 505

Please see pages 90-91 for additional new products.

For competitive quotes on larger quantities, contact www.safcglobal.com.

SAFC

Aldrichimica ACTA

VOL. 38, NO. 3 • 2005

Aldrich Chemical Co., Inc.
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-231-8327
SAFC™ 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 414-438-3850
24-Hour Emergency 414-438-3850
Web site sigma-aldrich.com
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201, USA

Subscriptions

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
P.O. Box 355
Milwaukee, WI 53201-9358
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the inside back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich family. © 2005 Sigma-Aldrich Co.

"PLEASE BOTHER US."



Joe Porwoll, President

Dr. Annegret Stark of the University of Jena (Friedrich-Schiller-Universität Jena), Germany, kindly suggested that we add high-purity ionic liquids (HPIL) to our existing extensive list of ionic liquids. Dr. Stark stated, "In our laboratories, HPILs serve as a benchmark in investigating impurity-dependent effects on reactions." The low water and halogen content of HPILs renders them more effective as solvents in transition-metal-catalyzed reactions.

(1) Seddon, K. R.; Stark, A. *Green Chem.* **2002**, *4*, 119. (2) Dyson, P. J.; Ellis, D. J.; Parker, D. G.; Welton, T. *Chem Commun.* **1999**, 25. (3) *Ionic Liquids in Synthesis*; Wasserscheid, P., Welton, T., Eds.; Wiley-VCH, 2003; p 26. (4) Stark, A.; Ajam, M.; Green, M.; Raubenheimer, H. G.; Ranwell, A.; Ondruschka, B. *Adv. Synth. Catal.*, submitted for publication, 2005.



R = Bu, Et; X = Cl, PF₆, BF₄

04129	1-Butyl-3-methylimidazolium chloride, puriss., dry, ≥99.0% (AT)	5 g	\$53.60
		25 g	214.30
55509	1-Butyl-3-methylimidazolium chloride, dry, ≥99.0% (AT)	5 g	\$57.40
		25 g	229.30
18122	1-Butyl-3-methylimidazolium hexafluorophosphate, for catalysis, ≥98.5% (T)	5 g	114.80
		50 g	820.00
39931	1-Butyl-3-methylimidazolium tetrafluoroborate, for catalysis, ≥98.5% (T)	5 g	105.00
		50 g	315.00
39736	1-Ethyl-3-methylimidazolium tetrafluoroborate, for catalysis, ≥98.5% (T)	5 mL	262.00
		50 mL	720.50

Naturally, we made these useful solvents. It was no bother at all, just a pleasure to be able to help.


Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the inside back cover.

TABLE OF CONTENTS

Palladium-Catalyzed Alkenylation by the Negishi Coupling	71
<i>Ei-ichi Negishi,* Qian Hu, Zhihong Huang, Mingxing Qian, and Guangwei Wang, Purdue University</i>	
Enantiopure Sulfoxides and Sulfinamides: Recent Developments in Their Stereoselective Synthesis and Applications to Asymmetric Synthesis	93
<i>Chris H. Senanayake,* Dhileepkumar Krishnamurthy, Zhi-Hui Lu, Zhengxu Han, and Isabelle Gallou, Boehringer Ingelheim Pharmaceuticals, Inc</i>	

ABOUT OUR COVER

Valdemosa, Majorca: Thistles and Herbage on a Hillside (oil on canvas, 55.8 × 71.1 cm) was painted by the artist John Singer Sargent in 1908. Sargent was born in 1856 in Florence, Italy, to expatriate American parents. His first artistic training was in Rome, but he later attended the Accademia delle Belle Arti in Florence, and studied drawing at the École des Beaux-Arts and painting in the studio of the portrait painter Charles Carolus-Duran in Paris. In 1877, he began to exhibit in the Salons, the French government sanctioned art exhibitions. Sargent copied works by Diego Velázquez on a trip to Spain in 1879 and by Frans Hals in Belgium and Holland in 1880, an experience that had a great impact on his artistic development.



Photograph by David of Trustees, National Gallery of Art, Washington


By the end of the nineteenth century, Sargent had become the most sought-after portrait painter of his time, but he always felt constrained by the limitations of portrait painting. By the early twentieth century, his success finally allowed him to free himself almost completely from the painting of formal portraits. He made annual trips to Spain, Italy, Austria, and Switzerland, and it was on a trip to the Balearic island of Majorca that he painted this small picture. He had no interest in what he called "enormous views and huge skies", and chose to concentrate on a small patch of vegetation growing in the earth of a hillside. This painting is not strictly a realistic image, but one which, with its strong formal contrasts, bright colors, and seemingly spontaneous execution, achieves an extraordinary intensity of expression. Confronted with such an immediate response to nature that seems almost spiritual in its intensity, it is easy to understand Sargent's disdain for the artificiality of the academic portrait painting that had dominated his career for most of his life.

This painting was acquired by the National Gallery of Art, Washington, DC, through the Avalon Fund and by Gift of Virginia Bailey Brown.

The Negishi Coupling Catalyzed by Palladium on Polymer Supports

Supported palladium catalysts are widely used in the Suzuki, Heck, and Sonogashira cross-coupling reactions. However, no examples of their use in the Negishi coupling have been reported in the literature. There are several advantages to using supported catalysts in organic synthesis. These

include reagent stability, suitability for automation, ease of workup, recyclability, and lower Pd contamination in the final product. Herein, we describe the application of four commercially available polymer-supported palladium reagents as catalysts in the Negishi coupling.



Typical Experimental Procedure

The palladium catalyst (0.01 mmol Pd) is charged into the reaction vessel. 3-Iodobenzotrifluoride (144 μL , 1 mmol) is then introduced, followed by addition of a THF solution of 4-methylphenylzinc iodide (0.5 M, 3 mL, 1.5 mmol). The resulting mixture is stirred at 50 °C for 18 h, cooled, and then filtered. The resin is washed with THF (2 \times 3 mL), the THF filtrates combined and evaporated. The evaporation residue is dissolved in a minimum

amount of THF and filtered through a silica gel pad to remove any residual zinc compounds. The pad is rinsed with ether, and the combined ether filtrates evaporated. The crude product thus obtained is purified by flash chromatography on silica gel (column size 1.5 \times 2.5 cm) using hexane as eluent. The purified product, 4-methyl-3'-trifluoromethylbiphenyl, is isolated as a colorless oil. (See the table below for yields.)

Cat. No.	Catalyst Name	Price	Product Yield	Byproduct Yield ^a
596930-1G 596930-5G	Dichlorobis(triphenylphosphine)palladium(II), polymer-bound	\$29.00 147.50	83%	4%
654167-5G 654167-25G	Diacetoxybis(triphenylphosphine)palladium(II), polymer-bound	280.00 800.00	84%	5%
646555-1G 646555-5G	Bis[(diphenylphosphanyl)methyl]aminepalladium(II) dichloride, polymer-bound	23.50 110.00	86%	4%
654159-1G 654159-5G	Bis[(diphenylphosphanyl)methyl]aminepalladium(II) diacetate, polymer-bound	187.50 700.00	81%	5%

^a 4,4'-Dimethylbiphenyl.

For more information, please contact your local Sigma-Aldrich office or visit us at sigma-aldrich.com.
For competitive quotes on larger quantities, please contact www.safcglobal.com.

sigma-aldrich.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
ALDRICH • BOX 355 • MILWAUKEE • WISCONSIN • USA

ALDRICH[®]
Advancing Science

Sigma-Aldrich Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA, S.A.
Tel: 54 11 4556 1472
Fax: 54 11 4552 1698

Australia

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: 612 9841 0555
Fax: 612 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: 43 1 605 81 10
Fax: 43 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA.
Free Tel: 0800-14747
Free Fax: 0800-14745
Tel: 03 899 13 01
Fax: 03 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Tel: 55 11 3732-3100
Fax: 55 11 3733-5151

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 800-565-1400
Free Fax: 800-265-3858
Tel: 905-829-9500
Fax: 905-829-9292

China

SIGMA-ALDRICH CHINA INC
Tel: 86-21-6386 2766
Fax: 86-21-6386 3966

Czech Republic

SIGMA-ALDRICH s.r.o.
Tel: +420 246 003 200
Fax: +420 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: 43 56 59 10
Fax: 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND
Tel: (09) 350 9250
Fax: (09) 350 9255

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Tel appel gratuit: 0800 211 408
Fax appel gratuit: 0800 031 052

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800-51 55 000
Free Fax: 0800-649 00 00

Greece

SIGMA-ALDRICH (O.M.) LTD
Tel: 30 210 9948010
Fax: 30 210 9943831

Hungary

SIGMA-ALDRICH Kft
Tel: 06-1-235-9054
Fax: 06-1-269-6470
Ingyenes zöld telefon: 06-80-355-355
Ingyenes zöld fax: 06-80-344-344

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: 91-80-5112-7272
New Delhi: 91-11-5165 4255
Mumbai: 91-22-2570 2364
Hyderabad: 91-40-5584 5488
Fax
Bangalore: 91-80-5112-7473
New Delhi: 91-11-5165 4266
Mumbai: 91-22-2579 7589
Hyderabad: 91-40-5584 5466

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: 353 1 4041900
Fax: 353 1 4041910

Israel

SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1-800-70-2222
Tel: 08-948-4100
Fax: 08-948-4200

Italy

SIGMA-ALDRICH S.r.l.
Telefono: 02 33417310
Fax: 02 38010737
Numero Verde: 800-827018

Japan

SIGMA-ALDRICH JAPAN K.K.
Tokyo Tel: 03 5796 7300
Tokyo Fax: 03 5796 7315

Korea

SIGMA-ALDRICH KOREA
Tel: 031-329-9000
Fax: 031-329-9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: 603-56353321
Fax: 603-56354116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01-800-007-5300
Free Fax: 01-800-712-9920

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Tel Gratis: 0800-0229088
Fax Gratis: 0800-0229089
Tel: 078-6205411
Fax: 078-6205421

New Zealand

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: 61 2 9841 0500
Fax: 61 2 9841 0500

Norway

SIGMA-ALDRICH NORWAY AS
Tel: 23 17 60 60
Fax: 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: +4861 829 01 00
Fax: +4861 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202180
Free Fax: 800 202178
Tel: 21 9242555
Fax: 21 9242610

Russia

SIGMA-ALDRICH RUS, LLC
Tel: +7 (095) 975-1917/3321
Fax: +7 (095) 975-4792

Singapore

SIGMA-ALDRICH PTE. LTD
Tel: 65-67791200
Fax: 65-67791822

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD.
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: 27 11 979 1188
Fax: 27 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA S.A.
Free Tel: 900 101376
Free Fax: 900 102028
Tel: 91 661 99 77
Fax: 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: 020-350510
Fax: 020-352522
Outside Sweden Tel: +46 8 7424200
Outside Sweden Fax: +46 8 7424243

Switzerland

SIGMA-ALDRICH CHEMIE LLC
Swiss Free Call: 0800 80 00 80
Tel: +41 81 755 2828
Fax: +41 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717181
Free Fax: 0800 378785
Tel: 01747 833000
Fax: 01747 833313
SAFC (UK): 01202 712305

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-free: 800-325-3010
Call Collect: 314-771-5750
Toll-Free Fax: 800-325-5052
Tel: 314-771-5765
Fax: 314-771-5757

Internet

sigma-aldrich.com

Introducing science beyond the ordinary

Process Development and Production-Scale Raw Materials in

Returnable Containers

A few of the chemicals that SAFC™ offers in large quantities are:

- Allyl Alcohol
- Allyl Chloride
- Chlorosulfonic Acid
- Dimethyl Sulfate
- Epichlorohydrin
- Solvents, Anhydrous
- Solvents, HPLC-Grade

Packaged under an inert nitrogen atmosphere in a variety of cylinders and portable tanks in 18- to 1,000-L sizes.

A safe, easy, and cost-effective way to purchase and transfer high-hazard, air-sensitive, and high-purity chemicals to your process vessels.

SAFC

For more information, please call Dan Zagrodnik at 800-213-1206 (USA) or 314-888-3888.

www.safcglobal.com

Member of the Sigma-Aldrich Group

Inspired by a culture of creativity



SAFC

SAFC is a leading provider of specialty chemicals and reagents for research and development.

SAFC is a leading provider of specialty chemicals and reagents for research and development. We are committed to providing the highest quality products and services to our customers. Our products are used in a wide range of applications, including pharmaceuticals, biotechnology, and materials science.

SAFC is a leading provider of specialty chemicals and reagents for research and development. We are committed to providing the highest quality products and services to our customers. Our products are used in a wide range of applications, including pharmaceuticals, biotechnology, and materials science.

SAFC™—your chance to shine.

www.safglobal.com

Member of the Sigma-Aldrich Group



ALDRICH®

P.O. Box 14508
St. Louis, MO 63178
USA

Return Service Requested

PSRT STD
U.S. POSTAGE PAID
MADISON, WI
PERMIT NO. 2223

ALDRICH CONGRATULATES THE 2006 ACS AWARD WINNERS

Aldrichimica ACTA

VOL. 39, NO. 1 • 2006



**Organic Synthesis with Light-Fluorous Reagents,
Reactants, Catalysts, and Scavengers**

**Synthetic Applications of Buchwald's Phosphines in
Palladium-Catalyzed Aromatic-Bond-Forming Reactions**

sigma-aldrich.com



ALDRICH[®]

Advancing Science


New Products from Aldrich R&D

Phenylpropionic Acid Derivatives

3-(4-Chlorophenyl)propionic acid, 97%

656151

[2019-34-3]
C₉H₉ClO₂




1 g \$20.80
5 g 52.00

3-(2,4-Dichlorophenyl)propionic acid, 97%

656186

[55144-92-8]
C₉H₇Cl₂O₂




1 g \$18.70
5 g 62.40

3-(3,4-Methylenedioxyphenyl)propionic acid, 97%

657565

[2815-95-4]
C₉H₉O₃




1 g \$24.50
5 g 85.00

3-(3-Methylphenyl)propionic acid, 97%

656178

[3751-48-2]
C₁₀H₁₂O₂




1 g \$22.90
5 g 62.40

3-(4-Methylthiophenyl)propionic acid, 97%

656135

C₁₀H₁₂O₂S



1 g \$20.80
5 g 52.00

Versatile building blocks employed in the synthesis of several biologically active compounds, including γ -secretase inhibitors^{1,2} and antagonists of human vanilloid receptor 1,^{3,4} adenosine A₁ receptor,⁵ and dopamine D₄ receptor.⁶


(1) Owens, A. P. et al. *Bioorg. Med. Chem. Lett.* **2003**, *13*, 4143. (2) Churcher, I. *Bioorg. Med. Chem. Lett.* **2003**, *13*, 179. (3) Jetter, M. C. et al. *Bioorg. Med. Chem. Lett.* **2004**, *14*, 3053. (4) McDonnell, M. E. et al. *Bioorg. Med. Chem. Lett.* **2004**, *14*, 531. (5) Van Calenbergh, S. et al. *J. Med. Chem.* **2002**, *45*, 1845. (6) Unangst, P. C. et al. *J. Med. Chem.* **1997**, *40*, 4026

Organometallic Reagents

(1,5-Cyclooctadiene)(1,3,5-cyclooctatriene)ruthenium, Ru(cod)(cot)

654418

[127382-91-6]
C₈H₁₂Ru



1 g \$165.00


Highly selective catalyst for amine alkylations, [2+2] cycloadditions, and enyne generation;¹ and key precursor for the preparation of stabilized Ru nanoparticles.²

(1) (a) Mitsudo, T. et al. *J. Am. Chem. Soc.* **1999**, *121*, 1839. (b) Watanabe, Y. et al. *J. Org. Chem.* **1996**, *61*, 4214. (2) (a) Pelzer, K. et al. *Chem. Mater.* **2004**, *16*, 4937. (b) Pan, C. et al. *J. Am. Chem. Soc.* **2001**, *123*, 7584

Dicyclohexyl(2-methylphenyl)phosphine, 95%

651885

[173593-25-4]
C₁₄H₂₂P




1 g \$26.50
10 g 145.00

Bulky ligand in catalyst systems utilized in coupling reactions.

Isopropylmagnesium chloride–lithium chloride complex solution, 1.0 M in tetrahydrofuran

656984

[807329-97-1]
C₃H₇Cl₂LiMg



10 mL \$18.00
100 mL 60.00

Performs halogen–magnesium exchange better than its counterparts, isopropylmagnesium chloride and diisopropylmagnesium


Ren, H. et al. *Org. Lett.* **2004**, *6*, 4215

1,2-Amino Alcohols

N-Boc-(*R*)-(+)-2-amino-1-butanol, 96%

660108

C₉H₁₉NO₃




1 g \$24.00
5 g 80.00

N-Boc-(*S*)-(-)-2-amino-1-butanol, 96%

660116

C₉H₁₉NO₃




1 g \$24.00
5 g 80.00

N-Boc-DL-2-amino-1-butanol, 97%

657786

[138373-86-1]
C₉H₁₉NO




1 g \$18.70
5 g 62.40

N-Boc-1-amino-1-cyclopentanemethanol, 97%

657689

C₁₁H₂₁NO₃




1 g \$65.50
5 g 218.50

N-Boc-DL-2-amino-1-hexanol, 90%

657794

C₁₁H₂₃NO₃




1 g \$40.00
10 g 275.00

N-Boc-(2*S*,3*S*)-(-)-2-amino-3-methyl-1-pentanol, 96%

660205

C₁₁H₂₃NO₃




2 g \$46.00
10 g 152.00

N-Boc-DL-2-amino-1-pentanol, 97%

657808

[179684-02-7]
C₁₀H₂₁NO₃




1 g \$31.20
5 g 104.00

N-Boc-2-amino-2-methyl-1-propanol, 97%

657778

[102520-97-8]
C₉H₁₉NO₃




1 g \$23.40
5 g 78.00

N-Boc-DL-2-amino-1-propanol, 95%

657816

[147252-84-4]
C₈H₁₇NO₃




1 g \$26.50
5 g 88.40

N-Boc-(1*R*,2*R*)-(-)-2-amino-1-(4-nitrophenyl)-1,3-propanediol

660213

C₁₄H₂₀N₂O₆



1 g \$24.00
5 g 80.00

Versatile building blocks for the synthesis of aldehydes,¹ carboxylic acids,² aziridines,³ azides,⁴ and oxathiazolidines.⁵

(1) Tietze, L. F.; Burkhardt, O. *Synthesis* **1994**, 1331. (2) (a) Kokotos, G. et al. *J. Med. Chem.* **2004**, *47*, 3615. (b) Sasaki, N. A. *Tetrahedron Lett.* **1987**, *28*, 6069. (3) (a) Wessig, P.; Schwarz, J. *Synlett* **1997**, 893. (b) Braga, A. L. et al. *Org. Lett.* **2003**, *5*, 2635. (4) Benalil, A. et al. *Tetrahedron* **1991**, *47*, 8177. (5) Posakony, J. J. et al. *J. Org. Chem.* **2002**, *67*, 5164

Aldrichimica ACTA

VOL. 39, NO. 1 • 2006

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-231-8327
SAFC™ 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 414-438-3850
24-Hour Emergency 414-438-3850
Web Site sigma-aldrich.com
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201, USA

Subscriptions

To request your FREE subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)

Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 355
Milwaukee, WI 53201-9358

Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the inside back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich family. © 2006 Sigma-Aldrich Co.


"PLEASE BOTHER US."

Joe Porwoll


Joe Porwoll, President
Aldrich Chemical Co., Inc

Professor Larry Overman of the University of California, Irvine, kindly suggested that we offer these cobalt(II) oxazoline palladacycle enantiomers. These acetate-bridged dimers (COP-OAc) effectively promote the asymmetric rearrangement of allylic trichloroacetimidates to allylic trichloroacetamides.¹ This grants ready access to valuable allylic amines of high enantiomeric purity, following the removal of the trichloroacetyl group under basic, acidic, or reductive conditions

(1) K.irsch, S. F.; Overman, L. E.; Watson, M. P. *J. Org. Chem.* **2004**, *69*, 8101



661716



661708

661716	(S)-(+)-COP-OAc Catalyst, 95%	250 mg	\$75.00
		1 g	215.00
661708	(R)-(-)-COP-OAc Catalyst, 95%	250 mg	\$35.00
		1 g	120.00

Naturally, we made these useful catalysts. It was no bother at all, just a pleasure to be able to help

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the inside back cover.

TABLE OF CONTENTS

Organic Synthesis with Light-Fluorous Reagents, Reactants, Catalysts, and Scavengers	3
Dennis P. Curran, University of Pittsburgh	
Synthetic Applications of Buchwald's Phosphines in Palladium-Catalyzed Aromatic-Bond-Forming Reactions	17
Christelle C. Mauger* and Gérard A. Mignani, Rhodia Recherches et Technologies	

ABOUT OUR COVER

Family Group (oil on canvas, 182.8 × 213.3 cm) was painted by the American painter William Glackens in 1910–1911. Glackens was one of a number of artists whose works reflected the rapid changes that were occurring in America in the early years of the twentieth century. America was evolving from a predominantly agricultural society into an industrial power, people were moving from the country to the cities, and massive immigration was rapidly increasing the population. These artists strove to document the realism of everyday life, painting scenes in fashionable cafés and restaurants, immigrant life on New York City's Lower East Side, theatergoers emerging from popular shows on Broadway, and views in Central Park and the city's streets. Rejecting any self-conscious aestheticism or romanticism, eight of them joined together in 1908 to stage their own independent exhibition, repudiating the prevailing academicism of the time. The culmination of this movement occurred in 1913 with *The Armory Show*, a huge exhibition of modern American and European art in New York, in which Glackens's *Family Group* was first shown.

The painting records a visit to the artist's Fifth Avenue apartment by Grace Morgan, a family friend who had recently returned to New York from France. Mrs. Edith Glackens leans on her sister Irene's chair, while Ira Glackens, the artist's son, stands between his mother and a table on which their visitor rests her elbow. The light streaming in from the window in the background illuminates the diverse colors and patterns in the room, and the casual arrangement of the figures and their closeness to the viewer underscore the naturalness of the artist's approach. Even the fact that the bottom edge of the picture cuts off the diagonal of the carpet implies that the carpet extends into the space occupied by the viewer, further enhancing the immediacy and realism of Glackens's representation of a singular, but otherwise ordinary event in the life of his family

This painting is a gift of Mr. and Mrs. Ira Glackens to the National Gallery of Art, Washington, DC.




Family Group, by William Glackens, National Gallery of Art, Washington, DC

Fluorinating Reagents from Sigma-Aldrich

The importance of selectively fluorinating compounds in medicinal chemistry, biology, and organic synthesis is well appreciated and provides a major impetus to the discovery of new and mild fluorinating agents that can operate safely and efficiently. Elemental fluorine and many electrophilic fluorinating agents have been used in synthesis; however, most of these fluorinating agents are highly aggressive, unstable, and require special equipment and care for safe handling. Sigma-Aldrich is pleased to offer the following alternatives, which lack these drawbacks

NEW 4-Iodotoluene Difluoride (Tol-IF₂)


4-Iodotoluene difluoride (Tol-IF₂) is easy to handle and less toxic than many fluorinating agents. Selective monofluorination of β-keto esters, β-keto amides, and β-diketones takes place under mild conditions without the use of HF-amine complexes.¹ A new methodology for the synthesis of fluorinated cyclic ethers was recently reported, which utilized Tol-IF₂ to achieve a fluorinative ring-expansion of four-, five-, and six-membered rings.² When one equivalent of Tol-IF₂ is reacted with phenylsulfanylated esters, the α-fluoro sulfide results through a fluoro-Pummerer reaction.³ When phenylsulfanylated lactams were treated with two equivalents of Tol-IF₂, the lactams were fluorinated in the α and β positions, resulting in diastereomeric difluorides.⁴



Diethylaminosulfur Trifluoride (DAST)

DAST has been regularly employed in selective fluorinations of alcohols, alkenols, carbohydrates, ketones, sulfides, epoxides, thioethers, and cyanohydrins. In addition, some novel organic cyclizations are possible when DAST is employed as a reagent.⁵ 1,2,2-Trifluorostyrenes can be synthesized by fluorination of the parent α-(trifluoromethyl)phenylethanol with DAST, followed by dehydrohalogenation with lithium bis(trimethylsilyl)amide (LHMDS). This method leads to the trifluorostyrene without requiring a palladium-catalyzed coupling.⁶

4-Iodotoluene difluoride



p-Tolyliododifluoride

[371-11-9] C₇H₇I₂F₂ FW 256.03

mp 102–106 °C

ρ₄: R: 34 S: 26-36/37/39-45

651117-1G 1 g \$32.80

651117-5G 5 g 112.50

(Diethylamino)sulfur trifluoride



DAST

[38078-09-0] (C₂H₅)₂NSF₃ FW 161.19

Feiser **6**, 183; **10**, 142; **12**, 183; **13**, 110; **16**, 128

bp 30–32 °C / 3 mm Hg

density 1.22 g/mL at 25 °C

ρ₄: R: 11-14-20/21/22-34 S: 16-26-36/37/39-45

235253-1G 1 g \$20.30


235253-5G 5 g 62.90

235253-25G 25 g 194.50

235253-125G 125 g 606.00

235253-250G 250 g 936.00

Selectfluor™ fluorinating reagent



F-TEDA

[740681-55-6] C₇H₁₄B₂ClF₉N₂ FW 354.26

>95% in F+ active


mp 260 °C

ρ₄: R: 22-36/37/38-41 S: 26-36/37/39

439479-5G 5 g \$20.50

439479-25G 25 g 70.80

Ishikawa's Reagent



N,N-Diethyl-1,1,2,3,3,3-hexafluoropropylamine

[309-88-6] CF₃CF₂CHFN(C₂H₅)₂ FW 223.16

bp 56–57 °C n_D²⁰ 1.3460


density 1.230 g/mL at 25 °C

ρ₄: R: 10-34 S: 16-26-36/37/39-45

564990-25G 25 g \$63.70

564990-100G 100 g 217.00

The one-pot synthesis of α,α-difluoroamides via direct fluorination was recently reported using DAST as the fluorinating reagent. Decreasing the molar ratio of DAST to substrate resulted in the formation of the corresponding α-keto amide.⁷




Selectfluor™ Fluorinating Reagent (F-TEDA)

Selectfluor™ fluorinating reagent ([1-chloromethyl-4-fluoro-1,4-diazoniabicyclo[2.2.2]octane bis(tetrafluoroborate), or F-TEDA) is a user friendly, mild, air- and moisture-stable, non-volatile reagent for electrophilic fluorination. Selectfluor™ is capable of introducing fluorine into organic substrates in one step, with a remarkably broad scope of reactivity, often with excellent regioselectivity.⁸ For example, allylic fluorides can be prepared using Selectfluor™ via a sequential cross-metathesis–electrophilic fluorodesilylation route. This route avoids the formation of byproducts that result from allylic transposition, which is observed when nucleophilic displacement or ring-opening reaction with DAST is attempted.⁹



Ishikawa's Reagent

Ishikawa's reagent [(N,N-diethyl-1,1,2,3,3,3-hexafluoropropylamine) has also demonstrated its utility in fluorination reactions. In a recent example, Ishikawa's reagent has been used in a clean and reproducible fluorination protocol on an intermediate in the total synthesis of 26-fluoropothione B.¹⁰ Ishikawa's reagent also displays the ability to insert a fluoro(trifluoromethyl)methylene moiety into unsaturated alcohols.¹¹



References: (1) Yoshida, M. et al. *Arkivoc* [Online] **2003**(vi), 36. (2) Inagaki, T. et al. *Tetrahedron Lett* **2003**, *44*, 4117. (3) Motherwell, W. B. et al. *J. Chem. Soc., Perkin Trans. 1* **2002**, 2809. (4) Greaney, M. F. et al. *Tetrahedron Lett* **2001**, *42*, 8523. (5) For a review, see Singh, R. P.; Shreeve, J. M. *Synthesis* **2002**, 2561. (6) Anilkumar, R.; Burton, D. J. *Tetrahedron Lett* **2003**, *44*, 6661. (7) Singh, R. P.; Shreeve, J. M. *J. Org. Chem.* **2003**, *68*, 6063. (8) For a review, see Singh, R. P.; Shreeve, J. M. *Acc. Chem. Res.* **2004**, *37*, 31. (9) Thibaudeau, S.; Gouverneur, V. *Org. Lett.* **2003**, *5*, 4891. (10) Koch, G. et al. *Synlett* **2004**, 693. (11) Ogu, K.-i. et al. *Tetrahedron Lett* **1998**, *39*, 305

Sigma-Aldrich Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA, S.A.
Tel: 54 11 4556 1472
Fax: 54 11 4552 1698

Australia

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: 612 9841 0555
Fax: 612 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: 43 1 605 81 10
Fax: 43 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA.
Free Tel: 0800-14747
Free Fax: 0800-14745
Tel: 03 899 13 01
Fax: 03 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Tel: 55 11 3732-3100
Fax: 55 11 3733-5151

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 800-565-1400
Free Fax: 800-265-3858
Tel: 905-829-9500
Fax: 905-829-9292

China

SIGMA-ALDRICH CHINA INC.
Tel: 86-21-6386 2766
Fax: 86-21-6386 3966

Czech Republic

SIGMA-ALDRICH s.r.o.
Tel: +420 246 003 200
Fax: +420 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: 43 56 59 10
Fax: 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND
Tel: (09) 350 9250
Fax: (09) 350 92555

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Tel appel gratuit: 0800 211 408
Fax appel gratuit: 0800 031 052

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800-51 55 000
Free Fax: 0800-649 00 00

Greece

SIGMA-ALDRICH (O.M.) LTD
Tel: 30 210 9948010
Fax: 30 210 9943831

Hungary

SIGMA-ALDRICH Kft
Tel: 06-1-235-9054
Fax: 06-1-269-6470
Ingyenes zöld telefon: 06-80-355-355
Ingyenes zöld fax: 06-80-344-344

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: 91-80-5112-7272
New Delhi: 91-11-5165 4255
Mumbai: 91-22-2570 2364
Hyderabad: 91-40-5584 5488
Fax
Bangalore: 91-80-5112-7473
New Delhi: 91-11-5165 4266
Mumbai: 91-22-2579 7589
Hyderabad: 91-40-5584 5466

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: 353 1 4041900
Fax: 353 1 4041910

Israel

SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1-800-70-2222
Tel: 08-948-4100
Fax: 08-948-4200

Italy

SIGMA-ALDRICH S.r.l.
Telefono: 02 33417310
Fax: 02 38010737
Numero Verde: 800-827018

Japan

SIGMA-ALDRICH JAPAN K.K.
Tokyo Tel: 03 5796 7300
Tokyo Fax: 03 5796 7315

Korea

SIGMA-ALDRICH KOREA
Tel: 031-329-9000
Fax: 031-329-9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: 603-56353321
Fax: 603-56354116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01-800-007-5300
Free Fax: 01-800-712-9920

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Tel Gratis: 0800-0229088
Fax Gratis: 0800-0229089
Tel: 078-6205411
Fax: 078-6205421

New Zealand

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: 61 2 9841 0500
Fax: 61 2 9841 0500

Norway

SIGMA-ALDRICH NORWAY AS
Tel: 23 17 60 60
Fax: 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: 061 829 01 00
Fax: 061 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202180
Free Fax: 800 202178
Tel: 21 9242555
Fax: 21 9242610

Russia

SIGMA-ALDRICH RU\$, LLC
Tel: +7 (095) 621-5828/6037
Fax: +7 (095) 975-4792

Singapore

SIGMA-ALDRICH PTE. LTD
Tel: 65-67791200
Fax: 65-67791822

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD.
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: 27 11 979 1188
Fax: 27 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA S.A.
Free Tel: 900 101376
Free Fax: 900 102028
Tel: 91 661 99 77
Fax: 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: 020-350510
Fax: 020-352522
Outside Sweden Tel: +46 8 7424200
Outside Sweden Fax: +46 8 7424243

Switzerland

SIGMA-ALDRICH CHEMIE LLC
Swiss Free Call: 0800 80 00 80
Tel: +41 81 755 2828
Fax: +41 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717181
Free Fax: 0800 378785
Tel: 01747 833000
Fax: 01747 833313
SAFC (UK): 01202 712305

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-free: 800-325-3010
Call Collect: 314-771-5750
Toll-Free Fax: 800-325-5052
Tel: 314-771-5765
Fax: 314-771-5757

Internet

sigma-aldrich.com

Introducing science beyond the ordinary

Process Development and Production-Scale Raw Materials in Returnable Containers

A few of the chemicals that SAFC™ offers in large quantities are:

- Alkyl Alcohol
- Alkyl Chloride
- Chloroalkanes
- Chlorosulfonic Acid
- Dimethyl Sulfate
- Epichlorohydrin
- Solvents, Anhydrous
- Solvents, HPLC-Grade

Packaged under an inert nitrogen atmosphere in a variety of cylinders and portable tanks in 18- to 1,000-L sizes.

A safe, easy, and cost-effective way to purchase and transfer high-hazard, air-sensitive, and high-purity chemicals to your process vessels.

SAFC

For more information, please call Dan Zagrodzki
at 800-293-1288 (USA) or 414-438-3800.

Inspired by a culture of creativity



SAFC

Inspiring Science™

SAFC is a leading provider of specialty chemicals. From the earliest stages of research and development through commercialization, we bring our unique brand of brilliance to the complex multi-step processes that bring some of the world's brightest organic, biochemical, and specialty ideas to life.

With rapidly advancing technologies, extensive manufacturing capabilities, and dedicated professionals, SAFC™ takes pride in getting groundbreaking products to market faster. Through a culture of innovation and excellence, we aim to inspire fellow scientists to create a brighter future for all.

SAFC™—your chance to shine.

www.safcglobal.com

Member of the Sigma-Aldrich Group



ALDRICH®

P.O. Box 14508
St. Louis, MO 63178
USA

Return Service Requested

PRSR STD
U.S. POSTAGE PAID
MADISON, WI
PERMIT NO. 2223

ORGANIC CHEMISTRY ENABLING MOLECULAR ELECTRONICS

Aldrichimica ACTA

VOL. 39, NO. 2 • 2006



**Chiral, Poly(Rare-Earth Metal) Complexes
in Asymmetric Catalysis**

**Organic Synthesis and Device Testing
for Molecular Electronics**

sigma-aldrich.com



ALDRICH®


Advancing Science

New Products from Aldrich R&D

Aldrich Is Pleased to Offer the Latest Catalysis Technology

Du Bois C-H Amination Catalyst


The Du Bois group at Stanford University has utilized $Rh_2(esp)_2$ to facilitate both inter- and intramolecular N-insertion in a range of benzylic, secondary, and tertiary C-H bonds, giving the corresponding aminated products in high yields.^{1,2}



(1) (a) Espino, C. G. et al. *J. Am. Chem. Soc.* **2004**, 126, 15378. (b) Fiori, K. W. et al. *Angew. Chem., Int. Ed.* **2004**, 43, 4349. (2) (a) Dauban, P.; Dodd, R. H. *Synlett* **2003**, 1571. (b) Muller, P.; Frut, C. *Chem. Rev.* **2003**, 103, 2905. (c) Diaz-Requejo, M. M. et al. *J. Am. Chem. Soc.* **2003**, 125, 12078


Bis[rhodium($\alpha,\alpha,\alpha',\alpha'$ -tetramethyl-1,3-benzenedipropionic acid)], 96%

$Rh_2(esp)_2$ 100 mg \$35.00
662623 500 mg 125.00
 [819050-89-0]
 $C_{32}H_{40}O_8Rh_2$
 FW: 758.47




Bis(tert-butylcarbonyloxy)iodobenzene, 97%

662283 5 g \$15.00
 [57357-20-7] 25 g 40.00
 $C_{16}H_{23}IO_4$ 100 g 110.00
 FW: 406.26




2,2,2-Trichloroethoxysulfonamide, 97%

663727 1 g \$25.00
 [69226-51-3] 10 g 154.00
 $C_2H_4Cl_3NO_2S$
 FW: 228.48



Hayashi Asymmetric Conjugate-Addition Catalyst and Precursors


The dimeric catalysts and catalyst precursors developed by Hayashi demonstrate impressive levels of enantiocontrol in the conjugate-addition reactions of both acyclic and cyclic enones of varying electronic character.¹⁻⁵



(1) Hayashi, T. et al. *J. Am. Chem. Soc.* **2002**, 124, 5052. (2) Hayashi, T.; Yamasaki, K. *Chem. Rev.* **2003**, 103, 2829. (3) Takaya, Y. et al. *J. Am. Chem. Soc.* **1998**, 120, 5579. (4) Takaya, Y. et al. *Tetrahedron: Asymmetry* **1999**, 10, 4047. (5) Takaya, Y. et al. *Tetrahedron Lett.* **1999**, 40, 6957


Hydroxy(cyclooctadiene)rhodium(I) dimer, 95%

661023 250 mg \$28.00
 [73468-85-6] 1 g 80.00
 $C_{16}H_{26}O_2Rh_2$
 FW: 456.19




Methoxy(cyclooctadiene)rhodium(I) dimer

661058 250 mg \$28.00
 [12148-72-0] 1 g 80.00
 $C_{18}H_{30}O_2Rh_2$
 FW: 484.24




Hydroxy[(S)-BINAP]rhodium(I) dimer, 90%

661007 1 g \$75.00
 [434314-10-0]
 $C_{88}H_{66}O_2P_4Rh_2$
 FW: 1485.17



Overman Asymmetric Allylic Rearrangement Catalysts


Overman and others have successfully utilized Co-based COP catalysts in the preparation of chiral amides and esters via cyclic rearrangement of allylic trichloroacetimidates. These rearrangements occur in high yield with excellent enantiocontrol.^{1,2}



(1) (a) Kirsch, S. F. et al. *J. Org. Chem.* **2004**, 69, 8101. (b) Anderson, C. E.; Overman, L. E. *J. Am. Chem. Soc.* **2003**, 125, 12412. (c) Kirsch, S. F.; Overman, L. E. *J. Am. Chem. Soc.* **2005**, 127, 2866. (2) Kwon, T. W. et al. *J. Org. Chem.* **1992**, 57, 6169.


(R)-(-)-COP-OAc Catalyst, 95%

661708 250 mg \$35.00
 [849592-74-1] 1 g 120.00
 $C_{82}H_{72}O_6N_2Co_2Pd_2$
 FW: 1512.17




(S)-(+)-COP-OAc Catalyst, 95%

661716 250 mg \$35.00
 [222400-03-5] 1 g 120.00
 $C_{82}H_{72}O_6N_2Co_2Pd_2$
 FW: 1512.17




(R)-(-)-COP-Cl Catalyst

661791 100 mg \$45.00
 [612065-00-6] 500 mg 245.00
 $C_{78}H_{66}Cl_2Co_2N_2O_2Pd_2$
 FW: 1464.98



(S)-(+)-COP-Cl Catalyst

646636 250 mg \$45.90
 [612065-01-7] 1 g 128.00
 $C_{78}H_{66}Cl_2Co_2N_2O_2Pd_2$
 FW: 1464.98



Please see pages 44-45 for additional new products. For competitive quotes on larger quantities, contact www.safcglobal.com

sigma-aldrich.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
 ALDRICH • BOX 355 • MILWAUKEE • WISCONSIN • USA

ALDRICH
 Advancing Science

Aldrichimica ACTA

VOL. 39, NO. 2 • 2006

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
 6000 N. Teutonia Ave.
 Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
 FAX 800-325-5052 (USA)
 or 414-438-2199
 Mail P.O. Box 2060
 Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
 Technical Service 800-231-8327
 SAFC™ 800-244-1173
 Custom Synthesis 800-244-1173
 Flavors & Fragrances 800-227-4563
 International 414-438-3850
 24-Hour Emergency 414-438-3850
 Web Site sigma-aldrich.com
 Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
 P.O. Box 355, Milwaukee, WI 53201, USA

Subscriptions

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
 Mail: **Attn: Mailroom**
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
 P.O. Box 355
 Milwaukee, WI 53201-9358
 Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the inside back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.


All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2006 Sigma-Aldrich Co.

"PLEASE BOTHER US."

Professor John G. Ekerdt of the University of Texas at Austin and his student, Wyatt Winkenwerder, kindly suggested that we offer (1,5-cyclooctadiene)(1,3,5-cyclooctatriene)ruthenium, or Ru(cod)(cot). This complex is used in the preparation of monodisperse ruthenium nanoparticles for catalysis,^{1,2} as well as a highly selective catalyst for amine alkylations, [2+2] cycloadditions,³ and enyne generation

(1) Pelzer, K.; Philippot, K.; Chaudret, B. *Z. Phys. Chem.* **2003**, *217*, 1539. (2) Hulea, V.; Brunel, D.; Galarneau, A.; Philippot, K.; Chaudret, B.; Kooymann, P. J.; Fajula, F. *Microporous Mesoporous Mater.* **2005**, *79*, 185. (3) Mitsudo, T.; Suzuki, T.; Zhang, S.-W.; Imai, D.; Fujita, K.; Manabe, T.; Shiotsuki, M.; Watanabe, Y.; Wada, K.; Kondo, T. *J. Am. Chem. Soc.* **1999**, *121*, 1839



654418-1G (1,5-Cyclooctadiene)(1,3,5-cyclooctatriene)- 1 g \$165.00
ruthenium, Ru(cod)(cot)

Naturally, we made this versatile reagent. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the inside back cover.

TABLE OF CONTENTS

Chiral, Poly(Rare-Earth Metal) Complexes in Asymmetric Catalysis	31
<i>Masakatsu Shibasaki,* Motomu Kanai, and Shigeaki Matsunaga</i> The University of Tokyo	
Organic Synthesis and Device Testing for Molecular Electronics	47
<i>Dustin K. James and James M. Tour,* Rice University</i>	

ABOUT OUR COVER

View from Vaekero near Christiania (oil on canvas, 60.5 × 96.5 cm) was painted by the Norwegian romantic painter Johan Christian Dahl in 1827. Dahl studied in Dresden and was directly influenced by his teacher and friend, the German painter Casper David Friedrich. Dahl's paintings also show his strong interest in the work of seventeenth century Dutch landscape painters such as Jacob van Ruisdael.



Photograph © Board of Trustees, National Gallery of Art, Washington

Dahl visited Christiania, present-day Oslo, in the summer of 1826. The following winter in Dresden, Dahl painted *View from Vaekero near Christiania* from memory for the Hamburger Kunstverein artists' cooperative. In this moody and melancholy nocturne, Dahl invites the viewer to imagine a romantic moonlit evening complete with sand, sea, and sky. His use of successive bands of light and dark clouds against a pink-and-blue backdrop shows an alluring distance, possibly unattainable. Harmoniously cascading hills, which meld into an illuminated sea, may also suggest adventure. Ethereal light and drying fishnets seem to envelop the mysterious, solitary couple, who stand in the center foreground contemplating the quixotic setting. True to his romantic spirit, Dahl presents a thought-provoking and poignant scene, allowing us to do what paintings should make us do—dream

This painting was purchased for the National Gallery of Art by the Patrons' Permanent Fund.


NEW Shibasaki Catalysts: La, Y, Gd, and Sm Trisisopropoxides

Rare-Earth Metals Used in Diversity-Oriented Organic Transformations

Product Highlights

- Dramatically enhance selectivities by varying the nature of the rare-earth (RE) metal and the ratio of catalyst to reaction partners.
- In most cases, the metal complexes are insensitive to oxygen after preparation.
- Bifunctional: Can perform effectively as both a Brønsted base and a Lewis acid.
- The catalysts can be recovered and recycled without loss of selectivities.
- RE catalyst systems can effectively facilitate a broad range of organic reactions.

Shibasaki and co-workers have developed rare-earth (RE) metal catalysts, utilized in conjunction with a variety of chiral ligands, to effect asymmetric transformations ranging from the formation of quaternary chiral centers to the epoxidation of unsaturated substrates. The Shibasaki research group has published extensively in the field of RE-metal catalysis and has optimized reaction conditions to afford high selectivities in C–C and C–O bond-forming reactions. Sigma-Aldrich is pleased to offer an array of RE-metal pre-catalysts that can be paired with our growing line of chiral ligands to accelerate your research discoveries.



References: (a) Masumoto, S. et al. *J. Am. Chem. Soc.* **2003**, *125*, 5634. (b) Kim, Y. S. et al. *J. Am. Chem. Soc.* **2000**, *122*, 6506. (c) Kakei, H. et al. *J. Am. Chem. Soc.* **2005**, *127*, 8962. (d) Nemoto, T. et al. *J. Am. Chem. Soc.* **2002**, *124*, 14544. (e) Sasai, H. et al. *J. Am. Chem. Soc.* **1993**, *115*, 10372. (f) Mita, T. et al. *J. Am. Chem. Soc.* **2005**, *127*, 11252. (g) Gröger, H. et al. *J. Am. Chem. Soc.* **1998**, *120*, 3089. (h) Yoshikawa, N. et al. *J. Am. Chem. Soc.* **2001**, *123*, 2466. (i) Shibasaki, M. et al. *Chem. Rev.* **2002**, *102*, 2187. (j) Yabu, K. et al. *J. Am. Chem. Soc.* **2001**, *123*, 9908. (k) Nemoto, T. et al. *J. Am. Chem. Soc.* **2001**, *123*, 2725

Lanthanum(III) isopropoxide

[19446-52-7]

C₉H₂₁LaO₃

FW: 316.17

665193-500MG

500 mg

\$45.00

665193-3G

3 g

185.00

Yttrium(III) isopropoxide

[2172-12-5]

C₉H₂₁YO₃

FW: 266.17

665916-500MG

500 mg

\$45.00

665916-3G

3 g

199.00

NEW

Gadolinium(III) isopropoxide

[14532-05-9]

C₉H₂₁GdO₃

FW: 334.51

663948-500MG

500 mg

\$75.00

663948-3G

3 g

365.00

NEW

NEW

Samarium(III) isopropoxide

[3504-40-3]

C₉H₂₁SmO₃

FW: 327.62

410136-1G

1 g

\$26.70

sigma-aldrich.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
ALDRICH • BOX 355 • MILWAUKEE • WISCONSIN • USA

ALDRICH
Advancing Science

Sigma-Aldrich Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA, S.A.
Tel: 54 11 4556 1472
Fax: 54 11 4552 1698

Australia

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: 612 9841 0555
Fax: 612 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: 43 1 605 81 10
Fax: 43 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA.
Free Tel: 0800-14747
Free Fax: 0800-14745
Tel: 03 899 13 01
Fax: 03 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Tel: 55 11 3732-3100
Fax: 55 11 3733-5151

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 800-565-1400
Free Fax: 800-265-3858
Tel: 905-829-9500
Fax: 905-829-9292

China

SIGMA-ALDRICH CHINA INC.
Tel: 86-21-6386 2766
Fax: 86-21-6386 3966

Czech Republic

SIGMA-ALDRICH s.r.o.
Tel: +420 246 003 200
Fax: +420 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: 43 56 59 10
Fax: 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND
Tel: (09) 350 9250
Fax: (09) 350 9255

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Tel appel gratuit: 0800 211 408
Fax appel gratuit: 0800 031 052

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800-51 55 000
Free Fax: 0800-649 00 00

Greece

SIGMA-ALDRICH (O.M.) LTD
Tel: 30 210 9948010
Fax: 30 210 9943831

Hungary

SIGMA-ALDRICH Kft
Tel: 06-1-235-9054
Fax: 06-1-269-6470
Ingyenes zöld telefon: 06-80-355-355
Ingyenes zöld fax: 06-80-344-344

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: 91-80-4112-7272
New Delhi: 91-11-4165 4255
Mumbai: 91-22-2570 2364
Hyderabad: 91-40-5584 5488
Fax
Bangalore: 91-80-4112-7473
New Delhi: 91-11-4165 4266
Mumbai: 91-22-2579 7589
Hyderabad: 91-40-5584 5466

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: 353 1 4041900
Fax: 353 1 4041910

Israel

SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1-800-70-2222
Tel: 08-948-4100
Fax: 08-948-4200

Italy

SIGMA-ALDRICH S.r.l.
Telefono: 02 33417310
Fax: 02 38010737
Numero Verde: 800-827018

Japan

SIGMA-ALDRICH JAPAN K.K.
Tokyo Tel: 03 5796 7300
Tokyo Fax: 03 5796 7315

Korea

SIGMA-ALDRICH KOREA
Tel: 031-329-9000
Fax: 031-329-9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: 603-56353321
Fax: 603-56354116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01-800-007-5300
Free Fax: 01-800-712-9920

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Tel Gratis: 0800-0229088
Fax Gratis: 0800-0229089
Tel: 078-6205411
Fax: 078-6205421

New Zealand

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: 61 2 9841 0500
Fax: 61 2 9841 0500

Norway

SIGMA-ALDRICH NORWAY AS
Tel: 23 17 60 60
Fax: 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: 061 829 01 00
Fax: 061 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202180
Free Fax: 800 202178
Tel: 21 9242555
Fax: 21 9242610

Russia

SIGMA-ALDRICH RUS, LLC
Tel: +7 (095) 621-5828/6037
Fax: +7 (095) 975-4792

Singapore

SIGMA-ALDRICH PTE. LTD
Tel: 65-67791200
Fax: 65-67791822

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: 27 11 979 1188
Fax: 27 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA S.A.
Free Tel: 900 101376
Free Fax: 900 102028
Tel: 91 661 99 77
Fax: 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: 020-350510
Fax: 020-352522
Outside Sweden Tel: +46 8 7424200
Outside Sweden Fax: +46 8 7424243

Switzerland

SIGMA-ALDRICH CHEMIE GmbH
Swiss Free Call: 0800 80 00 80
Tel: +41 81 755 2828
Fax: +41 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717181
Free Fax: 0800 378785
Tel: 01747 833000
Fax: 01747 833313
SAFC (UK): 01202 712305

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-free: 800-325-3010
Call Collect: 314-771-5750
Toll-Free Fax: 800-325-5052
Tel: 314-771-5765
Fax: 314-771-5757

Internet

sigma-aldrich.com

Sigma-Aldrich Career Opportunities

As a leading Life Science and High Technology company, we are always looking for talented individuals to join our team. At Sigma-Aldrich we value the contributions of our employees, and recognize the impact they have on our success. We strive to foster creativity and innovation, and encourage professional development.

We offer a highly motivational and rewarding work environment with an attractive salary and excellent benefits, including:

Medical	Pension
Dental	401(k)
Flexible Spending Accounts	Life Insurance
Tuition Reimbursement	

Accelerating Customers' Success through Leadership in Life Science, High Technology and Service

Over one million scientists and technologists use our biochemical and organic chemical products in scientific and genomic research, biotechnology, pharmaceutical development, disease diagnosis and manufacturing. Sigma-Aldrich operates in 35 countries, and has over 6,800 employees providing excellent service worldwide.

Learn more about our career opportunities by visiting
our award-winning Web site at sigma-aldrich.com/careers.

sigma-aldrich.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
SIGMA-ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI 63178 • USA

 SIGMA-ALDRICH

Sigma-Aldrich Corporation is an equal opportunity employer

Be First in Science with the Just-Released 2006–2007 Catalog from Sigma

2006–2007

Biochemicals, Reagents & Kits

for Life Science Research

INNOVATIVE TECHNOLOGIES

WORLD'S LARGEST SELECTION

OVER 3,000 NEW PRODUCTS

WORLD-CLASS SERVICE

LEADER IN WEB CONTENT

COMMITMENT TO QUALITY
PRODUCTS AND SYSTEMS

Available Now!

To order your copy, visit us online at
sigma-aldrich.com/catalog
and find out why it's the catalog
scientists reach for first.

sigma-aldrich.com

 **SIGMA**[®]



ALDRICH[®]

P.O. Box 14508
St. Louis, MO 63178
USA

Return Service Requested

PSRT STD
U.S. POSTAGE PAID
MADISON, WI
PERMIT NO. 2223

ASYMMETRIC SYNTHESIS ENABLED BY METAL-FREE CATALYSIS

Aldrichimica ACTA

VOL. 39, NO. 3 • 2006



Enzymes in Organic Synthesis

Modern Strategies in Organic Catalysis

sigma-aldrich.com



ALDRICH[®]

Advancing Science


New Products from Aldrich R&D

Aldrich Is Pleased to Offer the Latest Catalysis Technology


PEPPSI™ Cross-Coupling Catalyst

The PEPPSI™ catalyst offers efficient mediation of C–C and C–N bond-forming processes and robust air stability at a competitive cost structure. This powerful new catalyst successfully promotes Negishi, Suzuki, and Kumada couplings as well as Buchwald–Hartwig aminations in excellent yields.¹ These achievements make it attractive for widespread application in the global research and fine chemicals market.

Negishi Coupling




Suzuki Coupling



Kumada Coupling



Buchwald–Hartwig Amination



(1) (a) O' Brien, C. J. et al. *Eur. J. Org. Chem.* **2006**, in press. (b) O' Brien, C. J. et al. *Eur. J. Org. Chem.* **2006**, in press. (c) O' Brien, C. J. et al. *Manuscript in preparation* (d) O' Brien, C. J. et al. *Manuscript in preparation*


[1,3-Bis(2,6-diisopropylphenyl)imidazol-2-ylidene](3-chloropyridyl)palladium(II) dichloride, 98%

PEPPSI™-IPr Catalyst

669032

C₃₂H₄₀Cl₃N₃Pd

FW: 679.46




1 g \$60.00
5 g 270.00

Please see pages 76–77 for additional new products. For competitive quotes on larger quantities, contact www.safcglobal.com

Heriot-Watt Carbonylation Catalysts

These palladium catalysts, sold under exclusive license from Heriot-Watt University (Scotland),¹ efficiently carbonylate benzyl halides in methanol under low pressures and temperatures.²



(1) Manufactured by Aldrich under exclusive license from Heriot-Watt University, PCT/GB2005/002738 patent pending. (2) (a) Jones, R. V. H. et al. *Tetrahedron Lett.* **2005**, 46, 8695. (b) Lindsell, W. E. et al. *Organometallics* **2005**, 24, 1119.


Bromo[2-(hydroxy-κO)methyl]phenylmethyl-κC(triphenylphosphine)palladium(II)

666327

[849417-33-0]

C₂₁H₁₅BrOPPd

FW: 569.77




250 mg \$30.00
1 g 95.00

2-[Bis(triphenylphosphine)palladium(II) bromide]benzyl alcohol, 96%

665932

C₄₃H₃₇BrOP₂Pd


FW: 818.02



250 mg \$30.00
1 g 95.00

Trost Bis-ProPhenol Ligands

The combination of these phenolic ligands and diethylzinc generates catalyst systems that are capable of effecting a variety of asymmetric reactions. This methodology has been applied to asymmetric aldol condensations,¹ Mannich-type reactions,² and Henry reactions.³



(1) (a) Trost, B. M.; Ito, H. *J. Am. Chem. Soc.* **2000**, 122, 12003. (b) Trost, B. M.; Yeh, V. S. C. *Org. Lett.* **2002**, 4, 3513. (2) Trost, B. M.; Terrell, L. R. *J. Am. Chem. Soc.* **2003**, 125, 338. (3) Trost, B. M. et al. *Org. Lett.* **2002**, 4, 2621.


(R,R)-(-)-2,6-Bis[2-(hydroxydiphenylmethyl)-1-pyrrolidinylmethyl]-4-methylphenol, 95%

Trost (R,R)-Bis-ProPhenol Ligand

667625

C₄₃H₄₆N₂O₃

FW: 638.84



1 g \$20.50
5 g 81.20


(S,S)-(+)-2,6-Bis[2-(hydroxydiphenylmethyl)-1-pyrrolidinylmethyl]-4-methylphenol, 95%

Trost (S,S)-Bis-ProPhenol Ligand

668370

C₄₃H₄₆N₂O₃

FW: 638.84



1 g \$20.50
5 g 81.20

sigma-aldrich.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
ALDRICH • BOX 355 • MILWAUKEE • WISCONSIN • USA

ALDRICH
Advancing Science

PEPPSI™ patent pending

Aldrichimica ACTA

VOL. 39, NO. 3 • 2006

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-231-8327
SAFC™ 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 414-438-3850
24-Hour Emergency 414-438-3850
Web Site sigma-aldrich.com
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201, USA

Subscriptions

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 355
Milwaukee, WI 53201-9358
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the inside back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.


All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2006 Sigma-Aldrich Co.

"PLEASE BOTHER US."

Professor Barry Trost of Stanford University kindly suggested that we offer these enantiomeric Bis-ProPhenol ligands. Upon treatment of a Bis-ProPhenol ligand with diethylzinc, a versatile catalyst is formed, which is capable of performing a variety of reactions in an asymmetric fashion. Aldol, imine-aldol, and nitro-aldol condensations all proceed with impressive levels of stereocontrol.¹ Additionally, the catalyst system promotes the asymmetric alkynylation of unsaturated aldehydes.²

(1) See the references on the facing page. (2) Trost, B. M. et al. *J. Am. Chem. Soc.* **2006**, *128*, 8



667625	(R,R)-(-)-2,6-Bis[2-(hydroxydiphenylmethyl)-1-pyrrolidinyl-methyl]-4-methylphenol, 95% Trost (R,R)-Bis-ProPhenol Ligand	1 g \$20.50
		5 g 81.20
668370	(S,S)-(+)-2,6-Bis[2-(hydroxydiphenylmethyl)-1-pyrrolidinyl-methyl]-4-methylphenol, 95% Trost (S,S)-Bis-ProPhenol Ligand	1 g \$20.50
		5 g 81.20

Naturally, we made these useful ligands. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the inside back cover.

TABLE OF CONTENTS

Enzymes in Organic Synthesis: Aldolase-Mediated Synthesis of Iminocyclitols and Novel Heterocycles	63
<i>Lisa J. Whalen and Chi-Huey Wong,* The Scripps Research Institute</i>	
Modern Strategies in Organic Catalysis: The Advent and Development of Iminium Activation	79
<i>Gérald Lelais and David W. C. MacMillan,* California Institute of Technology</i>	

ABOUT OUR COVER

View of the Tiber near Perugia (oil on canvas, 98.0 × 161.5 cm) was painted in Italy between 1872 and 1874 by George Inness, the American landscape artist. Inness was born in 1825 in Newburgh, New York, and was raised in New York City and New Jersey. Although he had little formal artistic training, he developed his painting style through his association with artists from the American Hudson River School and his frequent visits to Europe, where 17th-century old masters and French Barbizon School landscape painters influenced him. Emanuel Swedenborg, the 18th-century Swedish scientist and theologian who stressed that the spiritual world was as much a reality as the material world, also inspired Inness to express his personal vision of spiritual harmonies in nature.



Photograph © Board of Trustees, National Gallery of Art, Washington

Subtle, yet dramatic atmospheric effects of harmonious color, light and shade, and mood and meaning are displayed in this painting. The dynamic diagonal created by the middle ground landscape reinforces the suggestion of depth and recession and separates the foreground from the background. A sense of the human proportion is understood by looking at the figures in the near ground. This poetic and spiritual view of the Perugia area is remarkably rendered with a sense of peace and calm.

Purchased for the National Gallery of Art, Washington, DC, through the Ailsa Mellon Bruce Fund.

Biocatalytic Aldol Reactions in Organic Synthesis


The construction of C–C bonds with complete stereochemical control of the reaction course is a key goal in organic synthesis. The desired relative and absolute configurations of the newly formed stereogenic centers can be achieved by using chiral starting materials and chiral auxiliaries. Catalytic approaches towards asymmetric aldol reactions have been a special focus in recent years.¹ Aldolases and catalytic Aldolase antibodies represent valuable tools for the biocatalytic asymmetric C–C-bond formation.

Sigma-Aldrich offers a comprehensive range of chiral auxiliaries, organocatalysts, organometallic catalysts, ligands, and biocatalysts for efficient asymmetric aldol reactions. Please visit our Web site at sigma-aldrich.com and take a look at the Sigma-Aldrich Enzyme Explorer for more biocatalysts.

(1) Machajewski, T. D.; Wong, C.-H. *Angew. Chem., Int. Ed.* **2000**, *39*, 1352.

FDP Aldolase


E.C. 4.1.2.13



Cat. No.	Name	Pack Sizes
94864	Fructose-1,6-bisphosphate Aldolase from <i>Staphylococcus carnosus</i>	5 mg
05518	Aldolase from rabbit muscle	10 mg 50 mg
A2714	Aldolase from rabbit muscle	100 units 200 units 500 units
A9329	Aldolase from spinach	25 units

DERA Aldolase


E.C. 4.1.2.4



Cat. No.	Name	Pack Sizes
91252	2-Deoxy-D-ribose 5-phosphate Aldolase, E. coli K12, recombinant from <i>Escherichia coli</i>	25 mg 100 mg
41228	2-Deoxyribose-5-phosphate Aldolase from <i>Lactobacillus plantarum</i>	10 mg


KHG Aldolase

E.C. 4.1.3.16



Cat. No.	Name	Pack Sizes
59892	4-Hydroxy-2-oxoglutarate Aldolase from <i>Escherichia coli</i>	10 mg

Aldolase Antibody 38C2



Cat. No.	Name	Pack Sizes
479950	Aldolase Antibody 38C2, murine catalytic monoclonal antibody, contains phosphate buffer salts (pbs)	10 mg
481572	Aldolase Antibody 38C2, murine catalytic monoclonal antibody	10 mg

Other Aldol Reaction Biocatalysts

Cat. No.	Name	E.C. Number	Pack Sizes
47153	<i>N</i> -Acetylneuraminic acid Aldolase from <i>Escherichia coli</i>	4.1.3.3	25 mg 100 mg
A6680	<i>N</i> -Acetylneuraminic acid Aldolase from <i>Escherichia coli</i>	4.1.3.3	25 units
96586	L-Threonine Aldolase from <i>Pseudomonas putida</i>	4.1.2.5	10 mg 50 mg
67891	3-Deoxy-D-manno-octulosonate Aldolase from <i>Escherichia coli</i>	4.1.2.23	10 mg
527858	Aldolase Antibody 84G3, murine catalytic monoclonal antibody		10 mg

For more information, please visit us at sigma-aldrich.com.

sigma-aldrich.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
SIGMA-ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI 63178 • USA



SIGMA-ALDRICH

Sigma-Aldrich Worldwide Locations

- Argentina**
SIGMA-ALDRICH DE ARGENTINA, S.A.
Tel: 54 11 4556 1472
Fax: 54 11 4552 1698
- Australia**
SIGMA-ALDRICH PTY., LIMITED
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: 612 9841 0555
Fax: 612 9841 0500
- Austria**
SIGMA-ALDRICH HANDELS GmbH
Tel: 43 1 605 81 10
Fax: 43 1 605 81 20
- Belgium**
SIGMA-ALDRICH NV/SA.
Free Tel: 0800-14747
Free Fax: 0800-14745
Tel: 03 899 13 01
Fax: 03 899 13 11
- Brazil**
SIGMA-ALDRICH BRASIL LTDA.
Tel: 55 11 3732-3100
Fax: 55 11 3733-5151
- Canada**
SIGMA-ALDRICH CANADA LTD.
Free Tel: 800-565-1400
Free Fax: 800-265-3858
Tel: 905-829-9500
Fax: 905-829-9292
- China**
SIGMA-ALDRICH (SHANGHAI)
TRADING CO. LTD.
Tel: +86-21-61415566
Fax: +86-21-61415567
- Czech Republic**
SIGMA-ALDRICH S.R.O.
Tel: +420 246 003 200
Fax: +420 246 003 291
- Denmark**
SIGMA-ALDRICH DENMARK A/S
Tel: 43 56 59 10
Fax: 43 56 59 05
- Finland**
SIGMA-ALDRICH FINLAND
Tel: (09) 350 9250
Fax: (09) 350 92555
- France**
SIGMA-ALDRICH CHIMIE S.à.r.l.
Tel appel gratuit: 0800 211 408
Fax appel gratuit: 0800 031 052
- Germany**
SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800-51 55 000
Free Fax: 0800-649 00 00
- Greece**
SIGMA-ALDRICH (O.M.) LTD
Tel: 30 210 9948010
Fax: 30 210 9943831
- Hungary**
SIGMA-ALDRICH Kft
Tel: 06-1-235-9054
Fax: 06-1-269-6470
Ingyenes zöld telefon: 06-80-355-355
Ingyenes zöld fax: 06-80-344-344
- India**
SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: 91-80-4112-7272
New Delhi: 91-11-4165-4255
Mumbai: 91-22-2570-2364
Hyderabad: 91-40-5584-5488
Fax
Bangalore: 91-80-4112-7473
New Delhi: 91-11-4165-4266
Mumbai: 91-22-2579-7589
Hyderabad: 91-40-5584-5466
- Ireland**
SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: 353 1 4041900
Fax: 353 1 4041910
- Israel**
SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1-800-70-2222
Tel: 08-948-4100
Fax: 08-948-4200
- Italy**
SIGMA-ALDRICH S.r.l.
Telefono: 02 33417310
Fax: 02 38010737
Numero Verde: 800-827018
- Japan**
SIGMA-ALDRICH JAPAN K.K.
Tokyo Tel: 03 5796 7300
Tokyo Fax: 03 5796 7315
- Korea**
SIGMA-ALDRICH KOREA
Tel: 031-329-9000
Fax: 031-329-9090
- Malaysia**
SIGMA-ALDRICH (M) SDN. BHD
Tel: 603-56353321
Fax: 603-56354116
- Mexico**
SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01-800-007-5300
Free Fax: 01-800-712-9920
- The Netherlands**
SIGMA-ALDRICH CHEMIE BV
Tel Gratis: 0800-0229088
Fax Gratis: 0800-0229089
Tel: 078-6205411
Fax: 078-6205421
- New Zealand**
SIGMA-ALDRICH PTY., LIMITED
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: 61 2 9841 0500
Fax: 61 2 9841 0500
- Norway**
SIGMA-ALDRICH NORWAY AS
Tel: 23 17 60 60
Fax: 23 17 60 50
- Poland**
SIGMA-ALDRICH Sp. z o.o.
Tel: 061 829 01 00
Fax: 061 829 01 20
- Portugal**
SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202180
Free Fax: 800 202178
Tel: 21 9242555
Fax: 21 9242610
- Russia**
SIGMA-ALDRICH RUS, LLC
Tel: +7 (495) 621-5828/5923/6037
Fax: +7 (495) 621-5828/5923/6037
- Singapore**
SIGMA-ALDRICH PTE. LTD.
Tel: 65-67791200
Fax: 65-67791822
- South Africa**
SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD.
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: 27 11 979 1188
Fax: 27 11 979 1119
- Spain**
SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 900 101376
Free Fax: 900 102028
Tel: 91 661 99 77
Fax: 91 661 96 42
- Sweden**
SIGMA-ALDRICH SWEDEN AB
Tel: 020-350510
Fax: 020-352522
Outside Sweden Tel: +46 8 7424200
Outside Sweden Fax: +46 8 7424243
- Switzerland**
SIGMA-ALDRICH CHEMIE GmbH
Swiss Free Call: 0800 80 00 80
Tel: +41 81 755 2828
Fax: +41 81 755 2815
- United Kingdom**
SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717181
Free Fax: 0800 378785
Tel: 01747 833000
Fax: 01747 833313
SAFC (UK): 01202 712305
- United States**
SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-Free: (800) 325-3010
Call Collect: (314) 771-5750
Toll-Free Fax: (800) 325-5052
Tel: (314) 771-5765
Fax: (314) 771-5757
- Internet**
sigma-aldrich.com

Sigma-Aldrich Career Opportunities

As a leading Life Science and High Technology company, we are always looking for talented individuals to join our team. At Sigma-Aldrich we value the contributions of our employees, and recognize the impact they have on our success. We strive to foster creativity and innovation, and encourage professional development.

We offer a highly motivational and rewarding work environment with an attractive salary and excellent benefits, including:

Medical	Pension
Dental	401(k)
Flexible Spending Accounts	Life Insurance
Tuition Reimbursement	

Accelerating Customers' Success through Leadership in Life Science, High Technology and Service

Over one million scientists and technologists use our biochemical and organic chemical products in scientific and genomic research, biotechnology, pharmaceutical development, disease diagnosis and manufacturing. Sigma-Aldrich operates in 35 countries, and has over 6,800 employees providing excellent service worldwide.

Learn more about our career opportunities by visiting
our award-winning Web site at sigma-aldrich.com/careers.

sigma-aldrich.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
SIGMA-ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI 63178 • USA

Sigma-Aldrich Corporation is an equal opportunity employer



SIGMA-ALDRICH

Introducing *Material Matters*

FREE

Material Matters
Clean. Strong. Durable. Performance.

**Polymerization
for Advanced Applications**

Vol. 1 No. 1



Aldrichimica ACTA and *ChemFiles* have long been respected publications designed to inspire advancement in chemistry. Now, Sigma-Aldrich is proud to launch a **technical guide** dedicated to researchers with a focus in **materials science and technology**—*Material Matters*. Each issue will feature technical reviews, application notes, and highlight new products and services.

Material Matters is a quarterly publication that is yours **FREE** just by completing our online form.

Subscribe for **FREE** at
sigma-aldrich.com/mm

ALDRICH

P.O. Box 14508
St. Louis, MO 63178
USA

Return Service Requested

Presorted Standard
U.S. Postage
PAID
Des Moines, IA
Permit No. 589

PEPPSI™: A NEW GENERATION OF AIR-STABLE Pd PRECATALYSTS

Aldrichimica **ACTA**

VOL. 39, NO. 4 • 2006



Pd-N-Heterocyclic Carbene (NHC) Catalysts
for Cross-Coupling Reactions

sigma-aldrich.com



ALDRICH®
Advancing Science

Aldrichimica ACTA

VOL. 39, NO. 4 • 2006

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Service

Customer Inquiries 800-325-3010
Technical Service 800-231-8327
SAFC™ 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 414-438-3850
24-Hour Emergency 414-438-3850
Web Site sigma-aldrich.com
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201, USA

Subscriptions

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 355
Milwaukee, WI 53201-9358
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the inside back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.


All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2006 Sigma-Aldrich Co.

"PLEASE BOTHER US."

Professor Andrew Whiting of the University of Durham, U.K., kindly suggested that we offer vinylboronic acid 2-methyl-2,4-pentanediol ester. This reagent functions as a vinyl dianion equivalent, since it can react at either end of the double bond under Suzuki–Miyaura or Heck coupling conditions. This compound also exhibits enhanced stability and lower volatility than related vinylboronate esters.^{1–3}

(1) Lightfoot, A. P.; Maw, G.; Thirsk, C.; Twiddle, S. J. R.; Whiting, A. *Tetrahedron Lett.* **2003**, *44*, 7645
(2) Lightfoot, A. P.; Twiddle, S. J. R.; Whiting, A. *Org. Biomol. Chem.* **2005**, *3*, 3167. (3) Lightfoot, A. P.; Twiddle, S. J. R.; Whiting, A. *Synlett* **2005**, 529



673641 Vinylboronic acid 2-methyl-2,4-pentanediol ester, 95% 1 g \$30.20
(4,4,6-Trimethyl-2-vinyl-1,3,2-dioxaborinane) 10 g 211.50

Naturally, we made this useful reagent. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the inside back cover.

TABLE OF CONTENTS

Pd–N-Heterocyclic Carbene (NHC) Catalysts for Cross-Coupling Reactions 97
*Eric Assen B. Kantchev, Christopher J. O'Brien, and Michael G. Organ,** York University, Toronto, Canada

ABOUT OUR COVER

No one will be surprised to read that the painting on our cover, *Ships Riding on the Seine at Rouen* (1872/1873, oil on canvas, 37.8 × 46.6 cm), was painted by the archetypal Impressionist, Claude Monet (French, 1840–1926). His distinctive style is recognizable, fresh, and still sought after over 130 years after his radical and expressive brushwork shocked the art world in Paris in the 1860s and 1870s.

Monet's devotion to the ideals of Impressionism lasted throughout his long and prolific life. It is fitting that the word Impressionism, which defines the movement and style, was coined by an art critic reviewing Monet's painting *Impression: Sunrise*, 1872 (Musée Marmottan, Paris).

This delightful, small painting displays large and small sailboats comfortably resting on the river Seine on a sunny, warm, and idyllic day. Monet's choice of primarily soft blues and greens creates an alluring calm and relaxed atmosphere. He characteristically captures the illusion of moving water by suggesting with broad brushstrokes that the landscape, sky, clouds, and boats are reflected in the gently rolling water.



Photograph © Board of Trustees, National Gallery of Art, Washington

This painting is in the Ailsa Mellon Bruce Collection of Small French Paintings at the National Gallery of Art, Washington, DC.

Relax and Enjoy Your Research!

Access 40,000+ products with one click.

700+
Boronic Acids
& Derivatives

170+
Pd Catalysts

250+
Privileged
Ligands

6000+
Heterocycles

The success of your research relies heavily on having the best products for your experiments. Tired of browsing through various catalogs and databases for the research products you need? Look no further! Sigma-Aldrich introduces the new **Chem Product Central™**, featuring the industry-leading range of 40,000+ products including **cutting-edge catalysts and ligands, novel and classical synthetic reagents, and diverse libraries of building blocks**. Unique to Chem Product Central™, all chemical products are categorized based on compound classes and application areas for easy browsing. In addition, new products are added on a weekly basis, providing you with the most innovative products for your research.

Check out the new Chem Product Central™ at
sigma-aldrich.com/gochem.

Built for Chemists
by Chemists.™

sigma-aldrich.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
SIGMA-ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI 63178 • USA

Chem Product Central and Built for Chemists by Chemists. are trademarks of Sigma-Aldrich Biotechnology, L.P.



SIGMA-ALDRICH

Sigma-Aldrich Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA, S.A.
Tel: 54 11 4556 1472
Fax: 54 11 4552 1698

Australia

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: 612 9841 0555
Fax: 612 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: 43 1 605 81 10
Fax: 43 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA.
Free Tel: 0800-14747
Free Fax: 0800-14745
Tel: 03 899 13 01
Fax: 03 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Tel: 55 11 3732-3100
Fax: 55 11 3733-5151

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 800-565-1400
Free Fax: 800-265-3858
Tel: 905-829-9500
Fax: 905-829-9292

China

SIGMA-ALDRICH (SHANGHAI)
TRADING CO. LTD.
Tel: +86-21-61415566
Fax: +86-21-61415567

Czech Republic

SIGMA-ALDRICH S.R.O.
Tel: +420 246 003 200
Fax: +420 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: 43 56 59 10
Fax: 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND
Tel: (09) 350 9250
Fax: (09) 350 9255

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Tel appel gratuit: 0800 211 408
Fax appel gratuit: 0800 031 052

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800-51 55 000
Free Fax: 0800-649 00 00

Greece

SIGMA-ALDRICH (O.M.) LTD
Tel: 30 210 9948010
Fax: 30 210 9943831

Hungary

SIGMA-ALDRICH Kft
Tel: 06-1-235-9054
Fax: 06-1-269-6470
Ingyenes zöld telefon: 06-80-355-355
Ingyenes zöld fax: 06-80-344-344

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: 91-80-6621-9600
New Delhi: 91-11-4165-4255
Mumbai: 91-22-2570-2364
Hyderabad: 91-40-6684-5488
Fax
Bangalore: 91-80-6621-9650
New Delhi: 91-11-4165-4266
Mumbai: 91-22-2579-7589
Hyderabad: 91-40-6684-5466

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: 353 1 4041900
Fax: 353 1 4041910

Israel

SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1-800-70-2222
Tel: 08-948-4100
Fax: 08-948-4200

Italy

SIGMA-ALDRICH S.r.l.
Telefono: 02 33417310
Fax: 02 38010737
Numero Verde: 800-827018

Japan

SIGMA-ALDRICH JAPAN K.K.
Tokyo Tel: 03 5796 7300
Tokyo Fax: 03 5796 7315

Korea

SIGMA-ALDRICH KOREA
Tel: 031-329-9000
Fax: 031-329-9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: 603-56353321
Fax: 603-56354116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01-800-007-5300
Free Fax: 01-800-712-9920

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Tel Gratis: 0800-0229088
Fax Gratis: 0800-0229089
Tel: 078-6205411
Fax: 078-6205421

New Zealand

SIGMA-ALDRICH PTY., LIMITED
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: 61 2 9841 0500
Fax: 61 2 9841 0500

Norway

SIGMA-ALDRICH NORWAY AS
Tel: 23 17 60 60
Fax: 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: 061 829 01 00
Fax: 061 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202180
Free Fax: 800 202178
Tel: 21 9242555
Fax: 21 9242610

Russia

SIGMA-ALDRICH RUS, LLC
Tel: +7 (495) 621-5828/5923/6037
Fax: +7 (495) 621-5828/5923/6037

Singapore

SIGMA-ALDRICH PTE. LTD.
Tel: 65-67791200
Fax: 65-67791822

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD.
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: 27 11 979 1188
Fax: 27 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 900 101376
Free Fax: 900 102028
Tel: 91 661 99 77
Fax: 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: 020-350510
Fax: 020-352522
Outside Sweden Tel: +46 8 7424200
Outside Sweden Fax: +46 8 7424243

Switzerland

SIGMA-ALDRICH CHEMIE GmbH
Swiss Free Call: 0800 80 00 80
Tel: +41 81 755 2828
Fax: +41 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717181
Free Fax: 0800 378785
Tel: 01747 833000
Fax: 01747 833313
SAFC (UK): 01202 712305

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-Free: (800) 325-3010
Call Collect: (314) 771-5750
Toll-Free Fax: (800) 325-5052
Tel: (314) 771-5765
Fax: (314) 771-5757

Internet

sigma-aldrich.com

Sigma-Aldrich Career Opportunities

As a leading Life Science and High Technology company, we are always looking for talented individuals to join our team. At Sigma-Aldrich we value the contributions of our employees, and recognize the impact they have on our success. We strive to foster creativity and innovation, and encourage professional development.

We offer a highly motivational and rewarding work environment with an attractive salary and excellent benefits, including:

Medical	Pension
Dental	401(k)
Flexible Spending Accounts	Life Insurance
Tuition Reimbursement	

Accelerating Customers' Success through Leadership in Life Science, High Technology and Service

Over one million scientists and technologists use our biochemical and organic chemical products in scientific and genomic research, biotechnology, pharmaceutical development, disease diagnosis and manufacturing. Sigma-Aldrich operates in 35 countries, and has over 7,200 employees providing excellent service worldwide.

Learn more about our career opportunities by visiting
our award-winning Web site at sigma-aldrich.com/careers.

sigma-aldrich.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
SIGMA-ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI 63178 • USA

SIGMA-ALDRICH

Sigma-Aldrich Corporation is an equal opportunity employer.

Looking For More Free Updates in Chemistry from Sigma-Aldrich?

ChemFiles

ChemFiles are monthly new-product application guides in chemical synthesis from Sigma-Aldrich.



Each issue introduces new products and services related to a research area of interest.

To subscribe to ChemFiles or download the corresponding

PDFs, visit our Web site at sigma-aldrich.com/chemfiles.

ChemNews


ChemNews are monthly chemistry e-newsletters featuring the latest products, literature, and news in chemical synthesis from Sigma-Aldrich.

To view or subscribe to ChemNews, visit sigma-aldrich.com/chemnews.

Cheminars™

Learn about the latest innovative chemical synthesis techniques and products with Cheminars, Sigma-Aldrich's Web-based chemistry seminars. These seminars are convenient to navigate, highly interactive, and can be accessed directly via your desktop browser.


To check out the latest Cheminars, please visit sigma-aldrich.com/cheminars.



ChemBlogs

ChemBlogs is designed to be an open forum for the global chemical community, with posts written by Sigma-Aldrich personnel and invited posts from leaders in academic and industry. We will highlight innovative and exciting new developments in catalysis as found in the broad scientific literature and in the greater free press.

To read Sigma-Aldrich's ChemBlogs, please visit chemblogs.com.



ALDRICH®

P.O. Box 14508
St. Louis, MO 63178
USA

PRESORTED
STANDARD
U.S. POSTAGE
PAID
SIGMA-ALDRICH
CORPORATION

Return Service Requested

CELEBRATING 40 YEARS OF THE ALDRICHIMICA ACTA

Aldrichimica ACTA

VOL. 40, NO. 1 • 2007



sigma-aldrich.com

Catalytic Azide–Alkyne Cycloaddition
Selecting and Driving Monolayer Structures


 **ALDRICH**[®]
Advancing Science

New Products from Aldrich R&D

Aldrich Is Pleased to Offer the Latest Asymmetric Allylation Reagents

Stable *Ipc*₂B(allyl) Solutions

The asymmetric allylboration of aldehydes is an extremely important method for the preparation of homoallylic alcohols, as demonstrated in numerous complex natural product syntheses.¹ Although a variety of reagents have been developed to execute this reaction, the most broadly adopted have been Brown's α -pinene-derived (-) and (+)-*B*-allyl-diisopinocampheylboranes.² Typically, the allylboration reagents are generated in situ and used immediately



R	Enantiomeric excess, ee (%)		
	Pentane	Dioxane	Water
Ph	94*	92	90
2-BrC ₆ H ₄	92	90	89
Et	92	90	

*96% ee at 78 °C

at -78 to -100 °C. We recently developed new, salt-free, *Ipc*₂B(allyl) reagents as stable solutions in pentane or dioxane.³ Under refrigerator storage, no appreciable decreases in selectivity are observed after several months. They exhibit superb reactivity, and reactions can be performed at 5 °C, even in water.⁴ We are excited to introduce these convenient allylboration reagents to the research community

(1) (a) Fürstner, A. et al. *Angew. Chem., Int. Ed.* **2006**, *45*, 5506. (b) Smith, A. B., II et al. *J. Am. Chem. Soc.* **2003**, *125*, 350. (2) For a review of allylboration reagents, see Ramachandran, P. V. *Aldrichimica Acta* **2002**, *35*, 23. (3) Patent pending. (4) Josyula, K. V. B. et al., 2006 National ACS meeting, San Francisco, MEDI 545

(+)-*Ipc*₂B(allyl) solution, 1 M in pentane ((+)-*B*-Allyl-diisopinocampheylborane)

678503	5 mL	\$36.50
[106356-53-0]	25 mL	128.00
C ₂₃ H ₃₉ B		
FW: 326.37		



(-)-*Ipc*₂B(allyl) solution, 1 M in pentane ((-)-*B*-Allyl-diisopinocampheylborane)

678538	5 mL	\$41.60
[85116-38-7]	25 mL	145.50
C ₂₃ H ₃₉ B		
FW: 326.37		



(+)-*Ipc*₂B(allyl) solution, 1 M in dioxane

678511	5 mL	\$57.30
[106356-53-0]	25 mL	199.50
C ₂₃ H ₃₉ B		
FW: 326.37		





(-)-*Ipc*₂B(allyl) solution, 1 M in dioxane

678546	5 mL	\$55.00
[85116-38-7]	25 mL	192.50
C ₂₃ H ₃₉ B		
FW: 326.37		



Soderquist Borabicyclodecanes


Treatment of these air-stable crystalline pseudoephedrine (PE) borinic acid complexes with the appropriate Grignard reagent is an excellent method for accessing non-racemic homoallylic or homopropargylic alcohols and α -allyl carbinols via asymmetric allyl-, propargyl-, and allenylboration reactions.¹⁻³ Allylboration reactions occur rapidly (<3 h), can be performed in the presence of Mg salts, and the degree of asymmetric induction is only minimally affected by temperature. Moreover, the reagent is easily recycled in good yield. The reagents also perform well in asymmetric crotylboration,¹ allylboration of imines,⁴ and the alkynylboration of *N*-acylimines to give chiral *N*-propargylamides.⁵



(1) Burgos, C. H. et al. *J. Am. Chem. Soc.* **2005**, *127*, 8044. (2) Lai, C.; Soderquist, J. A. *Org. Lett.* **2005**, *7*, 799. (3) Hernandez, E.; Soderquist, J. A. *Org. Lett.* **2005**, *7*, 5397. (4) Hernandez, E. et al. *Pure Appl. Chem.* **2006**, *78*, 1389. (5) Gonzales, A. Z. et al. *Org. Lett.* **2006**, *8*, 3331


B-Methoxy-10-trimethylsilyl-9-borabicyclo[3.3.2]decane, 95%

676667	5 g	\$84.00
C ₁₃ H ₂₇ BOSi	25 g	336.50
FW: 238.25		




(+)-9-(1*S*,2*S*-Pseudoephedrinyl)-(10*R*)-(trimethylsilyl)-9-borabicyclo[3.3.2]decane

676675	1 g	\$48.00
[848618-13-3]	5 g	168.50
C ₂₂ H ₃₈ BNOSi		
FW: 371.44		



(-)-9-(1*R*,2*R*-Pseudoephedrinyl)-(10*S*)-(trimethylsilyl)-9-borabicyclo[3.3.2]decane

676683	1 g	\$55.70
C ₂₂ H ₃₈ BNOSi	5 g	195.00
FW: 371.44		



Aldrichimica ACTA

VOL. 40, NO. 1 • 2007

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer Inquiries 800-325-3010
Technical Service 800-231-8327
SAFC™ 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 414-438-3850
24-Hour Emergency 414-438-3850
Web Site sigma-aldrich.com
Email aldrich@sial.com

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201, USA

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 355
Milwaukee, WI 53201-9358
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the inside back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All Prices Subject to Change Without Notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2007 Sigma-Aldrich Co.

"PLEASE RATHER US,"

Christopher P. Johnson, III, of the process development department at Boehringer Ingelheim Chemicals, Inc., kindly suggested that we offer boron trichloride as a solution in toluene, in addition to our other solutions of this reagent. Boron trichloride is a versatile Lewis acid and has been used extensively in the selective cleavage of alkyl aryl ethers¹ as well as a variety of C-C and C-O bond-forming reactions.²

(1) Bhatt, M. V.; Kulkarni, S. U. *Synthesis* **1983**, 249. (2) For recent examples, see: (a) Zhang, L.; Zhang, J. Y. *J. Comb. Chem.* **2006**, *8*, 361. (b) Kabalka, G. W.; Yao, M.-L.; Borella, S. J. *Am. Chem. Soc.* **2006**, *128*, 11320. (c) Bellur, E.; Langer, P. J. *Org. Chem.* **2005**, *70*, 3819



678732 Boron trichloride solution, 1.0 M in toluene 100 mL \$30.75
800 mL 138.25

Naturally, we made this useful reagent. It was no bother at all, just a pleasure to be able to help

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the inside back cover.

TABLE OF CONTENTS

What a Journey It Has Been! 3 <i>Sharbil J. Firsan, Sigma-Aldrich Corporation</i>
Catalytic Azide-Alkyne Cycloaddition: Reactivity and Applications 7 <i>Peng Wu and Valery V. Fokin, * The Scripps Research Institute</i>
Selecting and Driving Monolayer Structures through Tailored Intermolecular Interactions 21 <i>Thomas J. Mullen, Arrelaine A. Dameron, Anne M. Andrews, and Paul S. Weiss, * The Pennsylvania State University</i>

ABOUT THE COVER

Winter in the Country (oil on canvas, 66.7 × 92.1 cm) was painted by the American landscape artist, George Henry Durrie (1820–1863), in Connecticut in the winter of 1859. Durrie, the son of a stationer, married a choirmaster's daughter and lived in New Haven, Connecticut, most of his short life. Although he studied painting with a portraitist, landscapes, which first appeared as backgrounds for his portraits, became his primary focus. Durrie did many snow scenes and advertised them in New Haven newspapers, noting, "no collection ... is complete without ... [a] winter scene." After his death, he was immortalized when the lithographic firm of Currier and Ives successfully reproduced ten of his scenes.



© Board of Trustees, National Gallery of Art, Washington

In *Winter in the Country*, Durrie illustrates a number of his most famous compositional conventions. Snow blankets the cold, shadowy landscape, while the ominous sky, reminiscent of Dutch landscapes, covers most of the background. The composition is at eye level, balanced by the centered barn and complementary dark leaf-less trees. Farm workers and animals in the barnyard bring a sense of scale and reality and suggest a waning day. The viewer is invited to explore this everyday 19th-century New England winter scene, pleased to find curling smoke from the farmhouse in the left middle ground, signifying a welcome and warm refuge from the cold, hard day.


This painting was acquired by the National Gallery of Art, Washington, DC, through the Avalon Fund.

New Products for Cross-Coupling

Highly Robust Oxime Palladacycles

The air- and water-stable palladium complexes shown here are extremely active catalysts for C–C and C–X bond-forming processes, ranging from the Heck to the Sonogashira to the Suzuki reaction.¹ These oxime palladacycles represent novel and powerful tools to mediate an array of industrially useful transformations. Furthermore, the versatility of this family of Pd compounds and their electronic and steric features should ensure their wide application in research and manufacturing.


(1) (a) Botella, L.; Najera, C. *J. Org. Chem.* **2005**, *70*, 4360. (b) Alonso, D. A. et al. *Org. Lett.* **2000**, *2*, 1823.



Hindered Ferrocenyl(Dialkyl)phosphine Ligand (Q-Phos)

Developed by the Hartwig group, Q-Phos is a sterically hindered ferrocenylphosphine ligand that has demonstrated broad applicability in a variety of Pd-catalyzed C–C, C–N, and C–O bond-forming reactions including: amination and etherification of aryl chlorides, Suzuki coupling and, most recently, the α arylation of zinc amide enolates.¹


(1) (a) Hama, T. et al. *J. Am. Chem. Soc.* **2006**, *128*, 4976. (b) Kataoka, N. et al. *J. Org. Chem.* **2002**, *67*, 5553



Bulky Alkylphosphine Ligand for Amination Reactions

Shaughnessy and co-workers have utilized this neopentyl-substituted phosphine with a palladium source to form a highly active catalyst for the amination of a diverse set of aryl halides.¹ The increased steric bulk of this ligand, when compared to tri(*tert*-butyl)phosphine, enhances the effectiveness of this system and allows for the formation of the tetra-*ortho*-substituted product shown here.

(1) Hill, L. L. et al. *J. Org. Chem.* **2006**, *71*, 5117.



For more information on these and other new products from Sigma-Aldrich, visit sigma-aldrich.com/synthesis.

What a Journey It Has Been!





Sharbil J. Firsan
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA
Email: sfirsan@sial.com

This January, the *Aldrichimica Acta* begins its 40th year of uninterrupted publication. Reaching this milestone is a genuine measure of the high regard in which chemists worldwide have held it for four decades. Their enduring interest has been sustained by the high-quality scientific information they find within its pages—information contributed by their peers from around the world. Without the loyalty and commitment of these two pillars of the *Acta*, its authors and readers, I would not be sharing this significant anniversary with you today. Thank you both! Your interest and support over the years have fostered stimulating scientific exchanges, provided many answers, and inspired further quests. I am enthusiastic about the prospect of continuing this journey with a new generation of *Aldrichimica Acta* authors and readers.

While some readers may be quite familiar with the *Acta*, others have only recently been introduced to it. As we celebrate reaching this milestone, and in the interest of these readers, let me briefly reflect on the past four decades of the *Acta*.

The *Aldrichimica Acta* publishes review articles on innovative chemistry research that are written by leading experts. It is an ideal vehicle for disseminating the accumulated scientific knowledge, research findings, and insights of these experts to chemistry practitioners in academic, industrial, governmental, and nonprofit laboratories around the globe. In this regard, it serves an important function for both our readers and our company. The key to its popularity and longevity has been its unique format and look, which successfully combine high-quality scientific and technical information with an eye-catching cover featuring a beautiful painting and accompanying commentary, and with the latest useful product offerings from the publisher, Sigma-Aldrich Corporation. All of this at no cost to over 130,000 readers worldwide! While its appearance has undergone several modifications over the years, the premise the *Acta* was based on has not changed.

Aldrichimica



Now for a bit of history. The *Aldrichimica Acta* evolved from the *Kardindex Sheets* that Aldrich used to distribute to its best customers in the 1960s. Dr. Alfred R. Bader, the President of Aldrich Chemical Company at the time, came up with the name, in analogy to the names of such well-known chemistry journals as *Helvetica Chimica Acta* and *Inorganica Chimica Acta*. The name, of course, is a great choice because it combines in its first half the name of the company (“Aldri” from Aldrich) with the type of business it is in (“chimica” for chemical). The second half of the name, *Acta*, is derived from the Latin “*Ācta*”, which means “Proceedings”, and was intended to convey the idea that it is a record of scientific proceedings both at Aldrich and the chemistry community at large. The name *Aldrichimica Acta* has been in use since 1967, and is currently highly recognized worldwide.

The *Aldrichimica Acta* is treated on a par with other popular chemistry journals and its content is covered by such prominent chemical-abstract publishers as Chemical Abstracts Service[®] and Thomson ISI[®]. Collections of *Acta* issues are in practically every academic and industrial science library in the U.S. and many others abroad. In the field of organic chemistry, the Science Edition of Thomson ISI[®]'s *Journal Citation Reports*[®] has ranked the *Acta* #1 five times since 1998, out of 45–58 journals compared by Impact Factor. This is no coincidence, since some of the most prestigious names in chemistry have published in the *Acta* over the years, including six Nobel laureates: Derek H. R. Barton (1969), Vladimir Prelog (1975), Herbert C. Brown (1979), Charles J. Pedersen (1987), George A. Olah (1994), and K. Barry Sharpless (2001). On this occasion, I would like to state unequivocally that we owe them, as well as all the authors that have published in the *Acta* since 1967 (cf. table), a debt of gratitude for their outstanding contributions. Without these contributions, we would not be marking this significant milestone today.

I would be remiss if I didn't also acknowledge my debt and gratitude to all the editors that have preceded me, to the dedicated staff past and present, and to the members of Sigma-Aldrich management, past and present, who have believed in the vital role that the *Acta* plays and have shown unwavering support for it through the ups and downs of the often-turbulent world of modern business. I would also like to thank Sean Battles of Sigma-Aldrich for his stimulating questions and comments relating to this write-up.

Mindful of the great legacy of the past 40 years, we strive everyday not only to keep the *Acta* changing with the times, but also to maintain and enhance its appeal, relevance, and usefulness to our readers and to our company.

Here is to another 40 magnificent years
of the *Aldrichimica Acta*!

All Contributors to the *Aldrichimica Acta* from 1967 to the Present*

Professor Scott E. Denmark

University of Illinois

Abdulla, Riaz F.
 Abeles, Robert H.
 Abell, Andrew D.
 Adams, Richard D.
 Ager, David J.
 Alper, Howard
 Amery, Gordon
 Andrews, Anne M.
 Babu, G.
 Bader, Alfred R.
 Baker, Gregory R.
 Balakumar, R.
 Balding, Peter
 Barluenga, Jose
 Barrientos-Astigarraga, Rafael E.
 Barton, Derek H. R.
 Bashir-Hashemi, A.
 Behrens, Carl H.
 Belyakov, Sergei A.
 Beschke, Helmut
 Biel, John H.
 Bilodeau, M. T.
 Black, T. Howard
 Block, Eric
 Bosnich, Brice
 Bowman, Sheryl A.
 Bradley, M.
 Brewster, James H.
 Bristow, Garry S.
 Brown, Herbert C.
 Brown, Tom
 Bruice, Thomas C.
 Bryce, Martin R.
 Buth, William F.
 Butler, Douglas N.
 Campbell, James B., Jr.
 Campbell, John R.
 Cardillo, Giuliana

Chignell, Colin F.
 Cho, Byung Tae
 Choy, William
 Clark, Anthony M.
 Clive, Derrick L. J.
 Comasseto, Joao V.
 Cook, James M.
 Cori, Tom
 Cornelis, Andre
 Cornforth, John Warcup
 Coward, James K.
 Cox, James R., Jr.
 Crich, David
 Crockett, Geoffrey C.
 Crosby, Guy A.
 Cundy, Darren J.
 Curran, Dennis P.
 Daly, John
 Dameron, Arrelaine A.
 Danishefsky, Samuel J.
 Daub, Guido H.
 Daub, G. William
 Davenport, Derek A.
 Davies, Huw M. L.
 Davies, Stephen G.
 Davis, J. David
 Dawson, D. J.
 De Keukeleire, Denis
 Denmark, Scott E.
 Dervan, Peter B.
 Deslongchamps, Pierre
 Dhar, Raj K.
 Dolphin, David
 Donaldson, William A.
 Dondoni, Alessandro
 Doyle, Gerald
 Doyle, Michael P.
 Durst, H. Dupont



Professor Gary A. Molander

University of Pennsylvania

Emanuel, Ralph N.
 Epton, Roger
 Escalante, Jaime
 Evans, C. Anderson
 Evans, David A.
 Fergason, James L.
 Ferosie, Ishak
 Figueroa, Ruth
 Firsan, Sharbil J.
 Flynn, Daniel L.
 Fokin, Valery V.
 Friess, S. L.
 Frost, John W.
 Fry, Albert J.
 Fu, Xiaoyong
 Fukuyama, Tohru
 Furlan, Ricardo L. E.
 Gallou, Isabelle
 Gentilucci, Luca
 Gibson, David T.
 Gilchrist, Thomas L.
 Gill, G. B.
 Gill, Steven K.
 Gladysz, J. A.
 Gokel, George W.
 Gonzalez, David
 Goralski, Christian T.
 Gray, Harry B.
 Green, Floyd J.
 Greenlee, Kenneth W.
 Grieco, Paul A.
 Griffith, William P.
 Griffiths, David W.
 Gupton, John T.
 Gutsche, C. David
 Ha, Hyun-Joon
 Hall, H. K., Jr.
 Han, Zhengxu
 Hanessian, Stephen
 Hanson, Paul R.
 Harned, Andrew M.
 Harvey, David R.
 Harwood, Lawrence M.
 Hase, Tapio A.
 Hayes, Brittany L.
 Heaney, Harry
 Heathcock, Clayton H.
 Herpel, Russell H.
 Hewitt, Christopher D.
 Higbee, W. Edward
 Hoffman, Frances
 Hopps, Harvey B.



Professor Dennis P. Curran

University of Pittsburgh

Hormann, Robert E.
 Hossain, M. Mahmud
 Hu, Qian
 Huang, Zhihong
 Hudlicky, Tomas
 Humphries, R. E.
 Humphries, Ken C.
 Ingold, Keith U.
 Ireland, Robert E.
 James, Dustin K.
 Janda, Kim D.
 Johnson, Carl R.
 Jones, J. Bryan
 Jones, R. Alan
 Juaristi, Eusebio
 Kabalka, George W.
 Kanai, Motomu
 Kantchev, Eric Assen B.
 Kanth, Josyula V. B.
 Katritzky, Alan R.
 Kiesow, L. A.
 Kime, Katherine A.
 Kisanga, Philip B.
 Kishi, Yoshito
 Klitsner, Marvin E.
 Klundt, Irwin L.
 Knowles, Jeremy R.
 Koganti, Vara Prasad
 Koppel, Henry
 Kornblum, Nathan
 Koser, Gerald F.
 Koskimies, Jorma K.
 Kramer, Gary W.
 Krishnamurthy, Dhileepkumar
 Krishnamurthy, S.
 Kulkarni, Yashwant S.
 Lane, Clinton F.
 Laszlo, Pierre
 Lauret, Christelle

Sigma-Aldrich Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA S.A.
Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

SIGMA-ALDRICH PTY LTD
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA.
Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

China

SIGMA-ALDRICH (SHANGHAI)
TRADING CO. LTD
Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

SIGMA-ALDRICH S.R.O.
Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: (+45) 43 56 59 10
Fax: (+45) 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND OY
Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 00
Fax: (+33) 474 95 68 08

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Greece

SIGMA-ALDRICH (O.M.) LTD.
Tel: (+30) 210 994 8010
Fax: (+30) 210 994 3831

Hungary

SIGMA-ALDRICH Kft
Ingyenes zold telefon: 06 80 355 355
Ingyenes zold fax: 06 80 344 344
Tel: (+36) 1 235 9055
Fax: (+36) 1 235 9050

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: (+91) 80 6621 9600
New Delhi: (+91) 11 4165 4255
Mumbai: (+91) 22 2570 2364
Hyderabad: (+91) 40 6684 5488
Fax
Bangalore: (+91) 80 6621 9650
New Delhi: (+91) 11 4165 4266
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 6684 5466

Ireland

SIGMA-ALDRICH IRELAND LTD
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 1 404 1900
Fax: (+353) 1 404 1910

Israel

SIGMA-ALDRICH ISRAEL LTD
Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

SIGMA-ALDRICH S.r.l.
Numero Verde: 800 827018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

SIGMA-ALDRICH JAPAN K.K.
Tokyo Tel: (+81) 3 5796 7300
Tokyo Fax: (+81) 3 5796 7315

Korea

SIGMA-ALDRICH KOREA
Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: 52 722 276 1600
Fax: 52 722 276 1601

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

SIGMA-ALDRICH NEW ZEALAND LTD
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

SIGMA-ALDRICH NORWAY AS
Tel: (+47) 23 17 60 60
Fax: (+47) 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

SIGMA-ALDRICH RUS, LLC
Tel: +7 (495) 621 6037
Fax: +7 (495) 621 5923

Singapore

SIGMA-ALDRICH PTE. LTD
Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD.
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA, S.A
Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD
Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313
SAFC (UK) Free Tel: 0800 71 71 17

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Call Collect: (+1) 314 771 5750
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Internet

sigma-aldrich.com

Sigma-Aldrich Career Opportunities

As a leading Life Science and High Technology company, we are always looking for talented individuals to join our team. At Sigma-Aldrich we value the contributions of our employees, and recognize the impact they have on our success. We strive to foster creativity and innovation, and encourage professional development.

Our biochemical and organic chemical products and kits are used in scientific and genomic research, biotechnology, pharmaceutical development, the diagnosis of disease, and as key components in pharmaceutical and other high technology manufacturing. We have customers in life science companies, university and government institutions, hospitals, and in industry

UNLEASH YOUR TALENTS

Learn more about our career opportunities by visiting
our award-winning Web site at sigma-aldrich.com/careers.

sigma-aldrich.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
SIGMA-ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI 63178 • USA

Sigma-Aldrich Corporation is an equal opportunity employer



SIGMA-ALDRICH

It's Here!

NEW Aldrich Handbook



- ✓ 35,000 Chemical Listings
- ✓ 10,000 Structures
- ✓ 2,000 NEW Products
- ✓ Enhanced Application Index

Order your FREE copy today.

Contact your local office,

or in the US call **1-800-325-3010**

or visit sigma-aldrich.com/handbook23

RESEARCH IS EASIER WHEN YOU USE THE RIGHT TOOLS!

Aldrich Handbook of Fine Chemicals

- ✓ 35,000 Chemical Listings
- ✓ 10,000 Structures
- ✓ 2,000 NEW Products
- ✓ 6,100 Citations
- ✓ 3,000 Application Notes
- ✓ Enhanced Application Index
- ✓ Advanced Product Tables

Now in a convenient Dual-Volume Set that includes:

Sigma-Aldrich Labware Catalog

- ✓ 17,000 Complementary Products
- ✓ Convenient Procurement
- ✓ Customer Preferred Brands
- ✓ Special Offers
- ✓ Unprecedented Savings

Set includes the
**NEW Sigma-Aldrich
Labware Catalog**



sigma-aldrich.com

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE
SIGMA-ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI 63178 • USA



SIGMA-ALDRICH

...taloo Cover: The Alchemist (ca. 1937) by Newell Conner's Wyeth. Courtesy of the Chemical Heritage Foundation Image Archive.



ALDRICH®

P.O. Box 14508
St. Louis, MO 63178
USA

Return Service Requested

PRESORTED
STANDARD
U.S. POSTAGE
PAID
SIGMA-ALDRICH
CORPORATION

11-14508-01
011

ALDRICH CONGRATULATES THE 2007 ACS AWARD WINNERS

Aldrichimica **ACTA**

VOL. 40, NO. 2 • 2007



sigma-aldrich.com

Recent Advances in Intermolecular
Direct Arylation Reactions

Evolution and Applications of Second-Generation
Ruthenium Olefin Metathesis Catalysts


 **ALDRICH**[®]
Advancing Science

New Products from Aldrich R&D

Available in limited quantities. Contact your local distributor for details.

Optimization of Mitsunobu Reagents

The Mitsunobu reaction is one of the most extensively used coupling reactions in organic synthesis and typically employs azodicarboxylate reagents such as DEAD or DIAD. However, these reagents have drawbacks such as low room-temperature stability and difficulty in removing the hydrazine byproducts. Professor Bruce Lipshutz and co-workers have developed an attractive alternative to the existing reagents: di(4-chlorobenzyl) azodicarboxylate (DCAD). DCAD is a stable solid that




DCAD: 92%
DEAD: 94%
DIAD: 89%

has an activity comparable to those of DEAD and DIAD in typical Mitsunobu reactions such as substitutions, esterifications, and etherifications. However, unlike the standard reagents, the hydrazine byproduct can be removed by simple precipitation directly from the reaction mixture, and is easily recycled in high yield to regenerate DCAD.

Lipshutz, B. H. et al. *Org. Lett.* **2006**, *8*, 5069.


680850

$\text{C}_{16}\text{H}_{12}\text{Cl}_2\text{N}_2\text{O}_4$
FW: 367.18



1 g \$24.00
10 g 99.50


Because of the ease of preparation of *meso*-diols, synthetic methods that can desymmetrize these substrates are critically important. Professors Marc Snapper and Amir Hoveyda at Boston College have reported the first practical enantioselective silylation of *meso*-1,2- and 1,3-diols relying on an amino acid derived organocatalyst. The reactions do not require the rigorous exclusion of air or moisture, and the catalyst can be nearly quantitatively recovered by an aqueous wash. This catalyst greatly increases the efficiency with which optically enriched molecules can be prepared.



Zhao, Y. et al. *Nature* **2006**, *443*, 67

680826


$\text{C}_{17}\text{H}_{21}\text{N}_3\text{O}$
FW: 308.46



1 g \$90.00

Trichloroacetimidate Reagents


Trichloroacetimidates are useful reagents for protection of alcohols as their allyl and benzyl ethers. We are delighted to offer two new reagents, allyl 2,2,2-trichloroacetimidate and 4-methoxybenzyl 2,2,2-trichloroacetimidate, that have been extensively employed in organic synthesis. These reagents are particularly attractive in applications where base-sensitive functional groups are present that would not tolerate the standard alkoxide alkylation method of alcohol protection.



Clark, J. S. et al. *Tetrahedron* **2006**, *62*, 73.

678414


$\text{C}_5\text{H}_6\text{Cl}_3\text{NO}$
FW: 202.47



5 g \$75.00


679585

$\text{C}_{10}\text{H}_{10}\text{Cl}_3\text{NO}_2$
FW: 282.55



5 g \$50.00
25 g 198.00

Cyclopropyl groups are found in a variety of natural products and are increasingly incorporated into pharmaceuticals such as the broad-spectrum antibiotic ciprofloxacin. Both the Charette¹ and Deng² groups have reported success in the cross-coupling of potassium cyclopropyltrifluoroborates with aryl bromides in the presence of common palladium catalysts. The trifluoroborate salts exhibit enhanced stability and more certain stoichiometry relative to their boronic acid counterparts. However, like boronic acids, post-reaction byproducts are easily removed. We are pleased to add this useful reagent to our ever-growing arsenal of organoboron compounds.



(1) Charette, A. B. et al. *Synlett* **2005**, *11*, 1779. (2) Fang, G.-H. et al. *Org. Lett.* **2004**, *6*, 357

662984

$\text{C}_3\text{H}_5\text{BF}_3\text{K}$
FW: 147.98



1 g \$54.00
5 g 216.00

Aldrichimica ACTA

VOL. 40, NO. 2 • 2007

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer Inquiries	800-325-3010
Technical Service	800-231-8327
SAFC™	800-244-1173
Custom Synthesis	800-244-1173
Flavors & Fragrances	800-227-4563
International	414-438-3850
24-Hour Emergency	414-438-3850
Web Site	sigma-aldrich.com
Email	aldrich@sial.com

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201, USA

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 355
Milwaukee, WI 53201-9358
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the inside back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com.


Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2007 Sigma-Aldrich Co.

Professor Gregory B. Dudley of Florida State University kindly suggested that we make 2-benzyloxy-1-methylpyridinium triflate. This crystalline, neutral, and stable organic salt is an excellent reagent for the protection of an alcohol as a benzyl ether under mild conditions. Reaction with this reagent can be performed under near-neutral pH, unlike other benzylation protocols, which require strongly acidic or basic reaction media (e.g., the use of benzyl trichloroacetimidate or benzyl halides).¹

(1) Poon, K. W. C.; Dudley, G. B. *J. Org. Chem.* **2006**, *71*, 3923. (2) Poon, K. W. C.; House, S. E.; Dudley, G. B. *Synlett* **2005**, 3142



679674 2-Benzyloxy-1-methylpyridinium triflate 1 g \$75.00
5 g 250.00

Naturally, we made this useful reagent. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the inside back cover.

Recent Advances in Intermolecular Direct Arylation Reactions 35

Louis-Charles Campeau, David R. Stuart, and Keith Fagnou,* University of Ottawa

Evolution and Applications of Second-Generation Ruthenium Olefin Metathesis Catalysts 45

Yann Schrodi* and Richard L. Pederson, Matera, Inc

Panoramic Landscape with Hunters (oil on canvas, 105 x 135 cm) was painted in the mid-1660s by Philips Koninck (1619–1688), one of the great Baroque landscape artists of the Golden Age of Dutch Art (ca. 1600–1680). Although a contemporary of Rembrandt, Koninck is not believed to have studied with him. However, Koninck knew the master and some of his pupils and was certainly familiar with Rembrandt's paintings, which had some influence on him.





Photograph © Alfrede Bader

This painting illustrates Koninck's method of bringing together details of real-life scenes to create fictional but convincing sweeping landscapes featuring streams, fields, abundant flora, and rural dwellings. The translucent colors of the sky, the receding diagonal lines, and the horizontal striations denoting successive planes that recede into the distance add to the great allure of this landscape.

This painting is in the private collection of Isabel and Alfred Bader. Dr. Bader is a perennial "chemist collector" and a former Aldrich and Sigma-Aldrich president.


Reagents for Direct Arylation

Pd-catalyzed cross-coupling of organometallic nucleophiles with aryl halides has become the most commonly used method for biaryl synthesis. However, the range of biaryls that can be prepared is limited to those organometallic reagents that are commercially available or easily made. Nitrogen-containing heterocyclic organometallic reagents are often difficult to prepare and success of their coupling reactions can be sporadic. Professor Keith Fagnou and co-workers at the University of Ottawa have developed a novel method for biaryl synthesis by the direct arylation of heterocyclic *N*-oxides.¹⁻³ Yields are typically very good, and the oxide residue is easily reduced to give the free azine or diazine.



References

(1) Leclerc, J.-P.; Fagnou, K. *Angew. Chem., Int. Ed.* **2006**, *45*, 7781. (2) Campeau, L.-C.; Rousseau, S.; Fagnou, K. *J. Am. Chem. Soc.* **2005**, *127*, 18020. (3) Campeau, L.-C.; Stuart, D.R.; Fagnou, K. *Aldrichimica Acta* **2007**, *40*, 35.



For more information on these and other new products from Sigma-Aldrich, visit sigmaaldrich.com/go/chem.

Sigma-Aldrich Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA S.A.
Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

SIGMA-ALDRICH PTY LTD
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA
Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

China

SIGMA-ALDRICH (SHANGHAI)
TRADING CO. LTD.
Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

SIGMA-ALDRICH S.R.O.
Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: (+45) 43 56 59 10
Fax: (+45) 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND OY
Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 00
Fax: (+33) 474 95 68 08

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Greece

SIGMA-ALDRICH (O.M.) LTD.
Tel: (+30) 210 994 8010
Fax: (+30) 210 994 3831

Hungary

SIGMA-ALDRICH Kft
Ingyenes zöld telefon: 06 80 355 355
Ingyenes zöld fax: 06 80 344 344
Tel: (+36) 1 235 9055
Fax: (+36) 1 235 9050

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: (+91) 80 6621 9600
New Delhi: (+91) 11 4165 4255
Mumbai: (+91) 22 2570 2364
Hyderabad: (+91) 40 6684 5488
Fax
Bangalore: (+91) 80 6621 9650
New Delhi: (+91) 11 4165 4266
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 6684 5466

Ireland

SIGMA-ALDRICH IRELAND LTD
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 1 404 1900
Fax: (+353) 1 404 1910

Israel

SIGMA-ALDRICH ISRAEL LTD
Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

SIGMA-ALDRICH S.r.l.
Numero Verde: 800 827018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

SIGMA-ALDRICH JAPAN K.K.
Tokyo Tel: (+81) 3 5796 7300
Tokyo Fax: (+81) 3 5796 7315

Korea

SIGMA-ALDRICH KOREA
Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: 52 722 276 1600
Fax: 52 722 276 1601

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

SIGMA-ALDRICH NEW ZEALAND LTD
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

SIGMA-ALDRICH NORWAY AS
Tel: (+47) 23 17 60 60
Fax: (+47) 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

SIGMA-ALDRICH RUS, LLC
Tel: +7 (495) 621 6037
Fax: +7 (495) 621 5923

Singapore

SIGMA-ALDRICH PTE. LTD
Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD
Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313
SAFC (UK) Free Tel: 0800 71 71 17

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Call Collect: (+1) 314 771 5750
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Internet

sigma-aldrich.com

Reserve your copy of the 2007-2008 SUPELCO® Catalog

The 2007-2008 SUPELCO® Catalog features:

- Thousands of "Chromatography Products for Analysis & Purification"
- Hundreds of NEW applications

The catalog is scheduled to ship in April 2007.

Visit sigma-aldrich.com/supelco-catalog to reserve your copy today.




Looking for More Free Updates in Chemistry from Sigma-Aldrich?

ChemFiles

ChemFiles are monthly new-product application guides in chemical synthesis from Sigma-Aldrich. Each issue introduces new products and services related to a research area of interest.

To subscribe to *ChemFiles* or download the corresponding PDFs, visit our Web site at sigma-aldrich.com/chemfiles.




Cheminars™

Learn about the latest innovative chemical synthesis technologies and products with *Cheminars™*, Sigma-Aldrich's

Web-based chemistry seminars. These seminars are convenient to navigate, highly interactive, and can be accessed directly via your desktop browser.


To check out the latest *Cheminars™*, please visit sigma-aldrich.com/cheminars.



ChemNews

ChemNews are monthly chemistry e-newsletters featuring the latest products, literature, and news in chemical synthesis from Sigma-Aldrich.

To view or subscribe to *ChemNews*, visit sigma-aldrich.com/chemnews.




ChemBlogs

ChemBlogs is designed to be an open forum for the global chemical community, with posts written by Sigma-Aldrich personnel and invited posts from leaders in academia and industry. We will highlight

innovative and exciting new developments in catalysis as found in the broad scientific literature and in the greater free press.

To read Sigma-Aldrich's *ChemBlogs*, please visit chemblogs.com.



Cheminars is a trademark of Sigma-Aldrich Biotechnology, L.P.



ALDRICH®

P.O. Box 14508
St. Louis, MO 63178
USA

PRESORTED
STANDARD
U.S. POSTAGE
PAID
SIGMA-ALDRICH
CORPORATION

Return Service Requested

THE GROWING IMPACT OF ASYMMETRIC CATALYSIS

Aldrichimica ACTA

VOL. 40, NO. 3 • 2007



Palladium-Catalyzed Dynamic Kinetic Asymmetric
Allylic Alkylation with the DPPBA Ligands

Development and Applications of C_2 -Symmetric,
Chiral, Phase-Transfer Catalysts


SIGMA-ALDRICH™

New Products from Aldrich R&D

Sigma-Aldrich Is Pleased to Offer Cutting-Edge Tools for Organic Synthesis

DABAL-Me₃

Trimethylaluminum is a versatile methylation reagent in organic synthesis. However, because of its pyrophoric nature, it cannot be handled in open air. Developed by the Woodward group (University of Nottingham, U.K.), DABAL-Me₃ is a free-flowing solid adduct of trimethylaluminum and DABCO® that can be manipulated without the need for an inert atmosphere.¹ This bench-stable reagent has been employed in numerous reactions including methylations of aldehydes and imines,^{1,2} methylation of aryl and vinyl halides,³ conjugate additions to enones,⁴ and amide-bond formation.⁵ In the presence of the appropriate chiral ligand and catalyst, many of these reactions can be performed asymmetrically



(1) Woodward, S. *Synlett* **2007**, 1490. (2) Mata, Y. et al. *J. Org. Chem.* **2006**, *71*, 8159. (3) Cooper, T. et al. *Adv. Synth. Catal.* **2006**, *348*, 686. (4) Alexakis, A. et al. *Chem. Commun.* **2005**, 2843. (5) Novak, A. et al. *Tetrahedron Lett.* **2006**, *47*, 5767


Bis(trimethylaluminum)-1,4-diazabicyclo[2.2.2]octane adduct (DABAL-Me₃)

682101	1 g	\$20.90
[137203-34-0]	5 g	71.60
C ₁₄ H ₂₄ Al ₂ N ₂		
FW: 256.34		



White Catalyst for Allylic C-H Oxidation


Professor Christina White's group (University of Illinois) recently reported selective allylic C-H oxidation reactions catalyzed by a Pd(II)-bis-sulfoxide system that furnishes branched allylic esters from α -olefins and carboxylic acids.¹ These reactions can be performed in an inter- or intramolecular fashion, the latter being capable of yielding highly functionalized, large-ring macrolactone products.² Finally, the catalyst system allows for a one-pot sequential allylic oxidation-C-H arylation to afford the *E* arylated allylic ester from the corresponding olefin, carboxylic acid, and arylboronic acid.³



(1) Chen, M. S. et al. *J. Am. Chem. Soc.* **2005**, *127*, 6970. (2) Fraunhofer, K. J. et al. *J. Am. Chem. Soc.* **2006**, *128*, 9032. (3) Delcamp, J. H.; White, M. C. *J. Am. Chem. Soc.* **2006**, *128*, 15076


White Catalyst

684821	250 mg	\$27.50
[858971-43-4]	1 g	90.00
C ₁₈ H ₁₈ O ₂ PdS ₂		
FW: 502.90		



TarB-NO₂ Reducing Reagents


In conjunction with NaBH₄, Singaram's chiral TarB-NO₂ boronic esters rapidly reduce prochiral ketones to optically active secondary alcohols with enantiomeric excesses as high as 99%.¹⁻³ The reagents cleanly reduce aromatic ketones with high enantioselectivity and, in many cases, aliphatic ketones can be reduced with a similar degree of selectivity. Typically, TarB-NO₂ reagents perform as well as, or better than, existing hydride asymmetric reduction methods such as those employing DIP-Chloride™ or the CBS reagents.



(1) Kim, J.; Singaram, B. *Tetrahedron Lett.* **2006**, *47*, 3901. (2) Kim, J. et al. *Org. Process Res. Dev.* **2006**, *10*, 949. (3) Cordes, D. B. et al. *Eur. J. Org. Chem.* **2005**, 5289.


3-Nitrophenylboronic acid D-tartaric acid ester, 1 M in THF (D-TarB-NO₂)

682748	5 mL	\$54.30
C ₁₀ H ₈ BNO ₈	25 mL	187.50
FW: 280.98		




3-Nitrophenylboronic acid L-tartaric acid ester, 1 M in THF (L-TarB-NO₂)

682713	5 mL	\$54.30
[467443-01-2]	25 mL	187.50
C ₁₀ H ₈ BNO ₈		
FW: 280.98		



N-tert-Butylbenzenesulfenamide


In the presence of NCS, *N*-tert-butylbenzenesulfenamide catalyzes the selective oxidation of a variety of primary and secondary alcohols to the corresponding aldehydes and ketones in high yield and under mild conditions.^{1,2} The catalytic oxidation tolerates various functional groups including silyl ethers, epoxides, urethanes, esters, and olefins. The reaction is particularly useful for the preparation of labile or easily epimerized aldehydes.



(1) Mukaiyama, T. *Angew. Chem., Int. Ed.* **2004**, *43*, 5590. (2) Matsuo, J.-i. et al. *Tetrahedron* **2003**, *59*, 6739

N-tert-Butylbenzenesulfenamide, 97%

681792	1 g	\$65.90
[19117-31-8]	5 g	263.50
C ₁₀ H ₁₅ NS		
FW: 181.30		



DABCO is a registered trademark of Air Products and Chemicals, Inc. DIP-Chloride is a trademark of Sigma-Aldrich Biotechnology, L.P.

Aldrichimica ACTA

VOL. 40, NO. 3 • 2007

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
 6000 N. Teutonia Ave.
 Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
 FAX 800-325-5052 (USA)
 or 414-438-2199
 Mail P.O. Box 2060
 Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
 Technical Service 800-231-8327
 SAFC™ 800-244-1173
 Custom Synthesis 800-244-1173
 Flavors & Fragrances 800-227-4563
 International 414-438-3850
 24-Hour Emergency 414-438-3850
 Web Site sigma-aldrich.com
 Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
 P.O. Box 355, Milwaukee, WI 53201, USA

Subscriptions

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
 Mail: **Attn: Mailroom**
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 355
Milwaukee, WI 53201-9358
 Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the inside back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2007 Sigma-Aldrich Co.



"PLEASE BOTHER US."

Joe Porwoll

Joe Porwoll, President
 Aldrich Chemical Co., Inc

Professor Carsten Bolm of RWTH Aachen University, kindly suggested that we make 2-(trimethylsilyl)ethanesulfonyl chloride (SES-Cl). This reagent is employed to protect an amine in the form of its sulfonamide. In contrast to the harsh conditions sometimes needed to deprotect tosyl-protected amines, the SES group is readily cleaved under mild conditions using a fluoride ion source, regenerating the parent amine along with volatile byproducts. We have also prepared SES-NH₂, a useful reagent for the introduction of a protected nitrogen atom into a substrate.^{1,2}

(1) Weinreb, S. M. et al. *Tetrahedron Lett.* **1986**, 27, 2099. (2) Ribière, P. et al. *Chem. Rev.* **2006**, 106, 2249



681334	2-(Trimethylsilyl)ethanesulfonyl chloride (SES-Cl)	1 g \$89.40 5 g 330.50
681326	2-(Trimethylsilyl)ethanesulfonamide (SES-NH ₂)	1 g \$95.60

Naturally, we made these useful reagents. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the inside back cover.

TABLE OF CONTENTS

Palladium-Catalyzed Dynamic Kinetic Asymmetric Allylic Alkylation with the DPPBA Ligands	59
<i>Barry M. Trost* and Daniel R. Fandrick, Stanford University</i>	
Development and Applications of C₂-Symmetric, Chiral, Phase-Transfer Catalysts	77
<i>Takashi Ooi and Keiji Maruoka,* Kyoto University</i>	

ABOUT OUR COVER

Oarsmen at Chatou (oil on canvas, 81.2 × 100.2 cm) was painted in 1879 by the French impressionist painter, Pierre Auguste Renoir (1841–1919), on the river Seine west of Paris. His use of light fresh colors in this painting and throughout his career was the result of his love of paintings from the Rococo period and of his training in a porcelain factory as a young man.

Rowing was the foremost attraction at Chatou. The man in this boat—wearing the typical costume of a short jacket and a straw hat—may be the artist's brother, Edmond. The man standing on the bank, similarly attired, is probably the painter Gustave Caillebotte, a devoted rowing enthusiast and a friend of Renoir. The woman is most likely Aline Charigot, who was his favorite model and later became his wife.

The painting captures the brilliance of sun and water, summer and youth. In the water, strong blues and white alternate. Their shimmering intensity is enhanced by the equally strong presence of orange in the boat's reflection and the scarlet accent of Aline's bow. Renoir has put into practice the principle of simultaneous contrast: colors are perceived stronger when juxtaposed with their opposites—orange with blue, for example, or green with red. The silky texture of Renoir's feathery brushstrokes mirrors the languid and leisurely scene.

This painting is a gift of Sam A. Lewisohn to the National Gallery of Art, Washington, DC.




Photograph © Board of Trustees, National Gallery of Art, Washington


Trost Ligands for Asymmetric Allylic Alkylation

Asymmetric allylic alkylation is a versatile catalytic reaction allowing access to a diversity of chiral molecules. This transformation converts both enantiomers of the substrate into the same enantiomer of the product, allowing theoretical yields of 100% of one enantiomer. Professor Trost developed a series of ligands based on diphenylphosphinobenzoic acid (DPPBA) and used them with a variety of palladium complexes for the asymmetric allylic alkylation. These ligands perform with a high degree of enantioselectivity and high yields.

DACH-Phenyl Trost Ligands



Ovaa, H. et al. *Chem. Commun.* **2000**, 1501




692808


(*R,R*)-DACH-Phenyl Trost Ligand
692808

(*S,S*)-DACH-Phenyl Trost Ligand
692794

DACH-Naphthyl Trost Ligands



Trost, B. M.; Schroeder, G. M. *J. Org. Chem.* **2000**, 65, 1569




692778


(*R,R*)-DACH-Naphthyl Trost Ligand
692778

(*S,S*)-DACH-Naphthyl Trost Ligand
692786

DACH-Pyridyl Trost Ligands



Trost, B. M.; Hachiya, I. *J. Am. Chem. Soc.* **1998**, 120, 1104.



692743

(*R,R*)-DACH-Pyridyl Trost Ligand
692751

(*S,S*)-DACH-Pyridyl Trost Ligand
692743

For more information, see Professor Trost's review in this issue.

Sold in collaboration with DowPharmaSM for research purposes only. US Patent 5739396 applies

Sigma-Aldrich Worldwide Locations

- Argentina**
SIGMA-ALDRICH DE ARGENTINA S.A.
Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698
- Australia**
SIGMA-ALDRICH PTY LTD.
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500
- Austria**
SIGMA-ALDRICH HANDELS GmbH
Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20
- Belgium**
SIGMA-ALDRICH NV/SA.
Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11
- Brazil**
SIGMA-ALDRICH BRASIL LTDA.
Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895
- Canada**
SIGMA-ALDRICH CANADA LTD.
Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292
- China**
SIGMA-ALDRICH (SHANGHAI)
TRADING CO. LTD.
Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567
- Czech Republic**
SIGMA-ALDRICH S.R.O.
Tel: (+420) 246 003 200
Fax: (+420) 246 003 291
- Denmark**
SIGMA-ALDRICH DENMARK A/S
Tel: (+45) 43 56 59 10
Fax: (+45) 43 56 59 05
- Finland**
SIGMA-ALDRICH FINLAND OY
Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555
- France**
SIGMA-ALDRICH CHIMIE S.à.r.l.
Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 00
Fax: (+33) 474 95 68 08
- Germany**
SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160
- Greece**
SIGMA-ALDRICH (O.M.) LTD.
Tel: (+30) 210 994 8010
Fax: (+30) 210 994 3831
- Hungary**
SIGMA-ALDRICH Kft
Ingyenes zöld telefon: 06 80 355 355
Ingyenes zöld fax: 06 80 344 344
Tel: (+36) 1 235 9055
Fax: (+36) 1 235 9050
- India**
SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: (+91) 80 6621 9600
New Delhi: (+91) 11 4165 4255
Mumbai: (+91) 22 2570 2364
Hyderabad: (+91) 40 6684 5488
Fax
Bangalore: (+91) 80 6621 9650
New Delhi: (+91) 11 4165 4266
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 6684 5466
- Ireland**
SIGMA-ALDRICH IRELAND LTD
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 1 404 1900
Fax: (+353) 1 404 1910
- Israel**
SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200
- Italy**
SIGMA-ALDRICH S.r.l.
Numero Verde: 800 827018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737
- Japan**
SIGMA-ALDRICH JAPAN K.K.
Tokyo Tel: (+81) 3 5796 7300
Tokyo Fax: (+81) 3 5796 7315
- Korea**
SIGMA-ALDRICH KOREA
Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090
- Malaysia**
SIGMA-ALDRICH (M) SDN. BHD
Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116
- Mexico**
SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: 52 722 276 1600
Fax: 52 722 276 1601
- The Netherlands**
SIGMA-ALDRICH CHEMIE BV
Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421
- New Zealand**
SIGMA-ALDRICH NEW ZEALAND LTD.
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500
- Norway**
SIGMA-ALDRICH NORWAY AS
Tel: (+47) 23 17 60 60
Fax: (+47) 23 17 60 50
- Poland**
SIGMA-ALDRICH Sp. z o.o.
Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20
- Portugal**
SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610
- Russia**
SIGMA-ALDRICH RUS, LLC
Tel: +7 (495) 621 6037
Fax: +7 (495) 621 5923
- Singapore**
SIGMA-ALDRICH PTE. LTD.
Tel: (+65) 6779 1200
Fax: (+65) 6779 1822
- South Africa**
SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119
- Spain**
SIGMA-ALDRICH QUÍMICA, S.A
Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42
- Sweden**
SIGMA-ALDRICH SWEDEN AB
Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243
- Switzerland**
SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815
- United Kingdom**
SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313
SAFC (UK) Free Tel: 0800 71 71 17
- United States**
SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Call Collect: (+1) 314 771 5750
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757
- Internet**
sigma-aldrich.com

Sigma-Aldrich Career Opportunities

As a leading Life Science and High Technology company, we are always looking for talented individuals to join our team. At Sigma-Aldrich we value the contributions of our employees, and recognize the impact they have on our success. We strive to foster creativity and innovation, and encourage professional development.

Our biochemical and organic chemical products and kits are used in scientific and genomic research, biotechnology, pharmaceutical development, the diagnosis of disease, and as key components in pharmaceutical and other high technology manufacturing. We have customers in life science companies, university and government institutions, hospitals, and in industry.

UNLEASH YOUR TALENTS

Learn more about our career opportunities by visiting
our award-winning Web site at sigma-aldrich.com/careers.

Sigma-Aldrich Corporation is an equal opportunity employer

sigma-aldrich.com

SIGMA-ALDRICH™

SciBook™ Select


Best Books for the Best Minds
sigma-aldrich.com/books

The Sigma-Aldrich Library
is your guide to finding new and
best-selling chemistry books.

Whether your interest includes Drug Discovery,
Chemical Synthesis, Materials Science, or a wide range of
other areas of interest, we can help you find the right book.

View table of contents, search, browse, or order
from our entire library at sigma-aldrich.com/books.

Chemical Synthesis



New
for
2007


Z730203 - \$54.90
**Palladium in Heterocyclic Chemistry
A Guide for the Synthetic Chemist,
2nd Edition**
J. J. Li and G. W. Gribble, Elsevier, 2007,
658 pp. Softcover




Z705284 - \$200.00
**Handbook of Chiral Chemicals,
2nd Edition**
D. Ager, Ed., CRC Press, 2006, 664 pp.
Hardcover.




Z706043 - \$295.00
Chiral Analysis
K. W. Busch and M. A. Busch, Eds.,
Elsevier, 2006, 720 pp. Hardcover.



Z730181 - \$115.00
The Pilot Plant Real Book, 2nd Edition
F. X. McConville, FXM Engineering, 2006,
320 pp. Softcover



Z703354 - \$110.00
**Side Reactions in Organic Synthesis:
A Guide to Successful Synthesis Design**
F. Z. Dörwald, Wiley-VCH, 2005, 389 pp. Softcover



Z704113 - \$195.00
**Asymmetric Organocatalysis -
From Biomimetic Concepts to
Applications in Asymmetric Synthesis**
A. Berkessel and H. Gröger, Wiley-VCH, 2005, 454
pp. Hardcover

SciBook Select is a trademark of Sigma-Aldrich Biotechnology, L.P.

sigma-aldrich.com

SIGMA-ALDRICH™



P.O. Box 14508
St. Louis, MO 63178
USA

Change Service Requested

PRESORTED
STANDARD
U.S. POSTAGE
PAID
SIGMA-ALDRICH
CORPORATION


RR
02472-503200
0087

ALLENES AND IONIC LIQUIDS IN ASYMMETRIC SYNTHESIS

Aldrichimica ACTA

VOL. 40, NO. 4 • 2007

 **ALDRICH**
Chemistry



Recent Advances in the Chemistry of Allenes
Chiral Imidazolium Ionic Liquids: Their Synthesis and
Influence on the Outcome of Organic Reactions


SIGMA-ALDRICH™

New Products from Aldrich R&D

Aldrich Is Pleased to Offer Cutting-Edge Reagents for Organic Synthesis

Deoxygenation Reagent

Developed by the Movassaghi group at MIT, *N*-isopropylidene-*N'*-2-nitrobenzenesulfonyl hydrazine (IPNBSh) is a useful reagent for the mild deoxygenation of allylic and propargylic alcohols to give allylically transposed alkenes and allenes, respectively.¹ This reagent exhibits excellent reactivity in difficult reductive fragmentations, as demonstrated in the total syntheses of (-)-acylfulvene and (-)-irofulven.²



(1) Movassaghi, M.; Ahmad, O. K. *J. Org. Chem.* **2007**, *72*, 1838. (2) Movassaghi, M. et al. *Angew. Chem., Int. Ed.* **2006**, *45*, 5859.


N-Isopropylidene-*N'*-2-nitrobenzenesulfonyl hydrazine, 97% IPNBSh

687855

[6655-27-2]

C₉H₁₁N₃O₄S

FW: 257.27




1 g \$52.50

5 g 175.00

Bromination Reagent

Bromodimethylsulfonium bromide (BDMS) is an easy-to-handle and highly effective bromination reagent as well as a catalyst for various organic transformations.¹ BDMS has been employed in numerous reactions including the preparation of α -bromo- β -keto esters and α -bromo enones, from their corresponding keto esters and enones, respectively;² the conversion of epoxides and enamines to α -bromoketones;³ and electrophilic aromatic bromination.⁴ BDMS has also been employed in the synthesis of α -aminonitriles and homoallylic amines, and in Michael additions of amines to electron-deficient olefins.⁵



(1) Choudhury, L. H. *Synlett* **2006**, 1619. (2) (a) Khan, A. T. et al. *J. Org. Chem.* **2006**, *71*, 8961. (b) Chow, Y. L.; Bakker, B. H. *Can. J. Chem.* **1982**, *60*, 2268. (3) Olah, G. A. et al. *Tetrahedron Lett.* **1979**, *20*, 3653. (4) (a) Majetich, G. et al. *J. Org. Chem.* **1997**, *62*, 4321. (b) Megyeri, G.; Keve, T. *Synth. Commun.* **1989**, *19*, 3415. (5) (a) Das, B. et al. *Synthesis* **2006**, 1419. (b) Das, B. et al. *Tetrahedron Lett.* **2006**, *47*, 5041. (c) Khan, A. T. *Tetrahedron Lett.* **2007**, *48*, 3805

Bromodimethylsulfonium bromide, 95% BDMS

694142

[50450-21-0]

C₂H₆Br₂S




5 g \$60.00

25 g 217.50

Cumulated Ylide

The cumulated ylide (triphenylphosphoranylidene)ketene is a versatile two-carbon building block useful in preparing numerous classes of oxygen- and nitrogen-containing heterocycles. While typically not Wittig-active itself, this reagent reacts with a host of electrophiles to yield Wittig-active products that can participate in subsequent intra- or intermolecular olefination reactions.



(1) Schobert, R.; Jagusch, C. *Synthesis* **2005**, 2421. (2) Schobert, R. et al. *Synthesis* **2006**, 3902. (3) Boeckman, R. K., Jr. et al. *J. Am. Chem. Soc.* **2006**, *128*, 11032.

(Triphenylphosphoranylidene)ketene


Bestmann ylide

688185

[15596-07-3]

C₂₀H₁₅OP


FW: 302.31



1 g \$55.00

Mg(HMDS)₂

While lithium amides, such as LDA and LiHMDS, are the predominant bases of choice for the selective generation of enolates, magnesium amide bases have garnered recent attention due to their enhanced thermal stabilities and selectivity characteristics.¹ Magnesium bis(hexamethyldisilazide), Mg(HMDS)₂, has demonstrated efficacy in ketone-aldehyde aldol addition reactions² and in the regio- and stereoselective synthesis of silyl enol ethers.³



(1) He, X. et al. *J. Am. Chem. Soc.* **2006**, *128*, 13599. (2) Allan, J. F. et al. *Chem. Commun.* **1999**, 1325. (3) Bonafoux, D. et al. *J. Org. Chem.* **1996**, *61*, 5532.

Magnesium bis(hexamethyldisilazide)


Mg(HMDS)₂

692352

[857367-60-3]

C₁₂H₂₆MgN₂Si₄

FW: 345.07




5 g \$36.00

25 g 120.00

Reagent for Disulfide Synthesis

The reaction of alkyl halides or pseudohalides with the sulfur-based nucleophile sodium methanethiosulfonate (NaMTS) yields organic methanethiosulfonates, intermediates that are highly reactive towards sulfhydryl groups, yielding unsymmetrical disulfides. NaMTS has been extensively used in glycosylations, spin-labeling, and photoprobe chemistry.¹⁻³



(1) Grayson, E. J. et al. *J. Org. Chem.* **2005**, *70*, 9740. (2) Kálai, T. et al. *Synthesis* **2006**, 439. (3) Guo, L.-W. et al. *Bioconjugate Chem.* **2005**, *16*, 685.

Sodium methanethiosulfonate, 95% NaMTS

684538

[1950-85-2]

CH₃O₂NaS₂

FW: 134.15



1 g \$39.50

Aldrichimica ACTA

VOL. 40, NO. 4 • 2007

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
 6000 N. Teutonia Ave.
 Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
 FAX 800-325-5052 (USA)
 or 414-438-2199
 Mail P.O. Box 2060
 Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
 Technical Service 800-231-8327
 SAFC® 800-244-1173
 Custom Synthesis 800-244-1173
 Flavors & Fragrances 800-227-4563
 International 414-438-3850
 24-Hour Emergency 414-438-3850
 Web Site sigma-aldrich.com
 Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
 P.O. Box 355, Milwaukee, WI 53201, USA

Subscriptions

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
 Mail: Attn: Mailroom
 Aldrich Chemical Co., Inc.
 Sigma-Aldrich Corporation
 P.O. Box 355
 Milwaukee, WI 53201-9358
 Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the inside back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2007 Sigma-Aldrich Co.

"PLEASE BOTHER US."

Joe Porwoll, President
 Aldrich Chemical Co.

Professor Keith Fagnou of the University of Ottawa (Canada) kindly suggested that we make cesium pivalate. This versatile base has been used by several research groups in transition-metal-catalyzed direct arylations of indoles.¹ Cesium pivalate has also been employed as a base in the study of aryl-to-aryl palladium migration in Heck and Suzuki couplings of *o*-halobiaryls.²

(1) (a) Stuart, D. R.; Fagnou, K. *Science* **2007**, *316*, 1172. (b) Campeau, L.-C. et al. *Aldrichimica Acta* **2007**, *40*, 35. (c) Wang, X. et al. *J. Am. Chem. Soc.* **2005**, *127*, 4996. (2) Campo, M. A. et al. *J. Am. Chem. Soc.* **2007**, *129*, 6298



694037 Cesium pivalate, 98% 5 g \$40.00
 25 g 134.00

Naturally, we made this useful reagent. It was no bother at all, just a pleasure to be able to help

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the inside back cover.


TABLE OF CONTENTS

Recent Advances in the Chemistry of Allenes 91
 Shengming Ma, Shanghai Institute of Organic Chemistry

Chiral Imidazolium Ionic Liquids: Their Synthesis and Influence on the Outcome of Organic Reactions 107
 Allan D. Headley* and Bukuo Ni, Texas A&M University-Commerce

ABOUT OUR COVER

In the last decades of Monet's life, his prized water garden became his most important—and eventually his only—subject. Monet began to paint his lily pond and our cover, *The Japanese Footbridge* (oil on canvas, 81.3 × 101.6 cm), in his garden at Giverny in 1899. He constructed the water garden soon after he moved to Giverny with his family in 1893, petitioning local authorities to divert water from the nearby river. Monet remade the landscape with the same artifice he applied to his paintings—and then he used it, in turn, as his creative focus.



Photograph © Board of Trustees, National Gallery of Art, V.


When Monet exhibited his lily ponds, a number of critics mentioned his debt to Japanese art and the idea of the *hortus conclusus* (closed garden) of medieval images

Monet painted his garden from the same vantage point as our cover twelve times, focusing on the arching blue-green bridge and a microcosm of the water. He gave equal emphasis to the physical qualities of his painting materials and to the landscape motif he depicted. In this painting, the sky has already disappeared; the lush foliage rises all the way to the horizon; and the decorative arch of the bridge flattens the space. Floating lily pads and mirrored reflections assume equal stature, blurring distinctions between solid objects and transitory effects of light. Monet had always been interested in reflections, seeing their fragmented forms as the natural equivalence for his own broken brushwork. The artist—who, as a leader of the impressionists, had espoused the spontaneity of directly observed works that capture the fleeting effects of light and color—had in these later paintings subjected a nature he re-created to sustained, meditated scrutiny





This painting is a gift to the National Gallery of Art from Victoria Nebeker Coberly, in memory of her son, John W. Mudd, and Walter H. and Leonore Annenberg.

New Reactive Allenes







Allenes are becoming highly sought building blocks for their ability to react with many different classes of substrates. In addition, the allene moiety itself is present in many bioactive natural products and pharmaceutical agents. Recent work by Reissig¹ and Zhang² demonstrates the utility of lithiated allenyl ethers in the synthesis of various carbocycles and heterocycles. Other work by Suginome and co-workers reports the use of cyclohexyllallene in a palladium-catalyzed asymmetric silaboration.³ Sigma-Aldrich is pleased to add these and other allenenes to our expanding portfolio of reactive building blocks.



New Allenes



Other Allenes



References: (1) (a) Brasholz, M.; Reissig, H.-U. *Synlett* **2007**, 1294. (b) Brasholz, M.; Reissig, H.-U. *Angew. Chem., Int. Ed.* **2007**, 46, 1634. (c) Gwiazda, M.; Reissig, H.-U. *Synlett* **2006**, 1683. (2) Huang, X.; Zhang, L. *J. Am. Chem. Soc.* **2007**, 129, 6398. (3) Ohmura, T. et al. *J. Am. Chem. Soc.* **2006**, 128, 13682.

Sigma-Aldrich Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA S.A.
Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

SIGMA-ALDRICH PTY LTD.
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA.
Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

China

SIGMA-ALDRICH (SHANGHAI)
TRADING CO. LTD.
Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

SIGMA-ALDRICH spol. s r. o.
Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: (+45) 43 56 59 10
Fax: (+45) 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND OY
Tel: (+358) 9 350 9250
Free Fax: (+358) 9 350 92555

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 00
Fax: (+33) 474 95 68 08

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Greece

SIGMA-ALDRICH (O.M.) LTD.
Tel: (+30) 210 994 8010
Fax: (+30) 210 994 3831

Hungary

SIGMA-ALDRICH Kft
Ingyenes zöld telefon: 06 80 355 355
Ingyenes zöld fax: 06 80 344 344
Tel: (+36) 1 235 9055
Fax: (+36) 1 235 9050

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: (+91) 80 6621 9600
New Delhi: (+91) 11 4165 4255
Mumbai: (+91) 22 2570 2364
Hyderabad: (+91) 40 4015 5488
Fax
Bangalore: (+91) 80 6621 9650
New Delhi: (+91) 11 4165 4266
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 4015 5466

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 1 404 1900
Fax: (+353) 1 404 1910

Israel

SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

SIGMA-ALDRICH S.r.l.
Numero Verde: 800 827018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

SIGMA-ALDRICH JAPAN K.K.
Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

SIGMA-ALDRICH KOREA
Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A.
de C.V.
Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: 52 722 276 1600
Fax: 52 722 276 1601

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

SIGMA-ALDRICH NEW ZEALAND LTD.
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

SIGMA-ALDRICH NORWAY AS
Tel: (+47) 23 17 60 60
Fax: (+47) 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

SIGMA-ALDRICH RUS, LLC
Tel: +7 (495) 621 6037
+7 (495) 621 5828
Fax: +7 (495) 621 5923

Singapore

SIGMA-ALDRICH PTE. LTD.
Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD.
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313
SAFC (UK) Free Tel: 01202 712305

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Call Collect: (+1) 314 771 5750
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Internet

sigma-aldrich.com

Sigma-Aldrich Career Opportunities

As a leading Life Science and High Technology company, we are always looking for talented individuals to join our team. At Sigma-Aldrich we value the contributions of our employees, and recognize the impact they have on our success. We strive to foster creativity and innovation, and encourage professional development.

Our biochemical and organic chemical products and kits are used in scientific and genomic research, biotechnology, pharmaceutical development, the diagnosis of disease, and as key components in pharmaceutical and other high technology manufacturing. We have customers in life science companies, university and government institutions, hospitals, and in industry.

UNLEASH YOUR TALENTS

Learn more about our career opportunities by visiting
our award-winning Web site at sigma-aldrich.com/careers.

Sigma-Aldrich Corporation is an equal opportunity employer.

Simplify
Discovery Chemistry with ChemDose®



X-Phos

Pd(OAc)₂Pd(dppf)Cl₂·CH₂Cl₂

PEPPSI™-iPr

PdCl₂(PPh₃)₂Pd₂(dba)₃

HATU

(Sampling of available tablets)

Fast, Accurate, and Convenient Dosing of Catalysts and Reagents

ChemDose® is a novel High-Throughput Screening (HTS) technology allowing for the use of chemical reagents and catalysts in the form of tablets. Co-developed with Reaxa Ltd., the ~5 millimeter tablets are composed of a chemically inert magnesium aluminosilicate matrix, together with an absorbed reagent or catalyst. Upon exposure to solvents, the reagent or catalyst readily dissolves out, leaving behind an easily removed insoluble tablet.¹

- Convenient handling and dispensing of reagents and catalysts.
- Eliminates the tedious weighing process for milli- and micromolar chemical quantities.
- Simple reaction workup: inert tablet easily removed upon completion of reaction.
- Consistent chemical loadings of tablets, controlled release rates, microwave compatible.
- Virtually no "learning curve" when using ChemDose®.

For more information, visit sigma-aldrich.com/chemdose.

Reference: (1) Ruhland, T. et al. *J. Comb. Chem.* **2007**, *9*, 301.

PEPPSI is a trademark of Total Synthesis Ltd. (Toronto, Canada). ChemDose is a registered trademark of Reaxa Ltd



P.O. Box 14508
St. Louis, MO 63178
USA

PRESORTED
STANDARD
U.S. POSTAGE
PAID
SIGMA-ALDRICH
CORPORATION

Change Service Requested

IMPROVED CATALYSTS AND LIGANDS FOR ASYMMETRIC SYNTHESIS

Aldrichimica ACTA

VOL. 41, NO. 1 • 2008

 **ALDRICH**
Chemistry



Practical Organocatalysis with (*S*)- and
(*R*)-5-Pyrrolidin-2-yl-1*H*-tetrazoles

Aminophosphine Catalysts in Modern Asymmetric Synthesis


SIGMA-ALDRICH

New Products from Aldrich R&D

Aldrich Is Pleased to Offer Cutting-Edge Tools for Organic Synthesis

Togni Reagent for Electrophilic Trifluoromethylation

The direct transfer of a trifluoromethyl group usually requires harsh conditions that are often incompatible with more sensitive functionalities in a molecule. The Togni Reagent is an electrophilic reagent based on hypervalent iodine and is captured by a range of nucleophilic substrates under mild conditions. This reagent nicely complements the nucleophilic Ruppert's Reagent.



Eisenberger, P. et al. *Chem.—Eur. J.* **2006**, *12*, 2579


3,3-Dimethyl-1-(trifluoromethyl)-1,2-benziodoxole, 97%

696641 250 mg \$71.50
 [887144-97-0] 1 g 199.00
 C₉H₉F₃O
 FW: 330.09



Ligands for Aqueous Transfer Hydrogenation


When used in conjunction with [(Cp*IrCl₂)₂], the ligands *N*-tosylethylenediamine (Ts(en)) and *N*-(2-aminoethyl)-4-(trifluoromethyl)benzenesulfonamide (CF₃-Ts(en)) enable the facile and selective transfer hydrogenation of aldehydes with TOFs as high as 1.3 × 10⁵ h⁻¹. Furthermore, the reactions are carried out in aqueous media and exhibit very good chemoselectivity and functional-group tolerance. In cases where α,β-unsaturated aldehydes are employed, reduction occurs selectively on the formyl group. Aliphatic aldehydes are also readily converted when the substrate is added portionwise over the course of the reaction.



Wu, X. et al. *Angew. Chem., Int. Ed.* **2006**, *45*, 6718


N-Tosylethylenediamine, 97%

693855 5 g \$57.50
 [14316-16-6] C₉H₁₄N₂O₂S
 FW: 214.28




N-(2-Aminoethyl)-4-(trifluoromethyl)benzenesulfonamide, 97%

694967 500 mg \$49.00
 C₉H₉F₃N₂O₂S
 FW: 268.26



Thiocarbonyl Transfer Agent

1-(Methyldithiocarbonyl)imidazole is a stable, non-hazardous reagent that can replace thiophosgene, isothiocyanates, chlorothioformates, and other high-hazard reagents as a thiocarbonyl-transfer agent for the synthesis of dithiocarbonates,¹ dithiocarbamates, symmetrical and unsymmetrical thioureas,² and 2-thiohydantoins.³



(1) Sun, W. Y. et al. *Synlett* **1997**, 1279. (2) Mohanta, P. K. et al. *Tetrahedron* **2000**, *56*, 629. (3) Sundaram, G. S. M. et al. *Synlett* **2007**, 251


1-(Methyldithiocarbonyl)imidazole, 97%

694029 5 g \$74.50
 [74734-11-5] C₅H₈N₂S₂
 FW: 158.24



Air-Stable, Nucleophilic Alkylphosphine

1,3,5-Triaza-7-phosphaadamantane (PTA) is a convenient, efficient, and air-stable nucleophilic trialkylphosphine organocatalyst for the Baylis-Hillman reaction. Both aromatic and aliphatic aldehydes react with activated alkenes in the presence of 15–20 mol % of PTA to afford the corresponding adducts in fair-to-excellent yields. Furthermore, PTA displays activity that is superior to that of the structurally similar hexamethylenetetramine.



He, Z. et al. *Adv. Synth. Catal.* **2006**, *348*, 413

1,3,5-Triaza-7-phosphaadamantane

695467 500 mg \$25.00
 [53597-69-6] 2 g 90.00
 C₆H₁₂N₃P
 FW: 157.15



Aldrichimica ACTA

VOL. 41, NO. 1 • 2008

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-231-8327
SAFC® 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 414-438-3850
24-Hour Emergency 414-438-3850
Web Site sigma-aldrich.com
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201, USA

Subscriptions

To request your FREE subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 355
Milwaukee, WI 53201-9358
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the inside back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com/acta.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2008 Sigma-Aldrich Co.

"PLEASE BOTHER US."



Joe Porwoll, President
Aldrich Chemical Co., Inc.

Professor Matthew Clarke of the University of St. Andrews (U.K.) kindly suggested that we make 1,3,5,7-tetramethyl-6-phenyl-2,4,8-trioxo-6-phosphaadamantane. This phosphine is very stable to air and moisture and has the stereoelectronic properties of bulky phosphonites. When employed with rhodium complexes, this ligand shows high catalytic activity in the hydroformylation of various alkenes. High selectivities and conversions as high as 99% have been reported.

Clarke, M. L.; Roff, G. J. *Chem.—Eur. J.* **2006**, *12*, 7978



695459 1,3,5,7-Tetramethyl-6-phenyl-2,4,8-trioxo-6-phosphaadamantane, 97% 500 mg \$52.00
2 g 170.00

Naturally, we made this useful reagent. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the inside back cover.

TABLE OF CONTENTS

- Practical Organocatalysis with (S)- and (R)-5-Pyrrolidin-2-yl-1H-tetrazoles** 3
Deborah A. Longbottom, Vilijus Franckevičius, Sirirat Kumarn, Alexander J. Oelke, Veit Wascholowski, and Steven V. Ley; University of Cambridge
- Aminophosphine Catalysts in Modern Asymmetric Synthesis** 15
Dino Amoroso, Todd W. Graham, Rongwei Guo, Chi-Wing Tsang, and Kamaluddin Abdur-Rashid; Kanata Chemical Technologies, Inc

ABOUT OUR COVER

Landscape with Tobias and the Angel, with a View of Antwerp in the Background (oil on copper, 20.5 × 26.0 cm) was painted possibly around 1665 by Gillis Neyts, an enigmatic Flemish painter and engraver. Neyts (1623–1687) was born in Ghent, and spent a good part of his life in the city of Antwerp. He specialized in small, imaginary landscape scenes, which sometimes incorporated historical material or views of Flemish towns. His style approaches that of Lucas van Uden (1595–1672; Antwerp), who may have been his teacher.




Photograph © Alfred Bader

This small painting, with its soft and delicate handling, which was typical for Neyts, shows on the left just below the horizon a part of the skyline of the city of Antwerp. The spectacular form of the arching tree in the center frames the figures of two travelers (with walking sticks) in the foreground on the right. One of them appears to waive at the viewer, while the other—dressed in red and white and with wings rising from his shoulders—is identified as the Archangel Raphael accompanying young Tobias on his journey.


Neyts has painted here a fantasy landscape in which he transposes the ancient story of Tobias and the angel onto a contemporary setting, the outskirts of the 17th-century city of Antwerp. It would appear that Neyts's purpose is to help the viewer of that period identify more closely with the story.

This painting is in the private collection of Isabel and Alfred Bader.


Explore Organocatalysis




P0380




684341




688274



677183



671428




677019

Organocatalysis provides convenient new methods to construct complex chiral compounds with operational simplicity and without the need for metals. Excellent selectivities are observed in various organocatalyzed asymmetric transformations.

Sigma-Aldrich is pleased to offer an expanding portfolio of innovative organocatalysts to accelerate your research success.

To see a comprehensive listing of organocatalysts, please visit sigma-aldrich.com/organocat



Practical Organocatalysis with (*S*)- and (*R*)-5-Pyrrolidin-2-yl-1*H*-tetrazoles



Deborah A. Longbottom, Vilius Franckevičius, Sirirat Kumarn, Alexander J. Oelke, Veit Wascholowski, and Steven V. Ley*
 Department of Chemistry
 University of Cambridge
 Lensfield Road
 Cambridge CB2 1EW, U.K.
 Email: svl1000@cam.ac.uk

From left to right: Dr. Deborah A. Longbottom, Dr. Veit Wascholowski, Mr. Vilius Franckevičius, Prof. Steven V. Ley, Mr. Alexander J. Oelke, and Ms. Sirirat Kumarn.



Professor Ley (center) receiving the Sigma-Aldrich sponsored 2007 ACS Award for Creative Work in Synthetic Organic Chemistry. Pictured with Professor Ley are Dr. John Chan (left), Sigma-Aldrich Market Segment Manager, and Dr. Catherine T. Hunt (right), 2007 ACS President. Photo © Peter Cutts Photography, LLC.

Outline

1. Introduction
2. The Aldol Reaction
3. The Mannich Reaction
4. Conjugate Additions
5. α -Aminoxylation, α -Hydroxyamination, and α -Amination
6. One-Pot Reaction Processes
7. Conclusions
8. References and Notes

1. Introduction

The phenomenal renaissance of interest in organocatalysis has been fuelled by the ever-increasing repertoire of organocatalytic reactions of utility to the synthetic organic chemist. Innovative reactions appear weekly in research publications throughout the world, and these developments have spawned a search for newer and more effective catalysts to bring about a myriad of important chemical transformations.^{1,2}

As with any evolving scientific (sub)discipline, there exists a need to provide a range of tools to solve particular problems and stimulate the creation of new concepts. The properties, function, and mechanisms of action of the individual organocatalysts are of prime importance for they must tolerate a wide range of chemistries, functional groups, and reaction conditions and, ideally, be of broad synthetic utility.

The simple amino acids L- and D-proline (**1** and **2**, **Figure 1**) have been widely utilized in many organocatalytic reactions and are often considered the benchmark with respect to which other catalysts are evaluated. Nevertheless, their lack of solubility in certain solvents and sometimes slow turnover rates have caused concern and, therefore, led to the discovery of other related catalytic systems that overcome some of these drawbacks.

For example, (*S*)- and (*R*)-5-pyrrolidin-2-yl-1*H*-tetrazole (**3** and **4**) are isosteres of proline with similar pK_a 's but anticipated greater solubility and, hence, reactivity in more lipophilic organic solvents. They were originally synthesized for organocatalytic applications almost simultaneously by three groups, Yamamoto's,³ Arvidsson's,⁴ and ours,⁵ and have since proven very useful in a wide range of reactions.

Herein, we discuss the practical synthetic opportunities that have arisen through the development of these new catalytic species, which are shelf- and thermally stable, crystalline, and readily prepared on scale.⁶⁻⁸ Emphasis is given to reaction type, rather than to detailed mechanistic discussion, as this is still the subject of much study and debate. In each reaction table, several examples have been selected from the original publication(s) to represent the breadth in substrate substitution, electronic character, and general compatibility of functional groups.

2. The Aldol Reaction

The aldol reaction is one of the most important carbon-carbon bond-forming reactions and, therefore, the widespread interest in developing asymmetric variants of this transformation is not surprising. The direct asymmetric addition of unmodified ketones to aldehydes has been developed by Shibasaki's and Trost's groups by using heterobimetallic catalysts,⁹ whereas others have used more nature-inspired catalytic systems






Figure 1. L- and D-Proline and (S)- and (R)-5-Pyrrolidin-2-yl-1H-tetrazole.




Scheme 1. Transition States in the Pyrrolidinyltetrazole **3** Mediated Aldol Reaction.




Ref. 7

eq 1



Ref. 7

Scheme 2. Parasitic Consumption of the Catalyst in the Case of L-Proline (**1**), but not in the Case of Pyrrolidinyltetrazole **3**.



Ref. 3

eq 2

consisting of aldolase enzymes and catalytic antibodies.¹⁰ An organocatalytic approach that uses L-proline (**1**) as catalyst for an intramolecular aldol cyclization, known as the Hajos–Parrish–Eder–Sauer–Wiechert reaction, was reported around 35 years ago.¹¹ More recently, following List's important work,² a number of groups have confirmed that L-proline (**1**) can also mediate the analogous intermolecular aldol reaction of unmodified ketones and aldehydes.¹²

To date, there have been three publications that focus fully on the ability of pyrrolidinyltetrazole **3** to facilitate the intermolecular aldol reaction,^{3,4,7} and its usefulness has now been further demonstrated in an aldol reaction applied to natural product synthesis.¹³

As alluded to previously, the details of the mechanism of these reactions are generally the subject of much debate and discussion. However, two widely accepted transition state models for the aldol reaction catalyzed by **3** produce the same stereochemical outcome and involve an enamine intermediate reacting either via a coordinated Zimmerman–Traxler-type transition state, **5**, or via transition state **6** (Scheme 1).⁷

Hartikka and Arvidsson have shown that aliphatic aldehydes are generally less reactive than aromatic ones in the direct asymmetric aldol reaction between acetone (**8**) and a variety of aldehydes leading to β -hydroxy ketones **9** (eq 1).^{4,7} Nevertheless, the high catalytic activity of pyrrolidinyltetrazole **3** allowed even aliphatic aldehydes to be transformed into the corresponding chiral β -hydroxy ketones **9** with high enantioselectivities and fair yields in thirteen hours or less. The authors additionally proved that even a catalyst loading of 5% was still effective, though a longer reaction time was required.

It is interesting to note at this point that parasitic catalyst consumption¹⁴ is observed with L-proline (**1**) but not with **3**. Arvidsson carried out NMR studies using a mixture of **1** or **3** and 2,2-dimethylpropionaldehyde (**11**) and proved that, while L-proline (**1**) easily engages in parasitic formation of bicyclic oxazolidinone **10**, pyrrolidinyltetrazole **3** does not (Scheme 2).⁷ Consequently, in theory, this results in more of **3** being available to form the postulated enamine intermediate in the aldol reaction. The authors suggested that this could be the main reason for the increased reactivity of **3** compared to that of L-proline (**1**) in DMSO.¹⁵ However, this does not rule out the possibility that factors relating either to the increased solubility of **3** in DMSO or to alternative mechanistic pathways operating in other solvent systems may also be contributing to the observed enhancement in reactivity of **3**.

Yamamoto's work has focused on the formation of optically active 1,1,1-trichloro-2-alkanols (eq 2),³ previously demonstrated as being versatile tools for the preparation of variously functionalized compounds such as α -hydroxy and α -amino acids.¹⁶ The formation of 1,1,1-trichloro-2-alkanols by the asymmetric aldol reaction is challenging due to the propensity of the starting aldehydes to form hydrates. However, in Yamamoto's report, **3** displayed excellent catalytic efficiency and activity in the reaction of either chloral monohydrate (**14**) or chloral (**15**) and water with a variety of aliphatic and aromatic ketones.

Ward then proved that pyrrolidinyltetrazole **3** was also useful in a total synthesis program:¹³ serricornin (**21**), a sex pheromone produced by the female cigarette beetle *Lasioderma serricorne* F., was elegantly prepared in just seven steps and overall 31% yield from the readily available *racemic* aldehyde **17** with **18**, catalyzed by **3**, was the key step in the synthesis

and occurred with dynamic kinetic resolution to give adduct **19** with >98% ee. Aldol product **19** was then smoothly converted into diol **20** with excellent yield and diastereoselectivity. Because diastereomeric diols **22**, **23**, and **24** (Figure 2) are also readily prepared from **19**,¹⁷ it was proposed that this powerful strategy could be extended to afford stereoisomers of **21**, which could then be tested for biological activity.

3. The Mannich Reaction


The development of syntheses that provide enantiomerically pure α -amino acids has been the subject of generations of research by organic chemists. This has engendered an array of methodologies,¹⁸ which, not only allow for the stereoselective construction of naturally occurring amino acids, but also permit the rational design of optically active nonproteinogenic ones. These unnatural amino acids in particular have enjoyed increased popularity, mainly due to their incorporation into nonscissile peptide mimetics and peptide isosteres, known to exhibit reduced susceptibility to catabolism and thus increased bioavailability.¹⁹

In a similar way, chiral diamines are important building blocks for pharmaceuticals and are features that are frequently found in natural products.^{20,21} As synthetic tools, chiral diamines are also used extensively as chiral auxiliaries and catalysts.²² However, despite their significance, their asymmetric synthesis is not straightforward: they are most frequently synthesized from diols or aziridines²¹ or by addition of glycine ester enolates to imines.²³ The direct reductive coupling of imines has also been reported, but this approach is limited to the preparation of symmetrical vicinal diamines and results in relatively low stereoselectivity.²⁴

Thus, the organocatalytic synthesis of enantiomerically pure α -amino acids and diamines had so far represented a worthwhile challenge to organic chemists. Gratifyingly, pyrrolidinyltetrazole **3** has now been used to great effect in the synthesis of both classes of compound by employing the Mannich reaction as the key carbon-carbon-bond-forming step. In the synthesis of α -amino acid derivatives—which serves as an excellent comparison with previous work by Barbas (where L-proline (**1**) is the catalytic species)²⁵—our group has shown that **3** also effectively catalyzes this reaction (eq 3).^{5,26} Indeed, our method represents a very attractive alternative to Barbas's as the yields and stereoselectivities are comparable to those obtained with L-proline (**1**),²⁵ yet the reaction is carried out in solvents such as dichloromethane (avoiding dimethyl sulfoxide) and with catalyst loadings as low as 1%.²⁷

Following this report, Barbas showed that the organocatalytic asymmetric Mannich reaction of protected amino ketones with imines in the presence of **3** affords diamines with excellent yields and enantioselectivities of up to 99%.²⁸ The amino ketone protecting group controlled the regioselectivity of the reaction, providing access to chiral 1,2- and 1,4-diamines from azido and phthalimido ketones, respectively. Under optimized conditions, the three-component Mannich reaction of various combinations of azido ketones and aldehydes was investigated (eq 4).²⁸ All products were obtained regioselectively and with good diastereo- (syn:anti = 70:30 to 91:9) and enantioselectivities (82–99% ee, syn). A one-pot reduction–Boc-protection of Mannich product **31b** provided differentially protected 1,2-diamine **32** (eq 5), illustrating the potential utility of these compounds in further synthetic steps. The scope of this reaction seems very broad, and the azido ketone products, **31**, are in themselves interesting substrates for potential “click chemistry” based diversification.²⁹

The Mannich reaction of phthalimidoacetone (**33**), a phthaloyl-protected amino ketone, in *N*-methyl-2-pyrrolidone (NMP) as



Scheme 3. Enantioselective Total Synthesis of Serricornin (**21**) from Racemic Aldehyde **18**.






Figure 2. 19-Derived Diastereomeric Diols for the Synthesis of Serricornin Analogues.



^a Reaction carried out in acetone


Ref 5, 26

eq 3




Ref 28

eq 4



Ref 28


eq 5



35	R ¹	Yield	ee
a	CO ₂ Et	88%	91%
b	4-O ₂ NC ₆ H ₄	68%	97%
c	4-NCC ₆ H ₄	56%	95%
d	Ph	71%	77%


Ref. 28

eq 6



Ref. 28


Figure 3. Major Product of the Mannich Reaction of Phthalimidoacetone in the Presence of L-Proline (1).



39	R ¹ , R ²	Yield	Syn:Anti	ee
a	-(CH ₂) ₂ SCH ₂ -	62%	10:1	70%
b	-(CH ₂) ₂ OCH ₂ -	94%	6:1	54%
c	Me, Me	71%	10:1	32%
d	Me, H	72%	-	3%
e	Et, Me	68%	>19:1	65%

Ref. 26,30


eq 7



43	R ¹ , R ²	R ³	R ⁴	Yield	ee
a	-(CH ₂) ₃ -	Me	H	64%	91%
b	-(CH ₂) ₂ -	H	Me	62%	80%
c	Ph, Et	H	Me	78%	76%
d	thien-2-yl, Me	H	Me	61%	72%
e	CO ₂ Me, Me	H	Me	96%	82%
f	n-Pent, Me	H	Me	44%	58%

Ref. 32

eq 8



47	R ¹ , R ²	Yield	ee
a	-(CH ₂) ₃ -	87%	83%
b	Ph, Me	89%	84%
c	4-F ₃ CC ₆ H ₄ , Me	94%	75%
d	4-HOC ₆ H ₄ , Me	70%	64%
e	furan-2-yl, Me	69%	81%
f	thien-2-yl, Me	82%	84%

Ref. 33

eq 9

solvent exhibited excellent regioselectivity: Enamine formation was favored at the less hindered side of the carbonyl group and generated protected 1,4-diamines **35** in good yields and enantioselectivities (eq 6).²⁸ Interestingly, when the formation of **35a** was attempted using L-proline (**1**) as catalyst in NMP at room temperature, phthalimidoacetone (**33**) provided Mannich product **35a** only in trace amounts accompanying the formation of cycloadduct **36** (Figure 3) as the major product in 59% yield. It would thus appear that L-proline (**1**) is not a useful catalyst in this reaction.

To summarize, the pyrrolidinyltetrazole **3** mediated Mannich reaction provides efficient access to several highly important product types, namely chiral α -amino acids and diamines. It can be performed under environmentally favorable conditions without the requirement for inert atmosphere or dry solvents, and provides good-to-excellent yields and regio- and stereoselectivities.

4. Conjugate Additions

The organocatalytic conjugate addition of nucleophiles to nitroolefins^{26,30} and enones^{31,33} can also be mediated by **3** and leads to useful adducts such as γ -nitro ketones and 1,5-dicarbonyl compounds (eq 7-9).

In the former case (see eq 7),^{26,30} ketone-derived enamines add to electrophilic nitroalkenes to form Michael adducts **39**, which are useful synthetic precursors to other functionalities that are derived from the nitro group.³⁴ Interestingly, when compared to the results obtained with L-proline (**1**),³⁵ pyrrolidinyltetrazole **3** far outperformed it in terms of yield, enantioselectivity, reaction times, and stoichiometry. However, despite the fact that the results published were the best in the literature at that time, they still left some room for improvement, and it was only when the homoproline tetrazole derivative **50** (eq 10) was used that the yields and enantioselectivities moved to practically useful levels.^{36,37}

In the addition of nitroalkanes to enones, the same type of γ -nitro ketone adduct is formed but, due to the nature of the reaction, products with alternative structural features are obtained (see eq 8).^{31,32} In this case, pyrrolidinyltetrazole **3** proved to be a versatile catalyst for the asymmetric 1,4 addition of a variety of nitroalkanes to cyclic and acyclic enones, using *trans*-2,5-dimethylpiperazine (**42**) as a stoichiometric base additive. Using **3**, the reaction was also scalable, providing enantiomeric excesses of up to 98% in relatively short reaction times (1-3 days) and employing just two equivalents of the coupling nitroalkane.³⁸ Kinetic investigations, combined with the observed sensitivity of certain substrates to water, support the iminium catalysis mechanism.³²

The addition of malonates to enones (see eq 9)³³ leads to a variety of useful 1,5-dicarbonyl compounds. In the case of **3**,³⁹ only 1.5 equivalents of malonate is needed, and the reaction is readily scaled and practical to operate,³³ rendering the process potentially useful in a synthesis program. The utility of such addition products in synthesis has now been further proved by carrying out the decarboxylation of **47a** to the corresponding monomethyl ester (**53**, eq 11). While a loss in enantiomeric excess had been observed under various Krapcho conditions,^{40,41} it has now been shown that sodium hydroxide or porcine liver esterase (PLE) smoothly mediates the monohydrolysis of **47a**; subsequent decarboxylation provides the corresponding methyl ester, **53**, in excellent yield, with neither step resulting in any reduction in enantiomeric purity.⁴¹

Finally, it was thought that an extension of these conjugate addition methods might be useful in a new organocatalytic

asymmetric nitrocyclopropanation reaction. Cyclopropane-containing structures are compounds of interest within organic chemistry as they display a relatively large amount of stereochemical information over a small, rigid framework of just three carbon atoms. They serve as versatile synthetic intermediates in a variety of reactions⁴² and are widely distributed in a range of naturally occurring compounds⁴³ and peptidomimetics.⁴⁴ Consequently, their stereoselective preparation is a valuable goal and, to date, several methods have been developed towards this aim.⁴⁵


In particular, nitrocyclopropanes may be converted into a wide range of functional groups^{34,46} and can be prepared by a variety of methods.⁴⁷ Therefore, it was surprising that there were only two synthetic approaches,⁴⁸ prior to the one described below,⁴⁹ that detailed their enantioselective formation. Indeed, the novel organocatalytic method developed by our group provides a higher yield and enantioselection than is found in either: following optimization studies, 7-nitrobicyclo[4.1.0]heptanone **56** (eq 12) was provided in 80% yield and 77% ee, which was then easily improved to >98% ee upon a single recrystallization.⁴⁹ More recent experiments have indicated that, under further optimized reaction conditions, not only can this result be improved (87% yield, 90% ee),⁵⁰ but that the reaction is now generally applicable to a wider range of substrates such as aliphatic and aromatic enones, providing useful products in good yields and enantiomeric excesses.⁵⁰

Thus, these conjugate addition procedures can be extremely powerful, providing, not only the products of conjugate addition, but also of tandem reactions such as the nitrocyclopropane example given above. Many further applications of this concept can be envisaged and are currently being investigated in our laboratory.

5. α -Aminoxylation, α -Hydroxyamination, and α -Amination


The regio- and stereoselective replacement of hydrogen by oxygen or nitrogen results in a rapid increase in molecular complexity,⁵¹ and one can see that, with a nitrosobenzene electrophile⁵² under enamine catalysis, either the oxygen- (α -aminoxylation) or nitrogen-substituted (α -hydroxyamination) product might be observed. This can give rise to optically active α,α' -disubstituted oxy- or amino aldehydes.⁵³

The two major independent studies that have been carried out by Yamamoto and Kim, respectively, have shown that a reactivity pattern exists.^{6,54} When ketones and aldehydes with no branching at the α position are employed (eq 13),⁶ generally the preference is for α -aminoxylation; whereas if the substrate is α -branched, α -hydroxyamination is also observed, at least in the case of aldehydes (eq 14).⁵⁴ A plausible explanation for this inherent difference in reactivity is found when the reaction transition state is examined (Figure 4): if α branching is present, the clash between the α substituent and the phenyl group of nitrosobenzene in the usual transition state, **63**, will push the phenyl group into the pseudoequatorial position, **64**. This results in hydrogen bonding of the oxygen rather than the nitrogen atom with the tetrazole portion, thus changing the regiochemical outcome of the reaction. However, although the contrasting regioselectivity of this reaction is usually predictable, the selective formation of α -hydroxyamination products is not yet general: in order to introduce an α -nitrogen substituent, α -branched aldehydes must be utilized, and mixtures of O- and N-substituted products are usually observed.




Ref 36

eq 10




Ref 41

eq 11



Ref 49


eq 12



59	R ¹ , R ²	Yield	ee
a	-(CH ₂) ₄ -	94%	>99%
b	-CH ₂ O(CH ₂) ₂ -	87%	>99%
c	-CH ₂ O(CH ₂) ₂ (CH ₂) ₂ -	92%	99%
d	CH ₂ Ni(BnCO) ₂ (CH ₂) ₃ -	95%	>99%
e	H, Bn	67%	98%
f	H, <i>i</i> -Pr	65%	98%
g	H, <i>n</i> -Bu	69%	98%

Ref 6


eq 13



No.	R ¹	R ²	Yield	61:62	ee 61	ee 62
a	Bn	Me	96%	1.0:1	81%	37%
b	4-MeOC ₆ H ₄ CH ₂	Me	75%	1.7:1	90%	35%
c	4-BrC ₆ H ₄ CH ₂	Me	98%	1.4:1	86%	45%
d	BnOCH ₂	Me	89%	0.8:1	79%	5%
e	allyl	Me	91%	0.7:1	62%	27%
f	<i>n</i> -Bu	Et	55%	0.6:1	5%	2%


Ref 54

eq 14




Ref 6

Figure 4. Postulated Transition States Leading to α -Aminoxylation or α -Hydroxyamination.




Scheme 4. α -Amination of Aldehydes as a Step in the Total Synthesis of BIRT-377 (**68**).



Ref 58a


eq 15



74	R ¹ , R ²	R ³	R ⁴	Yield
a	-(CH ₂) ₄ -	H	H	60%
b	-(CH ₂) ₄ -	Me	H	39%
c	-(CH ₂) ₄ -	H	Me	65%
d	-(CH ₂) ₂ SCH ₂ -	H	H	51%
e	-(CH ₂) ₂ C(OCH ₂ CH ₂ O)CH ₂ -	H	H	50%
f	H, <i>n</i> -Pr	H	H	71%
g	H, <i>t</i> -Bu	H	H	82%
h	H, (CH ₂) ₂ CO ₂ Me	H	H	50%

Ref 59a,b

eq 16



78	R ¹ , R ²	R ³	Yield	ee
a	H, <i>i</i> -Pr	Et	89%	94%
b	H, <i>i</i> -Pr	<i>t</i> -Bu	81%	99%
c	H, <i>i</i> -Pr	Bn	63%	89%
d	H, <i>t</i> -Bu	Et	84%	99%
e	H, allyl	Et	69%	90%
f	H, (CH ₂) ₂ CO ₂ Me	Et	67%	69%
g	-(CH ₂) ₄ -	Et	52%	76%
h	-(CH ₂) ₂ SCH ₂ -	Et	40%	83%
i	-(CH ₂) ₂ C(OCH ₂ CH ₂ O)CH ₂ -	Et	50%	76%

Ref 59c, 60

eq 17

In their total synthesis of BIRT-377 (**68**),⁵⁵ Chowdari and Barbas have shown that, as a related reaction, pyrrolidinyltetrazole **3** mediated α -amination is possible with dibenzyl azodicarboxylate (**66**) as the nitrogen source, and that even a quaternary stereocenter can be formed (**Scheme 4**). The synthesis of quaternary amino acids through organocatalytic amination reactions is challenging, since the *cis* and *trans* enamines derived from α -branched aldehydes are energetically less distinct, as compared with their linear counterparts, and this can lead to low enantioselectivities.⁵⁵ The higher reactivity and enantioselectivity obtained with **3** relative to L-proline (**1**), in the reaction leading to **67**, was ascribed to the lower pK_a and increased steric bulk of the tetrazole relative to proline's carboxylic acid moiety. Indeed, the desired key intermediate **67** was formed in an excellent 95% yield and 80% ee (compared with those observed with L-proline (**1**): 5-day reaction time, 90% yield, and 44% ee). It was suggested that analogues of **67** could be accessed by simply changing the α,α' -disubstituted aldehyde and catalyst stereochemistry. This means much scope remains for investigations into this unexploited area of research.

6. One-Pot Reaction Processes⁵⁶

Presently, organic synthesis can be hampered by time-consuming and costly protecting-group strategies and (lengthy) purification procedures after each synthetic step. In order to circumvent these difficulties, the synthetic potential of multicomponent domino reactions has now been exploited in the efficient and stereoselective construction of complex molecules from simple precursors in a single process. These domino reactions often proceed with excellent stereoselectivities and are generally environmentally more appropriate. The efficiency of asymmetric domino reactions can easily be seen in the number of newly formed bonds, the number of new stereocenters, and the concomitant rapid increase in molecular complexity. In particular, domino reactions mediated by organocatalysts are of great utility as they are characterized by high efficiencies and are in a way biomimetic in origin: the same governing principles are often found in the biosynthesis of natural products.⁵⁷

This field has grown over the last few years and, often, the advantage of employing organocatalysts is their ability to promote several types of reactions through different activation modes. Pyrrolidinyltetrazole **3** has so far been useful in two major tandem reaction types, namely the enantioselective α -aminooxylation–Michael reaction (**eq 15**),⁵⁸ and the formation of chiral 1,2-oxazines (**eq 16**)⁵⁹ and their 3,6-dihydropyridazine congeners (**eq 17**).^{59c,60} In the former example (see **eq 15**), pyrrolidinyltetrazole **3** mediated a highly enantioselective synthesis of Diels–Alder nitroso adducts **71**,⁵⁸ and the results disclosed revealed opposite regioselectivities and increased stereochemical control over the more common and complementary Diels–Alder procedures used to make the same structural motif. In the latter tandem reaction type (see **eq 16**), a new method was developed for the synthesis of enantiomerically pure 1,2-oxazines **74** from achiral starting materials.⁵⁹ This procedure relies on initial α -aminooxylation of an enamine intermediate with nitrosobenzene (**58**), followed by nucleophilic attack on a vinylphosphonium salt **73**, and subsequent intramolecular Wittig reaction. This tandem process is a useful addition to the toolbox of the organic chemist: it is quite general and reliable, and methods for preparing this unit in an optically pure fashion are scarce.

In addition, an analogous method was published for the synthesis of a related heterocycle, the 3,6-dihydropyridazine unit

(78), from aldehydes and ketones.^{59c,60} In the case of aldehydes,⁶⁰ products were formed in generally good-to-excellent yields and enantioselectivities with a variety of nitrogen protecting groups. This method has now been extended to ketones,^{59c} greatly increasing its scope. Furthermore, the selective α amination of aldehydes with differentially protected azodicarboxylates (e.g., BocN=NTroc) has also been developed recently, serving as a useful platform for further selective derivatization of these products.⁶¹

Thus, it can be seen that early results on these one-pot reaction processes show that they can be very powerful, generating molecular complexity extremely quickly. We look forward with great excitement to further publications in this area.

7. Conclusions

In this short review, the variety of practical synthetic opportunities offered by the (*S*)- and (*R*)-pyrrolidinyltetrazole catalysts **3** and **4** has been illustrated. Their utility has been demonstrated beyond doubt: they are indeed worthy catalysts of a number of asymmetric organocatalytic processes and are undeniably useful additions to the rapidly developing armory of shelf-stable catalysts available to the synthetic organic chemist.

8. References and Notes

- For seminal publications in this area, see: (a) Ahrendt, K. A.; Borths, C. J.; MacMillan, D. W. C. *J. Am. Chem. Soc.* **2000**, *122*, 4243. (b) Jen, W. S.; Wiener, J. J. M.; MacMillan, D. W. C. *J. Am. Chem. Soc.* **2000**, *122*, 9874. (c) Bui, T.; Barbas, C. F., III *Tetrahedron Lett.* **2000**, *41*, 6951. (d) List, B. *J. Am. Chem. Soc.* **2000**, *122*, 9336. (e) Notz, W.; List, B. *J. Am. Chem. Soc.* **2000**, *122*, 7386. (f) For a comprehensive book concerning this area, see: Berkessel, A.; Gröger, H. *Asymmetric Organocatalysis: From Biomimetic Concepts to Applications in Asymmetric Synthesis*; Wiley-VCH: Weinheim, 2005. For a selection of leading reviews in the area, see: (g) List, B. *Synlett* **2001**, 1675. (h) List, B. *Tetrahedron* **2002**, *58*, 5573. (i) Jarvo, E. R.; Miller, S. J. *Tetrahedron* **2002**, *58*, 2481. (j) Notz, W.; Tanaka, F.; Barbas, C. F., III *Acc. Chem. Res.* **2004**, *37*, 580. (k) Dalko, P. I.; Moisan, L. *Angew. Chem., Int. Ed.* **2004**, *43*, 5138. (l) Janey, J. M. *Angew. Chem., Int. Ed.* **2005**, *44*, 4292. (m) Hayashi, Y. *J. Synth. Org. Chem. Jpn.* **2005**, *63*, 464. (n) Seayad, J.; List, B. *Org. Biomol. Chem.* **2005**, *3*, 719. (o) Limbach, M. *Chem. Biodiv.* **2006**, *3*, 119. (p) List, B. *Chem. Commun.* **2006**, 819. (q) Guillena, G.; Ramón, D. J. *Tetrahedron: Asymmetry* **2006**, *17*, 1465. (r) Gaunt, M. J.; Johansson, C. C. C.; McNally, A.; Vo, N. T. *Drug Discovery Today* **2007**, *12*, 8. See also the following special issues on organocatalysis: (s) *Asymmetric Organocatalysis*. Houk, K. N., List, B., Eds.; *Acc. Chem. Res.* **2004**, *37*, Issue 8 (August 2004); pp 487–631. (t) *Organic Catalysis. Adv. Synth. Catal.* **2004**, *346*, Issues 9–10 (August 2004); pp 1007–1249.
- List, B.; Lerner, R. A.; Barbas, C. F., III *J. Am. Chem. Soc.* **2000**, *122*, 2395.
- Torii, H.; Nakadai, M.; Ishihara, K.; Saito, S.; Yamamoto, H. *Angew. Chem., Int. Ed.* **2004**, *43*, 1983.
- Hartikka, A.; Arvidsson, P. I. *Tetrahedron: Asymmetry* **2004**, *15*, 1831.
- Cobb, A. J. A.; Shaw, D. M.; Ley, S. V. *Synlett* **2004**, 558.
- Momiyama, N.; Torii, H.; Saito, S.; Yamamoto, H. *Proc. Natl. Acad. Sci. U.S.A.* **2004**, *101*, 5374.
- Hartikka, A.; Arvidsson, P. I. *Eur. J. Org. Chem.* **2005**, 4287.
- (a) McManus, J. M.; Herbst, R. M. *J. Org. Chem.* **1959**, *24*, 1643. (b) Grzonka, Z.; Liberek, B. *Roczniki Chemii (Ann. Soc. Chim. Polonorum)* **1971**, *45*, 967. (c) Grzonka, Z.; Liberek, B. *Tetrahedron* **1971**, *27*, 1783. (d) Grzonka, Z.; Gwizdala, E.; Kofluk, T. *Pol. J. Chem.* **1978**, *52*, 1411. (e) Almquist, R. G.; Chao, W.-R.; Jennings-White, C. *J. Med. Chem.* **1985**, *28*, 1067. (f) Franckevičius, V.; Rahbek Knudsen, K.; Ladlow, M.; Longbottom, D. A.; Ley, S. V. *Synlett* **2006**, 889. (g) Aureggi, V.; Franckevičius, V.; Kitching, M. O.; Ley, S. V.; Longbottom, D. A.; Oelke, A. J.; Sedelmeier, G. *Org. Synth.* **2008**, *85*, 72. For an alternative method to form the tetrazole ring system, see: (h) Sedelmeier, G. Intl. Patent WO 2005/014602 A1, February 17, 2005. (i) Sedelmeier, G.; Aureggi, V. Intl. Patent WO 2007/009716, January 25, 2007.
- (a) Yoshikawa, N.; Yamada, Y. M. A.; Das, J.; Sasai, H.; Shibasaki, M. *J. Am. Chem. Soc.* **1999**, *121*, 4168. (b) Trost, B. M.; Ito, H. *J. Am. Chem. Soc.* **2000**, *122*, 12003.
- (a) Wagner, J.; Lerner, R. A.; Barbas, C. F., III *Science* **1995**, *270*, 1797. (b) Zhong, G.; Hoffmann, T.; Lerner, R. A.; Danishefsky, S.; Barbas, C. F., III *J. Am. Chem. Soc.* **1997**, *119*, 8131.
- (a) Eder, U.; Wiechert, R.; Sauer, G. German Patent DE 2014757.3, October 7, 1971. (b) Hajos, Z. G.; Parrish, D. R. German Patent DE 2102623 C2, July 29, 1971. (c) Eder, U.; Sauer, G.; Wiechert, R. *Angew. Chem., Int. Ed. Engl.* **1971**, *10*, 496. (d) Hajos, Z. G.; Parrish, D. R. *J. Org. Chem.* **1974**, *39*, 1615.
- (a) Sakthivel, K.; Notz, W.; Bui, T.; Barbas, C. F., III *J. Am. Chem. Soc.* **2001**, *123*, 5260. (b) Bogevig, A.; Kumaragurubaran, N.; Jorgensen, K. A. *Chem. Commun.* **2002**, 620. (c) Sekiguchi, Y.; Sasaoka, A.; Shimamoto, A.; Fujioka, S.; Kotsuki, H. *Synlett* **2003**, 1655. (d) Thayumanavan, R.; Tanaka, F.; Barbas, C. F., III *Org. Lett.* **2004**, *6*, 3541. (e) Suri, J. T.; Ramachary, D. B.; Barbas, C. F., III *Org. Lett.* **2005**, *7*, 1383. (f) Suri, J. T.; Mitsumori, S.; Albertshofer, K.; Tanaka, F.; Barbas, C. F., III *J. Org. Chem.* **2006**, *71*, 3822. (g) Grondal, C.; Enders, D. *Tetrahedron* **2006**, *62*, 329. (h) For the reaction of 1,2-diketones and ketones, see Samanta, S.; Zhao, C.-G. *Tetrahedron Lett.* **2006**, *47*, 3383. (i) Other amino acid tetrazole derivatives have also been examined and found to compare favourably with their unmodified amino acid counterparts in the intermolecular aldol reaction process: Córdova, A.; Zou, W.; Dziedzic, P.; Ibrahim, I.; Reyes, E.; Xu, Y. *Chem.—Eur. J.* **2006**, *12*, 5383.
- Ward, D. E.; Jheengut, V.; Beye, G. E. *J. Org. Chem.* **2006**, *71*, 8989.
- (a) Seebach, D.; Boes, M.; Naef, R.; Schweizer, W. B. *J. Am. Chem. Soc.* **1983**, *105*, 5390. (b) List, B.; Hoang, L.; Martin, H. J. *Proc. Natl. Acad. Sci. U.S.A.* **2004**, *101*, 5839. (c) For an alternative view, see Seebach, D.; Beck, A. K.; Badine, D. M.; Limbach, M.; Eschenmoser, A.; Treasurywala, A. M.; Hobi, R.; Prikoszovich, W.; Linder, B. *Helv. Chim. Acta* **2007**, *90*, 425. For selected publications discussing alternative aspects, see: (d) Iwamura, H.; Wells, D. H., Jr.; Mathew, S. P.; Klussmann, M.; Armstrong, A.; Blackmond, D. G. *J. Am. Chem. Soc.* **2004**, *126*, 16312. (e) Iwamura, H.; Mathew, S. P.; Blackmond, D. G. *J. Am. Chem. Soc.* **2004**, *126*, 11770.
- A Density Functional Theory (DFT) study has since supported these findings: Arnó, M.; Zaragoza, R. J.; Domingo, L. R. *Tetrahedron: Asymmetry* **2005**, *16*, 2764.
- (a) Kiehlmann, E.; Loo, P.-W.; Menon, B. C.; McGillivray, N. *Can. J. Chem.* **1971**, *49*, 2964. (b) Hatch, C. E., III; Baum, J. S.; Takashima, T.; Kondo, K. *J. Org. Chem.* **1980**, *45*, 3281. (c) Benner, J. P.; Gill, G. B.; Parrot, S. J.; Wallace, B. *J. Chem. Soc., Perkin Trans. 1* **1984**, 331. (d) Wyenberg, H.; Staring, E. G. *J. Chem. Soc., Chem. Commun.* **1984**, 1181. (e) Muljiani, Z.; Gadre, S. R.; Modak, S.; Pathan, N.; Mitra, R. B. *Tetrahedron: Asymmetry* **1991**, *2*, 239. (f) Song, C. E.; Lee, J. K.; Lee, S.

- H.; Lee, S. *Tetrahedron: Asymmetry* **1995**, *6*, 1063. (g) Donohoe, T. J.; Guyo, P. M. *J. Org. Chem.* **1996**, *61*, 7664. (h) Fujisawa, T.; Ito, T.; Fujimoto, K.; Shimizu, M.; Wynberg, H.; Staring, E. G. *J. Tetrahedron Lett.* **1997**, *38*, 1593. (i) Fujisawa, T.; Ito, T.; Nishiura, S.; Shimizu, M. *Tetrahedron Lett.* **1998**, *39*, 9735. (j) Corey, E. J.; Link, J. O.; Shao, Y. *Tetrahedron Lett.* **1992**, *33*, 3435. (k) Corey, E. J.; Link, J. O. *J. Am. Chem. Soc.* **1992**, *114*, 1906. (l) Corey, E. J.; Helal, C. J. *Tetrahedron Lett.* **1993**, *34*, 5227. (m) Tennyson, R. L.; Cortez, G. S.; Galicia, H. J.; Kreiman, C. R.; Thompson, C. M.; Romo, D. *Org. Lett.* **2002**, *4*, 533.
- (17) (a) Ward, D. E.; Sales, M.; Man, C. C.; Shen, J.; Sasmal, P. K.; Guo, C. *J. Org. Chem.* **2002**, *67*, 1618. (b) Ward, D. E.; Sales, M.; Sasmal, P. K. *J. Org. Chem.* **2004**, *69*, 4808.
- (18) (a) Williams, R. M.; Hendrix, J. A. *Chem. Rev.* **1992**, *92*, 889. (b) Williams, R. M. *Aldrichimica Acta* **1992**, *25*, 11. (c) Duthaler, R. O. *Tetrahedron* **1994**, *50*, 1539. (d) Arend, M. *Angew. Chem., Int. Ed.* **1999**, *38*, 2873.
- (19) (a) Juaristi, E.; Quintana, D.; Escalante, J. *Aldrichimica Acta* **1994**, *27*, 3. (b) Cole, D. C. *Tetrahedron* **1994**, *50*, 9517.
- (20) Reedijk, J. *Chem. Commun.* **1996**, 801.
- (21) Lucet, D.; Le Gall, T.; Mioskowski, C. *Angew. Chem., Int. Ed.* **1998**, *37*, 2580.
- (22) Whitesell, J. K. *Chem. Rev.* **1989**, *89*, 1581.
- (23) (a) Bernardi, L.; Gothelf, A. S.; Hazell, R. G.; Jørgensen, K. A. *J. Org. Chem.* **2003**, *68*, 2583. (b) Davis, F. A.; Deng, J. *Org. Lett.* **2004**, *6*, 2789. (c) Viso, A.; de la Pradilla, R. F.; López-Rodríguez, M. L.; García, A.; Flores, A.; Alonso, M. *J. Org. Chem.* **2004**, *69*, 1542. (d) Ooi, T.; Kameda, M.; Fujii, J.; Maruoka, K. *Org. Lett.* **2004**, *6*, 2397.
- (24) Annunziata, R.; Benaglia, M.; Caporale, M.; Raimondi, L. *Tetrahedron: Asymmetry* **2002**, *13*, 2727.
- (25) Córdova, A.; Notz, W.; Zhong, G.; Betancort, J. M.; Barbas, C. F., III *J. Am. Chem. Soc.* **2002**, *124*, 1842.
- (26) Cobb, A. J. A.; Shaw, D. M.; Longbottom, D. A.; Gold, J. B.; Ley, S. V. *Org. Biomol. Chem.* **2005**, *3*, 84.
- (27) Loadings of 5% were generally used for operational simplicity.
- (28) Chowdari, N. S.; Ahmad, M.; Albertshofer, K.; Tanaka, F.; Barbas, C. F., III *Org. Lett.* **2006**, *8*, 2839.
- (29) Demko, Z. P.; Sharpless, K. B. *Angew. Chem., Int. Ed.* **2002**, *41*, 2113.
- (30) Cobb, A. J. A.; Longbottom, D. A.; Shaw, D. M.; Ley, S. V. *Chem. Commun.* **2004**, 1808.
- (31) Mitchell, C. E. T.; Brenner, S. E.; Ley, S. V. *Chem. Commun.* **2005**, 5346.
- (32) Mitchell, C. E. T.; Brenner, S. E.; García-Fortanet, J.; Ley, S. V. *Org. Biomol. Chem.* **2006**, *4*, 2039.
- (33) Rahbek Knudsen, K.; Mitchell, C. E. T.; Ley, S. V. *Chem. Commun.* **2006**, 66.
- (34) Ballini, R.; Bosica, G.; Fiorini, D.; Palmieri, A.; Petrini, M. *Chem. Rev.* **2005**, *105*, 933.
- (35) (a) List, B.; Pojarliev, P.; Martin, H. *J. Org. Lett.* **2001**, *3*, 2423. (b) Enders, D.; Seki, A. *Synlett* **2002**, 26.
- (36) Mitchell, C. E. T.; Cobb, A. J. A.; Ley, S. V. *Synlett* **2005**, 611.
- (37) For recent related publications demonstrating the utility of alternative catalytic species, see: (a) Andrey, O.; Alexakis, A.; Tomassini, A.; Bernardinelli, G. *Adv. Synth. Catal.* **2004**, *346*, 1147. (b) Betancort, J. M.; Sakthivel, K.; Thayumanavan, R.; Tanaka, F.; Barbas, C. F., III *Synthesis* **2004**, 1509. (c) Ishii, T.; Fujioka, S.; Sekiguchi, Y.; Kotsuki, H. *J. Am. Chem. Soc.* **2004**, *126*, 9558. (d) Terakado, D.; Takano, M.; Oriyama, T. *Chem. Lett.* **2005**, *34*, 962. (e) Enders, D.; Chow, S. *Eur. J. Org. Chem.* **2006**, 4578. (f) Mase, N.; Watanabe, K.; Yoda, H.; Takabe, K.; Tanaka, F.; Barbas, C. F., III *J. Am. Chem. Soc.* **2006**, *128*, 4966. (g) Tsogoeva, S. B.; Wei, S. *Chem. Commun.* **2006**, 1451. (h) Xu, Y.; Zou, W.; Sundén, H.; Ibrahem, I.; Córdova, A. *Adv. Synth. Catal.* **2006**, *348*, 418. (i) Yalalov, D. A.; Tsogoeva, S. B.; Schmatz, S. *Adv. Synth. Catal.* **2006**, *348*, 826. (j) Zhu, R.; Zhang, D.; Wu, J.; Liu, C. *Tetrahedron: Asymmetry* **2006**, *17*, 1611.
- (38) For previous reports on this reaction, see: (a) Yamaguchi, M.; Shiraishi, T.; Igarashi, Y.; Hirama, M. *Tetrahedron Lett.* **1994**, *35*, 8233. (b) Yamaguchi, M.; Igarashi, Y.; Reddy, R. S.; Shiraishi, T.; Hirama, M. *Tetrahedron* **1997**, *53*, 11223. (c) Hanessian, S.; Pham, V. *Org. Lett.* **2000**, *2*, 2975. (d) Halland, N.; Hazell, R. G.; Jørgensen, K. A. *J. Org. Chem.* **2002**, *67*, 8331. For the latest improvements to this organocatalytic process, see: (e) Prieto, A.; Halland, N.; Jørgensen, K. A. *Org. Lett.* **2005**, *7*, 3897. (f) Hanessian, S.; Shao, Z.; Warrior, J. S. *Org. Lett.* **2006**, *8*, 4787.
- (39) For use of other catalytic systems in this reaction, see: (a) Kawara, A.; Taguchi, T. *Tetrahedron Lett.* **1994**, *35*, 8805. (b) Yamaguchi, M.; Shiraishi, T.; Hirama, M. *Angew. Chem., Int. Ed. Engl.* **1993**, *32*, 1176. (c) Yamaguchi, M.; Shiraishi, T.; Hirama, M. *J. Org. Chem.* **1996**, *61*, 3520. (d) Halland, N.; Aburel, P. S.; Jørgensen, K. A. *Angew. Chem., Int. Ed.* **2003**, *42*, 661. (e) Ooi, T.; Ohara, D.; Fukumoto, K.; Maruoka, K. *Org. Lett.* **2005**, *7*, 3195.
- (40) (a) Krapcho, A. P. *Synthesis* **1982**, 805. (b) Krapcho, A. P. *Synthesis* **1982**, 893.
- (41) Wascholowski, V.; Rahbek Knudsen, K.; Mitchell, C. E. T.; Ley, S. V. University of Cambridge, Cambridge, U.K. Unpublished work, 2007.
- (42) Reissig, H.-U.; Zimmer, R. *Chem. Rev.* **2003**, *103*, 1151.
- (43) (a) Corey, E. J.; Achiwa, K.; Katzenellenbogen, J. A. *J. Am. Chem. Soc.* **1969**, *91*, 4318. (b) Higgs, M. D.; Mulheirn, L. J. *Tetrahedron* **1981**, *37*, 4259. (c) Paul, V. J.; Fenical, W. *Science* **1983**, *221*, 747. (d) Kerr, R. G.; Baker, B. J. *Nat. Prod. Rep.* **1991**, *8*, 465. (e) Williams, R. M.; Fegley, G. J. *J. Am. Chem. Soc.* **1991**, *113*, 8796. (f) Donaldson, W. A. *Tetrahedron* **2001**, *57*, 8589. (g) Faust, R. *Angew. Chem., Int. Ed.* **2001**, *40*, 2251. (h) Chakraborty, T. K.; Reddy, V. R. *Tetrahedron Lett.* **2006**, *47*, 2099. (i) Yakambram, P.; Puranik, V. G.; Gurjar, M. K. *Tetrahedron Lett.* **2006**, *47*, 3781.
- (44) Reichelt, A.; Martin, S. F. *Acc. Chem. Res.* **2006**, *39*, 433.
- (45) For a selection of recent publications in this area, see: (a) Johansson, C. C. C.; Bremeyer, N.; Ley, S. V.; Owen, D. R.; Smith, S. C.; Gaunt, M. J. *Angew. Chem., Int. Ed.* **2006**, *45*, 6024. (b) Du, H.; Long, J.; Shi, Y. *Org. Lett.* **2006**, *8*, 2827. (c) Itagaki, M.; Masumoto, K.; Suenobu, K.; Yamamoto, Y. *Org. Proc. Res. Dev.* **2006**, *10*, 245. (d) McAllister, G. D.; Oswald, M. F.; Paxton, R. J.; Raw, S. A.; Taylor, R. J. K. *Tetrahedron* **2006**, *62*, 6681. (e) Mekonnen, A.; Carlson, R. *Synthesis* **2006**, 1657. (f) Werner, H.; Herrerias, C. I.; Glos, M.; Gissibl, A.; Fraile, J. M.; Pérez, I.; Mayoral, J. A.; Reiser, O. *Adv. Synth. Catal.* **2006**, *348*, 125. (g) Lebel, H.; Marcoux, J.-F.; Molinaro, C.; Charette, A. B. *Chem. Rev.* **2003**, *103*, 977.
- (46) Rosini, G.; Ballini, R. *Synthesis* **1988**, 833.
- (47) For a selection of publications in this area, see: (a) Kocór, M.; Kroszczyński, W. *Synthesis* **1976**, 813. (b) Russell, G. A.; Makosza, M.; Hershberger, J. J. *J. Org. Chem.* **1979**, *44*, 1195. (c) O'Bannon, P. E.; Dailey, W. P. *Tetrahedron* **1990**, *46*, 7341. (d) Yu, J.; Falck, J. R.; Mioskowski, C. *J. Org. Chem.* **1992**, *57*, 3757. (e) Kumaran, G.; Kulkarni, G. H. *Synthesis* **1995**, 1545. (f) Hübner, J.; Liebscher, J.; Pätzelt, M. *Tetrahedron* **2002**, *58*, 10485. (g) Wurz, R. P.; Charette, A. B. *J. Org. Chem.* **2004**, *69*, 1262.
- (48) (a) Arai, S.; Nakayama, K.; Ishida, T.; Shioiri, T. *Tetrahedron Lett.* **1999**, *40*, 4215. (b) McCooley, S. H.; McCabe, T.; Cannon, S. J. *J. Org. Chem.* **2006**, *71*, 7494.
- (49) Hansen, H. M.; Longbottom, D. A.; Ley, S. V. *Chem. Commun.* **2006**, 4838.
- (50) Wascholowski, V.; Hansen, H. M.; Longbottom, D. A.; Ley, S. V. *Synthesis* **2008**, in press.

- (51) Davis, F. A.; Chen, B.-C. *Chem. Rev.* **1992**, *92*, 919
- (52) For use of alternative oxidant species, see: (a) Sundén, H.; Engqvist, M.; Casas, J.; Ibrahim, I.; Córdova, A. *Angew. Chem., Int. Ed.* **2004**, *43*, 6532. (b) Engqvist, M.; Casas, J.; Sundén, H.; Ibrahim, I.; Córdova, A. *Tetrahedron Lett.* **2005**, *46*, 2053.
- (53) For prior work in this area using L-proline (**1**) as catalyst, see: (a) Brown, S. P.; Brochu, M. P.; Sinz, C. J.; MacMillan, D. W. C. *J. Am. Chem. Soc.* **2003**, *125*, 10808. (b) Zhong, G. *Angew. Chem., Int. Ed.* **2003**, *42*, 4247. (c) Hayashi, Y.; Yamaguchi, J.; Hibino, K.; Shoji, M. *Tetrahedron Lett.* **2003**, *44*, 8293. (d) Hayashi, Y.; Yamaguchi, J.; Sumiya, T.; Hibino, K.; Shoji, M. *J. Org. Chem.* **2004**, *69*, 5966. (e) Bøgevig, A.; Sundén, H.; Córdova, A. *Angew. Chem., Int. Ed.* **2004**, *43*, 1109.
- (54) Kim, S.-G.; Park, T.-H. *Tetrahedron Lett.* **2006**, *47*, 9067.
- (55) Chowdari, N. S.; Barbas, C. F., III *Org. Lett.* **2005**, *7*, 867.
- (56) For a recent minireview on asymmetric organocatalytic domino reactions, see Enders, D.; Grondal, C.; Hüttel, M. R. M. *Angew. Chem., Int. Ed.* **2007**, *46*, 1570.
- (57) (a) Khosla, C. *Chem. Rev.* **1997**, *97*, 2577. (b) Khosla, C.; Gokhale, R. S.; Jacobsen, J. R.; Cane, D. E. *Annu. Rev. Biochem.* **1999**, *68*, 219. (c) Katz, L. *Chem. Rev.* **1997**, *97*, 2557. (d) Staunton, J.; Weissman, K. J. *Nat. Prod. Rep.* **2001**, *18*, 380.
- (58) (a) Yamamoto, Y.; Momiyama, N.; Yamamoto, H. *J. Am. Chem. Soc.* **2004**, *126*, 5962. (b) Momiyama, N.; Yamamoto, Y.; Yamamoto, H. *J. Am. Chem. Soc.* **2007**, *129*, 1190.
- (59) (a) Kumarn, S.; Shaw, D. M.; Longbottom, D. A.; Ley, S. V. *Org. Lett.* **2005**, *7*, 4189. (b) Kumarn, S.; Shaw, D. M.; Ley, S. V. *Chem. Commun.* **2006**, 3211. (c) Kumarn, S.; Oelke, A. J.; Shaw, D. M.; Longbottom, D. A.; Ley, S. V. *Org. Biomol. Chem.* **2007**, *5*, 2678.
- (60) Oelke, A. J.; Kumarn, S.; Longbottom, D. A.; Ley, S. V. *Synlett* **2006**, 2548.
- (61) Oelke, A. J.; Knauer, S.; Ley, S. V. Towards the Synthesis of Chloptosin. Presented at the International Symposium on Advances in Synthetic and Medicinal Chemistry (ASMC 07), St. Petersburg, Russia, August 27–31, 2007; Poster No. 141

About the Authors


Deborah A. Longbottom received her undergraduate degree from the University of Durham in 1997 and, following a year working in the pharmaceutical industry, came to Cambridge to carry out her Ph.D. work under the guidance of Professor Steven Ley. In 2002, she joined Professor K. C. Nicolaou's research group at The Scripps Research Institute, San Diego, California, as a postdoctoral associate and returned to Professor Ley's Group early in 2004. Her research interests have encompassed both natural product synthesis (e.g., polyenoyltetramic acids and depsipeptides) and method development (e.g., novel uses of the Burgess reagent and organocatalytic methodologies). Currently, she is a senior research associate in Professor Ley's group, and concurrently holds teaching fellowships in both the Department of Chemistry and Trinity College, Cambridge.

Vilius Franckevičius was born in 1983 in Kaunas, Lithuania. He studied Natural Sciences at the University of Cambridge, where he undertook his final-year project on the development of new organocatalysts under the supervision of Professor Steven Ley, and subsequently obtained his M.Sc. degree in Natural Sciences (chemistry) in 2005 (Fitzwilliam College). He is currently a Ph.D. student in Ley's research group, where he is involved in the application of organocatalytic methodology in natural product synthesis.

Sirirat Kumarn was born in Sukhothai, Thailand. She received her undergraduate degrees in Natural Sciences (chemistry) in 2004 from St. Catharine's College, University of Cambridge. She is currently a Ph.D. student in Professor Ley's research group, where she is working on the development of an organocatalytic route to enantiopure 1,2-oxazines and its applications to natural product synthesis.

Alexander J. Oelke was born in 1980 in Reinbek, Germany. He studied chemistry at the University of Hamburg, where he obtained his Diploma in 2006 under the supervision of Professor Chris Meier, and in collaboration with Professor Steven Ley, for the development of an organocatalytic tandem procedure for the synthesis of chiral pyridazine derivatives. He is currently a Ph.D. student in Ley's group at the University of Cambridge, where he is involved in the application of organocatalytic methodology in natural product synthesis.

Veit Wascholowski was born in 1975 in Braunschweig, Germany. He studied chemistry at the University of Karlsruhe, Germany, and completed his Diploma in 2000. He received his Ph.D. degree in 2006 from the University of Leipzig, Germany, where he worked under the guidance of Professor Athanassios Giannis in the field of chemical biology, which involved the synthesis and biological evaluation of natural products and their analogues. In 2006, he joined Professor Ley's research group at the University of Cambridge as a postdoctoral research associate, where he is currently working on the development of new organocatalytic reactions and their application in the total synthesis of natural products.

Steven V. Ley received his Ph.D. degree from Loughborough University in 1972, after which he carried out postdoctoral research with Professor Leo Paquette at Ohio State University, and then with Professor Derek Barton at Imperial College, London. In 1975, he joined that Department as a lecturer and became Head of Department in 1989. In 1992, he moved to take up the 1702 BP Chair of Organic Chemistry at the University of Cambridge, and became a Fellow of Trinity College. He was elected to the Royal Society in 1990 and, between 2000 and 2002, was President of the Royal Society of Chemistry (RSC). In addition, Steve was made a Commander of the British Empire (CBE) early in 2002, and has been the recipient of many prizes and awards for his creative work and innovative solutions in the art of organic synthesis. Among the most recent of these are the Yamada-Koga Prize (2005), the Nagoya Gold Medal (2006), the ACS Award for Creative Work in Synthetic Organic Chemistry (2007), and the Paul Karrer Medal (2007).

2008 ACS Award Recipients

Aldrich, a proud sponsor of three ACS awards, congratulates the following recipients for their outstanding contributions to chemistry.

ACS Award for Creative Work in Synthetic Organic Chemistry
Professor Masakatsu Shibasaki
University of Tokyo

ACS Award in Inorganic Chemistry
Professor Kenneth N. Raymond
University of California, Berkeley

Herbert C. Brown Award for Creative Research in Synthetic Methods
Professor Eric N. Jacobsen
Harvard University


Congratulations to each and all!

Accelerate Catalysis


Buchwald Ligands

The Pd-catalyzed C–N bond formation has become an important synthetic reaction in the past 20 years. Buchwald and co-workers have been very active in synthesizing and developing a portfolio of phosphine ligands for this transformation and other cross-coupling reactions. The ligands chosen are based on a biaryl skeleton with a phosphorus moiety at the 2 position


of one aromatic ring and another moiety on the other aromatic ring. These ligands are very stable and active in a variety of cross-coupling reactions such as carbon–carbon, carbon–nitrogen, and carbon–oxygen coupling. Sigma-Aldrich is pleased to offer the following portfolio of Buchwald ligands.




(SPhos)
638072




(SPhos, water soluble)
677280




(RuPhos)
663131




(XPhos)
638064




(tBuXPhos)
638080




(Tetramethyl di-tBuXPhos)
675938




(MePhos)
695262



(tBuMePhos)
695211



(JohnPhos)
638439




(CyJohnPhos)
638099



(PhDavePhos)
695882



(tBuDavePhos)
695874



(DavePhos)
638021

Buchwald Ligands Kit I 659932

The kit contains:

638021	638072	638099
638064	638080	638439

Buchwald Ligands Kit II 687243

The kit contains:

638439	638080	638021
638099	663131	638064
675938	638072	677280

Buchwald Ligands Kit III 698903


The kit contains:

638072	675938	695882
677280	695262	695874
663131	695211	638021
638064	638439	
638080	638099	

HydraPhos Ligands

Hintermann and coworkers introduced a new set of ligands based on a pyridylphosphane backbone for the anti-Markovnikov hydration of terminal alkynes. When used with a ruthenium complex, high yields were reported for a variety of terminal alkynes.


Labonne, A. et al. *Org. Lett.* **2006**, 8, 5853



Dynamic Kinetic Resolution Catalysts


Dynamic Kinetic Resolution (DKR) catalysis is an essential methodology for the conversion of racemic substrates into single enantiomers. Kim et al. reported the (S)-selective DKR of a variety of alcohols by utilizing a combination of substilisin and an aminocyclopentadienyrruthenium complex. High yields and selectivities were observed for a variety of secondary alcohols.

Kim, M.-J. et al. *J. Am. Chem. Soc.* **2003**, 125, 11494



New Metal Precursors for Asymmetric Catalysis

Sigma-Aldrich is pleased to offer a comprehensive portfolio of rhodium and iridium BARF complexes for asymmetric transformations.



Presenting...

Sigma-Aldrich's *ChemBlogs*

An industry-first for open scientific discussion.




Visit chemblogs.com


Aminophosphine Ligands

Recently, there has been a growing interest in aminophosphine ligands for asymmetric synthesis. Researchers at Kanata Chemical Technologies, Inc., have synthesized several sets of aminophosphine ligands that show high reactivity and selectivity in a wide array of enantioselective reactions.


A growing area of application for aminophosphine ligands in asymmetric synthesis is in ruthenium-catalyzed hydrogenations. Chen et al. have described the use of ferrocenylaminophosphines in the ruthenium-catalyzed asymmetric hydrogenation of acetophenone. Using the precatalyst $[\text{Ru}(\text{C}_6\text{H}_5)_2\text{Cl}_2]$ and the ferrocenyl-based aminophosphine ligand, these researchers found that the hydrogenation proceeded efficiently with reasonable enantioselectivity. Aldrich is pleased to offer a portfolio of aminophosphine ligands and complexes.




Chen, W.; Mbafor, W.; Roberts, S. M.; Whittall, J. *Tetrahedron: Asymmetry* **2006**, *17*, 1161.




664979




664987




697583




697354




697362




696943




697370




697389




697397




697419



697427



43162



43163

Aminophosphine Catalysts in Modern Asymmetric Synthesis



Dr. Dino Amoroso



Dr. Todd W. Graham



Dr. Rongwei Guo

Dino Amoroso, Todd W. Graham, Rongwei Guo, Chi-Wing Tsang, and Kamaluddin Abdur-Rashid*

*Kanata Chemical Technologies, Inc.
MaRS Centre, South Tower
230-101 College Street
Toronto, ON M5G 1L7, Canada
Email: chemistry@kctchem.com*



Dr. Chi-Wing Tsang



Dr. K. Abdur-Rashid

Outline

1. Introduction
2. Ligand Synthesis
3. Hydrogenation
 - 3.1. Ruthenium Catalysts
 - 3.2. Rhodium Catalysts
 - 3.3. Iridium Catalysts
4. Allylic Alkylation
5. Hydroformylation
6. Conjugate Additions
 - 6.1. Asymmetric Michael Addition to Enones
 - 6.2. Asymmetric Addition of Organolithiums to Aldehydes
7. Cycloaddition Reactions
8. Conclusions
9. Acknowledgment
10. References

1. Introduction

The importance of ligand composition and structure in transition-metal-catalyzed asymmetric synthesis cannot be overstated. From the simplest lock-and-key model to the most complex transition state, the interaction between catalyst and substrate can be completely dictated by the chemical composition and spatial orientation of the supporting ligands. An excellent example of this is in the direct hydrogenation of ketones, aldehydes, and imines catalyzed by ruthenium complexes supported by phosphine or amine ligands.¹ In this process, the unsaturated substrate does not bind directly to the metal center, but rather interacts simultaneously with the Ru–H and N–H bonds of an amine-containing ligand (**Figure 1, Part A**). Often generically described as metal–ligand bifunctional catalysis, the importance of the ligand composition (e.g., incorporation of an

N–H bond) and structure (e.g., the spatial orientation of the N–H bond and the other spectator ligands that direct the approach of substrate and often govern stereochemistry) in this reaction manifold is clear.

A particular class of ligand, which is often involved in such ligand-dependent interactions, is chelating aminophosphines (**Figure 1, Part B**). The nature of the substituents at nitrogen, and the stereochemistry at phosphorus and in the ligand backbone, render this motif particularly versatile in catalysis. Because of the highly modular nature of this ligand type, it has found application in a broad range of asymmetric transformations and has become an invaluable tool for the preparation of chiral molecules. In this review, a subset of this ligand class, defined by restricting at least R¹ or R² to H, is considered. This class is of particular interest owing to the potential involvement of this functionality in catalysis. A further restriction that excludes ligands incorporating a direct P–N bond has also been imposed. However, in select instances—such as in the Rh-catalyzed hydrogenation where the direct P–N bond motif is almost exclusively employed, or in cases where ligands are readily derived from those which *are* included in the preceding subset—both restrictions have been overlooked. The synthesis of chelating aminophosphine ligands that results from imposing these two restrictions, and their application in asymmetric synthesis over the last 10 years, are reviewed.

2. Ligand Synthesis

The growth in popularity of aminophosphine ligands in asymmetric synthesis is due in part to the increasing number of convenient synthetic pathways leading to useful ligand sets. In recent years, several general routes have been described, which allow access to a broad range of versatile aminophosphines.

Amino acids constitute a convenient class of precursors to chiral aminophosphine ligands.^{2,3} Morimoto and Achiwa have described the use of L-valine (**1a**, **Scheme 1**) and other amino acids in the preparation of aminophosphine ligands of type **2**. These ligands have widespread applications in catalysis, particularly in hydrogenation, while other derivatives of L-valine have been exploited in palladium-catalyzed allylic transformations (vide infra).

Another common route to chiral aminophosphine ligands is through commercially available chiral amino alcohols such as ephedrine, norephedrine, and pseudoephedrine. Dahlenburg and Götz have reported the synthesis of chiral aminophosphines by the aziridination of amino alcohols.⁴ Ring-opening of the aziridines by nucleophilic attack with diphenylphosphine affords chiral ligands **3–5** (**Figure 2**). The appeal of this route is in its ability to dictate the stereochemistry of the ring-opened aminophosphine by controlling that of the aziridine (via a judicious choice of the aziridination protocol). Indeed the ring opening of aziridines is a convenient route to a range of chiral aminophosphines (**eq 1**).⁵ Yudin and co-workers employed secondary phosphines as nucleophiles to ring-open cyclohexene aziridines to cyclohexylaminophosphines,





Figure 1. (A) Interaction Between Ligands and Substrate in the Ruthenium-Catalyzed Hydrogenation of Ketones, Aldehydes, and Imines. (B) General Structure of a Chelating Aminophosphine Ligand.



Scheme 1. The Synthesis of Aminophosphines from Amino Acids.





Figure 2. Aminophosphines Derived from Commercially Available Amino Alcohols.



which were optically resolved with D-tartaric acid. It is worth noting that pyrazole derivatives were also prepared from the cyclohexylaminophosphines.

An alternate route employing amino alcohols has been disclosed by our group (**Scheme 2**)⁶ and others.⁷ The route involves the formation of cyclic sulfamidates from the corresponding N-protected sulfamidites. Treatment of the sulfamidate with a metal phosphide, followed by removal of sulfate with dilute acid and N-deprotection, yield the chiral or achiral aminophosphine. The route described by Hilmersson and co-workers differs in that they include an N-alkylation step prior to sulfamidite formation. This allows them to proceed without protecting the amine. Some representative examples of the ligands prepared by our method are shown in **Figure 3**.

We have also described a simple route to aminophosphines via haloalkylammonium salts (**Scheme 3**).⁶ Many haloalkylammonium salts are commercially available, although they can also be readily prepared from amino alcohols. The procedure involves neutralization of the salt and protection of the amine, followed by halide substitution with a metal phosphide. Hydrolysis then leads to the desired aminophosphine ligand. Representative examples are depicted in **Figure 4**.


While chiral amino alcohols and acids represent a convenient source of chirality for ligand construction, so too does the 1-ferrocenylalkyl fragment that has been exploited for the development of chiral ferrocenylaminophosphines by Boaz⁸ and Chen.⁹ Boaz prepared chiral ferrocenylaminophosphines of type **6**, which were subsequently derivatized into phosphinoferoceylaminophosphines and extensively studied in asymmetric synthesis (**Scheme 4**).⁸ Chen synthesized similar compounds by a modified procedure, which provides entry into a range of P-chiral ligands with nonidentical substituents on phosphorus. An attractive feature of all such ferrocenylaminophosphines is their remarkable stability toward air oxidation, as samples of material exposed to air for up to 3 years showed no loss in enantioselectivity or activity in Rh-catalyzed hydrogenations.^{8d}

A library of compounds aimed at elucidating ligand structural effects in the asymmetric transfer hydrogenation of prochiral ketones was described by Jubault and co-workers.¹⁰ This group employed two different routes to arrive at intermediate amides that were subsequently reduced and deprotected (**Scheme 5**). Both the coupling and reduction of the amides proceeded in moderate-to-high yields, while the removal of borane took place quantitatively.

A simple route to enantiopure β-aminophosphines through vinylphosphine oxides (**Scheme 6**) was employed by Maj et al.¹¹ While these workers described several variants of aminophosphine oxides and their use as ligands in transfer hydrogenation,^{11a} aminophosphines **7b** and **8b** were also prepared and tested vis-à-vis their oxidized precursors.^{11b}

Hii and co-workers have employed a similar methodology to prepare other N-alkylaminophosphine (and aminodiphosphine) ligands from vinylphosphine oxides, and reported on their use in ruthenium-catalyzed transfer hydrogenation.¹² These workers also described the development of a range of chiral aminophosphine ligands (**Figure 5**)^{12c,d} subsequent to their initial disclosure of several achiral variants.^{12a,b}

A strong base containing a guanidine functional group was utilized by Fu et al. to catalyze the phospho-Michael reaction between a number of diarylphosphine oxides and various aryl-substituted nitroalkenes. Subsequent reduction of the nitro and phosphine oxide groups led to the corresponding aminophosphines in good yields and >99% ee's (**Scheme 7**).¹³



Scheme 2. The Sulfamidate Route to Aminophosphines.





Figure 3. Representative Chiral Aminophosphines Prepared by the Sulfamidate Route.



Scheme 3. The Haloalkylammonium Route to Aminophosphines.






Figure 4. Representative Aminophosphines Prepared from Haloalkylammonium Salts.



Scheme 4. Aminophosphines Prepared from Ferrocene Derivatives.



Scheme 5. Aminophosphines Derived from Amides.





Scheme 6. β -Aminophosphines Derived from Vinylphosphine Oxides.

Figure 5. Chiral Aminophosphines Derived from Vinylphosphine Oxides.




Scheme 7. Synthesis of Aminophosphines Catalyzed by the Guanidine Functional Group.




Ref. 17

Scheme 8. Industrially Significant Alcohols Prepared by the Ruthenium-Aminophosphine-Catalyzed Hydrogenation.




Ref. 17

Figure 6. Representative Ketones and Imines that Have Been Reduced to the Corresponding Alcohols and Amines by the Ruthenium-Aminophosphine-Catalyzed Hydrogenation.




Ref. 18



Figure 7. Ruthenium-Aminophosphine Complexes Employed in the Hydrogenation and Transfer Hydrogenation of Ketones.



15	Time	Conv.	R
a	3 h	100%	67% (R)
b	3 h	88%	44% (R)
c	3 h	100%	79% (R)
d	4 h	4%	21% (R)
e	3 h	100%	67% (S)
f	3 h	56%	16% (S)



15	Ar	R	R'
a	Ph	Me	H
b	Ph	Ph	H
c	3,5-Me ₂ C ₆ H ₃	Me	H
d	Ph	Me	Bn



Ref. 9b

3. Hydrogenation

3.1. Ruthenium Catalysts

A heavily exploited application area for aminophosphine ligands in asymmetric synthesis is the ruthenium-catalyzed hydrogenation. This process is integral to the preparation of alcohols and amines that are useful in the flavor and fragrance, pharmaceutical, agrochemical, materials, and fine chemicals industries.¹⁴ Our group and others have reported extensively on the use of ruthenium aminophosphine complexes, and has studied the relationship between catalyst structure and enantioselectivity.^{3,15,16,17} The industrially relevant compounds (*E*)-2-ethyl-4-(2,2,3-trimethylcyclopent-3-en-1-yl) but-2-en-1-ol (**10**) and *cis*-4-*tert*-butylcyclohexanol (**11**) are two examples of the type of product that can be very efficiently produced from the corresponding aldehyde or ketone by using ruthenium aminophosphine catalysts (**Scheme 8**).¹⁷ With precatalysts of the type RuCl₂(aminophosphine)₂ or RuCl₂(diphosphine)(aminophosphine), substrate:catalyst ratios of 100,000:1 or greater are typical in the direct hydrogenation. This approach was applied to a broad range of ketones, aldehydes, and imines (**Figure 6**).¹⁷ The advantage of a ruthenium-catalyzed hydrogenation over a conventional stoichiometric reduction with a hydride-transfer reagent (e.g., alkali metal borohydrides or aluminum hydrides) is quite apparent: avoid the use of large quantities of expensive and difficult-to-handle materials that produce inorganic hydroxides in favor of utilizing catalytic amounts of robust materials. It should also be pointed out that these ruthenium catalysts selectively reduce only the ketone, aldehyde, or imine while leaving carbon-carbon double bonds intact.

Dahlenburg and Kühnlein have described the use of ruthenium-aminophosphine complexes **12–14** (**Figure 7**) in the transfer hydrogenation and, in the case of **14**, the direct hydrogenation of acetophenone.¹⁸ They observed an induction period in the transfer-hydrogenation experiments corresponding to the time needed for the active catalyst to arise from the precatalyst complexes. Variations in induction times were observed and found to correlate directly to the extent of steric shielding of the amino group. Catalyst **14** was effective in the direct hydrogenation of acetophenone to 1-phenylethanol as well. It was shown, through the use of isotopically labeled solvent, that, under the direct hydrogenation conditions, the source of hydrogen atoms in the alcohol product was indeed hydrogen gas.

Chen and co-workers have described the use of ferrocenylaminophosphines in the ruthenium-catalyzed asymmetric (direct) hydrogenation of 1-acetonaphthone (**eq 2**).^{9b} In the presence of precatalysts of type RuCl₂(**15**)₂, derived from the reaction of [RuCl₂(C₆H₆)₂] and **15**, the hydrogenation proceeded rapidly and with reasonable enantioselectivity when ligands **15a**, **15c**, and **15e** were used. When the steric bulk of the phosphine aryl substituent was largest (**15c**), enantioselectivity was highest. In contrast, increasing the steric bulk at or near nitrogen resulted in significantly diminished activity and selectivity (e.g., **15d** and **15f**). This was anticipated as increased steric bulk at nitrogen should limit the ability of the substrate to interact with the N-H bond. The carbon-centered chirality was found to dictate the sense of induction in the product.

Using aminophosphines derived from vinylphosphine oxides (see Scheme 6), a comparison of the catalytic activities between reduced and oxidized β-aminophosphines has been reported.¹¹ In the transfer hydrogenation of various aromatic ketones, the reduced aminophosphine ligands had activities comparable to those of the corresponding aminophosphine oxide ligands, but gave rise to consistently higher ee's. Hii and co-workers have also used similarly derived β-aminophosphines (and their oxides) in the

eq 2

transfer hydrogenation of ketones (eq 3).¹² They initially described the use of achiral ligands (**17a–d**) and reported that increasing the steric bulk of the N-substituent led to diminished activity. Upon examination of chiral ligands **9a–e** (see Figure 5), they found that incorporation of the alcohol functionality in the ligand resulted in dramatically improved enantioselectivity. They obtained a 79% ee (*R* form) in the hydrogenation of acetophenone using ligand **9e**, compared to 39% ee (*R* form) with **9b** as the next best of the five chiral ligands tested in this transformation. It is worth noting that Hii's group also discovered that the optimal metal-to-ligand ratio (when $[\text{RuCl}_2(p\text{-cymene})_2]$ was used as the Ru source) was 1:1, in contrast to the 2:1 ratio traditionally employed in Ru-catalyzed transfer hydrogenations.

The library of aminophosphines developed by Jubault and co-workers (see Scheme 5) was tested in the transfer hydrogenation of acetophenone, propiophenone, and isobutyrophenone.¹⁰ These researchers found that acetophenone and propiophenone were most effectively reduced when $\text{R}^1 = \text{Et}$ (*S*), $\text{R}^2 = \text{Ph}$ (*S*), and $\text{R}^3 = \text{Me}$. For isobutyrophenone, the best results were obtained when $\text{R}^1 = \text{H}$, $\text{R}^2 = \text{Ph}$, and $\text{R}^3 = \text{CH}_2\text{OBz}$. Jubault's group was also able to glean the impact of the nature and position of the substituents on this particular class of ligands. For instance, the introduction of a chiral center adjacent to phosphorus has a dramatic effect on enantioselectivity, while the sense of induction is most strongly governed by the chiral center adjacent to the amine.

3.2. Rhodium Catalysts

When combined with rhodium, the versatile chiral phosphine–aminophosphine and phosphine–phosphoramidate ligands, derived by combining phosphorus and nitrogen functional groups in one molecular structure, represent a particularly important class of catalyst.

In 2002, Boaz first introduced the hybrid phosphine–aminophosphine ligands **18–22** (Figure 8), based on a chiral ferrocenylethylamine backbone, for the rhodium-catalyzed asymmetric hydrogenation of olefins.⁸ Ligands **18** and **19** showed excellent enantioselectivities in the hydrogenation of dehydro- α -amino acid and itaconic acid derivatives, while **22** showed very good enantioselectivity in the hydrogenation of α -keto esters (Scheme 9).^{8c} Ligand (S_C, R_{Fe})-**19** was successfully utilized for the preparation of single-enantiomer 2-naphthylalanine derivatives on multikilogram scale by Eastman Chemical Company (Scheme 10).¹⁹ The catalyst system $\text{Rh}-(S_C, R_{Fe})$ -**19** showed high activity and enantioselectivity, as the hydrogenation product **23** was obtained in 96% isolated yield and >99% ee after one crystallization from toluene–heptane. Further elaboration of **23** led to the final product in high yield and >99.5% ee.¹⁹

While Rh -BoPhoz catalysts show high activity and enantioselectivity in the hydrogenation of dehydro- α -amino acid derivatives, they result in only moderate enantioselectivity in the hydrogenation of α -aryl enamides (~80% ee). Chan introduced fluorinated phosphinoferochenylaminophosphine ligands **24** and **25**, which showed excellent enantioselectivities in the hydrogenation of dehydro- α -amino acid derivatives ($\leq 99.7\%$ ee) and α -aryl enamide substrates ($\leq 99.7\%$ ee) (eq 4).²⁰ A significant feature of this catalyst system is that the rhodium complexes are exceptionally air- and moisture-stable, even when dissolved in an organic solvent.²⁰

Hu, Zheng, and co-workers recently introduced a modified BoPhoz ligand with three chiral centers, **26**, for the highly enantioselective ($\leq 97\%$ ee) rhodium-catalyzed hydrogenation of γ -phthalimido-substituted α, β -unsaturated carboxylic acid esters (Scheme 11).²¹ The catalyst system Rh -**26** was successfully applied to the synthesis of the optically active pharmaceuticals (*R*)-baclofen,







Figure 8. BoPhoz Ligands Derived from Phosphinoferochenylethylamine.



Scheme 9. Hydrogenations with Rh -BoPhoz Ligands.



Scheme 10. Multikilogram Hydrogenation Process with the $\text{Rh}-(S_C, R_{Fe})$ -**19** Catalyst System.



Scheme 11. Application of Rh-26 Catalyst System to the Preparation of (R)-Baclofen and (R)-Rolipram.





Figure 9. Phosphine-Phosphoramidite Ligands.



Scheme 12. Rhodium-33 Hydrogenation of Dehydro- β -amino Acid Derivatives.

which is widely used as an antispasmodic agent, and (R)-rolipram, which is employed as an antidepressant and anti-inflammatory agent.

In 2004, Hu and Zheng employed a new set of phosphine-phosphoramidite ligands, **27–32** (Figure 9), as alternatives to BoPhoz ligands for the Rh-catalyzed hydrogenation of olefins.²² The enantioselectivity with the Rh-**27** catalyst system in the hydrogenation of *N*-(1-phenylethyl)acetamide was greater than 99% ee (S/C = 5,000/1). Rh-**27** also led to >99% ee for the hydrogenation of dimethyl itaconate and (*Z*)-acetamidocinnamate at low catalyst loadings (S/C = 10,000). The rhodium-catalyzed hydrogenation reactions employing these ligands were carried out under ambient conditions, taking no precautions to exclude air or moisture, with no loss in activity or enantioselectivity.²²

While **27** showed high enantioselectivities in the rhodium-catalyzed hydrogenation of enamides, itaconates, and dehydro- α -amino acid derivatives, very low enantioselectivities were observed for the rhodium-catalyzed hydrogenation of the *Z* and *E* isomers of dehydro- β -amino acid derivatives (20–60% ee's). After expanding the ligand scope, Zheng's group found that **33**, bearing a proton instead of a methyl group on nitrogen, showed excellent enantioselectivity in the rhodium-catalyzed hydrogenation of *Z* and *E* β -aryl- or β -alkyl- β -(acylamino)acrylates, leading to the two products with opposite configurations (Scheme 12).²³

Leitner, Faraone, and co-workers introduced chiral phosphine-phosphoramidite ligands **34** and **35** (*n*-Bu-QuinaPhos) for the rhodium-catalyzed hydrogenation of dimethyl itaconate and methyl 2-acetamidoacrylate (eq 5).²⁴ Excellent activity and enantioselectivity were observed for the reaction. Moreover, the rhodium-catalyzed hydroformylation of styrene with *n*-Bu-QuinaPhos gave rise to high regioselectivity and moderate enantioselectivity.²⁴

Me-AniPhos (**36**) and IndolPhos ligands **37** and **38** were recently reported by the groups of Kostas and Reek (Figure 10).^{25,26} These chiral phosphine-phosphoramidite ligands, derived from achiral aminophosphine ligands, also show high enantioselectivity in the rhodium-catalyzed hydrogenation of enamides and itaconate derivatives ($\leq 97\%$ ee for enamides and 98% ee for itaconates). Owing to the simple structure and the wide range of available aminophosphine precursors, **36–38** represent a highly versatile ligand class.

The novel phosphine-phosphoramidite ligands **39–42** (PEAPhos), derived from chiral α -phenylethylamine, and **43**, derived from 1,2,3,4-tetrahydro-1-naphthylamine, were disclosed by Zheng and co-workers.²⁷ The Rh-**39** catalyst system showed excellent enantioselectivity (>99% ee) in the hydrogenation of olefins (Scheme 13).²⁷ Ligand **43** was also successfully applied to the synthesis of α -hydroxyphosphonic acid derivatives by the rhodium-catalyzed hydrogenation (a significant achievement) of β -substituted α -acyloxyphosphonates. A greater than 99% ee was achieved for a range of substrates bearing β -aryl, β -alkoxy, and β -alkyl substituents (eq 6).²⁷

3.3. Iridium Catalysts

The class of aminophosphine ligands discussed so far has found only limited application in iridium-catalyzed hydrogenations. Dahlenburg and collaborators have employed aminophosphine ligands in the iridium-catalyzed hydrogenation of unsaturated substrates.^{4,28} They described a series of chiral and achiral aminophosphine-chelated iridium(I) complexes prepared by treating $[\text{Ir}(\text{cod})_2]\text{BF}_4$ with the β -aminophosphine or by treating $\text{Ph}_2\text{PCH}_2\text{CMe}_2\text{N}(\text{Li})\text{H}$ and 2-(Ph_2P) $\text{C}_6\text{H}_4\text{N}(\text{Li})\text{Me}$ with $[\text{Ir}(\text{cod})(\mu\text{-Cl})_2]$ to give the neutral alkyl and aryl amido compounds. When

combined with an alkali or amine base in methanol, all of the iridium complexes acted as catalysts for the direct hydrogenation of alkyl aryl ketones to the corresponding 1-phenylalkanols. The reactions, carried out at 25–50 °C and 10–50 bar of hydrogen, occurred with modest-to-good enantioselectivities (20–75% ee).

4. Allylic Alkylation


The palladium-catalyzed allylic alkylation has emerged as a powerful carbon–carbon-bond-forming reaction, and is now widely used in organic synthesis. The reaction is believed to proceed by nucleophilic addition to either C-1 or C-3 of a coordinated η^3 -allyl ligand (Scheme 14).^{2,29} The asymmetric version of this reaction has become quite popular, and aminophosphine ligands may provide a distinct advantage over symmetrical analogues as alkylation tends to occur at the position that is trans to the more strongly π -acidic PR_2 group.^{2,29} The enantioselective C–C-bond-formation step occurs via the major diastereomer of the equilibrating diastereomeric π -allyl intermediates.

Achiwa and co-workers have reported the synthesis of the chiral amidine ligand VALAP (**44**) from L-valine by condensation of aminophosphine **2a** with $\text{Me}_2\text{NCH}(\text{OMe})_2$ (eq 7).^{2a} VALAP has been utilized in the Pd-catalyzed asymmetric allylation of 1,3-diphenyl-2-propenyl acetate and pivalate (eq 8) with dimethyl malonate in the presence of BSA and LiOAc, affording excellent yields and up to 95% ee's. Loadings of $[\text{Pd}(\eta^3\text{-C}_3\text{H}_5)\text{Cl}]_2$ as low as 0.01 mol % still allowed for reasonable reaction times.

Morimoto modified the VALAP ligand (and the *tert*-butyl leucine analogue) via reaction of **44** with pyrrolidine and piperidine, or reaction of **2a** with para-substituted aromatic aldehydes (Scheme 15).^{2d} An examination of the effect of ligands **45c–h** on the allylic alkylation reaction showed a clear electronic effect wherein electron-donating substituents in the para position resulted in higher yields and ee's (eq 9).^{2d,30} This effect is most dramatic when comparing R = CF_3 (entry 3) and CH_3 (entry 4) which have a similar steric profile, yet the presence of the CH_3 group resulted in a marked improvement in both yield and ee. With the strongly electron-donating substituent, NMe_2 , both the catalytic activity and enantioselectivity are higher still than those obtained with the less electron-donating substituents. Indeed, use of this substituent allowed the $[\text{Pd}(\eta^3\text{-C}_3\text{H}_5)\text{Cl}]_2$ loading to be reduced to 0.005 mol % while still retaining excellent reactivity and leading to only a slight decrease in selectivity.


Saitoh et al. have also investigated the allylation reaction with silyl acetals and ketals **46a–d** and found that $[\text{Pd}(\eta^3\text{-C}_3\text{H}_5)\text{Cl}]_2$ –VALAP and related systems exhibit low-to-moderate activities with moderate-to-high enantioselectivities: $\leq 93\%$ ee using ligand **44** with acetal **46d** (eq 10).^{2d} When the analogous reaction with $\text{RR}'\text{C}=\text{C}(\text{OMe})(\text{OM})$ (M = Li, NR_4) as the nucleophile was examined, a low enantioselective induction was observed.

Yudin's group has employed iminophosphine ligands of type **47** (eq 11) in the palladium-catalyzed allylation.³¹ The $[\text{Pd}(\eta^3\text{-C}_3\text{H}_5)\text{Cl}]_2$ -**47** catalyzed allylation of 1,3-diphenyl-2-propenyl acetate in the presence of BSA and diethyl malonate was explored in order to determine the efficiency of the new chiral ligands for asymmetric induction. In the presence of aminocyclohexylphosphines, the precursors to ligands **47**, the reaction resulted in low yields and low enantioselectivities. When the catalytic reaction was carried out in the presence of the iminophosphine ligands, more favorable results were obtained. The yield and asymmetric induction for **47a** and **47b** were similar ($\leq 89\%$ yield and 87% ee), indicating that the ortho-methoxy fragment had little effect, whereas an electron-withdrawing




Ref 24a

eq 5




Ref 25,26

Figure 10. Me-AnilaPhos and IndolPhos.




Ref 2/a

Scheme 13. Hydrogenation of Olefins with Rh-PEAPhos.




Ref 27b

eq 6




Ref 2a,29

Scheme 14. Palladium-Catalyzed Allylic Alkylation.




Ref 2a

eq 7




Ref 2a

eq 8



Ref 2d

Scheme 15. Modified VALAP and Leucine Ligands.



Ref 2d

eq 9

substituent in this position, as in **47d**, had a deleterious effect on the reaction (60% yield and 51% ee). The bulky anthryl group, **47e**, greatly enhanced the reaction rate (complete in ca. 5 min), but resulted in a large decrease in enantioselectivity (21% ee).

Zheng and co-workers have recently reported on ferrocenylaminophosphine ligands that are capable of producing high yields and excellent asymmetric induction in the catalytic alkylation of 1,3-diphenyl-2-propenyl pivalate with dimethyl malonate.³² The ligands are prepared in one step from aminophosphines or phosphinoacetates and chloropyrimidines, chlorotriazines, or aminopyridines (**Scheme 16**). Their report indicates that increasing the number of nitrogen atoms in the ligand dramatically increases both the catalytic activity and enantioselectivity. For example, substituting the NMe₂ group with MeN(2-Py) results in an increase in enantioselectivity from 48% to 81% ee. When the pyrimidine-substituted ligand **51b** is used, an ee of 93% is obtained. The triazine-substituted ligand **50b** results in an enantioselectivity of 98% ee.

Gong, Mi, and co-workers have disclosed a series of amino-phosphinite ligands, **53–54** (**Figure 11**), that give good-to-excellent asymmetric induction in the Pd-catalyzed allylation of 1,3-diphenyl-2-propenyl acetate with dimethyl acetate.³³ The chiral ligands were prepared in one step by the reaction of aminoethanols with chlorodiphenylphosphine. These workers found that ligands with an NHR fragment (**54a–c**) gave higher ee's than ligands with an NMeR group (**53a–d**). The authors indicated that the N–H group was essential for optimal catalyst activity and selectivity, and proposed that the selectivity was a result of substrate interaction with the NH group.

5. Hydroformylation

While the hydroformylation of olefins employing rhodium catalysts represents an area of significant interest,³⁴ few recent reports have focused on the use of aminophosphine ligands bearing an NH group. Despite the relative scarcity of information, much is understood about the role and efficacy of such ligands in this process.

In a report by Andrieu and co-workers, diastereomeric trifunctional diaminophosphine ligands were derived from bidentate aminophosphine ligands by nucleophilic addition of a phosphinoalkyl carbanion (generated by lithiation) onto an imine (**Scheme 17**).³⁵ Both the bifunctional precursors and the derived trifunctional ligands were tested in the hydroformylation of styrene. Andrieu's group found that, while there was no impact on the isomer ratio, a substantial increase in activity was observed as a result of variation in the ligand set. An approximate threefold increase in activity using **55** or **56** relative to its precursor ligand suggested a dependence on the proximity and/or basicity of the dangling amine functionality. In subsequent studies,³⁶ it was determined that under catalytically relevant conditions, the aminophosphine ligand binds in a monodentate fashion through phosphorus while the amine functionality remains uncoordinated. The role of Bronsted base was proposed for the uncoordinated amine, which could assist in either the heterolytic splitting of dihydrogen or in the reductive elimination of HCl. Either scenario leads to an ammonium functionality in the dangling ligand. Based on a series of experiments designed to elucidate the mechanism, the authors proposed that a key step in rhodium-catalyzed hydroformylations employing aminophosphine ligands involves Rh–acyl racemization. This occurs via interaction of the acyl intermediate with the ammonium functionality of the dangling ligand (eq 12).

6. Conjugate Additions

6.1. Asymmetric Michael Addition to Enones

Ruthenium complexes of aminophosphines catalyze the asymmetric Michael addition reaction.^{15b} A range of such complexes (**Figure 12**) containing borohydride ligands were employed in the addition of dimethylmalonate to 2-cyclohexenone and, in a tandem process, were subsequently used in the asymmetric hydrogenation of the Michael adduct to the alcohol (**Scheme 18**).

The results of the Michael addition reaction clearly showed that catalyst activity and enantioselectivity were sensitive to solvent and ligand structure, respectively. While the activity of all of the catalysts employed in the Michael addition was insensitive to ligand structure, their sensitivity to solvent was pronounced. Enantioselectivity also displayed a strong dependence on solvent as a clear preference for aprotic solvents emerged (strongly coordinating acetonitrile also displayed deleterious effects on enantioselectivity). A pronounced favorable effect of ligand rigidity on enantioselectivity was observed, with the more rigid BINAP-supported catalysts (**61** and **62**) affording the highest enantiomeric excess ($\leq 97\%$). Furthermore, while (*R*)-BINAP (**61**) provided the *R* product, (*S*)-BINAP (**62**) gave the *S* isomer and the highest ee. In the subsequent hydrogenation, excellent diastereoselectivity was observed as a 30:1 *trans*:*cis* ratio was achieved for the alcohol product.

Zhang and co-workers reported on the use of larger-bite-angle aminophosphines in the copper-catalyzed addition of diethylzinc to enones.³⁷ The performance of the chiral binaphthyl ligands **63** and **64** (**Figure 13**) was evaluated in the conjugate addition of diethylzinc to 2-cyclohexenone and several acyclic enones, including chalcone and substituted chalcones, as well as an entirely aliphatic acyclic enone. In the case of 2-cyclohexenone, Zhang's group found that nonpolar solvents favored higher conversions and enantioselectivities over coordinating solvents. Mixtures of solvents such as toluene–dichloroethane were also effective. Removal of dissociated CH_3CN from the copper precursor $[\text{Cu}(\text{CH}_3\text{CN})_4]\text{BF}_4$ was important in realizing higher conversions and ee's and improved enantioselectivity was also gained from decreasing the temperature. $[\text{Cu}(\text{OTf})_2 \cdot \text{C}_6\text{H}_6]$ was the preferred copper precursor, as it allowed for room-temperature reactions albeit with diminished selectivity. Enantioselectivity was dependent on the ligand:metal ratio and was highest with a ligand:metal ratio of 5:1, but only slightly better than with a ratio of 2.5:1. The methyl-substituted ligand, **64**, gave modest improvements in ee over the unsubstituted analogue **63**. Ligand **64** provided much higher ee's in the alkylation of acyclic enones. The mixed solvent system toluene–dichloroethane was optimal with respect to yield and enantioselection, likely owing to the improved solubility of the substrates. This system proved to be competent in the asymmetric addition of diethylzinc to the entirely aliphatic acyclic enone as well.

A pair of amidoarylferrocenyldiphenylphosphines have also found application in the copper-catalyzed asymmetric addition of diethylzinc to enones (**Scheme 19**).³⁸ Johannsen and co-workers reported that the alkylated product from the addition of diethylzinc to *trans*-chalcone was obtained in reasonable yields and modest enantioselectivities. A strong dependence on solvent was observed as the highest yield (95%) and ee (58%) were realized in toluene, whereas halogenated solvents resulted in a dramatic reduction in both yield and ee. The better performance by ligand (*S*)-**66** in this alkylation prompted the authors to investigate the asymmetric addition of diethylzinc to the more challenging substrate diethyl ethylenemalonate. In this case, a conversion of >98% was obtained with moderate enantioselectivity (55%). No dependence on the ligand–metal ratio was observed in either of the enone addition reactions.

Adding to the diversity of scaffolds of aminophosphine ligands for conjugate additions, a series of carbohydrate-based aminophosphines were tested by Diéguez and his team in the copper-catalyzed addition of diethylzinc to 2-cyclohexenone.³⁹ The furanoside-supported aminophosphines (**Figure 14**) showed good activity, with phosphoramidite **70** being best in this regard (TOF >1200). In these systems, dichloromethane was the preferred solvent giving the highest conversions and enantiomeric excesses ($\leq 63\%$). The optimal temperature was 0 °C, as either increasing or decreasing the temperature resulted in diminished selectivity. Replacing the *tert*-butyl group at the para position of the biphenyl moiety with a methoxy group resulted in a decreased enantioselectivity. The aminomethyl substituent that gave both the greatest enantioselectivity and TOF was the phenylaminomethyl group. The sense of enantioselectivity was also influenced by the aminomethyl substituent. The more sterically demanding *tert*-butylamino group of ligand **67** gave preferentially the *R* product, while the less demanding isopropylamino and phenylamino substituents of **68** and **69**, respectively, provided the *S* isomer. Ligand **72**, having the opposite configuration at C-3 of the furanoside ring to that of ligand **69**, showed similar activity to **69**, however the enantioselectivity was dramatically reduced (only 8% ee). As mentioned above, phosphoramidite **70** gave the highest reaction rate, but the corresponding enantioselectivity was lower than that of **68** and **69**.

6.2. Asymmetric Addition of Organolithiums to Aldehydes

The asymmetric addition of organolithium reagents to aldehydes is a recent entry into the repertoire of transformations in which aminophosphine ligands play an important role.⁷ A series of aminophosphine ligands have been employed in the addition of *n*-butyllithium to benzaldehyde (**eq 13**).⁴⁰ A comparison of the aminophosphines with the corresponding ether and thioether ligands showed that the aminophosphines gave comparable high yields of the alcohol in consistently (and sometimes substantially) higher enantiomeric excess ($\leq 98\%$ ee).

7. Cycloaddition Reactions

While a broad range of both metals and ligand scaffolds have been employed in selective cycloaddition reactions,⁴¹ only recently has the potential utility of aminophosphines bearing NH groups come to light. The enantioselective addition of dimethyl maleate to iminoester **73a** is efficiently catalyzed by silver acetate in the presence of ferrocenyl-based aminophosphines (**eq 14**).⁴² The significance of this account lies in the fact that incorporation of H, rather than alkyl or aryl substituents, on nitrogen leads to the opposite absolute configuration of the product pyrrolidine **74a** (compare **75a** with **75b**). The ability of the H substituent to participate in substrate–ligand hydrogen bonding is implicated in the observed results. Increasing the steric bulk at phosphorus leads to improved enantioselectivity and the same reversal of configuration (**75c** vs **75d**). The mixed, NHMe-containing ligand, **75e**, gives dramatically reduced enantioselectivity (19%). Lowering the temperature to –25 °C results in greater selectivity than when the temperature is equal to 0 °C. A broad range of iminoesters and dipolarophiles were tested, and successful reversal of absolute configuration was maintained (**eq 15**).⁴²

8. Conclusions

Aminophosphines are a highly versatile class of ligands for asymmetric synthesis. While their applications in metal-catalyzed hydrogenations predominate, involvement of these ligands in other


asymmetric processes continues to gather interest. Simplicity and diversity in structure and preparation, coupled with a unique ability to become an integral component in chiral syntheses, guarantee continued development. With the current level of understanding of the role that these ligands assume in catalysis, researchers can apply them in catalytic transformations based on the nature of the substrate of interest. That is, in catalytic transformations where hydrogen-bonding interactions may play a role, aminophosphine ligands should be leading candidates in the ligand selection process.

9. Acknowledgment

We would like to acknowledge all of our colleagues at Kanata Chemical Technologies, Inc., for their efforts and support.


10. References

- (1) (a) Noyori, R.; Ohkuma, T. *Angew. Chem., Int. Ed.* **2001**, *40*, 40. (b) Noyori, R.; Yamakawa, M.; Hashiguchi, S. *J. Org. Chem.* **2001**, *66*, 7931. (c) Sandoval, C. A.; Ohkuma, T.; Muñiz, K.; Noyori, R. *J. Am. Chem. Soc.* **2003**, *125*, 13490. (d) Abdur-Rashid, K.; Faatz, M.; Lough, A. J.; Morris, R. H. *J. Am. Chem. Soc.* **2001**, *123*, 7473. (e) Abdur-Rashid, K.; Clapham, S. E.; Hadzovic, A.; Harvey, J. N.; Lough, A. J.; Morris, R. H. *J. Am. Chem. Soc.* **2002**, *124*, 15104. (f) Hartmann, R.; Chen, P. *Angew. Chem., Int. Ed.* **2001**, *40*, 3581.
- (2) (a) Saitoh, A.; Morimoto, T.; Achiwa, K. *Tetrahedron: Asymmetry* **1997**, *8*, 3567. (b) Saitoh, A.; Misawa, M.; Morimoto, T. *Tetrahedron: Asymmetry* **1999**, *10*, 1025. (c) Saitoh, A.; Uda, T.; Morimoto, T. *Tetrahedron: Asymmetry* **1999**, *10*, 4501. (d) Saitoh, A.; Achiwa, K.; Tanaka, K.; Morimoto, T. *J. Org. Chem.* **2000**, *65*, 4227.
- (3) Abdur-Rashid, K.; Guo, R.; Lough, A. J.; Morris, R. H.; Song, D. *Adv. Synth. Catal.* **2005**, *347*, 571.
- (4) Dahlenburg, L.; Götz, R. *J. Organomet. Chem.* **2001**, *619*, 88.
- (5) Caiazzo, A.; Dalili, S.; Yudin, A. K. *Org. Lett.* **2002**, *4*, 2597.
- (6) Guo, R.; Chen, X.; Jia, W.; Abdur-Rashid, K. U.S. Patent Appl. 60 942,699, 2007.
- (7) Rönnholm, P.; Södergren, M.; Hilmersson, G. *Org. Lett.* **2007**, *9*, 3781.
- (8) (a) Boaz, N. W.; Debenham, S. D. U.S. Patent 6,590,115 B2, July 8, 2003. (b) Boaz, N. W. U.S. Patent Appl. 6,906,213 B1, July 14, 2005. (c) Boaz, N. W.; Debenham, S. D.; Mackenzie, E. B.; Large, S. E. *Org. Lett.* **2002**, *4*, 2421. (d) Boaz, N. W.; Mackenzie, E. B.; Debenham, S. D.; Large, S. E.; Ponasik, J. A., Jr. *J. Org. Chem.* **2005**, *70*, 1872. (e) Boaz, N. W.; Ponasik, J. A., Jr.; Large, S. E. *Tetrahedron: Asymmetry* **2005**, *16*, 2063.
- (9) (a) Chen, W.; Mbafor, W.; Roberts, S. M.; Whittall, J. *J. Am. Chem. Soc.* **2006**, *128*, 3922. (b) Chen, W.; Mbafor, W.; Roberts, S. M.; Whittall, J. *Tetrahedron: Asymmetry* **2006**, *17*, 1161.
- (10) (a) Léautey, M.; Jubault, P.; Pannecoucke, X.; Quirion, J.-C. *Eur. J. Org. Chem.* **2003**, 3761. (b) Léautey, M.; Castellet-Deliencourt, G.; Jubault, P.; Pannecoucke, X.; Quirion, J.-C. *Tetrahedron Lett.* **2002**, *43*, 9237.
- (11) (a) Maj, A. M.; Pietrusiewicz, K. M.; Suisse, I.; Agbossou, F.; Mortreux, A. *Tetrahedron: Asymmetry* **1999**, *10*, 831. (b) Maj, A. M.; Pietrusiewicz, K. M.; Suisse, I.; Agbossou, F.; Mortreux, A. *J. Organomet. Chem.* **2001**, *626*, 157.
- (12) (a) Rahman, M. S.; Steed, J. W.; Hii, K. K. *Synthesis* **2000**, 1320. (b) Rahman, M. S.; Prince, P. D.; Steed, J. W.; Hii, K. K. *Organometallics* **2002**, *21*, 4927. (c) Rahman, M. S.; Oliana, M.; Hii, K. K. *Tetrahedron: Asymmetry* **2004**, *15*, 1835. (d) Oliana, M.; King, F.; Horton, P. N.; Hursthouse, M. B.; Hii, K. K. *J. Org. Chem.* **2006**, *71*, 2472.
- (13) Fu, X.; Jiang, Z.; Tan, C.-H. *Chem. Commun.* **2007**, 5058.
- (14) (a) Saudan, L. A. *Acc. Chem. Res.* **2007**, *40*, 1309. (b) Guiry, P. J.; Saunders, C. P. *Adv. Synth. Catal.* **2004**, *346*, 497. (c) Ohkuma, T.




Ref. 39

Figure 14. Furanose-Supported Aminophosphines for Conjugate Additions.




Ref. 40

eq 13



Ref. 42

eq 14



Ref. 42

eq 15

- Ooka, H.; Hashiguchi, S.; Ikariya, T.; Noyori, R. *J. Am. Chem. Soc.* **1995**, *117*, 2675. (d) Doucet, H.; Ohkuma, T.; Murata, K.; Yokozawa, T.; Kozawa, M.; Katayama, E.; England, A. F.; Ikariya, T.; Noyori, R. *Angew. Chem., Int. Ed.* **1998**, *37*, 1703.
- (15) (a) Guo, R.; Lough, A. J.; Morris, R. H.; Song, D. *Organometallics* **2004**, *23*, 5524. (b) Guo, R.; Morris, R. H.; Song, D. *J. Am. Chem. Soc.* **2005**, *127*, 516.
- (16) Chen, X.; Jia, W.; Guo, R.; Abdur-Rashid, K. U.S. Patent Appl. 60 948,231, 2007.
- (17) (a) Rautenstrauch, V.; Challand, R.; Churlaud, R.; Morris, R. H.; Abdur-Rashid, K.; Brazi, E.; Mimoun, H. Int. Patent Appl. WO 02/22526 A2, March 21, 2002. (b) Chen, X.; Guo, R.; Abdur-Rashid, K. U.S. Patent Appl. 60/948,238, 2007.
- (18) Dahlenburg, L.; Kühnlein, C. *J. Organomet. Chem.* **2005**, *690*, 1.
- (19) Boaz, N. W.; Large, S. E.; Ponasik, J. A., Jr.; Moore, M. K.; Barnette, T.; Nottingham, W. D. *Org. Process Rev. Dev.* **2005**, *9*, 472.
- (20) Li, X.; Jia, X.; Xu, L.; Kok, S. H. L.; Yip, C. W.; Chan, A. S. C. *Adv. Synth. Catal.* **2005**, *347*, 1904.
- (21) Deng, J.; Duan, Z.-C.; Huang, J.-D.; Hu, X.-P.; Wang, D.-Y.; Yu, S.-B.; Xu, X.-F.; Zheng, Z. *Org. Lett.* **2007**, *9*, 4825.
- (22) Hu, X.-P.; Zheng, Z. *Org. Lett.* **2004**, *6*, 3585.
- (23) Hu, X.-P.; Zheng, Z. *Org. Lett.* **2005**, *7*, 419.
- (24) (a) Franció, G.; Faraone, F.; Leitner, W. *Angew. Chem., Int. Ed.* **2000**, *39*, 1428. (b) Burk, S.; Franció, G.; Leitner, W. *Chem. Commun.* **2005**, 3460.
- (25) Vallianatou, K. A.; Kostas, I. D.; Holz, J.; Börner, A. *Tetrahedron Lett.* **2006**, *47*, 7947.
- (26) Wassenaar, J.; Reek, J. N. H. *Dalton Trans.* **2007**, 3750.
- (27) (a) Huang, J.-D.; Hu, X.-P.; Duan, Z.-C.; Zeng, Q.-H.; Yu, S.-B.; Deng, J.; Wang, D.-Y.; Zheng, Z. *Org. Lett.* **2006**, *8*, 4367. (b) Wang, D.-Y.; Hu, X.-P.; Huang, J.-D.; Deng, J.; Yu, S.-B.; Duan, Z.-C.; Xu, X.-F.; Zheng, Z. *Angew. Chem., Int. Ed.* **2007**, *46*, 7810.
- (28) (a) Dahlenburg, L.; Götz, R. *Eur. J. Inorg. Chem.* **2004**, 888. (b) Dahlenburg, L.; Götz, R. *Inorg. Chem. Commun.* **2003**, *6*, 443.
- (29) (a) Sprinz, J.; Kiefer, M.; Helmchen, G.; Reggelin, M.; Huttner, G.; Walter, O.; Zsolnai, L. *Tetrahedron Lett.* **1994**, *35*, 1523. (b) Dawson, G. J.; Williams, J. M. J.; Coote, S. J. *Tetrahedron: Asymmetry* **1995**, *6*, 2535. (c) Allen, J. V.; Coote, S. J.; Dawson, G. J.; Frost, C. G.; Martin, C. J.; Williams, J. M. J. *J. Chem. Soc., Perkin Trans. 1* **1994**, 2065. (d) Chelucci, G.; Cabras, M. A. *Tetrahedron: Asymmetry* **1996**, *7*, 965.
- (30) Saitoh, A.; Misawa, M.; Morimoto, T. *Synlett* **1999**, 483.
- (31) Dalili, S.; Caiazzo, A.; Yudin, A. K. *J. Organomet. Chem.* **2004**, *689*, 3604.
- (32) Hu, X.-P.; Chen, H.-L.; Zheng, Z. *Adv. Synth. Catal.* **2005**, *347*, 541.
- (33) Chen, G.; Li, X.; Zhang, H.; Gong, L.; Mi, A.; Cui, X.; Jiang, Y.; Choi, M. C. K.; Chan, A. S. C. *Tetrahedron: Asymmetry* **2002**, *13*, 809.
- (34) *Rhodium Catalyzed Hydroformylation*; Van Leeuwen, P. W. N. M., Claver, C., Eds.; Catalysis by Metal Complexes Series; Kluwer: Dordrecht, 2000; Vol. 22.
- (35) Andrieu, J.; Richard, P.; Camus, J.-M.; Poli, R. *Inorg. Chem.* **2002**, *41*, 3876.
- (36) (a) Andrieu, J.; Camus, J.-M.; Richard, P.; Poli, R.; Gonsalvi, L.; Vizza, F.; Peruzzini, M. *Eur. J. Inorg. Chem.* **2006**, 51. (b) Andrieu, J.; Camus, J.-M.; Balan, C.; Poli, R. *Eur. J. Inorg. Chem.* **2006**, 62.
- (37) Hu, X.; Chen, H.; Zhang, X. *Angew. Chem., Int. Ed.* **1999**, *38*, 3518.
- (38) Jensen, J. F.; Sotefte, I.; Sørensen, H. O.; Johannsen, M. *J. Org. Chem.* **2003**, *68*, 1258.
- (39) Diéguez, M.; Ruiz, A.; Claver, C. *Tetrahedron: Asymmetry* **2001**, *12*, 2861.
- (40) Granander, J.; Eriksson, J.; Hilmersson, G. *Tetrahedron: Asymmetry* **2006**, *17*, 2021.
- (41) Carmona, D.; Lamata, M. P.; Oro, L. A. *Coord. Chem. Rev.* **2000**, *200–202*, 717.
- (42) Zeng, W.; Chen, G.-Y.; Zhou, Y.-G.; Li, Y.-X. *J. Am. Chem. Soc.* **2007**, *129*, 750.

DABCO is a registered trademark of Air Products and Chemicals, Inc.


About the Authors

Dino Amoroso received his B.Sc. degree in chemistry in 1997 from McMaster University. In 2002, he received his Ph.D. degree from the University of Ottawa under the supervision of Professor Deryn Fogg. His graduate studies focused on diversifying ligand scaffolds employed in ruthenium-catalyzed olefin metathesis reactions to affect stereochemical control. After graduation, he moved to industry where he has developed transition-metal catalysts for a range of transformations including olefin polymerization, C–X bond formation, and hydrogenation. He is currently a Senior Research Scientist with Kanata Chemical Technologies (KCT) in Toronto.

Todd W. Graham received his Ph.D. degree in 1999 from the University of Alberta under the supervision of Professor Martin Cowie, studying the synthesis and reactivity of early-late heterobimetallic transition-metal complexes incorporating bifunctional cyclopentadienylalkylidene phosphine ligands. He then joined the group headed by Peter Maitlis and Michael Turner at the University of Sheffield, studying C–C bond forming reactions related to Fischer–Tropsch chemistry. Next, he joined Professor Douglas Stephan's group at the University of Windsor, where he prepared and studied the reactivity of low-valent titanium phosphinimide ($R_3P=N$) complexes. He then moved to Professor Arthur Carty's group at the National Research Council of Canada in Ottawa, where he examined the chemistry of electrophilic aminophosphinidene complexes ($L_nM=PNR_2$), which are phosphorus analogues of Fischer carbenes. He is currently a Research Scientist at KCT in Toronto, Canada.

Rongwei Guo received his Ph.D. degree in 2002 from Hong Kong's Polytechnic University under the supervision of Professor Albert S. C. Chan. His thesis research focused on the synthesis of novel chiral ligands and their applications in asymmetric catalysis. In 2003, he joined Professor Morris's group at the University of Toronto, where he worked on the enantioselective hydrogenation of C=O and C=N double bonds and the formation of C–C bonds. Since 2005, he has been employed by KCT in Toronto, Canada, where he is currently a Senior Research Scientist.

Chi-Wing Tsang received his Ph.D. degree in 2000 in the field of inorganic clusters from the Chinese University of Hong Kong under the direction of Professor Zuowei Xie. He then joined the research group of Professor Derek Gates at the University of British Columbia as a postdoctoral fellow studying inorganic polymers. He later moved to Ottawa to take up the position of Visiting Fellow at the National Research Council of Canada in the field of metal-containing biodegradable polymers. He is currently a Research Scientist at KCT in Toronto.

Kamaluddin Abdur-Rashid received his Ph.D. degree in 1994 at the University of the West Indies, Mona Campus, Jamaica, under the supervision of Professor Tara Dasgupta. He was a research associate from 1998 to 2002 in Professor Bob Morris's group at the University of Toronto, where he spearheaded the group's quest into pure and applied catalysis research. His discoveries led to the development of new classes of organometallic catalysts and their applications in organic synthesis, including industrial use. In 2004, he founded Kanata Chemical Technologies, Inc., an R&D company that is dedicated to the development and application of innovative catalyst technologies and processes. 

Accelerate Chiral Separation with ChiroSolv[®] Kits



The ChiroSolv Kit Advantage

Sigma-Aldrich[®] is pleased to partner with ChiroSolve, Inc., to offer a series of ready-to-use, disposable kits for chiral resolution of both **solid and liquid** racemates. These kits allow scientists to quickly screen calibrated quantities of resolving agents and solvents against a racemate to find the optimum combination, as well as optimize reaction conditions in order to separate a racemic mixture into its constituent enantiomers. The kits' high throughput format allows scientists to identify **within 24 hours** the optimum resolution process that might otherwise take over 2 months.

ChiroSolv Resolving Kits for Liquid Racemates

681431	Acid Series 1
681423	Acid Series 2
681415	Acid Series 3
699217	Acid Series 4
681407	Base Series 1
681393	Base Series 2
681377	Base Series 3
699241	Base Series 4

ChiroSolv Resolving Kits for Solid Racemates

698881	Acid Series 1
699527	Acid Series 2
698873	Acid Series 3
699225	Acid Series 4
699233	Base Series 1
698938	Base Series 2
698946	Base Series 3
698954	Base Series 4

- Ready to use kits — only your racemate is required
- Accurate and consistent results
- Tremendous time and money savings

For more information visit
sigma-aldrich.com/chirosolv

ChiroSolv is a registered trademark of ChiroSolve, Inc.
Sigma-Aldrich is a registered trademark of Sigma-Aldrich Biotechnology, L.P., and Sigma-Aldrich Co

New Labware Products for Chemistry



Pigment Free, Filler Free Sleeve Stopper Septa

These sleeve stopper septa are manufactured from pure natural rubber without fillers or pigments for super tactility and resealability. Less additives in the rubber lowers the potential for contamination during use.

Joint Size	Cat. No.
10/30	Z566136
14/20	Z566144
15/25	Z566152
24/40	Z566160

Precision Seal® Rubber Septa Caps

The unique design provides a penetration point for cannulation and a dual seal inside the tube and on the outer sleeve.

Description	Red	White
For 5-mm NMR tubes and ampoules	Z554014	Z553891
For 7-mm o.d. tubes and ampoules	Z565784	Z565768
For 10-mm o.d. tubes and ampoules	Z565792	Z565776

To learn more about our new labware products or to order, please visit sigma-aldrich.com/labware

Precision Seal is a registered trademark of Sigma-Aldrich Biotechnology, L.P., and Sigma-Aldrich Co



Sigma-Aldrich Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA S.A.
Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

SIGMA-ALDRICH PTY LTD.
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: (+43) 1 605 81 90
Fax: (+43) 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA.
Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

China

SIGMA-ALDRICH (SHANGHAI)
TRADING CO. LTD.
Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

SIGMA-ALDRICH spol. s r. o.
Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: (+45) 43 56 59 10
Fax: (+45) 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND OY
Tel: (+358) 9 350 9250
Free Fax: 9 350 92555

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 00
Fax: (+33) 474 95 68 08

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Greece

SIGMA-ALDRICH (O.M.) LTD.
Tel: (+30) 210 994 8010
Fax: (+30) 210 994 3831

Hungary

SIGMA-ALDRICH Kft
Ingyenes telefon: 06 80 355 355
Ingyenes zöld fax: 06 80 344 344
Tel: (+36) 1 235 9055
Fax: (+36) 1 235 9050

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: (+91) 80 6621 9600
New Delhi: (+91) 11 4165 4255
Mumbai: (+91) 22 2570 2364
Hyderabad: (+91) 40 4015 5488
Fax
Bangalore: (+91) 80 6621 9650
New Delhi: (+91) 11 4165 4266
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 4015 5466

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 1 404 1900
Fax: (+353) 1 404 1910

Israel

SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

SIGMA-ALDRICH S.r.l.
Numero Verde: 800 827018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

SIGMA-ALDRICH JAPAN K.K.
Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

SIGMA-ALDRICH KOREA
Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

SIGMA-ALDRICH QUIMICA, S.A. de C.V.
Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: 52 722 276 1600
Fax: 52 722 276 1601

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

SIGMA-ALDRICH NEW ZEALAND LTD.
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

SIGMA-ALDRICH NORWAY AS
Tel: (+47) 23 17 60 60
Fax: (+47) 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

SIGMA-ALDRICH RUS, LLC
Tel: +7 (495) 621 6037
+7 (495) 621 5828
Fax: +7 (495) 621 5923

Singapore

SIGMA-ALDRICH PTE. LTD.
Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD
Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313
SAFC (UK) Free Tel: 01202 712305

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Call Collect: (+1) 314 771 5750
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Internet

sigma-aldrich.com

Sigma-Aldrich Career Opportunities

As a leading Life Science and High Technology company, we are always looking for talented individuals to join our team. At Sigma-Aldrich we value the contributions of our employees, and recognize the impact they have on our success. We strive to foster creativity and innovation, and encourage professional development.

Our biochemical and organic chemical products and kits are used in scientific and genomic research, biotechnology, pharmaceutical development, the diagnosis of disease, and as key components in pharmaceutical and other high technology manufacturing. We have customers in life science companies, university and government institutions, hospitals, and in industry.

UNLEASH YOUR TALENTS

Learn more about our career opportunities by visiting
our award-winning Web site at sigma-aldrich.com/careers.

Sigma-Aldrich Corporation is an equal opportunity employer.

Sigma-Aldrich®

New Lab Start-Up Program



- Are you starting a new lab?
- Are you moving to a new location?
- Have you received your first research grant?

If you have answered "yes" to any of the above questions, let us help you get your new lab set up in an easy, effective, and economical way with the **Sigma-Aldrich® New Lab Start-Up Program**.

Get started now at
[sigma-aldrich.com/
newlabacta](http://sigma-aldrich.com/newlabacta)

sigma-aldrich.com

SIGMA-ALDRICH®



P.O. Box 14508
St. Louis, MO 63178
USA

PRESORTED
STANDARD
U.S. POSTAGE
PAID
SIGMA-ALDRICH
CORPORATION

Change Service Requested

IMPROVED CATALYSTS AND LIGANDS FOR ASYMMETRIC SYNTHESIS

Aldrichimica ACTA

VOL. 41, NO. 1 • 2008

 **ALDRICH**
Chemistry



Practical Organocatalysis with (*S*)- and
(*R*)-5-Pyrrolidin-2-yl-1*H*-tetrazoles

Aminophosphine Catalysts in Modern Asymmetric Synthesis


SIGMA-ALDRICH

New Products from Aldrich R&D

Aldrich Is Pleased to Offer Cutting-Edge Tools for Organic Synthesis

Togni Reagent for Electrophilic Trifluoromethylation

The direct transfer of a trifluoromethyl group usually requires harsh conditions that are often incompatible with more sensitive functionalities in a molecule. The Togni Reagent is an electrophilic reagent based on hypervalent iodine and is captured by a range of nucleophilic substrates under mild conditions. This reagent nicely complements the nucleophilic Ruppert's Reagent.



Eisenberger, P. et al. *Chem.—Eur. J.* **2006**, *12*, 2579


3,3-Dimethyl-1-(trifluoromethyl)-1,2-benziodoxole, 97%

696641 250 mg \$71.50
 [887144-97-0] 1 g 199.00
 $C_{10}H_{10}F_3IO$
 FW: 330.09



Ligands for Aqueous Transfer Hydrogenation


When used in conjunction with $[(Cp^*IrCl_2)_2]$, the ligands *N*-tosylethylenediamine (Ts(en)) and *N*-(2-aminoethyl)-4-(trifluoromethyl)benzenesulfonamide (CF_3 -Ts(en)) enable the facile and selective transfer hydrogenation of aldehydes with TOFs as high as $1.3 \times 10^5 h^{-1}$. Furthermore, the reactions are carried out in aqueous media and exhibit very good chemoselectivity and functional-group tolerance. In cases where α,β -unsaturated aldehydes are employed, reduction occurs selectively on the formyl group. Aliphatic aldehydes are also readily converted when the substrate is added portionwise over the course of the reaction



Wu, X. et al. *Angew. Chem., Int. Ed.* **2006**, *45*, 6718


N-Tosylethylenediamine, 97%

693855 5 g \$57.50
 [14316-16-6]
 $C_9H_{14}N_2O_2S$
 FW: 214.28




N-(2-Aminoethyl)-4-(trifluoromethyl)benzenesulfonamide, 97%

694967 500 mg \$49.00
 $C_9H_{11}F_3N_2O_2S$
 FW: 268.26



Thiocarbonyl Transfer Agent

1-(Methyldithiocarbonyl)imidazole is a stable, non-hazardous reagent that can replace thiophosgene, isothiocyanates, chlorothioformates, and other high-hazard reagents as a thiocarbonyl-transfer agent for the synthesis of dithiocarbonates,¹ dithiocarbamates, symmetrical and unsymmetrical thioureas,² and 2-thiohydantoin.³



(1) Sun, W. Y. et al. *Synlett* **1997**, 1279. (2) Mohanta, P. K. et al. *Tetrahedron* **2000**, *56*, 629. (3) Sundaram, G. S. M. et al. *Synlett* **2007**, 251


1-(Methyldithiocarbonyl)imidazole, 97%

694029 5 g \$74.50
 [74734-11-5]
 $C_4H_6N_2S$
 FW: 158.24



Air-Stable, Nucleophilic Alkylphosphine

1,3,5-Triaza-7-phosphaadamantane (PTA) is a convenient, efficient, and air-stable nucleophilic trialkylphosphine organocatalyst for the Baylis-Hillman reaction. Both aromatic and aliphatic aldehydes react with activated alkenes in the presence of 15–20 mol % of PTA to afford the corresponding adducts in fair-to-excellent yields. Furthermore, PTA displays activity that is superior to that of the structurally similar hexamethylenetetramine.



He, Z. et al. *Adv. Synth. Catal.* **2006**, *348*, 413

1,3,5-Triaza-7-phosphaadamantane

695467 500 mg \$25.00
 [53597-69-6]
 $C_6H_{10}N_3P$
 FW: 157.15



Aldrichimica ACTA

VOL. 41, NO. 1 • 2008

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-231-8327
SAFC® 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 414-438-3850
24-Hour Emergency 414-438-3850
Web Site sigma-aldrich.com
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201, USA

Subscriptions

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 355
Milwaukee, WI 53201-9358
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the inside back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com/acta.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2008 Sigma-Aldrich Co.

"PLEASE BOTHER US."



Joe Porwoll, President
Aldrich Chemical Co., Inc.

Professor Matthew Clarke of the University of St. Andrews (U.K.) kindly suggested that we make 1,3,5,7-tetramethyl-6-phenyl-2,4,8-trioxa-6-phosphaadamantane. This phosphine is very stable to air and moisture and has the stereoelectronic properties of bulky phosphonites. When employed with rhodium complexes, this ligand shows high catalytic activity in the hydroformylation of various alkenes. High selectivities and conversions as high as 99% have been reported.

Clarke, M. L.; Roff, G. J. *Chem.—Eur. J.* **2006**, *12*, 79/8



695459 1,3,5,7-Tetramethyl-6-phenyl-2,4,8-trioxa-6-phosphaadamantane, 97% 500 mg \$52.00
2 g 170.00

Naturally, we made this useful reagent. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the inside back cover.

TABLE OF CONTENTS

- Practical Organocatalysis with (S)- and (R)-5-Pyrrolidin-2-yl-1H-tetrazoles** 3
Deborah A. Longbottom, Vilijus Franckevičius, Sirirat Kumarn, Alexander J. Oelke, Veit Waschowski, and Steven V. Ley, University of Cambridge*
- Aminophosphine Catalysts in Modern Asymmetric Synthesis** 15
Dino Amoroso, Todd W. Graham, Rongwei Guo, Chi-Wing Tsang, and Kamaluddin Abdur-Rashid,* Kanata Chemical Technologies, Inc*

ABOUT OUR COVER

Landscape with Tobias and the Angel, with a View of Antwerp in the Background (oil on copper, 20.5 × 26.0 cm) was painted possibly around 1665 by Gillis Neyts, an enigmatic Flemish painter and engraver. Neyts (1623–1687) was born in Ghent, and spent a good part of his life in the city of Antwerp. He specialized in small, imaginary landscape scenes, which sometimes incorporated historical material or views of Flemish towns. His style approaches that of Lucas van Uden (1595–1672; Antwerp), who may have been his teacher.




Photograph © Alfred Bader

This small painting, with its soft and delicate handling, which was typical for Neyts, shows on the left just below the horizon a part of the skyline of the city of Antwerp. The spectacular form of the arching tree in the center frames the figures of two travelers (with walking sticks) in the foreground on the right. One of them appears to waive at the viewer, while the other—dressed in red and white and with wings rising from his shoulders—is identified as the Archangel Raphael accompanying young Tobias on his journey.


Neyts has painted here a fantasy landscape in which he transposes the ancient story of Tobias and the angel onto a contemporary setting, the outskirts of the 17th-century city of Antwerp. It would appear that Neyts's purpose is to help the viewer of that period identify more closely with the story.

This painting is in the private collection of Isabel and Alfred Bader.


Explore Organocatalysis




P0380




684341




688274



677183




671428



677019

Organocatalysis provides convenient new methods to construct complex chiral compounds with operational simplicity and without the need for metals. Excellent selectivities are observed in various organocatalyzed asymmetric transformations.

Sigma-Aldrich is pleased to offer an expanding portfolio of innovative organocatalysts to accelerate your research success.



To see a comprehensive listing of organocatalysts, please visit sigma-aldrich.com/organocat

Sigma-Aldrich Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA S.A.
Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

SIGMA-ALDRICH PTY LTD.
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: (+43) 1 605 81 90
Fax: (+43) 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA.
Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

SIGMA-ALDRICH CANADA LTD
Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

China

SIGMA-ALDRICH (SHANGHAI)
TRADING CO. LTD.
Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

SIGMA-ALDRICH spol. s r. o.
Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: (+45) 43 56 59 10
Fax: (+45) 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND OY
Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 00
Fax: (+33) 474 95 68 08

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Greece

SIGMA-ALDRICH (O.M.) LTD.
Tel: (+30) 210 994 8010
Fax: (+30) 210 994 3831

Hungary

SIGMA-ALDRICH Kft
Ingyenes zold telefon: 06 80 355 355
Ingyenes zold fax: 06 80 344 344
Tel: (+36) 1 235 9055
Fax: (+36) 1 235 9050

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: (+91) 80 6621 9600
New Delhi: (+91) 11 4165 4255
Mumbai: (+91) 22 2570 2364
Hyderabad: (+91) 40 4015 5488
Fax
Bangalore: (+91) 80 6621 9650
New Delhi: (+91) 11 4165 4266
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 4015 5466

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 1 404 1900
Fax: (+353) 1 404 1910

Israel

SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

SIGMA-ALDRICH S.r.l.
Numero Verde: 800 827018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

SIGMA-ALDRICH JAPAN K.K.
Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

SIGMA-ALDRICH KOREA
Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

SIGMA-ALDRICH QUIMICA, S.A. de C.V.
Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: 52 722 276 1600
Fax: 52 722 276 1601

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

SIGMA-ALDRICH NEW ZEALAND LTD.
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

SIGMA-ALDRICH NORWAY AS
Tel: (+47) 23 17 60 60
Fax: (+47) 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

SIGMA-ALDRICH RUS, LLC
Tel: +7 (495) 621 6037
+7 (495) 621 5828
Fax: +7 (495) 621 5923

Singapore

SIGMA-ALDRICH PTE. LTD.
Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD.
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

SIGMA-ALDRICH QUIMICA, S.A.
Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313
SAFC (UK) Free Tel: 01202 712305

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Call Collect: (+1) 314 771 5750
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Internet

sigma-aldrich.com

Sigma-Aldrich Career Opportunities

As a leading Life Science and High Technology company, we are always looking for talented individuals to join our team. At Sigma-Aldrich we value the contributions of our employees, and recognize the impact they have on our success. We strive to foster creativity and innovation, and encourage professional development.

Our biochemical and organic chemical products and kits are used in scientific and genomic research, biotechnology, pharmaceutical development, the diagnosis of disease, and as key components in pharmaceutical and other high technology manufacturing. We have customers in life science companies, university and government institutions, hospitals, and in industry.

UNLEASH YOUR TALENTS

Learn more about our career opportunities by visiting
our award-winning Web site at sigma-aldrich.com/careers.

Sigma-Aldrich Corporation is an equal opportunity employer.

Sigma-Aldrich®

New Lab Start-Up Program



- Are you starting a new lab?
- Are you moving to a new location?
- Have you received your first research grant?

If you have answered "yes" to any of the above questions, let us help you get your new lab set up in an easy, effective, and economical way with the **Sigma-Aldrich® New Lab Start-Up Program**.

Get started now at
[sigma-aldrich.com/
newlabacta](http://sigma-aldrich.com/newlabacta)

sigma-aldrich.com

SIGMA-ALDRICH



P.O. Box 14508
St. Louis, MO 63178
USA

PRESORTED
STANDARD
U.S. POSTAGE
PAID
SIGMA-ALDRICH
CORPORATION

Change Service Requested

ASYMMETRIC CATALYSIS WITH TRIP AND NHC-Cu(I)

Aldrichimica ACTA

VOL. 41, NO. 2 • 2008

 **ALDRICH**
Chemistry



TRIP—A Powerful Brønsted Acid Catalyst for Asymmetric Synthesis

**N-Heterocyclic Carbene–Copper Complexes:
Synthesis and Applications in Catalysis**


SIGMA-ALDRICH

New Products from Aldrich R&D

Aldrich Is Pleased to Offer Cutting-Edge Tools for Organic Synthesis

Reagent for the Synthesis of Diazoacetates

Diazoacetates are extremely useful reagents for a variety of reactions, such as cyclopropanations and C-H insertions. *N,N'*-Bis-(*p*-toluenesulfonyl)hydrazine is a stable, crystalline reagent that offers a useful alternative to tosyl azide for the production of diazoacetates. A variety of diazoacetates can be synthesized from the corresponding bromoacetates by treatment with *N,N'*-bis(*p*-toluenesulfonyl)hydrazine in moderate-to-high yields in relatively short reaction times.



Toma, T. et al. *Org. Lett.* **2007**, 9, 3195


N,N'-Bis(*p*-toluenesulfonyl)hydrazine

700525 1 g \$39.50

[14062-05-6]


C₁₄H₁₆N₂O₄S₂

FW: 340.42



Zn(tmp)₂ for the Selective Functionalization of C-H bonds

The selective functionalization of C-H bonds by neutral organozincs has typically been limited to fairly acidic carbon acids due to the slow kinetics of most dialkyl- and diarylzinc reagents. Bis(2,2,6,6-tetramethylpiperidinyl)zinc, Zn(tmp)₂, is an effective and versatile base for the mild, selective deprotonation of a broad range of substrates. The resultant organozincs can be conveniently coupled with aryl bromides using typical Pd-catalyzed cross-coupling methods




Hlavinka, M. L.; Hagadorn, J. R. *Organometallics* **2007**, 26, 4105.

Bis(2,2,6,6-tetramethylpiperidinyl)zinc, 0.5 M in toluene

697486 50 mL \$165.00


C₁₈H₃₆N₂Zn

FW: 345.88



7-Azaindoles

The derivatives of 7-azaindole (1*H*-pyrrolo[2,3-*b*]pyridine) have emerged as an important and promising class of compounds in agrochemistry and medicinal chemistry. Recent years have seen an exponential increase in reports about pharmaceutically active 7-azaindole derivatives in the medicinal chemistry literature. One recent report details the practical synthesis of the key pharmaceutical intermediate 2-[(7-azaindol-4-yl)methylamino]-5-fluoronicotinic acid from the 7-azaindole *N*-oxide MCBA salt, **693588**.



Wang, X. et al. *J. Org. Chem.* **2006**, 71, 4021.


7-Azaindole *N*-oxide 3-chlorobenzoate, 95%

693588 1 g \$50.00

[611197-49-0]

C₁₄H₁₁ClN₂O₃

FW: 290.71




7-Azaindole *N*-oxide hemihydrate, 97%

696838 1 g \$62.50

[55052-24-9]

C₇H₆N₂O

FW: 134.14



4-Chloro-7-azaindole, 97%

696218 250 mg \$49.00

[55052-28-3]

C₇H₅ClN₂

FW: 152.58



3-Nitro-7-azaindole

699268 1 g \$49.50

[23709-47-9]

C₇H₅N₃O₂

FW: 163.13



1-Boc-7-azaindole-3-carboxaldehyde

696811 1 g \$76.50

[144657-66-9]

C₁₃H₁₄N₂O₃

FW: 246.26



3-(Dimethylaminomethyl)-7-azaindole, 97%

696846 1 g \$64.00

[5654-92-2]

C₁₀H₁₃N₃

FW: 175.23




1-Boc-3-[(dimethylamino)methyl]-7-azaindole, 97%

699209 1 g \$87.50

[144657-65-8]

C₁₅H₂₁N₃O₂

FW: 275.35



Aldrichimica ACTA

VOL. 41, NO. 2 • 2008

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-231-8327
SAFC® 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 414-438-3850
24-Hour Emergency 414-438-3850
Web Site sigma-aldrich.com
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201, USA

Subscriptions

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 355
Milwaukee, WI 53201-9358
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the inside back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com/acta.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2008 Sigma-Aldrich Co.

"PLEASE BOTHER US."




Joe Porwoll

Joe Porwoll, President
Aldrich Chemical Co., Inc.

Professor Antonio M. Echavarren from the Institute of Chemical Research of Catalonia (ICIQ) (Tarragona, Spain) kindly suggested that we make (acetonitrile)[(2-biphenyl)di-tert-butylphosphine]gold(I) hexafluoroantimonate. This crystalline gold complex is air-stable and can be handled under ordinary bench-top conditions. The cationic complex is capable of catalyzing various unique intramolecular cyclization reactions.¹⁻³

(1) Nieto-Oberhuber, C. et al. *J. Am. Chem. Soc.* **2008**, *130*, 269. (2) Jiménez-Núñez, E. et al. *Angew. Chem., Int. Ed.* **2006**, *45*, 5452. (3) Ferrer, C.; Echavarren, A. M. *Angew. Chem., Int. Ed.* **2006**, *45*, 1105



697575 (Acetonitrile)[(2-biphenyl)di-tert-butylphosphine]gold(I) hexafluoroantimonate 250 mg \$65.00
1 g 235.00

Naturally, we made this useful reagent. It was no bother at all, just a pleasure to be able to help

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the inside back cover.

TABLE OF CONTENTS

TRIP—A Powerful Brønsted Acid Catalyst for Asymmetric Synthesis 31
Gareth Adair, Santanu Mukherjee, and Benjamin List, Max-Planck-Institut für Kohlenforschung

N-Heterocyclic Carbene—Copper Complexes: Synthesis and Applications in Catalysis 43
Silvia Díez-González and Steven P. Nolan, Institute of Chemical Research of Catalonia (ICIQ)

ABOUT OUR COVER

Camille Pissarro (1830–1903), one of the creators of the impressionist style, painted our cover, *Orchard in Bloom, Louveciennes* (oil on canvas, 45.1 × 54.9 cm) in 1872. Early in his career, Pissarro designated himself a pupil of Corot and, in this painting, Pissarro's broad method of composing and choice of a tranquil rural setting inhabited by a few small peasant figures still recall the older artist




Photograph © Board of Trustees, National Gallery of Art, Washington, DC

Pissarro completed this work shortly after he had returned to his home in Louveciennes after fleeing France during the Franco-Prussian War and Paris Commune. (Born in the Virgin Islands, then a possession of Denmark, Pissarro was a Danish citizen.) During the war, his house had been used by Prussian troops, and many of the canvases he left there were destroyed. He must have viewed the freshly plowed earth, like the spring blossoms that bring life to the dormant landscape, as a signal of renewed hope for his adopted country and for his career


When the idea arose for a group exhibition of work by the artists who would come to be called impressionists, Pissarro was among the earliest and most enthusiastic supporters. Pissarro drafted the group's written statement of purpose and was the only artist to participate in all eight impressionist exhibitions. This painting was one of five he showed at the first exhibition in 1874

This painting is part of the Ailsa Mellon Bruce Collection at the National Gallery of Art, Washington, DC.


Chiral Phosphoric Acids




248932




674745




674605



689890



688320



695718

BINOL-derived chiral phosphoric acids are useful organocatalysts for a range of transformations: reductive aminations and Pictet–Spengler, aza-Diels–Alder, addition, and cascade reactions are all effected with excellent enantioselectivities.

Sigma-Aldrich is pleased to offer an expanding portfolio of innovative organocatalysts to accelerate your research success.



To view a comprehensive listing of organocatalysts, please visit sigma-aldrich.com/organocat

Sigma-Aldrich Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA S.A.
Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

SIGMA-ALDRICH PTY LTD.
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA.
Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

China

SIGMA-ALDRICH (SHANGHAI)
TRADING CO. LTD.
Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

SIGMA-ALDRICH spol. s r. o.
Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: (+45) 43 56 59 10
Fax: (+45) 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND OY
Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 00
Fax: (+33) 474 95 68 08

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Greece

SIGMA-ALDRICH (O.M.) LTD.
Tel: (+30) 210 994 8010
Fax: (+30) 210 994 3831

Hungary

SIGMA-ALDRICH Kft
Ingyenes zöld telefon: 06 80 355 355
Ingyenes zöld fax: 06 80 344 344
Tel: (+36) 1 235 9055
Fax: (+36) 1 235 9050

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: (+91) 80 6621 9600
New Delhi: (+91) 11 4165 4255
Mumbai: (+91) 22 2570 2364
Hyderabad: (+91) 40 4015 5488
Fax
Bangalore: (+91) 80 6621 9650
New Delhi: (+91) 11 4165 4266
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 4015 5466

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 1 404 1900
Fax: (+353) 1 404 1910

Israel

SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

SIGMA-ALDRICH S.r.l.
Numero Verde: 800 827018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

SIGMA-ALDRICH JAPAN K.K.
Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

SIGMA-ALDRICH KOREA
Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: 52 722 276 1600
Fax: 52 722 276 1601

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

SIGMA-ALDRICH NEW ZEALAND LTD.
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

SIGMA-ALDRICH NORWAY AS
Tel: (+47) 23 17 60 60
Fax: (+47) 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

SIGMA-ALDRICH RUS, LLC
Tel: +7 (495) 621 6037
+7 (495) 621 5828
Fax: +7 (495) 621 5923

Singapore

SIGMA-ALDRICH PTE. LTD.
Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD.
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD
Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313
SAFC (UK) Free Tel: 01202 712305

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Call Collect: (+1) 314 771 5750
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Internet

sigma-aldrich.com

Sigma-Aldrich Career Opportunities

As a leading Life Science and High Technology company, we are always looking for talented individuals to join our team. At Sigma-Aldrich we value the contributions of our employees, and recognize the impact they have on our success. We strive to foster creativity and innovation, and encourage professional development.

Our biochemical and organic chemical products and kits are used in scientific and genomic research, biotechnology, pharmaceutical development, the diagnosis of disease, and as key components in pharmaceutical and other high technology manufacturing. We have customers in life science companies, university and government institutions, hospitals, and in industry.

UNLEASH YOUR TALENTS

Learn more about our career opportunities by visiting
our award-winning Web site at sigma-aldrich.com/careers

Sigma-Aldrich Corporation is an equal opportunity employer.

Green Tidbit

Cyclopentyl Methyl Ether (CPME)

Alternative to Tetrahydrofuran, 2-Methyltetrahydrofuran, *tert*-Butyl Methyl Ether (MTBE), 1,4-Dioxane, and other Ether Solvents.

CPME provides a green alternative for those looking to improve their chemical process. CPME not only reduces energy waste, but also improves laboratory safety due to its unique composition which resists peroxide formation.

Cyclopentyl Methyl Ether, contains 50 ppm BHT	Cat. No.
Anhydrous, $\geq 99.9\%$	675970
ReagentPlus [®] , $\geq 99.90\%$	675989

Features and Benefits

- More stable than THF when it comes to forming peroxides
- High boiling point (106 °C)
- Narrow explosion range (1.1–9.9% by vol.)
- Stable under acidic and basic conditions
- Forms azeotropes with water
- Conventional drying is unnecessary for general organometallic reactions

To learn more about CPME's unique resistance to peroxide formations, visit us at sigma-aldrich.com/greensolvents

ReagentPlus is a registered trademark of Sigma-Aldrich Biotechnology, L.P., and Sigma-Aldrich Co

sigma-aldrich.com

SIGMA-ALDRICH



P.O. Box 14508
St. Louis, MO 63178
USA

PRESORTED
STANDARD
U.S. POSTAGE
PAID
SIGMA-ALDRICH
CORPORATION

Change Service Requested

GREEN, MILD, AND VERSATILE SYNTHETIC METHODS

Aldrichimica ACTA

VOL. 41, NO. 3 • 2008

 **ALDRICH**
Chemistry



Transition-Metal-Catalyzed Cross-Couplings
Going Green: in *Water* at Room Temperature

Preparation of Chiral Diamines by the
Diaza-Cope Rearrangement (DCR)


SIGMA-ALDRICH

New Products from Aldrich R&D

Aldrich Is Pleased to Offer Cutting-Edge Tools for Organic Synthesis

Oxidation Reagent

The selective oxidation of alcohols under conditions which are environmentally friendly and/or do not require toxic and hazardous chemicals is a well-precedented problem in organic synthesis. More recently, conditions have been developed which use a catalyst system comprised of DABCO[®]-CuCl in the presence of TEMPO and molecular oxygen as the oxidant. This mild catalytic oxidation generates water as the only byproduct




Mannam, S. et al. *Adv. Synth. Catal.* **2007**, 349, 2253

DABCO[®]-CuCl Complex

703141 1 g \$49.00
 $\text{C}_7\text{H}_{14}\text{N}_6\text{Cu}$
 FW: 211.17

Reagent for the Preparation of Organozincs

$\text{Zn}(\text{OMe})_2$ is an exceptional precursor for the formation of reactive and selective organozinc reagents under salt-free conditions. Historically, the difficulty in preparing organozinc reagents that has prevented their widespread use is due to their highly pyrophoric nature as well as the generation of byproducts. The use of additives, which remove these byproducts, presents other limitations. Therefore, a method to control the solubility of the byproducts so that they can be removed by simple filtration or centrifugation has been developed, and relies on the use of $\text{Zn}(\text{OMe})_2$ as a precursor to the organozinc reagent to be formed in solution. Following formation of the active organozinc compound, filtration or centrifugation is followed by the desired organic transformation. Examples of such a transformation include the catalytic enantioselective addition of organozincs to imines, conjugate addition, addition to aldehydes (example below), and addition to β -nitrostyrene.




Cote, A.; Charette, A. B. *J. Am. Chem. Soc.* **2008**, 130, 2771

Zinc Methoxide

702684 1 g \$75.00
 $\text{C}_2\text{H}_6\text{O}_2\text{Zn}$ 5 g 249.00
 FW: 127.46

Reagent for the PMB-protection of Alcohols

Protection with the *para*-methoxybenzyl protecting group (PMB) has been an ongoing challenge in organic chemistry due to the limitations with regard to the reagents that can be used. Common reagents typically require acidic or basic reaction media and present problems that relate to their long-term storage. The lepidine ether below (2-(4-methoxybenzyloxy)-4-methylquinoline) was developed to address some of these limitations and, in combination with MeOTf, affords an active reagent *in situ*. A plethora of alcohols are readily protected under neutral reaction conditions. Additionally, 2-(4-methoxybenzyloxy)-4-methylquinoline is stable, and byproducts generated by its use are easily removed through aqueous workup or chromatography.




Nwoye, E. O.; Dudley, G. B. *Chem. Commun.* **2007**, 1436

2-(4-Methoxybenzyloxy)-4-methylquinoline

701440 1 g \$39.00
 [937184-70-8] 5 g 130.00
 $\text{C}_{18}\text{H}_{17}\text{NO}_2$
 FW: 279.33

Reagent for Selective Methylhydrazine Addition

N-Boc-*N*-methylhydrazine is a useful reagent that can be employed when reaction at the less nucleophilic nitrogen of methylhydrazine is desired. Without protection of the methylated nitrogen with a Boc group, reaction would occur also at the methylated nitrogen. As demonstrated below, reaction of the chlorotriazine provides the corresponding, properly hydrazine-functionalized arene.



Kelly, T. R. et al. *J. Am. Chem. Soc.* **2006**, 128, 5646

1-Boc-1-methylhydrazine

699101 5 g \$56.50
 [21075-83-2] 25 g 187.50
 $\text{C}_7\text{H}_{12}\text{N}_2\text{O}_2$
 FW: 146.19

Aldrichimica ACTA

VOL. 41, NO. 3 • 2008

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries	800-325-3010
Technical Service	800-231-8327
SAFC®	800-244-1173
Custom Synthesis	800-244-1173
Flavors & Fragrances	800-227-4563
International	414-438-3850
24-Hour Emergency	414-438-3850
Web Site	sigma-aldrich.com
Email	aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201, USA

Subscriptions

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 355
Milwaukee, WI 53201-9358
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the inside back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com/acta.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.


Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2008 Sigma-Aldrich Co.

"PLEASE BOTHER US."



Joe Porwoll, President
Aldrich Chemical Co. Inc.

Professor Michael Krische from the University of Texas at Austin kindly suggested that we make Ir and Rh BARF (BARF = {3,5-(CF₃)₂C₆H₃}₂B⁻) salts. These compounds with loosely coordinating properties catalyze various transformations, including hydrogenation and reductive coupling, that otherwise do not proceed effectively.



(1) Ngai, M.-Y.; Barchuk, A.; Krische, M. J. *J. Am. Chem. Soc.* **2007**, *129*, 280. (2) Ngai, M.-Y.; Kong, J.-R.; Krische, M. J. *J. Org. Chem.* **2007**, *72*, 1063

693774 Bis(cyclooctadiene)iridium(I) tetrakis(3,5-bis(trifluoromethyl)phenyl)borate
500 mg \$116.00
2 g 435.00

692573 Bis(1,5-cyclooctadiene)rhodium(I)tetrakis[bis(3,5-trifluoromethyl)phenyl]borate
100 mg \$45.90
500 mg 183.00

Naturally, we made these useful reagents. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the inside back cover.

TABLE OF CONTENTS

Transition-Metal-Catalyzed Cross-Couplings Going Green: in Water at Room Temperature	59
<i>Bruce H. Lipshutz¹ and Subir Ghorai, University of California, Santa Barbara</i>	
Preparation of Chiral Diamines by the Diaza-Cope Rearrangement (DCR)	77
<i>Hyunwoo Kim, Soon Mog So, and Jik Chin,² University of Toronto; B. Moon Kim, Seoul National University</i>	

ABOUT OUR COVER

Looking at *The Bridge at Argenteuil* (oil on canvas, 60 × 79.7 cm) from a distance of ten feet or so, Claude Monet's brushstrokes blend to yield a convincing view of the river Seine and the pleasure boats that drew tourists to Argenteuil. Up close, however, each dab of paint is distinct, and the scene dissolves into a mosaic of paint—brilliant, unblended tones of blue, red, green, and yellow. In the water, quick, fluid skips of the brush mimic the lapping surface. In the trees, thicker paint is applied with denser, stubbier strokes. The figure in the sailboat is only a ghostly wash of dusty blue, and the women rowing nearby are indicated by mere shorthand.



Photograph © Board of Trustees, National Gallery of Art, Washington



In the early years of impressionism, Monet, Renoir, and others strove to capture the fleeting effects of light and atmosphere on the landscape and to transcribe directly and quickly their sensory experience of it. Monet advised his students, "When you go out to paint, try to forget what objects you have before you, a tree, a house, a field or whatever. Merely think here is a little square of blue, here an oblong of pink, here a streak of yellow, and paint it just as it looks to you, the exact color and shape, until it gives your own naive impression of the scene before you."

In this early work (1874), Monet (1840–1926) captures a warm, sunny, idyllic day—a motif he used often and for which he became famous. Today, Monet's characteristic style and distinctive brushstroke are still fresh, recognizable, and most popular.

This painting is part of the Collection of Mr. and Mrs. Paul Mellon at the National Gallery of Art, Washington, DC.

PTS—New Amphiphile for Metathesis and Cross-Coupling in Water

Recently introduced by Professor Bruce Lipshutz of UC, Santa Barbara, polyoxyethanyl α -tocopheryl sebacate (PTS) is a nonionic amphiphile that is proving to be a versatile "solubilizer" for organic molecules in water.¹ Lipophilic substrates and catalysts can efficiently enter micelles formed by PTS in water, leading to cross-coupling reactions at room temperature without the need for a co-solvent.²



Polyoxyethanyl α -tocopheryl sebacate, 15 wt. % in H₂O
698717

10 mL \$47.50

(1) Sold under license from Zymes, LLC. (2) (a) Lipshutz, B. H. et al. *Org. Lett.* **2008**, *10*, 1325. (b) Lipshutz, B. H. et al. *Adv. Synth. Catal.* **2008**, *350*, 953. (c) Lipshutz, B. H. et al. *Org. Lett.* **2008**, *10*, 1333. (d) Lipshutz, B. H.; Taft, B. R. *Org. Lett.* **2008**, *10*, 1329. (e) Lipshutz, B. H.; Ghorai, S. *Aldrichimica Acta* **2008**, *41*, in press. (f) Lipshutz, B. H. et al. *Org. Lett.* **2008**, *10*, ASAP.

Sigma-Aldrich® Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA S.A.
Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

SIGMA-ALDRICH PTY LTD.
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA.
Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

China

SIGMA-ALDRICH (SHANGHAI)
TRADING CO. LTD.
Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

SIGMA-ALDRICH spol. s r. o.
Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: (+45) 43 56 59 10
Fax: (+45) 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND OY
Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 00
Fax: (+33) 474 95 68 08

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Greece

SIGMA-ALDRICH (O.M.) LTD.
Tel: (+30) 210 994 8010
Fax: (+30) 210 994 3831

Hungary

SIGMA-ALDRICH Kft
Ingyenes zöld telefon: 06 80 355 355
Ingyenes zöld fax: 06 80 344 344
Tel: (+36) 1 235 9055
Fax: (+36) 1 235 9050

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: (+91) 80 6621 9600
New Delhi: (+91) 11 4165 4255
Mumbai: (+91) 22 2570 2364
Hyderabad: (+91) 40 4015 5488
Fax
Bangalore: (+91) 80 6621 9650
New Delhi: (+91) 11 4165 4266
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 4015 5466

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 1 404 1900
Fax: (+353) 1 404 1910

Israel

SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

SIGMA-ALDRICH S.r.l.
Numero Verde: 800 827018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

SIGMA-ALDRICH JAPAN K.K.
Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

SIGMA-ALDRICH KOREA
Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD.
Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: 52 722 276 1600
Fax: 52 722 276 1601

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

SIGMA-ALDRICH NEW ZEALAND LTD.
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

SIGMA-ALDRICH NORWAY AS
Tel: (+47) 23 17 60 60
Fax: (+47) 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

SIGMA-ALDRICH RUS, LLC
Tel: +7 (495) 621 6037
+7 (495) 621 5828
Fax: +7 (495) 621 5923

Singapore

SIGMA-ALDRICH PTE. LTD.
Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD.
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313
SAFC (UK) Free Tel: 01202 712305

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Call Collect: (+1) 314 771 5750
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Internet

sigma-aldrich.com



Mixed Sources

Product group from well-managed
forests, controlled sources and
recycled wood or fiber
www.fsc.org Cert no SGS-COC-003338
© 1996 Forest Stewardship Council



Sigma-Aldrich Career Opportunities

As a leading Life Science and High Technology company, we are always looking for talented individuals to join our team. At Sigma-Aldrich we value the contributions of our employees, and recognize the impact they have on our success. We strive to foster creativity and innovation, and encourage professional development.

Our biochemical and organic chemical products and kits are used in scientific and genomic research, biotechnology, pharmaceutical development, the diagnosis of disease, and as key components in pharmaceutical and other high technology manufacturing. We have customers in life science companies, university and government institutions, hospitals, and in industry.

UNLEASH YOUR TALENTS

Learn more about our career opportunities by visiting
our award-winning Website at sigma-aldrich.com/careers

Sigma-Aldrich Corporation is an equal opportunity employer.

Aldrichimica Acta Archival Collection on CD

40th Anniversary Edition (1968–2007)


Search, view, and print—at your convenience, at home or at work—review articles written by some of the world's leading chemists on some of the most compelling chemistry topics.



Request your complimentary copy of the CD at sigma-aldrich.com/acta40cd

To receive the *Aldrichimica Acta* on a regular basis, subscribe for FREE at sigma-aldrich.com/acta

Aldrichimica Acta
Growing Impact over the Past Six Years*



*Impact Factor data published by Thomson Scientific's *Journal Citation Reports*®

sigma-aldrich.com

SIGMA-ALDRICH



P.O. Box 14508
St. Louis, MO 63178
USA

PRESORTED
STANDARD
U.S. POSTAGE
PAID
SIGMA-ALDRICH
CORPORATION

Change Service Requested

CARBONYL COMPOUNDS: STILL CENTRAL TO ORGANIC SYNTHESIS

Aldrichimica ACTA

VOL. 41, NO. 4 • 2008

ALDRICH
Chemistry



Formation of C–C Bonds via Catalytic Hydrogenation and Transfer Hydrogenation: Vinylation, Allylation, and Enolate Addition

Amino Carbonyl Compounds in Organic Synthesis


SIGMA-ALDRICH

New Products from Aldrich R&D

Aldrich Is Pleased to Offer Cutting-Edge Tools for Organic Synthesis

Reagents for the Bromination of Alcohols

There are various methods for the conversion of alcohols to bromides; however, commonly employed methods either use or generate toxic HBr gas. The use of hexabromoacetone ($\text{Br}_3\text{CCOBr}_3$) and ethyl tribromoacetate ($\text{Br}_3\text{CCO}_2\text{Et}$) as less toxic, milder bromination reagents has recently been reported. Both reagents provide the desired alkyl bromide in excellent yield.



Tongkate, P. et al. *Tetrahedron Lett.* **2008**, 49, 1146

Ethyl tribromoacetate, 97%


704679	1 g	\$34.50
[599-99-5]	5 g	114.50
$\text{C}_4\text{H}_5\text{Br}_3\text{O}_2$		
FW: 324.79		

1,1,1,3,3,3-Hexabromoacetone, 97%

702404	5 g	\$61.50
[23162-64-3]	25 g	205.00
$\text{C}_3\text{Br}_6\text{O}$		
FW: 531.46		


New Boronic Acid Surrogates for Iterative Cross-Coupling

Professor Martin Burke and co-workers at the University of Illinois (Urbana-Champaign) have recently disclosed a technology employing boronic acid surrogates (termed "MIDA boronates") for use in iterative Suzuki cross-coupling reactions. The air-stable, chromatography-compatible, and easily deprotected boron building blocks permit difficult couplings through the attenuation of transmetalation via pyramidalization of the boron atom. The chemistry has been applied to the preparation of polyenyl MIDA boronates, for which the boronic acid counterpart is unstable. This subsequently led to the efficient synthesis of the left half of amphotericin B.




Lee, S. J. et al. *J. Am. Chem. Soc.* **2008**, 130, 466

trans-2-Bromovinylboronic acid MIDA ester


703478		500 mg	\$120.00
$\text{C}_7\text{H}_9\text{BBrNO}_4$		1 g	200.00
FW: 261.87			

New Aldehydes from Aldrich R&D


1-(2-Tetrahydropyryl)-1H-pyrazole-5-carboxaldehyde

699365		1 g	\$89.50
[957483-88-4]			
$\text{C}_8\text{H}_{12}\text{N}_2\text{O}_2$			
FW: 180.20			


3-Methylpyridine-2-carboxaldehyde, 97%

699071		1 g	\$82.00
[55589-47-4]			
$\text{C}_7\text{H}_9\text{NO}$			
FW: 121.14			


5-Hexylthiophene-2-carboxaldehyde, 97%

699187		1 g	\$89.50
[100943-46-2]			
$\text{C}_{11}\text{H}_{16}\text{OS}$			
FW: 196.31			

4-Oxazolecarboxaldehyde, 97%

697915		250 mg	\$84.40
[118994-84-6]		1 g	237.50
$\text{C}_4\text{H}_3\text{NO}_2$			
FW: 97.07			

4-Bromothiazole-2-carboxaldehyde, 96%

699284		1 g	\$47.50
[167366-05-4]			
$\text{C}_4\text{H}_2\text{BrNOS}$			
FW: 192.03			

Aldrichimica ACTA

VOL. 41, NO. 4 • 2008

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
 6000 N. Teutonia Ave.
 Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
 FAX 800-325-5052 (USA)
 or 414-438-2199
 Mail P.O. Box 2060
 Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
 Technical Service 800-231-8327
 SAFC® 800-244-1173
 Custom Synthesis 800-244-1173
 Flavors & Fragrances 800-227-4563
 International 414-438-3850
 24-Hour Emergency 414-438-3850
 Web Site sigma-aldrich.com
 Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
 P.O. Box 355, Milwaukee, WI 53201, USA

Subscriptions

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
 Mail: Attn: Mailroom
 Aldrich Chemical Co., Inc.
 Sigma-Aldrich Corporation
 P.O. Box 355
 Milwaukee, WI 53201-9358
 Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the inside back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com/acta.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2008 Sigma-Aldrich Co.

"PLEASE BOTHER US."




Joe Porwoll

Joe Porwoll, President
 Aldrich Chemical Co., Inc

Professor Hisashi Yamamoto of The University of Chicago kindly suggested that we make bis(hydroxamic acid) based ligands, which, in combination with VO(Oi-Pr)₃, generate highly active catalysts for the asymmetric epoxidation of allylic alcohols. Good-to-excellent yields and enantioselectivities of up to 97% ee have been reported.

Zhang, W. et al. *Angew. Chem., Int. Ed.* **2005**, *44*, 4389.



700592 (1R,2R)-N,N'-Dihydroxy-N,N'-bis(diphenylacetyl)-1,2-cyclohexanediamine, 97% (R)-CBHA-DPA	50 mg	\$ 115.00
700576 (1S,2S)-N,N'-Dihydroxy-N,N'-bis(diphenylacetyl)-1,2-cyclohexanediamine, 97% (S)-CBHA-DPA	50 mg	\$ 115.00

Naturally, we made these useful ligands. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the inside back cover.

TABLE OF CONTENTS

Formation of C-C Bonds via Catalytic Hydrogenation and Transfer Hydrogenation: Vinylation, Allylation, and Enolate Addition	95
<i>Ryan L. Patman, John F. Bower, In Su Kim, and Michael J. Krische,* University of Texas at Austin</i>	
Amino Carbonyl Compounds in Organic Synthesis	109
<i>Sivaraj Baktharaman, Ryan Hill, and Andrei K. Yudin,* University of Toronto</i>	

ABOUT OUR COVER

Pont Neuf, Paris (oil on canvas, 75.3 × 93.7 cm), was painted by the French Impressionist painter, Pierre Auguste Renoir (1841–1919), on a gloriously sunny and warm summer's day in 1872. While Renoir's paintings of women and children are better known, his landscapes resonate with a vigor and freshness new to the art scene at the time.

In this painting, Renoir uses a bright, light palette to emphasize an intense midday sun, but he deliberately suppresses incidental detail and clarity, displaying the basic tenets central to the development of Impressionism. Although Renoir presents us with an impressionistic view of the scene, the bridge and buildings in the background are accurate enough to be identifiable, then and now.

Renoir occupied an upper floor of a café on the left bank of the Seine to depict this famous view of the ninth bridge. Pont Neuf connects the Île de la Cité with the rest of Paris. Edmond Renoir, the artist's younger brother and novice journalist in 1872, later recounted in an interview that he helped his brother by periodically delaying a Parisian on the bridge long enough for the artist to record their appearance. Renoir captures Edmond, walking stick in hand, wearing a light-colored straw hat and slacks and a dark jacket in at least two locations. Can you find him?




Photograph © Board of Trustees, National Gallery of Art, Washington

This painting is part of the Ailsa Mellon Bruce Collection at the National Gallery of Art, Washington, DC.

Metal Complexes and Ligands for Enantioselective Reductive Coupling


Asymmetric hydrogenation is one of the most utilized reactions to induce chirality in a molecule. It is currently widely used in industry. Krische and co-workers have developed a new type of transformation based on the enantioselective reductive C–C bond formation mediated by hydrogen. Utilizing a rhodium-, iridium-, or ruthenium-based complex with a variety of ligands, Krische and co-workers demonstrated the potency of this reaction for the reductive coupling of conjugated enones, dienes, imines, enynes, and carbonyls. Aldrich is offering a series of complexes and ligands for enantioselective reductive coupling.

Enantioselective Imine Vinylation




Ngai, M.-Y. et al. *J. Am. Chem. Soc.* **2007**, 129, 12644. Skucas, E. et al. *J. Am. Chem. Soc.* **2007**, 129, 7242.

Enantioselective Reductive Coupling of Alkynes with Glyoxalates



Hong, Y.-T. et al. *Org. Lett.* **2007**, 9, 3745. Cho, C.-W.; Krische, M. *J. Org. Lett.* **2006**, 8, 3873




693774



692573




683159



275131




334995




(*R*)-(+)-MeO-BIPHEP **29510**

(*S*)-(-)-MeO-BIPHEP **29511**





(*R*)-(-)-Cl,MeO-BIPHEP **76854**

(*S*)-(+)-Cl,MeO-BIPHEP **96738**



(*R*)-(3,5-*t*-Bu-4-MeOPh)-MeO-BIPHEP **29512**

(*S*)-(3,5-*t*-Bu-4-MeOPh)-MeO-BIPHEP **29513**



(*R*)-Xylyl-WALPHOS **65683**

(*S*)-Xylyl-WALPHOS **65684**

New Aldrich® Handbook!

Over 6,000 Innovative New Products



Request your copy of the new 2009-2010
Aldrich Handbook set today.

2009-2010 Aldrich Handbook

- 10,000 chemical structures
- 8,500 updated literature citations
- Extensive chemical & physical data

Labware Catalog

- Innovative new products
- Improved product images
- New technical information index

To request your
complimentary Aldrich
Handbook set, visit
sigma-aldrich.com/aldrichcat



P.O. Box 14508
St. Louis, MO 63178
USA

PRESORTED
STANDARD
U.S. POSTAGE
PAID
SIGMA-ALDRICH
CORPORATION

Change Service Requested

Sigma-Aldrich® Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA S.A.
Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

SIGMA-ALDRICH PTY LTD.
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

SIGMA-ALDRICH NV/S.A.
Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

China

SIGMA-ALDRICH (SHANGHAI)
TRADING CO. LTD.
Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

SIGMA-ALDRICH spol. s r. o.
Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: (+45) 43 56 59 10
Fax: (+45) 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND OY
Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 00
Fax: (+33) 474 95 68 08

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Greece

SIGMA-ALDRICH (O.M.) LTD.
Tel: (+30) 210 994 8010
Fax: (+30) 210 994 3831

Hungary

SIGMA-ALDRICH Kft
Ingyenes telefonszám: 06 80 355 355
Ingyenes fax szám: 06 80 344 344
Tel: (+36) 1 235 9055
Fax: (+36) 1 235 9050

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: (+91) 80 6621 9600
New Delhi: (+91) 11 4358 8000
Mumbai: (+91) 22 2570 2364
Hyderabad: (+91) 40 4015 5488
Fax
Bangalore: (+91) 80 6621 9650
New Delhi: (+91) 11 4358 8001
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 4015 5466

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: +353 (0) 402 20370
Fax: +353 (0) 402 20375

Israel

SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

SIGMA-ALDRICH S.r.l.
Numero Verde: 800 827018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

SIGMA-ALDRICH JAPAN K.K.
Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

SIGMA-ALDRICH KOREA
Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: 52 722 276 1600
Fax: 52 722 276 1601

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

SIGMA-ALDRICH NEW ZEALAND LTD.
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

SIGMA-ALDRICH NORWAY AS
Tel: (+47) 23 17 60 60
Fax: (+47) 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

SIGMA-ALDRICH RUS, LLC
Tel: +7 (495) 621 6037
+7 (495) 621 5828
Fax: +7 (495) 621 5923

Singapore

SIGMA-ALDRICH PTE. LTD.
Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

Slovakia

SIGMA-ALDRICH spol. s r. o.
Tel: (+421) 255 571 562
Fax: (+421) 255 571 564

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD.
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313
SAFC (UK) Tel: 01202 712305

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Call Collect: (+1) 314 771 5750
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Vietnam

SIGMA-ALDRICH PTE LTD. VN R.O.
Tel: (848) 3516 2810
Fax: (848) 6258 4238

Internet

sigma-aldrich.com

World Headquarters

3050 Spruce St., St. Louis, MO 63103
(314) 771-5765
sigma-aldrich.com

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052

Technical Service (800) 325-5832 • sigma-aldrich.com/techservice

Development/Bulk Manufacturing Inquiries SAFC® (800) 244-1173

©2009 Sigma-Aldrich Co. All rights reserved. SIGMA, SAFC, SAFC®, SIGMA-ALDRICH, ALDRICH, FLUKA, and SUPELCO are trademarks belonging to Sigma-Aldrich Co. and its affiliate Sigma-Aldrich Biotechnology, L.P. Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in the Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.



Mixed Sources

Product group from well-managed
forests, controlled sources and
recycled wood or fiber
www.fsc.org Cert no. SGS-COC-003338
© 1996 Forest Stewardship Council



Accelerating Customers'
Success through Innovation and
Leadership in Life Science,
High Technology and Service

DEDICATED TO DR. ALFRED BADER ON THE OCCASION OF HIS 85TH BIRTHDAY

Aldrichimica ACTA

VOL. 42, NO. 1 • 2009

 **ALDRICH**
Chemistry



The Super Silyl Group in Diastereoselective Aldol
and Cascade Reactions

Iterative Cross-Coupling with MIDA Boronates: towards a
General Strategy for Small-Molecule Synthesis

SIGMA-ALDRICH

Editorial

"Please Call Me Alfred"

On April 28 of this year, Dr. Alfred Bader, undeniably the world's best-known Chemist Collector, celebrated his 85th birthday. Alfred's amazing life story has been covered extensively in books, magazines, and lectures, by Alfred himself and by others, and need not be repeated here. Furthermore, most of our readers are undoubtedly aware of Alfred's strong connections, past and present, not only to Sigma-Aldrich, but also to the *Aldrichimica Acta*, which he has showered with his attention for many years. Many of Alfred's paintings have graced the covers of the *Acta*, including this issue, which is featuring one of Alfred's favorite paintings.

We honor and thank Alfred for his outstanding contributions to Sigma-Aldrich and to the worlds of chemistry, business, and art. To what he calls the "ABC" (Art, Bible, Chemistry) of his life, one should add a "D" for Donating. Alfred's philanthropic activities are considerable and ongoing, and cover a wide range of causes that are near and dear to Alfred's heart.

The magnitude of Alfred's philanthropic efforts was made possible by the considerable financial rewards that he has reaped from two of his lifelong passions. The first is a spectacularly successful chemical business (Aldrich and Sigma-Aldrich) that he helped found and successfully managed for years. The second is his passion for collecting art works, particularly of Dutch and Flemish Masters, as well as rare stamps. His collection of about two hundred such paintings has been donated by Alfred to the Agnes Etherington Art

Centre of Queen's University in Kingston, Ontario (Canada). This gift is one of several sizeable ones that he has made to Queen's in gratitude for the education he received there in the 1940s.

Lesser known, but not less important to Alfred, is his Bible scholarship, particularly of Old Testament themes, which he has studied all his life and taught for a good many years. Another lesser known trait of Alfred is his modesty and unassuming lifestyle, which used to lead many a new employee of Aldrich, in the days when Alfred was company president, to think, when running into him for the first time, that he was just another employee. After meeting Alfred for the first time and calling him Dr. Bader, the editor of this publication was gently chided for calling him Dr. Bader, rather than Alfred, which is what he insists on being called even by people who don't know him that well.



Dr. Alfred Bader in 2005

Photo Courtesy of Alfred Bader

On behalf of all Sigma-Aldrich employees, past and present, we wish Alfred a very happy 85th birthday and many more in years to come.

SAFC® Launches New Vendor Audit Services


On March 19, 2009, SAFC® announced details of its new Vendor Audit Services, designed to extend audit reach for customers that need to secure global supply chain integrity through enhanced traceability and quality confirmation for critical raw materials and key intermediates. SAFC's Vendor Audit Services are an extension of in-house auditing programs, and provide pharmaceutical and biotechnology firms with a team of highly trained auditors that have extensive regulatory compliance expertise and are on the ground in local regions, including China

Standard audits available include:

- ISO 9001:2008 (previously ISO 9001:2000)
- ISO 13485:2003
- 21 CFR 211
- 21 CFR 820
- ICH Q7A
- GAMP 5
- GLP
- 21CFR58

SAFC can also perform audits to customer specifications upon request.

SAFC's Typical Audit Services Timeline



Please contact us at safcglobal.com/auditservices for more information.

Aldrichimica ACTA

VOL. 42, NO. 1 • 2009

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-231-8327
SAFC® 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 414-438-3850
24-Hour Emergency 414-438-3850
Web Site sigma-aldrich.com
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201, USA

Subscriptions

To request your FREE subscription to the
Aldrichimica Acta, please contact us by:

Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 2988
Milwaukee, WI 53201-2988
Email: sams-usa@sial.com

International customers, please contact your local
Sigma-Aldrich office. For worldwide contact infor-
mation, please see the back cover.

The *Aldrichimica Acta* is also available on the
Internet at sigma-aldrich.com/acta.

Aldrich brand products are sold through Sigma-
Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its
products conform to the information contained in
this and other Sigma-Aldrich publications. Purchaser
must determine the suitability of the product for its
particular use. See reverse side of invoice or packing
slip for additional terms and conditions of sale.

All prices listed in this publication are subject to
change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication
of Aldrich. Aldrich is a member of the Sigma-Aldrich
Group. © 2009 Sigma-Aldrich Co.


"PLEASE BOTHER US."



Joe Porwoil

Joe Porwoil, President
Aldrich Chemical Co., Inc.

Professor Stephen Buchwald of the Massachusetts Institute of Technology kindly suggested that we make the single-component palladacycle precatalysts of SPhos and XPhos. These complexes simplify the use of these Buchwald ligands in a pre-defined metal-to-ligand ratio. The precatalysts are air- and moisture-stable, and can achieve high yields in C-N cross-coupling reactions using catalyst loadings as low as 0.1 mol % and 10-minute reaction times.



Biscoe, M. R. et al. *J. Am. Chem. Soc.* 2008, 130, 6686

704946 SPhos-palladium(II)phenethylamine chloride (1:1 MTBE solvate)	250 mg	\$45.00
	1 g	150.00
704954 XPhos-palladium(II)phenethylamine chloride (1:1 MTBE solvate)	250 mg	\$69.00
	1 g	190.00

Naturally, we made these useful catalyst precursors. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the back cover.

TABLE OF CONTENTS

The Super Silyl Group in Diastereoselective Aldol and Cascade Reactions	3
<i>Matthew B. Boxer, Brian J. Albert, and Hisashi Yamamoto,* University of Chicago</i>	
Iterative Cross-Coupling with MIDA Boronates: towards a General Strategy for Small-Molecule Synthesis	17
<i>Eric P. Gillis and Martin D. Burke,* University of Illinois at Urbana-Champaign</i>	

ABOUT OUR COVER


Two Squirrels (oil on panel, 33.0 × 40.6 cm) was painted around 1616 and is attributed to the Flemish painter Jan Brueghel (or Bruegel) the Elder (1568–1625), also known as "Velvet" Brueghel and "Flower" Brueghel. He was the second son of Pieter Brueghel the Elder and trained with his older brother, Pieter Brueghel the Younger, in the family workshop in Antwerp.

The painting, which reflects the painter's focus later on in life on painting flowers and animals, depicts two sprightly squirrels clutching a gnarly, twisting branch set against an empty sky. This simple, engaging, and humorous depiction of a pair of lively small animals parallels the practice, at that time, of hanging portraits of husband and wife. Its high level of finish suggests that it probably was executed on commission or for sale in the market. The luminous effect created by the thin pigment layers and the fine detail are reminiscent of the style of the Brueghel family.



Detail from *Two Squirrels* Photograph © Alfred Bader

This painting is part of the Bader Collection of Dutch and Flemish Paintings at the Agnes Etherington Art Centre of Queen's University, Kingston, ON, Canada.




New Products from Aldrich R&D

Aldrich is Pleased to Offer Cutting-Edge Tools for Organic Synthesis

Double-Allylation Reagents

Reagents bearing multifunctional handles are of interest for the preparation of complex molecules. Professor Hall and coworkers developed a new multifunctional reagent that provides high diastereo- and enantiocontrol in a number of reactions, such as in the nucleophilic addition to aldehydes.




Peng, F.; Hall, D. G. *J. Am. Chem. Soc.* **2007**, *129*, 3070

(+)-Allylboronic acid pinanediol ester, 97%

694584

C₁₃H₂₁BO₂
FW: 220.12




1 g \$34.60
5 g 112.00

(+)-Vinylboronic acid pinanediol ester, 95%

691615


C₁₂H₁₉BO₂
FW: 206.09



1 g \$66.00
5 g 218.50

Substrate for Nickel-Catalyzed Negishi Coupling

Chiral building blocks are of the utmost importance in the synthesis of more complex molecules. Professor Fu and coworkers devised the first catalytic enantioselective cross-coupling of secondary α -bromo amides with organozinc reagents. This new method proved to be highly selective for the coupling of unfunctionalized and functionalized organozincs with good yields




Fischer, C.; Fu, G. C. *J. Am. Chem. Soc.* **2005**, *127*, 4594

N-Benzyl-2-bromo-N-phenylbutanamide, 97%

704911


[851073-30-8]
C₁₇H₁₉BrNO
FW: 332.23



5 g \$65.00

Copper(I) Fluoride Complex for Aldol Reaction

Chiral tertiary alcohols are important building blocks for the synthesis of more complex molecules such as biologically active compounds or potential drugs. Shibasaki and coworkers developed a new copper fluoride catalyzed aldol reaction of ketones using ketene silyl acetals. Various aromatic ketones were screened and led to the desired aldol products in good yields and high selectivity.



Oisaki, K.; Zhao, D.; Kanai, M.; Shibasaki, M. *J. Am. Chem. Soc.* **2006**, *128*, 7164. Taniaphos is a registered trademark of OMG AG and Co.

Fluorotris(triphenylphosphine)copper(I), 95% by HMR

706000


C₅₄H₄₅CuFP₃
FW: 869.40



250 mg \$39.00
1 g 95.00

Copper Chloride-Bis(lithium chloride) Solution for Transmetalation

Cross-coupling reactions are essential tools for chemists. In particular, the amination of aromatic halides has become a method highly relied upon to prepare aryl amines. Knochel and coworkers developed a new general amination procedure using amidocuprates. This method proved to be very versatile resulting in good yields of the amine products. This new method is a good complement to the Pd-catalyzed amination reactions.



Del Amo, V.; Dubbaka, S. R.; Krasovskiy, A.; Knochel, P. *Angew. Chem., Int. Ed.* **2006**, *45*, 7838.

Copper(I) chloride-bis(lithium chloride) complex, 1 M in tetrahydrofuran

701742

CuCl·2(LiCl)
FW: 183.79




50 mL \$95.00

The Super Silyl Group in Diastereoselective Aldol and Cascade Reactions



Dr. Matthew B. Boxer



Dr. Brian J. Albert

Matthew B. Boxer, Brian J. Albert, and
Hisashi Yamamoto*
Department of Chemistry
The University of Chicago, GHJ 409
5735 S. Ellis Avenue
Chicago, IL 60637, USA
Email: yamamoto@uchicago.edu



ACS Photo © Peter Cutts Photography, LLC

Professor Hisashi Yamamoto (second from right) receiving the Sigma-Aldrich sponsored 2009 ACS Award for Creative Work in Synthetic Organic Chemistry. Pictured with Professor Yamamoto are Dr. Thomas H. Lane (right), 2009 ACS President, and Dr. Joseph S. Francisco (left), 2009 ACS President-Elect. Presenting the award on behalf of Sigma-Aldrich is Dr. Mark Redlich (second from left), Product Line Manager, Chemical Synthesis.

Outline

1. Introduction
2. Installation of the Super Silyl Group
3. [2 + 2] Cyclizations
4. The Mukaiyama Aldol Reaction
 - 4.1. α -Substituted Enol Ethers
 - 4.2. Acetaldehyde-Derived Enol Ethers
 - 4.3. Ketone-Derived Super Silyl Enol Ethers (SEEs)
5. Sequential Aldol–Aldol Reactions
 - 5.1. Cascade Reactions of Aldehyde-Derived SEEs
 - 5.2. Cascade Reactions of Ketone- and Aldehyde-Derived SEEs
6. Other Sequential Reactions of SEEs
 - 6.1. Silyl Dinucleophile Reagents
 - 6.2. Sequential Aldol–Carbanion Addition Reactions
 - 6.3. The Super Silyl Group in Four-Component Reactions
7. Protodesilylation–Self-Repair Mechanism
8. Transition-State Calculations
9. Stability and Cleavage of Super Silyl Ethers
10. Conclusions and Outlook
11. Acknowledgements
12. References and Notes

1. Introduction

Silyl groups constitute a very important and distinct class of protective groups and serve as active participants in many thermal, and acid- and base-catalyzed reactions,¹ such as the Mukaiyama aldol, Mannich, Hosomi–Sakurai, and many cyclization reactions.²

Their unique and important reactivity profile is manifested in several ways: (i) The silyl ethers, silyl enol ethers (SEEs), and allylsilanes of a large variety of substrates can be prepared under mild conditions.^{1,3} (ii) A wide stability range exists for the many commercially available silyl-group-containing reagents.¹ Specifically, smaller silyl groups, such as TMS and TBS, can be cleaved under certain acidic or basic conditions, while larger silyl groups, like TIPS and TPS, are stable to those same conditions. This difference in reactivity allows for the development of selective protection–deprotection strategies. (iii) The particularly low electronegativity of silicon (1.8 vs 2.5 for carbon on the Pauling electronegativity scale) permits the stabilization of positive charges on silicon, while its large size (compared to carbon) stabilizes negative charges through polarization.^{2,4} These properties, combined with their aforementioned thermal and chemical stabilities, have made silyl groups remarkably useful in many synthetic reactions. Lastly, one of the most distinct features of silyl groups is their facile and selective cleavage by the fluoride ion—conditions that generally do not affect the rest of the organic molecule.¹

One silyl group that has been occasionally employed as a protective group, but rarely outside of radical reactions, is tris(trimethylsilyl)silyl $\{[(\text{CH}_3)_3\text{Si}]_3\text{Si}; \text{TTMSS}\}$, which is also known in the literature as the hypersilyl, sisy, or super silyl group. In this review, we will use the term super silyl that was coined by Hans Bock in 1993.⁵ Our group has found that the super silyl group exhibits peerless reactivity in a variety of diastereoselective C–C bond-forming reactions,⁶ which is the main topic of this article, as it has not been reviewed to date.

While many literature reports of the application of the super silyl group in chemistry have been published, the majority (61%) of these reports involve its use in radical reactions (**Figure 1**).⁷ Tris(trimethylsilyl)silane (TTMSSH) won Fluka's Reagent of the Year Prize in 1990 for its use in radical reactions as an alternative to the toxic (*n*-Bu)₃SnH. By substituting the three alkyl groups of a silane with three TMS groups, a Si-H bond weakening of 11 kcal/mol is observed ($\Delta_r H^\circ$: TESH = 90 kcal/mol vs TTMSSH = 79 kcal/mol).⁸ TTMSSH has been employed in hydrosilations of alkenes and alkynes,⁹⁻¹¹ radical cyclization






Figure 1. Pie Chart Showing the Breakdown by Application of Literature References Citing the Tris(trimethylsilyl)silyl (TTMSS) Group. (Ref. 7)



Scheme 1. Synthesis of Super Silyl Ethers under Standard Conditions.



Scheme 2. Radical Hydrosilylation Reactions. (Ref. 12)

reactions,⁹ the reduction of acid chlorides,^{9c} the hydrosilylation of carbonyl groups,^{9a,12} and the reduction of carbon-halogen bonds.^{9b}

The second most commonly reported application of the super silyl group deals with its complexation with transition metals and main group elements (14%).¹³ Its extremely bulky size and electron-donating ability stabilize metals and main group elements in various oxidation states. In 2004, Kornev published a review on the tris(trimethylsilyl)silyl group, with emphasis on its complexation with transition metals and main group elements and on its application in radical reactions.^{13c}

Surprisingly, the employment of the super silyl group in nonradical, standard acid-base, catalytic, and thermal reactions comprises only around 3% of the total literature.^{13c,14,15} Much of this 3% is devoted to papers describing its photochemical and strong-acid- or strong-base-induced rearrangements, which are rarely accompanied by C-C-bond formation. A small portion of the 3% is devoted to its application as an alcohol protecting group.¹⁵ Recently, the super silyl group has been incorporated into organic catalysts to increase their steric requirement and selectivity.¹⁶

The major structural difference between the super silyl group and the more typical silyl groups, which all contain only Si-C and Si-O bonds, is the presence of the Si-Si bonds. A number of experimental and theoretical calculations attribute many of the unique properties of this silyl group to its Si-Si bonds,¹⁷ which render it UV active (allowing for straightforward TLC analysis) and impart other distinct electronic properties. These unique characteristics are demonstrated in the thoroughly studied polysilanes and oligosilanes, where the electrons are actually delocalized along the Si-Si σ bonds.¹⁷ In these cases, photolysis reactions arise by promoting an electron from a σ orbital to the σ^* orbital, which often takes place with remarkably low excitation energy (320 nm).

Recently, the unique reactivity of the super silyl group has led to its successful application in a variety of C-C-bond-forming reactions: [2 + 2] cyclizations, Mukaiyama aldol synthesis, and various sequential reactions.⁶

2. Installation of the Super Silyl Group

The protection of alcohols with the super silyl group can be easily achieved under standard conditions (**Scheme 1**).^{7b,15,18} Brook and co-workers have demonstrated that protection of various primary and secondary alcohols can be accomplished using 1 equivalent of the alcohol, 1 equivalent of super silyl chloride, and 1.2 equivalents of DMAP.¹⁵ We have found that unhindered primary alcohols can be efficiently protected by employing super silyl chloride and triethylamine in either THF or CH₂Cl₂.^{7b} We have also been able to protect β -hydroxy esters using super silyl chloride and imidazole in DMF.^{7b} Super silyl triflate has also been utilized in conjunction with triethylamine to selectively protect carbohydrates.¹⁸

Because of the ability of TTMSSH to participate in hydrosilylation reactions, the radical reaction conditions shown in **Scheme 2** have been utilized to prepare super silyl ethers in one step from carbonyl compounds.¹²

The initial synthesis of the first super silyl enol ether (super SEE) utilized an *n*-BuLi induced THF fragmentation to generate the lithium enolate of acetaldehyde (**Scheme 3**).¹⁹ The synthesis of super SEEs evolved over time to use a metal-halogen exchange of AgOTf and TTMSSCl to generate the silyl triflate, which was then employed under soft enolization conditions to form **1**. The synthesis of TTMSSOTf was further modified

and took advantage of a triflic acid–super silane reaction that liberated H₂ gas and generated the silyl triflate in situ, which was again utilized to prepare **1**.^{6a-c}

3. [2 + 2] Cyclizations

The cyclobutane ring is found in various natural and unnatural bioactive molecules, and is commonly found in synthetic intermediates as well.²⁰ There is only a limited number of reports of efficient syntheses of this 4-membered ring by non-photochemical [2 + 2] cycloadditions,^{20–22} in large part because one of the synthetic routes, the concerted cyclization, is disallowed by the Woodward–Hoffmann rules.²³ Thus, both thermal and ground-state catalytic reactions must proceed stepwise, which forms the basis for the proposed Michael aldol mechanism (Scheme 4).^{6a,22a} This mechanism invokes a zwitterionic Michael addition intermediate, **A**, which may be present long enough to give rise to undesirable side-products. The use of aldehyde enol equivalents for this type of reaction had been previously unattainable, presumably due to this fact.

Takasu, Ihara, and co-workers have published numerous reports on formal [2 + 2] cyclizations of SEEs derived from ketones.²² Their original report took advantage of in situ SEE formation with a TMSI–HMDS system and the use of a chiral auxiliary (8-phenylmenthyl ester).^{22a} They later employed preformed TBS and TIPS enol ethers in cyclizations with hexafluoroisopropyl acrylate catalyzed by EtAlCl₂,^{22c,d} which worked well for forming cyclobutane-containing products including 5,4-, 6,4-, 7,4-, and 8,4-fused bicyclic compounds. A few years later, they reported the use of a strong Bronsted acid, Tf₃NH, to effect the [2 + 2] cyclization of SEEs and methyl acrylate.^{22c}


The same research group also investigated formal [2 + 2] cyclizations of α,β -unsaturated esters with acetaldehyde-derived SEEs, and reported that TBS^{22d} and TIPS SEEs^{22c} failed to give any desired products (eq 1). To prevent side reactions potentially stemming from high-energy intermediates such as **A**, it was envisaged that a bulky super silyl group, such as **1**, could stabilize and “shield” **A** from undesired reactions and decomposition. Fortunately, using **1**, hexafluoroisopropyl acrylate, and EtAlCl₂ as catalyst gave the desired cyclobutane in 45% yield with low diastereoselectivity.^{6a} Interestingly, the use of the pentamethyldisilyl (PMDS) derived enol ether **2**, containing one Si–Si bond, also afforded the cyclobutane adduct, albeit in 7% yield. The effect of the R group of the ester was investigated next: aliphatic esters gave no product, whereas phenyl acrylate gave the best yield and diastereoselectivity.

Various Lewis acids were screened, but only the EtAlCl₂-based catalyst gave acceptable results, while TiCl₄, GaCl₃, SnCl₄, AgNTf₂, TMSOTf, and HNTf₂ all gave <10% of the desired product. It has been observed that unwanted transfer of silyl groups during attempted asymmetric aldol synthesis could be prevented by using a bulky Lewis acid with the triflimide counteranion.²⁴ The use of bulky catalysts based on the methylaluminum bis(2,6-diphenylphenoxide) (MAPH)²⁵ scaffold were also investigated. The best result was obtained with bis(2,6-diphenylphenoxide) aluminum triflimide (BDAT).


Both aldehyde- and ketone-derived super SEEs were shown to succeed in this system. The reaction proceeded smoothly in all cases at –40 °C with a 3 mol % catalyst loading and gave high yields with high trans:cis ratios (eq 2).^{6a} Acetaldehyde SEE afforded cyclobutane **4** with a >99:1:0:0 dr by use of the chiral *trans*-(1*R*,2*S*)-2-phenylcyclohexanol derived ester (entry 2). SEE **5** led to *gem*-dimethylcyclobutane **6**, with excellent stereoselectivity (entry 3), while super SEE **7** gave cyclobutane **8**, which

contains a chiral quaternary carbon (entry 4). Interestingly, the highest yield (94%) and excellent selectivity were observed for cyclohexanecarboxaldehyde SEE **9** (entry 5). Three contiguous stereocenters were formed with excellent diastereoselectivity from the reaction of *E*-**11** with phenyl acrylate (entry 6).


The significance of these super SEEs was clearly demonstrated by their success in the [2 + 2] cyclizations. That is, these examples represented the first formal [2 + 2] condensations between an acetaldehyde SEE and an acrylic




Scheme 3. Evolution of the Synthesis of Super Silyl Enol Ether **1**. (Ref. 6a-c 19b)




Scheme 4. Proposed Michael Aldol Mechanism for the Formal [2 + 2] Cyclization. (Ref. 6a)



R	R'	Yield	Trans:Cis
TBS	a	0%	----
TIPS	(CF ₃) ₂ CH	0%	----
PMDS (2)	(CF ₃) ₂ CH	7%	1.5:1
TTMSS (1)	(CF ₃) ₂ CH	45%	2.0:1
TTMSS (1)	Et	0%	----
TTMSS (1)	Ph	58%	5.0:1



eq 1 (Ref. 6a,7b,22d)



1, 5, 7, 9, 11

3, 4, 6, 8, 10, 12

No	R ¹	R ² ,R ³	Yield	Trans:Cis
3	H	H,H	84%	10:1
4	H	H,H ^a	81%	>99:1:0:0
6	H	Me,Me	91%	>99:1
8	Me	H,H	81%	25:1
10	H	(CH ₂) ₅	94%	>99:1
12	H	Me,H ^b	93%	>99:1

^a *trans*-(1*R*,2*S*)-2-Phenylcyclohexyl acrylate was used. ^b *E:Z* (of starting silyl enol ether **11**) = 92:8

eq 2 (Ref. 6a)


ester. The necessity of the Si–Si bonds for these reactions was a first indication of the distinctiveness of the super silyl group in Lewis acid catalyzed reactions.

4. The Mukaiyama Aldol Reaction


The aldol reaction has emerged as a regular method for forming β -hydroxy carbonyl and/or 1,3-diol motifs typically seen in polyketides.²⁶ Of the various examples, a great deal utilize an ester, thioester, or ketone enolate (or enolate equivalent) as the nucleophile to circumvent problems associated with the aldehyde crossed aldol reaction. Frequently, the resulting products are reduced to the desired hydroxy-protected aldehydes through an additional one- or two-step procedure. Since Mukaiyama's seminal report on the titanium-catalyzed reaction of SEEs with aldehydes and ketones,²⁷ the Mukaiyama aldol reaction has developed into one of the most powerful and best-known synthetic reactions.²⁸ With a large number of reviews and books published in this area, this section will focus mainly on the aldehyde crossed aldol reaction.

4.1. α -Substituted Enol Ethers


Although the Mukaiyama aldol synthesis is one of the most powerful variants of the aldol reaction, the aldehyde crossed aldol reaction has only been realized in a few cases, many of which are limited in scope.^{29–31} In 1974, Mukaiyama followed his archetypal report with two examples of the isobutyraldehyde SEE crossed aldol reaction.^{29b} This system took advantage of



R	Catalyst	Yield	Syn:Anti
n-Hep	HNTf ₂	82%	80:20
n-Hep	TTMSSNTf ₂	85%	79:21
Cy	HNTf ₂	72%	85:15
Cy	TTMSSNTf ₂	71%	82:18
t-Bu	HNTf ₂	78%	95:5
t-Bu	TTMSSNTf ₂	79%	95:5
(S)-PhCH(Me)	HNTf ₂	84%	>95:5:0 ^a
(S)-PhCH(Me)	TTMSSNTf ₂	87%	>95:5:0 ^a



eq 3 (Ref. 6b)



R	Catalyst	Yield	Syn:Anti
n-Hep	HNTf ₂	87%	...
n-Hep	TTMSSNTf ₂	85%	...
Cy	HNTf ₂	89%	...
Cy	TTMSSNTf ₂	86%	...
t-Bu	HNTf ₂	90%	...
t-Bu	TTMSSNTf ₂	91%	...
(E,E)-MeCH=CHCH=CH	HNTf ₂	78%	...
(E,E)-MeCH=CHCH=CH	TTMSSNTf ₂	75%	...
Ph	HNTf ₂	83%	...
Ph	TTMSSNTf ₂	87%	...
(S)-PhCH(Me)	HNTf ₂	86%	>95:5
(S)-PhCH(Me)	TTMSSNTf ₂	85%	>95:5
(2R)-EiCH(Me)	HNTf ₂	93%	86:14
(2R)-EiCH(Me)	TTMSSNTf ₂	90%	85:15
(2S)-MeCH(OTIPS)CH ₂	HNTf ₂	88%	85:15
(2S)-MeCH(OTIPS)CH ₂	TTMSSNTf ₂	89%	88:12

eq 4 (Ref. 6b,7b)

the fact that neopentyl aldehydes were formed after the first aldol reaction and, thus, further aldol reactions were retarded.

Later use of α -substituted aldehydes took advantage of in situ enolate formation. In 1980, Heathcock reported the use of the in situ formed lithium enolate of propanal, which gave low selectivity when reacted with benzaldehyde.^{30a} In 1983, Kato and Mukaiyama reacted an in situ generated tin enolate of isobutyraldehyde with a variety of aryl and alkyl aldehydes.^{29c} Mahrwald's^{30b,d} and Oshima's^{30e,f} groups utilized in situ formed titanium enolates for reaction with a variety of aldehydes.

In 2001, Denmark reported that trichloro SEEs of propanal and heptanal successfully reacted with aromatic and aliphatic aldehydes in an enantioselective fashion catalyzed by a chiral phosphoramidate base.^{30g} The diastereoselectivity was controlled by the enol geometry, with the *Z* enols giving high syn selectivity and the *E* enols giving high anti selectivity.

After MacMillan's seminal report on the enantioselective proline-catalyzed aldehyde crossed aldol reaction,^{31a} an explosion of variants of this reaction emerged. Publications from Jørgensen, Córdova, and Barbas all employed proline-based catalysts to effect the direct aldehyde crossed aldol reaction of α -alkyl-substituted aldehydes.³¹

In 2006, our group published a diastereoselective reaction that worked very well with propionaldehyde-derived super SEE (**Z**)-**11**.^{6b} High yields and good syn:anti ratios were obtained using HNTf₂ as the precatalyst. This Brønsted acid is termed the precatalyst due to the fact that the use of (TTMSS)NTf₂ (0.05 mol %) as the catalyst in all reactions shown in equations 3 and 4 led to results identical to those obtained using HNTf₂, implying that the silyl triflimide is likely the true catalyst (see Section 7 for the proposed protodesilylation mechanism). Aliphatic and branched aldehydes successfully underwent this reaction, generating aldolates in moderate-to-high diastereoselectivities (eq 3).^{6b} The use of (*S*)-2-phenylpropanal exhibited high Felkin control in conjunction with syn selectivity, providing three adjacent stereocenters. Importantly, this provides a complementary method to the anti selectivity obtained by MacMillan.^{31a}

4.2. Acetaldehyde-Derived Enol Ethers

Interestingly, the direct installation of acetaldehyde had not been described in any broad sense before 2006. In 1958, Leech and co-workers reported the condensation of an alkyl vinyl ether with an aldehyde catalyzed by BF₃.^{29a} Paterson et al. later utilized the TBS enol ether of acetaldehyde in a reaction with a highly electrophilic oxonium ion in the presence of super stoichiometric Cl₂Ti(Oi-Pr)₂.³² There are a few reports on the use of the enzyme 2-deoxyribose-5-phosphate aldolase (DERA) for the aldol reaction of acetaldehyde with various aldehydes; however, the observed yields are rather low.³³

Denmark and Bui have reported a chiral-phosphoramidate-catalyzed enantioselective aldol reaction of the TMS-SEE of acetaldehyde with aldehydes.^{30b} High enantioselectivities and good yields were obtained for a variety of aryl aldehydes, but the products were isolated as the dimethyl acetals after addition of MeOH followed by NaHCO₃ to the reaction mixture.

In 2006, our group reported a broad, highly diastereoselective, aldehyde crossed aldol reaction of acetaldehyde super SEEs (eq 4).^{6b,7b} The use of HNTf₂ as the precatalyst gave consistently high yields of the aldol products with aliphatic, branched, aryl, and even $\alpha,\beta,\gamma,\delta$ -unsaturated aldehydes. Moreover, (*S*)-2-phenylpropanal was tested and showed extremely high Felkin selectivity. Somewhat amazingly, good selectivity was

obtained in the reaction with 2-methylbutanal, demonstrating the ability of this reaction to differentiate between methyl and ethyl groups. Syn selectivity was observed for the reaction of **1** with a β -siloxy aldehyde, in contrast to previous studies by Evans and co-workers, wherein anti selectivity was observed for the open-transition-state Mukaiyama aldol additions to β -alkoxy aldehydes.³⁴ The syn selectivity observed in our system is proposed to arise from the size of the TIPS and super silyl groups, whereas the study reported by Evans dealt with much smaller β -alkoxy aldehydes. A more detailed discussion of these observations, including DFT calculations, is presented in Section 8.

Recently, acetaldehyde has been introduced through the Mannich reaction. A singular example was given in the Yb(OTf)₃ catalyzed reaction of vinyloxytrimethylsilane.³⁵ List later published a paper on the direct use of acetaldehyde under proline catalysis.³⁶ While yields were relatively low, enantiomer ratios were very high for a range of substrates.

4.3. Ketone-Derived Super Silyl Enol Ethers (SEEs)

Ketone-derived super SEEs were also examined in simple Mukaiyama aldol reactions.⁶¹ The acetone-derived super SEE gave exceedingly high Felkin selectivity with 2-methylbutanal and 2-phenylpropionaldehyde (Scheme 5). The reaction of cyclohexanone super SEE with isobutyraldehyde gave the aldol product in high yield and unprecedented high anti selectivity for this type of Mukaiyama aldol reaction. This is in stark contrast to the TBS and TMS SEEs of cyclohexanone, which have been reported to give little-to-no selectivity in aldol reactions with a variety of catalysts.³⁷

5. Sequential Aldol–Aldol Reactions

For many obvious and necessary reasons, there is considerable interest in economical and environmentally friendly reactions.^{38,39} In this vein, one-pot, sequential and multicomponent reactions have emerged as an important means for accomplishing some of these goals. As Tietze and Beifuss wrote in a review on sequential reactions,³⁹ “if we compare our synthetic performance to date with that of Nature, then we must recognize that Nature is not just highly selective, but also very efficient, often employing sequential transformations. By this we understand a series of reactions steps in which several bonds are formed or broken, without the isolation of any intermediates”. They later go on to say, “the quantity of solvents and eluents required in comparison with stepwise processes is considerably reduced. Sequential reactions should, therefore, be more frequently included in future synthetic planning”. The use of sequential reactions not only saves bulk materials, but also time and labor.


Aldol–aldol reactions that incorporate more than three equivalents of starting aldehyde in the product have been described in discussions of aldolase-catalyzed reactions.^{33,40} As previously mentioned, these systems suffer from low yields and limited applicability. MacMillan expanded his proline-catalyzed aldehyde crossed aldol reaction to include a two-step synthesis of enantiopure hexoses.^{41,42} This elegantly designed system took advantage of the enantioselective aldol dimerization of α -siloxy- and α -benzyloxy aldehydes under proline catalysis. Various α -heteroatom-substituted TMS enol ethers were then added under TiCl₄ catalysis to give the cyclized hexoses in diastereo- and enantiomerically enriched forms, thereby providing an excellent means for protective-

group control in ¹³C-labelling experiments. Leighton and co-workers later used strained silacycles for a one-pot, ketone enol, aldol–aldol reaction, which provided hemiketal products containing two chiral quaternary carbon stereocenters.⁴³


5.1. Cascade Reactions of Aldehyde-Derived SEEs

Our group described a cascade aldol reaction, later termed sequential aldol–aldol (SAA) reaction.^{6b} Interestingly, this reaction stops at the 2:1 adduct (SEE:starting aldehyde) stage even when an excess of super SEE is employed. It is believed that after the first addition and silyl transfer, the steric encumbrance of the super silyl group kinetically slows down the addition of a second equivalent of SEE to a rate that does not compete with the rate of the first addition. When all of the aldehyde starting material has been consumed, a second addition occurs giving the 2:1 adducts with high diastereoselectivity. After the second aldol reaction, the substrate has two super siloxy groups (at the β and δ positions), and a third aldol reaction would require the coordination of a third super silyl group (as the catalyst), which should significantly slow the third aldol reaction due to steric hindrance.

The SAA reaction succeeded in generating a variety of protected β,δ -dihydroxy aldehydes in good yields and selectivities (eq 5).^{6b} Pival- and cyclohexyl aldehydes showed comparable selectivities. Octyl aldehyde led to a slightly lower selectivity, while (*S*)-2-phenylpropanal resulted in high selectivity for the all-syn isomer. The β -TIPSoxy aldehyde afforded the all-syn protected β,δ,ζ -tris-siloxy aldehyde, and



Scheme 5. Mukaiyama Aldol Reactions of Ketone-Derived Super SEEs. (Ref 61)




^a Major product has the opposite configuration at the β and δ carbons

eq 5 (Ref 6b)


α -benzyloxypropanal formed a β,δ,ϵ -trihydroxy aldehyde adduct—consistent with a mechanism involving chelation-controlled first aldol addition followed by a syn-selective second aldol addition.

5.2. Cascade Reactions of Ketone- and Aldehyde-Derived SEEs

The SAA reaction was extended to ketone- and aldehyde-derived super SEEs (Scheme 6),^{6d} which requires, in this case, the addition of



Scheme 6. SAA and SA-Hetero-Diels-Alder Reactions of Super SEEs. (Ref 6d)



^a Reactions run at room temperature.

Scheme 7. Sequential Aldol-Grignard Addition with Super SEE 1. (Ref 6d)

first the aldehyde and then the ketone SEEs due to their competitive reactivity with the initial aldehyde. Following the standard aldol reaction of **1** with pivalaldehyde, addition of the super SEE of cyclohexanone (**13**) to the reaction mixture gave the SAA adduct in excellent diastereoselectivity. A similar reaction sequence employing 2-phenylpropionaldehyde in the first step generated a product with three new stereocenters highly stereoselectively. Substituting acetophenone super SEE (**14**) in the second step of the same reaction cascade similarly yielded the product with high diastereoselectivity. Using a hetero-Diels-Alder reaction with a siloxy-diene, developed by Kozmin and Rawal,⁴⁴ as the second step in the cascade formed a dihydropyranone with comparable diastereoselectivity.

6. Other Sequential Reactions of SEEs

The 1,3-diol substructure is very common in many important medicinal compounds.²⁸ The majority of syntheses of molecules containing this motif require a multistep protocol to access the stereodefined diol. Significant contributions by Leighton's group and ours have employed SEEs in tandem and sequential reactions to access the 1,3-diol motif in one step.^{6b-f,43,45}

6.1. Silyl Dinucleophile Reagents

In 1999, Berrisford's group published a very interesting report of a silicon-tethered dinucleophile reagent that undergoes a sequential aldol-allylation reaction.⁴⁶ When this unique reagent was reacted with dimethyl acetals in the presence of $\text{BF}_3 \cdot \text{OEt}_2$, the aldol-allylation product was obtained in reasonable yield, but with rather low diastereoselectivity.

Leighton and co-workers later developed a number of strained silacycles that are very efficient at tandem aldol-allylation, aldol-crotylation, and aldol-aldol reactions.^{43,45} These reagents rely on an increase in Lewis acidity of the silicon center that arises from the strained 5-membered ring, thus enabling the aldol-allylation and aldol-crotylation to proceed without the need for a catalyst. Moderate yields and reasonable selectivities were obtained for various combinations of enol and allyl or crotyl moieties. The strained silacycle approach was also extended to include the use of ketone enol equivalents for the synthesis of tertiary carbinols.⁴³ In this latter case, the sequential aldol-allylation, aldol-methallylation, aldol-crotylation, and aldol-aldol reactions proceeded in good yields and selectivities producing a variety of relatively complex polyketide-like products in one pot.

6.2. Sequential Aldol-Carbanion Addition Reactions

Due to the very low catalyst loading (0.05–0.10 mol % of HNTf₂) utilized in the aldol reactions and the high diastereoselectivity obtained in the SAA reactions, it was anticipated that the addition of Grignards would be tolerated and would proceed with high stereoselectivity in generating secondary and tertiary carbinols, which still remains a challenging task in organic chemistry.⁴⁷ Following the aldol reaction of pivalaldehyde and **1**, allylmagnesium bromide was added to give the aldol-allylation product in 90:10 syn:anti selectivity in one pot (Scheme 7).^{6d} Chiral starting aldehydes 2-phenylpropionaldehyde and 3-TBSoxybutanal also led to good selectivities for the corresponding sequential SA-Grignard reaction product. The use of vinyl-, alkenyl-, and allylmagnesium halides generated synthetically useful allylic, propargylic, and homoallylic alcohols with high selectivity. Worthy of note is that enantioenriched (*S*)-2-phenylpropanal (98.3% ee) led to the major diastereomeric product in 97% ee, indicating no racemization during the acid-catalyzed aldol step.^{6d}

Polyhalomethanes (PHaMs) are small inexpensive molecules that are typically used as solvents (i.e., CH_2Cl_2 and CHCl_3). Halogens, of course, have a diverse reactivity profile and can render compounds electrophilic or nucleophilic.⁴⁸ Medicinal chemistry has revealed that fluorine- and chlorine-containing compounds have enhanced properties in biological settings as well as in crop management.⁴⁹ Despite considerable evidence of the importance of halogen-containing compounds, the stereocontrolled introduction of the polyhalomethyl group is not widely achievable.⁵⁰ To our knowledge, there exist only a few examples of the diastereoselective introduction of such groups.^{50c,h}


One of the most straightforward pathways to these PHaMs is the nucleophilic addition of polyhalomethylolithiums (PHaMLi's) to aldehydes (Scheme 8).^{6c} PHaMLi's are best generated by deprotonation of PHaMs with bulky lithium amides such as LiTMP (TMP = 2,2,6,6-tetramethylpiperidiny). A potential difficulty with this approach is dealing with the instability of such carbenoid-type species even at low temperatures.^{50,51} However, it has been shown that when a solution of an aldehyde or ketone, in the presence of excess PHaMs, is treated with a bulky lithium amide, the kinetically generated lithium carbenoid species reacts with the aldehyde or ketone before side reactions and decomposition occur.^{50a}

Thus, this reaction is a natural fit for our SA–nucleophile addition reaction sequence, since the requisite acid catalyst for the first step is present at only 0.05 mol % and the reaction proceeds in nonprotic solvents. Using our standard aldol reaction of **1** and 2-phenylpropionaldehyde, a variety of solvents and temperatures were screened for the diastereoselective sequential addition of dibromomethylolithium (Scheme 9).^{6c} 1,2-Dichloroethane was utilized as solvent for the aldol reaction since CH_2Cl_2 is competitively deprotonated under the subsequent PHaMLi-generating reaction conditions. For PHaMLi generation and addition reactions, employing THF as solvent at -100°C was optimal for producing the syn diols. The aldol solvent could also be utilized as the PHaMLi precursor. The highest selectivity was obtained with the largest such anion, tribromomethide, giving the product in high yield and high syn selectivity. Diiodomethane is also successfully deprotonated under these conditions and adds with high selectivity to give α -diiodomethylcarbinols in good yields.


A mixed α -polyhalomethylcarbinol has also been synthesized by Kuroboshi's group using a slightly different protocol.⁵² Following the standard aldol reaction, CFBr_3 was added, the solution diluted with 2:1 THF– Et_2O , cooled to -130°C , and Br_2FCl prepared in situ by lithium–bromine exchange with $n\text{-BuLi}$ (Scheme 10).^{6c} The α -dibromofluorocarbonol was produced in 55% yield with good selectivity, and was converted into the Z α -haloenol ester by treatment with acetic anhydride and then CrCl_2 in refluxing THF.⁵³ Moreover, a disubstituted Z fluoroalkene was prepared by an SA–Wittig-type olefination sequence.^{6c} This sequential reaction succeeded when the aldol reaction was followed by addition of the in situ prepared $(n\text{-Bu})_3\text{P}\text{-CF}\text{-P}(n\text{-Bu})_2\text{Cl}$.⁵⁴ After stirring for 12 h, 10% NaOH was added and stirred for an additional 12 h, inducing hydrolysis of the vinylphosphonium salt, thereby generating the Z fluoroalkene, in moderate yield and high selectivity.

With a clear indication that ketone-derived super SEEs were satisfactory in the basic Mukaiyama aldol reaction, their use in SA–Grignard addition reactions was also demonstrated.⁶¹ While a plethora of literature reports have been published on the diastereoselectivity of additions to β -oxygenated aldehydes,^{34,55}

significantly fewer reports can be found for the corresponding simple β -oxygenated ketones (not including hydrogenation–reduction reactions).^{47b,c,56} The majority of the reports that do exist, involve a β -hydroxy ketone and are proposed to undergo cyclic, 6-membered-ring transition states involving a Lewis acid catalyst or a metal from the organometallic species. While there is a report concerning syn selectivity for the methyl and butyl additions to β -TBSOxy-protected ketones,^{47b} few examples of this type of diastereoselective reaction exist. This likely arises from




Scheme 8. PHaMLi Generation and in Situ Reaction with Aldehydes. (Ref 6c)




R	R'	Solvent	$\text{CH}_2\text{-nX}_n$	Yield	dr
Ph	Me	$\text{ClCH}_2\text{CH}_2\text{Cl}$	CH_2Br_2	72%	85:5:5:5
Et	Me	CH_2Cl_2 ^a	CH_2Cl_2	72%	70:18:10:2
Bn	H	CHCl_3	CHCl_3	70%	90:10
Ph	Me	$\text{ClCH}_2\text{CH}_2\text{Cl}$	CHBr_3	70%	90:7:<2:<2
b	Me	$\text{ClCH}_2\text{CH}_2\text{Cl}$	CH_2I_2	76%	82:15:<2:<2
Ph	Me	$\text{ClCH}_2\text{CH}_2\text{Cl}$	CH_2I_2	76%	82:15:<2:<2

^a At -78°C . ^b R =


Scheme 9. Sequential Aldol–Polyhalomethylithium Addition Reactions. (Ref 6c)




Scheme 10. Synthesis of Mixed α -Polyhalomethyl Carbinols and Fluoroalkenes. (Ref 6c)



Scheme 11. Distinct Diastereomer Formation by Simple Substrate Choice. (Ref. 6f)



eq 6 (Ref. 6f)



Scheme 12. Total Synthesis of (+)-Cryptocarya Diacetate. (Ref. 6d)

two main factors: (i) ketones typically show lower selectivity and reactivity than aldehydes in many stereoselective reactions, and (ii) stereoselectivity induced by chirality at the β carbon is often lower than that induced by chirality at the α carbon.

Using the previously established, simple, one-pot SA–Grignard addition reaction protocol, acetone super SEE and 2-phenylpropionaldehyde underwent the aldol reaction initiated by 0.1 mol % of HNTf₂. PhMgBr was subsequently added leading to *anti*-**15** in good yield and excellent diastereoselectivity (**Scheme 11**).^{6f} Interestingly, the *anti* product was the major diastereomer, and this sense of stereoinduction results from nucleophilic attack on the π face of the carbonyl *opposite* to that which is seen for SAA and SA–Grignard additions utilizing **1**. DFT calculations were carried out to investigate the nature of the transition state (TS) in these reactions; the results will be discussed in Section 8. Aware that the choice of substrate could lead to predictable, distinct diastereomer formation, *syn*-**15** was prepared with high selectivity by simply employing acetophenone super SEE in the first (SA) step, and adding MeMgBr in the second step. Extending this idea to generate a tertiary carbinol with similarly sized substituents, *anti*-**16** was obtained from acetophenone super SEE and *p*-fluorophenylmagnesium bromide in a similar reaction sequence. Next, using 4'-fluoroacetophenone super SEE and phenylmagnesium bromide sequentially led to the expected isomer *syn*-**16**. These examples clearly demonstrate the TS control exhibited by the super siloxy substituent in these open-chain β -super-siloxy ketone intermediates, as well as the ability to generate the desired diastereomers by the judicious choice of the SEE and the Grignard reagent.^{6f}

The generality of this one-pot sequence was demonstrated by the success of a variety of ketone super SEEs, aldehydes, and Grignard reagents (eq 6).^{6f} The reaction scope was quite broad, giving products such as **17**, which contains three contiguous stereocenters. The use of vinyl and alkyl Grignards worked quite well, giving products **18–25** with good selectivity. The formation of **24** and **23** again showcases the super silyl group's powerful control of diastereoselection by first differentiating methyl and ethyl groups to give a large excess of the Felkin isomer, and then by stereoselectively controlling the Grignard addition via its presence in the β position. Product **21**, containing the valuable pyridine moiety, was generated by the in situ preparation of the heteroaryl Grignard through Knochel's powerful *i*-PrMgCl–2,6-dibromopyridine exchange reaction.⁵⁷

6.3. The Super Silyl Group in Four-Component Reactions

A four-component, one-pot reaction sequence was employed for the extremely concise synthesis of cryptocarya diacetate.^{6d} This compound is isolated from the bark of the South African plant, *Cryptocarya latifolia*, which is used for medicinal purposes.^{58,59} The synthesis was initiated with a one-pot, SA–Grignard addition–acylation sequence, which generated a diene in 63% yield along with a 24% yield of minor diastereomers (**Scheme 12**).^{6d} The use of Grubbs's second-generation catalyst in the ring-closing metathesis step gave a dihydropyranone, which was treated with HF–pyridine and then excess Ac₂O–pyridine to furnish cryptocarya diacetate. This total synthesis was accomplished in only three laboratory steps and 32% overall yield.

The success achieved with SAA and SA–Grignard addition sequences, led our group to combine these methods in a four-component, SAA–Grignard addition protocol.^{6f} The aldol reaction of **1** and cyclohexanecarboxaldehyde was followed by reaction

with a second ketone SEE and subsequent addition of MeMgCl, which gave the four-component product in moderate yield and high diastereoselectivity (**Scheme 13**).⁶¹ The same protocol was employed with **1**, isobutyraldehyde, 4'-fluoroacetophenone super SEE, and PhMgBr to give the product in moderate yield and selectivity. These four-component, one-pot reaction sequences truly showcase the power of this approach by producing relatively complex chiral architecture from simple starting materials.

7. Protodesilylation–Self-Repair Mechanism

Due to the prevalence of protodesilylation reactions of allyl silanes and SEEs under Brønsted acid conditions, we postulated that the reaction of super SEEs and aldehydes might be proceeding via this protodesilylation mechanism (**Scheme 14**).^{6b,60} In this case, HNTf₂ would actually be the precatalyst and (TTMSS)NTf₂ the true catalyst. This was supported by the use of (TTMSS)NTf₂ (generated from the exchange reaction of TTMSS–Cl and AgNTf₂) as the catalyst, which led to essentially identical results to those obtained using triflimide (see eq 3 and 4). This fact, in combination with the tolerance of an extremely low catalyst loading (S:C = 2000:1), led us to propose a self-repairing catalyst system, wherein the silyltriflimide can be generated and regenerated even in the presence of water or other protic Lewis bases.^{6b}

8. Transition-State Calculations


Much of the exceptional diastereoselectivity and control associated with the TTMSS group can likely be traced to its large steric size.^{2,61} The super silyl group is extraordinarily bulky, and has been stated to shield molecular skeletons with a “H₂C-skin.”²² It has been reported that this group has a local steric influence comparable to that of the *t*-Bu group,^{61a} and is among the strongest electron donors to π systems, lone-pair centers, and molecular cations.^{61b} The important work done by Evans's group on the selectivity of additions to β -oxygenated aldehydes, led us to initially propose TS **B** as the reason for the syn selectivity. This TS does not suffer the unfavorable steric interaction that is present in conformation **C** between the Lewis acid coordinated oxygen and the R group (**Figure 2**).^{6b,34}

This proposal was corroborated through DFT calculations at the B3LYP/6-31+G(d,p)//B3LYP/6-31G(d) level, which led to a TS structure particularly similar to **B**.⁶¹ Informative results were also gleaned regarding the reversal of selectivity observed for nucleophilic addition to β -super siloxy aldehydes versus that to β -super siloxy ketones. The calculations showed that the vinyl Grignard addition to the β -super siloxy aldehyde favored the syn pathway by 0.3 kcal/mol in the TS leading to the experimentally observed syn isomer. The calculations for the β -super siloxy methyl ketone gave a 2.6 kcal/mol preference in the TS for formation of the observed anti isomer. Significant steric repulsion between the methyl group of the ketone and the β -isopropyl group in the syn TS is a major reason for the divergence of TS energies. Akin to the work by Evans, the ketone and aldehyde transition states indicate a conformational preference that minimizes destabilizing electrostatic β -C–O and C=O dipole interactions.³⁴ In the case of the β -super siloxy aldehyde, the methyl group is replaced with a hydrogen, resulting in the oxygen being the larger atom (hydrogen vs oxygen in aldehyde and methyl vs oxygen in ketone). This leads to a preferable passing through the syn TS. A key feature that is revealed in these calculations is that the super silyl group creates a large “umbrella”-like structure under which the rest of the molecule aligns, and which restricts the conformational freedom of the remaining portion of the molecule.


The stereochemical outcome is then largely determined by the interaction of the carbonyl and its substituent (Me for ketone and H for aldehyde) with the medium-sized β -group (*i*-Pr). This is in contrast to typical open-chain ketones and aldehydes, which have much more freedom of rotation due to the absence of the “umbrella” effect, which is why we believe we see such high selectivities for these β -super siloxy carbonyl addition reactions.

9. Stability and Cleavage of Super Silyl Ethers

As mentioned in the introduction, silyl groups occupy a distinct place in protective-group chemistry.¹ Their ease of preparation, combined with the range of stabilities associated with the various commercially available silyl-group-containing reagents, allow for tailor-made syntheses, in which a distinct deprotection step can be planned well in advance. Significant work by Brook's group has shown that the super silyl group is a unique, photolabile protecting group that is stable to a range of typical synthetic conditions, such as Grignards, Wittig reagents, and oxidation conditions (Jones reagent).¹⁵ Through our research, we additionally found



Scheme 13. Four-Component SAA–Grignard Addition Reactions. (Ref 6f)



Scheme 14. Proposed Protodesilylation–Self-Repair Mechanism. (Ref 6b)




Figure 2. Original Transition-State Proposal for an Aldol Addition of **1** to a β -Super Siloxy Aldehyde. (Ref 6b)

that super silyl ethers are stable to (i) reducing agents (NaBH₄, DIBAL-H, and L-Selectride®); (ii) oxidizing reagents (SO₃-pyridine-DMSO, OsO₄, and Dess-Martin periodinane); and (iii) the Tebbe reagent and organometallic reagents derived from Ce, Mn, and Cu. Interestingly, the super silyl group is also stable towards CsF and KF-18-crown-6, but is cleaved with (*n*-Bu)₄NF in under 1 min. Super silyl ethers show limited stability in the presence of *n*-BuLi or LAH for prolonged periods of time.

Perhaps the most appealing aspect of the reactivity of the super silyl group is its photochemical lability. In 1997, Brook reported that the super silyl group could be cleaved in methanolic dichloromethane upon irradiation with UV light (eq 7).^{7b,15a} The absorption of super silane and the related super silyl ether derivatives at 254 nm was exploited to effect this deprotection. Moreover, the typical silyl group, TBS, wasn't deprotected under these conditions (254 nm irradiation in quartz cell). We decided to further explore this reactivity profile by irradiation of substrates containing both the super silyl group and other typical silyl groups.^{7b} 1,4-Butanediol was used as the starting material for the preparation of various disilylated compounds. The experimental setup was designed to test an extremely simple and practical application of this selective photochemical deprotection. While the use of a quartz round-bottom flask was necessary, the deprotection was carried out by irradiation with a common UV lamp designed for analysis of typical fluorescent TLC plates. All substrates tested gave high yields for the selective deprotection of the super silyl group.

10. Conclusions and Outlook

The super silyl group is superior in achieving many of the promises of one-pot sequential and multicomponent reactions. The reactivity profile of the super silyl group in a variety of C-C-bond-forming reactions is quite broad and useful. The super silyl group imposes significant steric bulk and possesses unique electronic properties that have enabled it to outperform commonly employed silyl groups in typical as well as atypical


reactions. In particular, its use in sequential reactions has assigned the super silyl group a distinctive place in efficient one-pot transformations. Thus far, super SEEs have succeeded in [2 + 2] cyclization reactions, Mukaiyama aldehyde crossed aldol reactions, SAA reactions, SA-hetero-Diels-Alder reactions, SA-Grignard addition reactions, and SAA-Grignard addition reactions. While a number of SEEs have succeeded in these cases, we believe that an array of sequentially added acidic and/or basic reagents and substrates can be combined with the SA system to generate complex chiral architecture in simple and efficient one-pot protocols. We intend to include these sequential reactions in future synthetic planning whenever we encounter the 1,3-diol- and 1,3,5-triol motifs (Figure 3).

11. Acknowledgements

We thank the University of Chicago and Merck Co. for their generous financial support of this work.

12. References and Notes

- (1) (a) Greene, W. T.; Wuts, P. G. M. In *Protective Groups in Organic Synthesis*, 3rd ed.; Greene, W. T., Wuts, P. G. M., Eds.; Wiley: New York, 1999. (b) White, J. D.; Carter, R. G. *Silyl Ethers*. In *Compounds of Group 15 (As, Sb, Bi) and Silicon Compounds*; Fleming, I., Ed.; Science of Synthesis: Houben-Weyl Methods of Molecular Transformations Series; Georg Thieme Verlag: New York, 2001; Vol. 4, pp 371-412. (c) Larson, G. L. In *The Chemistry of Organic Silicon Compounds*; Patai, S., Rappoport, Z., Eds.; The Chemistry of Functional Groups Series; Patai, S., Series Ed.; Wiley-Interscience: New York, 1989; Vol. 1, pp 763-808. (d) Weber, W. P. In *Reactivity and Structure Concepts in Organic Chemistry*; Springer: New York, 1983; Vol. 14. (e) Miura, K.; Hosomi, A. In *Main Group Metals in Organic Synthesis*; Yamamoto, H., Oshima, K., Eds.; Wiley-VCH: Weinheim, Germany, 2004; Vol. 2, pp 409-592.
- (2) Colvin, E. W. In *Silicon in Organic Synthesis*; Butterworths Monographs in Chemistry and Chemical Engineering Series. Butterworths: London, 1981.
- (3) (a) Kobayashi, S.; Manabe, K.; Ishitani, H.; Matsuo, J.-I. In *Compounds of Group 15 (As, Sb, Bi) and Silicon Compounds*; Fleming, I., Ed.; Science of Synthesis: Houben-Weyl Methods of Molecular Transformations Series; Georg Thieme Verlag: New York, 2001; Vol. 4, pp 317-369. (b) Birkofer, L.; Stuhl, O. In *The Chemistry of Organic Silicon Compounds*; Patai, S., Rappoport, Z., Eds.; The Chemistry of Functional Groups Series; Patai, S., Series Ed.; Wiley-Interscience: New York, 1989; Vol. 1, pp 655-761.
- (4) Corey, J. Y. In *The Chemistry of Organic Silicon Compounds*; Patai, S., Rappoport, Z., Eds.; The Chemistry of Functional Groups Series; Patai, S., Series Ed.; Wiley-Interscience: New York, 1989; Vol. 1, pp 1-56.
- (5) The term "super silyl" was coined by Hans Bock: Bock, H.; Meuret, J.; Ruppert, K. *Angew. Chem., Int. Ed. Engl.* **1993**, *32*, 414.
- (6) (a) Boxer, M. B.; Yamamoto, H. *Org. Lett.* **2005**, *7*, 3127. (b) Boxer, M. B.; Yamamoto, H. *J. Am. Chem. Soc.* **2006**, *128*, 48. (c) Boxer, M. B.; Yamamoto, H. *Nature Protocols* **2006**, *1*, 2434. (d) Boxer, M. B.; Yamamoto, H. *J. Am. Chem. Soc.* **2007**, *129*, 2762. (e) Boxer, M. B.; Yamamoto, H. *Org. Lett.* **2008**, *10*, 453. (f) Boxer, M. B.; Akakura, M.; Yamamoto, H. *J. Am. Chem. Soc.* **2008**, *130*, 1580.
- (7) (a) A comprehensive search of the literature utilizing SciFinder Scholar® and employing "tris(trimethylsilyl)silane", "tris(trimethylsilyl)silyl", "super silyl", "hypersilyl", and "silyl" as keywords was used to tabulate all references. These references



eq 7 (Ref. 7b,15a)




Figure 3. Anticipated Retrosynthetic Disconnections (All Would Be One-Pot Reactions).

- were then analyzed and classified into the categories shown in Figure 1. (b) Boxer, M. B. Ph.D. Thesis, University of Chicago, Chicago, IL, 2008.
- (8) The bond dissociation energy of X–H in TMS(Me)₂Si–H is 85 kcal mol, while that in (*n*-Bu)₂Sn–H is 74 kcal mol: Kanabus-Kaminska, J. M.; Hawari, J. A.; Griller, D.; Chatgililoglu, C. *J. Am. Chem. Soc.* **1987**, *109*, 5267 and reference 9 therein.
- (9) Select articles: (a) Alberti, A.; Chatgililoglu, C. *Tetrahedron* **1990**, *46*, 3963. (b) Ballestri, M.; Chatgililoglu, C.; Clark, K. B.; Griller, D.; Giese, B.; Kopping, B. *J. Org. Chem.* **1991**, *56*, 678. (c) Chatgililoglu, C. *Acc. Chem. Res.* **1992**, *25*, 188. (d) Kopping, B.; Chatgililoglu, C.; Zehnder, M.; Giese, B. *J. Org. Chem.* **1992**, *57*, 3994. (e) Ballestri, M.; Chatgililoglu, C.; Cardi, N.; Sommazzi, A. *Tetrahedron Lett.* **1992**, *33*, 1787. (f) Chatgililoglu, C.; Ferreri, C.; Gimisis, T. In *The Chemistry of Organic Silicon Compounds*; Rappoport, Z., Apeloig, Y., Eds.; The Chemistry of Functional Groups Series; Patai, S., Series Ed.; Wiley-Interscience: New York, 1998; Vol. 2, Parts 1–3, pp 1539–1579. (g) Postigo, A.; Kopsov, S.; Ferreri, C.; Chatgililoglu, C. *Org. Lett.* **2007**, *9*, 5159.
- (10) Miura, K.; Oshima, K.; Utimoto, K. *Bull. Chem. Soc. Jpn.* **1993**, *66*, 2356.
- (11) Liu, Y.; Yamazaki, S.; Yamabe, S. *J. Org. Chem.* **2005**, *70*, 556.
- (12) Kulicke, K. J.; Giese, B. *Synlett* **1990**, 91.
- (13) Select articles and reviews: (a) Elsner, F. H.; Woo, H.-G.; Tilley, T. D. *J. Am. Chem. Soc.* **1988**, *110*, 313. (b) Heyn, R. H.; Tilley, T. D. *Inorg. Chim. Acta* **2002**, *341*, 91. (c) Kornev, A. V. *Russ. Chem. Rev. (Engl. Transl.)* **2004**, *73*, 1065. (d) Yu, X.; Morton, L. A.; Xue, Z.-L. *Organometallics* **2004**, *23*, 2210. (e) Lerner, H.-W. *Coord. Chem. Rev.* **2005**, *249*, 781.
- (14) Select examples: (a) Naka, A.; Ishikawa, M.; Matsui, S.; Ohshita, J.; Kunai, A. *Organometallics* **1996**, *15*, 5759. (b) Ohshita, J.; Nekoda, E.; Masaoka, S.; Ishikawa, M. *J. Organomet. Chem.* **1997**, *544*, 49. (c) Berry, M. B.; Griffiths, R. J.; Howard, J. A. K.; Leech, M. A.; Steel, P. G.; Yufit, D. S. *J. Chem. Soc., Perkin Trans. 1* **1999**, 3645. (d) Naka, A.; Ohshita, J.; Kunai, A.; Lee, M. E.; Ishikawa, M. *J. Organomet. Chem.* **1999**, *574*, 50. (e) Ohshita, J.; Yoshimoto, K.; Iida, T.; Kunai, A. *J. Am. Chem. Soc.* **2001**, *123*, 8400.
- (15) (a) Brook, M. A.; Gottardo, C.; Balduzzi, S.; Mohamed, M. *Tetrahedron Lett.* **1997**, *38*, 6997. (b) Brook, M. A.; Balduzzi, S.; Mohamed, M.; Gottardo, C. *Tetrahedron* **1999**, *55*, 10027.
- (16) (a) Kosugi, Y.; Akakura, M.; Ishihara, K. *Tetrahedron* **2007**, *63*, 6191. (b) Zhang, D.; Yuan, C. *Tetrahedron* **2008**, *64*, 2480.
- (17) (a) Gilman, H.; Atwell, W. H.; Schwabke, G. L. *J. Organomet. Chem.* **1964**, *2*, 369. (b) Ramsey, B. G. *J. Organomet. Chem.* **1974**, *67*, C67. (c) West, R. *Pure Appl. Chem.* **1982**, *54*, 1041. (d) Sakurai, H. *Pure Appl. Chem.* **1987**, *59*, 1637. (e) Shimizu, M.; Oda, K.; Bando, T.; Hiyama, T. *Chem. Lett.* **2006**, *35*, 1022.
- (18) Flemming, A.; Mamat, C.; Köckerling, M.; Krempner, C.; Miethchen, R. *Synthesis* **2006**, 2685.
- (19) (a) Bates, R. B.; Kroposki, L. M.; Potter, D. E. *J. Org. Chem.* **1972**, *37*, 560. (b) Jung, M. E.; Blum, R. B. *Tetrahedron Lett.* **1977**, *18*, 3791.
- (20) (a) Nemoto, H.; Fukumoto, K. *Synlett* **1997**, 863. (b) Lee-Ruff, E.; Mladenova, G. *Chem. Rev.* **2003**, *103*, 1449.
- (21) Select examples: (a) Hayashi, Y.; Narasaka, K. *Chem. Lett.* **1989**, *18*, 793. (b) Ichikawa, Y.; Narita, A.; Shiozawa, A.; Hayashi, Y.; Narasaka, K. *J. Chem. Soc., Chem. Commun.* **1989**, 1919. (c) Narasaka, K.; Hayashi, Y.; Shimadzu, H.; Niihata, S. *J. Am. Chem. Soc.* **1992**, *114*, 8869. (d) Sweis, R. F.; Schramm, M. P.; Kozmin, S. A. *J. Am. Chem. Soc.* **2004**, *126*, 7442. (e) Ito, H.; Hasegawa, M.; Takenaka, Y.; Kobayashi, T.; Iguchi, K. *J. Am. Chem. Soc.* **2004**, *126*, 4520. (f) Avenozza, A.; Busto, J. H.; Canal, N.; Peregrina, J. M. *J. Org. Chem.* **2005**, *70*, 330.
- (22) (a) Takasu, K.; Ueno, M.; Ihara, M. *Tetrahedron Lett.* **2000**, *41*, 2145. (b) Takasu, K.; Ueno, M.; Ihara, M. *J. Org. Chem.* **2001**, *66*, 4667. (c) Takasu, K.; Ueno, M.; Inanaga, K.; Ihara, M. *J. Org. Chem.* **2004**, *69*, 517. (d) Takasu, K.; Nagao, S.; Ueno, M.; Ihara, M. *Tetrahedron* **2004**, *60*, 2071. (e) Inanaga, K.; Takasu, K.; Ihara, M. *J. Am. Chem. Soc.* **2005**, *127*, 3668.
- (23) For a seminal report on the Woodward–Hoffmann symmetry rules, see Woodward, R. B.; Hoffmann, R. In *The Conservation of Orbital Symmetry*; Verlag Chemie: Weinheim, Germany, 1970.
- (24) Hiraiwa, Y. Ph.D. Thesis, University of Chicago, Chicago, IL, 2003.
- (25) Marx, A.; Yamamoto, H. *Angew. Chem., Int. Ed.* **2000**, *39*, 178.
- (26) (a) Rychnovsky, S. D. *Chem. Rev.* **1995**, *95*, 2021. (b) Koskinen, A. M. P.; Karisalmi, K. *Chem. Soc. Rev.* **2005**, *34*, 677.
- (27) Mukaiyama, T.; Narasaka, K.; Banno, K. *Chem. Lett.* **1973**, *2*, 1011.
- (28) (a) Mukaiyama, T. *Org. React.* **1982**, *28*, 203. (b) Heathcock, C. H. In *Additions to C–X π -Bonds, Part 2*; Heathcock, C. H., Ed.; Comprehensive Organic Synthesis Series; Trost, B. M., Fleming, I., Series Eds.; Pergamon: Oxford, 1992; Vol. 2, pp 133–319. (c) Bach, T. *Angew. Chem., Int. Ed. Engl.* **1994**, *33*, 417. (d) Nelson, S. G. *Tetrahedron: Asymmetry* **1998**, *9*, 357. (e) Carreira, E. M. In *Comprehensive Asymmetric Catalysis*; Jacobsen, E. N., Pfaltz, A., Yamamoto, H., Eds.; Springer: Berlin, 1999; Vol. 3, pp 996–1066. (f) Mahrwald, R. *Chem. Rev.* **1999**, *99*, 1095.
- (29) Original reports on nonstereoselective variants of the reaction: (a) Hoaglin, R. I.; Kubler, D. G.; Leech, R. E. *J. Am. Chem. Soc.* **1958**, *80*, 3069. (b) Mukaiyama, T.; Banno, K.; Narasaka, K. *J. Am. Chem. Soc.* **1974**, *96*, 7503. (c) Kato, J.; Mukaiyama, T. *Chem. Lett.* **1983**, *12*, 1727. (d) Kohler, B. A. B. *Synth. Commun.* **1985**, *15*, 39.
- (30) Reports on stereoselective variants of the reaction: (a) Heathcock, C. H.; Buse, C. T.; Kleschick, W. A.; Pirrung, M. C.; Sohn, J. E.; Lampe, J. J. *J. Org. Chem.* **1980**, *45*, 1066. (b) Mahrwald, R.; Costisella, B.; Gündogan, B. *Tetrahedron Lett.* **1997**, *38*, 4543. (c) Mahrwald, R.; Gündogan, B. *Chem. Commun.* **1998**, 2273. (d) Mahrwald, R.; Costisella, B.; Gündogan, B. *Synthesis* **1998**, 262. (e) Maeda, K.; Shinokubo, H.; Oshima, K. *J. Org. Chem.* **1998**, *63*, 4558. (f) Uehira, S.; Han, Z.; Shinokubo, H.; Oshima, K. *Org. Lett.* **1999**, *1*, 1383. (g) Denmark, S. E.; Ghosh, S. K. *Angew. Chem., Int. Ed.* **2001**, *40*, 4759. (h) Denmark, S. E.; Bui, T. *J. Org. Chem.* **2005**, *70*, 10190.
- (31) Proline-catalyzed reactions: (a) Northrup, A. B.; MacMillan, D. W. C. *J. Am. Chem. Soc.* **2002**, *124*, 6798. (b) Bøgevig, A.; Kumaragurubaran, N.; Jørgensen, K. A. *Chem. Commun.* **2002**, 620. (c) Chowdari, N. S.; Ramachary, D. B.; Córdova, A.; Barbas, C. F., III. *Tetrahedron Lett.* **2002**, *43*, 9591. (d) Córdova, A. *Tetrahedron Lett.* **2004**, *45*, 3949. (e) Mase, N.; Tanaka, F.; Barbas, C. F., III. *Angew. Chem., Int. Ed.* **2004**, *43*, 2420. (f) Northrup, A. B.; Mangion, I. K.; Hetteche, F.; MacMillan, D. W. C. *Angew. Chem., Int. Ed.* **2004**, *43*, 2152. (g) Casas, J.; Engqvist, M.; Ibrahim, I.; Kaynak, B.; Córdova, A. *Angew. Chem., Int. Ed.* **2005**, *44*, 1343.
- (32) Paterson, I.; Smith, J. D.; Ward, R. A. *Tetrahedron* **1995**, *51*, 9413.
- (33) (a) Barbas, C. F., III; Wang, Y.-F.; Wong, C.-H. *J. Am. Chem. Soc.* **1990**, *112*, 2013. (b) Chen, L.; Dumas, D. P.; Wong, C.-H. *J. Am. Chem. Soc.* **1992**, *114*, 741. (c) Gijssen, H. J. M.; Wong, C.-H. *J. Am. Chem. Soc.* **1994**, *116*, 8422.
- (34) (a) Evans, D. A.; Duffy, J. L.; Dart, M. J. *Tetrahedron Lett.* **1994**, *35*, 8537. (b) Evans, D. A.; Dart, M. J.; Duffy, J. L.; Yang, M. G. *J. Am. Chem. Soc.* **1996**, *118*, 4322. (c) Evans, D. A.; Cee, V. J.; Siska, S. J. *J. Am. Chem. Soc.* **2006**, *128*, 9433.
- (35) Di Fabio, R.; Alvaro, G.; Bertani, B.; Donati, D.; Pizzi, D. M.; Gentile, G.; Pentassuglia, G.; Giacobbe, S.; Spada, S.; Ratti, E.; Corsi, M.; Quartaroli, M.; Barnaby, R. J.; Vitulli, G. *Bioorg. Med. Chem. Lett.* **2007**, *17*, 1176.
- (36) Yang, J. W.; Chandler, C.; Stadler, M.; Kampen, D.; List, B. *Nature* **2008**, *452*, 453.
- (37) (a) Ishihara, K.; Kondo, S.; Yamamoto, H. *J. Org. Chem.* **2000**, *65*, 9125. (b) Kobayashi, S.; Murakami, M.; Mukaiyama, T. *Chem. Lett.*

- 1985, 14, 1535. (c) Jankowska, J.; Mlynarski, J. *J. Org. Chem.* **2006**, 71, 1317.
- (38) (a) Ho, T.-L. *Tandem Organic Reactions*; Wiley: New York, 1992. (b) Hall, N. *Science* **1994**, 266, 32.
- (39) Tietze, L. F.; Beifuss, U. *Angew. Chem., Int. Ed. Engl.* **1993**, 32, 131.
- (40) (a) Gijzen, H. J. M.; Wong, C.-H. *J. Am. Chem. Soc.* **1995**, 117, 2947. (b) Gijzen, H. J. M.; Wong, C.-H. *J. Am. Chem. Soc.* **1995**, 117, 7585. (c) Heine, A.; DeSantis, G.; Luz, J. G.; Mitchell, M.; Wong, C.-H.; Wilson, I. A. *Science* **2001**, 294, 369. (d) Liu, J.; Wong, C.-H. *Angew. Chem., Int. Ed.* **2002**, 41, 1404.
- (41) Northrup, A. B.; MacMillan, D. W. C. *Science* **2004**, 305, 1752.
- (42) Around the same time, Córdova published a similar proline-catalyzed hexose synthesis; see reference 31g.
- (43) Wang, X.; Meng, Q.; Perl, N. R.; Xu, Y.; Leighton, J. L. *J. Am. Chem. Soc.* **2005**, 127, 12806.
- (44) Kozmin, S. A.; Rawal, V. H. *J. Org. Chem.* **1997**, 62, 5252.
- (45) (a) Wang, X.; Meng, Q.; Nation, A. J.; Leighton, J. L. *J. Am. Chem. Soc.* **2002**, 124, 10672. (b) Zacuto, M. J.; O' Malley, S. L.; Leighton, J. L. *Tetrahedron* **2003**, 59, 8889.
- (46) Frost, L. M.; Smith, J. D.; Berrisford, D. J. *Tetrahedron Lett.* **1999**, 40, 2183.
- (47) Select examples and reviews: (a) Lynch, J. E.; Eliel, E. L. *J. Am. Chem. Soc.* **1984**, 106, 2943. (b) Ukaji, Y.; Kanda, H.; Yamamoto, K.; Fujisawa, T. *Chem. Lett.* **1990**, 19, 597. (c) Weber, B.; Seebach, D. *Angew. Chem., Int. Ed. Engl.* **1992**, 31, 84. (d) Weber, B.; Seebach, D. *Tetrahedron* **1994**, 50, 6117. (e) Garcia Ruano, J. L.; Tito, A.; Culebras, R. *Tetrahedron* **1996**, 52, 2177. (f) Corey, E. J.; Guzman-Perez, A. *Angew. Chem., Int. Ed.* **1998**, 37, 388. (g) Yus, M.; Ramón, D. J. *Recent Res. Dev. Org. Chem.* **2002**, 6, 297. (h) Garcia, C.; Martin, V. S. *Curr. Org. Chem.* **2006**, 10, 1849. (i) Ren, H.; Dunet, G.; Mayer, P.; Knochel, P. *J. Am. Chem. Soc.* **2007**, 129, 5376. (j) Roy, S.; Sharma, A.; Dhotare, B.; Vichare, P.; Chattopadhyay, A.; Chattopadhyay, S. *Synthesis* **2007**, 1082.
- (48) (a) Heck, R. F. In *Additions to and Substitutions at C-C π -Bonds*; Semmelhack, M. F., Ed.; Comprehensive Organic Synthesis Series; Trost, B. M., Fleming, I., Series Eds.; Pergamon: Oxford, U.K., 1991; Vol. 4, pp 851-856. (b) Urch, C. J. In *Synthesis: Carbon with One Heteroatom Attached by a Single Bond*; Ley, S. V., Ed.; Comprehensive Organic Functional Group Transformations Series; Katritzky, A. R., Meth-Cohn, O., Rees, C. W., Series Eds.; Pergamon: Oxford, U.K., 1995; Vol. 2, pp 605-633. (c) *Handbook of Organopalladium Chemistry for Organic Synthesis*; Negishi, E., de Meijere, A., Eds.; Wiley: New York, 2002; Vol. 1 and 2.
- (49) (a) Naumann, K. *J. Prakt. Chem.* **1999**, 341, 417. (b) Isanbor, C.; O'Hagan, D. *J. Fluorine Chem.* **2006**, 127, 303. (c) Kirk, K. L. *J. Fluorine Chem.* **2006**, 127, 1013. (d) Lézé, M.-P.; Le Borgne, M.; Pinson, P.; Paluszczak, A.; Duflos, M.; Le Baut, G.; Hartmann, R. *W. Bioorg. Med. Chem. Lett.* **2006**, 16, 1134.
- (50) (a) Taguchi, H.; Yamamoto, H.; Nozaki, H. *J. Am. Chem. Soc.* **1974**, 96, 3010. (b) Villieras, J.; Perriot, P.; Normant, J. F. *Synthesis* **1979**, 968. (c) Mukaiyama, T.; Yamaguchi, M.; Kato, J. *Chem. Lett.* **1981**, 10, 1505. (d) Benner, J. P.; Gill, G. B.; Parrott, S. J.; Wallace, B. *J. Chem. Soc., Perkin Trans. 1* **1984**, 331. (e) Tabuchi, T.; Inanaga, J.; Yamaguchi, M. *Tetrahedron Lett.* **1986**, 27, 3891. (f) Teager, D. S.; Ward, H. D.; Murray, R. K., Jr. *J. Org. Chem.* **1993**, 58, 5493. (g) Hu, C.-M.; Tu, M.-H. *J. Fluorine Chem.* **1994**, 67, 9. (h) Concellón, J. M.; Bernad, P. L.; Pérez-Andrés, J. A. *Tetrahedron Lett.* **1998**, 39, 1409.
- (51) (a) Köbrich, G.; Akhtar, A.; Ansari, F.; Breckoff, W. E.; Büttner, H.; Drischel, W.; Fischer, R. H.; Flory, K.; Fröhlich, H.; Goyert, W.; Heinemann, H.; Hornke, I.; Merkle, H. R.; Trapp, H.; Zündorf, W. *Angew. Chem., Int. Ed. Engl.* **1967**, 6, 41. (b) Hoeg, D. F.; Lusk, D. I. *J. Am. Chem. Soc.* **1964**, 86, 928. (c) Kirmse, W. *Carbene Chemistry*, 2nd ed.; Academic Press: New York, 1971.
- (52) Kuroboshi, M.; Yamada, N.; Takebe, Y.; Hiya, T. *Synlett* **1995**, 987.
- (53) Bejot, R.; Tisserand, S.; Reddy, L. M.; Barma, D. K.; Baati, R.; Falck, J. R.; Mioskowski, C. *Angew. Chem., Int. Ed.* **2005**, 44, 2008.
- (54) Cox, D. G.; Gurusamy, N.; Burton, D. J. *J. Am. Chem. Soc.* **1985**, 107, 2811.
- (55) (a) Kiyooka, S.; Heathcock, C. H. *Tetrahedron Lett.* **1983**, 24, 4765. (b) Reetz, M. T.; Jung, A. *J. Am. Chem. Soc.* **1983**, 105, 4833. (c) Reetz, M. T.; Kessler, K.; Jung, A. *Tetrahedron Lett.* **1984**, 25, 729. (d) Reetz, M. T. *Angew. Chem., Int. Ed. Engl.* **1984**, 23, 556. (e) Heathcock, C. H. In *Asymmetric Organic Reactions*; Morrison, J. D., Ed.; Academic Press: New York, 1984; Vol. 3. (f) Davis, A. P.; Plunkett, S. J.; Muir, J. E. *Chem. Commun.* **1998**, 1797.
- (56) (a) Leitereg, T. J.; Cram, D. J. *J. Am. Chem. Soc.* **1968**, 90, 4019. (b) Ikeda, N.; Omori, K.; Yamamoto, H. *Tetrahedron Lett.* **1986**, 27, 1175. (c) Batra, M. S.; Aguilar, F. J.; Brunet, E. *Tetrahedron* **1994**, 50, 8169.
- (57) Krasovskiy, A.; Knochel, P. *Angew. Chem., Int. Ed.* **2004**, 43, 3333.
- (58) Isolation: Drewes, S. E.; Schlapelo, B. M.; Horn, M. M.; Scott-Shaw, R.; Sandor, P. *Phytochemistry* **1995**, 38, 1427.
- (59) Syntheses: (a) Jørgensen, K. B.; Suenaga, T.; Nakata, T. *Tetrahedron Lett.* **1999**, 40, 8855. (b) Hunter, T. J.; O'Doherty, G. A. *Org. Lett.* **2001**, 3, 2777. (c) Smith, C. M.; O'Doherty, G. A. *Org. Lett.* **2003**, 5, 1959. (d) Krishna, P. R.; Reddy, V. V. R. *Tetrahedron Lett.* **2005**, 46, 3905.
- (60) (a) Hollis, T. K.; Bosnich, B. *J. Am. Chem. Soc.* **1995**, 117, 4570. (b) Ishihara, K.; Hiraiwa, Y.; Yamamoto, H. *Synlett* **2001**, 1851. (c) Mathieu, B.; Ghosez, L. *Tetrahedron* **2002**, 58, 8219. (d) Ishihara, K.; Hiraiwa, Y.; Yamamoto, H. *Chem. Commun.* **2002**, 1564. (e) Hiraiwa, Y.; Ishihara, K.; Yamamoto, H. *Eur. J. Org. Chem.* **2006**, 1837.
- (61) (a) Frey, J.; Schottland, E.; Rappoport, Z.; Bravo-Zhivotovskii, D.; Nakash, M.; Botoshansky, M.; Kafory, M.; Apeloig, Y. *J. Chem. Soc., Perkin Trans. 2* **1994**, 2555. (b) Bock, H.; Meuret, J.; Baur, R.; Ruppert, K. *J. Organomet. Chem.* **1993**, 446, 113.

Trademarks: SciFinder Scholar® (American Chemical Society Corp.), L-Selectride® (Sigma-Aldrich Biotechnology, L.P., and Sigma-Aldrich Co.).


Keywords: super silyl group; silyl enol ethers; aldol reaction; sequential reactions; stereoselective reactions.

About the Authors

Matthew B. Boxer was born in 1980 in Peterborough, New Hampshire. He obtained a B.S. degree in chemistry in 2003 from the University of Vermont, working on DNA quadruplex intercalators under the supervision of Professor A. Paul Krapcho. In 2008, he earned a Ph.D. degree in chemistry from the University of Chicago, where he investigated the chemistry of the super silyl group under the direction of Professor Hisashi Yamamoto. He is currently working with Dr. Craig Thomas as a postdoctoral fellow at NIH.

Brian J. Albert was born in 1980 in Park Ridge, Illinois. He received a B.S. degree in chemistry in 2002 from the University of Illinois at Urbana-Champaign under the supervision of Professor J. A. Katzenellenbogen. In 2007, he obtained a Ph.D. degree in chemistry from the University of Pittsburgh under the direction of Professor K. Koide. He then moved to the University of Chicago, where he began postdoctoral studies on the chemistry of the super silyl group under Professor Yamamoto.

Hisashi Yamamoto was born in 1943 in Kobe, Japan. He received his bachelor's degree from Kyoto University and his Ph.D. degree from Harvard University under the mentorship of Professors H. Nozaki and E. J. Corey, respectively. His first academic position was as an assistant professor and lecturer at Kyoto University and, in 1977, he was appointed Associate Professor of Chemistry at the University of Hawaii. In 1980, he moved to Nagoya University, where he became Professor in 1983. In 2002, he moved to the University of Chicago as Arthur Holly Compton Distinguished Professor. His honors include: the Prelog Medal (1993), the Chemical Society of Japan Award (1995), the Max-Tishler Prize (1998), Le Grand Prix de la Fondation Maison de la Chimie (2002), National Prize of Purple Medal (Japan, 2002), Yamada Prize (2004), Tetrahedron Prize (2006), The

Karl-Ziegler Professorship (2006), The Japan Academy Prize (2007), Honorary Member of the Chemical Society of Japan (2008). He has also been named the 2009 recipient of the ACS Award for Creative Work in Synthetic Organic Chemistry. He has more than 500 original publications, 120 reviews and books, and 50 patents. He is on the board of editors or international advisory boards of over 20 international journals, and has given 154 plenary or invited lectures and 49 honorary lectures. His current interests are primarily the development of new synthetic reactions in the field of acid catalysis including designer Lewis acids, designer Brønsted acids, and a combination of these two acid systems. Recently, he has also become interested in a new field of asymmetric oxidations and metal catalyst design based on cis-beta configurations. 

Matthew B. Boxer, Brian J. Albert, and Hisashi Yamamoto




ALDRICH
Chemistry

Discover

Super Silyl Protecting Groups

Tris(trimethylsilyl)silane (**360716**) won the Fluka Reagent of the Year award in 1990 for its use as an alternative to tributylstannane (Bu_3SnH) in radical reactions. Currently, the tris(trimethylsilyl)silyl (TTMSS, or supersilyl) group is being utilized in carbon-carbon-bond forming reactions, where the researcher can take advantage of its unique reactivity that arises from its electronic properties and

large steric bulk. Professor Hisashi Yamamoto has employed supersilyl enol ethers in a variety of efficient one-pot sequential reactions. Sigma-Aldrich is pleased to offer acetaldehyde supersilyl enol ether (**698180**) and other tris(trimethylsilyl)silyl reagents that allow you to use the supersilyl protecting group in your research.



Boxer, M. B.; Yamamoto, H. J. *Am. Chem. Soc.* **2006**, 128, 48



360716



409405




698180

For more information on the applications of the supersilyl group in carbon-carbon-bond forming reactions please see the review by Boxer, Albert, and Yamamoto in this issue.

For product-related inquiries, please contact our Technical Services department at 800-231-8327 (USA) or aldrich@sial.com.

sigma-aldrich.com

SIGMA-ALDRICH



MIDA Boronates


For Suzuki–Miyaura Cross-Couplings

The Suzuki–Miyaura cross-coupling reaction is one of the most important and highly utilized reactions in organic chemistry, with applications in polymer science as well as in the fine chemicals and pharmaceutical industries. However, some classes of boronic acids are exceptionally unstable and susceptible to decomposition, which renders them inefficient in coupling reactions or makes long-term storage difficult. These limitations also make iterative cross-couplings challenging. Recently, Burke and coworkers have developed a method to allow such iterative couplings under mild conditions by attenuating the reactivity of boronic acids through complexation with *N*-methyliminodiacetic acid (MIDA). The MIDA ligand can be cleaved under mild conditions to liberate the corresponding boronic acid. MIDA boronates are easily handled, indefinitely bench-top stable under air, compatible with chromatography, unreactive under standard anhydrous cross-coupling conditions, even at temperatures up to 80 °C, and are stable to harsh reagents such as triflic acid and Jones reagent.

Palladium-catalyzed cross-coupling reactions are ideal methods for the synthesis of polyenes because of their stereospecificity and mildness. However, polyenylboronic acids are very unstable and therefore difficult to employ in the synthesis of polyenes by the Suzuki–Miyaura cross-coupling. In an exemplary demonstration of the stability and efficiency of MIDA boronates in iterative cross-couplings, Burke and coworkers utilized a common alkenyl, *trans*-2-bromovinylboronic acid MIDA ester (BB1, **703478**), to create a series of polyenyl building blocks. The MIDA boronate terminus is inert to Suzuki, Stille, and Heck couplings, yielding butadienyl MIDA boronates. Sonogashira and Negishi couplings, as well as Miyaura borylations also proved effective and yielded versatile bis-metallated lynchpin-type reagents.

References:

- (1) Gillis, E. P.; Burke, M. D. *J. Am. Chem. Soc.* **2007**, *129*, 6716
- (2) Lee, S. J.; Gray, K. C.; Paek, J. S.; Burke, M. D. *J. Am. Chem. Soc.* **2008**, *130*, 466
- (3) Gillis, E. P.; Burke, M. D. *J. Am. Chem. Soc.* **2008**, *130*, 14084
- (4) Uno, B. E.; Gillis, E. P.; Burke, M. D. *Tetrahedron* **2009**, *65*, 3130



To view all the MIDA boronates offered by Sigma-Aldrich, see page 28 of this issue or visit sigma-aldrich.com/mida.

Iterative Cross-Coupling with MIDA Boronates: towards a General Strategy for Small-Molecule Synthesis



Mr. Eric P. Gillis



Professor Martin D. Burke

Eric P. Gillis and Martin D. Burke*
 Department of Chemistry
 University of Illinois at Urbana-Champaign
 454 Roger Adams Laboratory, Box 52-5
 600 South Mathews Avenue
 Urbana, IL 61801, USA
 Email: burke@scs.uiuc.edu

Outline

1. Introduction
2. Synthesis of MIDA Boronates
3. Physical Properties of MIDA Boronates
4. Iterative Cross-Coupling (ICC) with Halogenated MIDA Boronates
5. Multistep Synthesis of Complex Boronic Acids from Simple MIDA Boronates
6. Conclusions and Prospects
7. Acknowledgements
8. References and Notes

1. Introduction

Many organic molecules are inherently modular in their constitution. With respect to the molecules found in living systems, this modularity is a direct consequence of the fact that nearly all biosynthetic systems are based on the iterative coupling of bifunctional building blocks. For example, polypeptides are built from amino acids, oligonucleotides are derived from nucleotide monomers, and oligosaccharides are stitched together from individual sugar units. Interestingly, most *small-molecule* natural products are similarly constructed by the iterative coupling of bifunctional building blocks: e.g., polyketides from malonyl-CoA or methylmalonyl-CoA units, nonribosomal peptides from amino acids, polyterpenes from isopentenyl pyrophosphate or dimethylallyl pyrophosphate, and fatty acids from malonyl-CoA.¹ Similarly, many man-made pharmaceuticals are also highly modular because they are constructed by using different reactions to assemble collections of small building blocks, typically cyclic and heterocyclic fragments and their associated appendages. Thus, modularity is a remarkably general feature of many of the molecules that are targeted for synthesis in the laboratory.

Despite this common modularity, the strategies utilized for making polypeptides, oligonucleotides, and oligosaccharides are very different from those typically used to prepare small molecules. Specifically, all of the former classes of compounds are almost always constructed via iterative coupling of suitably protected forms of their constituent monomers.² Organic polymers can be similarly prepared.³ Due to the powerfully

simple nature of this iterative coupling approach, these processes are now increasingly performed in a fully automated fashion.^{2,3} With peptides and oligonucleotides, the advanced development of such automation has made it possible for even nonchemists to routinely prepare these types of compound for a wide range of applications.

In stark contrast, it is typical for a synthetic chemist to develop a unique, customized strategy for each small molecule that is targeted for preparation in the laboratory. As a result, the synthesis of small molecules remains a relatively complex, unsystematized, and inflexible process that is practiced almost exclusively by highly trained specialists. Driven by the hypothesis that the inherent modularity in small molecules remains largely underutilized, we have established a research program that aims to develop a unified strategy for the construction of these compounds by the iterative coupling of bifunctional building blocks.⁴⁻⁷ Specifically, we have targeted the development of building blocks representing substructures that appear frequently in natural products and man-made pharmaceuticals and the chemistry that will enable their precise union via iterative, metal-mediated, cross-coupling reactions. In the idealized form of this "Iterative Cross-Coupling" (ICC) approach, building blocks having all of the required functional groups preinstalled in the correct oxidation state and with the desired stereochemical relationships are iteratively united using only stereospecific cross-coupling reactions (**Figure 1**). In addition to being simple, efficient, and potentially amenable to automation, the modularity of this approach makes it inherently well-suited for generating diverse collections of compounds simply by substituting modified building blocks into the same synthesis pathway. It is anticipated that the advanced development of this ICC strategy will substantially enable the laboratory synthesis of a wide range of natural products, pharmaceuticals, and organic materials, and may even extend the power of small-molecule synthesis to the nonchemist.

As described in this review, *N*-methyliminodiacetic acid (MIDA) boronates^{8,9} represent a highly promising platform for this type of synthesis strategy. These building blocks are remarkably convenient to prepare, analyze, purify, and store, and many are now commercially available. The MIDA boronate functional





Figure 1. Analogous Strategies for the Synthesis of Peptides and Small Molecules.



Scheme 1. Examples of Known Methods for the Synthesis of MIDA Boronates.

group is also stable to anhydrous cross-coupling conditions, but is easily hydrolyzed with mild aqueous base, thereby enabling the controlled ICC of B-protected "haloboronic acids".^{4–6} In addition, MIDA boronates are remarkably stable to a wide range of common reaction conditions and chromatography, which makes possible the facile preparation of complex borane building blocks from simple MIDA boronates via multistep synthesis.^{6,7} This review aims to enable the effective utilization of this platform and the ICC strategy to promote the simple, efficient, and flexible construction of small molecules.

2. Synthesis of MIDA Boronates

N-Methyliminodiacetic acid, $\text{MeN}(\text{CH}_2\text{CO}_2\text{H})_2$ (MIDA, **1**),¹⁰ is nontoxic, biodegradable,¹¹ and commercially available. It can also be conveniently, efficiently, and inexpensively synthesized on a large scale¹² from the commodity chemical iminodiacetic acid.¹³

Presently, four different methods for the synthesis of MIDA boronates are known (**Scheme 1**).^{4–9} Many boronic acids can be easily transformed into the corresponding MIDA boronates simply by condensation with MIDA under Dean-Stark conditions (Scheme 1, reaction 1).^{4,6,8,9,12} The removal of water by a variety of alternative techniques (e.g., molecular sieves, azeotropic drying with CH_3CN , etc.) can also promote full conversion to the MIDA boronate product. Typically, this condensation process requires heating to at least 40 °C, and the use of DMSO as a co-solvent is required to partially dissolve the MIDA ligand.

We have also developed several methods that enable the preparation of MIDA boronates without the intermediacy of a boronic acid. Alkenyl MIDA boronates can be synthesized via bromoboration of an alkyne to form the corresponding dibromoborane followed by trapping with MIDA in the presence of 2,6-lutidine (Scheme 1, reaction 2).⁵ Alternatively, a one-pot procedure has been developed in which organotrimethylsilanes can be converted directly into MIDA boronates via transmetalation with BBr_3 , followed by trapping with the disodium salt of MIDA (Na_2MIDA , Scheme 1, reaction 3).⁷ This approach was employed in the efficient synthesis of vinyl MIDA boronate (**9**), for which condensation of MIDA with the related vinylboronic acid or vinylboronate species failed. Lastly, a variety of olefins can be transformed directly into alkenyl MIDA boronates via cross-metathesis with **9** (Scheme 1, reaction 4).⁷ This approach is notable for its generality, efficiency, and mildness. Moreover, in contrast to previous reports involving the use of vinyl or propenyl pinacol boronic esters,^{14,15} cross-metathesis with vinyl MIDA boronate yields only the *E* isomer, and the products are uniformly compatible with silica gel chromatography (vide infra).

3. Physical Properties of MIDA Boronates

MIDA boronates possess a number of highly enabling physical properties that make them useful as a platform for ICC and as convenient alternatives to boronic acids for a wide range of other applications. These properties are remarkably general, i.e., aryl, heteroaryl, alkenyl, and alkyl MIDA boronates all behave similarly. Specifically, MIDA boronates are monomeric, free-flowing, crystalline solids which are stable to storage on the bench top in closed containers under air. MIDA boronates are also universally compatible with silica gel chromatography, allowing convenient reaction monitoring by TLC and facile product isolation and purification.^{4–7} If the goal is to separate different MIDA boronates of similar polarity, a ternary eluent of hexanes,

ethyl acetate, and up to 10% methanol is most effective. We have found that under these conditions even diastereomeric mixtures of MIDA boronates can be resolved. For the purification of nonpolar MIDA boronates, hexanes–ethyl acetate is an effective eluent. Additionally, acetic acid is generally compatible as a co-eluent in most solvent mixtures. Dichloromethane–methanol is a useful eluent for TLC analysis, but some decomposition of the MIDA boronates can occur if this eluent is used for preparative chromatography. Similarly, MIDA boronates should not be left to stand in solutions containing alcohols for more than an hour.

This compatibility with silica gel chromatography and the facility with which MIDA boronates can be formed from the corresponding boronic acids make it possible to utilize MIDA boronate formation and purification as a powerful tool to isolate high-purity boronic acids from crude mixtures that contain many nonboronic acid byproducts. Specifically, we have found that adding a small excess of MIDA to a crude mixture containing a boronic acid and performing a Dean–Stark complexation lead to the formation of the corresponding MIDA boronate while the other impurities remain largely unchanged. If the boronic acid is desired, a simple hydrolysis of the purified MIDA boronate with 1 M aqueous NaOH, followed by extraction of the boronic acid into an organic solvent, typically yields the pure boronic acid.

MIDA boronates are also easily purified by crystallization. A generally effective strategy is to dissolve the crude MIDA boronate in a minimum volume of acetone at 23 °C and then slowly add Et₂O to promote crystallization. The crystallization is complete when the addition of Et₂O no longer clouds the solution. Alternative crystallization solvents include MeCN–Et₂O and EtOAc–Et₂O. X-ray quality crystals are conveniently prepared via vapor diffusion of petroleum ether into an acetone solution of the MIDA boronate: The acetone solution of the MIDA boronate is placed in a small vial (7 mL or 15 mL; the solution is 2–5 mm in height), and this vial is placed in a closed jar containing petroleum ether (about 1 cm in height). Crystals form upon standing at room temperature overnight.

Another important property of MIDA boronates is their solubility in many organic solvents. Reactions are typically performed using THF, dioxane, dichloromethane, DMF, toluene, DMSO, acetonitrile, acetone, or 1,2-dichloroethane. Prolonged exposure of MIDA boronates to aqueous conditions or alcoholic solvents leads eventually to hydrolysis of the MIDA ligand, and this effect is accelerated with heating or in the presence of base. However, water or alcoholic solvents have been successfully employed as co-solvents in some reactions with MIDA boronates.⁹ Furthermore, MIDA boronates are generally compatible with aqueous extractions employing water, brine, aqueous acids (e.g., aq HCl or NH₄Cl), and even some oxidative or reductive aqueous solutions (e.g., aq H₂O₂ at pH < 6, or aq Na₂S₂O₃). Remarkably, even saturated aqueous NaHCO₃ is tolerated in the absence of alcoholic solvents. Aqueous extractions are typically performed using EtOAc or CH₂Cl₂ as the organic phase. For highly polar MIDA boronates, solvent mixtures of EtOAc–acetone (1:1) or THF–Et₂O (1:1) are convenient. As described below, despite this widespread stability, MIDA boronates are easily hydrolyzed to yield the corresponding boronic acids using very mild aqueous basic reagents at 23 °C.

Interestingly, in contrast to MIDA boronates such as **3**, similarly pyramidalized *N*-methyldiethanolamine adducts such as **12** (Figure 2) are not stable to silica gel. As described below, again in contrast to the MIDA boronates, *N*-methyldiethanolamine adducts are also reactive under cross-coupling and many other common reaction conditions.^{4,6,16} The remarkable (and in many cases unique) stability of MIDA boronates to storage under air, chromatography,

aqueous workups, as well as cross-coupling and many other reaction conditions is tentatively attributed to the unique conformational rigidity of the fused bicyclic [*N*-methyliminodiacetate-*O,O',N'*] borane framework. Specifically, as shown in Figure 2, variable-temperature NMR experiments with a DMSO-*d*₆ solution of **3** reveal no coalescence of the diastereotopic methylene protons of the MIDA backbone, even at 150 °C.^{6,17} In contrast, the same experiment with **12** reveals coalescence of the diastereotopic methylene protons of the diethanolamine backbone over a temperature range of 23 to 60 °C, suggesting that this complex is highly dynamic.^{6,17} Although the fundamental underpinnings of these striking differences in conformational rigidity remain to be elucidated, these studies suggest that, uniquely in MIDA boronates, the potentially reactive boron *p* orbital and nitrogen lone pair are kinetically inaccessible, even at elevated temperatures. This kinetic stability may be responsible for many of the unique physical properties of MIDA boronates.

4. Iterative Cross-Coupling (ICC) with Halogenated MIDA Boronates

The now routinely automated process of iterative peptide coupling^{2a} represents an inspiring benchmark for a potentially general strategy for making small molecules in the laboratory. It is interesting to note that peptides are quite complex in structure, having many different functional groups with varied oxidation states and a large number of stereogenic centers. However, the synthesis of many peptides is now very simple, involving the use of a single reaction to iteratively assemble a collection of amino acid building blocks having all of the required functional groups and stereochemistry preinstalled.

With the goal of developing an analogous process for the laboratory construction of small molecules, we decided to focus on the Suzuki–Miyaura reaction and the ICC of bifunctional “haloboronic acids” (see Figure 1).⁶ To avoid random oligomerization of a haloboronic acid under cross-coupling




Figure 2. Variable-Temperature NMR Studies in DMSO-*d*₆ with MIDA Boronate and *N*-Methyldiethanolamine Adducts That Demonstrate the Unique and Remarkable Conformational Rigidity of the MIDA Boronate Framework. (Ref. 6)

conditions, it is necessary to reversibly attenuate the reactivity of one end of this type of bifunctional reagent, in analogy to the use of a protective group to control the reactivity of the amine terminus of an amino acid.¹⁸ Toward this goal, we chose to focus on controlling the reactivity of the boronic acid functional group.

It is hypothesized that a vacant and Lewis acidic boron p orbital is required for transmetalation of a boronic acid under Suzuki–Miyaura cross-coupling conditions (**Figure 3a**).¹⁹ Consistent with this, complexation with electron-donating, Lewis basic ligands is known to attenuate the reactivity of boronic acids towards cross-coupling (**Figure 3b**).^{19a} For example, pinacol boronic esters can be less reactive towards cross-coupling than the corresponding boronic acids.²⁰ This reactivity attenuation can be attributed to the decreased Lewis acidity of the boron p orbital as a result of

conjugation with the lone pairs of the ligand heteroatoms.^{19a} This same approach has been utilized with a variety of other divalent heteroatomic ligands.²¹ There is an inherent limitation, however, that precludes the general utilization of this approach for complex small-molecule synthesis. Specifically, conjugation between the heteroatom lone pairs and the boron p orbital produces relatively strong boron–heteroatom bonds, creating a high energy barrier for bond cleavage. Moreover, the equilibrium between the boronic acid and the corresponding boronic ester typically lies strongly towards the latter, thereby disfavoring hydrolysis. As a result, cleaving this type of ligand to regenerate the boronic acid typically requires harsh conditions^{19–21} and/or additional reagents to destroy the divalent ligand after it has been cleaved.²² These types of conditions can be problematic in the context of complex small-molecule synthesis.

Recognizing the inherent limitations of this approach, we focused on an alternative strategy (**Figure 3c**).⁴ Specifically, given that the boron p orbital is predicted to be critical for the transmetalation of a boronic acid, we hypothesized that removing this p orbital through rehybridization of the boron atom from sp² to sp³ via complexation with a trivalent heteroatomic ligand would eliminate its reactivity towards cross-coupling. Further increasing our interest in this approach, it is known that boron–heteroatom bonds in tetrahedral adducts are weaker than those in their tricoordinate counterparts.²³ For example, the pyramidalization of trimethyl borate via complexation with ammonia weakens the boron–oxygen bonds by about 10–12 kcal/mol.²⁴ Thus, we felt that it might be possible to find relatively mild conditions that could hydrolyze this type of pyramidalized boronate and regenerate the reactive boronic acid. After surveying a series of trivalent heteroatomic ligands, we discovered that MIDA boronates embody all of these expectations and represent a powerful platform for ICC chemistry.

In a competition experiment between *p*-(*n*-butyl)phenylboronic acid (**13**) and *p*-tolyl MIDA boronate (**14a**) under Buchwald-type²⁵ anhydrous Suzuki–Miyaura cross-coupling conditions with *p*-bromoanisaldehyde, we observed a 24:1 ratio of products **15** and **16** (**eq 1**), consistent with a strong preference for cross-coupling of the sp²-hybridized boronic acid.⁴ Interestingly, a wide range of non-aryl substituents were tolerated on the nitrogen atom. The diethanolamine adduct, **14d**, lacking the carbonyl units of MIDA, was as reactive as boronic acid **13**. As described above, this difference in reactivity between **14a** and **14d** is attributed to differences in the conformational rigidity of these two complexes.

Encouraged by these results, we set out to prepare a series of bifunctional B-protected haloboronic acids and explore their capacity to undergo selective cross-coupling at the halide terminus. The efficient synthesis of aryl, heteroaryl, alkenyl, and alkyl derivatives was achieved via simple condensation of the corresponding boronic acids with MIDA under Dean–Stark conditions (**eq 2**).⁴

As shown in **Scheme 2**,⁴ this B-protection strategy is remarkably general, with the same ligand similarly protecting aryl, heteroaryl, alkenyl, and alkyl haloboronic acids, thereby enabling the highly selective coupling of the halide terminus of building blocks **18a–f**. Moreover, consistent with our initial hypothesis, the MIDA boronate products **19a–f** can all be hydrolyzed under mild aqueous basic conditions (1 N NaOH(aq), THF, 23 °C, 10 min) to generate the corresponding free boronic acids **20a–f**.

Polyenes are especially challenging synthetic targets because of the sensitivity of this framework to light, oxygen, and acid. It is also critical to control the stereochemistry of each double bond. The ICC approach is particularly well-suited to preparing these types of compound due to the mild and stereospecific nature of the





Figure 3. Under Suzuki–Miyaura Cross-Coupling Reaction Conditions: (a) Transmetalation Requires a Vacant and Lewis Acidic Boron p Orbital. (b) Strongly Electron-Donating Divalent Ligands Can Attenuate Boronic Acid Reactivity, but Typically Require Relatively Harsh Conditions for Cleavage. (c) The Reactivity of a Boronic Acid Can Be Reversibly Attenuated via Pyramidalization with a Trivalent Heteroatomic Ligand. (Ref. 4)



^a By HPLC; average of three runs.

eq 1 (Ref. 4)


metal-mediated cross-couplings. Given the prevalence of alkenyl and polyenyl subunits in both natural products and pharmaceutical targets, we developed a collection of bifunctional building blocks specifically designed to enable polyene synthesis via ICC.⁵

As described in Scheme 1, *trans*-(2-bromovinyl) MIDA boronate (**7**) can be prepared via bromoboration of acetylene²⁶ followed by complexation with MIDA in the presence of 2,6-lutidine. An alternative and more convenient procedure involves transmetalation of 1-bromo-2-(trimethylsilyl)ethylene with BBr_3 ²⁷ followed by trapping with Na_2MIDA .⁷ Bifunctional olefin **7** is a remarkably versatile cross-coupling partner (Scheme 3).⁵ Specifically, Suzuki–Miyaura, Stille, and Heck couplings are all achieved at the bromide terminus without perturbing the MIDA boronate. A series of bismetallated lynchpin-type reagents are also created via Sonogashira coupling with trimethylsilylacetylene, Miyaura borylation with pinacolodiborane (**25**), or a triply metal-selective (Zn vs Sn and B) Negishi coupling with bismetallated olefin **27**.

A generally useful strategy involves the boron-selective coupling of differentially ligated diboron reagents.^{5,28} In the first example of such a reaction, **26** was selectively coupled with *trans*-1-chloro-2-iodoethylene at the sp^2 -hybridized pinacol boronic ester terminus to generate chlorodiene MIDA boronate **29** (Scheme 4).⁵ A better-precedented Sn vs B coupling²⁹ between bismetallated diene **28** and *trans*-1-chloro-2-iodoethylene generated chlorotrienyl MIDA boronate **30**.⁵


The olefin cross-metathesis route to MIDA boronates is remarkably tolerant of a wide range of functional groups, including halogens. Thus, this method is also well-suited for preparing various haloalkenyl MIDA boronates, as shown for a series of bromostyrene derivatives (eq 3).⁷

The capacity to prepare and selectively couple bifunctional halo MIDA boronates enables one to envision the synthesis of natural products or pharmaceuticals by using only a single reaction iteratively to bring together a collection of pre-assembled building blocks. This strategy was first realized with the total synthesis of ratanahine,⁴ a complex neolignan isolated from the *Ratanhia radix* by Arnone and co-workers in 1990.³⁰ This natural product was retrosynthetically fragmented using recursive Suzuki–Miyaura transformations to generate four simpler building blocks, **35–38** (Figure 4). There were several challenges associated with this plan that were expected to test the limits of the MIDA-based




^a Yields refer to *analytically pure* materials obtained conveniently via silica gel chromatography

eq 2 (Ref. 4)



¹ B-Deprotection was also achieved by treatment with saturated aq. NaHCO₃–MeOH, 23 °C, 6 h (85%) ² X-Phos was used as the phosphine ligand

Scheme 2. Selective Coupling and Mild Deprotection of Halogenated Aryl, Heteroaryl, Alkenyl, and Alkyl MIDA Boronates. (Ref. 4)




Scheme 3. Bifunctional Halogenated MIDA Boronate Building Blocks Such as **7** Can Be Rapidly and Selectively Elaborated at the Halogen Terminus. (Ref. 5)

ICC methodology. First, couplings of alkenylboronic acids tend to be less efficient than those of their aryl counterparts, making the selective coupling between **35** and aryl MIDA boronate **36** unsecured. In addition, 2-substituted heterocyclic boronic acids such as **36** are notoriously unstable and difficult to purify, store, and cross-couple.³¹ Finally, the coupling of highly deactivated bromoaryl MIDA boronate **37** was expected to demand more forcing reaction conditions that would test the limits of stability of the MIDA boronate functionality.


Despite these challenges, the total synthesis of ratanhine was achieved via ICC as shown in **Scheme 5**.⁴ Specifically, selective coupling between propenylboronic acid (**35**) and 5-bromobenzofuran-2-yl MIDA boronate (**36**) proceeded smoothly to generate substituted benzofuran MIDA boronate **39**. Remarkably, while the corresponding benzofuranboronic acid decomposed over the course of several days, MIDA boronate **39** was stored on the bench top under air without noticeable decomposition for more than 6 months. This MIDA boronate was hydrolyzed under mild conditions, and the resulting boronic

acid was immediately utilized in a cross-coupling reaction with bromoaryl MIDA boronate **37**. As expected, this coupling required increased temperature (80 °C in a sealed tube) and an extended reaction time (28 h). Remarkably, the MIDA boronate functional group was stable to these forcing conditions, yielding the highly conjugated MIDA boronate product **40**. A final sequence of boronic acid deprotection and coupling with alkenyl bromide **38** and MOM-ether deprotection completed the first total synthesis of ratanhine. More importantly, to the best of our knowledge, this represents the first total synthesis of any natural product in which a single reaction was utilized iteratively to assemble all of the required building blocks.

This ICC strategy is also highly effective in the synthesis of polyene natural products.⁵ Specifically, *all-trans*-retinal³² was prepared simply via ICC of boronic acid **42**, bromoalkenyl MIDA boronate **7**, and alkenyl bromide **44** (**Scheme 6**).⁵ In a similar vein, β -parinaric acid³³ was prepared by ICC of butenylboronic acid (**46**), chlorodiethyl MIDA boronate **29**, and alkenyl iodide **48** (**Scheme 7**).⁵ Finally, despite the fact that



Scheme 4. The Preparation of Bifunctional Halogenated MIDA Boronates for Use in Polyene Synthesis. (Ref. 5)



eq 3 (Ref. 7)






Figure 4. Retrosynthetic Analysis of the Natural Product Ratanhine Using ICC Identified Four Building Blocks, of Which **36** and **37** Are Bifunctional Halogenated MIDA Boronates. (Ref. 4)



Scheme 5. Application of ICC in the First Total Synthesis of Ratanhine. (Ref. 4)



Scheme 6. Application of ICC in the Synthesis of *All-trans*-Retinal. (Ref. 5)

polyenylboronic acids can be very unstable,³⁴ the notoriously challenging heptaene framework of the polyene natural product amphotericin B was prepared using only the Suzuki–Miyaura reaction to assemble a collection of bifunctional haloalkenyl MIDA boronates (**Scheme 8**).⁵

As demonstrated by these examples, the ICC approach has significant potential to enable the simple, efficient, and flexible construction of small molecules.

5. Multistep Synthesis of Complex Boronic Acids from Simple MIDA Boronates


To avoid a general incompatibility with synthetic reagents, it is typically necessary to introduce the boronic acid functional group just prior to its utilization in a cross-coupling or other type of reaction. However, most of the methods that are available for achieving this have poor functional-group tolerance. Collectively, these limitations can render the synthesis of complex boronic acids very challenging. This can sometimes preclude the use of boronic acids in complex-molecule synthesis, and represents a potential bottleneck for the development of a truly general ICC-based approach.

Some sterically bulky boronic esters are known to be more tolerant of synthetic reagents;³⁵ however, removing these ligands to generate a targeted boronic acid usually requires harsh conditions that are generally incompatible with sensitive building blocks. Trifluoroborate salts represent very useful surrogates for boronic acids,³⁶ and Molander and co-workers have powerfully demonstrated that the trifluoroborate functional group is compatible with many synthetic reagents.³⁷ These features have provided novel access to many new organoborane building blocks. However, the incompatibility of trifluoroborate salts with chromatography can limit the utilization of these reagents in multistep synthesis, which is often necessary for accessing structurally and/or stereochemically complex building blocks.


Overcoming these limitations, we have recently found that the MIDA boronate functional group is stable to a wide range of common synthetic reagents, presumably due to the lack of a reactive boron *p* orbital.⁹ Combined with the general compatibility of MIDA boronates with chromatography and the capacity to release the corresponding boronic acids under very mild conditions, this stability enables the first reliable approach for the multistep synthesis of complex boronic acids from simple organoborane starting materials.

Specifically, *p*-(hydroxymethyl)phenyl MIDA boronate (**3**) can be smoothly oxidized under Swern conditions to generate the corresponding benzaldehyde (**Scheme 9**).⁶ Remarkably, this MIDA boronate is also stable to the very strongly acidic and oxidizing Jones conditions (H₂SO₄–CrO₃). This latter stability is highly unique; i.e., under these same conditions, the corresponding boronic acid (**56a**), pinacolboronic ester (**56b**), 1,8-diaminonaphthalene adduct (**56c**), trifluoroborate salt (**56d**), and *N*-methyldiethanolamine boronate^{16b,38} (**56e**) all decomposed. Similar to that which we observed under cross-coupling conditions, the remarkable difference in reactivity between the MIDA and diethanolamine boronates is likely related to the differences in conformational flexibility of the two complexes (see Figure 2).


This unique compatibility with strong acid and oxidants suggested that MIDA boronates could be stable to a wide range of reaction conditions. In fact, even triflic acid (p*K*_a–14) was tolerated, enabling the *p*-methoxybenzylation of **3** and the reversal of this transformation with DDQ (**Scheme 10**).⁹



Scheme 7. Application of ICC in the Synthesis of β -Parinaric Acid. (Ref. 5)



Scheme 8. Application of ICC in the Synthesis of the Heptaene Framework of Amphotericin B. (Ref. 5)




Scheme 9. MIDA Boronates, Such as **3**, Are Uniquely Stable to the Strongly Acidic and Oxidizing Jones Conditions, Whereas Boronic Acids and Boronates, Such as **56a–e**, Are Not. (Ref. 6)


Similarly, silylation and desilylation were well-tolerated, as was the transformation into the corresponding benzyl iodide **59** with $\text{PPh}_3\text{-I}_2$.

This latter reaction suggested compatibility with soft nucleophiles. In this vein, benzaldehyde **54** was successfully utilized in a series of carbon–carbon-bond-forming reactions including the Evans aldol, Horner–Wadsworth–Emmons, and Takai olefination reactions. Reductive amination and aldehyde reduction were also well-tolerated (Scheme 11).⁶

Whereas cyclopropanation of vinyl boranes^{39,40} and epoxidation of 1,2-disubstituted alkenyltrifluoroborate salts^{37b} are known, the versatility and broad compatibility of vinyl MIDA boronate (**9**) as a starting material has also been demonstrated (Scheme 12).⁷ Specifically, cyclopropanation of **9** produced cyclopropyl MIDA boronate (**64**) in excellent yield. Remarkably, epoxidation of this olefin with *m*CPBA was also well-tolerated, and even this epoxide, **65**, was stable to column chromatography and bench top storage under air. Boronate **9** was successfully engaged in the Heck reaction⁴¹



Scheme 10. MIDA Boronates Are Stable to a Number of Useful Reagents in Organic Synthesis. (Ref. 6)




Scheme 11. MIDA Boronate Building Blocks Can Be Elaborated Through a Number of Common Organic Reactions. (Ref. 6)


to yield styrenyl derivative **66**. Similarly, the White catalyst⁴² promoted an efficient oxidative Heck-type⁴³ reaction⁴⁴ to yield **67**. As described previously (see eq 3), **9** is also an excellent substrate for olefin cross-metathesis^{44,45} (analogous to *tert*-butylethylene), yielding (*E*)-octenyl MIDA boronate (**68**) as a chromatographically and air-stable crystalline solid and a single stereoisomer (see Scheme 12). Fortunately, this approach has proven to be quite general, and represents a very useful method for preparing a range of (*E*)-alkenyl MIDA boronates (eq 4).⁷

This broad compatibility of the MIDA boronate functional group with a wide range of reagents can enable the transformation of simple MIDA boronates into otherwise difficult-to-access complex boronic acids for use in a variety of synthesis applications. These include structurally complex B-protected haloboric acids for use in ICC.

To explore the enabling potential of this approach, we targeted the total synthesis of the natural product crocacin C.^{6,45} As shown in Figure 5, this molecule was retrosynthetically



Scheme 12. Vinyl MIDA Boronate (**9**) Is an Exceptionally Versatile Building Block. (Ref. 7)



eq 4 (Ref. 7)

fragmented via recursive cross-coupling transforms into known building blocks **72** and **74** as well as the novel, complex iodoalkenyl MIDA boronate **73**. The preparation of the latter represented a significant challenge that we hypothesized could be overcome via multistep synthesis starting with simple MIDA boronate **75** (Scheme 13).⁶

In practice, a Paterson aldol reaction between **75** and **76** followed by diastereoselective reduction of the resulting β -hydroxy ketone yielded diol **77**. Importantly, the small amounts of diastereomeric byproducts that are typically generated in these types of transformations were readily removed by taking advantage of the compatibility of the MIDA boronate functional group with silica gel chromatography. A subsequent sequence of permethylation with Meerwein's salt, oxidative cleavage of the PMB ether, oxidation of the resulting primary alcohol, and Takai olefination yielded the targeted, complex halogenated MIDA boronate **73**. Importantly, **73**, **75**, and all intermediates were compatible with chromatography and storage on the bench top under air. With B-protected haloboronic acid **73** in hand, the synthesis of (+)-crocacin C was readily achieved via ICC.⁶

6. Conclusions and Prospects

As described herein, the inherent modularity found in many of the small molecules targeted for synthesis in the laboratory stands to be more effectively harnessed via the ICC approach. Analogous to the synthesis of peptides, oligonucleotides, and oligosaccharides, this strategy has the potential to enable the preparation of a wide range of small molecules by the simple, iterative union of pre-assembled, bifunctional building blocks. Due to their ease of synthesis, purification, characterization, and storage; their capacity for reversibly attenuated reactivity under cross-coupling conditions; and their compatibility with a wide range of common synthetic reagents; MIDA boronates represent a powerful platform for the development of this type of synthesis strategy. Moreover, it was recently discovered that, under novel "slow-release cross-coupling" conditions, MIDA boronates can serve as highly effective surrogates for even notoriously unstable boronic acids, such as 2-heterocyclic (including 2-pyridyl), vinyl, and cyclopropyl derivatives.⁴⁶ This remarkably general approach has transformed a wide range of unstable boronic acids into air-stable and highly effective cross-coupling partners.

Looking forward, the ever-expanding scope of the Suzuki–Miyaura coupling suggests that the potential generality of this ICC approach could be substantial. Particularly critical to realizing this potential will be finding a way to form $\text{Csp}^3\text{–Csp}^2$

and even $\text{Csp}^3\text{–Csp}^1$ bonds with the same efficiency that is now routinely achieved stereospecifically with $\text{Csp}^2\text{–Csp}^3$ linkages. The discovery of additional methods to prepare MIDA boronates that do not proceed through the intermediacy of a difficult-to-access and/or unstable boronic acid will also be vital. Moreover, to realize the ultimate goal of developing a machine with the capacity for fully automated ICC, it will be important to further develop cross-coupling conditions that are maximally general²⁵ (to avoid the requirement for ad hoc optimization of conditions for each combination of coupling partners) and amenable to translation into the solid phase or some other form of iterative-synthesis-enabling technology. While these challenges are admittedly considerable, we are convinced that they each can be solved. Achieving these goals could have a substantial impact on the synthesis of small molecules in the laboratory, and may ultimately even extend the power of this discovery engine to the nonchemist.

7. Acknowledgements

Different aspects of this research were supported by NIH (GM 080436), NSF (CAREER 0747778), Bristol-Myers Squibb, Sigma-Aldrich, and the University of Illinois. EPG is a Seemon Pines Graduate Fellow. MDB is a Dreyfus New Faculty Awardee, a Beckman and Amgen Young Investigator, and a Sloan Research Fellow.

8. References and Notes

- (1) (a) Fischbach, M. A.; Walsh, C. T. *Chem. Rev.* **2006**, *106*, 3468. (b) Garret, R. H.; Grisham, C. M. *Biochemistry*; Saunders College Publishing: Fort Worth, TX, 1995.





Figure 5. Retrosynthetic Analysis of the Natural Product (+)-Crocacin C Identified the Relatively Complex Halogenated MIDA Boronate Building Block **73**. (Ref 6)



Scheme 13. The Broad Chemical Stability of the MIDA Boronate Functional Group Is Taken Advantage of in the Synthesis of (+)-Crocacin C. (Ref 6)

- (2) (a) Merrifield, R. B. *Angew Chem., Int. Ed. Engl.* **1985**, *24*, 799. (b) Caruthers, M. H. *Science* **1985**, *230*, 281. (c) Plante, O. J.; Palmacci, E. R.; Seeberger, P. H. *Science* **2001**, *291*, 1523.
- (3) (a) Zhang, J.; Moore, J. S.; Xu, Z.; Aguirre, R. A. *J. Am. Chem. Soc.* **1992**, *114*, 2273. (b) Young, J. K.; Nelson, J. C.; Moore, J. S. *J. Am. Chem. Soc.* **1994**, *116*, 10841. (c) Feuerbacher, N.; Vögtle, F. In *Dendrimers*; Vögtle, F., Ed.; Topics in Current Chemistry Series; Springer: Berlin, 1998; Vol. 197, pp 1–18.
- (4) Gillis, E. P.; Burke, M. D. *J. Am. Chem. Soc.* **2007**, *129*, 6716.
- (5) Lee, S. J.; Gray, K. C.; Paek, J. S.; Burke, M. D. *J. Am. Chem. Soc.* **2008**, *130*, 466.
- (6) Gillis, E. P.; Burke, M. D. *J. Am. Chem. Soc.* **2008**, *130*, 14084.
- (7) Uno, B. E.; Gillis, E. P.; Burke, M. D. *Tetrahedron* **2009**, *65*, 3130.
- (8) Mancilla, T.; Contreras, R.; Wrackmeyer, B. *J. Organomet. Chem.* **1986**, *307*, 1.
- (9) Mancilla, T.; Zamudio-Rivera, L. S.; Beltrán, H. I.; Santillan, R.; Farfán, N. *Arkivoc* **2005** (vi), 366.
- (10) (a) Berchet, G. *J. Org. Synth.* **1938**, *18*, 56. (b) Childs, A. F.; Goldsworthy, L. J.; Harding, G. F.; King, F. E.; Nineham, A. W.; Norris, W. L.; Plant, S. G. P.; Selton, B.; Tompsett, A. L. L. *J. Chem. Soc.* **1948**, 2174. (c) Chase, B. H.; Downes, A. M. *J. Chem. Soc.* **1953**, 3874.
- (11) Warren, C. B.; Malec, E. J. *Science* **1972**, *176*, 277.
- (12) Ballmer, S. G.; Gillis, E. P.; Burke, M. D. *Org. Synth.*, submitted for publication, 2008.
- (13) Tens of thousands of metric tons of iminodiacetic acid are produced annually worldwide: Yangong, Z. *China Chemical Reporter* **2005**, *16* (8), 16.
- (14) (a) Blackwell, H. E.; O'Leary, D. J.; Chatterjee, A. K.; Washenfelder, R. A.; Bussmann, D. A.; Grubbs, R. H. *J. Am. Chem. Soc.* **2000**, *122*, 58. (b) Morrill, C.; Grubbs, R. H. *J. Org. Chem.* **2003**, *68*, 6031. (c) Njardarson, J. T.; Biswas, K.; Danishefsky, S. J. *Chem. Commun.* **2002**, 2759. (d) Nicolaou, K. C.; Li, A.; Edmonds, D. J. *Angew. Chem., Int. Ed.* **2006**, *45*, 7086. (e) Jankowska, M.; Pietraszuk, C.; Marciniak, B.; Zaidlewicz, M. *Synlett* **2006**, 1695. (f) Dakas, P.-Y.; Barluenga, S.; Totzke, F.; Zirrgebel, U.; Winssinger, N. *Angew. Chem., Int. Ed.* **2007**, *46*, 6899. (g) Wu, Y.; Gao, J. *Org. Lett.* **2008**, *10*, 1533.
- (15) Ring-closing-metathesis reactions: (a) Ashe, A. J., III; Fang, X.; Kampf, J. W. *Organometallics* **2000**, *19*, 4935. (b) Ashe, A. J., III; Yang, H.; Fang, X.; Kampf, J. W. *Organometallics* **2002**, *21*, 4578.
- (16) For some important exceptions, see: (a) Ebdrup, S.; Jacobsen, P.; Farrington, A. D.; Vedso, P. *Bioorg. Med. Chem.* **2005**, *13*, 2305. (b) Vedso, P.; Olesen, P. H.; Hoeg-Jensen, T. *Synlett* **2004**, 892. (c) Kirihata, M.; Morimoto, T.; Mizuta, T.; Ichimoto, I. *Biosci. Biotechnol. Biochem.* **1995**, *59*, 2317. (d) Kirihara, M.; Morimoto, T.; Ichimoto, I. *Biosci. Biotechnol. Biochem.* **1993**, *57*, 1940. (e) Yamamoto, Y.; Seko, T.; Rong, F.; Nemoto, H. *Tetrahedron Lett.* **1989**, *30*, 7191.
- (17) These results are consistent with previous studies: (a) References 8 and 9. (b) Psarras, T. G.; Zimmerman, H. K., Jr.; Rasiel, Y.; Weidmann, H. *Liebigs Ann. Chem.* **1962**, *655*, 48. (c) Contreras, R.; Garcia, C.; Mancilla, T.; Wrackmeyer, B. *J. Organomet. Chem.* **1983**, *246*, 213.
- (18) (a) Bergmann, M.; Zervas, L. *Ber. Dtsch. Chem. Ges.* **1932**, *65*, 1192. (b) Carpino, L. A.; Han, G. Y. *J. Am. Chem. Soc.* **1970**, *92*, 5748. (c) Tarbell, D. S.; Yamamoto, Y.; Pope, B. M. *Proc. Natl. Acad. Sci. U.S.A.* **1972**, *69*, 730.
- (19) (a) Matos, K.; Soderquist, J. A. *J. Org. Chem.* **1998**, *63*, 461. (b) Miyaura, N. *J. Organomet. Chem.* **2002**, *653*, 54. (c) Braga, A. A. C.; Morgon, N. H.; Ujaque, G.; Maseras, F. J. *J. Am. Chem. Soc.* **2005**, *127*, 9298.
- (20) Holmes, D.; Chotana, G. A.; Maleczka, R. E., Jr.; Smith, M. R., III. *Org. Lett.* **2006**, *8*, 1407.
- (21) (a) Deng, X.; Mayeux, A.; Cai, C. *J. Org. Chem.* **2002**, *67*, 5279. (b) Hohn, E.; Pietruszka, J. *Adv. Synth. Catal.* **2004**, *346*, 863. (c) Noguchi, H.; Hojo, K.; Suginome, M. *J. Am. Chem. Soc.* **2007**, *129*, 758.
- (22) (a) Nakamura, H.; Fujiwara, M.; Yamamoto, Y. *Bull. Chem. Soc. Jpn.* **2000**, *73*, 231. (b) Falck, J. R.; Bondlela, M.; Venkataraman, S. K.; Srinivas, D. *J. Org. Chem.* **2001**, *66*, 7148. (c) Murphy, J. M.; Tzschucke, C. C.; Hartwig, J. F. *Org. Lett.* **2007**, *9*, 757.
- (23) Hall, D. G. In *Boronic Acids*; Hall, D. G., Ed.; Wiley-VCH: Weinheim, 2005; Chapter 1, pp 1–99.
- (24) Matteson, D. S. *Stereodirected Synthesis with Organoboranes*; Reactivity and Structure Series; Rees, C. W., Lehn, J. M., Zahradnik, R., Hafner, K., Series Eds.; Springer-Verlag: New York, 1995; Vol. 32, pp 1–20.
- (25) (a) Barder, T. E.; Walker, S. D.; Martinelli, J. R.; Buchwald, S. L. *J. Am. Chem. Soc.* **2005**, *127*, 4685. (b) Billingsley, K.; Buchwald, S. L. *J. Am. Chem. Soc.* **2007**, *129*, 3358. (c) Martin, R.; Buchwald, S. L. *Acc. Chem. Res.* **2008**, *41*, 1461.
- (26) Hyuga, S.; Chiba, Y.; Yamashina, N.; Hara, S.; Suzuki, A. *Chem. Lett.* **1987**, 1757.
- (27) Singleton, D. A.; Leung, S.-W. *J. Organomet. Chem.* **1997**, *544*, 157.
- (28) For some recent additional examples, see: (a) Molander, G. A.; Sandrock, D. L. *J. Am. Chem. Soc.* **2008**, *130*, 15792. (b) Noguchi, H.; Shioda, T.; Chou, C.-M.; Suginome, M. *Org. Lett.* **2008**, *10*, 377.
- (29) (a) Coleman, R. S.; Walczak, M. C. *Org. Lett.* **2005**, *7*, 2289. (b) Coleman, R. S.; Lu, X.; Modolo, I. *J. Am. Chem. Soc.* **2007**, *129*, 3826.
- (30) Arnone, A.; Modugno, V. D.; Nasini, G.; de Pava, O. V. *Gazz. Chim. Ital.* **1990**, *120*, 397.
- (31) Tyrrell, E.; Brookes, P. *Synthesis* **2004**, 469.
- (32) Previous cross-coupling-based syntheses of retinoids: (a) Negishi, E.; Owczarczyk, Z. *Tetrahedron Lett.* **1991**, *32*, 6683. (b) Uenishi, J.; Kawahama, R.; Yonemitsu, O.; Wada, A.; Ito, M. *Angew. Chem., Int. Ed.* **1998**, *37*, 320. (c) Uenishi, J.; Matsui, K.; Wada, A. *Tetrahedron Lett.* **2003**, *44*, 3093.
- (33) Some prior syntheses of β -parinaric acid: (a) Kuklev, D. V.; Smith, W. L. *Chem. Phys. Lipids* **2004**, *131*, 215. (b) Goerger, M. M.; Hudson, B. S. *J. Org. Chem.* **1988**, *53*, 3148. (c) Hayashi, T.; Oishi, T. *Chem. Lett.* **1985**, 413.
- (34) (a) Roush, W. R.; Brown, B. B. *J. Am. Chem. Soc.* **1993**, *115*, 2268. (b) Torrado, A.; Iglesias, B.; López, S.; de Lera, A. R. *Tetrahedron* **1995**, *51*, 2435.
- (35) (a) Matteson, D. S.; Kandil, A. A.; Soundararajan, R. *J. Am. Chem. Soc.* **1990**, *112*, 3964. (b) Matteson, D. S.; Man, H.-W.; Ho, O. C. *J. Am. Chem. Soc.* **1996**, *118*, 4560. (c) Luthle, J. E. A.; Pietruszka, J. *J. Org. Chem.* **2000**, *65*, 9194. (d) Hilt, G.; Bolze, P. *Synthesis* **2005**, 2091. (e) Jin, B.; Liu, Q.; Sulikowski, G. A. *Tetrahedron* **2005**, *61*, 401. (f) Gopalathnam, A.; Nelson, S. G. *Org. Lett.* **2006**, *8*, 7.
- (36) (a) Molander, G. A.; Ellis, N. *Acc. Chem. Res.* **2007**, *40*, 275. (b) Darses, S.; Genet, J.-P. *Chem. Rev.* **2008**, *108*, 288.
- (37) (a) Molander, G. A.; Petrillo, D. E. *J. Am. Chem. Soc.* **2006**, *128*, 9634. (b) Molander, G. A.; Ribagorda, M. *J. Am. Chem. Soc.* **2003**, *125*, 11148. (c) Molander, G. A.; Petrillo, D. E. *Org. Lett.* **2008**, *10*, 1795. (d) Molander, G. A.; Gormisky, P. E.; Sandrock, D. L. *J. Org. Chem.* **2008**, *73*, 2052. (e) Pucheault, M.; Darses, S.; Genet, J.-P. *J. Am. Chem. Soc.* **2004**, *126*, 15356.
- (38) Gravel, M.; Thompson, K. A.; Zak, M.; Bérubé, C.; Hall, D. G. *J. Org. Chem.* **2002**, *67*, 3.
- (39) Fontani, P.; Carboni, B.; Vaultier, M.; Carrié, R. *Tetrahedron Lett.* **1989**, *30*, 4815.
- (40) Fontani, P.; Carboni, B.; Vaultier, M.; Maas, G. *Synthesis* **1991**, 605.
- (41) Heck reactions with vinyl boranes: (a) Hunt, A. R.; Stewart, S. K.; Whiting, A. *Tetrahedron Lett.* **1993**, *34*, 3599. (b) Stewart, S. K.;


- Whiting, A. *J. Organomet. Chem.* **1994**, 482, 293. (c) Lightfoot, A. P.; Maw, G.; Thirsk, C.; Twiddle, S. J. R.; Whiting, A. *Tetrahedron Lett.* **2003**, 44, 7645. (d) Stewart, S. K.; Whiting, A. *Tetrahedron Lett.* **1995**, 36, 3925. (e) Hénaff, N.; Whiting, A. *Org. Lett.* **1999**, 1, 1137. (f) Hénaff, N.; Whiting, A. *Tetrahedron* **2000**, 56, 5193. (g) Lightfoot, A. P.; Twiddle, S. J. R.; Whiting, A. *Org. Biomol. Chem.* **2005**, 3, 3167. (h) Batsanov, A. S.; Knowles, J. P.; Whiting, A. *J. Org. Chem.* **2007**, 72, 2525. (i) Itami, K.; Tonogaki, K.; Nokami, T.; Ohashi, Y.; Yoshida, J. *Angew. Chem., Int. Ed.* **2006**, 45, 2404. (j) Itami, K.; Tonogaki, K.; Ohashi, Y.; Yoshida, J. *Org. Lett.* **2004**, 6, 4093.
- (42) (a) Chen, M. S.; White, M. C. *J. Am. Chem. Soc.* **2004**, 126, 1346. (b) Chen, M. S.; Prabakaran, N.; Labenz, N. A.; White, M. C. *J. Am. Chem. Soc.* **2005**, 127, 6970. (c) Fraunhoffer, K. J.; Prabakaran, N.; Sirois, L. E.; White, M. C. *J. Am. Chem. Soc.* **2006**, 128, 9032. (d) Fraunhoffer, K. J.; White, M. C. *J. Am. Chem. Soc.* **2007**, 129, 7274. (e) Reed, S. A.; White, M. C. *J. Am. Chem. Soc.* **2008**, 130, 3316. (f) Young, A. J.; White, M. C. *J. Am. Chem. Soc.* **2008**, 130, 14090.
- (43) (a) Cho, C. S.; Uemura, S. *J. Organomet. Chem.* **1994**, 465, 85. (b) Du, X.; Suguro, M.; Hirabayashi, K.; Mori, A.; Nishikata, T.; Hagiwara, N.; Kawata, K.; Okeda, T.; Wang, H. F.; Fugami, K.; Kosugi, M. *Org. Lett.* **2001**, 3, 3313. (c) Yoo, K. S.; Yoon, C. H.; Jung, K. W. *J. Am. Chem. Soc.* **2006**, 128, 16384. (d) Lindh, J.; Enquist, P.-A.; Pilotti, Å.; Nilsson, P.; Larhed, M. *J. Org. Chem.* **2007**, 72, 7957.
- (44) (a) Delcamp, J. H.; White, M. C. *J. Am. Chem. Soc.* **2006**, 128, 15076. (b) Delcamp, J. H.; Brucks, A. P.; White, M. C. *J. Am. Chem. Soc.* **2008**, 130, 11270.
- (45) (a) Kunze, B.; Jansen, R.; Höfle, G.; Reichenbach, H. *J. Antibiot* **1994**, 47, 881. Previous syntheses: (b) Dias, L. C.; de Oliveira, L. G. *Org. Lett.* **2001**, 3, 3951. (c) Chakraborty, T. K.; Jayaprakash, S. *Tetrahedron Lett.* **2001**, 42, 497. (d) Feutrill, J. T.; Lilly, M. J.; Rizzacasa, M. A. *Org. Lett.* **2000**, 2, 3365. (e) Sirasani, G.; Paul, T.; Andrade, R. B. *J. Org. Chem.* **2008**, 73, 6386.
- (46) (a) Knapp, D. M.; Gillis, E. P.; Burke, M. D. *J. Am. Chem. Soc.* [Online early access]. DOI: 10.1021/ja901416p. Published Online:

Apr 30, 2009. <http://pubs.acs.org/doi/abs/10.1021/ja901416p> (accessed May 14, 2009). (b) Borman, S. *Chem. Eng. News* **2009**, 87 (19), 8.

Keywords: MIDA boronate; iterative cross-coupling; automated synthesis; *N*-methyliminodiacetic acid; bifunctional molecules.

About the Authors

Eric Gillis grew up in Portland, ME, and received his undergraduate education at Grinnell College in Grinnell, IA. While at Grinnell, he worked with Professor T. Andrew Mobley on the synthesis and characterization of tungsten–tin complexes. In 2005, he began his doctoral studies under the direction of Professor Martin Burke at the University of Illinois at Urbana-Champaign. His current research focuses on the development of MIDA boronate esters as a platform to enable the automated synthesis of small molecules.

Marty Burke grew up in Manchester, MD, and graduated in 1998 with a B.A. degree in chemistry from Johns Hopkins University. He then moved to Harvard Medical School as a Ph.D. and M.D. student in the Health Sciences and Technology program. From 1999 to 2003, Marty completed his thesis research in organic synthesis under the direction of Professor Stuart Schreiber, and graduated from medical school in 2005. That same year, Marty began his independent career as an assistant professor in the Department of Chemistry at the University of Illinois at Urbana-Champaign. His research focuses on the synthesis and study of small molecules that perform protein-like functions. To enable these studies, Marty's group is pioneering a new synthesis strategy, dubbed *Iterative Cross-Coupling (ICC)*, that aims to make the process of complex-small-molecule synthesis as simple, efficient, and flexible as possible. The group is now harnessing the power of this new chemistry to systematically dissect the structure–function relationships that underlie the protein-like activities of a variety of prototypical small molecules, including the channel-forming polyene natural product amphotericin B. 

Purification of Laboratory Chemicals, 6th Edition


Description

This new edition keeps engineers, scientists, chemists, biochemists and students up-to-date with the processes for the purification of chemical reagents, and guides readers on critical safety and hazards for the safe handling of chemicals and processes. The 6th edition is updated and provides expanded coverage of the latest chemical products and processing techniques, safety, and hazards. The book has been reorganized and is now fully indexed by CAS Registry Numbers®. Compounds are now grouped to make navigation easier, and literature references for all substances and techniques have been added; ambiguous alternate names and cross references have been removed.

Table of Contents

Common Physical Techniques Used in Purification; Safety and Hazards; Distillation Techniques; Recrystallisation Techniques; Chromatography; Chemical Methods Used in Purification; Purification of Organic Compounds; Purification of Inorganic and Metal-Organic Compounds; Purification of Biochemicals and Related Products; Developing Techniques.


Sigma-Aldrich Cat. No. Z567930



By Wilfred L. F. Armarego and Christina L. L. Chai, Elsevier, 2009, 760 pp, ISBN: 978-1-85617-567-8

For more information, please visit us at sigma-aldrich.com/books.

CAS Registry Number is a registered trademark of American Chemical Society Corp



MIDA Boronates from Sigma-Aldrich



Greener Alternatives

Check out our new Web portal to find out about our Greener Products and Programs



Unlock a smarter, greener laboratory

Visit our new Greener Alternatives Web portal at sigma-aldrich.com/green to learn how you can minimize your lab's environmental impact—without compromising the quality of your data or research. Come see how Sigma-Aldrich is meeting the research community's needs for products and services which can significantly lessen environmental impact. Put our commitment to environmental sustainability to work in your lab.



Greener Products

Browse greener alternatives




Programs and Services

Discover solutions designed to reduce the environmental impact of your lab



Learning Center

Resources and links to information on environmental sustainability



For more information on our Greener Alternatives, Programs and Services, please visit us at

sigma-aldrich.com/green

Change Service Requested

Sigma-Aldrich® Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA S.A.
Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

SIGMA-ALDRICH PTY LTD.
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

SIGMA-ALDRICH NV/S.A.
Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

Chile

SIGMA-ALDRICH
QUIMICA LIMITADA
Tel: (+56) 2 495 7395
Fax: (+56) 2 495 7396

China

SIGMA-ALDRICH (SHANGHAI)
TRADING CO. LTD.
Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

SIGMA-ALDRICH spol. s r. o.
Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: (+45) 43 56 59 00
Fax: (+45) 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND OY
Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 00
Fax: (+33) 474 95 68 08

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Greece

SIGMA-ALDRICH (O.M.) LTD.
Tel: (+30) 210 994 8010
Fax: (+30) 210 994 3831

Hungary

SIGMA-ALDRICH Kft
Ingyenes telefonszám: 06 80 355 355
Ingyenes fax szám: 06 80 344 344
Tel: (+36) 1 235 9055
Fax: (+36) 1 235 9050

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: (+91) 80 6621 9400
New Delhi: (+91) 11 4358 8000
Mumbai: (+91) 22 2570 2364
Hyderabad: (+91) 40 4015 5488
Kolkata: (+91) 33 4013 8003
Fax
Bangalore: (+91) 80 6621 9650
New Delhi: (+91) 11 4358 8001
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 4015 5466
Kolkata: (+91) 33 4013 8016

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 402 20370
Fax: (+353) 402 20375

Israel

SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

SIGMA-ALDRICH S.r.l.
Numero Verde: 800 827018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

SIGMA-ALDRICH JAPAN K.K.
Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

SIGMA-ALDRICH KOREA
Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: (+52) 722 276 1600
Fax: (+52) 722 276 1601

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

SIGMA-ALDRICH NEW ZEALAND LTD.
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

SIGMA-ALDRICH NORWAY AS
Tel: (+47) 23 17 60 00
Fax: (+47) 23 17 60 10

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

SIGMA-ALDRICH RUS, LLC
Tel: (+7) 495 621 5579
Fax: (+7) 495 621 5923

Singapore

SIGMA-ALDRICH PTE. LTD.
Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

Slovakia

SIGMA-ALDRICH spol. s r. o.
Tel: (+421) 255 571 562
Fax: (+421) 255 571 564

South Africa

SIGMA-ALDRICH (PTY) LTD.
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313

United States

SIGMA-ALDRICH
Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Vietnam

SIGMA-ALDRICH PTE LTD. VN R.O.
Tel: (+84) 3516 2810
Fax: (+84) 6258 4238

Internet

sigma-aldrich.com

World Headquarters

3050 Spruce St., St. Louis, MO 63103
(314) 771-5765
sigma-aldrich.com

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052

Technical Service (800) 325-5832 • sigma-aldrich.com/techservice

Development/Custom Manufacturing Inquiries **SAPC**® (800) 244-1173

Accelerating Customers'
Success through Innovation and
Leadership in Life Science,
High Technology and Service

SIGMA, SAF, SAPC, SIGMA-ALDRICH, ALDRICH, FLUKA, and SUPELCO are trademarks belonging to Sigma-Aldrich Co. L.P. Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information provided on the label and in the literature. Purchaser must determine the suitability of the products for their particular use. Additional terms and conditions may apply.



Mixed Sources

Product group from well-managed
forests, controlled sources and
recycled wood or fiber
www.fsc.org Cert no. SGS-COC-003338
© 1996 Forest Stewardship Council



DEDICATED TO DR. ALFRED BADER ON THE OCCASION OF HIS 85TH BIRTHDAY

Aldrichimica ACTA

VOL. 42, NO. 1 • 2009

 **ALDRICH**
Chemistry



The Super Silyl Group in Diastereoselective Aldol
and Cascade Reactions

Iterative Cross-Coupling with MIDA Boronates: towards a
General Strategy for Small-Molecule Synthesis

SIGMA-ALDRICH

Editorial

"Please Call Me Alfred"

On April 28 of this year, Dr. Alfred Bader, undeniably the world's best-known Chemist Collector, celebrated his 85th birthday. Alfred's amazing life story has been covered extensively in books, magazines, and lectures, by Alfred himself and by others, and need not be repeated here. Furthermore, most of our readers are undoubtedly aware of Alfred's strong connections, past and present, not only to Sigma-Aldrich, but also to the *Aldrichimica Acta*, which he has showered with his attention for many years. Many of Alfred's paintings have graced the covers of the *Acta*, including this issue, which is featuring one of Alfred's favorite paintings.

We honor and thank Alfred for his outstanding contributions to Sigma-Aldrich and to the worlds of chemistry, business, and art. To what he calls the "ABC" (Art, Bible, Chemistry) of his life, one should add a "D" for Donating. Alfred's philanthropic activities are considerable and ongoing, and cover a wide range of causes that are near and dear to Alfred's heart.

The magnitude of Alfred's philanthropic efforts was made possible by the considerable financial rewards that he has reaped from two of his lifelong passions. The first is a spectacularly successful chemical business (Aldrich and Sigma-Aldrich) that he helped found and successfully managed for years. The second is his passion for collecting art works, particularly of Dutch and Flemish Masters, as well as rare stamps. His collection of about two hundred such paintings has been donated by Alfred to the Agnes Etherington Art

Centre of Queen's University in Kingston, Ontario (Canada). This gift is one of several sizeable ones that he has made to Queen's in gratitude for the education he received there in the 1940s.

Lesser known, but not less important to Alfred, is his Bible scholarship, particularly of Old Testament themes, which he has studied all his life and taught for a good many years. Another lesser known trait of Alfred is his modesty and unassuming lifestyle, which used to lead many a new employee of Aldrich, in the days when Alfred was company president, to think, when running into him for the first time, that he was just another employee. After meeting Alfred for the first time and calling him Dr. Bader, the editor of this publication was gently chided for calling him Dr. Bader, rather than Alfred, which is what he insists on being called even by people who don't know him that well.



Dr. Alfred Bader in 2005

Photo Courtesy of Alfred Bader

On behalf of all Sigma-Aldrich employees, past and present, we wish Alfred a very happy 85th birthday and many more in years to come.

SAFC® Launches New Vendor Audit Services


On March 19, 2009, SAFC® announced details of its new Vendor Audit Services, designed to extend audit reach for customers that need to secure global supply chain integrity through enhanced traceability and quality confirmation for critical raw materials and key intermediates. SAFC's Vendor Audit Services are an extension of in-house auditing programs, and provide pharmaceutical and biotechnology firms with a team of highly trained auditors that have extensive regulatory compliance expertise and are on the ground in local regions, including China.

Standard audits available include:

- ISO 9001:2008 (previously ISO 9001:2000)
- ISO 13485:2003
- 21 CFR 211
- 21 CFR 820
- ICH Q7A
- GAMP 5
- GLP
- 21CFR58

SAFC can also perform audits to customer specifications upon request.

SAFC's Typical Audit Services Timeline



Please contact us at safcglobal.com/auditservices for more information.

Aldrichimica ACTA

VOL. 42, NO. 1 • 2009

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-231-8327
SAFC® 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 414-438-3850
24-Hour Emergency 414-438-3850
Web Site sigma-aldrich.com
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 355, Milwaukee, WI 53201, USA

Subscriptions

To request your FREE subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 2988
Milwaukee, WI 53201-2988
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the back cover.

The *Aldrichimica Acta* is also available on the internet at sigma-aldrich.com/acta.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2009 Sigma-Aldrich Co.


"PLEASE BOTHER US."



Joe Porwoll

Joe Porwoll, President
Aldrich Chemical Co., Inc.

Professor Stephen Buchwald of the Massachusetts Institute of Technology kindly suggested that we make the single-component palladacycle precatalysts of SPhos and XPhos. These complexes simplify the use of these Buchwald ligands in a pre-defined metal-to-ligand ratio. The precatalysts are air- and moisture-stable, and can achieve high yields in C-N cross-coupling reactions using catalyst loadings as low as 0.1 mol % and 10-minute reaction times.



Biscoe, M. R. et al. *J. Am. Chem. Soc.* **2008**, *130*, 6686

704946 SPhos-palladium(II)phenethylamine chloride (1:1 MTBE solvate)

250 mg \$45.00
1 g 150.00

704954 XPhos-palladium(II)phenethylamine chloride (1:1 MTBE solvate)

250 mg \$69.00
1 g 190.00

Naturally, we made these useful catalyst precursors. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the back cover.

TABLE OF CONTENTS

The Super Silyl Group in Diastereoselective Aldol and Cascade Reactions <i>Matthew B. Boxer, Brian J. Albert, and Hisashi Yamamoto,* University of Chicago</i>	3
Iterative Cross-Coupling with MIDA Boronates: towards a General Strategy for Small-Molecule Synthesis <i>Eric P. Gillis and Martin D. Burke,* University of Illinois at Urbana-Champaign</i>	17

ABOUT OUR COVER


Two Squirrels (oil on panel, 33.0 × 40.6 cm) was painted around 1616 and is attributed to the Flemish painter Jan Brueghel (or Bruegel) the Elder (1568–1625), also known as "Velvet" Brueghel and "Flower" Brueghel. He was the second son of Pieter Brueghel the Elder and trained with his older brother, Pieter Brueghel the Younger, in the family workshop in Antwerp.

The painting, which reflects the painter's focus later on in life on painting flowers and animals, depicts two sprightly squirrels clutching a gnarly, twisting branch set against an empty sky. This simple, engaging, and humorous depiction of a pair of lively small animals parallels the practice, at that time, of hanging portraits of husband and wife. Its high level of finish suggests that it probably was executed on commission or for sale in the market. The luminous effect created by the thin pigment layers and the fine detail are reminiscent of the style of the Brueghel family.



Detail from *Two Squirrels*. Photograph © Alfred Bader

This painting is part of the Bader Collection of Dutch and Flemish Paintings at the Agnes Etherington Art Centre of Queen's University, Kingston, ON, Canada.




New Products from Aldrich R&D

Aldrich is Pleased to Offer Cutting-Edge Tools for Organic Synthesis

Double-Allylation Reagents

Reagents bearing multifunctional handles are of interest for the preparation of complex molecules. Professor Hall and coworkers developed a new multifunctional reagent that provides high diastereo- and enantiocontrol in a number of reactions, such as in the nucleophilic addition to aldehydes.




Peng, F.; Hall, D. G. *J. Am. Chem. Soc.* **2007**, *129*, 3070.

(+)-Allylboronic acid pinanediol ester, 97%

694584

C₁₃H₂₁BO₂

FW: 220.12



1 g \$34.60


5 g 112.00

(+)-Vinylboronic acid pinanediol ester, 95%

691615

C₁₂H₁₉BO₂

FW: 206.09




1 g \$66.00

5 g 218.50

Substrate for Nickel-Catalyzed Negishi Coupling

Chiral building blocks are of the utmost importance in the synthesis of more complex molecules. Professor Fu and coworkers devised the first catalytic enantioselective cross-coupling of secondary α -bromo amides with organozinc reagents. This new method proved to be highly selective for the coupling of unfunctionalized and functionalized organozincs with good yields.



Fischer, C.; Fu, G. C. *J. Am. Chem. Soc.* **2005**, *127*, 4594


N-Benzyl-2-bromo-N-phenylbutanamide, 97%

704911

[851073-30-8]

C₁₇H₁₈BrNO


FW: 332.23



5 g \$65.00

Copper(I) Fluoride Complex for Aldol Reaction

Chiral tertiary alcohols are important building blocks for the synthesis of more complex molecules such as biologically active compounds or potential drugs. Shibasaki and coworkers developed a new copper fluoride catalyzed aldol reaction of ketones using ketene silyl acetals. Various aromatic ketones were screened and led to the desired aldol products in good yields and high selectivity.



Oisaki, K.; Zhao, D.; Kanai, M.; Shibasaki, M. *J. Am. Chem. Soc.* **2006**, *128*, 7164. Taniaphos is a registered trademark of OMG AG and Co.

Fluorotris(triphenylphosphine)copper(I), 95% by HMR

706000

C₅₄H₄₅CuFP₃

FW: 869.40




250 mg \$39.00

1 g 95.00

Copper Chloride-Bis(lithium chloride) Solution for Transmetalation

Cross-coupling reactions are essential tools for chemists. In particular, the amination of aromatic halides has become a method highly relied upon to prepare aryl amines. Knochel and coworkers developed a new general amination procedure using amidocuprates. This method proved to be very versatile resulting in good yields of the amine products. This new method is a good complement to the Pd-catalyzed amination reactions.



Del Amo, V.; Dubbaka, S. R.; Krasovskiy, A.; Knochel, P. *Angew. Chem., Int. Ed.* **2006**, *45*, 7838

Copper(I) chloride-bis(lithium chloride) complex, 1 M in tetrahydrofuran

701742

CuCl·2(LiCl)

FW: 183.79



50 mL \$95.00

Greener Alternatives

Check out our new Web portal to find out about our Greener Products and Programs



Unlock a smarter, greener laboratory

Visit our new Greener Alternatives Web portal at sigma-aldrich.com/green to learn how you can minimize your lab's environmental impact—without compromising the quality of your data or research. Come see how Sigma-Aldrich is meeting the research community's needs for products and services which can significantly lessen environmental impact. Put our commitment to environmental sustainability to work in your lab.



Greener Products

Browse greener alternatives




Programs and Services

Discover solutions designed to reduce the environmental impact of your lab



Learning Center

Resources and links to information on environmental sustainability



For more information on our Greener Alternatives, Programs and Services, please visit us at sigma-aldrich.com/green



P.O. Box 14508
St. Louis, MO 63178
USA

PRESORTED
STANDARD
U.S. POSTAGE
PAID
SIGMA-ALDRICH
CORPORATION

Change Service Requested

Sigma-Aldrich® Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA S.A.
Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

SIGMA-ALDRICH PTY LTD.
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

SIGMA-ALDRICH NV/S.A.
Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

Chile

SIGMA-ALDRICH
QUIMICA LIMITADA
Tel: (+56) 2 495 7395
Fax: (+56) 2 495 7396

China

SIGMA-ALDRICH (SHANGHAI)
TRADING CO. LTD
Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

SIGMA-ALDRICH spol. s r. o.
Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: (+45) 43 56 59 00
Fax: (+45) 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND OY
Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 00
Fax: (+33) 474 95 68 08

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Greece

SIGMA-ALDRICH (O.M.) LTD.
Tel: (+30) 210 994 8010
Fax: (+30) 210 994 3831

Hungary

SIGMA-ALDRICH Kft
Ingyenes telefonszám: 06 80 355 355
Ingyenes fax szám: 06 80 344 344
Tel: (+36) 1 235 9055
Fax: (+36) 1 235 9050

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: (+91) 80 6621 9400
New Delhi: (+91) 11 4358 8000
Mumbai: (+91) 22 2570 2364
Hyderabad: (+91) 40 4015 5488
Kolkata: (+91) 33 4013 8003
Fax
Bangalore: (+91) 80 6621 9650
New Delhi: (+91) 11 4358 8001
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 4015 5466
Kolkata: (+91) 33 4013 8016

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 402 20370
Fax: (+353) 402 20375

Israel

SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

SIGMA-ALDRICH S.r.l.
Numero Verde: 800 827018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

SIGMA-ALDRICH JAPAN K.K.
Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

SIGMA-ALDRICH KOREA
Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: (+52) 722 276 1600
Fax: (+52) 722 276 1601

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

SIGMA-ALDRICH NEW ZEALAND LTD.
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

SIGMA-ALDRICH NORWAY AS
Tel: (+47) 23 17 60 00
Fax: (+47) 23 17 60 10

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

SIGMA-ALDRICH RUS, LLC
Tel: (+7) 495 621 5579
Fax: (+7) 495 621 5923

Singapore

SIGMA-ALDRICH PTE. LTD.
Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

Slovakia

SIGMA-ALDRICH spol. s r. o.
Tel: (+421) 255 571 562
Fax: (+421) 255 571 564

South Africa

SIGMA-ALDRICH (PTY) LTD.
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313

United States

SIGMA-ALDRICH
Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Vietnam

SIGMA-ALDRICH PTE LTD. VN R.O.
Tel: (+84) 3516 2810
Fax: (+84) 6258 4238

Internet

sigma-aldrich.com

World Headquarters

3050 Spruce St., St. Louis, MO 63103
314 771 5765
sigma-aldrich.com

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052

Technical Service (800) 325-5832 • sigma-aldrich.com/techservice

Development/Custom Manufacturing Inquiries **SAFC**® (800) 244-1173

©2009 Sigma-Aldrich Co. All rights reserved. SIGMA, SAF, SAFC, SIGMA-ALDRICH, ALDRICH, FLUKA, and SUPELCO are trademarks belonging to Sigma-Aldrich Co. and its affiliate Sigma-Aldrich Biotechnology, L.P. Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in the Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.



Mixed Sources

Product group from well-managed
forests, controlled sources and
recycled wood or fiber
www.fsc.org Cert. no. SGS-COC-003338
© 1996 Forest Stewardship Council



Accelerating Customers'
Success through Innovation and
Leadership in Life Science,
High Technology and Service

CHIRAL DIENE AND NHC LIGANDS FOR ASYMMETRIC CATALYSIS

Aldrichimica ACTA

VOL. 42, NO. 2 • 2009

 **ALDRICH**
Chemistry



Chiral Diene Ligands for Asymmetric Catalysis

Chiral, Chelating, Hydroxyalkyl and Hydroxyaryl N-Heterocyclic Carbenes:
Design, Synthesis, and Application in Copper-Catalyzed Asymmetric
Conjugate Addition (Cu-ACA)

SIGMA-ALDRICH




New Products from Aldrich R&D

Aldrich is pleased to offer cutting-edge tools for Organic Synthesis

Annulation Agent for Heterocyclic Synthesis

Professor Varinder Aggarwal recently reported an efficient annulation reaction to prepare a variety of 1,4-heterocyclic products in good-to-excellent yields. The mild reaction employs a bromoethylsulfonium salt in the presence of 1,2- and 1,3-amino alcohols, amino thiols, or diamines to give access to pharmacologically important morpholines, piperazines, thiomorpholines, benzoxazepines, and benzodiazepines



Yar, M. et al. *Org. Lett.* **2009**, *11*, 257.

(2-Bromoethyl)diphenylsulfonium trifluoromethanesulfonate

706345
C₁₅H₁₈BrF₃O₂S₂
FW: 443.30

1 g \$39.50
5 g 130.00

NEW Solutions of Common Reagents

Every researcher has experienced the frustration of using some lab reagent on a regular basis that is annoying to handle: it might be difficult to weigh out, is sticky, prone to static or nonflow, or flammable. Aldrich has designed a series of solutions in tetrahydrofuran to make many of these common reagents easier to measure, handle, and dispense.

Di-tert-butyl dicarbonate, 1 M in THF

436267
[24424-99-5]
C₁₀H₁₈O₅
FW: 218.25

100 mL \$41.50
800 mL \$330.00

Di-tert-butyl dicarbonate, 2 M in THF

703753
[24424-99-5]
C₁₀H₁₈O₅
FW: 218.25

100 mL \$51.50

Di-tert-butyl dicarbonate, 2 M in dichloromethane

703737
[24424-99-5]
C₁₀H₁₈O₅
FW: 218.25

100 mL \$51.50

Di-tert-butyl dicarbonate, 2 M in ethyl acetate

703745
[24424-99-5]
C₁₀H₁₈O₅
FW: 218.25

100 mL \$51.50

(S)-(-)-α,α-Diphenyl-2-pyrrolidinemethanol trimethylsilyl ether, 0.05 M in toluene

706442
[848821-58-9]
C₂₁H₂₇NOSi
FW: 325.52

25 mL \$49.50

(R)-(+)-α,α-Diphenyl-2-pyrrolidinemethanol trimethylsilyl ether, 0.05 M in toluene

706450
C₂₁H₂₇NOSi
FW: 325.52

25 mL \$49.50

1,1'-Carbonyldiimidazole, 0.4 M in dichloromethane

705284
[530-62-1]
C₇H₈N₄O
FW: 162.15

100 mL \$30.70

Triflic anhydride, 1 M in dichloromethane

704083
[598-23-6]
C₂F₆O₅
FW: 214.14

25 mL \$53.00
100 mL 145.50

Tributyltin hydride 1 M in cyclohexane

704091
[688-73-3]
C₁₁H₂₂Sn
FW: 291.06

10 mL \$28.50
50 mL 95.00

Iodine monochloride, 1 M in dichloromethane

291048
[7790-99-0]
ICl
FW: 162.36

100 mL \$41.80
800 mL 195.50

Acetyl chloride, 1 M in dichloromethane

708496
[75-36-5]
C₂H₃ClO
FW: 78.50


100 mL \$23.10
500 mL 77.00



Chiral Diene Technology


In the past 5 years, a number of chiral diene ligands have emerged for a variety of asymmetric transformations. Hayashi¹ and Carreira² pioneered this field by synthesizing chiral diene ligands that formed stable complexes with metals and exerted high catalytic activities as well as high enantioselectivities. To demonstrate the potency of these new ligands, Hayashi and co-workers synthesized

tolterodine, an antimuscarinic drug used to treat urinary incontinence, in a five-step sequence with an overall yield of 61%.³ The first step of the sequence is an asymmetric 1,4 addition of phenylboronic acid, using a chiral diene ligand with a rhodium complex, to give the corresponding phenylated product in 96% yield.



References:


- (1) Havasni, T. et al. *J Am Chem Soc* **2003**, 125, 11508
 (2) Fischer, C. et al. *J Am Chem Soc* **2004**, 126, 11508
 (3) Hayashi, T. et al. *Org. Lett* **2008**, 10, 589



679941



707406




707392




707163




707171




709492



669490



698725



698733

Aldrich is pleased to offer these chiral diene ligands for use in a variety of applications. For additional information, please visit sigma-aldrich.com/cpc



Easy Access to Freshly Prepared Anhydrous Solvents




Image shows four Pure Solv™ Micro systems side by side for processing more than one anhydrous solvent. Vacuum pump is not included

- Bench scale for laboratory use
- Assured solvent quality
- No electrical power
- Low-pressure operation
- Safe, nonsparking design

Pure Solv™ Micro Solvent Purification System Features:

- **Space-Saving:** unit (12 in. W × 13 in. D × 19 in. H) fits on bench top or in fume hood.
- **Economical:** disposable columns remove water from 100 L of raw solvent to low ppm range.
- **Integral Vacuum/Gas Purge Manifolds:** facilitate proper air-free dispensing techniques and include a size 24/40 stainless steel glassware adapter.
- **Low-Pressure Operation:** less than 5 psig of inert gas required.
- **4 L Raw Solvent Reservoir:** with leak-tight, clamp-style cap and Swagelok® quick-disconnect valves for easy handling and refilling.
- **Nonsparking Design:** systems are bonded and grounded to remove hazards associated with electrostatic discharge and ensure a safe operating environment.

Order Pure Solv Micro with the appropriate activated purification column. Solvent purification columns are disposable. Additional activated columns may be ordered when replacement is necessary.

For more information on Pure Solv Micro systems selection guide, please visit us at sigma-aldrich.com/puresolv

Pure Solv is a trademark of Innovative Technology, Inc. Swagelok is a registered trademark of Crawford Fitting Co



P.O. Box 14508
St. Louis, MO 63178
USA

PRESORTED
STANDARD
U.S. POSTAGE
PAID
SIGMA-ALDRICH
CORPORATION

Change Service Requested

Sigma-Aldrich® Worldwide Locations

Argentina

SIGMA-ALDRICH DE ARGENTINA S.A.
Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

SIGMA-ALDRICH PTY LTD
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

SIGMA-ALDRICH NV/S.A.
Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA
Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

Chile

SIGMA-ALDRICH
QUIMICA LIMITADA
Tel: (+56) 2 495 7395
Fax: (+56) 2 495 7396

China

SIGMA-ALDRICH (SHANGHAI)
TRADING CO. LTD.
Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

SIGMA-ALDRICH spol. s r. o.
Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: (+45) 43 56 59 00
Fax: (+45) 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND OY
Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 00
Fax: (+33) 474 95 68 08

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Greece

SIGMA-ALDRICH (O.M.) LTD.
Tel: (+30) 210 994 8010
Fax: (+30) 210 994 3831

Hungary

SIGMA-ALDRICH Kft
Ingyenes telefonszám: 06 80 355 355
Ingyenes fax szám: 06 80 344 344
Tel: (+36) 1 235 9055
Fax: (+36) 1 235 9050

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: (+91) 80 6621 9400
New Delhi: (+91) 11 4358 8000
Mumbai: (+91) 22 2570 2364
Hyderabad: (+91) 40 4015 5488
Kolkata: (+91) 33 4013 8003
Fax
Bangalore: (+91) 80 6621 9650
New Delhi: (+91) 11 4358 8001
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 4015 5466
Kolkata: (+91) 33 4013 8016

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 402 20370
Fax: (+353) 402 20375

Israel

SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

SIGMA-ALDRICH S.r.l.
Numero Verde: 800 827018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

SIGMA-ALDRICH JAPAN K.K.
Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

SIGMA-ALDRICH KOREA
Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: (+52) 722 276 1600
Fax: (+52) 722 276 1601

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

SIGMA-ALDRICH NEW ZEALAND LTD.
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

SIGMA-ALDRICH NORWAY AS
Tel: (+47) 23 17 60 00
Fax: (+47) 23 17 60 10

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

SIGMA-ALDRICH RUS, LLC
Tel: (+7) 495 621 5579
Fax: (+7) 495 621 5923

Singapore

SIGMA-ALDRICH PTE. LTD.
Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

Slovakia

SIGMA-ALDRICH spol. s r. o.
Tel: (+421) 255 571 562
Fax: (+421) 255 571 564

South Africa

SIGMA-ALDRICH (PTY) LTD.
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313

United States

SIGMA-ALDRICH
Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Vietnam

SIGMA-ALDRICH PTE LTD. VN R.O.
Tel: (+84) 3516 2810
Fax: (+84) 6258 4238

Internet

sigma-aldrich.com

World Headquarters

3050 Spruce St., St. Louis, MO 63103
(314) 771-5765
sigma-aldrich.com

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052

Technical Service (800) 325-5832 • sigma-aldrich.com/techservice

Development/Custom Manufacturing Inquiries **SAPC**® (800) 244-1173

Safety-Related Information sigma-aldrich.com/safetycenter

Accelerating Customers'
Success through Innovation and
Leadership in Life Science,
High Technology and Service

MA, SAPC, **SAPC**®, SIGMA-ALDRICH, ALDRICH, and SUPELCO are trademarks belonging to Sigma-Aldrich Co. and its affiliate companies. Sigma-Aldrich Co. warrants that its products conform to the specifications and standards of the International Organization of Standardization (ISO) and other applicable standards. © 1996 Forest Stewardship Council



Mixed Sources
Product group from well-managed
forests, controlled sources and
recycled wood or fiber
www.fsc.org Cert no. SGS-COC-003338
© 1996 Forest Stewardship Council



OPENING NEW FRONTIERS IN SYNTHESIS AND CATALYSIS

Aldrichimica ACTA

VOL. 42, NO. 3 • 2009

 **ALDRICH**
Chemistry



Discovering New Reactions with N-Heterocyclic Carbene Catalysis
Synthesis and Applications of Diorganozinc Reagents: Beyond Diethylzinc


SIGMA-ALDRICH

New Products from Aldrich R&D

Aldrich is Pleased to Offer Cutting-Edge Tools for Organic Synthesis


Phosphine for the Conversion of Azides into Diazo Compounds

Due to difficulties with their preparation, especially when sensitive functional groups are present, diazo compounds are often overlooked in synthesis despite their synthetic versatility. Myers and Raines have developed a mild method to convert an azido group with delicate functional groups in to a diazo compound by using the phosphine reagent shown below. Formally, this reaction is a reductive fragmentation of the azide, like the venerable Staudinger reaction, and is highly selective in most chemical environments




Myers, E. L.; Raines, R. T. *Angew. Chem., Int. Ed.* **2009**, *48*, 2359

N-Succinimidyl 3-(diphenylphosphino)propionate, 95%




Stable Precursor for Nazarov's Reagent

Nazarov's reagent is a commonly used annulating agent, but its synthesis is often fraught with poor yields or difficulties isolating the material. De Risi and coworkers recently reported a bench-stable powder that can be demasked under various conditions to generate Nazarov's reagent in situ



Benetti, S. et al. *Synlett* **2008**, 2609


Ethyl 5-[(4-methylphenyl)sulfonyl]-3-oxopentanoate, 95%




More NEW Solutions of Common Reagents

Every researcher has experienced the frustration of using on a regular basis some lab reagent that is annoying to handle: It might be difficult to weigh out, extremely volatile, prone to static, or noxious. Sigma-Aldrich has designed a series of solutions intended to make many of these common reagents easier to measure, handle, and dispense


2,2'-Azobis(2-methylpropionitrile) solution, 0.2 M in toluene




Iodine monochloride, 1 M in acetic acid




4-(Dimethylamino)pyridine solution, 0.5 M in ethyl acetate




4-(Dimethylamino)pyridine solution, 0.5 M in THF




1,8-Diazabicyclo[5.4.0]undec-7-ene solution, 1 M in ethyl acetate




1,8-Diazabicyclo[5.4.0]undec-7-ene solution, 1 M in THF



Acetaldehyde solution, 5 M in THF



Carbon disulfide, 5 M in THF



Aldrichimica ACTA

VOL. 42, NO. 3 • 2009

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-231-8327
SAFC® 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 414-438-3850
24-Hour Emergency 414-438-3850
Website sigma-aldrich.com
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 2988, Milwaukee, WI 53201, USA

Subscriptions

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 2988
Milwaukee, WI 53201-2988
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com/acta.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2009 Sigma-Aldrich Co.

"PLEASE BOTHER US."




Joe Porwoll

Joe Porwoll, President
Aldrich Chemical Co., Inc.

Dr. Biagetti Matteo from GlaxoSmithKline kindly suggested that we make lithium diisobutyl-t-butoxyaluminum hydride (LDBBA). LDBBA is an effective reducing agent for the conversion of esters into aldehydes without the problematic over-reduction or the inconvenience of requiring lower temperatures

Kim, M. S.; Choi, Y. M.; An, D. K. *Tetrahedron Lett.* **2007**, *48*, 5061



718386

718386 Lithium diisobutyl-t-butoxyaluminum hydride (LDBBA), 0.5 M in THF-hexanes

25 mL \$73.00
100 mL 201.00

Naturally, we made this useful reducing agent. It was no bother at all, just a pleasure to be able to help

Do you have a compound that you wish Aldrich could list, and that would help your research by saving you time and money? If so, please send us your suggestions and we will be glad to give it careful consideration. You can contact us in any one of the ways shown on this page and on the back cover.

TABLE OF CONTENTS

Discovering New Reactions with N-Heterocyclic Carbene Catalysis..... 55
Eric M. Phillips, Audrey Chan, and Karl A. Scheidt,* Northwestern University

Synthesis and Applications of Diorganozinc Reagents: Beyond Diethylzinc 71
Alexandre Lemire, Alexandre Côté, Marc K. Janes, and André B. Charette,* University of Montreal

ABOUT OUR COVER

Paul Cézanne (1839–1906), who painted *Landscape near Paris* (oil on canvas, 50.2 × 60 cm) around 1876, was born in Aix-en-Provence in 1839. His father, a prosperous businessman, decided that his only son should become a lawyer. Cézanne attended the Aix law school, but preferred classes at the Musée d'Aix (now the Musée Granet) and decided on a life as an artist instead.

In 1861, Cézanne, encouraged by his boyhood friend, the novelist Emile Zola, traveled to Paris. There, he frequented the Salon, studied the old masters, and copied Delacroix at the Louvre. He also forged friendships with many important artists, one of whom, Camille Pissarro, became a pivotal, lifelong influence on Cézanne.

Cézanne exhibited with the impressionists in 1874 and 1877. Our cover possibly painted in the company of Pissarro in Auvers-sur-Oise, a northwestern suburb of Paris, was completed sometime between the two exhibitions. The work clearly denotes Cézanne's departure from his early romantic and realist influences, and displays his enduring interest in *plein-air* painting. In this work, Cézanne placed an emphasis on the observation of nature and the rendering of light and atmospheric effects. He achieved structure by applying paint directly to the canvas, recording his response to the sensation of color.

This painting is part of the Chester Dale Collection at the National Gallery of Art, Washington, DC.



Photograph © Board of Trustees, National Gallery of Art, Washington




N-Heterocyclic Carbene Organocatalysts

Organocatalysis has been an active field of research with over 2,000 publications in the past 10 years. Interest in this relatively new field derives from the many advantages of organocatalysis such as easy experimental procedures, reduction of chemical waste, avoidance of metal contamination in the product, and cost saving from not using expensive metals for the catalysis.

Scheidt and coworkers have developed a series of N-heterocyclic carbene organocatalysts for various reactions. These catalysts are efficient, and the chiral ones have given rise to good enantioselectivities.


In 2005, Audrey Chan and Karl A. Scheidt reported the transformation of unsaturated aldehydes into saturated esters in good yields by using as low as 5 mol % of an N-heterocyclic carbene catalyst.









Reference:
Scheidt, K. A. *Org. Lett.* **2005**, *7*, 905

For ordering or more information, please visit
sigma-aldrich.com






Recently, Phillips et al. reported the utilization of a chiral triazole for the enantioselective addition of homoenolates to nitrones affording γ -amino esters. This reaction is very versatile and tolerates electron-rich and electron-poor groups on the aldehyde. Using 20 mol % of the chiral triazole at -25°C affords the desired products in good yields and selectivities.



Reference:
Phillips, E. M. et al. *J. Am. Chem. Soc.* **2008**, *130*, 2416



Phosphonics Oxidation Catalyst Kit
709034

Oxidation Catalyst	Substrate	Product	Re-Oxidant
	Allylic/Benzylic CH ₂ Allylic alcohols	Ketones Enones	<i>t</i> -BuOOH
	Allylic alcohols Sulfides	Epoxides Sulfoxides	<i>t</i> -BuOOH or NaBrO or H ₂ O
	Allylic CH ₂ or Benzylic CH	Ketones	<i>t</i> -BuOOH
	1° alcohols 2° alcohols Sulfides	Acids Ketones Sulfoxides	NaBrO or <i>t</i> -BuOOH
	Sulfides	Sulfoxides	NaBrO ₂ or <i>t</i> -BuOOH

For more information, please visit sigma-aldrich.com

Phosphonics is a trademark of Phosphonics Ltd



P.O. Box 14508
St. Louis, MO 63178
USA

PRESORTED
STANDARD
U.S. POSTAGE
PAID
SIGMA-ALDRICH
CORPORATION

Change Service Requested

Sigma-Aldrich® Worldwide Offices

Argentina

Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

Chile

Tel: (+56) 2 495 7395
Fax: (+56) 2 495 7396

China

Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

Tel: (+45) 43 56 59 00
Fax: (+45) 43 56 59 05

Finland

Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 88
Fax: (+33) 474 95 68 08

Germany

Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Hungary

Ingyenes telefonszám: 06 80 355 355
Ingyenes fax szám: 06 80 344 344
Tel: (+36) 1 235 9063
Fax: (+36) 1 269 6470

India

Telephone
Bangalore: (+91) 80 6621 9400
New Delhi: (+91) 11 4358 8000
Mumbai: (+91) 22 4087 2364
Hyderabad: (+91) 40 4015 5488
Kolkata: (+91) 33 4013 8000

Fax

Bangalore: (+91) 80 6621 9550
New Delhi: (+91) 11 4358 8001
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 4015 5466
Kolkata: (+91) 33 4013 8016

Ireland

Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 402 20370
Fax: (+353) 402 20375

Israel

Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

Free Tel: 800 827 018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: (+52) 722 276 1600
Fax: (+52) 722 276 1601

The Netherlands

Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

Tel: (+47) 23 17 60 00
Fax: (+47) 23 17 60 10

Poland

Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

Tel: (+7) 495 621 5828
Fax: (+7) 495 621 6037

Singapore

Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

Slovakia

Tel: (+421) 255 571 562
Fax: (+421) 255 571 564

South Africa

Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313

United States

Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Vietnam

Tel: (+84) 3516 2810
Fax: (+84) 6258 4238

Internet

sigma-aldrich.com

Accelerating Customers'
Success through Innovation and
Leadership in Life Science,
High Technology and Service

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052
Technical Service (800) 325-5832 • sigma-aldrich.com/techservice
Development/Custom Manufacturing Inquiries **SAPC**® (800) 244-1173
Safety-related Information sigma-aldrich.com/safetycenter

World Headquarters
3050 Spruce St.
St. Louis, MO 63103
(314) 771-5765
sigma-aldrich.com



SIGMA-ALDRICH

EXPLOITING RING STRAIN IN METAL-FREE SYNTHESIS

Aldrichimica ACTA

VOL. 43, NO. 1 • 2010

 **ALDRICH**
Chemistry



Powerful and Versatile Silicon Lewis Acids for Asymmetric Chemical Synthesis
Copper-Free Click Chemistry: Bioorthogonal Reagents for Tagging Azides


SIGMA-ALDRICH

New Products from Aldrich R&D

Aldrich is Pleased to Offer Cutting-Edge Tools for Organic Synthesis

Efficient Lewis Base for Mukaiyama Aldol Reaction with Ketones

There have been few reported methods for the synthesis of tertiary aldols, due to the low reactivity of the starting ketones and the rapid retro-aldol reaction, even at low temperatures. Ishihara and coworkers recently reported the Lewis base catalyzed Mukaiyama aldol reaction using a simple, mixed sodium phenoxide-phosphine oxide catalyst. The catalyst was effective with a wide variety of TMS enolates and ketones, providing the aldol adducts in generally good-to-excellent yields. Additionally, the catalyst was also effective in Mannich-type reactions when benzyl-, Boc-, or Cbz-protected aldimines were used as substrates




Hatano, M. et al. *Org. Lett.* **2007**, 9, 4527

Bis[1,2-bis(diphenylphosphine oxide)benzene] sodium phenoxide, 95%

708836

C₂₆H₂₀NaO₂P₂

FW: 1073.01



1 g

Indoles and Azaindoles

The indole subunit is a near-ubiquitous component of biologically active natural products, and the study of indoles has been a major focus of research for generations. Substituted indoles have frequently been referred to as privileged structures¹ since they are capable of binding to multiple receptors with high affinity, and thus have applications across a wide range of therapeutic areas. Due to this activity, it is not surprising that the indole ring system has become an important structural motif in many pharmaceutical agents.

The azaindole moiety differs only by the presence of an additional ring nitrogen and, thus, it exhibits excellent potential as a bioisostere of the indole ring system. Although considerably more rare in nature, azaindoles still constitute essential subunits in many pharmaceutically important compounds, and have been very valuable to synthetic and medicinal chemists. 7-Azaindoles are of particular interest because of their ability to mimic purine in its role as a hydrogen-bonding partner

(1) Horton, D. A. et al. *Chem. Rev.* **2003**, 103, 893 and references therein. (2) (a) Popowycz, F. et al. *Tetrahedron* **2007**, 63, 1031. (b) Popowycz, F. et al. *Tetrahedron* **2007**, 63, 8659. (c) Song, J. J. et al. *Chem. Soc. Rev.* **2007**, 36, 1120


4-Bromo-1-methylindole, 95%

720860

[590417-55-3]

C₉H₈BrN

FW: 210.07



500 mg


5-Bromo-1-methylindole, 97%

718300

[10075-52-2]

C₉H₈BrN

FW: 210.07




5 g

5-Bromo-7-methylindole, 97%

710822

C₉H₈BrN

FW: 210.07




1 g

6-Bromo-7-methylindole, 97%

707155

C₉H₈BrN

FW: 210.07



250 mg


6-Fluoro-7-methylindole, 97%

707147

[57817-10-4]

C₉H₈FN

FW: 149.16



250 mg


5,7-Dibromoindole, 97%

708798

[36132-08-8]

C₈H₆Br₂N

FW: 274.94




1 g

7-Fluoro-5-iodoindole, 97%

707058

C₈H₆FIN

FW: 261.03



500 mg


5-Methoxy-4-azaindole, 95%

707953

[17288-40-3]

C₈H₈N₂O

FW: 148.16



250 mg


6-Methoxy-7-azaindole, 97%

702838

[896722-53-5]

C₈H₈N₂O

FW: 148.16



250 mg

7-Azaindole-4-carbonitrile, 97%

706426

[344327-11-3]

C₈H₆N₂

FW: 143.15



250 mg

Aldrichimica ACTA

VOL. 43, NO. 1 • 2010

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries	800-325-3010
Technical Service	800-231-8327
SAFC®	800-244-1173
Custom Synthesis	800-244-1173
Flavors & Fragrances	800-227-4563
International	414-438-3850
24-Hour Emergency	414-438-3850
Website	sigma-aldrich.com
Email	aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 2988, Milwaukee, WI 53201, USA

Subscriptions

To request your FREE subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 2988
Milwaukee, WI 53201-2988
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com/acta.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2010 Sigma-Aldrich Co.


"PLEASE BOTHER US."



Joe Porwoll, President
Aldrich Chemical Co., Inc

Professor Phil Baran of The Scripps Research Institute kindly suggested that we make silver picolinate. Baran and coworkers recently utilized silver picolinate to effect a chemoselective oxidation of the guanidine ring en route to their notable total synthesis of palau'amine.¹ This approach avoids overoxidation and laborious protecting group chemistry. Silver picolinate was also employed in the total synthesis of (±)-axinellamines A and B

(1) Seiple, I. B.; Su, S.; Young, I. S.; Lewis, C. A.; Yamaguchi, J.; Baran, P. S. *Angew. Chem., Int. Ed.* **2010**, *49*, 1095. (2) O'Malley, D. P.; Yamaguchi, J.; Young, I. S.; Seiple, I. B.; Baran, P. S. *Angew. Chem., Int. Ed.* **2008**, *47*, 3581. (3) Burns, N. Z.; Baran, P. S.; Hoffmann, R. W. *Angew. Chem., Int. Ed.* **2009**, *48*, 2854



718157

718157 Silver(II) picolinate, 97%

5 g
25 g

Naturally, we made this useful oxidation reagent. It was no bother at all, just a pleasure to be able to help

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the back cover

TABLE OF CONTENTS

Powerful and Versatile Silicon Lewis Acids for Asymmetric Chemical Synthesis..... 3
James L. Leighton, Columbia University

Copper-Free Click Chemistry: Bioorthogonal Reagents for Tagging Azides 15
Jeremy M. Baskin and Carolyn R. Bertozzi, University of California, Berkeley

ABOUT OUR COVER

The Grand Canal (oil on canvas, 23.0 × 33.0 cm) is one of twenty works of art by the British artist Richard Parkes Bonington (1802-1828) in the collection of the National Gallery of Art, and was painted in 1815/1827. At the age of 16, Bonington moved to Paris with his parents and studied there at the French National Art School, the Ecole des Beaux-Arts. During his brief career as an artist, he adopted the ideals of the Romantic Movement, and was influenced by the French artists Antoine-Jean Gros and Eugène Delacroix and by his fellow Englishman, John Constable



Bonington was greatly admired for his exceptional ability to capture the effects of daylight and atmosphere with unerring assurance. In this work, the lovely play of sunlight on the building facades, the delicate reflections on the water, and the sweep of the clouds across the sky entice the eye and can be appreciated independently of the subject and its imitative appeal

In the nineteenth century, painted scenes of Venice were very popular with the "Grand Tour" set, as post cards or photographs are today. Bonington produced many studies of this most beautiful city on the Adriatic while visiting Italy in 1826. He turned his sketches into some of his finest paintings, like the one seen here, which further enhanced his flourishing reputation as a landscape painter in London and Paris


This painting is a gift of Roger and Victoria Sant to the National Gallery of Art, Washington, DC.




Leighton's Chiral Silane Reagents

The asymmetric allylation of carbonyl compounds remains one of the most important and fundamental addition reactions for the synthesis of optically active chiral building blocks.

In 2002, Leighton and co-workers developed strained silacycle compounds as versatile reagents for the practical enantioselective allylation of aldehydes.¹ A newly developed chiral auxiliary based on the cyclohexane-1,2-diamine scaffold successfully allylated a broad range of aldehydes highly enantioselectively.²




The development of practical enantioselective syntheses of chiral amines is of great importance to synthetic organic and medicinal chemists. In 2003, Leighton and co-workers successfully used a pseudoephedrine-derived, five-membered-ring, strained silacycle reagent for the enantioselective allylation of acylhydrazones.³



The reaction scope of these silacycles was extended to a practical method for the enantioselective synthesis of tertiary carbinamines based on the addition of this chiral allylsilane reagent to a structurally diverse array


of ketone-derived benzoylhydrazones.⁴ While many methods for the synthesis of quaternary α -amino acids have been published, far fewer reports have dealt with the synthesis of tertiary carbinamines. The free amines are easily accessed in good yields by treating the product hydrazides with SmI_2 .



The reagents are effective with other types of imine electrophiles, and may be derivatized efficiently using cross-metathesis reactions. In addition, the corresponding phenylsilanes are effective general Lewis acids for a variety of enantioselective (nonallylation) nucleophilic addition reactions with acylhydrazones and other imine electrophiles.

References:

- (1) Kinnaird, J. W. A.; Ng, P. Y.; Kubota, K.; Wang, X.; Leighton, J. L. *J. Am. Chem. Soc.* **2002**, 124, 7920. (2) Kubota, K.; Leighton, J. L. *Angew. Chem., Int. Ed.* **2003**, 42, 946. (3) Berger, R.; Rabbat, P. M. A.; Leighton, J. L. *J. Am. Chem. Soc.* **2003**, 125, 9596. (4) Berger, R.; Duff, K.; Leighton, J. L. *J. Am. Chem. Soc.* **2004**, 126, 5686



For more information on the applications of the Leighton silacycle reagents, please see Professor Leighton's review in this issue or visit sigma-aldrich.com

of soluble tetraazido star polymers (red) and bicyclooctyne-containing cross-linkers (blue) causes the formation of a 3D gel network (triazole linkages are shown in purple). The use of different star polymers and cross-linkers enables the modulation of gel properties (e.g., solubility, ability to be degraded by UV light, ability to be further functionalized post-gelation with imaging agents, drugs, or other probes). Incorporation of different cyclooctynes into the cross-linker (e.g., MOFO vs DIFO) allows for the tuning of gelation kinetics. See also reference 48.

- (48) DeForest, C. A.; Polizzotti, B. D.; Anseth, K. S. *Nat. Mater.* **2009**, *8*, 659.
- (49) Baskin, J. M.; Bertozzi, C. R. *QSAR Comb. Sci.* **2007**, *26*, 1211.
- (50) Blackman, M. L.; Royzen, M.; Fox, J. M. *J. Am. Chem. Soc.* **2008**, *130*, 13518.
- (51) Devaraj, N. K.; Weissleder, R.; Hilderbrand, S. A. *Bioconjugate Chem.* **2008**, *19*, 2297.
- (52) Pipkorn, R.; Waldeck, W.; Didinger, B.; Koch, M.; Mueller, G.; Wiessler, M.; Braun, K. *J. Pept. Sci.* **2009**, *15*, 235.
- (53) Devaraj, N. K.; Upadhyay, R.; Haun, J. B.; Hilderbrand, S. A.; Weissleder, R. *Angew. Chem., Int. Ed.* **2009**, *48*, 7013.
- (54) Gutmiedl, K.; Wirges, C. T.; Ehmke, V.; Carell, T. *Org. Lett.* **2009**, *11*, 2405.
- (55) Dondoni, A. *Angew. Chem., Int. Ed.* **2008**, *47*, 8995.
- (56) Song, W.; Wang, Y.; Qu, J.; Madden, M. M.; Lin, Q. *Angew. Chem., Int. Ed.* **2008**, *47*, 2832.
- (57) Song, W.; Wang, Y.; Qu, J.; Lin, Q. *J. Am. Chem. Soc.* **2008**, *130*, 9654.
- (58) Wang, Y.; Song, W.; Hu, W. J.; Lin, Q. *Angew. Chem., Int. Ed.* **2009**, *48*, 5330.

Trademarks: ALEXA FLUOR[®] (Molecular Probes, Inc.), Hoechst (Hoechst GmbH Corp.).

Keywords: click chemistry; [3 + 2] cycloaddition; bioorthogonal reagent; cyclooctynes; azides.

About the Authors

Jeremy M. Baskin was born in Montreal, Canada. He received his B.S. degree in chemistry in 2004 from the Massachusetts Institute of Technology (M.I.T.), with minors in biology and music. While at M.I.T., he performed research in the laboratories of Professor Stephen L. Buchwald and Professor Alice Y. Ting. In 2004, Jeremy began graduate studies in the laboratory of Professor Carolyn R. Bertozzi at the University of California,

Berkeley. As a graduate student, Jeremy's research focused on the development of fluorinated cyclooctyne reagents for copper-free click chemistry and their application to imaging glycans in vivo. He earned his Ph.D. degree in 2009 and is currently conducting postdoctoral research under the guidance of Professor Pietro De Camilli at the Yale School of Medicine.

Carolyn R. Bertozzi is the T. Z. and Irmgard Chu Distinguished Professor of Chemistry and Professor of Molecular and Cell Biology at UC Berkeley; an Investigator of the Howard Hughes Medical Institute; and Director of the Molecular Foundry, a DOE Nanoscale Science and Research Center at the Lawrence Berkeley National Laboratory. She received her undergraduate degree in chemistry from Harvard University in 1988 and her Ph.D. degree in chemistry from UC Berkeley in 1993. After postdoctoral work at UC San Francisco in the field of cellular immunology, she joined the UC Berkeley faculty in 1996.

Professor Bertozzi's research interests span the disciplines of chemistry and biology with an emphasis on studies of cell surface glycosylation pertinent to disease states. Her lab focuses on profiling changes in cell surface glycosylation associated with cancer, inflammation, and bacterial infection, and on exploiting this information for the development of diagnostic and therapeutic approaches. In addition, her group develops nanoscience-based technologies for probing cell function and for medical diagnostics.

Professor Bertozzi has been recognized with many honors and awards for both her research and teaching accomplishments. She is an elected member of the National Academy of Sciences, the American Academy of Arts and Sciences, and the German Academy of Sciences Leopoldina. She has been awarded the Whistler Award, the Ernst Schering Prize, a MacArthur Foundation Fellowship, the ACS Award in Pure Chemistry, a Presidential Early Career Award in Science and Engineering, and the Irving Sigal Young Investigator Award of the Protein Society, among many others. Her efforts in undergraduate education have earned her a UC Berkeley Distinguished Teaching Award and the Donald Sterling Noyce Prize for Excellence in Undergraduate Teaching. Professor Bertozzi participates in high-school outreach programs such as the Catalyst Program sponsored by the Camille and Henry Dreyfus Foundation, as well as programs that promote the participation of women in science. She was recently presented with the Li Ka Shing Award for Women in Science in recognition of these efforts. ☞

Of Chemistry, Art, and Other Matters...

Read all about them in the two fascinating autobiographies by Dr. Alfred Bader, cofounder and former president of Aldrich and Sigma-Aldrich and one of the most influential chemistry business leaders and art collectors of the past fifty some years.

Bader, A. *Adventures of a Chemist Collector*; Weidenfeld & Nicolson: London, U.K., 1995.

Bader, A. *Chemistry & Art: Further Adventures of a Chemist Collector*; Weidenfeld & Nicolson: London, U.K., 2008.

The two volumes are available from Amazon.com. Order your copies today.



New Catalysis Products from Aldrich R&D

Palladium(II) Acetate, Recrystallized

It has been shown that in some applications recrystallized palladium(II) acetate, Pd(OAc)₂, performs better than typical grades. White and Overman have both independently demonstrated that in some applications a particular grade of Pd(OAc)₂ is necessary. In particular, the preparation of the White catalyst as well as the preparation of the COP catalysts—[COP-OAc]₂ and [COP-Cl]₂ by Overman successfully employ and require recrystallized Pd(OAc)₂.




720070

References: Anderson, C. E.; Kirsch, S. F.; Overman, L. E.; Richards, C. J.; Watson, M. P. *J. Org. Synth.* **2007**, 84, 148

Cationic Palladium Complexes, [Pd(dppp)(PhCN)₂](BF₄)₂

Cationic palladium(II) complexes are utilized in a variety of reactions. [Pd(dppp)(PhCN)₂](BF₄)₂ catalyzes the hetero-Diels-Alder reaction of dienes with aldehydes. The reaction yields substituted 5,6-dihydro-2H-pyrans without the use of Lewis acids and is believed to proceed through a stepwise mechanism.






References: (1) Oi, S.; Kashiwagi, K.; Terada, E.; Oshii, K.; Inoue, Y. *Tetrahedron Lett.* **1996**, 37, 6351. (2) Davies, J. A.; Hartley, F. R.; Murray, S. G. *J. Chem. Soc., Dalton Trans.* **1980**, 2246

For more information on the range of catalysts we offer, visit aldrich.com/catalysis

Air-Stable Precatalysts for Amination

C–N bond-forming cross-coupling reactions typically require a palladium source along with associated ligands. Most Pd(0) sources are not air-stable, while the commonly employed air-stable Pd(0) source, Pd₂(dba)₃, contains associated ligands which could impede the reaction in some cases. Stable Pd(II) precursors require reduction under the reaction conditions. In either case, a ligand must be added to the reaction in order to lead to the active Pd species. Buchwald and coworkers recently reported the use of highly active air- and moisture-stable precatalysts, which, under the standard reaction conditions, form the active monoligated Pd species. These precatalysts are exceptionally efficient even under challenging conditions, such as coupling electron-poor anilines with deactivated aryl chlorides. The catalyst precursors also offer other advantages including low catalyst loadings and short reaction times.




Reference: Biscoe, M. R.; Fors, B. P.; Buchwald, S. L. *J. Am. Chem. Soc.* **2008**, 130, 6686

For more information on our Buchwald catalysts and ligands, visit aldrich.com/buchwald

Nonproprietary Catalysts for Cross-Coupling Reactions


The cross-coupling reaction of heteroaryl halides is of particular interest to the pharmaceutical industry since many biologically active compounds are accessed through use of the Suzuki–Miyaura reaction. However, the efficient coupling of boronic acids with five-membered-ring heteroaryl halides or six-membered-ring heteroaryl chlorides bearing heteroatom substituents has not been well-developed. Catalysts are thought to form inactive complexes with many of these types of substrates, and thus, they typically require high catalyst loadings in order to achieve good yields. Guram's group at Amgen has recently reported the development of an air-stable palladium complex, $(A^t\text{Phos})_2\text{PdCl}_2$, for Suzuki–Miyaura cross-coupling reactions. The catalyst is very effective at coupling a wide range of substrates with arylboronic acids, including amino-substituted 2-chloropyridines and five-membered-ring heteroaryl halides. The products are obtained in excellent yields and high turnover numbers (up to 10,000 TON) are typically achieved. A series of new $\text{PdCl}_2\{\text{PR}_2(\text{Ph-R}')\}_2$ catalysts were developed with various reactivities.



entry	Product	Yield (%)
1		93
2		98
3		95

References: (1) Singer, R. A. et al. *Tetrahedron Lett* **2006**, 47, 3727. (2) Singer, R. A. et al. *Synthesis* **2003**, 1727. (3) Guram, A. S. et al. *Org. Lett.* **2006**, 8, 1787

For more information on the range of catalysts we offer, visit aldrich.com/catalysis



Palladium(II) [1,3-bis(diphenylphosphino)propane]bis(benzonitrile) bis(tetrafluoroborate)	
696617	250 mg
[175079-12-6]	
C ₂₄ H ₂₀ F ₄ N ₂ P ₂	
FW: 898.71	
(XPhos) palladium(II) phenethylamine chloride	
704954	250 mg
C ₂₁ H ₂₆ ClNPPd	1 g
FW: 738.76	
(SPhos) palladium(II) phenethylamine chloride (1:1 MTBE solvate)	
704946	250 mg
C ₂₁ H ₂₆ ClNPPd	1 g
FW: 760.72	

RuPhos palladium(II) phenethylamine chloride	
707589	250 mg
C ₂₁ H ₂₆ ClNPPd	
FW: 738.76	
t-BuXPhos palladium(II) phenethylamine chloride	
708739	250 mg
C ₂₁ H ₂₆ ClNPPd	
FW: 686.69	
Bis(di-tert-butyl(4-dimethylaminophenyl)phosphine) dichloropalladium(II)	
678740	1 g
[887919-35-9]	5 g
C ₄₈ H ₆₄ N ₂ P ₂ Cl ₂	
FW: 708.07	
Bis(dicyclohexyl)(4-dimethylaminophenyl)phosphine]palladium(II) chloride	
692913	250 mg
C ₄₀ H ₆₄ Cl ₂ N ₂ P ₂	1 g
FW: 812.22	
Bis(dicyclohexyl(4-dimethylaminophenyl)phosphine) palladium(II) chloride	
692921	250 mg
C ₃₆ H ₅₄ Cl ₂ F ₂ P ₂	1 g
FW: 757.93	



P.O. Box 14508
St. Louis, MO 63178
USA

Sigma-Aldrich® Worldwide Offices

Argentina

Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

Chile

Tel: (+56) 2 495 7395
Fax: (+56) 2 495 7396

China

Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

Tel: (+45) 43 56 59 00
Fax: (+45) 43 56 59 05

Finland

Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 88
Fax: (+33) 474 95 68 08

Germany

Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Hungary

Ingyenes telefonszám: 06 80 355 355
Ingyenes fax szám: 06 80 344 344
Tel: (+36) 1 235 9063
Fax: (+36) 1 269 6470

India

Telephone
Bangalore: (+91) 80 6621 9400
New Delhi: (+91) 11 4358 8000
Mumbai: (+91) 22 4087 2364
Hyderabad: (+91) 40 4015 5488
Kolkata: (+91) 33 4013 8000

Fax

Bangalore: (+91) 80 6621 9550
New Delhi: (+91) 11 4358 8001
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 4015 5466
Kolkata: (+91) 33 4013 8016

Ireland

Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 402 20370
Fax: (+353) 402 20375

Israel

Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

Free Tel: 800 827 018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: (+52) 722 276 1600
Fax: (+52) 722 276 1601

The Netherlands

Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

Tel: (+47) 23 17 60 00
Fax: (+47) 23 17 60 10

Poland

Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

Tel: (+7) 495 621 5828
Fax: (+7) 495 621 6037

Singapore

Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

Slovakia

Tel: (+421) 255 571 562
Fax: (+421) 255 571 564

South Africa

Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313

United States

Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Vietnam

Tel: (+84) 3516 2810
Fax: (+84) 6258 4238

Internet

sigma-aldrich.com

*Accelerating Customers'
Success through Innovation and
Leadership in Life Science,
High Technology and Service*

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052
Technical Service (800) 325-5832 • sigma-aldrich.com/techservice
Development/Custom Manufacturing Inquiries **SAFC**® (800) 244-1173
Safety-related Information sigma-aldrich.com/safetycenter

World Headquarters
3050 Spruce St.
St. Louis, MO 63103
(314) 771-5765
sigma-aldrich.com

© 2004 Sigma-Aldrich. All rights reserved. SIGMA-ALDRICH, ALDRICH, and FSC are trademarks of Sigma-Aldrich. All other trademarks are the property of their respective owners. Sigma-Aldrich and its affiliate Sigma-Aldrich are not responsible for the content of any website or other information source. The information contained in this and other Sigma-Aldrich publications is for informational purposes only. The purchaser must determine the suitability of the products for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.

MF
73542 50845
1050



Mixed Sources

Product group from well-managed
forests, controlled sources and
recycled wood or fiber.
www.fsc.org Cert no. SGS-COC-003338
© 1996 Forest Stewardship Council



SIGMA-ALDRICH

THE COMING OF AGE OF GOLD CATALYSIS

Aldrichimica ACTA

VOL. 43, NO. 2 • 2010

 **ALDRICH**
Chemistry




Gold-Catalyzed Addition of X-H Bonds to C-C Multiple Bonds

Copper-, Silver-, and Gold-Catalyzed Migratory Cycloisomerizations

Leading to Heterocyclic Five-Membered Rings

Enantioselective Synthesis Based on Catalysis by Chiral Oxazaborolidinium Cations

SIGMA-ALDRICH



To advance your chemistry research

Aldrich ChemFiles 10.3 is available now!

Aldrich ChemFiles is a complimentary quarterly publication written by our experts from R&D and product management. It's specifically for research chemists who are looking to implement the latest innovative synthetic strategies using Aldrich Chemistry's cutting-edge products from:

- Building Blocks
- Synthetic Reagents
- Catalysis and Inorganic Reagents
- Chemical Biology
- Labware
- Asymmetric Synthesis
- Stockroom Reagents
- Organometallic Reagents
- Specialty Synthesis
- Chemistry Services

Each edition highlights common and useful applications of our new products, and shows where these and related products can be quickly found on the Aldrich Chemistry website.

Aldrich ChemFiles helps you keep informed of the new Aldrich Chemistry products that facilitate the latest research methodologies and trends, and allows you to access key starting materials and reagents more efficiently.

Advancing Chemistry. Aldrich ChemFiles.

ALDRICH
Chemistry



Aldrich ChemFiles
Newsletter Volume 10,
Number 3, available
now by visiting
aldrich.com/chemfiles

Aldrichimica ACTA

VOL. 43, NO. 2 • 2010

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries	800-325-3010
Technical Service	800-231-8327
SAFC®	800-244-1173
Custom Synthesis	800-244-1173
Flavors & Fragrances	800-227-4563
International	414-438-3850
24-Hour Emergency	414-438-3850
Website	sigma-aldrich.com
Email	aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 2988, Milwaukee, WI 53201, USA

Subscriptions

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 2988
Milwaukee, WI 53201-2988
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com/acta.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2010 Sigma-Aldrich Co.

"PLEASE BOTHER US."



Joe Porwoll, President
Aldrich Chemical Co., Inc

Professor Bruce Lipshutz from the University of California, Santa Barbara, kindly suggested that we make copper-in-charcoal. This useful heterogeneous catalyst promotes a variety of catalytic transformations such as asymmetric hydrosilylations, click reactions, and is also used in the preparation of diaryl ethers

(1) Lipshutz, B. H.; Frieman, B. A.; Tomaso, A. E., Jr. *Angew. Chem., Int. Ed.* **2006**, *45*, 1259
(2) Lipshutz, B. H.; Taft, B. R. *Angew. Chem., Int. Ed.* **2006**, *45*, 8235. (3) Lipshutz, B. H.; Unger, J. B.; Taft, B. R. *Org. Lett.* **2007**, *9*, 1089

Cu/C
709107

709107 Copper-in-charcoal, 3 wt. % 5 g \$71.80

Naturally, we made this useful heterogeneous catalyst. It was no bother at all, just a pleasure to be able to help

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestions. We're delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the back cover

TABLE OF CONTENTS

Gold-Catalyzed Addition of X-H Bonds to C-C Multiple Bonds	27
A. Stephen K. Hashmi* and Miriam Bührle, Ruprecht Karls-Universität Heidelberg	
Copper-, Silver-, and Gold-Catalyzed Migratory Cycloisomerizations Leading to Heterocyclic Five-Membered Rings	37
Alexander S. Dudnik, Natalia Chernvak, and J. G. de Vries, University of Illinois at Chicago	
Enantioselective Synthesis Based on Catalysis by Chiral Oxazaborolidinium Cations	49
Santanu Mukherjee, Indian Institute of Science, and E. J. Corey,* Harvard University	

ABOUT OUR COVER

Detail from *Travellers near the Edge of a Wood* (oil on panel, 51.5 × 67.2 cm) was painted ca. 1638 by Jacob van Geel (ca.1585–1638), the seventeenth century Dutch painter. Born in Middelburg, Zeeland Province, he became a member of the guild of St. Luke (Middelburg) at about age 30, and later joined the painters' guilds of Delft and Dordrecht, and seems to have led a private life dominated by debt and conflict

Having been influenced by the painting styles of Jan Brueghel the Elder (1568–1625) and Hercules Seghers (ca. 1589–1638), he specialized in painting fantasy landscapes often dominated by a large cluster of trees, as is the case here. In this complex composition, the viewer is drawn to at least three scenes featuring travellers at various stages of their journey. The bright daylight penetrates from the left side and illuminates the twisted and bulging trees that fill the composition and attract the viewer's eyes to the center of this work




This painting is part of the Bader Collection of Dutch and Flemish Paintings, whose future home will be the Agnes Etherington Art Centre of Queen's University, Kingston, ON, Canada.




Gold Catalysts from Sigma-Aldrich


Gold Catalysts




697575




254037




715050




679771




404217




288225




420727




696277




687510



702749



677922



298174




398470




481130



334049



677876



665142

Silver Salt Additives



668001



176427



227730



208361




227722



176435

For more information, please visit aldrich.com/gold



Heterocyclic MIDA Boronates from Aldrich

Aldrich has prepared a series of 2-, 3-, and 4-heterocyclic MIDA boronates, including the notoriously unstable 2-pyridyl boron anion equivalent 2-pyridinylboronic acid MIDA ester (**719390**).

(1) Gillis, E. P.; Burke, M. D. *Aldrichimica Acta* **2009**, 42, 17. (2) Knapp, D. M.; Gillis, E. P.; Burke, M. D. *J. Am. Chem. Soc.* **2009**, 131, 6961.

Benefits of MIDA Boronates:

- Bench- and air-stable
- Stable to anhydrous cross-couplings conditions
- Good solubility properties
- Monomeric composition
- Silica gel compatible
- Controlled release for in situ cross-coupling
- Straightforward surrogate to unstable boronates
- Simple, IP-included model





Figure 1. 2-Pyridinylboronic acid MIDA boronate, stable 2-pyridyl boron anion equivalent.

Aldrich has prepared a growing series of heterocyclic MIDA boronates, many of which are shown below.



For more information on MIDA boronates, please visit aldrich.com/mida



P.O. Box 14508
St. Louis, MO 63178
USA

PRESORTED
STANDARD
U.S. POSTAGE
PAID
SIGMA-ALDRICH
CORPORATION

Change Service Requested

Sigma-Aldrich® Worldwide Offices

Argentina

Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

Chile

Tel: (+56) 2 495 7395
Fax: (+56) 2 495 7396

China

Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

Tel: (+45) 43 56 59 00
Fax: (+45) 43 56 59 05

Finland

Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 88
Fax: (+33) 474 95 68 08

Germany

Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Hungary

Ingyenes telefonszám: 06 80 355 355
Ingyenes fax szám: 06 80 344 344
Tel: (+36) 1 235 9063
Fax: (+36) 1 269 6470

India

Telephone
Bangalore: (+91) 80 6621 9400
New Delhi: (+91) 11 4358 8000
Mumbai: (+91) 22 4087 2364
Hyderabad: (+91) 40 4015 5488
Kolkata: (+91) 33 4013 8000

Fax

Bangalore: (+91) 80 6621 9550
New Delhi: (+91) 11 4358 8001
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 4015 5466
Kolkata: (+91) 33 4013 8016

Ireland

Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 402 20370
Fax: (+353) 402 20375

Israel

Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

Free Tel: 800 827 018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

Free Tel: 800 202 180
Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: (+52) 722 276 1600
Fax: (+52) 722 276 1601

The Netherlands

Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

Tel: (+47) 23 17 60 00
Fax: (+47) 23 17 60 10

Poland

Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

Tel: (+7) 495 621 5828
Fax: (+7) 495 621 6037

Singapore

Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

Slovakia

Tel: (+421) 255 571 562
Fax: (+421) 255 571 564

South Africa

Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313

United States

Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Vietnam

Tel: (+84) 3516 2810
Fax: (+84) 6258 4238

Internet

sigma-aldrich.com

Enabling Science to
Improve the Quality of Life

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052
Technical Service (800) 325-5832 • sigma-aldrich.com/techservice
Development/Custom Manufacturing Inquiries **SAFC**® (800) 244-1173
Safety-related Information sigma-aldrich.com/safetycenter

World Headquarters
3050 Spruce St.
St. Louis, MO 63103
(314) 771-5765
sigma-aldrich.com



Mixed Sources

Product group from well-managed
forests, controlled sources and
recycled wood or fiber
www.fsc.org Cert no. SCS-COC-003338
© 1996 Forest Stewardship Council



SIGMA-ALDRICH

DEDICATED TO DR. JAI P. NAGARKATTI

Aldrichimica ACTA

VOL. 43, NO. 3 • 2010



Remembering Dr. Jai P. Nagarkatti (1947–2010)

Recent Advances in the Application of α -Phenylethylamine (α -PEA)
in the Preparation of Enantiopure Compounds

2-Iodoxybenzenesulfonic Acid (IBS) Catalyzed Oxidation of Alcohols

SIGMA-ALDRICH

Aldrich Congratulates the 2010 Winners of the Nobel Prize in Chemistry

Professors **Ei-ichi Negishi**, **Akira Suzuki**, and **Richard F. Heck** have been awarded the 2010 Nobel Prize in Chemistry by The Royal Swedish Academy of Sciences for their trailblazing contributions to the area of "palladium-catalyzed cross couplings in organic synthesis".



Photo courtesy of Purdue U/
Andrew Hancock

Distinguished Professor
Ei-ichi Negishi
Purdue University, U.S.A.



Photo courtesy of Hokkaido U

Distinguished Professor Emeritus
Akira Suzuki
Hokkaido University, Japan



Photo courtesy of U of Delaware

Professor Emeritus
Richard F. Heck
University of Delaware, U.S.A

Aldrich congratulates these distinguished chemistry pioneers on this achievement and thanks them for their lasting contributions to organic synthesis. Our company is proud to have had close collaborations with Professors Negishi and Suzuki, two former associates of another Nobel laureate and longtime Aldrich collaborator and member of the Board of Directors, the late Professor H. C. Brown.

The *Aldrichimica Acta* has had the distinct privilege of publishing high-impact review articles by a number of chemistry Nobel laureates, such as Professors Derek H. R. Barton, Herbert C. Brown, Elias James Corey, Ei-ichi Negishi, George A. Olah, Charles J. Pedersen, Vladimir Prelog, and K. Barry Sharpless.

The tradition of supporting future scientific leaders is strong at Sigma-Aldrich. Throughout our history, we have acknowledged scientific excellence through sponsorship of awards, symposia, and graduate-level research in the fields of chemistry, life science, and materials science. A number of previous winners of Sigma-Aldrich sponsored American Chemical Society awards (e.g., Professors George A. Olah, K. Barry Sharpless, and Ei-ichi Negishi) have gone on to receive the ultimate recognition, the Nobel Prize in Chemistry.

To access highly informative *Aldrichimica Acta* reviews written by Nobel laureates and other chemistry luminaries, visit aldrich.com/acta

Aldrichimica ACTA

VOL. 43, NO. 3 • 2010

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-231-8327
SAFC® 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 414-438-3850
24-Hour Emergency 414-438-3850
Website sigma-aldrich.com
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 2988, Milwaukee, WI 53201, USA

Subscriptions

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 2988
Milwaukee, WI 53201-2988
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com/acta.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2010 Sigma-Aldrich Co.

"PLEASE BOTHER US."




Joe Porwoll

Joe Porwoll, President
Aldrich Chemical Co., Inc.

A number of customers have suggested that we produce 2-pyridinylboronic acid MIDA ester. This isolable, chromatographically and benchtop-stable reagent serves as a viable 2-pyridylboron anion equivalent, and is used as a versatile cross-coupling partner.

(1) Knapp, D. M.; Gillis, E. P.; Burke, M. D. *J. Am. Chem. Soc.* **2009**, *131*, 6961



719390 2-Pyridinylboronic acid MIDA ester 1 g
5 g

Naturally, we made this useful boronic acid surrogate. It was no bother at all, just a pleasure to be able to help

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the back cover.

TABLE OF CONTENTS

Remembering Dr. Jai P. Nagarkatti (1947–2010)	63
<i>Sharbil J. Firsan, Sigma-Aldrich Corp</i>	
Recent Advances in the Application of α-Phenylethylamine (α-PEA) in the Preparation of Enantiopure Compounds	65
<i>Yamir Bandala and Eusebio Juaristi, Centro de Investigación y de Estudios Avanzados de Instituto Politécnico Nacional</i>	
2-Iodoxybenzenesulfonic Acid (IBS) Catalyzed Oxidation of Alcohols	83
<i>Muhammet Uyanik and Kazuaki Ishihara,* Nagoya University</i>	

ABOUT OUR COVER

Taos Pueblo Snow (oil on canvas, 40.6 x 70.0 cm) was painted in 2009 by the American painter and sculptor Rosie Sandifer (b. 1946) following a visit to the Pueblo. It depicts the centuries-old, Native American village with the same name. Taos Pueblo is located about two miles north of the town of Taos, New Mexico, and is one of the longest continually inhabited places in the U.S. It has been designated a World Heritage Site by UNESCO and a National Historic Landmark by the U.S. Government.



Photograph © Rosie Sandifer.

Equally gifted at painting and sculpting, Sandifer paints her impressions of what she observes in nature by simplifying to the essentials the effects of fleeting light on the subject, as evidenced by the white, brown, and turquoise that dominate in this painting. While her subject matter has included landscapes, figures, and animals, she has tended to focus on Western landscapes. She has participated in numerous solo, group, and juried exhibitions, and her paintings and sculptures are on display in a number of museums and public places in the U.S. Currently a resident of Santa Fe, NM, Sandifer received her extensive art education and training most notably at the Froman Painting School in Cloudcroft, NM, and the Art Students League in Stowe, VT, where she was influenced by Frank Mason, a classical realist painter and one of the most acclaimed modern American painters and teachers.

This painting is provided courtesy of the artist and is in her private collection.




New Oxetane Building Blocks

Oxetanes have historically been of modest interest to synthetic and medicinal chemists, perhaps with the natural product paclitaxel or TAXOL® being the best known example of an oxetane-containing substance. Presently, oxetanes are receiving greater attention as attractive modules for drug discovery, largely due to a series of reports from Rogers-Evans, Carreira, and coworkers. These reports have demonstrated the improved physico- and biochemical properties of a molecular scaffold when an oxetane unit replaces a *gem*-dimethyl unit¹ and the ability for an oxetane ring to function as a surrogate for a carbonyl group.^{2,1b} Another recent report has disclosed the use of 1,6-substituted azaspiro[3.3]heptanes containing an oxetane ring as alternatives to unstable 1,3-heteroatom-substituted cyclohexanes.³ In most cases, 3-oxetanone, **731536**, was the principal building block employed by the authors to install the oxetane unit (Scheme 1).


References: (1) (a) Wuitschik, G. et al. *Angew. Chem., Int Ed* **2006**, 45, 1736. (b) Wuitschik, G. et al. *J. Med. Chem.* **2010**, 53, 3227. (2) Wuitschik, G. et al. *Angew. Chem., Int Ed* **2008**, 47, 4512. (3) Burkhard, J. A. et al. *Org. Lett.* **2010**, 12, 1944

For a complete list of available oxetanes, please visit aldrich.com/oxetane



Scheme 1. Oxetane Derivatives Synthesized from 3-Oxetanone.

We are pleased to offer a wide selection of new oxetane building blocks for a variety of applications in synthetic and medicinal chemistry.



TAXOL is a registered trademark of Bristol-Myers Squibb Co

Remembering Dr. Jai P. Nagarkatti (1947–2010)

Dr. Sharbil J. Firsan¹
 Editor, *Aldrichimica Acta*
 Sigma-Aldrich Corp.


The *Aldrichimica Acta* has lost a dear friend, a lifelong staunch supporter, and a strong believer in its value and mission, both to the company and to the chemistry community. Dr. Jai P. Nagarkatti—the company's Chairman, CEO, and President—died suddenly on Saturday, November 13, 2010.^{2–4} While his Sigma-Aldrich “family” has been devastated by his unexpected departure and mourn his loss, this editorial will not be about the facts and figures of Jai's career and accomplishments, nor about Sigma-Aldrich's growth under his tenure. Suffice it to say that he gave over 34 years of dedicated service to Aldrich and Sigma-Aldrich, the only employer he ever worked for after receiving his Ph.D. degree from Texas A&M University, Commerce, TX. Not that facts and figures are not important, but the intention here is to offer fond recollections that shed some light on the kind of man Jai really was.

It was during my first interview with the company that I met Jai. He was the president of Aldrich at the time, and came across as energetic, eager, engaged, and passionate about Aldrich and the business and science of chemistry. He was keen to find out from me what new ideas I was bringing in and what I, as a customer, thought of Aldrich. He took notes during our meeting. During the interview, I could tell he was very proud of working for Aldrich.

After joining the company, I became well acquainted with Jai, who consistently showed strong interest in the *Acta*: He eagerly anticipated each issue, kept copies in his office, and took a strong interest in what was published and what was advertised in it. He recommended authors, emailed thank-you letters to them, and never missed a chance to offer this editor sincere words of appreciation and encouragement. In fact, one little story sums him up in this regard: I sometimes would be crossing St. Paul Avenue in Milwaukee going from one Aldrich building to another, while Jai would be crossing the same street in the opposite direction. He, invariably, would stop me in the middle of the street to ask me about the *Acta*—totally oblivious to fast-moving traffic and the fact that both of us could be run over by an inattentive driver.

I wholeheartedly concur with one of Jai's longtime associates who related to me that he was always amazed at how eminently reachable Jai was, even after becoming the Chairman, President, and CEO of Sigma-Aldrich, a \$2+ billion company, considering how incredibly busy and in demand he was. While he had administrative assistants, Jai never walled himself off; one could walk into his office anytime. He also did not hesitate to answer a phone call or an email himself, or walk into any area of the company and chat with anyone. When he wanted to talk to a subordinate, he preferred to walk over to his/her cube or office, rather than call him/her to come to his office. He did not hesitate,

as this editor can attest to, to pick up the phone and call someone in the company to thank them for their efforts and offer words of encouragement. He was incredibly humble and had a knack for conversing with equal ease with a custodian at Sigma-Aldrich, a Wall Street analyst, or a highly regarded university professor. He made sure to attend and serve at as many of the company functions for employees as he could. He almost was a fixture at our Company's annual picnic and other events. That was his way of staying connected to the employees. As Dr. Alfred Bader put it, “[Jai] was inspirational for employees who appreciated his dedication and warmth. Jai was so very human...”



Jai P. Nagarkatti in 2007 at Sigma-Aldrich's Headquarters in St. Louis, MO.

By keeping a full schedule and working hard to the point of being a workaholic, Jai led by example and brought an excitement, and a depth of perception, to projects and tasks that went beyond the customary business reasons and logic: If you met with him to discuss a project, you learned not only about the project, but also about the deeper significance of the project to the company, to customers, and even to society. This, of course, meant that you also left his office with a lot more work to do than when you walked in.

Being focused on his work and responsibilities did not mean that he lacked a sense of humor. In fact, it wasn't unusual to overhear him being lighthearted. One noteworthy moment that I remember occurred at an Aldrich employee award dinner in Milwaukee, when he was introducing a longtime, distinguished associate, who was well-known for being outspoken. Jai said that, “when [name omitted] reported to me, I was never sure who reported to whom.” Being as humble as he was, Jai was not afraid of making fun of his weaknesses, as when he picked up a set of golf clubs at a company outing and played a round of golf, never having swung a golf club before! Jai just wanted employees to see that he wasn't good at everything. ☺

Farewell my friend, Jai. Your enthusiasm, dedication, thoughtfulness, humility, and kindness will be forever missed.

Notes


- (1) I thank John Radke, Joe Porwoll, and Ali Ataei for sharing their recollections of Jai. Some of these recollections have been weaved into this article. I am also grateful to Ms. Melissa Jacobs of Sigma-Aldrich for reviewing this article and offering helpful suggestions.
- (2) A summary of Jai's career at Sigma-Aldrich is part of a press release issued by the company on November 14, 2010, announcing Jai's sudden death and the appointment of a successor. It can be accessed at [http://](http://investor.sigmaaldrich.com/releasedetail.cfm?ReleaseID=530331)

investor.sigmaaldrich.com/releasedetail.cfm?ReleaseID=530331.

- (3) The obituary released by Jai's family was published in the St. Louis Post-Dispatch of Monday, November 15, 2010, which can be accessed at <http://www.legacy.com/obituaries/stltoday/obituary.aspx?n=jai-nagarkatti&pid=146642609>.
- (4) The St. Louis Post-Dispatch ran a story on Jai by Jim Doyle on Tuesday, November 16, 2010, which can be accessed at http://www.stltoday.com/business/article_46d399bc-61e2-521a-9ff8-b7048bf6ab3h.html.



New Reductive Amination Reagent

PEMB: A Stable, Liquid Pyridine-Borane Complex



5-Ethyl-2-methylpyridine borane (PEMB)
725080











Select Substrate Scope



Advantageous Properties of PEMB:

- Excellent for reductive aminations
- Mild reducing agent for imines and oximes
- Reaction with protic solvents is very slow
- Soluble in aromatic hydrocarbons, alcohols, and ether solvents
- Can be used solvent-free for reductive aminations
- Chemically efficient: two of three hydrides are utilized

Examples

Aldehyde and Amine	Conditions	Product	% Yield in MeOH (% Yield Neat)
 PhNH ₂	$\xrightarrow[\text{MeOH, 25}^\circ\text{C}]{\text{PEMB, AcOH}}$		72 (80)
 Pr ₂ NH	$\xrightarrow[\text{MeOH, 25}^\circ\text{C}]{\text{PEMB, AcOH}}$		0 (96)
 PhNH ₂	$\xrightarrow[\text{MeOH, 25}^\circ\text{C}]{\text{PEMB}}$		92 (94)
 PhNH ₂	$\xrightarrow[\text{MeOH, 25}^\circ\text{C}]{\text{PEMB, AcOH}}$		92 (93)
 PhNH ₂	$\xrightarrow[\text{MeOH, 50}^\circ\text{C}]{\text{PEMB, AcOH}}$		74 (94)


For more examples and experimental details, please see: Burkhardt, E. R.; Coleridge, B. M. *Tetrahedron Lett.* **2008**, 49, 5152

For more information, please visit aldrich.com/pemb

For research quantities under 500 grams, please contact Aldrich Chemistry.

For quantities over 500 grams, please contact BASF.

Recent Advances in the Application of α -Phenylethylamine (α -PEA) in the Preparation of Enantiopure Compounds



Dr. Yamir Bandala



Professor Eusebio Juaristi

Yamir Bandala and Eusebio Juaristi
 Departamento de Química
 Centro de Investigación y de Estudios Avanzados
 del Instituto Politécnico Nacional
 Apartado Postal 14-740
 07000 México, D. F., México
 Email: vyamir@relaq.mx, juaristi@relaq.mx

Outline

1. Introduction
2. Incorporation of α -Phenylethylamine (α -PEA) in Chiral Ligands for Asymmetric Catalysis
 - 2.1. Enantioselective Reduction of Prochiral Unsaturated Double Bonds
 - 2.2. Enantioselective Addition of Diethylzinc to Aromatic Aldehydes
 - 2.3. α -PEA Derivatives as Ligands in Miscellaneous Enantioselective Reactions
3. α -Phenylethylamine as Resolving Agent
 - 3.1. Resolution via Diastereomeric Salt Formation
 - 3.2. Resolution via Covalent Attachment of α -PEA to Racemic Substrates
 - 3.3. Chiral Solvating Agents Containing the α -PEA Fragment
4. α -Phenylethylamine as Chiral Auxiliary
 - 4.1. Diastereoselective Reactions of Chiral Imines Incorporating the α -Phenylethyl Moiety
 - 4.2. Diastereoselective Reactions of Optically Active Enamines
 - 4.3. Diastereoselective Alkylation of α -PEA-Containing Enolates
 - 4.4. α -PEA in Diastereoselective Cyclization Reactions
 - 4.5. α -PEA as Chiral Auxiliary in Miscellaneous Procedures
5. α -Phenylethylamine as Chiral Reagent in the Stereodifferentiation of Prochiral Substrates
6. α -Phenylethylamine in Organocatalysts
7. Concluding Remarks
8. Acknowledgements
9. References and Notes

1. Introduction

Mostly as a consequence of the increased interest in the production of enantiomerically pure compounds by the pharmaceutical and agrochemical industries, the chemistry community has intensified its efforts to develop more efficient ways to obtain chiral products with the desired configuration. In this regard, α -phenylethylamine (α -PEA, 1-phenylethylamine, 1-phenylethylamine, α -methylbenzenemethanamine, α -methylbenzylamine, α -aminoethylbenzene, $\text{CH}_3\text{CH}(\text{Ph})\text{NH}_2$) has emerged as an

important chiral reagent principally because of its low price, accessibility in both enantiomeric forms, and facile introduction and removal in most asymmetric syntheses. Two reviews appeared in the late 1990s on the use of α -PEA in the preparation of enantioenriched compounds.¹ However, a significant number of new developments have been reported in the last decade, which warrants an update on the more significant advances in the area that were reported between 2000 and early 2010. This review focuses on the applications of α -PEA as a chiral ligand for asymmetric catalysis, as a chiral auxiliary in the resolution of racemates, and as a chiral reagent in the stereodifferentiation of prochiral substrates, as well as on its incorporation into novel organocatalysts.


2. Incorporation of α -Phenylethylamine (α -PEA) in Chiral Ligands for Asymmetric Catalysis

A central development in asymmetric catalysis has been the design and preparation of enantiopure catalysts containing an active metal and a regulating chiral organic ligand and their application in the synthesis of enantioenriched building blocks from prochiral substrates. Following the great success of complexes of chiral phosphine ligands and transition metals in asymmetric synthesis,² intense research activity has been directed towards the development of comparable ligands containing nonracemic amines. In particular, catalysts incorporating ligands with the α -PEA moiety have been quite effective in several asymmetric transformations.


2.1. Enantioselective Reduction of Prochiral Unsaturated Double Bonds

Hu, Zheng, and co-workers³ have described the use of chiral phosphine-aminophosphine ligand **1** (easily synthesized from (*S*)- α -PEA, *n*-BuLi, and Ph_2PCI), in combination with $\text{Rh}(\text{cod})_2\text{BF}_4$, for the enantioselective hydrogenation of the C–C double bond in various α -enol ester phosphonates and α -enamido phosphonates. In particular, ligand **1** was utilized in the successful preparation of biologically important α -hydroxy and α -amino phosphonates in 94–99% yields and 93–97% ee's (**Scheme 1**). In the same way, optically active ligand **2** was developed by Brauer et al.⁴ for the enantioselective hydrogenation of methyl acetamidocinnamate to the corresponding amido ester with 95% ee (see **Scheme 1**).


In a related development, Knochel and collaborators reported on the effectiveness of the nonracemic, ferrocene-based ligand **3** in



Scheme 1. Highly Enantioselective Hydrogenations of C=C Double Bonds in the Presence of Chiral Ligands Incorporating α -PEA.



Scheme 2. Chiral Ferrocene Derivative **3** in the Ru- and Rh-Catalyzed Hydrogenation of C=C and C=O Bonds. (Ref. 5)



Scheme 3. Highly Enantioselective Hydrogenations Catalyzed by 4-Rh(I) Complexes. (Ref. 6)


the enantioselective ruthenium- and rhodium-catalyzed hydrogen addition to prochiral substrates containing C=C or C=O bonds (**Scheme 2**).⁵

Similarly, Huang et al. prepared an interesting phosphine-phosphoramidite ligand, **4**, starting from (*S*)- α -PEA. Together with Rh(I), **4** efficiently catalyzes the reduction of α -dehydroamino esters, enamines, and dimethyl itaconate to afford the reduced products in $\geq 99\%$ ee's. These results show that the central chirality of the α -phenylethylamine moiety in **4** dictates the absolute configuration of the hydrogenation product, regardless of the axial chirality of the binaphthyl moiety (**Scheme 3**).⁶


The asymmetric transfer hydrogenation of prochiral ketones can be carried out successfully by employing α -PEA derivatives. For example, the enantioselective hydrosilylation and reduction of prochiral ketones to alcohols has been effected by using polymethylhydrosiloxane (PMHS) as the hydride source and a set of *all-S*, α -PEA-containing, chiral diamines, **5**, and thioureas, **6**. This procedure afforded the respective secondary *R* alcohol (**eq 1**).⁷ Along similar lines, Adolfsson and co-workers showed that the use of α -amino-acid-containing chiral thioamide ligands, *i*-PrOH as a hydride source, and Ru(II)- or Rh(III)-based catalyst complexes, led to several secondary alcohols derived from acetophenones in high yields and good enantioselectivities.⁸

2.2. Enantioselective Addition of Diethylzinc to Aromatic Aldehydes

Asymmetric organozinc addition to prochiral aldehydes and ketones allows the synthesis of chiral alcohols, which are ubiquitous features in the structures of natural products and pharmaceuticals. The reaction between diethylzinc (Et_2Zn) and benzaldehyde has become a benchmark reaction for testing newly designed ligands for catalytic enantioselective synthesis.⁹ In this context, Wilken et al. have reported the use of several nonracemic tridentate azetidines in the enantiocontrolled catalytic addition of diethylzinc to several aromatic aldehydes, and demonstrated that the *all-S* chiral ligand, **7**, efficiently catalyzes this reaction in good yields and ee's (**eq 2**).¹⁰ In a related development, Guofu Zi and collaborators prepared



eq 1 (Ref. 7)



eq 2 (Ref. 10)

various chiral azetidines by reaction of (*R*)- or (*S*)- α -PEA, benzaldehyde, and acetoxyacetyl chloride. These researchers observed that the resulting *R,R,R* azetidine ligand led to the *S* alcohols, while the *R,S,S* diastereomer gave rise to the *R* alcohols in 73–94% yields and 57–95% ee's.¹¹

Moreover, Wang's group¹² has described the application of optically active ligand **8**—prepared by a one-pot condensation of the corresponding aldehyde with 2-naphthol and (*S*)- α -PEA—in the enantioselective addition of diethylzinc to aldehydes. Most interesting are the high ee values obtained at room temperature (>99% ee's for the *R* alcohols). Furthermore, Mikami and collaborators¹³ have demonstrated that ligand **9** efficiently catalyzes diethylzinc addition to several aromatic aldehydes in good yields and high ee's (eq 3).

Based on the efficient opening of cyclohexene oxide by α -PEA,¹⁴ Mastranzo et al. prepared a series of chiral ligands derived from *trans*- β -aminocyclohexanol. The versatility of this approach resulted in the preparation of a large number of optically active ligands that were successfully tested in the benchmark reaction of Et₂Zn with benzaldehyde (eq 4).¹⁵

2.3. α -PEA Derivatives as Ligands in Miscellaneous Enantioselective Reactions



Recently, Isleyen and Dogan described the use of chiral ferrocenyl aziridinylmethanol **14**, which is prepared in three steps from an acryloylferrocene precursor and (*S*)- α -PEA. In particular, these authors showed that **14** can be employed in the catalytic, enantioselective conjugate addition of diethylzinc to enones to give β -ethylated ketones in up to 80% ee's, with the *R* enantiomer being favored in the majority of cases (eq 5).¹⁶

Another case of successful application of α -PEA derivatives as ligands in asymmetric synthesis is the palladium-catalyzed enantioselective allylic alkylation with dimethyl malonate in the presence of chiral ligands **15** and **16**. Ding and co-workers prepared chiral ligand **15** by the asymmetric aminoalkylation of 2-naphthol with (*R*)- α -PEA and benzaldehyde. Ligand **15** efficiently catalyzes the asymmetric substitution of the acetate group in 1,3-diphenylprop-2-en-1-yl substrates in 41–96% yields and with ee's as high as 70% in favor of the *S* product (Scheme 4).^{17a} A few years later, Huang et al. employed a non-racemic imine, **16**, synthesized from (*S*)- α -PEA, to catalyze a similar reaction, achieving better yields and enantioselectivities (see Scheme 4).^{17b} Furthermore, Tsogoeva and collaborators have recently obtained modest enantioselectivities in the allylation of aldimines and the reduction of ketimines with trichlorosilane by utilizing chiral proline–formamide derivatives.¹⁸



In a very recent elegant work, Alonso et al. examined a set of chiral phosphoramidite–gold(I) complexes, incorporating the bis(α -PEA) moiety, as catalysts in the intramolecular, enantioselective [4 + 2] cycloaddition reaction of allenes and dienes. In one example, Au–phosphoramidite complex **17** catalyzed the enantioselective allene and diene [4 + 2] cycloaddition over the [4 + 3] cycloaddition, giving rise to enantiomeric ratios greater than 95:5 (eq 6).¹⁹

Chiral, cyclic hydroxynaphthylphosphonodiamide **18**—synthesized from (*S*)- α -PEA, 1,2-dibromoethane, and *O*-aryl phosphorodichloridates—has been employed as a chiral ligand in the asymmetric cyanation of aromatic aldehydes catalyzed by Ti(*i*-PrO)₄. The corresponding cyanohydrins were obtained in good yields and ee's (eq 7).²⁰



A remarkable application of a chiral ligand containing α -PEA was reported by Mikami's group, who prepared nonracemic ligand **19** from bis(*S*)- α -PEA, PCl₃, and 2,2'-methylenebis(4-






eq 3




eq 4 (Ref 14)




eq 5




Scheme 4. Enantioselective Palladium-Catalyzed Allylic Alkylation of Allyl Acetates with Dimethyl Malonate in the Presence of **15** and **16**. (Ref 17)




eq 6 (Ref. 19)



eq 7 (Ref. 20)



Scheme 5. Asymmetric Cu-Catalyzed Conjugate Additions to Nitroalkenes and Nitroacrylates. (Ref. 21)

Scheme 6. Asymmetric Oxidations Catalyzed by Complexes of Vanadium and Chiral Ligands Incorporating α -PEA. (Ref. 21)

methylphenol) for use in copper-catalyzed conjugate additions of organometallics to nitroalkenes and nitroacrylates. In the case of copper-catalyzed addition of diethylzinc to nitroalkenes, metal complexes containing ligand **19** showed high catalytic activity and very good enantioinduction (>99%, 91-99% ee). Similarly, the conjugate addition of trimethylaluminum to nitroacrylate proceeded in high yield (>99%) and enantioselectivity (93% ee) (Scheme 5).²¹ The product of this latter reaction was hydrogenated with palladium-on-charcoal to give the corresponding β^2 -alanine ethyl ester.²¹

On the other hand, the asymmetric oxidation of prochiral olefins and thioethers, catalyzed by chiral vanadium complexes, has emerged as an important alternative to other standard oxidation procedures. Traber et al.²² have described the application of a series of chiral hydroxamic acids as ligands for the vanadium-catalyzed asymmetric epoxidation of allylic alcohols. A noteworthy example is the preparation of the *S,S* epoxy alcohol in 77% yield and 95% ee by the asymmetric epoxidation of 2,3-diphenylallyl alcohol catalyzed by a vanadium complex with chiral ligand **20** (Scheme 6). In the same context, Rehder and co-workers²³ evaluated the catalytic properties of a series of chiral oxovanadium(V) complexes containing a chiral tridentate amino-bis(alcoholate) in the asymmetric oxidation of prochiral sulfides by organic hydroperoxides. (*S*)-Methyl *p*-tolyl sulfoxide was obtained with modest ee in the presence of catalyst **21** and cumyl hydroperoxide (CHP) as the oxidant (see Scheme 6).

Millward et al.²⁴ employed an optically active cyclopentadienylamido titanocene complex in the enantioselective ethylaluminum of allylbenzene and styrene, which proceeded with modest yields and enantioselectivities. Bianchini and Lee²⁵ carried out the cyclopropanation of styrene with ethyl diazoacetate catalyzed by a chiral ruthenium 2,6-bis(imino)pyridyl complex in good yield but modest enantioselectivity.

3. α -Phenylethylamine as a Resolving Agent

The resolution of racemates through the formation of ionic or covalent diastereomeric derivatives is widely used in academic research as well as in industry, in particular for the manufacture of pharmaceuticals on an industrial scale.

3.1. Resolution via Diastereomeric Salt Formation

A variety of chiral racemic acids—mandelic,²⁶ organophosphorus,²⁷ acetic and benzoic acid derivatives,^{28,29} and others³⁰—have been resolved by separation of the diastereomeric salts that they form with α -PEA. Ma and co-workers³¹ developed an elegant synthesis of enantiopure β -halobutenolides (*R*)- and (*S*)-**24** that relies on the resolution of racemic 2-methyl-4-phenyl-2,3-butadienoic acid with (*R*)- or (*S*)- α -PEA (Scheme 7).

Similarly, Kato and collaborators employed (*S*)-*N*-benzyl- α -PEA, ((*S*)-**26**) to form a mixture of diastereomeric salts with the key chiral carboxylic acid precursor to the cardioprotective drug (*S*)-**27** (Scheme 8).³² Fractional crystallization of the diastereomeric salts and further elaboration of the precursor led to enantiopure (*S*)-**27**.

Another successful application of α -PEA as resolving agent in the pharmaceuticals area was disclosed by Trung et al.³³ In this work, fractional crystallization of the diastereomeric salts, prepared from (*R*)- α -PEA and racemic ibuprofen, was a key step in the preparation of (*R*)-ibuprofen, an effective analgesic and anti-inflammatory agent.

3.2. Resolution via Covalent Attachment of α -PEA to Racemic Substrates

Very recently, Aitken's group³⁴ carried out the preparation of the four possible stereoisomers of 2-aminocyclobutanecarboxylic acid starting with the condensation of the cis racemate with (*R*)- α -PEA to give a mixture of the corresponding diastereoisomeric amides. The amides were separated by column chromatography, (*R*)- α -PEA cleaved, and a controlled epimerization to the more stable trans isomers carried out, leading to the isolation of all four stereoisomers in enantiomerically pure form.

The 1,3-dipolar cycloaddition of a proline-derived, cyclic nitrene with acrylamide, followed by a reductive cleavage-cyclization domino process, afforded the racemic trans hydroxy ester **28**. Reaction of **28** with (*R*)- α -PEA led to two diastereomeric intermediates, which were separated and subjected to reductive cleavage of the α -PEA moiety to give enantiomerically pure, bicyclic pyrrolizidinones **30** (Scheme 9).³⁵ Pyrrolizidinones **30** were subsequently coupled to α -amino acids to give pseudotripeptides that induce β -turn peptidomimetic foldamers.


Recently, Morales-Rios and co-workers further confirmed the efficiency of α -PEA in the enantioselective synthesis of natural products, by carrying out the preparation of debromoflustramine B (**31**). This approach involved reacting the racemic γ -lactone precursor with (*S*)- α -PEA and chromatographic separation of the resulting diastereomeric *N*-(α -phenylethyl)lactam intermediates. Subsequent hydrolytic cleavage of the α -PEA moiety regenerated the lactone; this was followed by formation of the *N*-methyl lactam with MeNH₂ in methanol and reduction of the amide carbonyl to afford both enantiomers of **31** (Scheme 10).³⁶

An efficient resolution of chiral, C₂-symmetric biphenols, e.g., **32**, was described by Delogu et al., who employed the α -PEA-containing reagent (*S*)-**33** to derivatize **32** into diastereomeric phosphorothioamidates **34**. Separation of the phosphorothioamidate intermediates by fractional crystallization and reductive cleavage with LiAlH₄ afforded biphenols (*M*)-**32** and (*P*)-**32** (Scheme 11).³⁷ A similar approach was employed by Rozenberg and collaborators in the resolution of 15-hydroxy[2.2]paracyclophane-4-carbaldehyde.³⁸ In another interesting development, Zakrzewski and co-workers reported the resolution of a racemic diphosphaferrrocenecarboxylic acid by

chromatographic separation of the corresponding diastereomeric amides formed by reaction of the acid with (*S*)- α -PEA.³⁹


3.3. Chiral Solvating Agents Containing the α -PEA Fragment

One advantage of chiral recognition is that it can be easily evaluated spectroscopically, for example by NMR using chiral discriminating reagents. A variety of chiral solvating agents based on α -PEA are known; they efficiently differentiate the enantiomers of several types of molecules such as chiral Kemp's acid derivatives,⁴⁰ iminoboronate derivatives of chiral alcohols,⁴¹ chiral phospholene oxides,⁴² chiral unsaturated ethers,^{41b,43} and chiral carboxylates.⁴⁴ In a relevant example, our group^{44a} has demonstrated that simple chiral thiourea (*S,S*)-**35** (Scheme 12) is an efficient receptor for chiral carboxylates, such as α -amino




^a The yield of the desired *S* acid was improved to ~70% after recovering the undesirable *R* isomer, racemizing it, and subjecting it to the resolution steps outlined here.

Scheme 8. Key Enantiomeric Separation in the Synthesis of the Cardioprotective Drug (*S*)-27.




Scheme 7. Resolution of *rac*-2-Methyl-4-phenyl-2,3-butadienoic Acid (*rac*-**22**), and Subsequent Formation of Enantioenriched β -Halobutenolides. (Ref. 31a)

Scheme 9. Synthesis of Enantiopure **30** by Covalent Attachment of Enantiopure α -PEA to Racemic Substrates. (Ref. 35a)



Scheme 10. Synthesis of Enantiopure Debromoflustramine B, (R,S)- and (S,R)-31. (Ref. 36a)



Scheme 11. Resolution of C₂-Symmetric Biphenols 32. (Ref. 37a)

and α -hydroxy acids. The diastereomeric complexes obtained from such complexation give rise to distinguishable signals in ¹H NMR spectra, which can be used for determining the optical purity. Furthermore, thiourea (S,S)-35 has proved useful for assigning the absolute configuration of the carboxylic acids; indeed, it was observed that the *R* enantiomer of the carboxylate exhibits a C α -H chemical shift at higher frequency.


4. α -Phenylethylamine as Chiral Auxiliary

The ready availability of α -PEA and its facile incorporation in, and removal from, organic molecules have contributed to making α -PEA into an excellent auxiliary in asymmetric synthesis. Thus, α -PEA-controlled reactions have become an important tool in the preparation of enantioenriched products.

4.1. Diastereoselective Reactions of Chiral Imines Incorporating the α -Phenylethyl Moiety

The condensation of aldehydes and ketones with α -PEA provides the corresponding chiral imines, in which the difference in size between substituents at the stereogenic center helps to differentiate the diastereotopic faces at the prochiral C=N bond, especially in those cases where one conformation of the chiral auxiliary predominates in the transition state. This working hypothesis has been exploited in several methods, which are described below.

The apparent addition of hydride from the side of the less bulky group in the reduction of α -PEA imines has allowed the synthesis of several types of enantioenriched product such as ferrocenes,⁴⁵ amines,⁴⁶ glycosyl amino acids,⁴⁷ amino esters,⁴⁸ and spiro compounds.⁴⁹ Indeed, an excellent procedure for the preparation of chiral primary amines from both cyclic and acyclic dialkyl or aryl-alkyl prochiral ketones was described by Nugent and co-workers.^{46c} The key step in this reaction consists of the asymmetric reductive amination of the ketone with (*R*)- or (*S*)- α -PEA in the presence of a Lewis acid and a conventional hydrogenation catalyst. The correct combination of hydrogenation catalyst, solvent, and temperature is essential for achieving high reaction yields and high diastereomeric excesses. In one example, the reduction of ethyl or hexyl methyl ketone afforded the corresponding amines in good yields and diastereoselectivities (**Scheme 13**, Part (a)). When Raney[®] Ni was employed as the hydrogenation catalyst in the reduction of ethyl methyl ketone, the product was formed with a 74% de, whereas the use of rhodium-on-charcoal led to a significantly



Scheme 12. Thiourea (S,S)-35 and its Proposed Coordination in the Stereodifferentiation of Chiral α -Amino Carboxylates. (Ref. 44a)


lower diastereoselectivity (*de* = 36%, **Scheme 13**, Part (b)). Finally, the α -phenylethylimine derivative of hexyl methyl ketone was reduced with 72% diastereomeric excess in THF as solvent, and 49% *de* in DMSO (**Scheme 13**, Part (c)).

A variety of nucleophiles also add stereoselectively to α -PEA-containing imines. When the nucleophile is cyanide, this leads to enantioenriched non-natural amino acids (**Figure 1**).⁵⁰ An illustrative procedure for the enantioselective preparation of carbocyclic α -amino acids was described by Frahm's group.⁵¹ Chiral ketimines were obtained by condensation of racemic 2-alkylcyclopentanones and (*S*)- α -PEA as chiral auxiliary. In the key stereodifferentiating step, cyanide addition provided diastereomeric mixtures of nitriles, whose composition is dramatically influenced by the nature of the solvent, the temperature, and the size of the substituents. Hydrolysis of the nitriles with conc. H₂SO₄ yielded diastereomeric mixtures of carboxamides, which were separated, hydrogenolyzed, and hydrolyzed to yield the pure stereoisomers of 1-amino-2-methoxycyclopentanecarboxylic acid.

In a similar fashion, the addition of allylmetal nucleophiles to imines provides a valuable route to homoallylic amines, which can undergo any number of transformations at the C=C bond of the allylic fragment.⁵² Gálvez and co-workers⁵³ reported a good application of this approach in the synthesis of non-natural α -amino acid **37**. Thus, chiral *N*-phenylethylimines were prepared from (*R*)-2,3-di-*O*-benzylglyceraldehyde and (*S*)- or (*R*)- α -PEA and treated with allylmagnesium bromide or allyl-9-borabicyclo[3.3.1]nonane (allyl-9-BBN). Conversion of the resulting homoallylic amine into the corresponding syn or anti *N*-Boc-aminodiol was conveniently performed by hydrogenolysis with Pd(OH)₂ on charcoal in the presence of (Boc)₂O. Subsequent treatment with an excess of sodium periodate in the presence of ruthenium trichloride followed by acid hydrolysis, gave α -amino acid (*R*)- or (*S*)-**37**, depending on the nature of the starting material (**Scheme 14**).

4.2. Diastereoselective Reactions of Optically Active Enamines

The enantioselective reduction of enamines, is a useful tool for the synthesis of biologically important molecules such as β -amino



Scheme 13. The Asymmetric Reductive Amination of Prochiral Ketones Evaluated under a Variety of Conditions. (Ref 46e.)

acids, as demonstrated by Wright and collaborators.⁵⁴ These workers prepared diverse cyclic β -amino acids by amination of cyclic β -keto esters with (*R*)- α -PEA in the presence of acetic acid. The enamines were reduced with NaBH₄CN in acetic acid, the resulting mixtures of diastereomers purified by chromatography, separated by recrystallization, and the chiral auxiliary cleaved by hydrogenation over Pd/C. Final saponification of the methyl ester and Fmoc protection of the amino group provided cyclic β -amino acids, e.g. (*S,S*)-**38**, that are suitable for peptide synthesis (**Scheme 15**).⁵⁴ Other remarkable examples that follow this strategy were reported by Lhommet's group⁵⁵ and by Jona







Figure 1. Select Intermediates and α -Amino Acids Obtained by the Diastereoselective Addition of Cyanide to Chiral Imines Containing the α -PEA Moiety.




Scheme 14. The Diastereoselective Allylation of *N*-(α -Phenylethyl)imines in the Enantioselective Preparation of α -Amino Acids.



Scheme 15. Preparation of Cyclic β -Amino Acids (*S,S*)-**38** by Reduction of Enamines. (Ref. 55)



Scheme 16. Michael Addition to Methyl Crotonate of Chiral Imines Derived from 2-Methylcyclopentanone and α -PEA. (Ref. 56)



Scheme 17. α -PEA as Chiral Auxiliary in the Desymmetrization and Alkylation or Acylation of an α -PEA-Containing 2-Piperidinone Enolate. (Ref. 59)

and co-workers⁵⁶ in the synthesis of chiral enamino lactones and dihydrobenzimidazolepiperidine-1,3-dicarboxylate, respectively, in good yields and ee's.

On the other hand, Michael addition of enamines to electron-deficient olefins has emerged as a simple and efficient tool for the stereocontrolled synthesis of quaternary stereocenters, as demonstrated in several recent reports.⁵⁷ A relevant example is the report by Dumas and co-workers⁵⁸ of the stereoselective addition to methyl crotonate of chiral imines derived from 2-methylcyclopentanone and enantiopure α -PEA to give, in the case of the *R* imine, product **39**. Hydrolysis of **39** afforded keto ester **40** in 70% yield, $\geq 98\%$ de, and 96% ee (**Scheme 16**).


4.3. Diastereoselective Alkylation of α -PEA-Containing Enolates

The diastereoselective alkylation of enolates containing the α -phenylethylamino moiety has been employed in the synthesis of open-chain or cyclic enantioenriched compounds, notably piperidinones,⁵⁹ β -amino acids,⁶⁰ diamino dicarboxylic acids,⁶¹ disubstituted phosphonopropanamides⁶² or alkylcyclobutanones.⁶³ An illustrative example was described by Liu and co-workers, who prepared nonracemic disubstituted piperidinones and piperidines by alkylation or acylation of the enolate derived from α -PEA-containing 2-piperidinone (**Scheme 17**).⁵⁹

Our group has successfully employed C₂-symmetric bis(α -PEA) and Evans-type hexahydrobenzoxazolidinone chiral auxiliaries in the asymmetric synthesis of β -amino acids, such as (*R*)- β -homoDopa, (*S*)- β -homophenylalanine, (*S*)- β -homovaline, (*S*)- β -homoleucine, and (*S*)- β -homotryptophan.⁶⁰ Thus, chiral amide **41** is conveniently metallated with lithium hexamethyldisilazide (LiHMDS) whose diastereotopic faces in the resulting enolate allow the stereoselective addition of benzyl bromide in 90% de. In this instance, hydrogenolysis of the benzylated product and hydrolysis under microwave (μW) irradiation result in the formation of (*S*)- β -homophenylalanine with 98% ee (**Scheme 18**).^{60b}

4.4. α -PEA in Diastereoselective Cyclization Reactions

Piperidine, pyrrolidine, and morpholine derivatives have received considerable attention because of their applications in organic, materials, and pharmaceutical chemistry. The use of α -PEA has permitted the preparation of several promising molecules in this class, such as tetrahydroisoquinolines (used in the treatment of Parkinson's disease);⁶⁴ benzodiazepines (sedatives);⁶⁵ oxazolidines,



Scheme 18. Enantioselective Synthesis of β^2 -Homophenylalanine. (Ref. 60b)

oxazolidinones, and lactams (antibiotics);⁶⁶ piperidine and pyrrolidine derivatives;⁶⁷ phthalides (enzyme inhibitors);⁶⁸ noranabasamine alkaloids (anti-inflammatory); and others.⁶⁹

In this regard, Yamaguchi and collaborators⁷⁰ reported an interesting N-heterocyclization of primary amines with diols catalyzed by an iridium complex. When the reaction of (*R*)- α -PEA and 1-phenyl-1,5-pentanediol was carried out in the presence of a catalytic amount of [Cp*IrCl₂]₂, a diastereoisomeric mixture of phenylpiperidines formed in 76% yield, 92% de, and good enantiomeric excess. Reduction of this mixture with Pd/C gave (*S*)-2-phenylpiperidine **42** in 96% yield and 78% ee (**Scheme 19**).

The potential of α -PEA as chiral auxiliary in diastereoselective Diels–Alder reactions was successfully demonstrated in recent years.⁷¹ In a relevant example, Badorrey et al.^{71a} employed a chiral imine dienophile prepared from (*S*)-2,3-di-*O*-benzylglyceraldehyde and (*R*)- α -PEA. In the presence of Danishefsky's diene and ZnI₂, this chiral imine underwent a hetero-Diels–Alder reaction leading to **43** (**Scheme 20**), which is a valuable intermediate in the synthesis of palinavir—a potent inhibitor of the human immunodeficiency virus (HIV-1 and HIV-2).


4.5. α -PEA as Chiral Auxiliary in Miscellaneous Procedures

In addition to the processes described above, the applications of α -PEA as chiral auxiliary include several other kinds of reactions such as diastereoselective carbonylation of propargylic amines,⁷² aldimine coupling of imines with methoxyketene methyltrimethylsilyl acetal in the synthesis of β -amino- α -hydroxy acids,⁷³ preparation of α,β -unsaturated amides from 2-phosphonamides via the Horner–Wadsworth–Emmons reaction,⁷⁴ preparation of imidazole derivatives by the thio-Ugi reaction,⁷⁵ synthesis of azabicyclooctane carboxylic acids,⁷⁶ α -pyridylation of amines through urea coupling, stereoselective lithiation and rearrangement,⁷⁷ development of amino vinyl cyclohexenes,⁷⁸ and synthesis of therapeutic agents such as DPC 961 (HIV non-nucleoside reverse transcriptase inhibitor).⁷⁹


5. α -Phenylethylamine as Chiral Reagent in the Stereodifferentiation of Prochiral Substrates

Conjugate addition of a nitrogen nucleophile to α,β -unsaturated carboxylic acid derivatives is one of the most useful and simplest methods for the formation of N–C bonds.⁸⁰ In particular, addition


of “chiral ammonia” nucleophiles to the acceptor unsaturated carbonyl group produces a new stereogenic center at the β position. The diastereoselectivity of this reaction depends on the nature of the starting ester, the “chiral ammonia” equivalent, and on the particular reaction conditions. Davies, the recognized pioneer and leader of this methodology, has demonstrated the versatility of this approach in the synthesis of a large number of chiral β -amino acids.⁸¹ In a striking example, Davies and co-workers⁸² achieved the diastereoselective conjugate addition of lithium *N*-allyl-(*S*)- α -PEA, (*S*)-**44**, to a wide range of α,β -unsaturated esters (**Scheme 21**). Subsequent ring-closing metathesis (RCM) afforded a variety of substituted cyclic β -amino esters in high de's. Final reduction and




Scheme 19. Diastereoselective Cyclization in the Enantioselective Synthesis of (*S*)-2-Phenylpiperidine, (*S*)-**42**.




Scheme 20. Hetero-Diels–Alder Reaction in the Synthesis of **43**, an Intermediate in the Synthesis of Palinavir. (Ref)



Scheme 21. The Synthesis of β -Amino Acids by Addition of “Chiral Ammonia” Equivalent (*S*)-**44** to α,β -Unsaturated Esters Followed by RCM. (Ref 82)




Scheme 22. Formation of β^2 -Amino Acids by Diastereoselective Michael Addition. (Ref. 8)



R	Yield	ee
Me	55%	41%
<i>t</i> -Bu	94%	32%
FCH ₂	61%	43%
F ₂ CH	66%	59%
F ₃ C	86%	92%
CF ₃ CF ₂	92%	83%
Me ₂ NCH ₂	92%	18%

eq 8



R ^a	Yield	ee
2-MeO	74%	>99%
3-MeO	66%	>99%
4-MeO	80%	>99%

^a Position numbers are with respect to starting material

eq 9 (Ref. 92)




Figure 2. Chiral, α -PEA-Containing Organocatalysts Employed in Asymmetric Aldol Reactions. (Ref. 94–97)

hydrolysis generated the corresponding cyclic β -amino acids; i.e., (*S*)-homoproline, (*S*)-homopipecolic acid, and carbocyclic (*S,S*)-*trans*-pentacin.

A vast number of important molecules have been synthesized by applying this methodology. For example, Michael's group⁸³ prepared the optically active bicyclic compound (–)-indolizidine; Coleman and collaborators⁸⁴ synthesized nonracemic dihydrobenzofuran β -amino acids, which are aspartic acid mimetic that are structurally related to benzodioxole systems; and Podlech⁸⁵ reported the preparation of enantiomerically pure β -amino acids. Moreover, Gellman and co-workers described the enantioselective Michael addition of chiral hydroxylamine (*S*)-**45** to α -alkylacrylates followed by cyclization to give a diastereomeric mixture of α -substituted isoxazolidinones. A set of α -substituted β^2 -amino acids were obtained after separation, hydrogenation, and final Fmoc protection (**Scheme 22**).⁸⁶ Furthermore, Herrera et al.⁸⁷ reported the preparation of chiral α -hydroxy- β -amino acid derivatives based on the diastereoselective addition of (*R*)- α -PEA to representative captodative olefins, a strategy that Blanchet and co-workers⁸⁸ used later in the enantioselective synthesis of (*S*)- β -proline. Lhomme and collaborators⁸⁹ carried out a similar intramolecular Michael addition of (*S*)- α -PEA to chloroacetylenic esters in acetonitrile to form chiral piperidine enamino esters.

Several examples have been described in which α -PEA-lithium amide has been employed for the selective removal of one enantiotopic proton in symmetrical carbonyl compounds.⁹⁰ In a particularly interesting example, Aoki and Koga⁹¹ reported an efficient procedure for the enantioselective deprotonation of 4-*tert*-butylcyclohexanone. These researchers utilized α -PEA-derived chiral lithium amides containing alkyl or fluoroalkyl substituents at the amide nitrogen. The resulting trapped enolates were obtained in yields as high as 94%, and enantiomeric excesses up to 92% (eq 8).⁹¹ In a related study, Clayden and co-workers⁹² reported a lithiation of chiral benzamides containing the α -PEA moiety, which triggered a stereospecific cyclization process (eq 9).

6. α -Phenylethylamine in Organocatalysts

The use of small organic molecules to accelerate chemical reactions, has been of great relevance in organic chemistry. The lack of sensitivity to moisture and oxygen, ready availability, low cost, and low toxicity confer high advantages to organocatalysts over metal catalysts in the production of drugs.⁹³ Several asymmetric organic reactions have made use of organocatalysts containing a chiral α -PEA fragment to induce chirality in the final products. In recent years, numerous reports have described the development of the asymmetric aldol reaction (**Figure 2**).^{93–97} In an excellent example, Chimni and Mahajan⁹⁴ employed protonated chiral prolinamide derivative **46** as asymmetric organocatalyst for the enantioselective aldol reaction in water with good yields and moderate enantioselectivity. Gryko and co-workers⁹⁵ designed organocatalyst **47** by replacing the amide functionality in **46** with the thioamide group in order to increase the acidity of the N–H bond and, as a consequence, to increase yields and enantiomeric excesses (up to 99% ee). Recently, Singh and collaborators⁹⁶ utilized chiral primary and tertiary diamine **48** to catalyze the *syn*-aldol (with unprotected hydroxyacetone) and *anti*-aldol (with cyclic ketones) reactions in aqueous media with high diastereo- and enantioselectivities. Recently, our group described the synthesis of chiral urea **49**, which incorporates (*R*)- or (*S*)- α -PEA, and its application as a Lewis base in the asymmetric aldol reaction.⁹⁷

In other relevant organocatalytic reactions, Tsoogoeva and co-workers⁹⁸ have developed a group of chiral thiourea

derivatives, e.g. **50**, as bifunctional catalysts for the nitro-Michael addition (addition of acetone to nitroolefins), achieving higher enantioselectivities than with proline derivatives. Similarly, Kelleher et al.⁹⁹ synthesized spiro lactam organocatalyst **51** for the asymmetric conjugate addition of aldehydes to nitroolefins in excellent yields, with good diastereoselectivity and enantioselectivity using low catalyst loadings (5 mol %). Our group demonstrated the potential of α -PEA derivatives such as **52**, in the enantioselective amination of α -substituted α -cyanoacetates with azodicarboxylates as electrophiles.¹⁰⁰ González-Olvera et al. described the preparation and use of chiral diazabicyclo[2.2.1]heptane **53** as organocatalyst for the enantioselective Biginelli reaction, giving rise to good yields and enantioselectivities.¹⁰¹ On the other hand, Zhao, Zhou, and collaborators employed a chiral cyclic β -amino alcohol, **54**, as co-catalyst (in combination with L-proline) in the enantioselective Baylis–Hillman reaction between *o*-nitrobenzaldehyde and methyl vinyl ketone, obtaining the corresponding keto alkenyl alcohol with good selectivity.¹⁰² In a preliminary study, Peris and Miller reported the desymmetrization of prochiral ketones through the Baeyer–Villiger oxidation by employing an α -PEA-containing peptide, **55**, to give the corresponding lactones with modest enantioselectivities.¹⁰³ Finally, Hansch et al.¹⁰⁴ prepared enantiomerically pure thiomorpholine **56**, which includes the α -PEA fragment, for application in sulfur ylide mediated asymmetric epoxidation of aldehydes with excellent yields, enantioselectivities, and diastereoselectivities (Figure 3).

7. Concluding Remarks

Because of its low price, accessibility in both enantiomeric forms, facile incorporation in starting materials, and easy removal from final products, α -phenylethylamine (α -PEA) continues to stand as an important chiral reagent in asymmetric synthesis. This is attested to by the great number of methodologies that are available for its application as chiral ligand, in the resolution of racemic mixtures, as chiral auxiliary, as chiral base, and in the development of new organocatalysts. The present review has described noteworthy recent developments in the chemistry of α -PEA, which may motivate chemists involved in enantioselective synthesis to take advantage of this versatile molecule.

8. Acknowledgements

The authors are grateful to the Consejo Nacional de Ciencia y Tecnología (CONACyT, México) for financial support. Y. Bandala is also indebted to CINVESTAV-IPN for a postdoctoral fellowship.

9. References and Notes

- (1) (a) Juaristi, E.; Escalante, J.; León-Romo, J. L.; Reyes, A. *Tetrahedron: Asymmetry* **1998**, *9*, 715. (b) Juaristi, E.; León-Romo, J. L.; Reyes, A.; Escalante, J. *Tetrahedron: Asymmetry* **1999**, *10*, 2441.
- (2) See for example: (a) Noyori, R.; Ohkuma, T. *Pure Appl. Chem.* **1999**, *71*, 1493. (b) Noyori, R.; Ohkuma, T. *Angew. Chem., Int. Ed.* **2001**, *40*, 40. (c) For a review on mono- and bidentate ligands based on arene–chromium complexes starting from (*R*)- or (*S*)- α -PEA and its derivatives, see Salzer, A. *Coord. Chem. Rev.* **2003**, *242*, 59.
- (3) (a) Wang, D.-Y.; Huang, J.-D.; Hu, X.-P.; Deng, J.; Yu, S.-B.; Duan, Z.-C.; Zheng, Z. *J. Org. Chem.* **2008**, *73*, 2011. (b) Zhou, X.-M.; Huang, J.-D.; Luo, L.-B.; Zhang, C.-L.; Zheng, Z.; Hu, X.-P. *Tetrahedron: Asymmetry* **2010**, *21*, 420.
- (4) Brauer, D. J.; Kottsieper, K. W.; Roßenbach, S.; Stelzer, O. *Eur. J. Inorg. Chem.* **2003**, 1748.




Figure 3. Remarkable, α -PEA-Containing Organocatalysts Utilized in Various Asymmetric Reactions. (Ref. 98–104)

- (5) Ireland, T.; Tappe, K.; Grossheimann, G.; Knochel, P. *Chem.—Eur. J.* **2002**, *8*, 843.
- (6) (a) Huang, J.-D.; Hu, X.-P.; Duan, Z.-C.; Zeng, Q.-H.; Yu, S.-B.; Deng, J.; Wang, D.-Y.; Zheng, Z. *Org. Lett.* **2006**, *8*, 4367. See also: (b) Kangying, L.; Zhenghong, Z.; Goufeng, Z.; Chuchi, T. *Heteroatom Chem.* **2003**, *14*, 546. (c) Boezio, A. A.; Charette, A. B. *J. Am. Chem. Soc.* **2003**, *125*, 1692.
- (7) (a) Mastranzo, V. M.; Quintero, L.; Anaya de Parrodi, C.; Juaristi, E.; Walsh, P. J. *Tetrahedron* **2004**, *60*, 1781. (b) Santacruz, E.; Huelgas, G.; Angulo, S. K.; Mastranzo, V. M.; Hernández-Ortega, S.; Aviña, J. A.; Juaristi, E.; Anaya de Parrodi, C.; Walsh, P. J. *Tetrahedron: Asymmetry* **2009**, *20*, 2788.
- (8) (a) Zaitsev, A. B.; Adolfsson, H. *Org. Lett.* **2006**, *8*, 5129. (b) Ahlford, K.; Ekström, J.; Zaitsev, A. B.; Ryberg, P.; Eriksson, L.; Adolfsson, H. *Chem.—Eur. J.* **2009**, *15*, 11197.
- (9) Melgar-Fernández, R.; González-Olvera, R.; Olivares-Romero, J. L.; González-López, V.; Romero-Ponce, L.; Ramírez-Zárate, M. R.; Demare, P.; Regla, I.; Juaristi, E. *Eur. J. Org. Chem.* **2008**, 655.
- (10) Wilken, J.; Erny, S.; Wassmann, S.; Martens, J. *Tetrahedron: Asymmetry* **2000**, *11*, 2143.
- (11) (a) Zhang, Z.; Li, M.; Zi, G. *Chirality* **2007**, *19*, 802. (b) Zhang, Z.; Bai, X.; Liu, R.; Zi, G. *Inorg. Chim. Acta* **2009**, *362*, 1687.
- (12) Liu, D.-X.; Zhang, L.-C.; Wang, Q.; Da, C.-S.; Xin, Z.-Q.; Wang, R.; Choi, M. C. K.; Chan, A. S. C. *Org. Lett.* **2001**, *3*, 2733.
- (13) Xu, L.; Shen, X.; Zhang, C.; Mikami, K. *Chirality* **2005**, *17*, 476.
- (14) Anaya de Parrodi, C.; Juaristi, E. *Synlett* **2006**, 2699.
- (15) Mastranzo, V. M.; Santacruz, E.; Huelgas, G.; Paz, E.; Sosa-Rivadeneira, M. V.; Bernès, S.; Juaristi, E.; Quintero, L.; Anaya de Parrodi, C. *Tetrahedron: Asymmetry* **2006**, *17*, 1663.
- (16) Isleyen, A.; Dogan, Ö. *Tetrahedron: Asymmetry* **2007**, *18*, 679.
- (17) (a) Wang, Y.; Li, X.; Ding, K. *Tetrahedron: Asymmetry* **2002**, *13*, 1291. (b) Huang, J.-D.; Hu, X.-P.; Yu, S.-B.; Deng, J.; Wang, D.-Y.; Duan, Z.-C.; Zheng, Z. *J. Mol. Catal. A: Chem.* **2007**, *270*, 127. See also: (c) Huang, J. D.; Hu, X.-P.; Zheng, Z. *Chin. Chem. Lett.* **2008**, *19*, 261. (d) Zhang, A.; RajanBabu, T. V. *J. Am. Chem. Soc.* **2006**, *128*, 5620.
- (18) Baudequin, C.; Chaturvedi, D.; Tsogoeva, S. B. *Eur. J. Org. Chem.* **2007**, 2623.
- (19) Alonso, I.; Trillo, B.; López, F.; Montserrat, S.; Ujaque, G.; Castedo, L.; Lledós, A.; Mascareñas, J. L. *J. Am. Chem. Soc.* **2009**, *131*, 13020.
- (20) (a) Yang, Z.; Zhou, Z.; He, K.; Wang, L.; Zhao, G.; Zhou, Q.; Tang, C. *Tetrahedron: Asymmetry* **2003**, *14*, 3937. (b) For another interesting enantioselective alkyl addition to carbonyl compounds,

- see Duguet, N.; Petit, S. M.; Marchand, P.; Harrison-Marchand, A.; Maddaluno, J. J. *Org. Chem.* **2008**, *73*, 5397.
- (21) (a) Wakabayashi, K.; Aikawa, K.; Kawauchi, S.; Mikami, K. *J. Am. Chem. Soc.* **2008**, *130*, 5012. (b) Also, see Rimkus, A.; Sewald, N. *Org. Lett.* **2003**, *5*, 79.
- (22) Traber, B.; Jung, Y.-G.; Park, T. K.; Hong, J.-I. *Bull. Korean Chem. Soc.* **2001**, *22*, 547.
- (23) (a) Santoni, G.; Licini, G.; Rehder, D. *Chem.—Eur. J.* **2003**, *9*, 4700. (b) Rehder, D.; Santoni, G.; Licini, G. M.; Schulzke, C.; Meier, B. *Coord. Chem. Rev.* **2003**, *237*, 53.
- (24) Millward, D. B.; Cole, A. P.; Waymouth, R. M. *Organometallics* **2000**, *19*, 1870.
- (25) Bianchini, C.; Lee, H. M. *Organometallics* **2000**, *19*, 1833.
- (26) (a) Dyer, U. C.; Shapland, P. D.; Tiffin, P. D. *Tetrahedron Lett.* **2001**, *42*, 1765. (b) He, Q.; Peng, Y.-F.; Rohani, S. *Chirality* **2010**, *22*, 16.
- (27) (a) Vitharana, D.; France, J. E.; Scarpetti, D.; Bonneville, G. W.; Majer, P.; Tsukamoto, T. *Tetrahedron: Asymmetry* **2002**, *13*, 1609. (b) Reddy, P. M.; Kovach, I. M. *Tetrahedron Lett.* **2002**, *43*, 4063. (c) Kobayashi, Y.; Maeda, J.; Saigo, K. *Tetrahedron: Asymmetry* **2006**, *17*, 1617. (d) Kaboudin, B.; Haghghat, H.; Yokomatsu, T. *Tetrahedron: Asymmetry* **2008**, *19*, 862.
- (28) Boussac, H.; Crassous, J.; Dutasta, J.-P.; Grosvalet, L.; Thozet, A. *Tetrahedron: Asymmetry* **2002**, *13*, 975.
- (29) (a) Szabó, D.; Nemes, A.; Kövesdi, I.; Farkas, V.; Hollósi, M.; Rábai, J. *J. Fluorine Chem.* **2006**, *127*, 1405. (b) Uchida, Y.; Uematsu, T.; Nakayama, Y.; Takahashi, H.; Tsue, H.; Tanaka, K.; Tamura, R. *Chirality* **2008**, *20*, 282. (c) Bolchi, C.; Pallavicini, M.; Fumagalli, L.; Marchini, N.; Moroni, B.; Rusconi, C.; Valoti, E. *Tetrahedron: Asymmetry* **2005**, *16*, 1639.
- (30) (a) Kato, Y.; Niiyama, K.; Jona, H.; Okada, S.; Akao, A.; Hiraga, S.; Tsuchiya, Y.; Tomimoto, K.; Mase, T. *Chem. Pharm. Bull.* **2002**, *50*, 1066. (b) Bálint, J.; Schindler, J.; Egri, G.; Hanusz, M.; Marthi, K.; Juvancz, Z.; Fogassy, E. *Tetrahedron: Asymmetry* **2004**, *15*, 3401. (c) Ashcroft, C. P.; Challenger, S.; Clifford, D.; Derrick, A. M.; Hajikarimian, Y.; Slucock, K.; Silk, T. V.; Thomson, N. M.; Williams, J. R. *J. Org. Process Res. Devel.* **2005**, *9*, 663. (d) Hu, Y.; Yu, S.-L.; Yang, Y.-J.; Zhu, J.; Deng, J.-G. *Chinese J. Chem.* **2006**, *24*, 795. (e) Rasmy, O. M.; Vaid, R. K.; Semo, M. J.; Chelius, E. C.; Robey, R. L.; Alt, C. A.; Rhodes, G. A.; Vicenzi, J. T. *Org. Process Res. Devel.* **2006**, *10*, 28. (f) Hu, B.; Prashad, M.; Har, D.; Prasad, K.; Repič, O.; Blacklock, T. J. *Org. Process Res. Devel.* **2007**, *11*, 90. (g) Kremsner, J.; Wallfisch, B. C.; Belaj, F.; Uray, G.; Kappe, C. O.; Wentrup, C.; Kollenz, G. *Eur. J. Org. Chem.* **2008**, 3382.
- (31) (a) Ma, S.; Wu, S. *Chem. Commun.* **2001**, 441. (b) See also Ma, S.; Shi, Z. *Chem. Commun.* **2002**, 540.
- (32) Kato, T.; Ozaki, T.; Tsuzuki, K.; Ohi, N. *Org. Process Res. Devel.* **2001**, *5*, 122.
- (33) Trung, T. Q.; Kim, J. M.; Kim, K. H. *Arch. Pharm. Res.* **2006**, *29*, 108.
- (34) Fernandes, C.; Pereira, E.; Faure, S.; Aitken, D. J. *J. Org. Chem.* **2009**, *74*, 3217.
- (35) (a) Cordero, F. M.; Pisaneschi, F.; Batista, K. M.; Valenza, S.; Machetti, F.; Brandi, A. *J. Org. Chem.* **2005**, *70*, 856. (b) Salvati, M.; Cordero, F. M.; Pisaneschi, F.; Bucelli, F.; Brandi, A. *Tetrahedron* **2005**, *61*, 8836. For additional examples of the applications of α -PEA in 1,3-dipolar cycloadditions, see: (c) Carmona, D.; Lamata, M. P.; Viguri, F.; Rodríguez, R.; Oro, L. A.; Lahoz, F. J.; Balana, A. I.; Tejero, T.; Merino, P. J. *Am. Chem. Soc.* **2005**, *127*, 13386. (d) Carmona, D.; Lamata, M. P.; Viguri, F.; Rodríguez, R.; Fischer, T.; Lahoz, F. J.; Dobrinovitch, I. T.; Oro, L. A. *Adv. Synth. Catal.* **2007**, *349*, 1751.
- (36) (a) Morales-Rios, M. S.; Rivera-Becerril, E.; Joseph-Nathan, P. *Tetrahedron: Asymmetry* **2005**, *16*, 2493. For related examples, see: (b) Paquette, L. A.; Boulet, S. L. *Synthesis* **2002**, 888. (c) Pawłowska, J.; Krawczyk, K. K.; Wojtasiewicz, K.; Maurin, J. K.; Czarnocki, Z. *Monatsh. Chem.* **2009**, *140*, 83. (d) Marchalin, Š.; Cvopová, K.; Križ, M.; Baran, P.; Oulyadi, H.; Daich, A. *J. Org. Chem.* **2004**, *69*, 4227.
- (37) (a) Delogu, G.; Fabbri, D.; Dettori, M. A.; Forni, A.; Casalone, G. *Tetrahedron: Asymmetry* **2000**, *11*, 4417. For the synthesis of α -PEA phosphorus amide derivatives, see: (b) Kolodyazhnyi, O. I.; Gryshkun, E. V.; Andrushko, N. V.; Freytag, M.; Jones, P. G.; Schmutzler, R. *Tetrahedron: Asymmetry* **2003**, *14*, 181. (c) Kolodyazhnyi, O. I.; Andrushko, N. V.; Grishkun, E. V. *Russ. J. Gen. Chem.* **2004**, *74*, 515.
- (38) (a) Rozenberg, V. I.; Antonov, D. Y.; Sergeeva, E. V.; Vorontsov, E. V.; Starikova, Z. A.; Fedyanin, I. V.; Schulz, C.; Hopf, H. *Eur. J. Org. Chem.* **2003**, 2056. (b) A for the resolution of C_3 -symmetrical phosphanes using α -PEA, see Wyatt, P.; Eley, H.; Charmant, J.; Daniel, B. J.; Kantacha, A. *Eur. J. Org. Chem.* **2003**, 4216.
- (39) (a) Klys, A.; Nazarski, R. B.; Zakrzewski, J. J. *Organomet. Chem.* **2001**, 627, 135. (b) See also Tebben, L.; Kehr, G.; Fröhlich, R.; Erker, G. *Synthesis* **2004**, 1971.
- (40) (a) Hirose, T.; Naito, K.; Shitara, H.; Nohira, H.; Baldwin, B. W. *Tetrahedron: Asymmetry* **2001**, *12*, 375. (b) Hirose, T.; Naito, K.; Nakahara, M.; Shitara, H.; Aoki, Y.; Nohira, H.; Baldwin, B. W. *J. Inclusion Phenom. Macrocycl. Chem.* **2002**, *43*, 87.
- (41) (a) Kelly, A. M.; Pérez-Fuertes, Y.; Arimori, S.; Bull, S. D.; James, T. D. *Org. Lett.* **2006**, *8*, 1971. (b) Uccello-Barretta, G.; Bernardini, R.; Lazzaroni, R.; Salvadori, P. *J. Organomet. Chem.* **2000**, 598, 174.
- (42) Pakulski, Z.; Demchuk, O. M.; Kwiatosz, R.; Osinski, P. W.; Świerczyńska, W.; Pietrusiewicz, K. M. *Tetrahedron: Asymmetry* **2003**, *14*, 1459.
- (43) Uccello-Barretta, G.; Bernardini, R.; Balzano, F.; Lazzaroni, R.; Salvadori, P. *J. Organomet. Chem.* **2000**, 605, 68.
- (44) (a) Hernández-Rodríguez, M.; Juaristi, E. *Tetrahedron* **2007**, *63*, 7673. (b) Wang, W.; Ma, F.; Shen, X.; Zhang, C. *Tetrahedron: Asymmetry* **2007**, *18*, 832.
- (45) Cayuela, E. M.; Xiao, L.; Sturm, T.; Manzano, B. R.; Jalon, F. A.; Weissensteiner, W. *Tetrahedron: Asymmetry* **2000**, *11*, 861.
- (46) (a) Gong, Y.; Kato, K. *Tetrahedron: Asymmetry* **2001**, *12*, 2121. (b) Storace, L.; Anzalone, L.; Confalone, P. N.; Davis, W. P.; Fortunak, J. M.; Giangiordano, M.; Haley, J. J., Jr.; Kamholz, K.; Li, H.-Y.; Ma, P.; Nugent, W. A.; Parsons, R. L., Jr.; Sheeran, P. J.; Silverman, C. E.; Waltermire, R. E.; Wood, C. C. *Org. Process Res. Devel.* **2002**, *6*, 54. (c) Saidi, M. R.; Azzizi, N. *Tetrahedron: Asymmetry* **2003**, *14*, 389. (d) Grishina, G. V.; Luk'yanenko, E. R.; Borisenko, A. A. *Russ. J. Org. Chem.* **2005**, *41*, 807. (e) Nugent, T. C.; Wakchaure, V. N.; Ghosh, A. K.; Mohanty, R. R. *Org. Lett.* **2005**, *7*, 4967. (f) Kanai, M.; Ueda, K.; Yasumoto, M.; Kuriyama, Y.; Inomiya, K.; Ootsuka, T.; Katsuhara, Y.; Higashiyama, K.; Ishii, A. *J. Fluorine Chem.* **2005**, *126*, 377. (g) Jaroskova, L.; van der Veken, L.; de Belser, P.; Diels, G.; de Groot, A.; Linders, J. T. M. *Tetrahedron Lett.* **2006**, *47*, 8063. (h) Nugent, T. C.; Ghosh, A. K.; Wakchaure, V. N.; Mohanty, R. R. *Adv. Synth. Catal.* **2006**, *348*, 1289. (i) Federsel, H.-J.; Hedberg, M.; Qvarnström, F. R.; Sjögren, M. P. T.; Tian, W. *Acc. Chem. Res.* **2007**, *40*, 1377. (j) Camps, P.; Muñoz-Torrero, D.; Rull, J.; Font-Bardia, M.; Solans, X. *Tetrahedron: Asymmetry* **2007**, *18*, 2947.
- (47) Grison, C.; Coutrot, F.; Coutrot, P. *Tetrahedron* **2001**, *57*, 6215.
- (48) (a) Lee, K.-D.; Suh, J.-M.; Park, J.-H.; Ha, H.-J.; Choi, H. G.; Park, C. S.; Chang, J. W.; Lee, W. K.; Dong, Y.; Yun, H. *Tetrahedron* **2001**, *57*, 8267. (b) Nugent, T. C.; Ghosh, A. K. *Eur. J. Org. Chem.* **2007**, 3863.
- (49) Kelleher, F.; Kelly, S. *Tetrahedron Lett.* **2006**, *47*, 5247.

- (50) (a) Jiang, B.; Yang, C.-G.; Gu, X.-H. *Tetrahedron Lett.* **2001**, *42*, 2545. (b) Jiang, B.; Huang, Z.-G. *Synthesis* **2005**, 2198. (c) Leclerc, E.; Mangency, P.; Henryon, V. *Tetrahedron: Asymmetry* **2000**, *11*, 3471. (d) Iovel, I.; Golomba, L.; Popelis, J.; Lukevics, E. *Appl. Organomet. Chem.* **2002**, *16*, 133. (e) Warmuth, R.; Munsch, T. E.; Stalker, R. A.; Li, B.; Beatty, A. *Tetrahedron* **2001**, *57*, 6383. (f) Schlosser, M.; Brügger, N.; Schmidt, W.; Amrhein, N. *Tetrahedron* **2004**, *60*, 7731. (g) Resnick, L.; Galante, R. J. *Tetrahedron: Asymmetry* **2006**, *17*, 846.
- (51) (a) Wede, J.; Volk, F.-J.; Frahm, A. W. *Tetrahedron: Asymmetry* **2000**, *11*, 3231. (b) Fondecar, K. P. P.; Volk, F.-J.; Khaliq-uz-Zaman, S. M.; Bisel, P.; Frahm, A. W. *Tetrahedron: Asymmetry* **2002**, *13*, 2241. (c) Bisel, P.; Fondecar, K. P.; Volk, F.-J.; Frahm, A. W. *Tetrahedron* **2004**, *60*, 10541. (d) Meyer, U.; Breitling, E.; Bisel, P.; Frahm, A. W. *Tetrahedron: Asymmetry* **2004**, *15*, 2029. (e) Meyer, U.; Bisel, P.; Bräuner-Osborne, H.; Madsen, U.; Höfner, G.; Wanner, K. T.; Frahm, A. W. *Chirality* **2005**, *17*, 99.
- (52) (a) Martelli, G.; Morri, S.; Savoia, D. *Tetrahedron* **2000**, *56*, 8367. (b) Alvaro, G.; Grepioni, F.; Grilli, S.; Maini, L.; Martelli, G.; Savoia, D. *Synthesis* **2000**, 581. (c) Periasamy, M.; Suresh, S.; Ganesan, S. S. *Tetrahedron: Asymmetry* **2006**, *17*, 1323. (d) Wakchaure, V. N.; Mohanty, R. R.; Shaikh, A. J.; Nugent, T. C. *Eur. J. Org. Chem.* **2007**, 959.
- (53) (a) Badorrey, R.; Cativiela, C.; Díaz-de-Villegas, M. D.; Díez, R.; Gálvez, J. A. *Eur. J. Org. Chem.* **2002**, 3763. (b) Díez, R.; Badorrey, R.; Díaz-de-Villegas, M. D.; Gálvez, J. A. *Eur. J. Org. Chem.* **2007**, 2114.
- (54) (a) Wright, K.; Crisma, M.; Toniolo, C.; Török, R.; Péter, A.; Wakselman, M.; Mazaleyrat, J.-P. *Tetrahedron Lett.* **2003**, *44*, 3381. (b) Wright, K.; Sarciaux, M.; de Castries, A.; Wakselman, M.; Mazaleyrat, J.-P.; Toffoletti, A.; Corvajá, C.; Crisma, M.; Peggion, C.; Formaggio, F.; Toniolo, C. *Eur. J. Org. Chem.* **2007**, 3133. (c) Wright, K.; de Castries, A.; Sarciaux, M.; Formaggio, F.; Toniolo, C.; Toffoletti, A.; Wakselman, M.; Mazaleyrat, J.-P. *Tetrahedron Lett.* **2005**, *46*, 5573. (d) See also Matsuo, J.; Okano, M.; Takeuchi, K.; Tanaka, H.; Ishibashi, H. *Tetrahedron: Asymmetry* **2007**, *18*, 1906.
- (55) David, O.; Bellec, C.; Fargeau-Bellassoued, M.-C.; Lhomme, G. *Heterocycles* **2001**, *55*, 1689.
- (56) Jona, H.; Shibata, J.; Asai, M.; Goto, Y.; Arai, S.; Nakajima, S.; Okamoto, O.; Kawamoto, H.; Iwasawa, Y. *Tetrahedron: Asymmetry* **2009**, *20*, 2439.
- (57) (a) Gassama, A.; d'Angelo, J.; Cavé, C.; Mahuteau, J.; Riche, C. *Eur. J. Org. Chem.* **2000**, 3165. (b) Schenato, R. A.; dos Santos, É. M.; Tenius, B. S. M.; Costa, P. R. R.; Caracelli, I.; Zukerman-Schpector, J. *Tetrahedron: Asymmetry* **2001**, *12*, 579. (c) Nour, M.; Tan, K.; Jankowski, R.; Cavé, C. *Tetrahedron: Asymmetry* **2001**, *12*, 765. (d) Jabin, I.; Revial, G.; Pfau, M.; Netchitaïlo, P. *Tetrahedron: Asymmetry* **2002**, *13*, 563. (e) Desmaële, D.; Delarue-Cochin, S.; Cavé, C.; d'Angelo, J.; Morgant, G. *Org. Lett.* **2004**, *6*, 2421. (f) Hendra, F.; Nour, M.; Baglin, I.; Morgant, G.; Cavé, C. *Tetrahedron: Asymmetry* **2004**, *15*, 1027. (g) Delarue-Cochin, S.; Pan, J.-J.; Dauteloup, A.; Hendra, F.; Angoh, R. G.; Joseph, D.; Stephens, P. J.; Cavé, C. *Tetrahedron: Asymmetry* **2007**, *18*, 685. (h) Delarue-Cochin, S.; Bahlaouan, B.; Hendra, F.; Ourévitch, M.; Joseph, D.; Morgant, G.; Cavé, C. *Tetrahedron: Asymmetry* **2007**, *18*, 759. (i) Huang, X.; Chen, X.; Chen, Y.; Zhang, A.; Li, X. *Tetrahedron: Asymmetry* **2008**, *19*, 2529.
- (58) Camara, C.; Joseph, D.; Dumas, F.; d'Angelo, J.; Chiaroni, A. *Tetrahedron Lett.* **2002**, *43*, 1445.
- (59) (a) Liu, L. T.; Hong, P.-C.; Huang, H.-L.; Chen, S.-F.; Wang, C.-L. J.; Wen, Y.-S. *Tetrahedron: Asymmetry* **2001**, *12*, 419. (b) See also Sosa-Rivadeneira, M.; Quintero, L.; Anaya de Parrodi, C.; Bernés, S.; Castellanos, E.; Juaristi, E. *ARKIVOC* **2003** (xi), 61.
- (60) (a) Avila-Ortiz, C. G.; Reyes-Rangel, G.; Juaristi, E. *Tetrahedron* **2005**, *61*, 8372. (b) Guzmán-Mejía, R.; Reyes-Rangel, G.; Juaristi, E. *Nature Protocols* **2007**, *2*, 2759. (c) Reyes-Rangel, G.; Jiménez-González, E.; Olivares-Romero, J. L.; Juaristi, E. *Tetrahedron: Asymmetry* **2008**, *19*, 2839. (d) For a recent review on the synthesis of β -amino acids, see Bandala, Y.; Juaristi, E. Recent Developments in the Synthesis of β -Amino Acids. In *Amino Acids, Peptides and Proteins in Organic Chemistry*; Hughes, A. B., Ed.; Origins and Synthesis of Amino Acids Series, Vol. 1; Wiley-VCH: Weinheim, 2009; pp 291–365.
- (61) (a) Paradisi, F.; Porzi, G.; Rinaldi, S.; Sandri, S. *Tetrahedron: Asymmetry* **2000**, *11*, 1259. (b) Paradisi, F.; Porzi, G.; Rinaldi, S.; Sandri, S. *Tetrahedron: Asymmetry* **2000**, *11*, 4617. (c) Piccinelli, F.; Porzi, G.; Sandri, M.; Sandri, S. *Tetrahedron: Asymmetry* **2003**, *14*, 393. (d) Balducci, D.; Contaldi, S.; Lazzari, I.; Porzi, G. *Tetrahedron: Asymmetry* **2009**, *20*, 1398.
- (62) Ordóñez, M.; Hernández-Fernández, E.; Rojas-Cabrera, H.; Labastida-Galván, V. *Tetrahedron: Asymmetry* **2008**, *19*, 2767.
- (63) Hazelard, D.; Fadel, A. *Tetrahedron: Asymmetry* **2005**, *16*, 2067.
- (64) Toda, J.; Matsumoto, S.; Saitoh, T.; Sano, T. *Chem. Pharm. Bull.* **2000**, *48*, 91.
- (65) Broggin, G.; Garanti, L.; Molteni, G.; Pilati, T. *Tetrahedron: Asymmetry* **2001**, *12*, 1201.
- (66) (a) David, O.; Vanucci-Bacqué, C.; Fargeau-Bellassoued, M. C.; Lhomme, G. *Heterocycles* **2004**, *62*, 839. (b) Sugiyama, S.; Arai, S.; Ishii, K. *Heterocycles* **2005**, *65*, 149. (c) Todorov, A. R.; Kurteva, V. B.; Bontchev, R. P.; Vassilev, N. G. *Tetrahedron* **2009**, *65*, 10339.
- (67) (a) Biala, J.; Czarnocki, Z.; Maurin, J. K. *Tetrahedron: Asymmetry* **2002**, *13*, 1021. (b) Pawlowska, J.; Czarnocki, Z.; Wojtasiewicz, K.; Maurin, J. K. *Tetrahedron: Asymmetry* **2003**, *14*, 3335. (c) Hashimoto, N.; Funatomi, T.; Misaki, T.; Tanabe, Y. *Tetrahedron* **2006**, *62*, 2214. (d) Hermet, J.-P. R.; Viterisi, A.; Wright, J. M.; McGrath, M. J.; O'Brien, P.; Whitwood, A. C.; Gilday, J. *Org. Biomol. Chem.* **2007**, *5*, 3614. (e) Fadel, A.; Lahrache, N. *J. Org. Chem.* **2007**, *72*, 1780.
- (68) Karnik, A. V.; Kamath, S. S. *Synthesis* **2008**, 1832.
- (69) (a) Miao, L.; DiMaggio, S. C.; Shu, H.; Trudell, M. L. *Org. Lett.* **2009**, *11*, 1579. (b) Mlostoń, G.; Mucha, P.; Urbaniak, K.; Broda, K. *Helv. Chim. Acta.* **2008**, *91*, 232.
- (70) (a) Fujita, K.; Fujii, T.; Yamaguchi, R. *Org. Lett.* **2004**, *6*, 3525. (b) Fujita, K.; Fujii, T.; Komatsubara, A.; Enoki, Y.; Yamaguchi, R. *Heterocycles* **2007**, *74*, 673. (c) For a similar procedure in the synthesis of azetidinecarboxylic acids, see Couty, F.; Evano, G.; Vargas-Sanchez, M.; Bouzas, G. *J. Org. Chem.* **2005**, *70*, 9028.
- (71) See for example: (a) Badorrey, R.; Cativiela, C.; Díaz-de-Villegas, M. D.; Gálvez, J. A. *Tetrahedron* **2002**, *58*, 341. (b) Ekegren, J. K.; Modin, S. A.; Alonso, D. A.; Andersson, P. G. *Tetrahedron: Asymmetry* **2002**, *13*, 447. (c) Adams, H.; Bawa, R. A.; McMillan, K. G.; Jones, S. *Tetrahedron: Asymmetry* **2007**, *18*, 1003. (d) Caballero, E.; Figueroa, J.; Puebla, P.; Peláez, R.; Tomé, F.; Medarde, M. *Eur. J. Org. Chem.* **2008**, 4004.
- (72) Rezaei, H.; Marek, I.; Normant, J. F. *Tetrahedron* **2001**, *57*, 2477.
- (73) Ha, H.-J.; Ahn, Y.-G.; Woo, J.-S.; Lee, G. S.; Lee, W. K. *Bull. Chem. Soc. Jpn.* **2001**, *74*, 1667.
- (74) (a) Hernández-Fernández, E.; Fernández-Zertuche, M.; Garcia-Barradas, O.; Muñoz-Muñoz, O.; Ordóñez, M. *Synlett* **2006**, 440. (b) For a recent review on the synthesis of aminophosphonic acids using α -PEA, see Ordóñez, M.; Rojas-Cabrera, H.; Cativiela, C. *Tetrahedron* **2009**, *65*, 17.
- (75) Gulevich, A. V.; Balenkova, E. S.; Nenajdenko, V. G. *J. Org. Chem.* **2007**, *72*, 7878.
- (76) Gelmi, M. L.; Cattaneo, C.; Pellegrino, S.; Clerici, F.; Montali, M.; Martini, C. *J. Org. Chem.* **2007**, *72*, 9811.

- (77) (a) Clayden, J.; Hennecke, U. *Org. Lett.* **2008**, *10*, 3567. See also: (b) Dem'yanovich, V. M.; Shishkina, I. N.; Borisenko, A. A. *Russ. J. Org. Chem.* **2000**, *36*, 1168. (c) Dem'yanovich, V. M.; Shishkina, I. N.; Gritsyuk, Y. A.; Potekhin, K. A.; Chesnova, A. V.; Ashkinadze, L. D. *Russ. J. Org. Chem.* **2005**, *41*, 556.
- (78) Jacquet, O.; Charnay, F.; Fiaud, J.-C.; Ollivier, J. *Synthesis* **2009**, 3047.
- (79) Magnus, N. A.; Confalone, P. N.; Storace, L. *Tetrahedron Lett.* **2000**, *41*, 3015.
- (80) For a recent review on the asymmetric aza-Michael addition, see Krishna, P. R.; Sreeshailam, A.; Srinivas, R. *Tetrahedron* **2009**, *65*, 9657.
- (81) (a) Bull, S. D.; Davies, S. G.; Kelly, P. M.; Gianotti, M.; Smith, A. D. *J. Chem. Soc., Perkin Trans. 1* **2001**, 3106. (b) Davies, S. G.; Iwamoto, K.; Smethurst, C. A. P.; Smith, A. D.; Rodriguez-Solla, H. *Synlett* **2002**, 1146. (c) Davies, S. G.; Diez, D.; Dominguez, S. H.; Garrido, N. M.; Kruchinin, D.; Price, P. D.; Smith, A. D. *Org. Biomol. Chem.* **2005**, *3*, 1284. (d) Davies, S. G.; Smith, A. D.; Price, P. D. *Tetrahedron: Asymmetry* **2005**, *16*, 2833. (e) Davies, S. G.; Garrido, N. M.; Kruchinin, D.; Ichihara, O.; Kotchie, L. J.; Price, P. D.; Mortimer, A. J. P.; Russell, A. J.; Smith, A. D. *Tetrahedron: Asymmetry* **2006**, *17*, 1793. (f) Davies, S. G.; Mulvaney, A. W.; Rusell, A. J.; Smith, A. D. *Tetrahedron: Asymmetry* **2007**, *18*, 1554. (g) Davies, S. G.; Fletcher, A. M.; Roberts, P. M.; Smith, A. D. *Tetrahedron* **2009**, *65*, 10192. (h) Davies, S. G.; Hughes, D. G.; Price, P. D.; Roberts, P. M.; Russell, A. J.; Smith, A. D.; Thomson, J. E.; Williams, O. M. H. *Synlett* **2010**, 567.
- (82) Chippindale, A. M.; Davies, S. G.; Iwamoto, K.; Parkin, R. M.; Smethurst, C. A. P.; Smith, A. D.; Rodriguez-Solla, H. *Tetrahedron* **2003**, *59*, 3253.
- (83) (a) Michael, J. P.; Gravestock, D. *J. Chem. Soc., Perkin Trans. 1* **2000**, 1919. (b) De Koning, C. B.; Michael, J. P.; Riley, D. L. *Heterocycles* **2009**, *79*, 935.
- (84) Coleman, P. J.; Hutchinson, J. H.; Hunt, C. A.; Lu, P.; Delaporte, E.; Rushmore, T. *Tetrahedron Lett.* **2000**, *41*, 5803.
- (85) Podlech, J. *Synth. Commun.* **2000**, *30*, 1779.
- (86) Lee, H.-S.; Park, J.-S.; Kim, B. M.; Gellman, S. H. *J. Org. Chem.* **2003**, *68*, 1575.
- (87) Herrera, R.; Jiménez-Vázquez, H. A.; Tamariz, J. *ARKIVOC* **2005** (vi), 233.
- (88) Blanchet, J.; Pouliquen, M.; Lasne, M.-C.; Rouden, J. *Tetrahedron Lett.* **2007**, *48*, 5727.
- (89) David, O.; Fargeau-Bellassoued, M.-C.; Lhomme, G. *Tetrahedron Lett.* **2002**, *43*, 3471.
- (90) (a) Hammerschmidt, F.; Hanbauer, M. *J. Org. Chem.* **2000**, *65*, 6121. (b) Majewski, M.; DeCaire, M.; Nowak, P.; Wang, F. *Can. J. Chem.* **2001**, *79*, 1792. (c) Gambacorta, A.; Tofani, D.; Lupattelli, P.; Tafi, A. *Tetrahedron Lett.* **2002**, *43*, 2195. (d) Varie, D. L.; Beck, C.; Borders, S. K.; Brady, M. D.; Cronin, J. S.; Ditsworth, T. K.; Hay, D. A.; Hoard, D. W.; Hoying, R. C.; Linder, R. J.; Miller, R. D.; Moher, E. D.; Remacle, J. R.; Rieck, J. A., III; Anderson, D. D.; Dodson, P. N.; Forst, M. B.; Pierson, D. A.; Turpin, J. A. *Org. Process Res. Devel.* **2007**, *11*, 546.
- (91) Aoki, K.; Koga, K. *Chem. Pharm. Bull.* **2000**, *48*, 571.
- (92) Clayden, J.; Knowles, F. E.; Menet, C. J. *Synlett* **2003**, 1701.
- (93) For recent reviews on organocatalysis, see: (a) Dalco, P. I.; Moisan, L. *Angew. Chem., Int. Ed.* **2004**, *43*, 5138. (b) List, B. *Chem. Commun.* **2006**, 819. (c) Tsogoeva, S. B. *Eur. J. Org. Chem.* **2007**, 1701. (d) Gaunt, M. J.; Johansson, C. C. C.; McNally, A.; Vo, N. T. *Drug Discovery Today* **2007**, *12*, 8. (e) *Chem. Rev.* **2007**, *107* (12), 5413–5883 (special issue on organocatalysis). (f) McMillan, D. W. C. *Nature* **2008**, *455*, 304. (g) Barbas, C. F., III. *Angew. Chem., Int. Ed.* **2008**, *47*, 42. (h) Melchiorre, P.; Marigo, M.; Carlone, A.; Bartoli, G. *Angew. Chem., Int. Ed.* **2008**, *47*, 6138. (i) Dondoni, A.; Massi, A. *Angew. Chem., Int. Ed.* **2008**, *47*, 4638. (j) Bertelsen, S.; Jørgensen, K. A. *Chem. Soc. Rev.* **2009**, *38*, 2178.
- (94) Chimni, S. S.; Mahajan, D. *Tetrahedron: Asymmetry* **2006**, *17*, 2108.
- (95) (a) Gryko, D.; Lipiński, R. *Adv. Synth. Catal.* **2005**, *347*, 1948. (b) Gryko, D.; Lipiński, R. *Eur. J. Org. Chem.* **2006**, 3864. (c) Gryko, D.; Zimmnicka, M.; Lipiński, R. *J. Org. Chem.* **2007**, *72*, 964.
- (96) Raj, M.; Parashari, G. S.; Singh, V. K. *Adv. Synth. Catal.* **2009**, *351*, 1284.
- (97) Hernández-Rodríguez, M.; Avila-Ortiz, C. G.; del Campo, J. M.; Hernández-Romero, D.; Rosales-Hoz, M. J.; Juaristi, E. *Aust. J. Chem.* **2008**, *61*, 364.
- (98) (a) Tsogoeva, S. B.; Yalalov, D. A.; Hateley, M. J.; Weckbecker, C.; Huthmacher, K. *Eur. J. Org. Chem.* **2005**, 4995. (b) Tsogoeva, S. B.; Wei, S. *Chem. Commun.* **2006**, 1451.
- (99) Kelleher, F.; Kelly, S.; Watts, J.; McKee, V. *Tetrahedron* **2010**, *66*, 3525.
- (100) Liu, Y.; Melgar-Fernández, R.; Juaristi, E. *J. Org. Chem.* **2007**, *72*, 1522.
- (101) (a) González-Olvera, R.; Demare, P.; Regla, I.; Juaristi, E. *ARKIVOC* **2008** (vi), 61. (b) See also Bandala, Y.; Melgar-Fernández, R.; Guzmán-Mejía, R.; Olivares-Romero, J. L.; Díaz-Sánchez, B. R.; González-Olvera, R.; Vargas-Caporali, J.; Juaristi, E. *J. Mex. Chem. Soc.* **2009**, *53*, 147.
- (102) Tang, H.; Zhao, G.; Zhou, Z.; Zhou, Q.; Tang, C. *Tetrahedron Lett.* **2006**, *47*, 5717.
- (103) Peris, G.; Miller, S. J. *Org. Lett.* **2008**, *10*, 3049.
- (104) Hansch, M.; Illa, O.; McGarrigle, E. M.; Aggarwal, V. K. *Chem.—Asian J.* **2008**, *3*, 1657.

Trademarks: Raney® (W. R. Grace & Co.).

Keywords: α -Phenylethylamine; chiral ligand; resolving agent; chiral auxiliary; chiral reagent; organocatalyst.



About the Authors

Yamir Bandala was born in 1979 in Perote, México. He studied chemistry at the Centro de Investigación y de Estudios Avanzados (CINVESTAV), México, where he completed his doctoral thesis on the "Synthesis and Conformational Analysis of Cyclic and Open-Chain β - and α/β -Peptides" at the beginning of 2009 in Prof. Juaristi's research group. After one year of postdoctoral work with Prof. Juaristi, he moved in 2010 to a second postdoctoral position in Professor Gerardo Corzo's group at the Instituto de Biotecnología (UNAM), México, where he is working on the development of new β -amino acids and α/β -peptides.

Eusebio Juaristi was born in 1950 in Querétaro, México. He studied chemistry with Prof. E. L. Eliel at the University of North Carolina, Chapel Hill, where he received his Ph.D. degree in 1977. Following postdoctoral stays at UC Berkeley (with A. Streitwieser) and the Syntex Diagnostics Division (Palo Alto, CA), he returned to Mexico where he is now Professor of Chemistry at the Centro de Investigación y de Estudios Avanzados del Instituto Politécnico Nacional. He has also served as Visiting Professor at ETH-Zurich (1985–1986 and 1992–1993) and UC Berkeley (1999–2000). In 1998, he received the National Medal of Science and, in 2006, he became a member of El Colegio Nacional (highest academic honor in México).✉

Chemrus[®] Disposable Filter Funnels

SIGMA-ALDRICH
Labware



Clockwise from top: Chemrus funnels, flask adapter and vial adapter.

For more information, see
sigma-aldrich.com

Convenient and inexpensive, these disposable filter funnels have polypropylene bodies with tapered stems that fit 7/10 joints without seals. The funnels are available in two styles, Buchner and Hirsch, and a choice of 10-micron polyethylene frit, Celite[®], or perforated plate for use with filter paper. Maximum use temperature is 110 °C. Order vacuum adapters, filter paper, vials, and 20–400 to 24–400 vial connector separately below.

Cat. No.	Style	Filter type	Capacity (mL)
Z679798	Buchner	PE frit, 10 micron	18
Z679801	Buchner	PE frit, 10 micron	40
Z679828	Buchner	PE frit, 10 micron	110
Z679852	Buchner	Celite [®] (0.5 g)	18
Z679860	Buchner	Celite [®] (1.5 g)	40
Z679836	Hirsch	Perforated plate	20
Z679844	Hirsch	Perforated plate	60

Accessories

Vacuum adapters

Glass body with PTFE stopcock. All adapters have size 7/10 top joint that fits the disposable filter funnels.

Cat. No.	Outer joint or thread
Z679879	14/20
Z679887	19/22
Z679895	24/25 (fits all 24 joints)
Z679909	29/26 (fits all 29 joints)
Z679917	20–400 vial thread

Filter paper

Cat. No.	For Hirsch funnel (mL)	Diam. (cm)
Z274844	20	1.5
Z274852	60	2.5


Vials, clear glass (vial only)

Cat. No.	
Z7379	40 mL, clear glass, O.D. × H 27.5 mm × 95 mm, thread 24–400

Threaded vial connector

Cat. No.	
Z503193	Thread 20–400 x 24–400, made of glass-filled nylon with inert PTFE inner seals


Chemrus is a registered trademark of Jianjian Cai and Xiaogao Liu, individuals. Celite is a registered trademark of Celite Corp.




BASF's ChiPros® Now Available from Aldrich

Over 75 BASF ChiPros chiral molecules are now available in quantities up to kilogram size, including chiral amines, alcohols, epoxides and acids.


ChiPros® Chiral Amines




726621




727024




727229




726907




726729




727296




726583




727113




727288




726850




726559




726915




726664



726966




726648




727164




727261




726494




726591



727156




726486




727172




727245



726826




727121




726702




726532




727105




726680



726524




727148




726710




726540




726931




726796



726516



727202




726974




726737



726893




726605




726818

ChiPros is a registered trademark of BASF

ChiPros Chiral Amines (Cont'd)



727199




726958




727083




726885




726613




726761




726842




727334




726923




727180




727091



726656




726869



726877




727342




727032


ChiPros Chiral Acids




727075




727067




727016




727318




727326




726990



727008




726982




727350


ChiPros Chiral Alcohols




726753




727237




726672




727059



726788



727210



726567

ChiPros Chiral Epoxides




726508




726834



726699



727253



726745

For more chiral building blocks, please visit aldrich.com/chiralbb

New Fluorinating Reagents



XtalFluor Salts: Convenient, Crystalline Deoxyfluorination Reagents

XtalFluor reagents are crystalline dialkylaminodifluorosulfonium tetrafluoroborate salts. They are useful for the deoxyfluorination of hydroxyl and carbonyl moieties when used in conjunction with a promoter such as DBU, Et₃N·3HF, or Et₃N·2HF.




Advantages of XtalFluor salts as deoxyfluorination reagents:

- Air-stable solids
- Greater thermal stability than DAST or Deoxo-Fluor®
- Broad substrate scope
- Predictable and high chemoselectivity

XtalFluor Reagents



Representative Scope



^a Promoters: DBU, Et₃N·3HF, or Et₃N·2HF

References: (1) Beaulieu, F. et al. *Org. Lett.* **2009**, *11*, 5050. (2) L'Heureux, A. et al. *J. Org. Chem.* **2010**, *75*, 3401

SIGMA-ALDRICH

Research quantities available through Sigma-Aldrich Co.

XtalFluor-E®: Aldrich Cat. No. 719439

XtalFluor-M®: Aldrich Cat. No. 719447

MANCHESTER
ORGANICS

Multi-kilogram quantities available through
Manchester Organics

XtalFluor-E®: Manchester Cat. No. J11026

XtalFluor-M®: Manchester Cat. No. K11027

For more information, visit aldrich.com/xtalfluors

XtalFluor-E and XtalFluor-M are registered trademarks of OmegaChem, Inc. Deoxo-Fluor is a registered trademark of Air Products & Chemicals, Inc.

2-Iodoxybenzenesulfonic Acid (IBS) Catalyzed Oxidation of Alcohols



Prof. Muhammet Uyanik

Prof. Kazuaki Ishihara

Muhammet Uyanik and Kazuaki Ishihara*
 Graduate School of Engineering
 Nagoya University
 Furo-cho, Chikusa
 Nagoya, 464-8603, Japan
 Email: ishihara@cc.nagoya-u.ac.jp

Outline

1. Introduction
2. Hypervalent Iodine Catalysts for Alcohol Oxidation
3. 2-Iodoxybenzenesulfonic Acid (IBS)
 - 3.1. Related Hypervalent Iodine Reagents
 - 3.1.1. HMBI and λ^3 -Iodane Derivatives
 - 3.1.2. IBS and λ^5 -Iodane Derivatives
 - 3.2. Discovery of IBS as a Catalyst for Alcohol Oxidation with Oxone[®]
 - 3.3. Large-Scale Oxidations
 - 3.4. Application to the Oxidative Rearrangement of Tertiary Allylic Alcohols
 - 3.5. Theoretical Calculations and Reaction Mechanism
4. Conclusion
5. Acknowledgement
6. References and Notes

1. Introduction

The oxidation of alcohols to the corresponding carbonyl compounds is one of the most important transformations in synthetic organic chemistry.¹ Over the past two decades, hypervalent iodine compounds have received a great deal of attention due to their mild and chemoselective oxidizing properties and to the fact that, unlike toxic heavy-metal-based reagents, they are environmentally benign.² We recently published a review on the oxidation of alcohols with stoichiometric or catalytic hypervalent iodine compounds.³ The present review focuses on the discovery and development of 2-iodoxybenzenesulfonic acid (IBS, **1**) and related hypervalent iodine compounds. It also covers advances in hypervalent iodine catalyzed alcohol oxidations and related transformations reported between 2005 and 2010.

2. Hypervalent Iodine Catalysts for Alcohol Oxidation

To date, various hypervalent iodine (λ^3 - and λ^5 -iodane) reagents have been developed as oxidants for alcohols.³ 2-Iodoxybenzoic acid (IBX, **2**) and Dess–Martin periodinane (DMP, **3**) are the best known and most commonly used hypervalent iodine reagents (Figure 1). Dess and Martin used **2** as a precursor for **3**.⁴ In the early studies, **2** was synthesized from 2-iodobenzoic acid (**4**) and potassium bromate (KBrO₃) in aqueous sulfuric acid.^{4a} Although **2** was later reported to be explosive,⁵ Dess and Martin subsequently speculated that some bromate or other impurity may

have been included in samples that were found to be explosive, but that IBX (**2**) itself should be non-explosive.^{4b} The first use of **2** for alcohol oxidation in DMSO was reported in 1994 by Frigerio and Santagostino.⁶ The simple, one-step preparation of **2** from **4** with Oxone[®], an environmentally safe oxidant, has made **2** a popular reagent (eq 1).⁷ To date, **2** has been employed as a powerful and selective oxidant that mediates a variety of transformations such as the oxidation of alcohols, phenols, and amines; the dehydrogenation of ketones, aldehydes, and N-heterocycles; and the oxidative cleavage of dithioacetals.^{2,3} Several research groups have attempted to improve on **2** by structurally modifying it, or by developing polymer-supported analogues.³ Additionally, Quideau's group reported a stabilized formulation of IBX, SIBX, containing benzoic acid and isophthalic acid, which offered some advantages, such as safety and ease of workup.⁸

In 2005 and 2006, Vinod⁹ and Giannis¹⁰ independently reported the oxidation of alcohols catalyzed by **2**—generated in situ from **4** or 2-iodosobenzoic acid (IBA)—in the presence of Oxone[®] as co-oxidant (Scheme 1). Vinod's group employed 20–40 mol % of **4** in a water–acetonitrile biphasic solvent system, in which primary and secondary alcohols were oxidized to carboxylic acids and ketones, respectively (eq 2).⁹ In contrast, Giannis's group utilized a water–ethyl acetate biphasic solvent system in the presence of 10 mol % each of **4** and tetrabutylammonium hydrogen sulfate [(*n*-Bu)₄NHSO₄] as a phase-transfer catalyst. Under these conditions, primary benzylic alcohols were oxidized to the corresponding aldehydes, which did not undergo further oxidation (eq 3).¹⁰ These two reports stated that it was not necessary to isolate beforehand the hypervalent iodine compounds, which are potentially shock-sensitive explosive oxidants.

Page and co-workers demonstrated that several primary and secondary alcohols can be oxidized to the respective aldehydes and ketones under reflux conditions in acetonitrile or dichloroethane in the presence of tetraphenylphosphonium monoperoxydisulfate (TPPP, Ph₄P⁺HSO₅⁻) and a catalytic amount of **4** (eq 4).¹¹ TPPP was derived from Oxone[®] by simple counterion exchange with tetraphenylphosphonium chloride. This catalytic system enables the oxidation of primary alcohols to the corresponding aldehydes without further oxidation to the carboxylic acids.

The selective oxidation of benzylic C–H bonds to the corresponding carbonyl functionalities has been achieved using a catalytic amount of **4** and Oxone[®] in aqueous acetonitrile (Scheme 2).¹² The authors hypothesized that the active









Figure 1. IBS (1), IBX (2), and DMP (3).






Scheme 1. The Oxidation of Alcohols Catalyzed by IBX (2) Generated in Situ. (Ref. 9,10)




hypervalent iodine oxidant generated in situ might not be IBX (2), but a soluble derivative of IBX, 5, incorporating a peroxy ligand (KHSO_5). This species is believed to oxidize a benzylic C-H bond via a single-electron-transfer (SET) mechanism.¹²

Yakura's group reported that *para*-alkoxyphenols and *para*-arylphenols are oxidized in excellent yields to the corresponding *para*-quinones and *para*-quinols, respectively, using catalytic amounts of 4-iodophenoxyacetic acid (6) with Oxone[®] as a co-oxidant in aqueous acetonitrile (eq 5).¹³

An efficient, catalytic aerobic oxidation of primary and secondary alcohols to the corresponding aldehydes and ketones has been carried out by Liu and co-workers by using a mixture of iodoxybenzene (PhIO_2 , 7, 1 mol %), Br_2 (2 mol %), and NaNO_2 (1 mol %) in water (Scheme 3).¹⁴ The proposed reaction mechanism includes three redox cycles. In the first redox cycle, 7 is the active species that oxidizes the alcohol to the corresponding carbonyl compound, and is reduced to dihydroxyiodobenzene ($\text{PhI}(\text{OH})_2$). In the second cycle, $\text{PhI}(\text{OH})_2$ is reoxidized to 7 with Br_2 , which is reduced to HBr . In the third and final cycle, the oxidation of NO with O_2 produces NO_2 , which reoxidizes HBr to Br_2 . However, we were not able to effect the oxidation of alcohols under Liu's

Scheme 2. The Oxidation of Benzylic C-H Bonds with Oxone[®] Catalyzed by 5 Generated in Situ. (Ref. 13)

$\text{R}^1 = \text{H, Me, Ph}$
 $\text{R}^2 = \text{Ph, Bz, aryl, arylalkyl, alkyl}$

Scheme 3. Liu's Aerobic Oxidation of Alcohols Catalyzed by PhIO_2 (7). (Ref. 14)

conditions,¹⁵ which led us to suggest, on the basis of several control experiments, that the actual oxidant of the alcohols in this case is Br₂ rather than PhIO₂.¹⁵


Moreover, our group disclosed that the simple and environmentally benign catalytic system consisting of HBr and NaNO₂ is very effective for the aerobic oxidation of alcohols (**Scheme 4**).¹⁵ Primary benzylic alcohols are selectively oxidized to the corresponding aldehydes in acetonitrile under a balloon pressure of O₂. Secondary alcohols are oxidized to the corresponding ketones in the presence of a small amount of water. Furthermore, the aerobic oxidation of alcohols can also be achieved under a balloon pressure of air instead of pure O₂ with the HBr–NaNO₂–HNO₃ catalytic system.¹⁵

Li and co-workers developed an effective system for the oxidation of alcohols under an atmosphere of oxygen, without the need for any additional solvent or transition-metal catalyst, by using catalytic amounts of phenyliodine diacetate (PIDA (**8**))–TEMPO–KNO₂ (**Scheme 5**).¹⁶


In 2009, Yusubov, Zagulyaeva, and Zhdankin reported an efficient tandem catalytic system, based on a Ru(V)–catalyzed reoxidation of iodosobenzene (PhIO, **9**) to PhIO₂ (**7**), for the oxidation of alcohols and hydrocarbons to carbonyl compounds by employing stoichiometric amounts of Oxone[®] at room temperature (**Scheme 6**).¹⁷

Very recently, Nemykin, Zhdankin, and co-workers reported the room temperature, Fe(III)–porphyrin-catalyzed oxygenation of the anthracene ring to the anthraquinone system in the presence of a substoichiometric amount of iodobenzene and an excess of Oxone[®] (**Scheme 7**).¹⁸ The proposed reaction mechanism includes two catalytic redox cycles. Accordingly, the active λ³-iodane species, generated in solution by treatment of iodobenzene with Oxone[®], is responsible for the oxidation of the Fe(III)–porphyrin to the oxo-Fe(IV)–porphyrin complex, which then acts as the actual oxygenating agent.

The same year, a catalytic, room-temperature oxidation of alcohols was effected with *meta*-chloroperbenzoic acid (*m*-CPBA) or potassium peroxodisulfate (K₂S₂O₈) in the presence of a catalytic amount of iodobenzene and *N*-hydroxyphthalimide (NHPI) or TEMPO in aqueous acetonitrile (**Scheme 8**).¹⁹ In both cases, iodobenzene was oxidized by *m*-CPBA or K₂S₂O₈ in situ to a λ³-iodane species, a reoxidant of NHPI or TEMPO to phthalimide-*N*-oxyl radical (PINO) or TEMPO, which then oxidizes the alcohols.



Scheme 4. Ishihara's Bromine-Catalyzed Aerobic Oxidation of Alcohols. (Ref. 15)




Scheme 4. Ishihara's Bromine-Catalyzed Aerobic Oxidation of Alcohols. (Ref. 15)

3. 2-Iodoxybenzenesulfonic Acid (IBS)


3.1. Related Hypervalent Iodine Reagents

3.1.1. HMBI and λ³-Iodane Derivatives



In 1993, Koser and co-workers first reported that the oxidation of 2-iodo-5-methylbenzenesulfonic acid (**10**) with peracetic acid leads to 1*H*-1-hydroxy-5-methyl-1,2,3-benziodoxathiole 3,3-dioxide (HMBI, **11**) (**Scheme 9**).²⁰ In 2007, Justik reported that the sodium salt (**10**·Na) could be directly oxidized to **11** through in situ protonation with a small amount of concentrated sulfuric acid (see **Scheme 9**).²¹ HMBI (**11**) was utilized to




Scheme 5. TEMPO–PIDA (**8**) Co-catalyzed Aerobic Oxidation of Alcohols. (Ref. 16)




Scheme 6. Ru(III)–PhIO₂ (**7**) Co-catalyzed Oxidations with Oxone[®].




Scheme 7. Fe(III)–Porphyrin–Iodine(III) Co-catalyzed Oxidation of the Anthracene Nucleus. (Ref. 18)




Scheme 8. Nitroxyl Radical-Iodine(III) Co-catalyzed Oxidation of Alcohols with *m*-CPBA or $K_2S_2O_8$. (Ref. 19)




Scheme 9. Koser's and Justik's Syntheses of HMBS. (Ref. 20,21)



eq 6 (Ref. 21)



eq 7 (Ref. 22)



eq 8 (Ref. 23)

convert alkanophenones to alkyl esters of 2-arylalkanoic acids by oxidative rearrangement (eq 6).²¹ The use of **11** offers several advantages over other similar oxidizing agents, including facile workup of the reaction mixture and recovery of the reduced iodine reagent.


In the same 1993 report, Koser also disclosed the synthesis of 1*H*-1-(1-alkynyl)-5-methyl-1,2,3-benziodoxathiole 3,3-dioxides (AMBIs, **12**) from **11** and the corresponding terminal alkynes in low-to-good yields (eq 7).²⁰ Such AMBIs were later employed by Ishiwa and Togo to prepare thiazoles from the corresponding thioamides (eq 8).²² The reduced iodo compound **10•K** was recovered quantitatively by simple filtration of the reaction mixture, and could be converted to **12** to be reused for the preparation of thiazoles.²²

In a subsequent report, Justik's group disclosed that the reaction of 1*H*-1-hydroxy-1,2,3-benziodoxathiole 3,3-dioxide (HBI, **13**) with arenes produces 2-[(aryliodonio)benzenesulfonates (**14**) in moderate-to-high yields (Scheme 10).²³ HBI (**13**) was prepared from sodium 2-iodobenzenesulfonate (**15•Na**) similarly to the way its methyl analogue, **11**, was prepared (see Scheme 9). The zwitterionic sulfonates, **14**, were employed in the synthesis of 2-sulfonyloxyarenes and as potential recyclable aryl-transfer reagents in transition-metal-catalyzed cross-coupling reactions (see Scheme 10).²³

3.1.2. IBS and λ^5 -Iodane Derivatives

Zhdankin and co-workers first reported in 2006 the preparation and full characterization of IBS (**1**), a λ^5 -iodane and a thia analogue of IBX (**2**).²⁴ IBS (**1**)—which is the cyclic tautomeric form of 1*H*-1-hydroxy-1,2,3-benziodoxathiole 1,3,3-trioxide—was synthesized by two different methods: (i) the direct oxidation of 2-iodobenzenesulfonic acid (**15**) with Oxone[®] (Method A) or (ii) by hydrolysis of the methyl ester of IBS, **16**. (Method B)²⁵ (Scheme 11). Method A results in a low-purity IBS, while Method B yields IBS of a high purity. IBS (**1**) is thermally unstable and highly reactive toward organic solvents such as acetonitrile, DMSO, and methanol.²⁴ Because of this instability (reductive decomposition to the corresponding λ^3 -iodane) and reactivity, these researchers did not investigate its oxidative ability.²⁴

The pseudocyclic λ^5 -iodanes, IBS ester (**16**)²⁵ and amides (**17**),²⁶ have been prepared by oxidation of the corresponding 2-iodobenzenesulfonic ester and 2-iodobenzenesulfonamide with dimethyldioxirane (Scheme 12).²⁶ The starting material



Scheme 10. Preparation and Application of 2-[(Aryl)iodonio]benzenesulfonates (**14**). (Ref. 23)


for these monovalent iodines, 2-iodobenzenesulfonyl chloride, was synthesized from **15**•Na (see Scheme 12).²⁵ IBS amides **17** can selectively oxidize benzyl alcohols to aldehydes.²⁶ While IBS esters **16** could not oxidize alcohols, they are useful for the oxidation of other organic functional groups, such as sulfides and secondary amines to the respective sulfoxides and imines (Scheme 13).²⁵

3.2. Discovery of IBS as a Catalyst for Alcohol Oxidation with Oxone®

To develop a more powerful hypervalent iodine catalyst for selective alcohol oxidation, we became interested in modifying the arene moiety of IBX (**2**) (Figure 2).²⁷ We discovered that electron-donating groups in the benzene ring, as in **18** and **19**, resulted in λ^5 -iodanes that were superior to **2** as catalysts for the oxidation of alcohols with Oxone® under nonaqueous conditions, even though Oxone® was almost insoluble in most organic solvents. In contrast, electron-withdrawing groups, as in **20**, resulted in λ^5 -iodanes with lower reactivity than **2**. The stability of these IBX catalysts was strongly influenced by solvents and the substituents; for instance, the 3-methyl- and 4,5-dimethoxy-substituted IBXs, **21** and **22**, decomposed under both aqueous and nonaqueous oxidation conditions.²⁷

We then turned our attention to IBS (**1**) for the purpose of developing an even more powerful catalyst for alcohol oxidation. Our reasoning was that the Lewis acidity of the iodine(V) atom in **1** would be higher than that in **2** due to the strong electron-withdrawing sulfonate group in **1**.^{27,28} Although **1** is not stable enough to isolate in pure form and examine in stoichiometric reactions, it can be prepared in situ in the presence of the alcohol substrate from **15**, **15**•Na, or **15**•K and Oxone®. We compared the

catalytic efficiencies of in situ generated **1** and **2** in the oxidation of 5-nonanol (Figure 3, **1** in red and **2** in blue).²⁷ In nonaqueous nitromethane and acetonitrile, **1** was superior to **2** (red square and circle vs blue open-square and open circle; in nitromethane, TOF $\geq 50 \text{ h}^{-1}$ vs 16 h^{-1} , respectively). In sharp contrast, **2** was superior to **1** in aqueous media (red diamond vs blue open diamond).



Scheme 13. Facile Oxidations with IBS Ester **16** and IBS Amide **17**. (Ref. 25,26)






Figure 2. Catalytic Efficiencies of Arene-Substituted IBXs. (Ref. 27)



Scheme 11. Zhdankin's First Reported Synthesis of IBS (**1**). (Ref. 24)



Scheme 12. Synthesis of Pseudocyclic λ^5 -Iodanes **16** and **17**. (Ref. 25)




Figure 3. The Catalytic Efficiencies, in Aqueous or Nonaqueous Media, of IBS (**1**) and IBX (**2**) Generated in Situ from **15**•Na and **4**, Respectively. (Ref. 27)

Moreover, **1** was sufficiently active in less polar, but more environmentally benign, ethyl acetate; while the **2**-catalyzed oxidation of alcohols was very slow in nonaqueous solvents such as ethyl acetate or acetonitrile. We investigated the substituent effect on the IBS-catalyzed oxidation of alcohols: IBS derivatives substituted with electron-donating groups, as in 5-Me-IBS (**23**)






Figure 4. IBS Derivatives Substituted with Electron-Donating Groups. (Ref. 27)



Scheme 14. The IBS-Catalyzed Selective Oxidation of **25** to **26** or **27**. (Ref. 27)


and 4,5-Me₂-IBS (**24**) (**Figure 4**), were superior to **1**, although the differences in the catalytic efficiencies were not large. The oxidation rates of reactions catalyzed by IBS were further accelerated by the use of powdered Oxone[®] due to its increased surface area. Additionally, vigorous stirring of the heterogeneous reaction mixture is essential for the efficient grinding of Oxone[®] and efficient alcohol oxidation. In a nonaqueous solvent system, the desired carbonyl products were obtained in nearly pure form by simple filtration of most wastes derived from Oxone[®] and washing with water to remove catalyst derivatives.²⁷

Structurally diverse primary and secondary alcohols were oxidized with **1** or **23** under optimized conditions (eq 9).²⁷ Not only primary, α,β-unsaturated alcohols—such as allylic, propargylic, and benzylic alcohols—but also aliphatic alcohols were selectively oxidized to the corresponding aldehydes and carboxylic acids in excellent yields by controlling the amount of Oxone[®] added in the presence of 15•Na or 10•K. In particular, Methods B and C were effective for the selective oxidation of acid-sensitive alcohols and primary aliphatic alcohols to the corresponding aldehydes, respectively. This protocol was also applied to the chemoselective oxidation of alcohols bearing several functional or protective groups such as silyloxy, benzyloxy, ketal, alkenyl, alkynyl, and halo groups (see eq 9).²⁷

In 2000, Nicolaou and co-workers reported the direct oxidative dehydrogenation, via a single-electron-transfer (SET) mechanism, of saturated alcohols and carbonyl compounds by using stoichiometric amounts of IBX (**2**) in DMSO.^{29,30} During the course of our research on the oxidation of cycloalkanols, we found that the selective oxidation of 4-*tert*-butylcyclohexanol (**25**) to 4-*tert*-butylcyclohexanone (**26**) and the subsequent oxidation of **26** to 4-*tert*-butyl-2-cyclohexenone (**27**) and 5-*tert*-butyl-2-oxocyclohexanone (**28**) proceeded in excellent yields by controlling the amounts of 15•Na and Oxone[®] (**Scheme 14**).²⁷ Five- and six-membered cycloalkanols were transformed into the corresponding enones in good yields (eq 10).²⁷

3.3. Large-Scale Oxidations

The selective oxidation of primary alcohols to the corresponding aldehydes or carboxylic acids would be a powerful tool in organic synthesis. The transition-metal- or nitroxyl-radical-catalyzed oxidation of alcohols to ketones or aldehydes has attracted a lot of attention because aqueous H₂O₂ or gaseous O₂ can be used as a stoichiometric oxidant.¹ However, it is technically difficult to control the amount of gaseous O₂ added as an oxidant. Moreover, while aqueous H₂O₂ and gaseous O₂ are often concentrated under evaporation and high pressure, respectively, to increase their reactivity, such treatments are dangerous because these materials can be explosive. Although aqueous H₂O₂ and gaseous O₂ are more atom-economical than Oxone[®], the latter offers several advantages for selective large-scale oxidations such as stability, simple handling, controllable addition, and nontoxic nature.



eq 10 (Ref. 27)

4-Bromobenzyl alcohol (6 g) was selectively oxidized to the corresponding aldehyde and carboxylic acid in excellent yield by controlling the amount of Oxone[®] added in the presence of 10•K (1 mol %): 0.65 equiv was used for selective oxidation to the aldehyde, while 1.3 equiv was employed for oxidation to 4-bromobenzoic acid (eq 11).³¹


3.4. Application to the Oxidative Rearrangement of Tertiary Allylic Alcohols

The oxidative rearrangement of tertiary allylic alcohols to β -disubstituted α,β -unsaturated ketones or aldehydes using oxochromium(VI)-based reagents (Collins reagent, PCC, PDC) has been widely used in synthetic organic chemistry.³² In 2004, Iwabuchi and co-workers reported that IBX (**2**) could be used instead of hazardous Cr(VI) for the oxidative rearrangement of tertiary allylic alcohols.³³ In 2008, the same research group reported that TEMPO-derived oxoammonium salts (TEMPO⁺ BF₄⁻ and TEMPO⁺ SbF₆⁻) were more effective as stoichiometric reagents for this transformation of acyclic tertiary allylic alcohols in acetonitrile.³⁴ In that same year, Iwabuchi³⁵ and Vatèle³⁶ independently reported the first *catalytic* oxidative rearrangement of tertiary allylic alcohols. Iwabuchi's group utilized catalytic amounts of TEMPO with NaIO₄-SiO₂ as co-oxidant in dichloromethane to convert several cyclic and acyclic tertiary allylic alcohols into the corresponding enones in good yields.³⁵ In contrast, Vatèle developed a Lewis acid [Bi(OTf)₃ or Re₂O₇] promoted oxidative rearrangement using catalytic amounts of TEMPO with PhIO (**9**) as co-oxidant.³⁶ As part of our continuing interest in the use of IBS/Oxone[®] catalytic oxidation systems in organic synthesis, we reported the development of an oxidative rearrangement of tertiary allylic alcohols to enones with powdered Oxone[®] promoted by catalytic quantities of 10•Na (Scheme 15).³⁷ 5-Me-IBS (**23**) is generated in situ and serves as the actual catalyst for the oxidation. Interestingly, **1** was less effective than **23** as catalyst in this case. The addition of inorganic bases, such as K₂CO₃, and a phase-transfer catalyst, such as tetra-*n*-butylammonium hydrogen sulfate [(*n*-Bu)₄NHSO₄], extended the substrate scope for oxidative rearrangement reactions (see Scheme 15). Cyclic and acyclic substrates gave the corresponding enones in moderate to high yields. Notably, sterically demanding steroid alcohol **29** was converted into the desired enone **30** in 69% yield (see Scheme 15). In contrast, Iwabuchi's group reported that **29** was rearranged to allylic alcohol **31**, which was not oxidized to **30** under TEMPO-mediated conditions, due to considerable steric hindrance.^{34,35} Our protocol should be recognized as a practical method for the oxidative rearrangement of tertiary alcohols, since it does not require any toxic, dangerous, or expensive reagents.³⁷



3.5. Theoretical Calculations and Reaction Mechanism

The oxidation of alcohols with IBS consists of two steps, which are essentially identical to those of the IBX oxidation proposed by Santagostino:³⁸ a fast pre-equilibrium step and a rate-determining disproportionation step. The catalytic cycle of **1**, which is prepared in situ from **15**, can be accomplished by regenerating **1** through the oxidation of **15** with Oxone[®] (Scheme 16).


Based on the theoretical study by Su and Goddard on the oxidation of alcohols with IBX (**2**),³⁹ we also determined that the twisting of alkoxyperiodinanes **32-A** to the intermediates **32-B** was the rate-limiting step (Scheme 17).²⁷ Interestingly, **32** has a much lower twisting barrier than **33** (**33-A**→**33-B**: 10.3 kcal/mol vs **32-A**→**32-B**: 6.5 kcal/mol).²⁷ The I(V)-OCO and I(V)-OSO₂ distances in **33-A** and **32-A** are correlated with the twisting barriers: 2.252 Å for **32-A** > 2.193 Å for **33-A**. Based



eq 11 (Ref. 31)



Scheme 15. IBS-Catalyzed Oxidative Rearrangement of Tertiary Allylic Alcohols. (Ref. 37)



Scheme 16. Proposed Catalytic Cycle for IBS- and IBX-Catalyzed Oxidations. (Ref. 27)

on theoretical calculations, we assumed that the relatively ionic character of the intramolecular hypervalent iodine–OSO₂ bond of IBS might lower the twisting barrier of the alkoxyperiodinane intermediate **32-A**. We are able to confirm that Goddard's *hypervalent twisting* would be the rate-determining step for the *stoichiometric oxidation* of alcohols with not only IBX (**2**) but also IBS (**1**). In contrast, it was reasonable to assume that the rate-determining step of *I(V)-catalyzed oxidations* might be the regeneration of I(V) species because: (i) the catalytic oxidation of alcohols was accelerated with powdered Oxone[®], and (ii) control experiments indicated that the 5-methyl substituent of 5-Me-IBX had no influence on the oxidation of alcohols, while the 5-methyl substituent of 5-Me-IBS (**23**) accelerated the oxidation from I(III) to I(V).²⁷

4. Conclusion


Over the past two decades, IBX (**2**) and other hypervalent iodine compounds have received considerable attention because of their mild and chemoselective oxidizing properties and because they are environmentally benign in contrast to the well-known but highly toxic heavy metal oxidants.² However, the stoichiometric use of **2** has been limited because it is potentially shock-sensitive and weakly soluble in common organic solvents.² Over the past five years, several research groups have reported alcohol oxidation reactions catalyzed by in situ generated IBX or related hypervalent iodines. These developments have been highlighted in this review.^{9–19} Although IBS (**1**) is not stable enough to isolate in pure form and examine in stoichiometric reactions, we have reported that IBS's **1** and **23** can be prepared in situ from pre-IBS **15**•Na and **10**•K and Oxone[®] in the presence of an alcohol substrate. We have also disclosed that **1** and **23** show greater catalytic activity than IBX (**2**) under nonaqueous conditions.²⁷ Consequently, we have developed a highly efficient and chemoselective oxidation of various alcohols to carbonyl compounds with powdered Oxone[®] in the presence of catalytic amounts of IBS. Our findings could lead to new possibilities in hypervalent-iodine-catalyzed oxidative transformations.

5. Acknowledgement

We gratefully acknowledge support from Nissan Chemical Industries, Ltd. Financial support for this project was partially provided by JSPS.KAKENHI (20245022), the GSC Project of METI, the Toray Science Foundation, Kyowa Hakko Chemical Co., Ltd. (The 1st Seeds Contest), and the Global COE Program of MEXT.

6. References and Notes

- (1) Tojo, G.; Fernández, M. *Oxidation of Primary Alcohols to Carboxylic Acids*; Basic Reactions in Organic Synthesis Series, Vol. 2; Tojo, G., Series Ed.; Springer: Berlin, 2007. (b) Tojo, G.; Fernández, M. *Oxidation of Alcohols to Aldehydes and Ketones*; Basic Reactions in Organic Synthesis Series, Vol. 1; Tojo, G., Series Ed.; Springer: Berlin, 2006. (c) *Modern Oxidation Methods*; Bäckwall, J.-E., Ed.; Wiley-VCH: New York, 2004. (d) Ley, S. V. In *Comprehensive Organic Synthesis*, 3rd ed.; Pergamon: Oxford, 1999; Vol. 7, Chap. 2, pp 251–327.
- (2) (a) *Hypervalent Iodine Chemistry*; Wirth, T., Ed.; Topics in Current Chemistry Series, Vol. 224; Springer: Berlin, 2003. (b) Wirth, T. In *Organic Synthesis Highlights V*; Schmalz, H.-G., Wirth, T., Eds.; Wiley-VCH: Weinheim, 2003, 144–150. (c) Zhdankin, V. V.; Stang, P. J. *Chem. Rev.* **2002**, *102*, 2523. (d) Tohma, H.; Kita, Y. *Adv. Synth. Catal.* **2004**, *346*, 111. (e) Richardson, R. D.; Wirth, T. *Angew. Chem., Int. Ed.* **2006**, *45*, 4402. (f) Ladziata, U.; Zhdankin, V. V. *ARKIVOC* **2006** (ix), 26. (g) Ochiai, M. *Chem. Rec.* **2007**, *7*, 12. (h) Ciufolini, M. A.; Braun, N. A.; Canesi, S.; Ousmer, M.; Chang, J.; Chai, D. *Synthesis* **2007**, 3759. (i) Ladziata, U.; Zhdankin, V. V. *Synlett* **2007**, 527. (j) Quideau, S.; Pouységu, L.; Deffieux, D. *Synlett* **2008**, 467. (k) Zhdankin, V. V.; Stang, P. J. *Chem. Rev.* **2008**, *108*, 5299. (l) Ochiai, M.; Miyamoto, K. *Eur. J. Org. Chem.* **2008**, 4229. (m) Dohi, T.; Kita, Y. *Chem. Commun.* **2009**, 2073. (n) Pouységu, L.; Deffieux, D.; Quideau, S. *Tetrahedron* **2010**, *66*, 2235.
- (3) Uyanik, M.; Ishihara, K. *Chem. Commun.* **2009**, 2086.
- (4) (a) Dess, D. B.; Martin, J. C. *J. Org. Chem.* **1983**, *48*, 4155. (b) Dess, D. B.; Martin, J. C. *J. Am. Chem. Soc.* **1991**, *113*, 7277.
- (5) Plumb, J. B.; Harper, D. J. *Chem. Eng. News.* **1990**, July 16, 3.
- (6) (a) Frigerio, M.; Santagostino, M. *Tetrahedron Lett.* **1994**, *35*, 8019. (b) Frigerio, M.; Santagostino, M.; Sputore, S.; Palmisano, G. *J. Org. Chem.* **1995**, *60*, 7272.
- (7) Frigerio, M.; Santagostino, M.; Sputore, S. *J. Org. Chem.* **1999**, *64*, 4537.
- (8) (a) Ozanne-Beaudenon, A.; Quideau, S. *Tetrahedron Lett.* **2006**, *47*, 5869. (b) Ozanne, A.; Pouységu, L.; Depernet, D.; François, B.; Quideau, S. *Org. Lett.* **2003**, *5*, 2903.
- (9) Thottumkara, A. P.; Bowsheer, M. S.; Vinod, T. K. *Org. Lett.* **2005**, *7*, 2933.
- (10) Schulze, A.; Giannis, A. *Synthesis* **2006**, 257.
- (11) Page, P. C. B.; Appleby, L. F.; Buckley, B. R.; Allin, S. M.; McKenzie, M. J. *Synlett* **2007**, 1565.
- (12) Ojha, L. R.; Kudugunti, S.; Maddukuri, P. P.; Kommareddy, A.; Gunna, M. R.; Dokuparthi, P.; Gottam, H. B.; Botha, K. K.; Parapati, D. R.; Vinod, T. K. *Synlett* **2009**, 117.
- (13) (a) Yakura, T.; Konishi, T. *Synlett* **2007**, 765. (b) Yakura, T.; Yamauchi, Y.; Tian, Y.; Omoto, M. *Chem. Pharm. Bull.* **2008**, *56*, 1632. (c) Yakura, T.; Tian, Y.; Yamauchi, Y.; Omoto, M.; Konishi, T. *Chem. Pharm. Bull.* **2009**, *57*, 252. (d) Yakura, T.; Omoto, M. *Chem. Pharm. Bull.* **2009**, *57*, 643.
- (14) Mu, R.; Liu, Z.; Yang, Z.; Liu, Z.; Wu, L.; Liu, Z.-L. *Adv. Synth. Catal.* **2005**, *347*, 1333.
- (15) Uyanik, M.; Fukatsu, R.; Ishihara, K. *Chem.—Asian J.* **2010**, *5*, 456.
- (16) Herreras, C. I.; Zhang, T. Y.; Li, C.-J. *Tetrahedron Lett.* **2006**, *47*, 13.



Scheme 17. Goddard's Hypervalent Twisting Mechanism of I(V)-Mediated Alcohol Oxidation. (Ref. 27,39)


- (17) Yusubov, M. S.; Zagulyaeva, A. A.; Zhdankin, V. V. *Chem.—Eur. J.* **2009**, *15*, 11091.
- (18) Yoshimura, A.; Neu, H. M.; Nemykin, V. N.; Zhdankin, V. V. *Adv. Synth. Catal.* **2010**, *352*, 1455.
- (19) (a) Zhu, C.; Ji, L.; Zhang, Q.; Wei, Y. *Can. J. Chem.* **2010**, *88*, 362. (b) Zhu, C.; Ji, L.; Wei, Y. *Monatsh. Chem.* **2010**, *141*, 327.
- (20) Koser, G. F.; Sun, G.; Porter, C. W.; Youngs, W. J. *J. Org. Chem.* **1993**, *58*, 7310.
- (21) Justik, M. W. *Tetrahedron Lett.* **2007**, *48*, 3003.
- (22) Ishiwata, Y.; Togo, H. *Synlett* **2008**, 2637.
- (23) Justik, M. W.; Protasiewicz, J. D.; Updegraff, J. B. *Tetrahedron Lett.* **2009**, *50*, 6072.
- (24) Kuposov, A. Y.; Litvinov, D. N.; Zhdankin, V. V.; Ferguson, M. J.; McDonald, R.; Tykwinski, R. R. *Eur. J. Org. Chem.* **2006**, 4791.
- (25) Zhdankin, V. V.; Goncharenko, R. N.; Litvinov, D. N.; Kuposov, A. Y. *ARKIVOC* **2005** (iv), 8.
- (26) Kuposov, A. Y.; Litvinov, D. N.; Zhdankin, V. V. *Tetrahedron Lett.* **2004**, *45*, 2719.
- (27) Uyanik, M.; Akakura, M.; Ishihara, K. *J. Am. Chem. Soc.* **2009**, *131*, 251.
- (28) Many examples exist of the increase in the oxidation power of hypervalent aryl- λ^3 -iodanes following activation by strong Lewis acids. For a review, see reference 2g.
- (29) (a) Nicolaou, K. C.; Zhong, Y.-L.; Baran, P. S. *J. Am. Chem. Soc.* **2000**, *122*, 7596. (b) Nicolaou, K. C.; Montagnon, T.; Baran, P. S.; Zhong, Y.-L. *J. Am. Chem. Soc.* **2002**, *124*, 2245.
- (30) For the synthesis of α,β -unsaturated carbonyl compounds, see: (a) Larock, R. C. *Comprehensive Organic Transformations*; Wiley: New York, 1999; pp 251–256. (b) Buckle, D. R.; Pinto, I. L. In *Comprehensive Organic Synthesis*; Trost, B. M., Ed.; Pergamon: Oxford, 1991; Vol. 7, pp 119–146.
- (31) Uyanik, M.; Ishihara, K. *Org. Synth.*, submitted for publication, 2010
- (32) (a) Luzzio, F. A. *Org. React.* **1998**, *53*, 1. (b) Wietzerbin, K.; Bernadou, J.; Meunier, B. *Eur. J. Inorg. Chem.* **2000**, 1391.
- (33) Shibuya, M.; Ito, S.; Takahashi, M.; Iwabuchi, Y. *Org. Lett.* **2004**, *6*, 4303.
- (34) Shibuya, M.; Tomizawa, M.; Iwabuchi, Y. *J. Org. Chem.* **2008**, *73*, 4750.
- (35) Shibuya, M.; Tomizawa, M.; Iwabuchi, Y. *Org. Lett.* **2008**, *10*, 4715.
- (36) Vatele, J.-M. *Synlett* **2008**, 1785.
- (37) Uyanik, M.; Fukatsu, R.; Ishihara, K. *Org. Lett.* **2009**, *11*, 3470.
- (38) De Munari, S.; Frigerio, M.; Santagostino, M. *J. Org. Chem.* **1996**, *61*, 9272.
- (39) Su, J. T.; Goddard, W. A., III. *J. Am. Chem. Soc.* **2005**, *127*, 14146.

Trademarks: Oxone® (E. I. du Pont de Nemours and Company, Inc.).

Keywords: λ^5 -Iodanes; IBS; Oxone®; alcohol oxidation; chemoselective oxidation; catalytic oxidation.


About the Authors

Muhammet Uyanik was born in 1981 in Samsun, Turkey, and received his Ph.D. degree in 2007 from Nagoya University under the direction of Professor Kazuaki Ishihara. He was appointed Assistant Professor at Nagoya University in 2007. In 2009, Dr. Uyanik received the Shionogi Award in Synthetic Organic Chemistry, Japan, and, in 2010, an Incentive Award for Young Scientists from the Tokai Branch of The Society of Synthetic Organic Chemistry, Japan. His research interests include oxidation reactions and asymmetric catalysis.

Kazuaki Ishihara was born in 1963 in Aichi, Japan, and received his Ph.D. degree in 1991 from Nagoya University under the direction of Professor Hisashi Yamamoto. He had the opportunity to work under the direction of Professor Clayton H. Heathcock at the University of California, Berkeley, as a visiting graduate student for three months in 1988. From 1989 to 1991, he was a Japan Society for the Promotion of Science (JSPS) Fellow under the Japanese Junior Scientists Program. Beginning in 1991, he spent 15 months carrying out postdoctoral studies with Professor E. J. Corey at Harvard University. In 1992, Dr. Ishihara returned to Japan to join Professor Hisashi Yamamoto's group at Nagoya University as an assistant professor. In 1997, he was promoted to the rank of associate professor and, in 2002, he was appointed to his current position as a full professor at Nagoya University. Dr. Ishihara received the Inoue Research Award for Young Scientists (1994); the Chemical Society of Japan Award for Young Chemists (1996); the Thieme Chemistry Journal Award (2001); the Green and Sustainable Chemistry Award from the Ministry of Education, Culture, Sports, Science and Technology (2003); the JSPS Prize (2005); the BCSJ Award (2005); the International Conference on Cutting-Edge Organic Chemistry in Asia Lectureship Award (2006); a Japan/UK GSC Symposium Lectureship (2007); the IBM Japan Science Prize (2007); and the Mukaiyama Award (2009). His research interests include asymmetric catalysis, biomimetic catalysis induced by artificial enzymes, dehydrative condensation catalysis towards green and sustainable chemistry, and acid-base combination chemistry. 

TPGS-750-M

A second-generation amphiphile for performing organometallic chemistry in water.



For more information, please see sigma-aldrich.com/tpgs750m


Boronic Acid Surrogates

The use boronic acids and esters for Suzuki–Miyaura cross-coupling reactions is an invaluable tool for the construction of complex molecules. However, complications associated with the purification and handling of these sensitive organometallic reagents have detracted from the use of this method. To address these issues, two new classes of stable boron anion equivalents have been developed and are now available from Aldrich: MIDA boronates and trifluoroborate salts.


New MIDA Boronates from Aldrich

Benefits of MIDA Boronates


- Bench-top and chromatographically stable
- Good solubility properties and monomeric in nature
- Controlled release for in situ cross-couplings




Reference: (1) Gillis, E. P.; Burke, M. D. *Aldrichimica Acta* **2009**, *42*, 17
 (2) Knapp, D. M.; Gillis, E. P.; Burke, M. D. *J. Am. Chem. Soc.* **2009**, *131*, 6961




719390




723053




723959




704563




700908




702269




703370




699845




701084




697443




717193




708828




701092



701106



701017




704547

New Trifluoroborate Salts from Aldrich

Benefits of Trifluoroborate Salts

- Air- and moisture-stable
- Tetracoordinated monomeric species
- Less prone to protodeborylation


Reference: (1) Molander, G. A.; Canturk, B.; Kennedy, L. E. *J. Org. Chem.* **2009**, *74*, 973. (2) Molander, G. A.; Figueroa, R. *Aldrichimica Acta* **2005**, *38*, 49




711098



711101




711136




729299




711144



592846



719420



576093




576107



723592



659746




563056




684937



723916



710083



710059

To see our full line of boronic acid surrogates, please visit aldrich.com/boron



P.O. Box 14508
St. Louis, MO 63178
USA

Prsrt Std.
U.S. Postage
PAID
St. Louis MO
Permit #750

Change Service Requested

#BYNWDXZ*****AUTO**3-DIGIT 532
#03552552W075G04T# P-12 P4
Alfred R. Bader
Researcher
924 E Juneau Ave Unit 622
Milwaukee, WI 53202-6827



Sigma-Aldrich* Worldwide Offices

Argentina

Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

Chile

Tel: (+56) 2 495 7395
Fax: (+56) 2 495 7396

China

Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

Tel: (+45) 43 56 59 00
Fax: (+45) 43 56 59 05

Finland

Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 88
Fax: (+33) 474 95 68 08

Germany

Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Hungary

Ingyenes telefonszám: 06 80 355 355
Ingyenes fax szám: 06 80 344 344
Tel: (+36) 1 235 9063
Fax: (+36) 1 269 6470

India

Telephone
Bangalore: (+91) 80 6621 9400
New Delhi: (+91) 11 4358 8000
Mumbai: (+91) 22 4087 2364
Hyderabad: (+91) 40 4015 5488
Kolkata: (+91) 33 4013 8000

Fax

Bangalore: (+91) 80 6621 9550
New Delhi: (+91) 11 4358 8001
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 4015 5466
Kolkata: (+91) 33 4013 8016

Ireland

Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 402 20370
Fax: (+353) 402 20375

Israel

Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

Free Tel: 800 827 018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: (+52) 722 276 1600
Fax: (+52) 722 276 1601

The Netherlands

Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

Tel: (+47) 23 17 60 00
Fax: (+47) 23 17 60 10

Poland

Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

Tel: (+7) 495 621 5828
Fax: (+7) 495 621 6037

Singapore

Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

Slovakia

Tel: (+421) 255 571 562
Fax: (+421) 255 571 564

South Africa

Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313

United States

Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Vietnam

Tel: (+84) 3516 2810
Fax: (+84) 6258 4238

Internet

sigma-aldrich.com

Enabling Science to
Improve the Quality of Life

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052
Technical Service (800) 325-5832 • sigma-aldrich.com/techservice
Development/Custom Manufacturing Inquiries **SAFC** (800) 244-1173
Safety-related Information sigma-aldrich.com/safetycenter

World Headquarters
3050 Spruce St.
St. Louis, MO 63103
(314) 771-5765
sigma-aldrich.com



SIGMA-ALDRICH

DEDICATED TO DR. JAI P. NAGARKATTI

Aldrichimica ACTA

VOL. 43, NO. 3 • 2010

 **ALDRICH**
Chemistry




Remembering Dr. Jai P. Nagarkatti (1947–2010)

Recent Advances in the Application of α -Phenylethylamine (α -PEA)
in the Preparation of Enantiopure Compounds

2-Iodoxybenzenesulfonic Acid (IBS) Catalyzed Oxidation of Alcohols

SIGMA-ALDRICH



Aldrich Congratulates the 2010 Winners of the Nobel Prize in Chemistry

Professors **Ei-ichi Negishi**, **Akira Suzuki**, and **Richard F. Heck** have been awarded the 2010 Nobel Prize in Chemistry by The Royal Swedish Academy of Sciences for their trailblazing contributions to the area of "palladium-catalyzed cross couplings in organic synthesis".



Photo courtesy of Purdue U/
Andrew Hancock

Distinguished Professor
Ei-ichi Negishi
Purdue University, U.S.A.



Photo courtesy of Hokkaido U.

Distinguished Professor Emeritus
Akira Suzuki
Hokkaido University, Japan



Photo courtesy of U. of Delaware

Professor Emeritus
Richard F. Heck
University of Delaware, U.S.A.

Aldrich congratulates these distinguished chemistry pioneers on this achievement and thanks them for their lasting contributions to organic synthesis. Our company is proud to have had close collaborations with Professors Negishi and Suzuki, two former associates of another Nobel laureate and longtime Aldrich collaborator and member of the Board of Directors, the late Professor H. C. Brown.

The *Aldrichimica Acta* has had the distinct privilege of publishing high-impact review articles by a number of chemistry Nobel laureates, such as Professors Derek H. R. Barton, Herbert C. Brown, Elias James Corey, Ei-ichi Negishi, George A. Olah, Charles J. Pedersen, Vladimir Prelog, and K. Barry Sharpless.

The tradition of supporting future scientific leaders is strong at Sigma-Aldrich. Throughout our history, we have acknowledged scientific excellence through sponsorship of awards, symposia, and graduate-level research in the fields of chemistry, life science, and materials science. A number of previous winners of Sigma-Aldrich sponsored American Chemical Society awards (e.g., Professors George A. Olah, K. Barry Sharpless, and Ei-ichi Negishi) have gone on to receive the ultimate recognition, the Nobel Prize in Chemistry.

To access highly informative *Aldrichimica Acta* reviews written by Nobel laureates and other chemistry luminaries, visit aldrich.com/acta

Aldrichimica ACTA

VOL. 43, NO. 3 • 2010

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-231-8327
SAFC® 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 414-438-3850
24-Hour Emergency 414-438-3850
Website sigma-aldrich.com
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 2988, Milwaukee, WI 53201, USA

Subscriptions

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 2988
Milwaukee, WI 53201-2988
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the back cover.

The *Aldrichimica Acta* is also available on the Internet at sigma-aldrich.com/acta.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices listed in this publication are subject to change without notice.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2010 Sigma-Aldrich Co.

"PLEASE BOTHER US."




Joe Porwoll

Joe Porwoll, President
Aldrich Chemical Co., Inc.

A number of customers have suggested that we produce 2-pyridinylboronic acid MIDA ester. This isolable, chromatographically and benchtop-stable reagent serves as a viable 2-pyridylboron anion equivalent, and is used as a versatile cross-coupling partner

(1) Knapp, D. M.; Gillis, E. P.; Burke, M. D. *J. Am. Chem. Soc.* **2009**, *131*, 6961



719390 2-Pyridinylboronic acid MIDA ester 1 g
5 g

Naturally, we made this useful boronic acid surrogate. It was no bother at all, just a pleasure to be able to help

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the back cover

TABLE OF CONTENTS

Remembering Dr. Jai P. Nagarkatti (1947–2010)	63
<i>Sharbil J. Firsan, Sigma-Aldrich Corp</i>	
Recent Advances in the Application of α-Phenylethylamine (α-PEA) in the Preparation of Enantiopure Compounds	65
<i>Yamir Bandala and Eusebio Juanisti, Centro de Investigación y de Estudios Avanzados del Instituto Politécnico Nacional</i>	
2-Iodoxybenzenesulfonic Acid (IBS) Catalyzed Oxidation of Alcohols	83
<i>Muhammet Uyanik and Kazuaki Ishihara,* Nagoya University</i>	

ABOUT OUR COVER

Taos Pueblo Snow (oil on canvas, 40.6 × 70.0 cm) was painted in 2009 by the American painter and sculptor Rosie Sandifer (b. 1946) following a visit to the Pueblo. It depicts the centuries-old, Native American village with the same name. Taos Pueblo is located about two miles north of the town of Taos, New Mexico, and is one of the longest continually inhabited places in the U.S.A. It has been designated a World Heritage Site by UNESCO and a National Historic Landmark by the U.S. Government



Photograph © Rosie Sandifer

Equally gifted at painting and sculpting, Sandifer paints her impressions of what she observes in nature by simplifying to the essentials the effects of fleeting light on the subject, as evidenced by the white, brown, and turquoise that dominate in this painting. While her subject matter has included landscapes, figures, and animals, she has tended to focus on Western landscapes. She has participated in numerous solo, group, and juried exhibitions, and her paintings and sculptures are on display in a number of museums and public places in the U.S. Currently a resident of Santa Fe, NM, Sandifer received her extensive art education and training most notably at the Froman Painting School in Cloudcroft, NM, and the Art Students League in Stowe, VT, where she was influenced by Frank Mason, a classical realist painter and one of the most acclaimed modern American painters and teachers.

This painting is provided courtesy of the artist and is in her private collection.




New Oxetane Building Blocks

Oxetanes have historically been of modest interest to synthetic and medicinal chemists, perhaps with the natural product paclitaxel or TAXOL® being the best known example of an oxetane-containing substance. Presently, oxetanes are receiving greater attention as attractive modules for drug discovery, largely due to a series of reports from Rogers-Evans, Carreira, and coworkers. These reports have demonstrated the improved physico- and biochemical properties of a molecular scaffold when an oxetane unit replaces a *gem*-dimethyl unit¹ and the ability for an oxetane ring to function as a surrogate for a carbonyl group.^{2,1b} Another recent report has disclosed the use of 1,6-substituted azaspiro[3.3]heptanes containing an oxetane ring as alternatives to unstable 1,3-heteroatom-substituted cyclohexanes.³ In most cases, 3-oxetanone, **731536**, was the principal building block employed by the authors to install the oxetane unit (**Scheme 1**).


References: (1) (a) Wuitschik, G. et al. *Angew. Chem., Int Ed.* **2006**, *45*, 7736. (b) Wuitschik, G. et al. *J. Med. Chem.* **2010**, *53*, 3227. (2) Wuitschik, G. et al. *Angew. Chem., Int Ed.* **2008**, *47*, 4512. (3) Burkhard, J. A. et al. *Org. Lett.* **2010**, *12*, 1944

For a complete list of available oxetanes, please visit aldrich.com/oxetane



Scheme 1. Oxetane Derivatives Synthesized from 3-Oxetanone


We are pleased to offer a wide selection of new oxetane building blocks for a variety of applications in synthetic and medicinal chemistry.



New Trifluoroborate Salts from Aldrich

Benefits of Trifluoroborate Salts

- Air- and moisture-stable
- Tetracoordinated monomeric species
- Less prone to protodeborylation




Reference: (1) Molander, G. A.; Canturk, B.; Kennedy, L. E. *J. Org. Chem.* **2009**, *74*, 973. (2) Molander, G. A.; Figueroa, R. *Aldrichimica Acta* **2005**, *38*, 49



711098



711101




711136




729299




711144



592846



719420



576093




576107



723592



659746




563056




684937



723916



710083



710059

To see our full line of boronic acid surrogates, please visit aldrich.com/boron

THE GROWING ROLE OF ORGANOCATALYSIS IN ASYMMETRIC SYNTHESIS

Aldrichimica **ACTA**

VOL. 44, NO. 1 • 2011

 **ALDRICH**
Chemistry



Asymmetric N-Heterocyclic Carbene (NHC) Catalyzed Acyl Anion Reactions

Recent Advances in the Asymmetric Catalytic Mannich Reaction

SIGMA-ALDRICH

When you need Privileged Structures

Add  Aldrich

Indoles and Indole Isosteres


Substituted indoles are often referred to as "privileged structures". They are capable of binding to multiple receptors with high affinity and have applications across a wide range of therapeutic areas.¹ Similarly, isosteres of the indole ring:

- Allow researchers to attenuate or amplify the activity of a target structure without altering steric bulk
- Produce significant biological activity²
- Comprise essential subunits in many pharmaceutically relevant compounds
- Mimic purines in its role as a hydrogen-bonding partner


When you need Indoles and Indole Isosteres, Add Aldrich to your research program.

Aldrich.com/bb


References: (1) Horton, D. A. et al. *Chem. Rev.* **2003**, *103*, 893 and references therein. (2) (a) Popowycz, F. et al. *Tetrahedron* **2007**, *63*, 1031. (b) Popowycz, F. et al. *Tetrahedron* **2007**, *63*, 8689. (c) Song, J. J. et al. *Chem. Soc. Rev.* **2007**, *36*, 1120. (d) Huang, W.-S. et al. *J. Med. Chem.* **2010**, *53*, 4701. (e) Buckley, G. M. et al. *Bioorg. Med. Chem. Lett.* **2008**, *18*, 3656. (f) Buckley, G. M. et al. *Bioorg. Med. Chem. Lett.* **2008**, *18*, 3291.




733040




723789




716529




724378




732168




707953



723770



721050



732141

Aldrichimica ACTA

VOL. 44, NO. 1 • 2011

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-231-8327
SAFC® 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 414-438-3850
24-Hour Emergency 414-438-3850
Website sigma-aldrich.com
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 2988, Milwaukee, WI 53201, USA
sharbil.firsan@sial.com

Subscriptions

To request your FREE subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 2988
Milwaukee, WI 53201-2988
Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the back cover.

The *Aldrichimica Acta* is also available at Aldrich.com/acta

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc., warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2011 Sigma-Aldrich Co.

"PLEASE BOTHER US."




John Radke

John Radke
Director of Marketing, Chemistry

Professor Alison Thompson of the Department of Chemistry at Dalhousie University recently suggested that we introduce *N*-Boc-pyrrole-2-boronic acid MIDA ester. This air- and moisture-stable reagent serves as a 2-heterocyclic boronic acid surrogate that is useful in Suzuki-Miyaura cross-coupling reactions.

Knapp, D. M.; Gillis, E. P.; Burke, M. D. *J. Am. Chem. Soc.* **2009**, *131*, 6961



733539 *N*-Boc-pyrrole-2-boronic acid MIDA ester, 95%

250 mg
1 g

Naturally, we made this useful heterocyclic reagent. It was no bother at all, just a pleasure to be able to help.

Do you have a compound that you wish Aldrich could list, and that would help you in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the back cover.

TABLE OF CONTENTS

Asymmetric N-Heterocyclic Carbene (NHC) Catalyzed Acyl Anion Reactions.....	3
<i>Harit U. Vora and Tomislav Rovis,* Colorado State University</i>	
Recent Advances in the Asymmetric Catalytic Mannich Reaction	15
<i>Sandro José Greco,* Valdemar Lacerda, Jr., and Reginaldo Bezerra dos Santos, Universidade Federal do Espírito Santo (Brasil)</i>	

ABOUT OUR COVER

River Landscape with a Resting Traveller

(oil on panel, 46 × 66.5 cm) was painted most likely in the mid-1650s by Jan Lievens (1607–1674), who was one of the most prominent painters of the Dutch Golden Age. He is usually mentioned in the same breath as Rembrandt, because the two artists shared much during their early years. Because of Rembrandt's greater achievement, Lievens was long assumed to have been his follower, but recently it has become clear that he was initially the leader. Although younger, he was a celebrated child prodigy who established himself years before his colleague.

Lievens's later success lay in following and interpreting the grand figurative style of Flemish masters such as Peter Paul Rubens and Anthony van Dyck, much favored among the elite. He could study such works while in Antwerp, where he also saw the contemplative landscapes by Rubens and other artists. These reflected the Flemish tradition of rich, lively, and decorative fantasy landscapes, in contrast to the more sober realism exercised by Dutch landscapists. In his own landscapes, Lievens tried a more fluid brush handling and rolling forms. In this late scene, he places a weary traveler reclining against a tree in a lush forest. He also betrays his penchant for arranging a screen of trees against a light background that pierces between the trunks to generate drama and rhythm. The result is slightly otherworldly, an invitation to escape from the urban confines and daily pressures, that would have appealed to many collectors in both Flanders and Holland.

This painting is part of the Bader Collection of Dutch and Flemish Paintings, whose future home will be the Agnes Etherington Art Centre of Queen's University, Kingston, ON, Canada.



Detail from *River Landscape with a Resting Traveller*
By *Flavio de M. M. Lacerda, Jr.*

When you need a stable liquid Reductive Amination Reagent

Add Aldrich

PEMB is a liquid Pyridine Borane Complex useful for reductive amination chemistry

Advantageous Properties of PEMB:


- Excellent for reductive aminations
- Mild reducing agent for imines and oximes
- Reaction with protic solvents is very slow
- Soluble in aromatic hydrocarbons, alcohols, and ether solvents
- Can be used solvent-free for reductive aminations
- Chemically efficient: two of three hydrides are utilized

Add Aldrich to your research program.

Aldrich.com/pemb



For research quantities under 500 grams, please contact Aldrich Chemistry.

For quantities over 500 grams, please contact BASF.



5-Ethyl-2-methylpyridine borane (PEMB)
725080

Select Substrate Scope



Examples

Aldehyde and Amine	Conditions	Product	% Yield in MeOH (% Yield Neat)
	PEMB, AcOH MeOH, 25 °C		72 (80)
	PEMB, AcOH MeOH, 25 °C		0 (96)
	PEMB MeOH, 25 °C		92 (94)
	PEMB, AcOH MeOH, 25 °C		92 (93)
	PEMB, AcOH MeOH, 50 °C		74 (94)

For more examples and experimental details, please see: Burkhardt, E. R.; Coleridge, B. M. *Tetrahedron Lett.* **2008**, 49, 5152.

Distillation Adapter for On-The-Fly Sampling

SIGMA-ALDRICH
Labware

Distillation is the most widely used bulk separation method in the laboratory as well as in industry. Beyond purification, it is widely utilized to characterize complex fluids (such as fuels) through measurement of the distillation curve, a plot of the boiling temperature against volume distilled. A common theme in both of these applications is the desire to understand the composition. In purification, the goal is to monitor the distillation progress; and in fluid characterization, one seeks to relate the composition to the temperature data.

Features and Benefits

A distillate sampling adapter (**Figure 1**), installed after a condenser or distillation column, can provide this important capability without the need for cumbersome, expensive, and often unreliable fraction collectors.¹⁻³ The flow of the distillate is focused to drop into a 0.05 mL "hammock" that is positioned directly below the flow path. The sampling port, equipped with a vacuum-tight valve, allows access to the hammock with a standard chromatographic syringe through a septum. To sample the distillate, one simply positions the chromatographic syringe, preferably equipped with a blunt-tipped needle, in the well of the hammock. It is a simple matter to withdraw samples as the distillation progresses. The sample can then be directly injected into the gas chromatograph or spectrometer, or injected into an autosampler vial for analysis later. Indeed, any analytical technique that is applicable for liquid samples ranging in volume from 1 to 50 microliters can be used to characterize the distillate.

The distillate sampling adapter is offered with several joint options, and is designed to work with existing modular distillation glassware. It is simply a matter of replacing the existing vacuum adapter with a distillate sampling adapter. No further modification to the setup is required.

For more information, visit Aldrich.com/labware




Figure 1. Distillate sampling adapter

Aldrich® Distillate Sampling Adapter, with Vacuum Connection and PTFE Valve

Joint Size	Cat. No.
14/20	Z569895
24/40	Z569909
29/32	Z569917
Replacement valve septa	33310-U
Septum inserter for valve	33311
Hamilton 7015NR syringe, 10 µL, 22s gauge blunt tip needle	58380-U

References:

- (1) Bruno, T. J.; Ott, L. S.; Smitn, B. L.; Lovestead, T. M. Complex fluid analysis with the advanced distillation curve approach. *Anal. Chem.* **2010**, *82*, 777
- (2) Bruno, T. J.; Ott, L. S.; Lovestead, T. M.; Huber, M. L. The composition-explicit distillation curve technique: relating chemical analysis and physical properties of complex fluids. *J. Chromatogr., A* **2010**, *1217*, 2703
- (3) Bruno, T. J.; Ott, L. S.; Lovestead, T. M.; Huber, M. L. Relating complex fluid composition and thermophysical properties with the advanced distillation curve approach. *Chem. Eng. Tech.* **2010**, *33*, 363



P.O. Box 14508
St. Louis, MO 63178
USA

Sigma-Aldrich® Worldwide Offices

Argentina

Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

Chile

Tel: (+56) 2 495 7395
Fax: (+56) 2 495 7396

People's Republic of China

Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

Tel: (+45) 43 56 59 00
Fax: (+45) 43 56 59 05

Finland

Tel: (+358) 9 350 9250
Fax: (+358) 9 350 9255

France

Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 88
Fax: (+33) 474 95 68 08

Germany

Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1169

Hungary

Ingyenes telefonszám: 06 80 355 355
Ingyenes fax szám: 06 80 344 344
Tel: (+36) 1 235 9055
Fax: (+36) 1 269 6470

India

Telephone

Bangalore: (+91) 80 6621 9400
New Delhi: (+91) 11 4358 8000
Mumbai: (+91) 22 4087 2364
Hyderabad: (+91) 40 4015 5488
Kolkata: (+91) 33 4013 8000

Fax

Bangalore: (+91) 80 6621 9550
New Delhi: (+91) 11 4358 8001
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 4015 5466
Kolkata: (+91) 33 4013 8016

Ireland

Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 402 20370
Fax: (+353) 402 20375

Israel

Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4222
Fax: (+972) 8 948 4200

Italy

Free Tel: 800 827 018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Luxembourg

Tel: (+32) 3 899 1301
Fax: (+32) 3 899 1311

Malaysia

Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: (+52) 722 276 1600
Fax: (+52) 722 276 1601

The Netherlands

Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

Tel: (+47) 23 17 60 00
Fax: (+47) 23 17 60 10

Poland

Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

Tel: (+7) 495 621 5828
Fax: (+7) 495 621 6037

Singapore

Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

Slovakia

Tel: (+421) 255 571 562
Fax: (+421) 255 571 564

South Africa

Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2511
Fax: (+41) 81 756 5449

Taiwan

SAFC Hitech
Tel: (+886) 7 695 5066
Fax: (+886) 7 695 5088

Thailand

Tel: (+66) 2 126 8141
Fax: (+66) 2 126 8080

United Kingdom

Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313

United States

Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Vietnam

Tel: (+84) 3516 2810
Fax: (+84) 6258 4238

Internet

sigma-aldrich.com

Enabling Science to
Improve the Quality of Life

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052
Technical Service (800) 325-5832 • sigma-aldrich.com/techservice
Development/Custom Manufacturing Inquiries **SAFC**® (800) 244-1173
Safety-related Information sigma-aldrich.com/safetycenter

World Headquarters
3050 Spruce St.
St. Louis, MO 63103
(314) 771-5765
sigma-aldrich.com



SIGMA-ALDRICH

ORGANIC SYNTHESIS FACILITATED BY HALOGENS AND TRANSITION METALS

Aldrichimica ACTA

VOL. 44, NO. 2 • 2011

 **ALDRICH**
Chemistry



Halonium-Induced Cyclization Reactions

The A³-Coupling (Aldehyde–Alkyne–Amine) Reaction:
A Versatile Method for the Preparation of Propargylamines

SIGMA-ALDRICH

Order your 2012–2014 Aldrich® Handbook.

Add  Aldrich

The new Aldrich Handbook contains the widest selection of chemistry and materials science products and is your resource for chemical structures, literature references, and extensive chemical and physical data. Our complimentary catalog includes new and innovative reagents and building blocks, plus a focused line of Labware products to support your chemistry needs.

The Aldrich Handbook's portfolio supports the research community with:

- More than 40,000 research chemicals
- Over 4,000 new products
- 10,000 chemical structures
- 8,500 updated literature citations
- Extensive chemical and physical data

For reliable, high-quality chemicals you can trust, add your free copy of the Aldrich Handbook to your laboratory by visiting:

Aldrich.com/catalogs



Aldrichimica ACTA

VOL. 44, NO. 2 • 2011

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)

FAX 800-325-5052 (USA)

or 414-438-2199

Mail P.O. Box 2060

Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010

Technical Service 800-231-8327

SAFC® 800-244-1173

Custom Synthesis 800-244-1173

Flavors & Fragrances 800-227-4563

International 414-438-3850

24-Hour Emergency 414-438-3850

Website sigma-aldrich.com

Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.

P.O. Box 2988, Milwaukee, WI 53201, USA

sharbil.firsan@sial.com

Subscriptions

To request your **FREE** subscription to the *Aldrichimica Acta*, please contact us by:

Phone: 800-325-3010 (USA)

Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 2988
Milwaukee, WI 53201-2988

Email: sams-usa@sial.com


International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see the back cover.

The *Aldrichimica Acta* is also available at Aldrich.com/acta

Aldrich brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product for its particular use. See product information on the Sigma-Aldrich website at www.sigma-aldrich.com and/or on the reverse side of invoice or packing slip for additional terms and conditions of sale.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group. © 2011 Sigma-Aldrich Co. LLC.


"PLEASE BOTHER US."



John Radke
Director of Marketing, Chemistry

Professor Steven Leighton at the University of Wisconsin-Madison recently suggested that we introduce Leighton's chiral crotylsilane reagents. These air-stable, crystalline solids may be stored at room temperature and are useful for the enantioselective crotylation of aldehydes and other carbonyl compounds

(a) Kim, H.; Ho, S.; Leighton, J. L. *J. Am. Chem. Soc.* **2011**, *133*, 6517. (b) Leighton, J. L. *Aldrichimica Acta* **2010**, *43*, 3



733199	(R,R)-1,3-Bis(4-bromobenzyl)-2-chlorooctahydro-2-(2Z)-crotyl-1H-1,3,2-benzodiazasilole, 97%	1 g
733180	(R,R)-1,3-Bis(4-bromobenzyl)-2-chlorooctahydro-2-(2E)-crotyl-1H-1,3,2-benzodiazasilole	250 mg 1 g
733075	(S,S)-1,3-Bis(4-bromobenzyl)-2-chlorooctahydro-2-(2Z)-crotyl-1H-1,3,2-benzodiazasilole, 97%	1 g 5 g
733083	(S,S)-1,3-Bis(4-bromobenzyl)-2-chlorooctahydro-2-(2E)-crotyl-1H-1,3,2-benzodiazasilole	250 mg 1 g

Naturally, we made these useful crotylation reagents available. It was no bother at all, just a pleasure to be able to help

Do you have a compound that you are looking for and that we don't carry? You in your research by saving you time and money? If so, please send us your suggestion; we will be delighted to give it careful consideration. You can contact us in any one of the ways shown on this page and on the back cover

TABLE OF CONTENTS

Halonium-Induced Cyclization Reactions..... 27

Scott A. Snyder, Daniel S. Treitler, and Alexandra F. Drake, Columbia University

The A³-Coupling (Aldehyde–Alkyne–Amine) Reaction: A Versatile Method for the Preparation of Propargylamines..... 43

Woo-Jin Yoo, Liana Zhao, and Chao-Jun Li, McGill University

ABOUT OUR COVER

Aelbert Jacobsz Cuyp (1620–1691) of the Golden Age of Dutch Art painted *The Maas at Dordrecht* (oil on canvas, 114.9 × 170.2 cm) around 1650. It is possible that Cuyp was commissioned to commemorate an event that may have occurred during the summer of 1646. At that time, an enormous fleet of ships carrying around thirty thousand soldiers was anchored at Dordrecht presumably for symbolic purposes rather than for specific military ones as peace was finally at hand. The Peace of Münster, which ended the Eighty Years' War with Spain, was signed only two years later, in 1648.



Detail from *The Maas at Dordrecht* Photograph

One of Cuyp's finest paintings, this work depicts a flurry of maritime activity, while the intricately painted pleyt (riverboat) and adjacent rowboat in the right foreground are clearly the painting's focal point. Not only are the boats authentically depicted, but they also contain numerous figures that have personality and purpose. Most of the ships on the river have their sails raised and flags flying as though they are about to embark. The early morning light floods the scene and creates striking patterns on the clouds, sails, and water, adding a dramatic character to the setting. The variety and massiveness of Cuyp's forms give this work a tangibility that few other marine painters could achieve.

This painting is part of the Andrew W. Mellon Collection at the National Gallery of Art, Washington, DC.

When you need a stable liquid Reductive Amination Reagent.

Add Aldrich

PEMB is a liquid Pyridine Borane Complex useful for reductive amination chemistry.

Advantageous Properties of PEMB:


- Excellent for reductive aminations
- Mild reducing agent for imines and oximes
- Reaction with protic solvents is very slow
- Soluble in aromatic hydrocarbons, alcohols and ether solvents
- Can be used solvent-free for reductive aminations
- Chemically efficient: two-of-three hydrides are utilized

Add Aldrich to your research program.

Aldrich.com/pemb


Contact Aldrich Chemistry for research quantities under 500 grams.

For quantities over 500 grams, contact BASF.



5-Ethyl-2-methylpyridine borane (PEMB)
725080

Select Substrate Scope



Examples

Aldehyde and Amine	Conditions	Product	% Yield in MeOH (% Yield Neat)
	PEMB, AcOH MeOH, 25 °C		72 (80)
	PEMB, AcOH MeOH, 25 °C		0 (96)
	PEMB MeOH, 25 °C		92 (94)
	PEMB, AcOH MeOH, 25 °C		92 (93)
	PEMB, AcOH MeOH, 50 °C		74 (94)

For more examples and experimental detail
Burkhardt, E. R.; Coleridge, B. M. *Tetrahedron Lett.* **2008**, 49, 5152.



Missing out on the latest research developments in Chemistry?

Add **Aldrich**

Aldrich ChemFiles, a quarterly Chemical Synthesis product guide written by our experts, enables you to advance your chemistry research more effectively by implementing the latest innovative synthetic strategies. We help keep you informed of new Aldrich Chemistry products which facilitate the latest research methodologies and trends, allowing you to access key starting materials and reagents more efficiently.

Each edition highlights useful applications of these new products, and shows where they can be quickly found on the Aldrich Chemistry website.

Request your complimentary copy today.

Aldrich.com/chemfiles



P.O. Box 14508
St. Louis, MO 63178
USA

Sigma-Aldrich® Worldwide Offices

Argentina

Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

Chile

Tel: (+56) 2 495 7395
Fax: (+56) 2 495 7396

People's Republic of China

Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

Tel: (+45) 43 56 59 00
Fax: (+45) 43 56 59 05

Finland

Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 88
Fax: (+33) 474 95 68 08

Germany

Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1169

Hungary

Tel: (+36) 1 235 9055
Fax: (+36) 1 235 9068

India

Telephone

Bangalore: (+91) 80 6621 9400
New Delhi: (+91) 11 4358 8000
Mumbai: (+91) 22 4087 2364
Pune: (+91) 20 4146 4700
Hyderabad: (+91) 40 3067 7450
Kolkata: (+91) 33 4013 8000

Fax

Bangalore: (+91) 80 6621 9550
New Delhi: (+91) 11 4358 8001
Mumbai: (+91) 22 2579 7589
Pune: (+91) 20 4146 4777
Hyderabad: (+91) 40 3067 7451
Kolkata: (+91) 33 4013 8016

Ireland

Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: +353 (0) 402 20370
Fax: +353 (0) 402 20375

Israel

Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4222
Fax: (+972) 8 948 4200

Italy

Free Tel: 800 827 018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Luxembourg

Tel: (+32) 3 899 1301
Fax: (+32) 3 899 1311

Malaysia

Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: (+52) 722 276 1600
Fax: (+52) 722 276 1601

The Netherlands

Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

Tel: (+47) 23 17 60 00
Fax: (+47) 23 17 60 10

Poland

Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

Tel: (+7) 495 621 5828
Fax: (+7) 495 621 6037

Singapore

Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

Slovakia

Tel: (+421) 255 571 562
Fax: (+421) 255 571 564

South Africa

Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2511
Fax: (+41) 81 756 5449

Thailand

Tel: (+66) 2 126 8141
Fax: (+66) 2 126 8080

United Kingdom

Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313

United States

Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Vietnam

Tel: (+84) 8 3516 2810
Fax: (+84) 8 6258 4238

Internet

sigma-aldrich.com

Enabling Science to
Improve the Quality of Life

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052
Technical Service (800) 325-5832 • sigma-aldrich.com/techservice
Development/Custom Manufacturing Inquiries **SAFC**® (800) 244-1173
Safety-related Information sigma-aldrich.com/safetycenter

World Headquarters
3050 Spruce St.
St. Louis, MO 63103
(314) 771-5765
sigma-aldrich.com



© 2005 Sigma-Aldrich. All rights reserved. Sigma-Aldrich, the Sigma-Aldrich logo, and the SAFC logo are trademarks of Sigma-Aldrich. FSC® is a registered trademark of the Forest Stewardship Council. All other trademarks are the property of their respective owners.

NTH
76536-512540
1081

SIGMA-ALDRICH

VOL. 44, NO. 3 • 2011

ALDRICH
Chemistry

Aldrichimica **ACTA**



Peroxide-Mediated Wacker Oxidations for Organic Synthesis

Organofluorine Chemistry: Deoxyfluorination Reagents
for C–F Bond Synthesis

SIGMA-ALDRICH

Don't forget to order
your 2012–2014
Aldrich Handbook.

Add  Aldrich

The new Aldrich Handbook contains the widest selection of Chemistry and Materials Science products and is your resource for chemical structures, literature references, and extensive chemical and physical data. Our complimentary catalog includes new and innovative reagents and building blocks, plus a focused line of Labware products to support your chemistry needs.

The Aldrich Handbook's portfolio supports the research community with:

- More than 40,000 research chemicals
- Over 4,000 new products
- 10,000 chemical structures
- 8,500 updated literature citations
- Extensive chemical and physical data

For reliable, high-quality chemicals you can trust, add your free copy of the Aldrich Handbook to your laboratory by visiting:

Aldrich.com/catalogs



Aldrich Chemical Co., Inc.
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-231-8327
SAFC® 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 414-438-3850
24-Hour Emergency 414-438-3850
Website sigma-aldrich.com
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 2988, Milwaukee, WI 53201, USA
sharbil.firsan@sial.com

Subscriptions

Request your FREE subscription to the
Aldrichimica Acta, through:


Web: Aldrich.com/acta
Email: sales-usa@sial.com
Phone: 800-325-3010 (USA)
Mail: Attn: Mailroom
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 2988
Milwaukee, WI 53201-2988

International customers, please contact your local Sigma
Aldrich office. For worldwide contact information, please see
the last numbered page of this issue.

The entire *Aldrichimica Acta* archive is also available free of
charge at Aldrich.com/acta.

Aldrich brand products are sold through Sigma-Aldrich, Inc.
Purchaser must determine the suitability of the product for
their particular use. See product information on the Sigma
Aldrich website at www.sigma-aldrich.com and/or on the
reverse side of invoice or packing slip for additional terms and
conditions of sale.

Aldrichimica Acta (ISSN 0002-5700) is a publication of
Aldrich. Aldrich is a member of the Sigma-Aldrich Group.
© 2011 Sigma-Aldrich LLC.




"PLEASE BOTHER US."


Dear Fellow Chemists,

Dr. Matthew Tudge at Merck & Co. recently suggested that we introduce air- and moisture-stable medicinal chemistry oriented palladium(II) precatalysts. These stable, crystalline precatalysts rapidly form highly active L1 palladium(0) catalysts for C-C and C-N couplings with loadings as low as 0.15 mol %


Kinzel, T., Zhang, Y., Buchwald, S. L. *J Am Chem Soc* **2010**, *132*, 14073



753246



753009



741825

753246	Chloro(2-dicyclohexylphosphino-2',6'-diisopropoxy-1,1'-biphenyl)-[2-(2'-amino-1,1'-biphenyl)]palladium(II)	250 mg 1 g 5 g
753009	Chloro(2-dicyclohexylphosphino-2',6'-dimethoxy-1,1'-biphenyl)-[2-(2'-amino-1,1'-biphenyl)]palladium(II)	250 mg 1 g 5 g
741825	Chloro(2-dicyclohexylphosphino-2',4',6'-triisopropyl-1,1'-biphenyl)-[2-(2'-amino-1,1'-biphenyl)]palladium(II)	250 mg 1 g 5 g

We welcome fresh product ideas from you. Do you need a compound that isn't listed in our catalog? Ask Aldrich! For over 60 years, your research needs and suggestions have shaped the Aldrich product offering. To submit your suggestion visit Aldrich.com/pleasebotherus

John Radke
Director of Marketing, Chemistry


TABLE OF CONTENTS

Peroxide-Mediated Wacker Oxidations for Organic Synthesis	55
<i>Brian W. Michel and Matthew S. Sigman,* University of Utah</i>	
Organofluorine Chemistry: Deoxyfluorination Reagents for C-F Bond Synthesis.....	65
<i>Nawaf Al-Maharik and David O'Hagan,* University of St Andrews</i>	

ABOUT OUR COVER

Seashore with Fishermen (oil on canvas, 101.9 × 127.6 cm) was painted around 1781/1782 by the British portrait and landscape painter, Thomas Gainsborough (1727–1788). While better known as a portraitist for the more lucrative portraits of royalty and nobility that he painted most of his life to support his family, he never lost his fondness for landscape painting.

Gainsborough painted this landscape in his later years, when he reportedly declared that he was tired of painting portraits. Depicting fishermen struggling against strong winds and waves to launch their boat into the water, he imparts a measure of spontaneity and sensibility to the scene. His personal style is reflected in the way he merges the figures with the scene behind them, and in his handling of paint for which he was much admired.



Photograph © Board of Trustees, National Gallery of Art, Washington


This painting is part of the Ailsa Mellon Bruce Collection at the National Gallery of Art, Washington, DC.

When you need selectivity.


Add Aldrich

Aldrich Chemistry is proud to offer Prof. Matt Sigman's Pd(quinox)Cl₂ (**713929**) for the selective oxidation of terminal alkenes to methyl ketones (**Scheme 1**).¹⁻³ The extreme selectivity of this catalyst is particularly apparent when allylic functionalized substrates are employed in these Wacker-type oxidations.


Quinox Products from Aldrich:




713910





713929




778737




731005




81%, 4 h




82%, 20 min




69%, 16 h




95%, 2.5 h



88%, 4 h



74%, 2.5 h



89%, 17 h

Scheme 1: Pd-quinox catalyzed Wacker oxidation of allylic functionalized terminal alkenes.

References: (1) Michel, B. W.; Camelio, A. M.; Cornell, C. N.; Sigman, M. S. *J. Am. Chem. Soc.* **2009**, 131, 6076. (2) Michel, B. W.; McCombs, J. R.; Winkler, A.; Sigman, M. S. *Angew. Chem., Int. Ed.* **2010**, 49, 7312. (3) Michel, B. W.; Steffens, L. D.; Sigman, M. S. *J. Am. Chem. Soc.* **2011**, 133, 8317.

For the latest chemistry innovations, Add Aldrich.

To view or order these quinox ligands and catalyst, visit
Aldrich.com/quinox



All-PTFE Filter Reactor Assemblies

For Use in Fluorine Chemistry Applications and with Strong Alkaline Compounds

SIGMA-ALDRICH
Labware

PTFE filter reactors from Sigma-Aldrich® can be used for fluorine chemistry and with strong alkaline compounds that can etch glass. The head, head joints, body, stirrer-bearing, bottom filter assembly with valve, and stirring shaft with agitator are all PTFE. Even the O-rings are PTFE encapsulated. The filter supports are polypropylene. The bottom filter assembly is either #50 or #80 Ace-Thred and threads in/out for easy capture of filtrate, or for cleaning. A Halar® coated stainless steel internal coil is included (in the 1 liter size) for heating or cooling the contents. IKA® RW 20 overhead stirrer recommended.



Product Description	Cat. No.
250 mL Reactor Assembly Complete	Z690368
Components	
Reactor body	Z690376
PTFE head (60 mm diam., 4 openings, 24/40)	Z690384
CAPFE O-ring (60 mm diam.)	Z690392
Quick release clamp (60 mm diam.)	Z690406
Bottom filter assembly (#50 screw-thread joint)	Z690414
Bottom valve (1/4 in. x 1/4 in.)	Z690422
PTFE adapter	Z690430
Polypropylene filter support (#50 screw-thread joint)	Z690449
Stirrer bearing (10 mm)	Z690457
PTFE coated SS shaft with agitator (10 mm diam.)	Z690465
500 mL Reactor Assembly Complete	Z690473
Components	
Reactor body	Z690481
<i>Other components same as above</i>	
1 L Reactor Assembly Complete	Z690503
Components	
Reactor body	Z690511
PTFE head (100 mm diam., 7 openings, 4 x 24/40, 3 x #7 screw-thread joint)	Z690538
CAPFE O-ring (100 mm diam.)	Z690546
Quick release clamp (100 mm diam.)	Z690554
Bottom filter assembly (#80 screw-thread joint)	Z690562
Bottom valve (1/4 in. x 1/4 in.)	Z690422
PTFE adapter	Z690430
Polypropylene filter support (#80 screw-thread joint)	Z690570
Stirrer bearing (10 mm)	Z690457
PTFE coated SS shaft with agitator (10 mm diam.)	Z690589
Halar coated heating coil	Z690597
IKA RW 20 Digital Overhead Stirrer	
110 V	Z645087
230 V	Z645095

Halar is a registered trademark of Solvay Solexis, Inc. IKA is a registered trademark of IKA Werke GmbH & Co. KG

SIGMA-ALDRICH



P.O. Box 14508
St. Louis, MO 63178
USA

Prsrt Std.
U.S. Postage
PAID
St. Louis MO
Permit #750

Change Service Requested

#BYNWDNZ*****AUTO**3-DIGIT 532
#03552552ACTA 44.3# P-13 P4
Alfred R. Bader
Researcher
924 E. Juneau Ave Unit 622
Milwaukee, WI 53202-6827



Professor Paul Knochel

Winner of the 2011 EROS Best Reagent Award



Aldrich Chemistry and John Wiley, proud sponsors of the annual EROS Best Reagent Award, congratulate this year's winner—Professor Paul Knochel. The award was created to honor outstanding contributors to the online edition of *Encyclopedia of Reagents for Organic Synthesis* (EROS and e-EROS).


About Prof. Paul Knochel



Currently at the Department of Chemistry of Ludwig Maximilian University in Munich, Germany, Professor Paul Knochel's research interests include the development of novel organometallic reagents and methods for use in organic synthesis, asymmetric catalysis, and natural product synthesis.

About the Reagent

Discovered in 2004, the reagent is known as isopropylmagnesium chloride–lithium chloride complex or TurboGrignard™ (Product No. 656984). Since then it has found a wide range of elegant applications in laboratory synthesis, and has been produced on a large scale for industrial applications. The article on the award-winning reagent by Paul Knochel and Andrei Gavryushin was published in EROS in October 2010.



Isopropylmagnesium chloride–lithium chloride complex
Product No. 656984

Add Aldrich

Advance your research. Add Aldrich.
Aldrich.com/turbogrignard

VOL. 45, NO. 1 • 2012

ALDRICH
Chemistry

Aldrichimica **ACTA**



"Designer"-Surfactant-Enabled Cross-Couplings in
Water at Room Temperature

The Growing Impact of Niobium in Organic
Synthesis and Catalysis



SIGMA-ALDRICH


For convenient trifluoromethylation of aryl iodides.

Add Aldrich

Hartwig Trifluoromethylations

Aldrich Chemistry would like to introduce Trifluoromethylator™, a highly efficient trifluoromethylation reagent from the laboratory of Professor John Hartwig and Catylix, Inc. Trifluoromethylator™ (Phen)Cu-CF₃ (Prod. No. L510009) is a convenient, thermally stable, single-component reagent for the trifluoromethylation of aryl iodides. It provides:

- Tolerance of electrophilic functional groups, including aldehydes, nitroarenes, ketones and esters
- High yields under mild conditions for electron-rich and electron-deficient iodoarenes, as well as sterically hindered iodoarenes
- Increased flexibility for carrying out ester hydrolysis in your synthesis




Solutions of Trifluoromethylator™ in toluene are safe to use at process scale.




To request custom and bulk quantities of Trifluoromethylator™, visit Aldrich.com/catylix

Trifluoromethylator is a trademark of Catylix, Inc.

General Reaction Equation



Representative Reaction Scope



Reference: Morimoto, H.; Tsuboqo, T.; Litvinas, N. D.; Hartwig, J. F. *Angew. Chem. Int. Ed.* 2011, 50, 3793.

Aldrichimica ACTA

VOL. 45, NO. 1 • 2012

Aldrich Chemical Co., Inc.
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone: 800-325-3010 (USA)
FAX: 800-325-5052 (USA)
or 414-438-2199
Mail: P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries	800-325-3010
Technical Service	800-231-8327
SAFC®	800-244-1175
Custom Synthesis	800-244-1171
Flavors & Fragrances	800-227-4563
International	414-438-3850
24-Hour Emergency	414-438-3850
Website	sigma-aldrich.com
Email	aldrich@stl.com


General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 2988, Milwaukee, WI 53201, USA
sharbil.firsan@stl.com

Subscriptions

Request your **FREE** subscription to the Aldrichimica Acta, through:

Web:	Aldrich.com/acta
Email:	sams-usa@stl.com
Phone:	800-325-3010 (USA)
Mail:	Attn: Mailroom Aldrich Chemical Co., Inc. Sigma-Aldrich Corporation P.O. Box 2988 Milwaukee, WI 53201-2988




"PLEASE BOTHER US."


Dear Fellow Chemists,

We have had many suggestions for the DalPhos ligands shown below. These were recently developed by Professor Mark Stradiotto and his group at Dalhousie University. Both are air-stable, chelating N,P ligands that enable the Pd-catalyzed cross-coupling of ammonia. Mor-DalPhos is also effective in coupling hydrazine and acetone, and in the gold-catalyzed hydroamination. The catalysts are used with low loadings and exhibit excellent functional group tolerance and chemoselectivity.

(a) Lundgren, R. J.; Peters, B. D.; Alsabeh, P. G.; Stradiotto, M. *Angew. Chem., Int. Ed.* **2010**, *49*, 4071
(b) Lundgren, R. J.; Stradiotto, M. *Angew. Chem., Int. Ed.* **2010**, *49*, 8686



751618



751413

751618	Di(1-adamantyl)-2-morpholinophenylphosphine (Mor-DalPhos)	250 mg 1 g
751413	Di(1-adamantyl)-2-dimethylaminophenylphosphine (Me-DalPhos)	250 mg 1 g

We welcome fresh product ideas from you. Do you need a compound that isn't listed in our catalog? Ask Aldrich! For over 60 years, your research needs and suggestions have shaped the Aldrich product offering. To submit your suggestion visit Aldrich.com/pleasebotherus

John Radke
Director of Marketing, Chemistry

TABLE OF CONTENTS


"Designer"-Surfactant-Enabled Cross-Couplings in Water at Room Temperature	3
<i>Bruce H. Lipshutz* and Subir Ghorai</i> , University of California, Santa Barbara, and Sigma-Aldrich Co., Sheboygan Falls, Wisconsin	
The Growing Impact of Niobium in Organic Synthesis and Catalysis	19
<i>Valdemar Lacerda, Jr.,* Deborah Araujo dos Santos, Luiz Carlos da Silva-Filho, Sandro José Greco, and Reginaldo Bezerra dos Santos</i> , Universidade Federal do Espírito Santo and Universidade Estadual Paulista (Brasil)	

ABOUT OUR COVER

Venice: The Dogana and San Giorgio Maggiore (oil on canvas, 91.5 × 122.0 cm) was completed in 1834 by Joseph Mallord William Turner (London, 1775–1851). A prolific British artist and an influential member of the Romantic Movement of the late eighteenth and first half of the nineteenth century, he elevated landscape painting to a level not achieved before and introduced several innovations to the genre. His talent manifested itself very early in life, and he was much admired as an artist in life and death, with institutions, works, and artistic prizes dedicated to preserving his legacy.

This painting depicts a bustling maritime scene in the Grand Canal of Venice, which Turner had visited several times. (Another Venetian cityscape, painted seven years earlier by another British Romantic Movement artist, R. P. Bonington, has graced the cover of *Aldrichimica Acta*, Vol. 43, No. 1.) Turner intended this work as a celebration of commerce, as symbolized by the statue of Fortune atop the Dogana (Customs building) in the foreground. While Turner was not concerned with a faithful depiction of the scene, as evidenced by the apparent widening of the canal and placement of San Giorgio's church, his skill as a draftsman and his mastery of perspective drawing are evident in his precise, linear rendering of the buildings and sharp angles. His color combinations freely handled layers of paint, and the sparkling water and light sky exemplify his lifelong preoccupation with the effects and significance of light in an awe-inspiring natural world.

This painting is part of the Widener Collection at the National Gallery of Art, Washington, DC.



From *Venice: The Dogana and San Giorgio Maggiore*
Photograph © Board of Trustees, National Gallery of Art, Wash.



The entire Aldrichimica Acta archive is also available free of charge at Aldrich.com/acta.

Aldrichimica Acta (ISSN 0002-5138) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group.

© 2012 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA, SAFC, SIGMA-ALDRICH, ALDRICH, and SUPPLIC are trademarks of Sigma-Aldrich Co. LLC, registered in the US and other countries. Add Aldrich is a trademark of Sigma-Aldrich Co. LLC. Aldrich brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.

Want to run transition-metal-catalyzed reactions in *water* at room temperature?

Add Aldrich

Congratulations to Professor Bruce Lipshutz, University of California, Santa Barbara, winner of the 2011 *Presidential Green Chemistry Challenge Award: Academic*, for the development of TPGS-750-M, an amphiphile that enables transition-metal-catalyzed reactions to be run in *water*, at room temperature.












- Applicable to many transition-metal-catalyzed reactions
- Scalable
- Highly water-soluble
- No frothing at workup
- Easily separated from reaction product(s) and recycled
- Amenable to high-throughput screening



For greener chemistry in your lab, Add Aldrich.
Aldrich.com/tpgs750m

TPGS-750-M



Cat. No. 733857



Lipshutz, B. H. et al. *J. Org. Chem.* 2011, 76, 4379.

Easily activated Palladium Precatalysts for facile C–N and C–C bond formations.

Use of XPhos-Pd-G2 in the Suzuki coupling of
2,6-difluoroboronic acids



References:

(1) Surry, D. S.; Buchwald, S. L. *Chem. Sci.* **2011**, 2, 27. (2) Martin, R.; Buchwald, S. L. *Acc. Chem. Res.* **2008**, 41, 1461. (3) Surry, D. S.; Buchwald, S. L. *Angew. Chem. Int. Ed.* **2008**, 47, 6338. (4) Kinzel, T.; Zhang, Y.; Buchwald, S. L. *J. Am. Chem. Soc.* **2010**, 132, 14073.




For inquiries about bulk quantities,
contact catalysis@sial.com.
Aldrich.com/buchwald

Add **Aldrich**

Buchwald palladium precatalysts are highly reactive, air-stable crystalline white solids that generate the active monoligated Pd(0) species in situ, providing a more efficient and easily accessible catalytic system.

G1 = 1st Generation Precatalysts (Aliphatic Amine)

G2 = 2nd Generation Precatalysts (2-Aminobiphenylamine)



Professor Bruce H. Lipshutz

Winner of the 2012 EROS Best Reagent Award



Aldrich Chemistry and John Wiley, proud sponsors of the annual EROS Best Reagent Award, congratulate this year's winner, Professor Bruce Lipshutz. The award was created to honor outstanding contributors to the online edition of *Encyclopedia of Reagents for Organic Synthesis* (EROS and e-EROS).



About Prof. Bruce Lipshutz




Bruce Lipshutz has been on the faculty at UC Santa Barbara for the past 33 years. From his early contributions in the form of reagents such as SEM-Cl and higher order cyanocuprates to heterogeneous catalysts in the form of nickel- and copper-in-charcoal, his research focus has been on providing technologies that are broadly applicable


to synthetic problems. More recently, he has turned his attention to "green" chemistry, with the goal of developing micellar catalysis as an environmentally benign way of carrying out transition-metal-catalyzed cross-coupling reactions in water at room temperature.

About the Reagent

Lipshutz and co-workers have developed a new method for the asymmetric hydrosilylation of various carbonyl compounds using a preformed DTBM-SEGPHOS[®]-CuH species in the presence of excess polymethylhydrosiloxane (PMHS), with high yields and ee's reported for the reaction products. The resulting chiral building blocks are valuable precursors for the synthesis of various drugs and naturally occurring molecules. (Lipshutz, B. H. et al. *Org. Lett.* 2006, 8, 2969.)



762458



762938

Add **Aldrich**



VOL. 45, NO. 3 • 2012

 **ALDRICH**
Chemistry

Aldrichimica **ACTA**



Recent Advances in the Buchwald–Hartwig Amination
Reaction Enabled by the Application of Sterically
Demanding Phosphine Ancillary Ligands

Transition-Metal-Mediated Fluorination,
Difluoromethylation, and Trifluoromethylation



SIGMA-ALDRICH

Want a Cleaner Coupling Substrate?

Add **Aldrich**

Aldrich® Chemistry is proud to offer over 1,000 iodinated building blocks for your research. Iodine-containing substrates are valuable building blocks for a diverse array of synthetic methodologies. They can participate in various cross-coupling reactions for the generation of carbon-carbon, carbon-nitrogen and carbon-oxygen bonds. Iodinated substrates are often preferred to their brominated analogues, because of their cleaner reactions. They are also handy precursors for organolithium, organozinc and organomagnesium reagents.




779601



713457




731560



729396




756334




728594




778540



728586




728268




739758




755389




ADE000435




751251




ADE000340




716286



707120



ADE000915



ADE001132



For a complete list of available building blocks, visit Aldrich.com/bbacta

Aldrich Chemical Co. LLC.
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-231-8327
SAFC® 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4561
International 114-771-5765
24-Hour Emergency 114-776-6555
Website sigma-aldrich.com
Email aldrich@si-al.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 2988, Milwaukee, WI 53201, USA
sharbil.firsan@si-al.com

Subscriptions

Request your FREE subscriptions to the
Aldrichimica Acta, through:

Web: aldrich.com/acta
Email: sams-wa@si-al.com
Phone: 800-325-3010 (USA)
Mail: Acta: Marketing Services
Sigma-Aldrich Corporation
1450 Spruce Street
St. Louis, MO 63103



The entire *Aldrichimica Acta* archive is also available free of charge at aldrich.com/acta.

Aldrichimica Acta (ISSN 0002-1700) is a publication of Aldrich. Aldrich is a member of the Sigma-Aldrich Group.

©2012 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA, SAFC, SIGMA-ALDRICH, ALDRICH, and SUPELCO are trademarks of Sigma-Aldrich Co. LLC, registered in the US and other countries. All Aldrich is a trademark of Sigma-Aldrich Co. LLC. Aldrich brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.




"PLEASE BOTHER US."

Dear Fellow Chemists,

Professor Michael Willis of the University of Oxford suggested we offer DABSO, a charge-transfer complex between DABCO® and SO₂, and a convenient replacement for gaseous SO₂. It is a bench-stable, colorless, crystalline solid that has found numerous applications in the synthesis of sulfamides, sulfonamides, N-aminosulfonamides, and sulfolenes

(a) Emmett, E. J.; Richards-Taylor, C. S.; Nguyen, B.; Garcia-Rubia, A.; Hayter, B. R.; Willis, M. C. *Org. Biomol. Chem.* **2012**, *10*, 4007. (b) Ye, S.; Wu, J. *Chem. Commun.* **2012**, *48*, 7753. (c) Woolven, H.; González-Rodríguez, C.; Marco, I.; Thompson, A. L.; Willis, M. C. *Org. Lett.* **2011**, *13*, 4876. (d) Nguyen, B.; Emmett, E. J.; Willis, M. C. *J. Am. Chem. Soc.* **2010**, *132*, 16372

DABCO is a registered trademark of Air Products and Chemicals, Inc.



742937 1,4-Diazabicyclo[2.2.2]octane bis(sulfur dioxide) adduct (DABSO), 1 g
≥95% (sulfur, elemental analysis) 5 g

We welcome fresh product ideas from you. Do you need a compound that isn't listed in our catalog? Ask Aldrich! For more than 60 years, your research needs and suggestions have shaped the Aldrich product offering. To submit your suggestion, visit Aldrich.com/pleasebotherus

John Radke
Director of Marketing, Chemistry

TABLE OF CONTENTS

Recent Advances in the Buchwald-Hartwig Amination Reaction Enabled by the Application of Sterically Demanding Phosphine Ancillary Ligands	59
<i>Rylan J. Lundgren and Mark Stradiotto,*</i> Dalhousie University	
Transition-Metal-Mediated Fluorination, Difluoromethylation, and Trifluoromethylation	67
<i>Zhuang Jin, Gerald B. Hammond,* and Bo Xu,*</i> University of Louisville	

ABOUT OUR COVER

Forest of Fontainebleau (oil on canvas, 175.6 × 242.6 cm) was painted in 1834 by the French artist Jean-Baptiste-Camille Corot (1796-1875). Corot was a prolific painter with a talent that ranged widely from mainly landscapes to portraits, nudes, etchings, and other forms. He even dabbled in photography later on in his life. Unlike a number of his famous contemporaries, his rise to prominence in the artistic world was slow: He did not take up painting until his mid-twenties, and it wasn't until his mid-forties that he began to get the recognition he deserved, first from fellow artists, and then from art critics, official France, collectors, and the public. His artistic training consisted mainly of apprenticeships with Neoclassicist painters Achille-Etna Michallon and Jean-Victor Bertin and three trips a few years apart to Italy where he honed his landscape painting skills. Corot lived his life as a simple and humble man who was a mentor and a philanthropist to many young and struggling artists, and his work is thought to foreshadow the impressionist movement



Detail from *Forest of Fontainebleau* Photograph

© Board of Trustees, National Gallery of Art, Washington

Forest of Fontainebleau is a good example of the genre of painting known as "historical" landscape, which Corot excelled at, and in which a historical or biblical theme is incorporated into the painting in order to "elevate" it in the eyes of art critics and gain public attention. In this case, the biblical connection is inferred from several clues* Corot inserted in the painting. Corot's strong naturalist strain, controlled paint strokes, and tendency to favor brown and black colors are reflected in this unpretentious and simply rendered composition depicting a rough forest that stands in sharp contrast to the serene young woman

This painting is part of the Chester Dale Collection at the National Gallery of Art, Washington, DC.

* Can you find the clues in the painting? What do the clues reveal about the young woman reading? To find out visit Aldrich.com/acta/aoc453



Dalphos Ligands for Cross-Coupling

Add Aldrich


A Dalphos ligand contains the bulky di(1-adamantyl)phosphino [P(1-Ad),] fragment. These chelating N,P ligands are useful for Pd-catalyzed C–N and C–C bond formation. The more reactive Mor-Dalphos improves the scope and utility of ammonia coupling at room temperature, and is also effective in the coupling of hydrazine and acetone.

- Monoarylation of ammonia, hydrazine and acetone
- Air-stable
- Good functional group tolerance
- Mild conditions


Pd-Catalyzed Monoarylation of Ammonia, Hydrazine and Acetone



751618
Mor-DalPhos



751413
Me-DalPhos



759198
New

References:

- (1) Lundgren, R. J.; Peters, B. D.; Alsabeh, P. G.; Stradiotto, M. *Angew. Chem. Int. Ed.* **2010**, *49*, 4071. (2) Lundgren, R. J.; Stradiotto, M. *Angew. Chem. Int. Ed.* **2010**, *49*, 8686. (3) Hesp, K. D.; Lundgren, R. J.; Stradiotto, M. *J. Am. Chem. Soc.* **2011**, *133*, 5194. (4) Lundgren, R. J.; Sapping-Kumankumah, A.; Stradiotto, M. *Chem.—Eur. J.* **2010**, *16*, 1983. (5) Hesp, K. D.; Stradiotto, M. *J. Am. Chem. Soc.* **2010**, *132*, 18026.



To learn more about
DalPhos ligands, visit
Aldrich.com/dalposacta

New Equipment for Chemical Synthesis

Add Aldrich

Aldrich® Chemical Dryer

Dries Samples Using Heat, Vacuum and Desiccant

- Digital thermocouple control
- Sample viewing port
- Block heater enclosed for safety
- Removable drying tube assembly




Cat. No.	Product Description
Aldrich Chemical Dryer (Order Glass Drying Tube Assembly Separately)	
Z683507	120 V US plug
Z683515	230 V Euro plug
Z683523	230 V UK plug
Glass Drying Tube Assembly	
Z222712	Complete, ST/NS 34/35 joints

For product and ordering information, visit
Aldrich.com/chemicaldryer

Aldrich® Kugelrohr Bulb-to-Bulb Vacuum Distillation Apparatus

Quickly Distills the Most Difficult Materials with Minimal Loss

- Distills up to 230 °C
- 25 mL to 2 L flask capacity
- High vacuum capability
- Safe, front loading design
- Precise digital controls
- Oscillating drive motor



Cat. No.	Product Description
Kugelrohr Distillation Apparatus (Order Glassware Separately)	
Z683477	120 V US plug
Z683485	230 V Euro plug
Z683493	230 V UK plug
Kugelrohr Glassware Sets	
Z748293	25 and 100 mL flasks, ST/NS 14/20 joints
Z748307	1 L flasks, ST/NS 24/40 joints
Z748315	2 L flasks, ST/NS 24/40 joints

For product and ordering information, visit
Aldrich.com/kugelrohr

Speed, Convenience, Service

ALDRICH MARKET SELECT

Aldrich® Market Select is a comprehensive solution for identification and procurement of readily available screening compounds and building blocks. The structure database and purchasing website include regularly updated catalog data from more than 60 of the most reliable suppliers around the world.


We make compound and building block procurement easy!

- Quote & Order Management
- Custom Formatting & Packaging
- International Shipping and Compliance
- Consolidate Invoices and Supplier Payment
- Post-procurement Follow-up Support


Purchase your chemistry with ease at
www.AldrichMarketSelect.com




SEARCH



SELECT



PURCHASE



VOL. 46, NO. 1 • 2013

ALDRICH
Chemistry

Aldrichimica **ACTA**



Polymer-Incarcerated Metals: Highly Reactive,
Recoverable, and Multifunctional Nanocluster Catalysts
for Organic Synthesis

Advances in Acylation Methodologies Enabled by
Oxyima-Based Reagents

SIGMA-ALDRICH


Baran Reagents – diversification with ease.

Add Aldrich


The Baran Group has teamed with Aldrich® Chemistry to provide a toolkit of diversification reagents which are:

- Bench-stable, free-flowing solids
- Highly efficient under open-flask conditions
- Compatible with unprotected late-stage intermediates
- Able to perform selective one-pot sequential alkylations


Baran Reagents allow for easy diversification of heterocycles as described in the recent *Nature* paper "Practical and Innate Carbon-Hydrogen Functionalization of Heterocycles" by Professor Baran and co-workers.




DFMS
767840




TFMS
771406



DFES-Na
L511404




IPS
L511161



TFES
L511234

Toolkit for Easy Diversification



Reference:
Fujiwara, Y.; Dixon, J. A.; O'Hara, F.; Funder, E. D.; Dixon, D. D.; Rodriguez, R. A.; Baxter, R. D.; Herle, B.; Sach, N.; Collins, M. R.; Ishihara, Y.; Baran, P. S. *Nature* 2012, 492, 95.



To find out more about the Baran reagents available from Aldrich Chemistry, visit Aldrich.com/baranreagents

Aldrichimica ACTA

VOL. 46, NO. 1 • 2013

Aldrich Chemical Co. LLC
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-325-5837
SAFC® 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 114-771-5765
24-Hour Emergency 114-776-6555
Website sigma-aldrich.com
Email aldrich@sigmaaldrich.com

General Correspondence

Editor: Sharbil J. Friesan, Ph.D.
P.O. Box 2060, Milwaukee, WI 53201, USA
sharbil.friesan@sigmaaldrich.com

Subscriptions

Request your FREE subscription to the
Aldrichimica Acta, through:

Web: Aldrich.com/acta
Email: sams-usa@sigmaaldrich.com
Phone: 800-325-3010 (USA)
Mail: Attn: Marketing Services
Sigma-Aldrich Corporation
3050 Spruce Street
St. Louis, MO 63118



The online Aldrichimica Acta website is available at www.sigmaaldrich.com/acta.

Aldrichimica Acta (ISSN 0031-5100) is a publication of Aldrich, a member of the Sigma-Aldrich Group.

©2013 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA, SAFC, SIGMA-ALDRICH, ALDRICH, and SURE/SEAL are trademarks of Sigma-Aldrich Co. LLC, registered in the US and other countries. Aldrich is a trademark of Sigma-Aldrich Co. LLC. Aldrich brand products are sold through Sigma-Aldrich, Inc. Purchases must demonstrate the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.



"PLEASE BOTHER US."

Dear Fellow Chemists,

Professor Thomas R. Hove of The University of Minnesota recently suggested that we introduce convenience packs of smaller Sure/Seal™ bottles for reagent solutions that are air- or moisture-sensitive. We have already added over 30 products that customers can purchase in convenience packs of 4 x 25 mL Sure/Seal™ bottles, and more are being added daily. If there are other products you would like to have in 4 x 25 mL convenience packs, please let us know at Aldrich.com/25mlsureseal; we look forward to hearing from you!

Sure/Seal is a trademark of Sigma-Aldrich Co. LLC.



We welcome fresh product ideas from you. Do you need a compound that isn't found on our website or in our catalog? Ask Aldrich! For more than 60 years, your research needs and suggestions have shaped the Aldrich product offering. To submit your suggestion, visit Aldrich.com/pleasebotherus.

Scott Batcheller
R&D and Marketing Manager
Chemical Synthesis

TABLE OF CONTENTS

Polymer-Incarcerated Metals: Highly Reactive, Recoverable, and Multifunctional Nanocluster Catalysts for Organic Synthesis 3
Shū Kobayashi and Hirovuki Mivamura*, The University of Tokyo

Advances in Acylation Methodologies Enabled by Oxyma-Based Reagents..... 21
Ramon Subirós-Funosas, Sherine N. Khattab, Lidia Nieto-Rodriguez, Ayman El-Faham*, and Fernando Albericio** IRB Barcelona, CIBER-BBN Barcelona, Alexandria University, King Saud University, University of Barcelona, and University of KwaZulu-Natal

ABOUT OUR COVER

As a young man, van Gogh worked in various jobs in Holland, England, France, and Belgium. Only at age 27, in 1880, did he start his formal art training in Brussels, and continued his training in Holland, where he discovered the work of Rubens, before moving permanently to France in 1886. In Paris, van Gogh studied with Cormon; became acquainted with Pissarro, Monet, and Bernard; and befriended Gauguin. Between 1880 and 1890, van Gogh created close to 900 paintings, before his career was cut short by his suicide at the age of 37. He lived a difficult life: unhappy romances, lack of commercial success, relying on family for financial support, and severe mental illness most of his adult life. It is ironic that an individual who would influence the art world so greatly after his death thought of his life as a terrible waste and of himself as an utter failure.



Farmhouse in Provence
National Gallery of Art, Washington

Van Gogh's early "Dutch period" painting style, with somber tones and dark colors dominating, reflected the artistic influences he came under while in Holland and his strong interest in drawing and painting figures. However, he quickly abandoned the dark colors in favor of brighter, more vibrant colors and bold brushstrokes after moving to Paris and getting introduced to the impressionist and post-impressionist styles. *Farmhouse in Provence* (oil on canvas, 46.1 x 60.9 cm), which was completed in June 1888 after van Gogh relocated to the South of France, is a good example of this later style and of his belief that color ought to be the primary medium of expression in a painting. His love of the various shades of yellow* and his strive for realism are evident in the rich golden tones that dominate this painting. Also evident are his quick and thick brushstrokes and the swirling cloud lines that foreshadow the more dramatic swirling-cloud pattern in his masterpiece of about a year later, *The Starry Night*. His use of juxtaposed complementary colors mirrors the technique used by impressionists to enhance the vividness of their works.

This painting is part of the Ailsa Mellon Bruce Collection at the National Gallery of Art, Washington, DC.

What could be behind van Gogh's affinity for the color yellow? To find out, visit Aldrich.com/acta461

Air-Stable Palladacycles.

Add Aldrich


Aldrich® Chemistry offers new second- and third-generation palladium precatalysts that improve reactivity in Pd-catalyzed cross-coupling reactions. These air-stable catalysts have been extensively applied in the synthesis of pharmaceuticals, natural products, polymers and new materials.

- Highly efficient
- Utilized under mild conditions
- Short reaction times
- Low catalyst loadings


Gram to ton quantities available.
For bulk quantity inquiries, contact: catalysis@sial.com




To learn more about our portfolio of Buchwald Precatalysts and Ligands, visit Aldrich.com/buchwaldguide




RuPhos-Pd-G3
763403




BrettPhos-Pd-G3
761605




XPhos-Pd-G3
763381




JackiePhos-Pd-G3
762830




CPhos-Pd-G3
763004




tBuXPhos-Pd-G3
762229




XantPhos-Pd-G3
763039




CPhos-Pd-G2
763012



P(o-Tol)₃-Pd-G2
759163




SPhos-Pd-G2
753009




P(Cy)₃-Pd-G2
756741




P(t-Bu)₃-Pd-G2
756482



XantPhos-Pd-G2
763047



cataCXium®-Pd-G2
761311



cataCXium®-Pd-G3
761435

Reference:
Bruno, N. C.; Tudge, M. T.; Buchwald, S. L. *Chem. Sci.* 2013, 4, 916

New Glass Filter Dryers.

Add **Aldrich**

Save valuable research time with filtration and drying in a single step

The Glass Filter Dryer (GFD™) is an innovative piece of laboratory or pilot plant equipment. This scaled-down version of a production filter dryer allows easy lab-scale filtration and drying with the same reliability and benefits as full-size filter dryer products from Powder Systems Ltd. and larger Neutech filters.

The GFD technology supersedes old style laboratory filtration and drying methods such as Buchner filtration and tray drying.

- Works under vacuum and pressure
- Heated jacket for drying
- Improved yield
- Easy cleaning
- Reduced product exposure
- Batch homogeneity
- Process visibility
- Minimum process disruption
- Removable filtration basket for easy collection of product



GFDs are shipped ready-to-use as shown with all glassware, supports and bi-directional stirring motors.

GFDs are available in three sizes

	Vessel Volume
MINI LAB Glass Filter Dryer	0.3 L
LAB Glass Filter Dryer	2 L
MAXI LAB Glass Filter Dryer (not shown)	10 L



For additional product and ordering information, visit sigma-aldrich.com/glass-filter-dryers




P.O. Box 14508
St. Louis, MO 63178
USA

Discover the latest developments in Chemistry research.

Add  Aldrich™

Get your FREE subscription to the *Aldrichimica Acta* today.

For over 45 years, the *Aldrichimica Acta* has been an international forum for the frontiers of chemical research. Articles, written by chemists from around the world, cover a variety of topics based on a synthetic theme involving organic, organometallic, bio-organic or inorganic chemistry. The *Aldrichimica Acta* has been ranked **number one by Impact Factor** in nine of the past ten years in the field of organic chemistry.*



To subscribe, visit
[Aldrich.com/chemistry/
subscriptions.html](http://Aldrich.com/chemistry/subscriptions.html)

*As ranked in the Journal Citation Reports®, Annual Science Editions
© 2011 Aldrichimica Acta, a trademark of Thomson Reuters (Scientific) Inc.



PHP
19881
1043

SIGMA-ALDRICH

VOL. 46, NO. 2 • 2013

ALDRICH
Chemistry

Aldrichimica **ACTA**



Aminoacyl Benzotriazolides: Versatile Reagents for the
Preparation of Peptides and Their Mimetics and Conjugates

Recent Advances in the Prins Cyclization

SIGMA-ALDRICH

Access the most diverse product portfolio



















Add Aldrich

Aldrich has partnered with numerous trusted suppliers to continually add rare and hard-to-find building blocks and synthetic reagents to our portfolio. For example, we just added more than 1,000 high-quality building blocks from Synthonix, BioBlocks, Indofine, and Otava. Browse these new additions or discover the power of our new, upgraded Structure Search that makes finding the right product easier.

- More than 40,000 quality building blocks and synthetic reagents in stock and ready for immediate delivery
- Easily search our entire portfolio using JME, MarvinSketch™, or ChemDraw® editors



For more information, visit
Aldrich.com/affiliateacta



Aldrich Chemical Co. LLC.
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone: 800-325-3010 (USA)
FAX: 800-325-5052 (USA)
or 414-438-2199
P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries: 800-325-3010
Technical Service: 800-325-5812
SAFC®: 800-244-1173
Custom Synthesis: 800-244-1173
Flavors & Fragrances: 800-227-4563
International: 314-771-5765
24-Hour Emergency: 314-776-6555
Website: sigma-aldrich.com
Email: aldrich@sigmaal.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 2060, Milwaukee, WI 53201, USA
sharbil.firsan@sigmaal.com

Subscriptions

Request your FREE subscription to the Aldrichimica Acta through:

Web: Aldrich.com/acta
Email: sams-usa@sigmaal.com
Phone: 800-325-3010 (USA)
Mail: Attn: Marketing Services
Sigma-Aldrich Corporation
1050 Spruce Street
St. Louis, MO 63103



The entire Aldrichimica Acta archive is available at Aldrich.com/acta.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich, a member of the Sigma-Aldrich Group.

©2013 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA, SAFC, SIGMA-ALDRICH, ALDRICH, and SUPRECO are trademarks of Sigma-Aldrich Co. LLC, registered in the US and other countries. Aldrich is a trademark of Sigma-Aldrich Co. LLC. Aldrich brand products are sold through Sigma-Aldrich, Inc. Purchases must determine the suitability of the product for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.





"PLEASE BOTHER US."

Dear Fellow Chemists.

Professor Meike Niggemann of RWTH Aachen University, Germany, kindly suggested that we offer Ca(NTf₂), a main-group-metal catalyst used to effect a number of synthetic transformations that typically require stoichiometric quantities of other Lewis acids. When utilized together with a catalytic quantity of Bu₄NPF₆, Ca(NTf₂), can effectively catalyze regio- and stereoselective Friedel-Crafts alkylation,¹ hydroarylations of alkenes,² and cyclopropanations

(1) Begouin, J.-M.; Capitta, F.; Wu, X.; Niggemann, M. *Org. Lett.* **2013**, *15*, 1370. (2) Niggemann, M.; Meel, M. *J. Angew. Chem., Int. Ed.* **2010**, *49*, 3684. (3) Niggemann, M.; Bisek, N. *Chem.—Eur. J.* **2010**, *16*, 11246. (4) Haven, T.; Kubik, G.; Haubenreisser, S.; Niggemann, M. *Angew. Chem., Int. Ed.* **2013**, *52*, 401.



753122	Calcium(II) bis(trifluoromethanesulfonimide), 95%	1 g
281026	Tetrabutylammonium hexafluorophosphate, 98%	5 g 25 g 100 g

We welcome your product ideas. Do you need a compound that is not featured on our website or in our catalog? Ask Aldrich! For more than 60 years, your research needs and suggestions have shaped the Aldrich product offering. Email us your suggestion at techserv@sigmaal.com

Josef Zihlmann
Vice President of Product Management

TABLE OF CONTENTS

Aminoacyl Benzotriazolides: Versatile Reagents for the Preparation of Peptides and Their Mimetics and Conjugates 43
*Siva S. Panda, C. Dennis Hall, Eric Scriven, and Alan R. Katritzky,** The University of Florida and King Abdulaziz University (Saudi Arabia)

Recent Advances in the Prins Cyclization 59
Sandro José Greco, Rodolfo Goetze Fiorot, Valdemar Lacerda, Jr., and Reginaldo Bezerra dos Santos,* Federal University of Espírito Santo (Brasil)

ABOUT OUR COVER

Winslow Homer completed *Breezing Up (A Fair Wind)* (oil on canvas, 61.5 x 97 cm) in 1876, the U.S. centennial year. One of the most celebrated American artists of the nineteenth and early part of the twentieth centuries, Homer had little formal training as a painter, and seems to have developed his plein-air style late in life. In *Breezing Up*, Homer depicts a sailor who excelled in several media and travelled widely, depicting a variety of themes, from battlegrounds of the Civil War to the surely scenes of recreation.



Detail from *Breezing Up (A Fair Wind)*
National Gallery of Art, Washington

An iconic American painting and one of Homer's most celebrated works, *A Fair Wind* received wide acclaim when it was first exhibited in 1876, and was chosen by the U.S. Postal Service in 1962 for a commemorative postage stamp honoring the artist. This scene touches upon a favorite theme of his, the struggles of man against powerful natural forces. Here the light sail boat is returning home with the day's catch and a seemingly relaxed crew, unbothered or perhaps accustomed to the choppy waters. Unlike some of his later seascapes, this has a warm feel and more optimistic message,* perhaps meant to make a statement about the future of the young country following the uncertain and dangerous years of the Civil War.

This painting is a gift of the W. L. and May T. Mellon Foundation to the National Gallery of Art, Washington, DC.


* Homer included a traditional symbol of hope in this painting. Can you guess what it is? To find out, visit Aldrich.com/acta462

To efficiently synthesize peptides and peptidomimetics


Add **Aldrich**[™]

Benzotriazole amino acids, or aminoacyl benzotriazolides, are versatile and efficient acylating reagents for synthesizing peptides and peptidomimetics. This methodology provides solid- and solution-phase techniques for accessing complex peptides and peptide conjugates and has been effective in many applications including:


- Polypeptidal benzotriazolides¹
- Peptidomimetics, such as aminoxypeptides², depsipeptides³ and heterocyclic peptidomimetics⁴
- Tagged peptides with fluorescent dyes as well as other labels⁵
- N-, O-, S-, and C-linked peptide conjugates⁶




737747
Z-Ala-Bt




737798
Fmoc-Leu-Bt




737623
Z-Val-Bt




737801
Fmoc-Gly-Bt




737763
Fmoc-Pro-Bt




737720
Z-Gly-Bt



737739
Z-Ile-Bt



737712
Z-Phe-Bt



737704
Z-Trp-Bt



For more information, visit
Aldrich.com/bztaacta

References:

- (1) Abdelmajid, A.; Tala, S. R.; Amine, M. S.; Katritzky, A. R. *Synthesis* 2011, 2995.
- (2) Katritzky, A. R.; Avani, T.; Tala, S. R. *J. Org. Chem.* 2009, 74, 8690.
- (3) Avan, I.; Tala, S. R.; Steel, P. J.; Katritzky, A. R. *J. Org. Chem.* 2011, 76, 4884.
- (4) Hansen, F. K.; Beagle, L. K.; Todadze, E.; Katritzky, A. R. *Heterocycles* 2012, 84, 515.
- (5) Katritzky, A. R.; Narindoshvili, T.; Angrish, P. *Synthesis* 2008, 2013.
- (6) Katritzky, A. R.; Angrish, P.; Todadze, E. *Synlett* 2009, 2392.



Aminoacyl Benzotriazolides: Versatile Reagents for the Preparation of Peptides and Their Mimetics and Conjugates



Dr. S. S. Panda



Dr. C. D. Hall



Dr. E. Scriven



Prof. A. R. Katritzky

Siva S. Panda,^a C. Dennis Hall,^a Eric Scriven,^a and Alan R. Katritzky^{*,a,b}

^a Center for Heterocyclic Compounds
Department of Chemistry
University of Florida
Gainesville, FL 32611-7200, USA
Email: katritzky@chem.ufl.edu

^b Chemistry Department
King Abdulaziz University
Jeddah 21589, Saudi Arabia

Keywords. benzotriazole; peptide synthesis; peptidomimetics; conjugates; methodology; coupling reagent.

Abstract. *N*-(Protected α -aminoacyl)benzotriazoles are efficient acylating reagents that offer many advantages in the preparation of peptides and their mimetics and conjugates. Advances in methodology, made possible by these novel reagents, have given rise to solution- and solid-phase preparative techniques for generating complex peptides and peptide conjugates, which are useful in the construction of diverse libraries of building blocks for medicinal chemistry.

Outline

1. Introduction
2. Synthesis of Aminoacylbenzotriazolides
 - 2.1. From Fmoc- or Cbz-Protected Amino Groups
 - 2.2. From Boc-Protected Amino Groups
3. Synthesis of Oligopeptidoyl Benzotriazolides
4. Application to the Synthesis of Peptides
 - 4.1. Natural Peptides and Isopeptides
 - 4.2. Difficult-to-Prepare Peptide Sequences
 - 4.3. Cyclic Dipeptides (2,5-Diketopiperazines) and Tripeptides
5. Synthesis of Peptidomimetics
 - 5.1. Aminoxy peptides
 - 5.2. Depsipeptides
 - 5.3. Azapeptides
 - 5.4. Hydrazinopeptides
 - 5.5. Heterocyclic Peptidomimetics
 - 5.6. Cyclic Peptidomimetics
6. Synthesis of Tagged Peptides and Peptidomimetics
 - 6.1. Fluorescent Labels
 - 6.1.1. Coumarin-Labeled Peptides and Peptidomimetics

- 6.1.2. 6-Chloro-2,3-naphthalimides and Water-Soluble Fluorescent Tags
- 6.2. Azo-Dye-Labeled Peptides
7. Peptide Conjugates
 - 7.1. Conjugates of Sugars
 - 7.2. Conjugates of Heterocycles
 - 7.3. Conjugates of Vitamins
 - 7.4. Conjugates of Pharmaceuticals
 - 7.5. Conjugates of Plant Hormones
 - 7.6. Aminoxy Acid Conjugates of Peptidomimetics and Hydrazino Acid Conjugates
8. Differential *N*-, *S*-, and *O*-Acylation
 - 8.1. *S*- and *O*-Acyl Isopeptides
 - 8.2. Ligation at a Distance
 - 8.2.1. *S* to *N* Acyl Migration
 - 8.2.2. *O* to *N* Acyl Migration
 - 8.2.3. *N* to *N* Acyl Migration
9. Conclusions and Comparison with Alternative Methodologies
 - 9.1. Carboxyl Group Activation by Isolation of an Intermediate
 - 9.2. Carboxyl Group Activation without Isolation of an Intermediate
10. References


1. Introduction

Azolides are compounds in which an azole residue is attached to an acyl group. It has long been known that the azole nucleus in azolides can behave as a leaving group, and this property was explored widely by Staab as early as 1960.^{1,2} Katritzky and co-workers have reported that, of the azoles, benzotriazole is a particularly versatile synthetic auxiliary with attractive properties, easily inserted into or removed from molecules and endowing them with a range of useful reactivities.³ Much of this work has been summarized⁴ and, from 1985 to 2012,


some 1,000 publications have appeared dealing with benzotriazole. We have since found that benzotriazolides (*N*-acylbenzotriazoles) have many advantageous properties relative to acid chlorides. More recently, the utility of benzotriazolides in peptide chemistry has achieved prominence, and reports of numerous applications have appeared as highlighted in this review.

2. Synthesis of Aminoacylbenzotriazoles


Acybenzotriazoles are versatile reagents in which the benzotriazol-1-yl (Bt) group serves as a surrogate for halogen and is easily displaced by nitrogen, sulfur, oxygen, or carbon nucleophiles.⁴ Acylbenzotriazoles,



Scheme 1. Preparation of Fmoc- or Cbz-Protected α -Aminoacylbenzotriazoles. (Ref. 5,9)



Scheme 2. Preparation of Boc-Protected α -Aminoacylbenzotriazoles from Acid-Sensitive or Boc-Protected Amino Acids. (Ref. 6,7)



eq 1 (Ref. 11)

however, offer numerous advantages over their halogen analogues since they are isolated in high yields, in crystalline form, are usually stable in air (and even to water for short periods at 20 °C), and are more reactive than the corresponding *N*-acylimidazoles. Thus, they can effect peptide coupling in H₂O-THF or H₂O-MeCN by using unprotected amino acids with the distinct advantage that, if the temperature is controlled, chirality within the component amino acids is preserved. Although synthesis of aminoacylbenzotriazoles requires protection of the α -amino function with either Boc, Fmoc, or Cbz (or by protonation); other common functional groups, especially OH but also SH, CONH₂, and CO₂H can be left unprotected.

2.1. From Fmoc- or Cbz-Protected Amino Groups

The first general method for the preparation of aminoacylbenzotriazoles involves the condensation of a protected amino acid with benzotriazole and thionyl chloride (1.0–1.2 equiv) in THF or DCM (Scheme 1, Part (a)).^{5–9} Excess benzotriazole is required to neutralize the two equivalents of HCl formed, and the method is not applicable to benzotriazolides containing the acid-sensitive Boc group. A wide range of amino acids may be used affording high yields. In addition, di-Bt derivatives are generated when aspartic acid, glutamic acid, or the S–S dimer of cysteine comprise the starting amino acids (Scheme 1, Parts (b) and (c)). The method is extremely versatile, and excess benzotriazole is easily removed by washing with either acid or base. Most significantly, chirality within the starting amino acids is retained in the majority of cases and the protected aminoacylbenzotriazoles may be stored at room temperature for many weeks.

2.2. From Boc-Protected Amino Groups

In cases where the amino acid is Boc-protected or sensitive to thionyl chloride, the sodium or trialkylammonium salt of the amino acid may be converted into the benzotriazolide by treatment with 1-(methanesulfonyl)benzotriazole (BtO₂SMe or BtMs) (Scheme 2).^{6,7} This clean, preparative method is enhanced in some cases by crystallization of the triazolide from water with concurrent removal of the water-soluble methanesulfonate byproduct.

The thionyl chloride method can often be employed to prepare *N*-protected α -aminoacylbenzotriazoles (61–99%) without protection of potentially reactive side chains such as aliphatic OH (Ser, Thr), aromatic OH (Tyr), thiol SH (Cys), indole NH (Trp), and amide NH (Asn or Gln).^{7–9} There are, however, amino acids (e.g., His, Glu, Lys, and Asp) with side chains that do require protection (e.g., with Ts, Bn, or Cbz groups) in order to generate good yields of the benzotriazolides by reaction with either BtH-SOCl₂ or BtMs.¹⁰

3. Synthesis of Oligopeptidoyl Benzotriazolides

All of the methods employed to prepare Fmoc-, Cbz-, and Boc-protected aminoacylbenzotriazoles may also be used to prepare *N*-protected benzotriazolides of dipeptides. Likewise, *N*-protected tri- through hexapeptides may be converted into the corresponding benzotriazolides (eq 1).¹¹ Each benzotriazolide may then couple with amino acids or peptides to form *N*-protected tetra- through heptapeptides (vide infra).

4. Application to the Synthesis of Peptides

4.1. Natural Peptides and Isopeptides

A wide range of amino acids including many of those with additional, unprotected functional groups (Ser, Tyr, Cys, Trp, Pro, Asp, Glu, Lys, and Arg) couple with *N*-protected aminoacylbenzotriazoles in aqueous acetonitrile (MeCN-H₂O, 7:3) at 20 °C to produce dipeptides in 47–98% yields (eq 2).^{9,12–17} Enantiopure dipeptides (LL and LD)

are obtained in high purity (>99% in most cases) without the use of chromatography, thus highlighting the significant utility of the method. However, protection of a carboxylic acid function as its benzyl ester is advantageous in some instances.

The synthesis of tripeptides is achieved by either a fragment-coupling procedure^{9,12,16} or by stepwise coupling.^{12,15} In the former, N-protected dipeptides are converted into their benzotriazole derivatives at -10 °C (in order to avoid racemization), which are then coupled with amino acids to form tripeptides (Scheme 3).^{15,18,19} Tri- and tetrapeptides are similarly prepared in good yields (74–94%) by a stepwise procedure, but usually at a lower temperature (-10 °C) again to prevent racemization.⁹


The first examples of penta-, hexa-, and heptapeptides prepared by the Bt technology were generated using microwave-assisted, solid-phase peptide synthesis (SPPS). In this technique, N-Fmoc- α -aminoacyl benzotriazolides were employed to attach Fmoc-protected α -aminoacyl groups to a Rink resin, which was then utilized for the synthesis of di- to heptapeptides in 20–68% yields.^{18,20,21} In contrast, and although useful, SPPS and enzymatic techniques often afford low yields,^{22,23} and are unsuitable for larger-scale preparations.

Recently, N-protected α -tri-, tetra-, and pentapeptidylbenzotriazoles were coupled with free amino acids, dipeptides, and tripeptides in aqueous acetonitrile at 0 °C to afford N-terminal-protected polypeptides in isolated yields of 61–92%.¹¹ This reaction has been applied to amino acids (or small peptides) containing free OH, SH, or indole NH groups, and proceeds with no detectable racemization (eq 3);¹¹ it applies equally well to the synthesis of isopeptides.²⁴

4.2. Difficult to Prepare Peptide Sequences

Peptaibols, a group of antibiotics isolated from soil fungi, contain hindered amino acids such as 2-methylalanine (Aib), 2-ethylalanine, and 2,2-diethylglycine and, hence, α -substituted and α,α -disubstituted amino acids are important in peptide and medicinal chemistry.²⁵ Likewise, N-substituted peptides are important since they are constituents of cyclosporins²⁶ and exhibit antibiotic,²⁷ anticancer,²⁸ and antiviral²⁹ activities. Benzotriazole technology has proved valuable in the synthesis of hindered dipeptides from either N-protected Aib benzotriazolides and amino acids (Scheme 4, Part (a)) or C-terminus Aib dipeptides from N-protected aminoacylbenzotriazoles (Scheme 4, Part (b)) in isolated yields of 67–92%.¹⁷

It is well known that certain peptide sequences, particularly those containing valine (e.g., H-Leu-Met-Val-Gly-Gly-Val-Ile-Ala-NH₂), are difficult to prepare, and the classical approaches are often characterized by low yields, aggregation, and β -sheet formation leading to racemization. The stepwise SPPS synthesis of difficult peptides utilizing N-protected aminoacylbenzotriazoles and microwave acceleration has been shown to facilitate amide-bond formation (22–37% yields, 89 to >99% purities) and to reduce aggregation.²¹




PG = Boc, Cbz, Fmoc

R¹ = Me, *i*-Pr, Bn, MeS(CH₂)₂, H₂N(O)C(CH₂)₂, (indol-3-yl)CH₂

R² = Me, *i*-Pr, Bn, HOCH₂, HSCH₂, MeS(CH₂)₂, H₂N(O)C(CH₂)₂, (indol-3-yl)CH₂

eq 2 (Ref. 9,12)




PG = Cbz, Fmoc; R¹ = Me, Bn, MeS(CH₂)₂, (indol-3-yl)CH₂

R² = Me, Bn, MeS(CH₂)₂, H₂N(O)C(CH₂)₂, (indol-3-yl)CH₂


R³ = Me, *s*-Bu, Bn, HOCH₂, HSCH₂, MeS(CH₂)₂, (indol-3-yl)CH₂

Scheme 3. Preparation of Tripeptides by Fragment Coupling. (Ref. 15,18,19)




Example of the synthesis of N-terminal-protected polypeptides

eq 3 (Ref. 11)



PG = Cbz, Fmoc; AA = L-Phe-OH, L-Trp-OH, L-Met-OH




L-Phe-Bt, L-Trp-Bt, L-Met-Bt, Gly-Bt


Scheme 4. Peptide Sequences with Sterically Hindered Amino Acid Residues from N-(Cbz- α -aminoacyl)benzotriazoles

4.3. Cyclic Dipeptides (2,5-Diketopiperazines) and Tripeptides

Solution- and solid-phase (Scheme 5, Part (a)) Staudinger-type cyclizations afford efficient routes to hetero-2,5-diketopiperazines from protected azido dipeptide thioesters under microwave irradiation.³⁰ However, attempts to synthesize cyclic tripeptides by this method resulted in the unprecedented cleavage of an amide group rather than a thioester to form 2,5-diketopiperazines again (Scheme 5, Part (b)),³⁰ an observation that highlights the stability of a six-membered ring over its nine-membered analogue. Seven- and eight-membered-ring cyclic dipeptides can, however, be prepared in moderate-to-good yields by a Staudinger-type ring closure of a series of azido dipeptide thioesters (Scheme 5, Part (c)).³¹ The work was extended to the solution-phase synthesis of a ten-membered-ring cyclic tripeptide in 48% yield (Scheme 6).³¹



Scheme 5. Attempted Syntheses of Cyclic Tripeptides, Leading Instead to Cyclic Dipeptides (2,5-Diketopiperazines). (Ref. 30,31)



Scheme 6. Solution-Phase Synthesis of a Ten-Membered-Ring Cyclic Tripeptide. (Ref. 31)

5. Synthesis of Peptidomimetics


Peptidomimetics are small, protein-like molecules designed to mimic natural peptides by replacement of an amide bond or other element of the natural peptide. Clinical applications of bioactive natural peptides, for instance as hormones or enzyme inhibitors, have been limited by their susceptibility to rapid hydrolysis by peptidases. The corresponding peptidomimetics are not subject to this limitation, and, consequently, have been designed to exhibit high affinity for specific receptors, good metabolic stability toward endogenous proteases, greater oral bioavailability, and longer duration of action. To meet the need for good synthetic approaches to these peptidomimetics, flexible, high-yield, enantiospecific benzotriazole-mediated synthetic routes to six different structural types of peptidomimetics have been developed. These syntheses employ microwave-assisted benzotriazole acylation as the key step.^{32–34}

5.1. Aminoxypeptides

α -Aminoxy acids [RCH(OH)CO₂H] are peptidomimetics that resist enzymatic degradation. Peptidomimetics as a class are of interest as bioisosteric α -amino acids and as analogues of β -amino acids.^{35–38} They have been used to prepare aminoxypeptides, α -aminoxy- α -hybrid dipeptides, and α -aminoxy- α,α -hybrid tripeptides from unprotected amino acids. This methodology has been extended to other aminoxy-hybrid dipeptides and tripeptides starting from protected aminoxyacyl benzotriazolides (Scheme 7).³⁸

5.2. Depsipeptides

Depsipeptides contain both amino acid and hydroxy acid units with amide and/or ester bonds. Natural depsipeptides exhibit antifungal, antimicrobial, and anti-inflammatory activities, and certain depsipeptides have been used in cancer treatment.³⁹ A comparative study of standard coupling agents used to produce depsipeptides revealed variable yields and often long coupling times.⁴⁰ *N*-Cbz(depsidipeptidoyl)benzotriazoles were found to be useful for coupling with amino acids (*N*-acylation)



Scheme 7. Preparation of Aminoxy and Other Peptides from α -Aminoxy Benzotriazolides. (Ref. 38)

and α -hydroxy acids (O-acylation) to give depsidi- and depsitriptides in good yields under mild conditions (Scheme 8).⁴¹ This approach was applied to the preparation of unprotected depsiptides, which were converted into unprotected depsiptides.


5.3. Azapeptides

Azapeptides are peptidomimetics in which the α -CH group of one or more amino acid residues is replaced by a nitrogen atom. They are of interest for the generation of receptor ligands, enzyme inhibitors, and clinically approved drugs,⁴² and those with electrophilic moieties act as inhibitors of cysteine proteases.^{43,44} *N*-(*N*-Pg- α -Azadipeptidoyl)-


benzotriazoles have been prepared from amino acid esters in four steps and 48–77% overall yields and utilized for the synthesis of *N*-Pg-azatriptides, *N*-Pg-azatetraptides, and hybrid azapeptides (Scheme 9).⁴⁵ This methodology proved valuable for the insertion of an aza-amino acid unit into an azatriptide chain for the synthesis of the previously unknown aza-analogues of the endogenous opioid peptide neurotransmitter Leu-enkephalin, found in animals and humans.⁴⁵

5.4. Hydrazinopeptides

The replacement of an α -amino acid unit by a β -amino acid unit is a well-known strategy in the search for pharmacologically active peptides,⁴⁶




Scheme 8. Depsiptides from Cbz(depsiideptidoyl)benzotriazoles. (Ref. 41)




Scheme 9. Preparation of Azapeptides from *N*-(*N*-Pg- α -Azadipeptidoyl)benzotriazoles. (Ref. 45)

and replacement of the C^α and/or the C^β atom in such β-amino acid building blocks by a hetero atom offers an attractive extension of the β-peptide concept.⁴⁷

An alternate pathway to chirally pure α-hydrazino acids⁴⁸ is based on microwave irradiation during the conversion of α-bromo acids into α-hydrazino acids by hydrazine hydrate. The corresponding benzotriazolides afford hydrazine hybrid peptides (Scheme 10).⁴⁹



Scheme 10. Preparation of Hydrazinodipeptides. (Ref. 45)



Scheme 11. Microwave-Assisted Acylation of (a) 3,5-Diamino-1,2,4-triazole and (b) the Exocyclic Amino Group of 3,5-Diamino-1,2,4-triazole. (Ref. 52)

5.5. Heterocyclic Peptidomimetics

1,2,4-Triazoles have been employed for the bioisosteric replacement of the amide bond,⁵⁰ since 1,2,4-triazoles exhibit a wide range of antifungal and antibacterial activities.⁵¹ 3,5-Diamino-1,2,4-triazole has been coupled to di- and tripeptides at either ring or exocyclic nitrogens, using the benzotriazole methodology, to give potential building blocks for the preparation of peptidomimetics (Scheme 11).⁵²

5.6. Cyclic Peptidomimetics

In 2009, the benzotriazole methodology was extended to achieve cysteine S-acylation under mild conditions.⁵³ This has now been employed to couple an *N*-acylbis(benzotriazole) with cysteine to give a bis(*S*-acylcysteine), which, upon treatment with another equivalent of *N*-acylbis(benzotriazole), affords cyclic peptide mimetics in 64–72% yields (Scheme 12).⁵⁴

6. Synthesis of Tagged Peptides and Peptidomimetics

6.1. Fluorescent Labels


Fluorescent labeling of biological systems is of great importance. There is also increasing interest in establishing methods for the incorporation of non-natural amino acids into proteins without suppression of binding ability.⁵⁵

6.1.1. Coumarin-Labeled Peptides and Peptidomimetics

Coumarins are constituents of many commercially important fluorescent dyes since they offer high-emission quantum yields, photostability, and good solubility in most solvents. Synthetic methods for the incorporation of coumarin into amino acids and peptides are now available. (Coumarin-3-ylcarbonyl)benzotriazoles react readily with a variety of aminoxy acids in aqueous acetonitrile at room temperature to form coumarin-labeled aminoxy acids.⁵⁶ Coumarin-labeled aminoxy acids have been prepared in a similar way⁵⁷ and coumarin-labeled aminoxy hybrid peptides have been obtained in two steps from the respective benzotriazoles (Scheme 13).⁵⁶ The 7-methoxycoumarin derivatives have quantum yields of 0.35–0.71 and may therefore be useful in peptide assays.

6.1.2. 6-Chloro-2,3-naphthalimides and Water-Soluble Fluorescent Tags

Organic fluorophores that contain a naphthalene nucleus are of interest since, on binding with a substrate, they often exhibit significant changes



Scheme 12. Preparation of Cyclic Peptidomimetics. (Ref. 54b)

in their fluorescence spectra, quantum yields, and lifetimes in different solvents. The benzotriazole methodology offers access to new 6-chloro-2,3-naphthalimide derivatives (Scheme 14).⁵⁸

6.2. Azo-Dye-Labeled Peptides

Azo-arene carboxylic acids are widely used as molecular switches in life sciences.⁵⁹⁻⁶¹ The coupling of an azo-dye carboxylic acid to a biological moiety has, in many cases, required harsh conditions and given poor yields.^{62,63} In contrast, *N*-(4-arylazobenzoyl)benzotriazole and glycine in DMF-water at 20 °C give the coupled product in 99% yield. Similarly, other amino acids undergo this facile coupling (Scheme 15).⁶⁴

7. Peptide Conjugates


Conjugates comprise peptides attached to another molecular skeleton, usually through either the carboxyl group (C-terminus) or the amino group (N-terminus) of an amino acid.

7.1. Conjugates of Sugars


Numerous naturally occurring carbohydrate conjugates link a sugar glycosidically to an α -amino acid unit of a peptide or protein. Considerable effort has been devoted to the synthesis of N- and O-linked glycopeptide conjugates utilizing both solution- and solid-phase methodologies.⁶⁵⁻⁶⁹ The benzotriazole methodology is advantageous for the solution-phase synthesis of chirally pure *O*-(α -aminoacyl)-⁷⁰ and *N*-(α -aminoacyl)sugar conjugates.⁷¹ A typical procedure utilizes DMAP catalysis and microwave irradiation to give 82–92% yields of O- or N-coupled products (Scheme 16, Part (a)).⁷¹ The same benzotriazole-based method also provides a convenient and efficient route to Cbz-protected tri- and tetrapeptide conjugates with sugars (Scheme 16, Part (b)).⁷²

7.2. Conjugates of Heterocycles


(α -Aminoacyl)amino-substituted heterocycles are of considerable importance as synthetic intermediates (e.g., for endomorphin-2 (EM-2) analogues),⁷³ and because of their diverse biological activity (e.g., as inhibitors of bacterial RND efflux pump⁷⁴⁻⁷⁶ and of tumor necrosis factor- α converting enzyme (TACE) GW 333377). Until recently, only a few reports had described the preparation of α -aminoacyl conjugates




Scheme 13. Preparation of Coumarin-Labeled Amino Hybrid Peptides (Ref. 56)



Scheme 14. Preparation of 6-Chloro-2,3-naphthalimide and Water-Soluble Fluorescent Taqs. (Ref. 58)



Scheme 15. Preparation of Azo-Dye-Labeled Amino Acids. (Ref. 64)




Scheme 16. Examples of (a) (α -Aminoacyl)sugar Conjugate and (b) Tetrapeptide Sugar Conjugate Synthesis. (Ref. 71)

of heterocycles by C- or N-acylation of heterocycles with amino acids. Kraus and co-workers⁷⁸ reported that "N-acylation of weakly nucleophilic heterocyclic amines by protected amino acid is not a straightforward reaction which could be achieved under any standard coupling conditions".

N-Aminoacyl benzotriazolides now enable the synthesis of chirally pure α -aminoacyl conjugates of heterocycles even from weakly nucleophilic heterocyclic amines by N-acylation in DMF under microwave irradiation.⁷⁹ The analogous C-acylation of lithiated methylpyridine or methylquinolone in THF at -20 °C for 1–3 h with *N*-aminoacyl benzotriazolides gave the corresponding α -aminoacyl C-linked conjugates (Scheme 17).⁸⁰ The convenient and efficient formation of Cbz-protected tri- (e.g., Z-L-Val-L-Phe-Gly-NH-(2-Pyr)) and tetrapeptide (e.g., Z-L-Phe-Gly-L-Leu-Gly-NH-(*N*-methylpiperazine)) conjugates with heterocyclic nuclei of biological importance succeeds under a variety of reaction conditions.^{72,80}

7.3. Conjugates of Vitamins

Recent approaches to enhance vitamin uptake include covalently bonding the vitamins to peptides. Water-soluble vitamins are usually



R = H, Me, *i*-Pr, Bn, *t*-Bu, *s*-Bu, MeS(CH₂)₂, (indol-3-yl)CH₂

Het = 2-thiazolyl, 2-(6-methoxybenzothiazolyl), 2-(1-benzylbenzimidazolyl), 5-(3-methoxy-1,2,4-thiadiazolyl), 2-pyridyl, 2-(4,6-dimethylpyridyl), 2-(4-methylpyridyl), 4-pyridyl, 2-quinolyl

Scheme 17. α -Aminoacyl Conjugates of Amino-Substituted Heterocycles (Ref. 79,80)




Figure 1. Peptide Conjugates of Vitamins Prepared by the Benzotriazole Methodology. (Ref. 84)



transported into cells by potocytosis.⁸¹ Zhang and McCormick⁸² have proposed the delivery of vitamin B6 by receptor-mediated transport in eukaryotic cells with the amine of a peptide–vitamin conjugate, which facilitates the cell uptake of peptide and transport into the cytosol.^{82,83} The benzotriazole methodology enables the efficient coupling of amino acids and peptides to vitamins, again utilizing microwave irradiation to shorten reaction times and minimize epimerization.⁸⁴

Niacin and biotin benzotriazolides couple with free amino acids, dipeptides, and tripeptides (NEt₃, MeCN–H₂O, μ w, 70 °C) to give the corresponding bioconjugates **1** and **2** in yields of 43–81% and 35–82%, respectively (Figure 1).⁸⁴ Amino acid and peptide conjugates of vitamin D3 (Figure 1, **3**) are obtained by O-acylation of cholecalciferol with Cbz-protected acylbenzotriazoles in the presence of DMAP in THF and under microwave irradiation (50 W, 70 °C) for 1–2 h.⁸⁴ Amino acid and peptide conjugates of α -tocopherol (Figure 1, **4**) are formed by O-acylation of α -tocopherol with Cbz-protected acylbenzotriazoles under microwave irradiation (20 W, 50 °C, 0.3 h) in anhydrous DMF in the presence of potassium carbonate.⁸⁴

7.4. Conjugates of Pharmaceuticals

Improving the efficacy of therapeutics, particularly through enhanced local delivery to a diseased cell, is an important topic in pharmaceutical R&D. These techniques are helping to solve traditional drug delivery challenges, such as poor cellular uptake and/or non-specific toxicity.

The utilization of prodrugs that temporarily mask the acidic groups of NSAIDs may increase uptake and reduce irritation caused by direct contact.^{85,86} Indomethacin, diclofenac,⁸⁷ ibuprofen, and naproxen,^{87–91} are among well-known NSAIDs that have been modified by linking to natural amino acids, as reported by numerous investigators. We have recently shown that ibuprofen and naproxen bioconjugates are readily



Scheme 18. Preparation of Peptide Conjugates of Pharmaceuticals. (Ref. 92,100)

prepared by reacting NSAID benzotriazolides with amino acids and dipeptides in aqueous acetonitrile–triethylamine at 20 °C (Scheme 18, Part (a)).⁹²

Prodrugs formed by linking quinolone acids with amino acid esters are more lipophilic than the parent drugs,^{93,94} and show enhanced in vivo antibacterial properties^{95,96} with pronounced therapeutic effects against *Pseudomonas aeruginosa*,⁹⁷ *Escherichia coli*,⁹⁸ *Staphylococcus aureus*,⁹⁸ and *Salmonella typhi*.⁹⁶ Additional wide-ranging biological activities include anti-allergic,⁹⁹ antihypertensive,⁹⁹ bronchodilation,⁹⁹ and binding to bovine serum albumin.⁹² Amino acid conjugates of quinolone antibiotics were prepared in 39–88% yields by coupling free amino acids with benzotriazole-activated oxolinic acid, nalidixic acid, cinoxacin, or flumequine (Scheme 18, Part (b)).^{100,101}

7.5. Conjugates of Plant Hormones

Indole-3-acetic acid (IAA), an indispensable plant hormone (auxin), also occurs naturally as “conjugates” linked to amino acids.¹⁰² Gene expression and cell division, elongation, and differentiation in plant tissue are all regulated by these indole-3-acetic acid auxins. Another endogenous plant hormone, indole-3-propionic acid (IPA), and its amino acid conjugates interact with serum albumin.¹⁰³ Plant hormone benzotriazolides—prepared in 86–90% yields by standard treatment of indole-3-acetic acid and indole-3-propionic acid—coupled with diverse amino acids to give conjugates in 40–70% yields (Scheme 19).¹⁰⁴

7.6. Aminoxy Acid Conjugates of Peptidomimetics and Hydrazino Acid Conjugates

The benzotriazole methodology enables convenient and efficient synthesis of novel aminoxy acid containing conjugates even at hindered nucleophilic centers in steroids, terpenes, sugars, and nucleosides.¹⁰⁵ The benzotriazolides of α -hydrazino acids were used to generate hydrazine acid conjugates through N-, O-, S-, and C-acylations in good yields (49–88%).⁴⁹


8. Differential N-, O-, and S-Acylation

Isopeptides are used for the detection and capture of ubiquitinating and de-ubiquitinating enzymes using activity-based protein profiling (ABPP).¹⁰⁶ The presence of an additional amino group in N-, O-, or S-acyl isopeptides generally increases their hydrophilicity, which is advantageous in effecting their purification by HPLC. The native peptides can then be generated from the corresponding N-, O-, or S-acyl isopeptides via an N to N,¹⁰⁷ O to N,¹⁰⁸ or S to N^{6,109,110} intramolecular acyl migration reaction. These findings have led to the synthesis of peptides containing difficult sequences.¹¹¹


8.1. S- and O-Acyl Isopeptides

S-Acyl isopeptides are usually less likely to aggregate in solution and therefore are easier to synthesize and purify relative to the corresponding native peptides. S-Acylation of protected cysteine-containing peptides was carried out in the presence of KHCO₃ at 20 °C in acetonitrile (Schemes 20, Part (a)).^{6,109,110} Selective S-acylation of cysteine was also carried out in acetonitrile–water mixture in the absence of base.


O-Acylation of protected serine and protected threonine with various N-Pg-(α -aminoacyl)benzotriazolides in the presence of diisopropylethylamine in acetonitrile at 20 °C for 12 h gave O-acylserine and O-acylthreonine dipeptides without racemization (Scheme 20, Part (b)). O-Acylisotyrosine tripeptides were also prepared in yields of 74–91% by reacting tyrosine-containing protected dipeptides with N-Pg-(α -aminoacyl)benzotriazole in the presence of DBU in acetonitrile at 20 °C for 12 h.²⁴



Scheme 19. Preparation of Conjugates of Plant Hormones. (Ref. 104)



Scheme 20. Synthesis of S- and O-acyl Isopeptides. (Ref. 6,9,24,110)



8.2. Ligation at a Distance

8.2.1. S to N Acyl Migration


The S to N acyl migration through various cyclic transition states was investigated by carrying out the ligation experiment on mono-isopeptides under microwave irradiation (50 W) at 50 °C for 1–3 h using 1 M NaH_2PO_4 – Na_2HPO_4 phosphate buffer to maintain pH 7.3 (eq 4). The rates and yields of long-range S to N acyl transfers were found to depend significantly on the size of the macrocyclic transition state (TS), with the rates qualitatively following the TS ring-size trend $5 > 10 > 11 > 14, 16, 17 > 12 > 13, 15, 19 > 18 >>> 9 > 8$.^{6,109,110}

8.2.2. O to N Acyl Migration


The chemical ligation of serine isopeptide through O to N acyl transfer via 8- and 11-membered-ring transition states occurs without the use of an auxiliary group (eq 5).¹⁰⁸ In contrast, threonine isopeptide failed to undergo acyl migration even under more basic conditions and longer reaction times.¹¹² Chemical ligation studies of tyrosine isopeptides (μw , 50 W, 50 °C, 3 h, using 1 M phosphate buffer and DMF–piperidine) via 12- to 19-membered-ring cyclic transition states showed that intramolecular O to N acyl transfer occurs with 12- to 14-membered-ring TS's under basic conditions and with 15- to 19-membered-ring TS's in aqueous media.¹⁰⁸

8.2.3. N to N Acyl Migration

Tryptophan isopeptides with α -, β -, or γ -amino acid units were synthesized, and the acyl migration from the indole nitrogen to the terminal NH_2 was studied under microwave irradiation. Intramolecular



eq 5 (Ref. 108)



eq 6 (Ref. 107)

acyl transfer through 10-, 11-, and 12-membered-ring transition states was favored over that through a 7-membered-ring TS, and acyl migration occurred more readily in basic, nonaqueous media relative to aqueous buffered conditions (eq 6).¹⁰⁷

9. Conclusions and Comparison with Alternative Methodologies

9.1. Carboxyl Group Activation by Isolation of an Intermediate

The most obvious method for activating the carboxyl group of an amino acid for peptide bond formation at room temperature or below is by forming the corresponding acid chloride.¹¹³ This type of activation has been carried out with chlorinating reagents such as pivaloyl chloride,¹¹⁴ phthaloyl dichloride,¹¹⁵ thionyl chloride,¹¹⁶ and oxalyl chloride.¹¹⁷ However, an amino acid chloride bearing an acid labile protecting group can easily racemize through the oxazolone, which limits the application of acid chlorides despite their high reactivity and low cost. Amino acid fluorides are less moisture-sensitive than acyl chlorides, but the fluorinating reagents are expensive and hazardous, and the peptide-forming reactions require purification by chromatography.¹¹⁸

The acyl azide method of peptide coupling was developed about 100 years ago. It is not attractive for routine use because it involves four distinct steps, including two stable intermediates that require purification.¹¹⁹ An additional side reaction that occurs at higher temperature is rearrangement of the acyl azide into the alkyl isocyanate, which can react with nucleophiles to yield a peptide urea that is difficult to remove from the product.¹²⁰ Recently, El-Faham and Albericio published a review on the use of different peptide coupling reagents including benzotriazoles.¹²¹

9.2. Carboxyl Group Activation without Isolation of an Intermediate

Besides acyl halides and acyl azides, other methods for peptide coupling include the use of various reagents, where the intermediates are not isolated. A traditional approach to form peptide bonds is the carbodiimide method, using dicyclohexylcarbodiimide (DCC). However, despite being compatible with solid-phase synthesis (SPS) that uses *tert*-butoxycarbonyl (Boc) chemistry, DCC is not compatible with the fluorenylmethoxycarbonyl (Fmoc) group. When DCC is utilized in solution, traces of the byproduct, DCU, are difficult to remove, even after passage through a chromatography column. Thus, DCC has been replaced by reagents such as diisopropylcarbodiimide (DIC), *N*-ethyl-*N'*-(3-dimethylaminopropyl)carbodiimide (EDC), and *N*-cyclohexyl-*N'*-isopropylcarbodiimide (CIC), all of which are relatively soluble in DCM and therefore more suitable for Fmoc-SPS. Additives such as 1-hydroxybenzotriazole (HOBT), 1-hydroxy-7-azabenzotriazole (HOAt), *N*-hydroxysuccinimide (HOSu), and others increase the efficiency of carbodiimide-mediated reactions and decrease the degree of racemization.¹²²

Phosphonium reagents were developed to avoid racemization and side reactions that can occur with carbodiimide reagents. Coste et al. introduced chloro- and bromotris(dimethylamino)phosphonium hexafluorophosphate (Cl₃P and Br₃P) as peptide-coupling reagents with noticeable racemization in the Young test.¹²³ HOBt may be used in combination with (benzotriazol-1-yloxy)tris(dimethylamino)phosphonium hexafluorophosphate (BOP) to suppress racemization. However, the intermediates formed in the coupling reaction are highly unstable and BOP is reported to be highly carcinogenic.¹²⁴

The search for better coupling reagents based on DCC led to carbonyl diimidazole (CDI).¹²⁵ Rapoport introduced the imidazolium

reagent 1,10-carbonylbis(3-methylimidazolium) triflate (CBMIT) by bis-methylating CDI with methyl triflate.¹²⁶ This reagent showed no sign of racemizing the amino acid residues in the presence of CuCl₂ or Cu(OTf)₂. However, CBMIT is moisture-sensitive and, due to its polarity, the method is restricted to polar solvents such as nitromethane.¹²⁷ The reactivity of these reagents also increases in the presence of additives like HOAt, HOBt, or DMAP.

Aminium or uronium reagents such as *N*-[(1*H*-benzotriazol-1-yl)(dimethylamino)methylene]-*N*-methylmethanaminium hexafluorophosphate *N*-oxide (HBTU) and *N*-[(dimethylamino)-1*H*-1,2,3-triazolo[4,5-*b*]pyridin-1-ylmethylene]-*N*-methylmethanaminium hexafluorophosphate *N*-oxide (HATU) all react directly with the amine moiety of the amino acid residue to give a guanidine side product, which terminates the peptide chain.¹²⁸

Benzotriazole offers an extremely useful alternative to all the above methods by affording a versatile range of coupling procedures under the mild conditions required to avoid racemization.

10. References

- (1) Staab, H. A.; Wendel, K. *Angew. Chem.* **1960**, *72*, 708
- (2) Staab, H. A.; Wendel, K. *Chem. Ber.* (presently *Eur. J. Inorg. Chem.*) **1960**, *93*, 2902.
- (3) (a) Katritzky, A. R.; Rachwal, S. *Chem. Rev.* **2010**, *110*, 1564. (b) Katritzky, A. R.; Rachwal, S. *Chem. Rev.* **2011**, *111*, 7063
- (4) Katritzky, A. R.; Xiangfu, L.; Yang, J. Z.; Denisko, O. V. *Chem. Rev.* **1998**, *98*, 409.
- (5) (a) Katritzky, A. R.; Angrish, P.; Todadze, E. *Synlett* **2009**, 2392. (b) Katritzky, A. R.; Singh, A.; Haase, D. N.; Yoshioka, M. *ARKIVOC* **2009**, viii, 47.
- (6) Panda, S. S.; El-Nachef, C.; Bajaj, K.; Al-Youbi, A. O.; Oliferenko, A.; Katritzky, A. R. *Chem. Bio. Drug Des.* **2012**, *80*, 821.
- (7) Katritzky, A. R.; Wang, M.; Yang, H.; Zhang, S.; Akhmedov, N. G. *ARKIVOC* **2002**, viii, 134.
- (8) Katritzky, A. R.; Vincek, A. S.; Suzuki, K. *ARKIVOC* **2005**, v, 116.
- (9) Katritzky, A. R.; Angrish, P.; Suzuki, K. *Synthesis* **2006**, 411.
- (10) Panda, S. S.; Ibrahim, M. A.; Küçükbay, H.; Meyers, M. J.; Sverdrup, F. M.; El-Feky, S. A.; Katritzky, A. R. *Chem. Bio. Drug Des.* **2013**, *81*, DOI: 10.1111/cbdd.12134.
- (11) Abdelmajeid, A.; Tala, S. R.; Amine, M. S.; Katritzky, A. R. *Synthesis* **2011**, 2995.
- (12) Katritzky, A. R.; Suzuki, K.; Singh, S. K. *Synthesis* **2004**, 2645.
- (13) Katritzky, A. R.; Meher, G.; Angrish, P. *Chem. Biol. Drug Des.* **2006**, *68*, 326.
- (14) Katritzky, A. R.; Todadze, E.; Shestopalov, A. A.; Cusido, J.; Angrish, P. *Chem. Biol. Drug Des.* **2006**, *68*, 42.
- (15) Katritzky, A. R.; Todadze, E.; Cusido, J.; Angrish, P.; Shestopalov, A. A. *Chem. Biol. Drug Des.* **2006**, *68*, 37.
- (16) Katritzky, A. R.; Meher, G.; Narindoshvili, T. *J. Org. Chem.* **2008**, *73*, 7153.
- (17) Katritzky, A. R.; Todadze, E.; Angrish, P.; Draghici, B. *J. Org. Chem.* **2007**, *72*, 5794.
- (18) Katritzky, A. R.; Yoshioka, M.; Narindoshvili, T.; Chung, A.; Khashab, N. M. *Chem. Biol. Drug Des.* **2008**, *72*, 182.
- (19) Ibrahim, T. S.; Tala, S. R.; El-Feky, S. A.; Abdel-Samii, Z. K.; Katritzky, A. R. *Chem. Biol. Drug Des.* **2012**, *80*, 194.
- (20) Katritzky, A. R.; Khashab, N. M.; Yoshioka, M.; Haase, D. N.; Wilson, K. R.; Johnson, J. V.; Chung, A.; Haskell-Luevano, C. *Chem. Biol. Drug Des.* **2007**, *70*, 465.
- (21) Katritzky, A. R.; Haase, D. N.; Johnson, J. V.; Chung, A. *J. Org. Chem.* **2009**, *74*, 2028.
- (22) Hjørringgaard, C. U.; Pedersen, J. M.; Vosegaard, T.; Nielson, N. C.; Skrydstrup, T. *J. Org. Chem.* **2009**, *74*, 1329.
- (23) Haddoub, R.; Dauner, M.; Stefanowicz, F. A.; Barattini, V.; Laurent, N.; Flitsch, S. L. *Org. Biomol. Chem.* **2009**, *7*, 665.
- (24) El Khatib, M.; Jauregui, L.; Tala, S. R.; Khelashvili, L.; Katritzky, A. R. *Med. Chem. Commun.* **2011**, *2*, 1087.
- (25) Wenschuh, H.; Beyermann, M.; Haber, H.; Seydel, J. K.; Krause, E.; Bienert, M.; Carpino, L. A.; El-Faham, A.; Albericio, F. *J. Org. Chem.* **1995**, *60*, 405.
- (26) Wenger, R. M. *Angew. Chem.* **1985**, *97*, 88.
- (27) Abe, H.; Sato, K.; Kato, T.; Izumiya, M. *Bull. Chem. Soc. Jpn.* **1976**, *49*, 3113.
- (28) Aoki, S.; Cao, L.; Matsui, K.; Rachmat, R.; Akiyama, S.; Kobayashi, M. *Tetrahedron* **2004**, *60*, 7053.
- (29) Printsevskaya, S. S.; Solovieva, S. E.; Olsufyeva, E. N.; Mirchink, E. P.; Isakova, E. B.; De Clerq, E.; Balzarini, J.; Preobrazhenskaya, M. N. *J. Med. Chem.* **2005**, *48*, 3885.
- (30) Beagle, L. K.; Hansen, F. K.; Monbaliu, J.-C. M.; DesRosiers, M. P.; Phillips, A. M.; Stevens, C. V.; Katritzky, A. R. *Synlett* **2012**, 2337.
- (31) Ha, K.; Monbaliu, J. C. M.; Williams, B. C.; Pillai, G. G.; Ocampo, C. E.; Zeller, M.; Stevens, C. V.; Katritzky, A. R. *Org. Biomol. Chem.* **2012**, *10*, 8055.
- (32) Lee, M.; Lee, J.; Baek, B.; Shin, I. *Synlett* **2003**, 325.
- (33) Duléry, V.; Renaudet, O.; Dumy, P. *Tetrahedron* **2007**, *63*, 11952
- (34) Yang, D.; Zhang, Y.-H.; Li, B.; Zhang, D.-W. *J. Org. Chem.* **2004**, *69*, 7577.
- (35) Zhang, Y.-H.; Song, K.; Zhu, N.-Y.; Yang, D. *Chem. Eur. J.* **2010**, *16*, 577.
- (36) Briggs, M. T.; Morley, J. S. *J. Chem. Soc., Perkin Trans. 1* **1979**, 2138.
- (37) Brückner, A. M.; Chakraborty, P.; Gellman, S. H.; Diederichsen, U. *Angew. Chem., Int. Ed.* **2003**, *42*, 4395.
- (38) Katritzky, A. R.; Avan, I.; Tala, S. R. *J. Org. Chem.* **2009**, *74*, 8690.
- (39) Chen, J.; Zhang, M.; Ju, W.; Waldmann, T. A. *Blood* **2009**, *113*, 1287.
- (40) Kuisle, O.; Quiñoá, E.; Riguera, R. *J. Org. Chem.* **1999**, *64*, 8063.
- (41) Avan, I.; Tala, S. R.; Steel, P. J.; Katritzky, A. R. *J. Org. Chem.* **2011**, *76*, 4884.
- (42) Proulx, C.; Sabatino, D.; Hopewell, R.; Spiegel, J.; Garcia Ramos, Y.; Lubell, W. D. *Future Med. Chem.* **2011**, *3*, 1139.
- (43) Otto, H.-H.; Schirmeister, T. *Chem. Rev.* **1997**, *97*, 133.
- (44) Martin, N. I.; Beeson, W. T.; Woodward, J. J.; Marletta, M. A. *J. Med. Chem.* **2008**, *51*, 924.
- (45) Abo-Dya, N. E.; Biswas, S.; Basak, A.; Avan, I.; Alamry, K. A.; Katritzky, A. R. *J. Org. Chem.* **2013**, *78*, 3541.
- (46) Stigers, K. D.; Soth, M. J.; Nowick, J. S. *Curr. Opin. Chem. Biol.* **1999**, *3*, 714.
- (47) Günther, R.; Hofmann, H.-J. *J. Am. Chem. Soc.* **2001**, *123*, 247.
- (48) Ronwin, E.; Roach, J. V.; Tucker, B. L. *Can. J. Chem.* **1968**, *46*, 3013.
- (49) Panda, S. S.; El-Nachef, C.; Bajaj, K.; Katritzky, A. R. *Eur. J. Org. Chem.* **2013**, 4156.
- (50) Boeglin, D.; Cantel, S.; Heitz, A.; Martinez, J.; Fehrentz, J.-A. *Org. Lett.* **2003**, *5*, 4465
- (51) Panda, S. S.; Malik, R.; Chand, M.; Jain, S. C. *Med. Chem. Res.* **2012**, *21*, 3750.
- (52) Hansen, F. K.; Beagle, L. K.; Todadze, E.; Katritzky, A. R. *Heterocycles* **2012**, *84*, 515.
- (53) Katritzky, A. R.; Tala, S. R.; Abo-Dya, N. E.; Gyanda, K.; El-Gendy, B. E.-D. M.; Abdel-Samii, Z. K.; Steel, P. J. *J. Org. Chem.* **2009**, *74*, 7165
- (54) (a) Berhanu, W. M.; Ibrahim, M. A.; Pillai, G. G.; Oliferenko, A. A.; Khelashvili, L.; Jabeen, F.; Mirza, B.; Ansari, F. L.; Ul-Haq, I.; El-

- Feky, S. A.; Katritzky, A. R. *Beil. J. Org. Chem.* **2012**, *8*, 1146. (b) Ibrahim, M. M. A.; Panda, S. S.; Nhon, L.; Hamed, A.; El-Feky, S. A.; Katritzky, A. R. *Synthesis* **2013**, 767.
- (55) Katritzky, A. R.; Narindoshvili, T. *Org. Biomol. Chem.* **2009**, *7*, 627.
- (56) Katritzky, A. R.; Abdelmajeid, A.; Tala, S. R.; Amine, M. S.; Steel, P. J. *Synthesis* **2011**, 83.
- (57) Katritzky, A. R.; Narindoshvili, T.; Angrish, P. *Synthesis* **2008**, 2013.
- (58) Katritzky, A. R.; Ozcan, S.; Todadze, E. *Org. Biomol. Chem.* **2010**, *8*, 1296.
- (59) Natansohn, A.; Rochon, P. *Chem. Rev.* **2002**, *102*, 4139.
- (60) Pieroni, O.; Fissi, A.; Angelini, N.; Lenci, F. *Acc. Chem. Res.* **2001**, *34*, 9.
- (61) Renner, C.; Moroder, L. *ChemBioChem* **2006**, *7*, 868.
- (62) Rose, T. M.; Prestwich, G. D. *Org. Lett.* **2006**, *8*, 2575.
- (63) Grossmann, T. N.; Seitz, O. *J. Am. Chem. Soc.* **2006**, *128*, 15596.
- (64) Katritzky, A. R.; Chen, Q.-Y.; Tala, S. R. *Org. Biomol. Chem.* **2008**, *6*, 2400.
- (65) Geurtsen, R.; Boons, G.-J. *Eur. J. Org. Chem.* **2002**, 1473.
- (66) Horvat, S.; Varga-Defterdarović, L.; Roščić, M.; Horvat, J. *Chem. Commun.* **1998**, 1663.
- (67) Burt, J.; Dean, T.; Warriner, S. *Chem. Commun.* **2004**, 454.
- (68) Tennant-Eyles, R. J.; Davis, B. G.; Fairbanks, A. J. *Tetrahedron: Asymmetry* **2003**, *14*, 1201.
- (69) Park, O.-J.; Jeon, G.-J.; Yang, J.-W. *Enzyme Microb. Technol.* **1999**, *25*, 455.
- (70) Katritzky, A. R.; Angrish, P.; Narindoshvili, T. *Bioconjugate Chem* **2007**, *18*, 994.
- (71) Katritzky, A. R.; Narindoshvili, T.; Draghici, B.; Angrish, P. *J. Org. Chem.* **2008**, *73*, 511.
- (72) Bajaj, K.; Panda, S. S.; El-Nachef, C.; Katritzky, A. R. *Chem. Biol. Drug Des.* **2012**, *80*, 17.
- (73) Fujita, Y.; Tsuda, Y.; Li, T.; Motoyama, T.; Takahashi, M.; Shimizu, Y.; Yokoi, T.; Sasaki, Y.; Ambo, A.; Kita, A.; Jinsmaa, Y.; Bryant, S. D.; Lazarus, L. H.; Okada, Y. *J. Med. Chem.* **2004**, *47*, 3591.
- (74) Watkins, W. J.; Landaverry, Y.; Léger, R.; Litman, R.; Renau, T. E.; Williams, N.; Yen, R.; Zhang, J. Z.; Chamberland, S.; Madsen, D.; Griffith, D.; Tembe, V.; Huie, K.; Dudley, M. N. *Bioorg. Med. Chem. Lett.* **2003**, *13*, 4241.
- (75) Renau, T. E.; Léger, R.; Filonova, L.; Flamme, E. M.; Wang, M.; Yen, R.; Madsen, D.; Griffith, D.; Chamberland, S.; Dudley, M. N.; Lee, V. J.; Lomovskaya, O.; Watkins, W. J.; Ohta, T.; Nakayama, K.; Ishida, Y. *Bioorg. Med. Chem. Lett.* **2003**, *13*, 2755.
- (76) Fuwa, H.; Okamura, Y.; Morohashi, Y.; Tomita, T.; Iwatsubo, T.; Kan, T.; Fukuyama, T.; Natsugari, H. *Tetrahedron Lett.* **2004**, *45*, 2323.
- (77) Rabinowitz, M. H.; Andrews, R. C.; Becherer, J. D.; Bickett, D. M.; Bubacz, D. G.; Conway, J. G.; Cowan, D. J.; Gaul, M.; Glennon, K.; Lambert, M. H.; Leesnitzer, M. A.; McDougald, D. L.; Moss, M. L.; Musso, D. L.; Rizzolio, M. C. *J. Med. Chem.* **2001**, *44*, 4252.
- (78) Quéléver, G.; Burlet, S.; Garino, C.; Pietrancosta, N.; Laras, Y.; Kraus, J.-L. *J. Comb. Chem.* **2004**, *6*, 695.
- (79) Katritzky, A. R.; El-Gendy, B. E.-D. M.; Todadze, E.; Abdel-Fattah, A. A. *J. Org. Chem.* **2008**, *73*, 5442.
- (80) Katritzky, A. R.; Bajaj, K.; Charpentier, M.; Zadeh, E. G. *J. Org. Chem.* **2010**, *75*, 3938.
- (81) Rothberg, K. G.; Ying, Y.; Kolhouse, J. F.; Kamen, B. A.; Anderson, R. G. W. *J. Cell Biol.* **1990**, *110*, 637.
- (82) Zhang, Z.; McCormick, D. B. *Proc. Natl. Acad. Sci. U.S.A.* **1991**, *88*, 10407.
- (83) Zhu, T.; Stein, S. *Bioconjugate Chem.* **1994**, *5*, 312.
- (84) El-Nachef, C.; Bajaj, K.; Koblick, J.; Katritzky, A. R. *Eur. J. Org. Chem.* **2012**, 4412.
- (85) Nielsen, N. M.; Bundgaard, H. *J. Pharm. Sci.* **1988**, *77*, 285.
- (86) Bundgaard, H.; Nielsen, N. M. *J. Med. Chem.* **1987**, *30*, 451.
- (87) Otis, M.-F.; Levesque, L.; Marceau, F.; Lacroix, J.; Gaudreault, R. C. *Inflammopharmacology* **1992**, *1*, 201.
- (88) Lane, M. E.; Levis, K.; McDonald, G. S. A.; Corrigan, O. I. *Int. J. Pharm.* **2006**, *312*, 196.
- (89) Kwapiszewski, W.; Krupinska, J.; Koziej, P.; Stephanus, R. *Acta Pol. Pharm.* **1985**, *42*, 545.
- (90) Shanbhag, V. R.; Crider, A. M.; Gokhale, R.; Harpalani, A.; Dick, R. M. *J. Pharm. Sci.* **1992**, *81*, 149.
- (91) Levit, G. L.; Anikina, L. V.; Vikharev, Y. B.; Demin, A. M.; Safin, V. A.; Matveeva, T. V.; Krasnov, V. P. *Pharm. Chem. J.* **2002**, *36*, 232.
- (92) Katritzky, A. R.; Jishkariani, D.; Narindoshvili, T. *Chem. Biol. Drug Des.* **2009**, *73*, 618.
- (93) Aboul-Fadl, T.; Fouad, E. A. *Pharmazie* **1996**, *51*, 30.
- (94) Kohli, D. V.; Uppadhyay, R. K.; Saraf, S. K.; Vishwakarma, K. K. *Pharmazie* **1992**, *47*, 57.
- (95) Sanchez, J. P. Eur. Patent EP 304 087 A2, 1989; *Chem. Abstr.* **1989**, *111*, 153777.
- (96) Leonova, T. S.; Padeiskaya, E. N.; Yashunskii, V. G. *Khim.-Farm. Zh.* **1987**, *21*, 692.
- (97) Pongracz, K.; Frank, J.; Kulcsar, G.; Medzihradzsky, K.; Nadasi, L. Hung. Patent 36483, 1985; *Chem. Abstr.* **1987**, *106*, 120225.
- (98) Jefferson, E. A.; Swayze, E. E.; Osgood, S. A.; Miyaji, A.; Risen, L. M.; Blyn, L. B. *Bioorg. Med. Chem. Lett.* **2003**, *13*, 1635.
- (99) Laruelle, C.; Lepant, M.; Raynier, B. Fr. Patent 2 564 832, 1985; *Chem. Abstr.* **1986**, *105*, 153540.
- (100) Katritzky, A. R.; Munawar, M. A.; Kovacs, J.; Khelashvili, L. *Org. Biomol. Chem.* **2009**, *7*, 2359.
- (101) Panda, S. S.; Bajaj, K.; Meyers, M. J.; Sverdrup, F. M.; Katritzky, A. R. *Org. Biomol. Chem.* **2012**, *10*, 8985.
- (102) *Biochemistry and Molecular Biology of Plant Hormones*; Hooykaas, P. J. J., Hall, M. A., Libbenga, K. R., Eds.; Elsevier: Amsterdam, 1999; Vol. 115.
- (103) Tomašić, A.; Bertoša, B.; Tomić, S.; Šoškić, M.; Magnus, V. *J. Chromatogr. A* **2007**, *1154*, 240.
- (104) Katritzky, A. R.; Khelashvili, L.; Munawar, M. A. *J. Org. Chem.* **2008**, *73*, 9171.
- (105) Draghici, B.; Hansen, F. K.; Buciumas, A.-M.; El-Gendy, B. E.-D. M.; Todadze, E.; Katritzky, A. R. *RSC Advances* **2011**, *1*, 602.
- (106) Sheppard, R. J. *Peptide Sci.* **2003**, *9*, 545.
- (107) Popov, V.; Panda, S. S.; Katritzky, A. R. *Org. Biomol. Chem.* **2013**, *11*, 1594.
- (108) Popov, V.; Panda, S. S.; Katritzky, A. R. *J. Org. Chem.* **2013**, *78*, DOI: 10.1021/jo4009468.
- (109) Ha, K.; Chahar, M.; Monbaliu, J.-C. M.; Todadze, E.; Hansen, F. K.; Olfierenko, A. A.; Ocampo, C. E.; Leino, D.; Lillioctch, A.; Stevens, C. V.; Katritzky, A. A. *J. Org. Chem.* **2012**, *77*, 2637.
- (110) Katritzky, A. R.; Tala, S. R.; Abo-Dya, N. E.; Ibrahim, T. S.; El-Feky, S. A.; Gyanda, K.; Pandya, K. M. *J. Org. Chem.* **2011**, *76*, 85.
- (111) Monbaliu, J.-C. M.; Katritzky, A. R. *Chem. Commun.* **2012**, *48*, 11601.
- (112) El Khatib, M.; Elagawany, M.; Jabeen, F.; Todadze, E.; Bol'shakov, O.; Olfierenko, A.; Khelashvili, L.; El-Feky, S. A.; Asiri, A.; Katritzky, A. R. *Org. Biomol. Chem.* **2012**, *10*, 4836.
- (113) *The Peptides, Methods of Peptide Synthesis*; Schroder, E., Lubke, K., Eds.; Academic Press: New York, 1965; p 77.
- (114) Weisz, I.; Roboz, J.; Bekesi, J. G. *Tetrahedron Lett.* **1996**, *37*, 563.
- (115) Rozov, L. A.; Rafalko, P. W.; Evans, S. M.; Brockunier, L.; Ramig, K. *J. Org. Chem.* **1995**, *60*, 1319.
- (116) Zhang, L.; Chung, J. C.; Costello, T. D.; Valvis, I.; Ma, P.; Kauffman,

- S.; Ward, R. *J. Org. Chem.* **1997**, *62*, 2466.
- (117) Senokuchi, K.; Nakai, H.; Nagao, Y.; Sakai, Y.; Katsube, N.; Kawamura, M. *Bioorg. Med. Chem.* **1998**, *6*, 441.
- (118) Hudlicky, M. *Org. React.* **1988**, *35*, 513.
- (119) Meienhofer, J. The Azide Method in Peptide Synthesis. In *The Peptides: Analysis, Synthesis, Biology*; Gross, E., Meienhofer, J., Eds.; Academic Press: New York, 1979; p 197.
- (120) Curtius, T. *Ber. Deutsch. Chem. Ges.* (presently *Eur. J. Inorg. Chem.*) **1902**, *35*, 3226.
- (121) (a) El-Faham, A.; Albericio, F. *Chem. Rev.* **2011**, *111*, 6557. (b) Subirós-Funosas, R.; Khattab, S. N.; Nieto-Rodríguez, L.; El-Faham, A.; Albericio, F. *Aldrichimica Acta* **2013**, *46*, 21.
- (122) Carpino, L. A. *J. Am. Chem. Soc.* **1993**, *115*, 4397.
- (123) Coste, J.; Dufour, M. N.; Pantaloni, A.; Castro, B. *Tetrahedron Lett.* **1990**, *31*, 669.
- (124) Le Nguyen, D.; Heitz, A.; Castro, B. *J. Chem. Soc., Perkin Trans. 1* **1987**, 1915.
- (125) Anderson, G. W.; Paul, R. *J. Am. Chem. Soc.* **1958**, *80*, 4423.
- (126) Saha, A. K.; Schultz, P.; Rapoport, H. *J. Am. Chem. Soc.* **1989**, *111*, 4856.
- (127) Gibson, F. S.; Rapoport, H. *J. Org. Chem.* **1995**, *60*, 2615.
- (128) Albericio, F.; Boffill, J. M.; El-Faham, A.; Kates, S. A. *J. Org. Chem.* **1998**, *63*, 9678.

About the Authors

Siva S. Panda was born in Orissa, India, and received his bachelor's degree (B.Pharm.) in 2002 from the Roland Institute of Pharmaceutical Sciences, Berhampur, India, and his master's degree (M.Pharm.) in 2005 from Manipal College of Pharmaceutical Sciences, Manipal, India. He then joined the Dabur Research Foundation, Sahibabad, India, as a research scientist for a period of one year and worked on the synthesis of Mitoxantrone, an anticancer drug. Siva obtained his Ph.D. degree in synthetic organic chemistry in 2010 under the supervision of Prof. Subhash C. Jain at the University of Delhi, Delhi, India. He is currently working as a postdoctoral associate with Prof. Alan R. Katritzky at the Center for Heterocyclic Compounds, University of Florida, Gainesville, FL, where his research involves the synthesis of novel peptides, peptide bioconjugates, peptidomimetics in solution, and ligation studies of various cyclic transition states in S to N, N to N, and O to N acyl migrations.

C. Dennis Hall retired from his academic position at King's College, London, in 1999, and joined Prof. Katritzky's research group at the University of Florida, where he serves as a group leader,

administrator for the online journal *ARKIVOC*, and co-organizer of the Florida Heterocyclic and Synthetic Conferences (Flohets). Since joining Katritzky's group, he has co-authored some 40 papers in the fields of heterocyclic chemistry, QSAR, insect control, and synthetic ion channels.

Eric Scriven is a native of Wales, U.K. After working at BISRA and Esso Ltd, he attended the University of Salford and graduated in 1965 with a degree in chemistry. He obtained his M.Sc. Degree from the University of Guelph, and his Ph.D. degree in 1969 from the University of East Anglia (with Prof. Katritzky). After postdoctoral years at the University of Alabama and University College, London, he was appointed Lecturer in Organic Chemistry at the University of Salford. There, his research interests centered on the reactivity of azides and nitrenes. He joined Reilly Industries in 1979, where he served as Director of R&D from 1991 to 2003. He is currently at the University of Florida. He has edited two books: *Azides and Nitrenes: Reactivity and Utility* (1984) and *Pyridines: from Lab to Production* (2013). He and Professor H. Suschitzky were founding editors of *Progress in Heterocyclic Chemistry* that has been published annually since 1989. He collaborated with Professors Katritzky and Rees as Editors-in-Chief of *Comprehensive Heterocyclic Chemistry II* (1996), and with Professors Katritzky, Ramsden, and Taylor on the third edition of the work. Currently, he is Publishing Editor of *ARKIVOC*, an online journal of organic chemistry free to readers and authors.

Alan R. Katritzky was educated at Oxford and Cambridge (Lecturer and Founder Fellow of Churchill College). Founder Dean of the School of Chemical Sciences at East Anglia from 1962, he transferred in 1980 as inaugural Kenan Professor to the University of Florida. His research in heterocyclic chemistry has covered inter alia *N*-oxides, benzotriazole methodology, electrophilic and nucleophilic substitutions, computational QSPR relationships, and peptide chemistry. He holds 14 honorary doctorates from 10 Eurasian countries and associate or foreign memberships of five national academies. He is Cavaliere Ufficiale (Italy) and Honorary Fellow of St. Catherine's College, Oxford, and of the Polish and Italian Chemical Societies. Over 1,000 graduate students and postdocs have trained in his group. He created the not-for-profit Arkat USA, Inc., which organizes the Flohets Conferences and publishes the open-access journal *ARKIVOC*. Contributions to the secondary literature include editing *Comprehensive Heterocyclic Chemistry* (40 volumes in 3 editions), *Advances in Heterocyclic Chemistry* (106 volumes), *Handbook of Heterocyclic Chemistry* (3rd edition, 2010), and *Heterocycles in Life and Society* (2nd edition, 2011). ☞

For a custom synthesis partner you can trust

Add Aldrich

While our extensive listing of catalysts sets us apart, our capabilities extend far beyond the pages of a catalog. The Catalysis R&D team at Sigma-Aldrich is ready to work with you to meet your scale-up needs.

Knowledge

- Tap into our experience with over 45,000 manufactured products
- Benefit from our in-house expertise preparing the newest catalysts

Service

- Utilize our dedicated R&D resources for rapid development needs
- Leave the management of intellectual property to Sigma-Aldrich Technology Services

Partnership

- Work closely with a project manager to track and update progress
- Rely on seamless support from development to commercial scale



For more information, contact
catalysis@sial.com



For Grubbs' Z-Selective Ring-Closing Metathesis Catalyst

Add **Aldrich**

As the *exclusive* Materia partner for Grubbs Metathesis Technology, we offer Grubbs Ru-based Z-Selective Metathesis Catalyst (Product No. 771082). This catalyst effectively promotes both Z-selective metathesis and cross-metathesis reactions, as well as Z-selective macrocyclic ring-closing metathesis—a long-standing challenge in synthesis.


For inquiries regarding this Metathesis portfolio:
catalysis@sial.com

Z-Selective Macrocyclic Ring-Closing Metathesis




771082
Grubbs' Z-Selective
Catalyst

Complementing Our *Exclusive* Offering of Materia Metathesis Catalysts



569747
Grubbs Catalyst
2nd Generation



569755
Hoveyda-Grubbs Catalyst
2nd Generation

Reference:
Marx, V. M.; Herbert, M. B.; Keitz, B. K.; Grubbs, R. H. *J. Am. Chem. Soc.* 2013, 135, 94



For a complete listing of
Grubbs catalysts, visit
Aldrich.com/grubbsacta


New Heller-Sarpong reagents for acid derivatization

Add Aldrich



Introducing the Heller-Sarpong Reagents, which provide esters and Weinreb amides in a single step from the corresponding carboxylic acids. This new portfolio provides:

- An operationally simple reaction protocol
- Broad functional group tolerance
- A safer alternative to diazoalkanes

Esterification




Weinreb Amide Synthesis




748838
WImC




748730
MImC




748781
PlmC




748811
AllyImC




748773
EtImC




748765
IPlmC



752886



761052



748846


References:

- (1) Heller, S. T.; Sarpong, R. *Org. Lett.* 2010, 12, 4572
 (2) Heller, S. T.; Fu, T.; Sarpong, R. *Org. Lett.* 2012, 14, 1970




Find out more about
Heller-Sarpong products:
[Aldrich.com/
imidazolecarbamatesacta](http://Aldrich.com/imidazolecarbamatesacta)

Recent Advances in the Prins Cyclization



Prof. Dr. S. J. Greco



Mr. R. G. Fiorot



Prof. Dr. V. Lacerda, Jr.



Prof. Dr. R. B. dos Santos

Sandro José Greco,* Rodolfo Goetze Fiorot, Valdemar Lacerda, Jr., and Reginaldo Bezerra dos Santos

Department of Chemistry
Federal University of Espírito Santo
Avenida Fernando Ferrari, 514
Goiabeiras, Vitória, ES 29075-910, Brasil
Email: sandro.greco@ufes.br

Keywords. Prins reaction; tetrahydropyrans; dihydropyrans; tetrahydrofurans; dioxanes; piperidines; azepines; lactones; spiro compounds; macrocycles; natural products.

Abstract. The Prins reaction is often a key step in the synthesis of various heterocyclic rings that are important structural components of many classes of biologically active compounds and natural products. This review presents and discusses recent significant applications of this important reaction, and offers insight into its mechanism and regio- and stereochemical outcomes.

Outline

1. Introduction
2. Mechanistic Considerations
 - 2.1. Regioselectivity for 5- vs 6-Membered Rings
 - 2.2. Stereoselectivity of Nucleophile Capture at C4 of the Tetrahydropyran Ring
 - 2.2.1. Alder's Model: Equatorial Selectivity
 - 2.2.2. Rychnovsky's Model: Axial Selectivity
 - 2.3. Diastereoselectivity of the Prins Cyclization
3. Recent Synthetic Applications
 - 3.1. Substituted 1,3-Dioxanes
 - 3.2. Spiro and Bicyclic Tetrahydropyrans (THPs)
 - 3.3. Nitrogen Heterocycles via the Aza-Prins Cyclization
 - 3.4. Synthesis of Furan Derivatives
 - 3.5. Synthesis of Halogenated Tetrahydropyrans: Halo-Prins Cyclization
 - 3.5.1. Fluorinated THPs
 - 3.5.2. Chlorinated and Brominated THPs
 - 3.5.3. Iodinated THPs
 - 3.6. Macrocyclization Involving the Prins Reaction
4. Conclusion
5. References


1. Introduction

The Prins reaction is often related to the Kriewitz reaction, an example of which is the reaction of β -pinene with paraformaldehyde to produce an unsaturated alcohol through a thermal ene rearrangement (Scheme 1, Part (a)).¹⁻³ When the reaction between an alkene and formaldehyde is conducted in the presence of an acid catalyst—such as the reaction of styrene with paraformaldehyde in the presence of aqueous sulfuric acid to give diol **1** (Scheme 1, Part (b))—it is called the Prins reaction. Numerous protic and Lewis acids are known to catalyze the reaction, and excellent reviews have been published on the early work.¹⁻³ The products generally obtained in this condensation are formed as complex mixtures of 1,3-dioxanes, 1,3-glycols, tetrahydropyrans, and allylic and homoallylic alcohols, with the composition of the mixture being dependent on the specific experimental conditions employed (Scheme 1, Part (c)). In the presence of water, the intermediate carbocation leads to the formation of 1,3-glycols and 1,3-dioxanes, while 3-alkyl-4-halotetrahydropyrans are obtained through the intermediacy of the homoallylic alcohols.


The Prins reaction plays a key role in the synthesis of such important product classes as dihydrofurans, dihydropyrans, piperidines, and oxabicyclo and spiro compounds. This review presents a survey and a discussion of pertinent and interesting recent developments relating to the stereochemical course and mechanism of the Prins reaction and to its advantageous application in organic synthesis.

2. Mechanistic Considerations


The Prins cyclization involves a homoallylic alcohol, an aldehyde, and a Lewis acid. The latter acts as catalyst and, depending on experimental conditions, it can also serve as a source of a nucleophilic anion. In the presently accepted mechanism, the reaction is initiated by complexation of the Lewis acid with the aldehyde, which activates the carbonyl carbon toward attack by the hydroxyl group of the alcohol, generating the hemiacetal intermediate **2**. Loss of the Lewis acid fragment from



Scheme 1. Condensation Reactions of Olefins with Paraformaldehyde and the Dependence of Product Distribution on the Reaction Conditions. (Ref. 1–3)



Scheme 2. General Mechanism of the Prins Cyclization. (Ref. 2)



Scheme 3. Regioselectivity in the Prins Cyclization. (Ref. 4)

the hemiacetal forms the key oxonium ion intermediate, **3**, which assumes the more stable chair conformation in which the substituents are pseudoequatorial. Subsequent 6-*endo* cyclization of **3** selectively leads to secondary tetrahydropyranyl carbocation **4**, which captures the halide to give rise to the 2,4,6-trisubstituted tetrahydropyran product (**Scheme 2**).²

2.1. Regioselectivity for 5- vs 6-Membered Rings

When the double bond geometry in the homoallylic alcohol is switched from *E* to *Z*, tetrahydrofurans can be formed in competition with tetrahydropyrans. This regioselectivity can be studied by examining the stereochemistry of intermediates present in the accepted mechanism of this reaction. Under the Prins cyclization conditions, the *Z* homoallylic alcohol reacts with the activated aldehyde to give rise to oxonium ion **3**. Two competing transition states can then be formed from **3**: six-membered-ring transition state **6** has a 1,3-diaxial interaction between H and the substituent R^1 , while five-membered-ring transition state **8** has greater torsional and angular strains. When the R^1 substituent is sufficiently large, an increase in the activation barrier of the process results, which slows the formation of tetrahydropyran product **7** in favor of tetrahydrofuran product **9** (**Scheme 3**).⁴

2.2. Stereoselectivity of Nucleophile Capture at C4 of the Tetrahydropyran Ring

2.2.1. Alder's Model: Equatorial Selectivity

On the basis of theoretical calculations employing Density Functional Theory (DFT), Alder and co-workers concluded that the all-*cis* 2,4,6-trisubstituted product of the Prins cyclization is favored by stabilization of the cationic intermediate **10** through hyperconjugation. When the hydrogen attached to the carbocation center is pseudoaxial, the empty p orbital of the positively charged carbon overlaps more efficiently with the coplanar σ_{C-C} and σ_{C-O} orbitals and with the orbital of the nonbonding electrons of oxygen. Nucleophilic attack thus occurs from the exo face (convex), leading to the 2,4,6-trisubstituted tetrahydropyran product with all three substituents in equatorial positions (**eq 1**).⁵

2.2.2. Rychnovsky's Model: Axial Selectivity

Rychnovsky investigated the capture of bromide and iodide at C4 of **10** by reacting α -acetoxy ester **11** with TMSBr, AcBr, HBr, or TMSI and lutidine in dichloromethane. High axial stereoselectivity at C4 was observed for the resulting Prins cyclization product **16** (**Scheme 4**).⁶ In contrast, when $SnBr_4$ was employed, the major product was the equatorial epimer **19**. In the proposed mechanism, some Lewis acids, such as TMSBr, AcBr, and HBr, act as donors of bromide by forming the intimately associated ion pairs **14** and **15**. The slightest movement (least motion pathway) in **15**, results in Br⁻ attacking C4 in the axial position (endo attack) to form **16**. When $SnBr_4$ is employed as the Lewis acid, the in situ formed $[SnBr_3]^+$ in ion pair **17** is less nucleophilic than bromide, allowing separation of the ion pair by the solvent. In the resulting intermediate, **18**, exo (convex) attack leads to the formation of product **19** with an equatorial bromine at C4.

2.3. Diastereoselectivity of the Prins Cyclization

Substituents at C2 and C6 of tetrahydropyrans formed by the Prins cyclization are preferentially *cis*.^{7a} Methodologies for forming the C2/C6 anti isomers are not well established. These isomers are present in some structures of natural products, such as the psymberins^{7b} and the apicularens.^{7c} Panek's group has succeeded in synthesizing enantiomerically pure *anti*-2,6-dihydropyrans by the Prins cyclization with the aid of TMSOTf.^{7c-e}

Loh and co-workers investigated steric and electronic effects in the Prins cyclization of homoallylic *anti*- α -hydroxy esters, leading to 4-chloro-2,6-disubstituted THPs. This study demonstrated that groups with high electron density in the pseudoaxial position stabilize the oxonium ion by inductive electronic effects, and favor a transition state that forms the *anti* isomer, **20a**; while steric effects favor the transition state leading to the *syn* isomer, **20b**. In both cases, equatorial attack of the nucleophile is preferred (Scheme 5).⁸



Alder's model for nucleophile capture at C4 in the Prins cyclization

eq 1 (Ref. 5)

3. Recent Synthetic Applications

Saturated six-membered-ring oxygen and sulfur heterocycles are features found in the structures of a variety of biologically important natural products such as polyether antibiotics, marine toxins, pheromones, and pharmaceutical agents. Tetrahydropyran is also the structural core of most carbohydrates, oligomers, and polymers, which play crucial roles in living organisms. It is therefore not surprising that considerable efforts have been expended toward developing facile and viable syntheses of tetrahydropyran-containing compounds.

3.1. Substituted 1,3-Dioxanes


When olefins are condensed with aldehydes in aqueous solutions of mineral acid catalysts, alkyldioxanes (cyclic formals or acetals of 1,3-butanediols) and 1,3-butanediols are formed. The distribution of these two products varies with the concentration of the solution of the acid catalyst and the reaction temperature. Amrute et al. studied the catalytic activity of $\text{MoO}_3/\text{SiO}_2$ (7 wt %) in the Prins cyclization of a series of olefins with paraformaldehyde (2 equiv) in 1,2-dichloroethane at 80 °C.⁹ 72–90% conversions and 96–100% selectivities were observed for the corresponding 4-alkyl- and 4-phenyl-substituted 1,3-dioxane products. Du and Tian synthesized 1,3-dioxanes in moderate-to-high yields from formalin (aqueous formaldehyde), styrene derivatives, and trifluoromethanesulfonic acid (TfOH).¹⁰ The use of organic acid as catalyst for the Prins reaction was unprecedented. It is worth noting that this approach avoids the use of organic solvents by conducting the reaction in water, which makes it a more environmentally friendly process.

Yang and co-workers explored the use of water-stable and recyclable Brønsted acidic ionic liquids as environmentally benign catalysts for the Prins cyclization. The effectiveness of these ionic liquids was compared in the model reaction of styrene with formaldehyde at 94–96 °C, whereby $[\text{BMIM}][\text{HSO}_4]$ was found to be the most effective catalyst (eq 2).¹¹ The 1,3-dioxane products were obtained in good yields, and the catalysts, after vacuum distillation at 80 °C, were recovered and reused in subsequent runs, thus reducing the risks to the environment by avoiding the use of organic solvents and enabling large-scale applications of the Prins cyclization.


3.2. Spiro and Bicyclic Tetrahydropyrans (THPs)

Gais and co-workers have reported a modular asymmetric approach to spiroketals, spiroethers, and oxabicycles that employs a spiro- or bicycloannulation of α -hydroxydihydropyrans. The synthesis included a stereoselective Ferrier-type O- and C-glycosidation, ring-closing metathesis, and stereoselective Prins cyclization as key steps. When α -hydroxydihydropyran **21** was treated with TiCl_4 in dichloromethane at –78 °C, the regioisomeric spiro ethers **22** and **23** resulted from a Prins cyclization (eq 3).¹²



Nakamura and co-workers have described a versatile method for the synthesis of 4-substituted 6-methyl-3-oxabicyclo[3.3.1]non-6-ene-1-methanol derivatives using a Prins-type cyclization reaction between aldehydes and O-protected or unprotected 4-methylcyclohex-




Scheme 4. Synthesis of 1,3-dioxanes by Prins cyclization




Scheme 5. Diastereoselectivity of the Prins Cyclization. (Ref. 8)




eq 2 (Ref. 11)




eq 3 (Ref. 12)



Six additional examples of aromatic aldehydes, 72–96%



eq 4 (Ref. 13)



eq 5 (Ref. 17)

3-ene-1,1-dimethanol (**24**). Under the optimized reaction conditions employing hafnium triflate, various aldehydes, including functionalized benzaldehydes and heteroaromatic aldehydes, afforded the cyclization products in high yields (eq 4).¹³ The zinc and lanthanum triflates form preferably spirodioxane **26**.

3.3. Nitrogen Heterocycles via the Aza-Prins Cyclization

The Prins cyclization of homoallylic amines (the aza-Prins cyclization) takes place in a fashion similar to that of homoallylic alcohols, whereby the nonbonding electrons on the nitrogen initiate the sequence of reaction steps by attacking the electrophilic site of the aldehyde activated by an acid catalyst. The key intermediate of the aza-Prins cyclization is an iminium ion, in analogy to the oxonium ion. Piperidines are commonly found subunits in many biologically relevant molecules including alkaloids, and are attractive structural scaffolds for drug discovery.¹⁴ Subba Reddy showed that the BF₃·OEt₂ catalyzed aza-Prins reaction of benzaldehyde and *N*-tosyl-3-butenamine (a homoallylic amine), in the presence of anisole as solvent and nucleophile, produces the *trans*-2,4-diarylpiperidine in 83% yield.^{15,16} This is the first report of the preparation of 4-arylpiperidines via an aza-Prins–Friedel–Crafts reaction sequence.

The coupling of *E*- (**27**) and *Z*-3-hexene-1,6-ditosylamides with various aldehydes, including cinnamaldehyde (**28**) in the presence of 10 mol % Sc(OTf)₃ in 1,2-dichloroethane gave the corresponding *trans*- and *cis*-fused saturated pyrrolopyridines **29** and **30**, respectively, in good yields by an intramolecular aza-Prins cyclization (eq 5).¹⁷ Other aromatic aldehydes such as benzaldehyde, *para*-anisaldehyde, and thiophene-2-carboxaldehyde were not effective substrates in the reaction. Ketones, such as cyclohexanone, failed to give the spirodiaz bicyclic product. In contrast, aliphatic aldehydes, such as isovaleraldehyde (70%, *trans*:*cis* = 95:5), cyclohexanecarboxaldehyde (73%, *trans*:*cis* = 95:5), and propionaldehyde (66%); participated well in this reaction.

Camara et al. synthesized azepines fused to a naphthoquinone moiety by an intramolecular aza-Prins cyclization starting with an amino derivative, **31**, of lapachol (Scheme 6).¹⁸ Products **32a** and **32b** were formed in 42% yield as a diastereomeric mixture, with a *trans*:*cis* ratio of about 7:3. The mechanism of C–C bond formation leading to **32a** and **32b** appears to resemble that of the intramolecular ene reaction between a carbonyl group and an alkene.¹⁹ The authors proposed that the formation of intermediates occurs through a Prins reaction, via nucleophilic attack of H₂O or MeOH at the isoprenyl double bond, possibly followed by a concerted attack onto the protonated carbonyl. The observed diastereoselectivity is possibly induced by steric hindrance of the 4-isopropyl and 3-hydroxyl groups, despite the fact that the resulting seven-membered ring is conformationally less restricted than the corresponding six-membered ring. This synthetic method is important, since there are very few publications on the synthesis of such heterocyclic systems, which are of great interest in the scientific community because of their pharmacological applications.^{20,21}

3-Azabicyclo[3.3.1]non-6-enes are structural motifs in many natural products and, with the proper choice of substituents, can serve as templates for complexity-generating transformations. Krasavin and co-workers have reported a facile synthesis of this ring system in a diastereomerically pure form by an aza-Prins cyclization involving a δ,ϵ -unsaturated imine and an equivalent of BF₃·OEt₂ under microwave irradiation at 180 °C for 1 h (Scheme 7).²² As in the mechanism proposed earlier by Overman for the aza-Prins cyclization employed in the total synthesis of (+)-nankakurines A and B,²³ it is believed that the pair of nonbonding electrons on nitrogen participate in a regioselective

intramolecular deprotonation of the hydrogen vicinal to the initially formed carbocation in intermediate **A**.

Subba Reddy and co-workers²⁴ have reported the synthesis of 2-aryl- and 2-alkyl-4-amidopiperidines in good yields and high selectivities by an aza-Prins–Ritter tandem reaction employing a slight excess of triflic acid in MeCN at 0 °C. It was observed that in the absence of triflic acid, no aza-Prins cyclization occurred even in refluxing acetonitrile. Other Brønsted acids—such as acetic acid, formic acid, and trifluoroacetic acid—were tested and, in all cases, the reaction proceeded rapidly at 0 °C; however, triflic acid provided the best conversion.

3.4. Synthesis of Furan Derivatives

Substituted dihydrofuranones (γ -butyrolactones or GBLs) are important intermediates in synthetic organic chemistry, and are commonly found as structural fragments in natural products, receptor ligands, and drug molecules.²⁵ Compounds containing a GBL moiety exhibit pharmacological effects some of which are muscarinic (pilocarpine) and antimuscarinic (Kaiser lactones) activities, convulsant (picrotoxin, β -substituted GBL) and anticonvulsant (α -alkyl-substituted GBLs) activities, and the ability to modulate quorum sensing.²⁶

Gao and Canney reported a novel and concise approach for the synthesis of structurally diverse, substituted 5-(2-hydroxyethyl)-3,3-dihydrofuran-2(3*H*)-ones. This method relies on a modified Prins reaction that employs a catalytic amount of H₂SO₄ in glacial HOAc. In the proposed mechanism, acetic acid captures the initially formed carbocation to provide a protected hydroxyl group. The remaining aliphatic alcohol is also protected by an intramolecular esterification, leading to a caprolactone intermediate. Sequential treatment of the resulting seven-membered-ring lactone with aqueous base and then acid affords the desired hydroxyethyl lactones in moderate-to-good yields (**Scheme 8**).²⁶

Two examples of the use of the Prins cyclization in the synthesis of furan structures fused with pyran rings (furopyrans) have recently been reported.^{27,28} This ring system is commonly found in natural products such as flavonoids, pterocarpan, and catechins.²⁹ In the first example, 3-hexene-1,6-diol (**36**) was reacted with 4-bromobenzaldehyde in an intramolecular Prins cyclization in the presence of 10% TsOH in 1,2-dichloroethane to give the cis-fused bicyclic product **37** in 72% yield (**Scheme 9**, Part (a)).²⁷ In the second example, a substituted 2,6-dioxabicyclo[3.2.1]octane was similarly prepared by a tandem acetalization–intramolecular Prins cyclization starting with 4-pentene-1,2-diol (**Scheme 9**, Part (b)).²⁸ The overall process is catalyzed by a combination of Sc(OTf)₃ and TsOH and leads to good yields, high selectivities, and faster reaction times.


3.5. Synthesis of Halogenated Tetrahydropyrans: Halo-Prins Cyclization

In analogy to the oxygen and nitrogen variants, the halo-Prins cyclization involves nucleophilic attack by halogen present in the reaction medium on the carbocation intermediate that arises from the oxonium ion formed after cyclization.


3.5.1. Fluorinated THPs

The introduction of fluorine atoms into organic molecules alters in important ways their biological activity, solubility, hydrophobicity, metabolism, and bulk properties.^{30–32} However, few methods for the synthesis of fluorinated pyranyl motifs are known and, of these, the ones that employ BF₃•Et₂O and Et₃NF•5HF as both Lewis acids and fluorine sources successfully achieve the Prins cyclization of homoallylic alcohols into fluorinated pyranyl motifs.^{33a–c} When BF₃•Et₂O is utilized

in stoichiometric quantities, it contributes fluoride ion to quench the intermediate carbocation, giving rise to the fluorinated products. As an example of this approach, O'Hagan and co-workers investigated the oxa- and aza-Prins reactions for the synthesis of 4-fluoropyrans and 4-fluoropiperidines starting from homoallylic alcohols and various aldehydes. The fluorinated THP products were obtained in good yields, but with only moderate diastereoselectivity.³⁴ This method was extended to the aza-Prins reaction that utilizes *N*-tosylhomoallyl amines to generate the corresponding 4-fluoropyrrolidines.




Scheme 6. Intramolecular Aza-Prins Cyclization of a Derivative of Lapachol (Ref. 18).




Scheme 7. Synthesis of Azabicyclo Compounds through an Aza-Prins Cyclization. (Ref. ...)

O'Hagan's group then investigated the Prins fluorination reactions under microwave conditions, and observed significantly reduced reaction times and higher conversions. However, there was a slight decrease in the diastereoselectivity of the reaction and, in some cases, an inversion of the diastereoselectivity. When a series of low-temperature ($-20\text{ }^{\circ}\text{C}$) experiments were carried out in an attempt to improve the diastereoselectivity, dr increased from $\sim 2:1$ to $10:1$ and yields remained good, but, not surprisingly, a significant increase in the reaction time was observed.

Prior to Loh and co-workers' recent disclosure,³⁵ all studies of the Prins fluorination reaction reported the almost exclusive formation of *cis*-2,6-disubstituted fluorinated di- or tetrahydropyrans. An efficient, highly diastereoselective synthesis of the *trans*-2,6-disubstituted



Scheme 8. Synthesis of Dihydrofuranones (Ref. 34)



Scheme 9. Use of Unsaturated 1,6-Diol **36** for the Preparation of Furopyrans (Ref. 27,28)

counterparts—useful in the development of new pharmaceuticals—would be highly desirable. Loh's group explored the Prins reaction of various allenic alcohols, e.g. **38**, with a variety of aldehydes using different Lewis acids (LAs) and fluorine sources.³⁵ Their research demonstrated that $\text{BF}_3 \cdot \text{Et}_2\text{O}$, acting both as an efficient Lewis acid and as a source of fluoride, gives the best results. The authors proposed a mechanism in which the Prins cyclization of the allenic alcohol takes place through a distorted chair transition state, in which a lone electron pair on the carbonyl oxygen of the ester group stabilizes the partial positive charge on the oxocarbenium carbon. This forces the carbonyl group to adopt an axial orientation, leading to the desired intermediate **40** and suppressing the generation of the undesirable intermediate **41**. In turn, intermediate **40** gives rise to the desired *trans*-2,6-disubstituted fluorinated dihydropyran **39**, selectively (Scheme 10).³⁵

Saikia and co-workers have reported that TiF_4 can efficiently be employed for the stereoselective synthesis of substituted all-*cis* 4-fluorotetrahydropyrans via the halo-Prins cyclization.³⁶ A variety of aliphatic and aromatic aldehydes were reacted with a number of homoallylic alcohols to give good yields and high diastereoselectivities of the corresponding THPs. Moreover, acyclic and cyclic ketones were subjected to the reaction and found to be less reactive (5–6 h vs 2.5–4 h for the aldehydes), giving only moderate yields (50–70% vs 80–92% for the aldehydes). Cyclic ketones afford spiro compounds, as illustrated by the reaction of cyclohexanone, which leads to spirocyclic compound **42** in 70% yield (eq 6).³⁶

1.5.2. Chlorinated and Brominated THPs

InCl_3 has been demonstrated to be an excellent Lewis acid for the insertion of a chlorine atom at the 4 position of THPs by the Prins cyclization. For example, the InCl_3 -promoted diastereoselective Prins reaction of **43** with benzaldehyde led to the pentasubstituted tetrahydropyran derivative **44** (essentially as a single product) in which five stereogenic centers (up-down-up-down-up) were controlled (Scheme 11, Part (a)).^{4,37,38} Subba Reddy and co-workers³⁹ reported another example of the effectiveness of InCl_3 , whereby the synthesis of *cis*-fused hexahydro-1*H*-furo[3,4-*c*]pyran scaffolds containing chlorine proceeded smoothly under mild conditions (Scheme 11, Part (b)).³⁹

FeCl_3 has been employed as an inexpensive, environmentally friendly, and stable Lewis acid to promote the halo-Prins cyclization of 3-buten-1-ol with several aldehydes.⁴⁰ The cyclization affords the corresponding *cis*-4-halo-2-alkyltetrahydropyrans in generally excellent yields.⁴⁰ The reaction works quite well with both aliphatic and aromatic aldehydes and, when FeBr_3 is employed, the 4-bromo-substituted analogue is formed. Liu and Loh have disclosed an efficient and highly stereoselective Prins cyclization leading to *cis*-2,6-dialkyl-3,4-dibromotetrahydropyrans from terminal vinyl bromides. This method employs InBr_3 as the Lewis acid and TMSBr as the source of bromide ion.⁴¹

Cascade reactions can be very powerful transformations in organic synthesis.⁴² The first examples of a Mukaiyama aldol–Prins (MAP) cascade cyclization reaction were reported by Rychnovsky's group, whereby a very reactive allylsilane served as the internal nucleophile in a rapid and clean Prins cyclization.⁴³ Rychnovsky and co-workers also described the use of simple alkene substrates in MAP cyclizations and the importance of selecting the appropriate Lewis acid to promote the reaction. The attraction of such a sequence is that it forms two new C–C bonds, a ring, and three new stereogenic centers. Initially, the Mukaiyama aldol addition and Prins cyclization with the simple alkene **45** was evaluated using previously optimized conditions.⁴³ The reaction of **45** with 2.5 equiv of dihydrocinnamaldehyde in the presence

of $\text{BF}_3 \cdot \text{OEt}_2$ and 2,6-di-*tert*-butylpyridine (2,6-DTBP) at -78°C led to the unexpected product **46** in 82% yield. The more powerful Lewis acids, TiCl_4 and TiBr_4 , did not produce 1,3-dioxane **46** and gave the best yields of **47**. TiBr_4 was particularly effective and gave adduct **47** in 72% yield (eq 7).⁴⁴

3.5.3. Iodinated THPs

The mild Lewis acidic nature of molecular iodine⁴⁵ has been exploited by Yadav and co-workers in the first direct and metal-catalyst-free Prins cyclization of homoallylic alcohols with aldehydes for the rapid synthesis of highly substituted iododihydro- and iodotetrahydropyrans in good yields and selectivities under neutral conditions.^{46,47} Other reagents such as LiI , KI , and NaI failed to produce the desired product. Aliphatic, simple aromatic, and moderately activated aldehydes gave higher yields of products than strongly activated or deactivated aldehydes. This same research group reported a simple and metal-catalyst-free Prins cyclization for the synthesis of highly substituted tetrahydropyrans from sugar-based homoallylic alcohols and aldehydes using molecular iodine under neutral conditions.⁴⁶


When silylated secondary homopropargylic alcohols were subjected to the same experimental protocol with aldehydes, highly substituted dihydropyrans were rapidly formed by a Prins cyclization (Scheme 12).⁴⁷ The reactions were completed in 3 hours or less, and the all-cis products were obtained in 72–90% yields. The cis selectivity presumably arises from an *E* oxocarbenium ion formed via a chair-like transition state. The optimal geometry of this oxocarbenium ion places the hydrogen atom at C4 in a pseudoaxial position, which favors equatorial attack of the nucleophile.

TMSI, generated in situ from TMSCl and NaI , has been employed as an iodide source in the Prins cyclization of homoallylic and homopropargylic alcohols with various ketones leading to 2,2-disubstituted 4-iodotetrahydropyrans, spirocyclic 4-iodotetrahydropyrans, and spirocyclic 4-iodo-5,6-dihydro-2H-pyrans.⁴⁸ In the presence of iodide, no trapping of the intermediate 4-tetrahydropyranyl carbocation by acetonitrile was detected.


The synthesis of 4-iodotetrahydropyrans by the Prins cyclization can also be performed in the presence of a combination of $\text{CeCl}_3 \cdot 7\text{H}_2\text{O}$ and LiI .⁴⁹ This approach can be applied to both aldehydes and ketones, but requires a higher temperature (reflux in dichloroethane), with the best results achieved with 1 equiv of CeCl_3 . In contrast, CAN and $\text{Ce}(\text{OTf})_3$ were not effective and led to hydroxylated side products. Gallium triiodide (35 mol %) has also been employed for the synthesis of 4-iodotetrahydropyrans at room temperature in 10–25 min and in 82–89% yields.⁵⁰

3.6. Macrocyclization Involving the Prins Reaction


Oxacyclic macrolidomers constitute an important class of natural products that possess a wide range of structural complexity and bioactivity,⁵¹ and are popular targets for synthetic chemists.⁵² Rychnovsky and co-workers introduced a new sequential dimerization–macrocyclization based on the Prins cyclization for forming symmetrical macrocycles.⁵³ This approach is illustrated by the optimized conditions in the example in equation 8. A variety of Lewis acids, including other Re and non-Re ones, were investigated and found to be inferior to $\text{O}_3\text{ReOSiPh}_3$. Other reaction parameters; such as temperature, concentration, and substrate scope; were also examined: Both acetals and aldehydes were found to be viable substrates, whereas increasing the reaction temperature to 40°C did shorten the reaction times but did not improve the yields. The usefulness of this strategy was demonstrated in a successful synthesis of a model for clavosolide A, a marine sponge metabolite.




Scheme 10. Fluorinated Dihydropyrans by the Prins Cyclization of Allenic Alcohols. (Ref. 35)



eq 6 (Ref. 36)



Scheme 11. Chlorinated THPs by the Prins Cyclization of Homoallylic Alcohols. (Ref. 4.39)




Six additional examples using TiBr₄:
53–82% of 47, dr = 1:1 to 1.3:1


	LA	X	46	47
BF ₃ ·OEt ₂	---	---	82%	nd
TiCl ₄	Cl	Cl	nd	62%
TiBr ₄	Br	Br	nd	72%

nd = not detected.

eq 7 (Ref. 44)



Scheme 12. Prins Cyclization of Silylated Secondary Homopropargylic Alcohols with Aldehydes under Mild, Metal-Catalyst-Free Conditions. (Ref. 47)



eq 8 (Ref. 53)

4. Conclusion

The Prins cyclization is an efficient reaction for the stereoselective synthesis of substituted tetrahydropyran rings, and has significantly advanced in the past few years, as demonstrated by the number of applications described in the literature. This reaction is, in some cases, the method of choice for the preparation of natural products and biologically active compounds that feature the tetrahydropyran moiety in their structures.

5. References

- Arundale, E.; Mikeska, L. A. *Chem. Rev.* **1952**, *51*, 505.
- Olier, C.; Kaafarani, M.; Gastaldi, S.; Bertrand, M. P. *Tetrahedron* **2010**, *66*, 413.
- (a) Hanschke, E. *Chem. Ber.* (presently *Eur. J. Inorg. Chem.*) **1955**, *88*, 1053. (b) Vasconcellos, M. L. A. A.; Miranda, L. S. M. *Quim. Nova* **2006**, *29*, 834.
- Yang, X.-F.; Mague, J. T.; Li, C.-J. *J. Org. Chem.* **2001**, *66*, 739.
- Alder, R. W.; Harvey, J. N.; Oakley, M. T. *J. Am. Chem. Soc.* **2002**, *124*, 4960.
- Jasti, R.; Vitale, J.; Rychnovsky, S. D. *J. Am. Chem. Soc.* **2004**, *126*, 9904.
- (a) Capim, S. L.; Carneiro, P. H. P.; Castro, P. C.; Barros, M. R. M.; Marinho, B. G.; Vasconcellos, M. L. A. A. *Eur. J. Med. Chem.* **2012**, *58*, 1. (b) Cichewicz, R. H.; Valeriote, F. A.; Crews, P. *Org. Lett.* **2004**, *6*, 1951. (c) Su, Q.; Panek, J. S. *J. Am. Chem. Soc.* **2004**, *126*, 2425. (d) Lowe, J. T.; Panek, J. S. *Org. Lett.* **2000**, *7*, 3231. (e) Huang, H.; Panek, J. S. *J. Am. Chem. Soc.* **2000**, *122*, 9836.
- Chan, K.-P.; Seow, A.-H.; Loh, T.-P. *Tetrahedron Lett.* **2007**, *48*, 37.
- Amrute, A. P.; Sahoo, S.; Bordoloi, A.; Hwang, Y. K.; Hwang, J.-S.; Halligudi, S. B. *Catal. Commun.* **2009**, *10*, 1404.
- Du, Y.; Tian, F. *Catal. Commun.* **2007**, *8*, 2012.
- Wang, W.; Shao, L.; Cheng, W.; Yang, J.; He, M. *Catal. Commun.* **2008**, *9*, 337.
- Lejkowski, M.; Banerjee, P.; Runsink, J.; Gais, H.-J. *Org. Lett.* **2008**, *10*, 2713.
- Nakamura, M.; Niiyama, K.; Yamakawa, T. *Tetrahedron Lett.* **2009**, *50*, 6462.
- Pastor, I. M.; Yus, M. *Curr. Org. Chem.* **2012**, *16*, 1277.
- Yadav, J. S.; Subba Reddy, B. V.; Ramesh, K.; Narayana Kumar, G. G. K. S.; Grée, R. *Tetrahedron Lett.* **2010**, *51*, 818.
- Yadav, J. S.; Subba Reddy, B. V.; Ramesh, K.; Narayana Kumar, G. G. K. S.; Grée, R. *Tetrahedron Lett.* **2010**, *51*, 1578.
- Subba Reddy, B. V.; Borkar, P.; Chakravarthy, P. P.; Yadav, J. S.; Grée, R. *Tetrahedron Lett.* **2010**, *51*, 3412.
- Camara, C. A.; Pinto, A. C.; Vargas, M. D.; Zukerman-Schpector, J. *Tetrahedron* **2002**, *58*, 6135.
- Achmatowicz, O.; Bialecka-Florjańczyk, E. *Tetrahedron* **1996**, *52*, 8827.
- Taylor, E. C.; Dowling, J. E. *Bioorg. Med. Chem. Lett.* **1997**, *7*, 453.
- Link, A.; Kunick, C. *J. Med. Chem.* **1998**, *41*, 1299.
- Parchinsky, V.; Shumsky, A.; Krasavin, M. *Tetrahedron Lett.* **2011**, *52*, 7157.
- Altman, R. A.; Nilsson, B. L.; Overman, L. E.; de Alaniz, J. R.; Rohde, J. M.; Taupin, V. *J. Org. Chem.* **2010**, *75*, 7519.
- Subba Reddy, B. V.; Ramesh, K.; Ganesh, A. V.; Narayana Kumar, G. G. K. S.; Yadav, J. S.; Grée, R. *Tetrahedron Lett.* **2011**, *52*, 495.
- (a) Wu, X.; Cao, W.; Zhang, H.; Chen, J.; Jiang, H.; Deng, H.; Shao, M.; Zhang, J.; Chen, H. *Tetrahedron* **2008**, *64*, 10331. (b) Ok, T.; Jeon, A.; Lee, J.; Lim, J. H.; Hong, C. S.; Lee, H.-S. *J. Org. Chem.* **2007**, *72*, 7390. (c) Mori, K. *Tetrahedron* **1983**, *39*, 3107.
- Gao, R.; Canney, D. J. *Tetrahedron Lett.* **2009**, *50*, 5914, and references therein.
- Yadav, J. S.; Chakravarthy, P. P.; Borkar, P.; Subba Reddy, B. V.; Sarma, A. V. S. *Tetrahedron Lett.* **2009**, *50*, 5998.
- Subba Reddy, B. V.; Narasimhulu, G.; Vikram Reddy, Y.; Chakravarthy, P. P.; Yadav, J. S.; Sridhar, B. *Tetrahedron Lett.* **2012**, *53*, 3100.
- Génisson, Y.; Tyler, P. C.; Young, R. N. *J. Am. Chem. Soc.* **1994**, *116*, 759.
- Welch, J. T. *Tetrahedron* **1987**, *43*, 3123.
- Purser, S.; Moore, P. R.; Swallow, S.; Gouverneur, V. *Chem. Soc. Rev.*



- 2008, 37, 320.
- (32) McCarthy, J. R.; Matthews, P. D.; Stemerick, D. M.; Huber, E. W.; Bey, P.; Lippert, B. J.; Snyder, R. D.; Sunkara, P. S. *J. Am. Chem. Soc.* **1991**, 113, 7439.
- (33) (a) Kishi, Y.; Nagura, H.; Inagi, S.; Fuchigami, T. *Chem. Commun.* **2008**, 3876. (b) Kataoka, K.; Ode, Y.; Matsumoto, M.; Nokami, J. *Tetrahedron* **2006**, 62, 2471. (c) Dobbs, A. P.; Pivnevi, L.; Penny, M. J.; Martinović, S.; Iley, J. N.; Stephenson, P. T. *Chem. Commun.* **2006**, 3134. (d) Al-Mutairi, E. H.; Crosby, S. R.; Darzi, J.; Harding, J. R.; Hughes, R. A.; King, C. D.; Simpson, T. J.; Smith, R. W.; Willis, C. L. *Chem. Commun.* **2001**, 835. (e) Jaber, J. J.; Mitsui, K.; Rychnovsky, S. D. *J. Org. Chem.* **2001**, 66, 4679.
- (34) Launay, G. G.; Slawin, A. M. Z.; O'Hagan, D. *Beil. J. Org. Chem.* **2010**, 6, No. 41.
- (35) Luo, H.-Q.; Hu, X.-H.; Loh, T.-P. *Tetrahedron Lett.* **2010**, 51, 1041.
- (36) Bondalapati, S.; Reddy, U. C.; Kundu, D. S.; Saikia, A. K. *J. Fluorine Chem.* **2010**, 131, 320.
- (37) Alexakis, A.; Jachiet, D. *Tetrahedron* **1989**, 45, 6197.
- (38) Miles, R. B.; Davis, C. E.; Coates, R. M. *J. Org. Chem.* **2006**, 71, 1493.
- (39) Subba Reddy, B. V.; Sreelatha, M.; Kishore, C.; Borkar, P.; Yadav, J. S. *Tetrahedron Lett.* **2012**, 53, 2748.
- (40) Miranda, P. O.; Díaz, D. D.; Padrón, J. I.; Bermejo, J.; Martín, V. S. *Org. Lett.* **2003**, 5, 1979.
- (41) Liu, F.; Loh, T.-P. *Org. Lett.* **2007**, 9, 2063.
- (42) (a) De Meijere, A.; von Zezschwitz, P.; Nüske, H.; Stulgies, B. *J. Organomet. Chem.* **2002**, 653, 129. (b) Malacria, M. *Chem. Rev.* **1996**, 96, 289. (c) Heumann, A.; Réglér, M. *Tetrahedron* **1996**, 52, 9289.
- (43) Kopecky, D. J.; Rychnovsky, S. D. *J. Am. Chem. Soc.* **2001**, 123, 8420.
- (44) Patterson, B.; Marumoto, S.; Rychnovsky, S. D. *Org. Lett.* **2003**, 5, 3163.
- (45) (a) Lin, X.-F.; Cui, S.-L.; Wang, Y.-G. *Tetrahedron Lett.* **2006**, 47, 4509. (b) Royer, L.; De, S. K.; Gibbs, R. A. *Tetrahedron Lett.* **2005**, 46, 4595. (c) Banik, B. K.; Fernandez, M.; Alvarez, C. *Tetrahedron Lett.* **2005**, 46, 2479. (d) Ko, S.; Sastry, M. N. V.; Lin, C.; Yao, C.-F. *Tetrahedron Lett.* **2005**, 46, 5771.
- (46) Yadav, J. S.; Subba Reddy, B. V.; Narayana Kumar, G. G. K. S.; Swamy, T. *Tetrahedron Lett.* **2007**, 48, 2205.
- (47) Yadav, J. S.; Thrimurtulu, N.; Anantha Lakshmi, K.; Prasad, A. R.; Subba Reddy, B. V. *Tetrahedron Lett.* **2010**, 51, 661.
- (48) Sabitha, G.; Bhaskar Reddy, K.; Bhikshapathi, M.; Yadav, J. S. *Tetrahedron Lett.* **2006**, 47, 2807.
- (49) Yadav, J. S.; Subba Reddy, B. V.; Narayana Kumar, G. G. K. S.; Reddy, G. M. *Chem. Lett.* **2007**, 36, 426.
- (50) Yadav, J. S.; Subba Reddy, B. V.; Gupta, M. K.; Biswas, S. K. *Synthesis* **2004**, 2711.
- (51) Norcross, R. D.; Paterson, I. *Chem. Rev.* **1995**, 95, 2041.
- (52) Kang, E. J.; Lee, E. *Chem. Rev.* **2005**, 105, 4348.


- (53) Gesinski, M. R.; Tadpetch, K.; Rychnovsky, S. D. *Org. Lett.* **2009**, 11, 5342.

About the Authors

Sandro José Greco was born in 1974 in Rio de Janeiro, RJ, Brazil. He received his B.Sc. degree in chemistry in 1997 and his M.Sc. and Ph.D. degrees in 2001 and 2005 from the Federal University Fluminense (Rio de Janeiro, Brazil), working under the guidance of Professor Sergio Pinheiro on studies of the use of terpenes and terpenoids in the enantioselective synthesis of potential anticholinergic agents, and on the synthesis of amino alcohol based, new chiral phase-transfer catalysts. In 2006, he joined Professor Maria D. Vargas's group at the Federal University Fluminense as a postdoctoral researcher to work on the synthesis and pharmacological evaluation of new anticancer drugs containing the ferrocenyl group. He is currently an associate professor of organic chemistry at the Federal University of Espírito Santo, with research interests in the design and synthesis of potential bioactive compounds and the development of new organocatalysts and chiral phase-transfer catalysts for asymmetric synthesis.

Rodolfo Goetze Fiorot was born in 1991 in Linhares, ES, Brazil. He is currently pursuing his graduate-level studies in chemistry under the guidance of Professor Sandro J. Greco at the Federal University of Espírito Santo. His research investigations in the Laboratory of Organic & Medicinal Synthesis center on the development of multicomponent organic reactions of naphthoquinones.

Valdemar Lacerda, Jr. was born in 1975 in Goiânia, GO, Brazil. He received a B.Sc. degree in chemistry in 1997 from the Federal University of Goiás, where he worked in the laboratory of Professor Pedro Henrique Ferri. He received his M.Sc. degree in 2000 and his Ph.D. degree in 2004 from São Paulo University (Ribeirão Preto), working with Professor Mauricio Gomes Constantino in organic synthesis and NMR studies. In 2004, he began working as a postdoctoral researcher at the NMR laboratory coordinated by Professor Gil Valdo José da Silva. In 2006, he joined the Department of Chemistry of the Federal University of Espírito Santo (ES State, Brazil) as an associate professor. His current research interests focus on organic synthesis, NMR studies, theoretical calculations, and petroleum studies. He has been Head of the Department of Chemistry since 2007, and is presently also a CNPq level 2 researcher.

Reginaldo Bezerra dos Santos was born in Matão, SP, Brazil, and obtained his B.Sc. degree in chemistry in 1986 from the Federal University of São Carlos (SP State, Brazil). He received his M.Sc. and Ph.D. degrees at the same University in 1990 and 1995, working under the supervision of Prof. U. Brocksmo in the field of organic synthesis. In 1991, he accepted a position at the Federal University of Espírito Santo (ES State, Brazil) as an assistant professor in the Department of Chemistry, and was promoted to Associate Professor in 1995. 

Essential equipment for chemical synthesis

Add Aldrich

Aldrich® Chemical Dryer

Dries Samples Using Heat, Vacuum and Desiccant

- Digital thermocouple control
- Sample viewing port
- Block heater enclosed for safety
- Removable drying tube assembly




Cat. No.	Product Description
Aldrich Chemical Dryer (Order Glass Drying Tube Assembly Separately)	
Z683507	120 V US plug
Z683515	230 V Euro plug
Z683523	230 V UK plug
Glass Drying Tube Assembly	
Z222712	Complete, ST/NS 34/35 joints

For product and ordering information, visit
Aldrich.com/chemicaldryer

Aldrich® Kugelrohr Bulb-to-Bulb Vacuum Distillation Apparatus

**Quickly Distills the Most Difficult Materials
with Minimal Loss**

- Distills up to 230 °C
- Safe, front loading design
- 25 mL to 2 L flask capacity
- Precise digital controls
- High vacuum capability
- Oscillating drive motor



Cat. No.	Product Description
Aldrich Kugelrohr Distillation Apparatus (Order Glassware Separately)	
Z683477	120 V US plug
Z683485	230 V Euro plug
Z683493	230 V UK plug
Kugelrohr Glassware Sets	
Z748293	25 and 100 mL flasks ST/NS 14/20 joints
Z748307	1 L flasks, ST/NS 24/40 joints
Z748315	2 L flasks, ST/NS 24/40 joints

For product and ordering information, visit
Aldrich.com/kugelrohr

Looking for free-flowing salts that are ready when you are?

Add Aldrich™

Introducing Redi-Dri™ salts, available exclusively from Sigma-Aldrich

There are many advantages to using the Redi-Dri product line for hygroscopic powders:

- **Reduce cost:** Minimize the need to break up solids, reducing labor costs and saving valuable lab time
- **Safer to use:** Reduce the risk of injury and contamination by eliminating the need to chip away at solids
- **Longer shelf life:** Redi-Dri remains free flowing, even after repeated exposures to air
- **Less waste:** Eliminates unusable material

See why Dr. Clump switched to Redi-Dri.



Work smarter, not harder. Try Redi-Dri today!

sigma-aldrich.com/redi-dri-video





P.O. Box 14508
St. Louis, MO 63178
USA

Prsrt Std.
U.S. Postage
PAID
St. Louis MO
Permit #750

Change Service Requested

#BYNWDXZ*****AUTO**3-DIGIT 532
#89117089ACTA462# P-10 P4
Alfred Bader
Researcher
924 E Juneau Ave Unit 622
Milwaukee, WI 53202-6827



Search Select Purchase

ALDRICH MARKET SELECT

- Over 8 million products
- 60+ suppliers
- One website

www.AldrichMarketSelect.com



For more information, visit
Aldrich.com/aldrichmarketselect



POU
80733
1085

SIGMA-ALDRICH

VOL. 46, NO. 3 • 2013

ALDRICH
Chemistry

Aldrichimica **ACTA**



Scientist-Led High-Throughput Experimentation (HTE) and
Its Utility in Academia and Industry

Cross-Linked Enzyme Aggregates (CLEAs) in
Organic Synthesis

SIGMA-ALDRICH

Air-stable precatalysts for Suzuki–Miyaura couplings

Add Aldrich

Introducing two new versatile nickel precatalysts (**Prod. Nos. L511706** and **L511714**) for Suzuki–Miyaura cross-couplings. These precatalysts, provided in collaboration with the Hartwig Group, effectively mediate cross-couplings with low catalyst loadings.

They also overcome two long-standing challenges in the area of transition-metal-catalyzed cross-couplings by suppressing protodeboronation and catalyst poisoning by heteroatoms in the substrates.

Other advantages:



- High air and moisture stability
- Versatile substrate compatibility
- No added ligand required



Aldrich has collaborated with Professor John F. Hartwig to bring you these innovative air-stable precatalysts.



For more information:
Aldrich.com/hartwig



Reference:
Ge, S.; Hartwig, J. F. *Angew. Chem. Int. Ed.* 2012, 51, 12837

Aldrich Chemical Co. LLC
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone: 800-325-3010 (USA)
FAX: 800-325-5052 (USA)
or 414-438-2199
Mail: P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries: 800-325-3010
Technical Service: 800-325-5832
SAFC®: 800-244-1171
Custom Synthesis: 800-244-1171
Flavors & Fragrances: 800-227-4563
International: 314-771-5765
24-Hour Emergency: 314-776-6555
Website: sigma-aldrich.com
Email: aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 2060, Milwaukee, WI 53201, USA
sfirsan@sial.com

Subscriptions

Request your FREE subscription to the *Aldrichimica Acta*, through

iPad® App: Aldrich.com/acta-app
Web: Aldrich.com/acta
Email: sales-usa@sial.com
Phone: 800-325-3010 (USA)
Mail: Attn: Marketing Services
Sigma-Aldrich Corporation
3050 Spruce Street
St. Louis, MO 63103



The entire Aldrichimica Acta archive is available at Aldrich.com/acta.

Aldrichimica Acta (ISSN 0002-3100) is a publication of Aldrich, a member of the Sigma-Aldrich Group.

©2013 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA, SAFC, SIGMA-ALDRICH, ALDRICH, and SUPERCO are trademarks of Sigma-Aldrich Co. LLC, registered in the US and other countries. Aldrich is a trademark of Sigma-Aldrich Co. LLC. Aldrich brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the products for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.

iPad is a registered trademark of Apple Inc.




"PLEASE BOTHER US."

Dear Fellow Chemists,

GreenCentre Canada of Kingston, Ontario, Canada, kindly suggested that we offer a technology developed in the lab of Prof. Dennis Hall at the University of Alberta, Canada. Boronic acids have been shown to be effective catalysts in the formation of amide bonds from coupling carboxylic acids with amines, with the only major byproduct being water. 5-Methoxy-2-iodophenylboronic acid (MIBA) is an optimal and recyclable catalyst for the coupling of a wide range of carboxylic acids and amines,

and facilitates amide bond formation at room temperature

Gernigon, N.; Al-Zoubi, R. M.; Hall, D. G. *J. Org. Chem.* **2012**, *77*, 8386



We welcome your product ideas. Do you need a compound that is not featured on our website or in our catalog? Ask Aldrich! For more than 60 years, your research needs and suggestions have shaped the Aldrich product offering. Email us your suggestion at techserv@sial.com

J. Zihlmann

Josef Zihlmann
Vice President of Product Management

TABLE OF CONTENTS

Scientist-Led High-Throughput Experimentation (HTE) and Its Utility in Academia and Industry..... 71
Jason R. Schmink, Ana Bellomo, and Simon Berritt, Bryn Mawr College and the University of Pennsylvania*

Cross-Linked Enzyme Aggregates (CLEAs) in Organic Synthesis 81
Roger A. Sheldon, Sander van Pelt, Seda Karbak-Aksu, Jo-Anne Rasmussen, in Emanuel H. A. Jansen, CLEA Technologies B.V.*

ABOUT OUR COVER

Born and raised in England, Thomas Cole (1801–1848) immigrated with his family to the USA in 1818. Following a series of unremarkable jobs as engraver (England and Philadelphia) and itinerant portrait painter (Ohio), he moved, at 22, to Philadelphia to study at the Pennsylvania Academy of the Fine Arts. While he was largely self-taught, he benefitted from two trips to Western Europe, where he traveled extensively, met prominent artists such as Turner and Constable, and studied the Old Masters, such as Raphael. A founder of the Hudson River school of painting, he was a big influence on American landscape painters of the mid-19th century among whom is his pupil, Frederic Edwin Church. The subjects for many of his canvases were his beloved Hudson River valley, Catskill Mountains (New York), and White Mountains of New Hampshire.



Detail from *A View of the Mountain Pass Called the Notch of the White Mountains (Crawford Notch)*. Courtesy: National Gallery of Art, Washington, DC

A View of the Mountain Pass Called the Notch of the White Mountains (Crawford Notch) (oil on canvas, 102 x 155.8 cm) is a wonderful example of Cole's desire to create a more sophisticated form of landscape painting that expresses a higher meaning.* This realistic and exquisitely detailed landscape, with its warm autumnal yellow and red colors and contrasting brightly lit and dark areas, embodies Cole's romance with America's unspoiled and beautiful wilderness and his dismay at its fast disappearance as a result of encroaching human activity (represented in the foreground by the cut trees, road, structures, and farming). By rendering the human figures small and on the move and the landscape majestic, he is reflecting on man's insignificance and fleeting existence vis-à-vis nature's timeless beauty and power.

This painting was acquired by the National Gallery of Art, Washington, DC, through the Andrew W. Mellon Fund.

* Cole has included elements in this painting that hint at what his message likely is. Can you guess what they are? To find out, visit Aldrich.com/acta463

PRODUCT HIGHLIGHT

Run Hydrogenations on Your Bench!

Q-Tube™ Gas Purging Set and Aldrich® Lecture-Bottle Station

Safe and Simple Systems for Handling Gaseous Reagents

- Hydrogenations
- Catalyst screenings
- Liquid screenings
- Reactions involving gaseous reactants

Would you like to test new reactions and conditions, improve the yield, or minimize byproducts and know the results in 30 minutes?

- Safe Q-Tube Pressure Tube Reactors (Prod. No. Z757993)

For additional Q-Tube products, visit Aldrich.com/qtube.

- Aldrich Lecture-Bottle Station for non-corrosive gases (Prod. No. Z562602)

For Aldrich lecture-bottle gases and additional dispensing stations, visit Aldrich.com/airchem.

For larger-scale reactors up to 1,850 mL, contact Jim Brien at jim.brien@sial.com



Z757993



Z562602

SERVICE HIGHLIGHT

Aldrich Market Select

A Reliable Chemistry Source for Your Library Enrichment Program

- Structural diversity available with over 8 million products provided through our relationship with over 60 global suppliers.
- Custom packaging options are available to meet your research needs.
- Eliminate busy work by letting us facilitate international shipping, compliance, and vendor invoice consolidation.

For more information, visit Aldrich.com/ams-libraries.



and Fine Chemical Industries; Dunn, P. J., Hii, K. K., Krische, M. J., Williams, M. T., Eds.; Wiley: Hoboken, NJ, 2013; pp 75–88.

- (97) Ghislieri, D.; Green, A. P.; Pontini, M.; Willies, S. C.; Rowles, L.; Frank, A.; Grogan, G.; Turner, N. J. *J. Am. Chem. Soc.* **2013**, *135*, 10863.
- (98) Rodriguez-Mata, M.; Frank, A.; Wells, E.; Leipold, F.; Turner, N. J.; Hart, S.; Turkenburg, J. P.; Grogan, G. *ChemBioChem* **2013**, *14*, 1372.
- (99) (a) Gotor-Fernandez, V.; Gotor, V. *Curr. Org. Chem.* **2006**, *10*, 1125. (b) Guranda, D. T.; van Langen, L. M.; van Rantwijk, F.; Sheldon, R. A.; Švedas, V. K. *Tetrahedron: Asymmetry* **2001**, *12*, 1645. (c) Ismail, H.; Madeira Lau, R.; van Langen, M.; van Rantwijk, F.; Švedas, V. K.; Sheldon, R. A. *Green Chem.* **2008**, *10*, 415. (d) Ismail, H.; Madeira Lau, R.; van Rantwijk, F.; Sheldon, R. A. *Adv. Synth. Catal.* **2008**, *350*, 1511. (e) Päiviö, M.; Perkiö, P.; Kanerva, L. T. *Tetrahedron: Asymmetry* **2012**, *23*, 230. (f) Van Rantwijk, F.; Sheldon, R. A. *Tetrahedron* **2004**, *60*, 501. (g) Ditrich, K. *Synthesis* **2008**, 2283.
- (100) (a) Parvulescu, A.; Janssens, J.; Vanderleyden, J.; De Vos, D. *Top. Catal.* **2010**, *53*(13–14), 931. (b) Lee, J. H.; Han, K.; Kim, M.-J.; Park, J. *Eur. J. Org. Chem.* **2010**, 999. (c) Ahn, Y.; Ko, S.-B.; Kim, M.-J.; Park, J. *Coord. Chem. Rev.* **2008**, *252*, 647. (d) Jiang, C.; Cheng, G. *Curr. Org. Chem.* **2013**, *17*, 1225.
- (101) Pohle, S.; Rasmussen, J.; Moody, T.; Brown, G.; Wharry, S.; Sheldon, R. A. CLEA Technologies B.V., Delft, The Netherlands. Unpublished work, 2012.
- (102) Konst, P. CLEA Technologies B.V., Delft, The Netherlands. Unpublished work, 2013.
- (103) López-Iglesias, M.; Busto, E.; Gotor, V.; Gotor-Fernández, V. *Adv. Synth. Catal.* **2011**, *353*, 2345.
- (104) Steunenbergh, P.; Sijm, M.; Zuilhof, H.; Sanders, J. P. M.; Scott, E. L.; Franssen, M. C. R. *J. Org. Chem.* **2013**, *78*, 3802.

Trademarks. Alcalase®, Lipolase®, Lipozyme®, Novozym®, and Ultraflo® (Novozymes A/S); DSM® (DSM IP Assets B.V.); Lyrica® (C.P. Pharmaceuticals International C.V.); Pfizer® (Pfizer Inc.); Plavix® (Sanofi SA); Relenza® (Glaxo Group Limited).

About the Authors


Roger Sheldon is a recognized authority on Green Chemistry and Catalysis and Professor Emeritus of Biocatalysis and Organic Chemistry

at Delft University of Technology. He is widely known for developing the concept of E factors for assessing the environmental footprint of chemical processes. He has a Ph.D. in organic chemistry from the University of Leicester (U.K.). Prior to joining Delft University in 1991, he was Vice President for Research and Development at DSM-Andeno from 1980 to 1990 and with Shell Research Amsterdam from 1969 to 1980. He is currently CEO of CLEA Technologies B.V., a biotech company.

Sander van Pelt studied Chemical Technology and Bioprocess Technology at Delft University of Technology. After obtaining a master's degree in 2005, during which he was an intern at Evonik and Codexis, he joined Prof. Sheldon's research group as a Ph.D. candidate. This research culminated in 2010 in a thesis entitled *The Application of Nitrile Hydratases in Organic Synthesis*. In 2010, he joined CLEA Technologies B.V., where he is currently the Chief Technical Officer.

Seda Kanbak-Aksu studied pharmacy at Hacettepe University, Ankara, Turkey. After her graduation with a B.S. degree as pharmacist, she continued her education in the Pharmaceutical Chemistry group of the same university where she obtained her M.Sc. degree. In 2004, Seda moved to the Netherlands, and, in 2005, she joined the Biocatalysis and Organic Chemistry group at Delft University of Technology, where she completed her Ph.D. studies. In 2009, she joined CLEA Technologies and is currently a Product Manager.

Jo-Anne Rasmussen obtained her Ph.D. degree at the Australian National University in Canberra, Australia, and then took a position at the Commonwealth Scientific and Industrial Research Organisation (CSIRO) in Melbourne, Australia. In 2006, she started working as a scientific researcher at the Technical University of Graz, Austria, and, in 2010, moved to CLEA Technologies in Delft, where she is currently a Senior Scientist.

Michiel Janssen received a master's degree in bioprocess engineering from the then Wageningen Agricultural University (now Wageningen University) and a Ph.D. degree in applied biocatalysis from Delft University of Technology. Before joining CLEA Technologies in 2006, he spent another two years as a postdoctoral researcher at the latter university and a half year at DSM Research (Geleen, The Netherlands). His fields of interest include biocatalytic process development, enzyme immobilization, and nitroxyl radical catalytic oxidations. 

PRODUCT HIGHLIGHT

Looking for Free-Flowing Salts That Are Ready When You Are?

Introducing Redi-Dri™ Salts, Available Exclusively from Sigma-Aldrich

There are many advantages to using the Redi-Dri product line for hygroscopic powders

- **Reduce cost:** Minimize the need to break up solids, reducing labor costs and saving valuable lab time
- **Safer to use:** Reduce the risk of injury and contamination by eliminating the need to chip away at solids
- **Longer shelf life:** Redi-Dri remains free flowing, even after repeated exposure to air
- **Less waste:** Eliminates unusable material

Work smarter, not harder. Try Redi-Dri today! Visit sigma-aldrich.com/redi-dri-video



To access *Aldrichimica Acta* on your iPad®

Add Aldrich™


Automatically Delivered FREE to
Your iPad®

- Timely, topical reviews by leaders in chemical synthesis
- Leading-edge synthetic methodologies, catalysts, and reagents
- Convenient links to web resources and other useful features

Current issues now available



For more information, visit
Aldrich.com/acta-app



VOL. 47, NO. 1 • 2014

ALDRICH
Chemistry

Aldrichimica **ACTA**

DEDICATED TO DR. ALFRED BADER ON THE OCCASION OF HIS 90TH BIRTHDAY



Asymmetric Catalysis Using Chiral, Enantiopure
Disulfonimides

Enabling Novel Photoredox Reactivity via
Photocatalyst Selection

SIGMA-ALDRICH

For the most affordable designer surfactant yet

Add Aldrich

Aldrich collaborated with Professor Bruce Lipshutz to bring you designer surfactants (PTS and TPGS-750-M), which allow for several useful synthetic reactions to be performed in water, rather than in common organic solvents.


We now offer an even more cost-effective designer surfactant, Nok (SPGS-550-M), as an aqueous solution for convenient reaction set up. Nok displays comparable reactivity and recyclability as a third-generation surfactant, complementing the growing portfolio of products that enable sustainable synthesis in water at room temperature.

Room Temperature Reactions with Nok


- Suzuki-Miyaura Cross-Coupling
- Negishi Cross-Coupling
- Olefin Metathesis
- Miyaura Borylation
- Sonogashira Cross-Coupling
- Heck Reaction
- Stille Cross-Coupling
- Amination Reaction



For other innovative Aldrich products from Professor Lipshutz, visit:
Aldrich.com/lipshutz



Comparable Reactivity



surfactant (2 wt % aq.)	conversion (determined by ¹ H NMR)
PTS (698717)	71%
TPGS-750-M (733857)	77%
Nok (776033)	76%

Sustainable Solvent



cycle	conversion (determined by ¹ H NMR)
1	>99%
2	>99%
3	>99%
4	>99%
5	>99%
6	>99%

References:

- (1) Klumpphu, P.; Lipshutz, B. H. *J. Org. Chem.* **2014**, *79*, 888.
- (2) Lipshutz, B. H.; Ghorai, S. *Aldrichimica Acta* **2012**, *45*, 3.

Aldrich Chemical Co. LLC.
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-325-5832
SAFC* 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 314-771-5765
24-Hour Emergency 314-776-6555
Website sigma-aldrich.com
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 2060, Milwaukee, WI 53201, USA
sharbil.firsan@sial.com

Subscriptions

Request your **FREE** subscription to the
Aldrichimica Acta, through:

iPad* App: Aldrich.com/acta-app
Web: Aldrich.com/acta
Email: sams-usa@sial.com
Phone: 800-325-3010 (USA)
Mail: Aldrich Marketing Services
Sigma-Aldrich Corporation
3050 Spruce Street
St. Louis, MO 63103



The entire *Aldrichimica Acta* archive is available at Aldrich.com/acta.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich, a member of the Sigma-Aldrich Group.

©2014 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA, SAFC, SIGMA-ALDRICH, ALDRICH, and SUPELCO are trademarks of Sigma-Aldrich Co. LLC, registered in the US and/or other countries. Add Aldrich is a trademark of Sigma-Aldrich Co. LLC. Aldrich brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the products for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.

iPad is a registered trademark of Apple Inc.



"PLEASE BOTHER US."

Dear Fellow Chemists,

Dr. Jason Mulder, Principal Scientist at Boehringer Ingelheim (Chemical Development, US) suggested that we reconsider the formulation of our LDA solution due to its murky appearance. Not only did we take Jason's suggestion and produce a new formulation in THF-hexanes (**774766**), but also greatly improved the appearance of our long-standing 2.0 M listing (**361798**)! Thanks to his suggestion, we are now proud to offer our customers two high-quality listings of LDA

High-quality LDA: syringe in the background contains the improved, clear 2.0 M solution (**361798**), syringe in the foreground contains the new LDA formulation in THF-hexanes (**774766**).



774766	Lithium diisopropylamide solution, 1.0 M in THF-hexanes	100 mL 4x25 mL
361798	Lithium diisopropylamide solution, 2.0 M in THF-heptane-ethylbenzene	100 mL 4x25 mL 4x100 mL

We welcome your product ideas. Do you need a compound that is not featured on our website or in our catalog? Ask Aldrich! For more than 60 years, your research needs and suggestions have shaped the Aldrich product offering. Email your suggestion to techserv@sial.com


Josef Zihlmann
Vice President of Product Management

TABLE OF CONTENTS

Asymmetric Catalysis Using Chiral, Enantiopure Disulfonimides	3
<i>Manuel van Gemmeren, Frank Lay, and Benjamin List,* Max-Planck-Institut für Kohlenforschung</i>	
Enabling Novel Photoredox Reactivity via Photocatalyst Selection	15
<i>James J. Douglas, John D. Nguyen, Kevin P. Cole, and Corey R. J. Stephenson,* University of Michigan and Lilly Research Laboratories</i>	

ABOUT OUR COVER

Moonlight on the Yare (oil on canvas, 98.4 x 125.7 cm) was painted ca. 1816/1817 by John Crome (1768–1821), and is considered one of his most important works. Crome had a limited formal education, and his artistic training consisted of a seven-year apprenticeship with a house and sign painter and a few lessons in London with Sir William Beechey. Save for a few short visits to London and his friendship with Robert Ladbrooke and John Opie, he had little or no contact with the great artists of his time. Crome supplemented his painting and sketching career by lifelong work teaching drawing in schools and to children of the wealthy.



Detail from *Moonlight on the Yare*. Photo courtesy of National Gallery of Art, Washington

Crome was a founding member and later president of the Norwich School of Painters. While he did not enjoy great or wide fame during his life, his works became much more appreciated after his death, and have earned him a spot among England's prominent landscape painters. As with most of his landscapes, which are noted for their exquisite and realistic rendering of trees and their fidelity to nature, this painting is a romantic depiction of a rustic and tranquil nocturnal scene with a full moon rising into the cloudy sky over the River Yare in the Norwich area where Crome spent all his life. It is reminiscent of landscapes by seventeenth- and eighteenth-century Dutch and English masters,* which Crome greatly admired and copied early in his life as a way to hone his painting skills.

This painting is part of the Paul Mellon Collection at the National Gallery of Art, Washington, DC.

* Which landscapes by seventeenth- and eighteenth-century Dutch and English masters does this painting remind you of? To find out, visit Aldrich.com/acta471

For the newest wave of chiral acid catalysts

Add **Aldrich**

Chiral phosphoric acids such as TRIP have demonstrated the great potential of Brønsted acid catalysts in applications ranging from reductive amination to Friedel–Crafts alkylation.






Recently, a new class of chiral acid has been reported with reactivity that sets it apart from its phosphoric acid predecessors. With a lower pK_a and similarly tunable sterics and electronics, the binaphthyl-derived chiral disulfonimide (DSI) catalyst system shows immediate applicability to the Mukaiyama aldol reaction, Friedel–Crafts alkylation of indoles, and imine allylation.



Brought to you in collaboration with Prof. Dr. Benjamin List



For more information, visit Aldrich.com/dsi




References:

- (1) Garcia-Garcia, P.; Lay, F.; Garcia-Garcia, P.; Rabalakos, C.; List, B. *Angew. Chem. Int. Ed.* **2009**, *48*, 4363
- (2) Chen, L.-Y.; He, H.; Chan, W.-H.; Lee, A. W. M. *J. Org. Chem.* **2011**, *76*, 7141
- (3) Gandhi, S.; List, B. *Angew. Chem. Int. Ed.* **2013**, *52*, 2573

Twin Cities campus, and worked in the laboratory of Prof. Richard Hsung. In 2007, he completed postdoctoral studies in the laboratory of Prof. K. C. Nicolaou at The Scripps Research Institute in San Diego, CA. In the same year, Kevin joined Eli Lilly and Company in Indianapolis, IN, and was promoted to Sr. Research Scientist in 2010. He is currently working as a process chemist in Lilly's continuous processing group.

Corey R. J. Stephenson received an Honours B.Sc. degree in 1998 from the University of Waterloo and a Ph.D. degree in 2005 from the University of Pittsburgh, where he worked under the direction of Prof.

Peter Wipf. After carrying out postdoctoral research in the laboratory of Prof. Erick M. Carreira at ETH Zürich, Switzerland, he joined the Department of Chemistry at Boston University in 2007 as Assistant Professor and co-Principal Investigator in the Center for Chemical Methodology and Library Development (CMLD-BU). He was granted tenure and promoted to Associate Professor in February 2013. As of July 2013, he has been Associate Professor at the University of Michigan, Ann Arbor, where his research group focuses on catalysis, natural product synthesis, and continuous flow chemistry. 

PRODUCT HIGHLIGHT

QianCap™ Safety Cap System with Gas Purging

For Screening Catalysts (15 mL to 1,850 mL)



A QianCap safety cap with automatic pressure release and reseal mechanism to prevent Ace type pressure bottles/tubes from exploding due to excessive internal pressure.

B QianCap with 0–300 ps pressure gauge designed for monitoring internal pressure in Ace type pressure bottles/tubes with internal threads (nos. 15 and 25).

C QianCap with gas purging set
When used with Aldrich® lecture-bottle stations, QianCap can be operated easily, safely, and accurately for reaction reproduction with a wide range of gases, including H₂, D₂, O₂, CO, NH₃, N₂, Ar, CO₂, etc.

For additional QianCap information, visit:
Aldrich.com/qiancap

For higher pressure purging systems, up to 200 psi, visit Q-Tube-Purging-35:
<http://store.qlabtech.com/qipu253.html>


Part Number	Description
A	
Z568724	QianCap safety cap system for no. 15 internal thread
Z568732	QianCap safety cap system for no. 25 internal thread
B	
Z740238	QianCap safety cap system with pressure gauge for no. 15 internal thread
Z740239	QianCap safety cap system with pressure gauge for no. 25 internal thread
C	
Z740240	QianCap gas purging set for no. 15 internal thread
Z740242	QianCap gas purging set for no. 25 internal thread
Z562599	Aldrich lecture bottle station for non-corrosive gases

Celebrating *Alfred Bader's* 90th Birthday

For more than 60 years, Dr. Alfred Bader has been a significant influence in the world of chemistry, business, art and philanthropy. As a tangible manifestation of his vision and dedication to chemistry, he co-founded Aldrich Chemical Company in 1951. Through his contributions, Dr. Bader helped make Sigma-Aldrich the leading Life Science and High Technology company it is today.

His love for cutting-edge research prompted Alfred to start the *Aldrichimica Acta* as an open-access forum for such research. Now in its 47th year of continuous production, the *Acta* has consistently been a top-ranked journal in organic chemistry.

On behalf of Sigma-Aldrich, we thank Dr. Bader for his contributions, as he has made a significant difference in the lives of numerous employees, students, scientists and business leaders around the globe.



Dr. Bader, 2005

*Happy Birthday,
Alfred!*



To wish Dr. Bader a Happy Birthday,
sign our electronic guest book:

Aldrich.com/acta471-bader



VOL. 47, NO. 2 • 2014

ALDRICH
Chemistry

Aldrichimica **ACTA**



Highly Enantioselective Hydroformylation of Alkenes by
Rhodium-Diazaphospholane Catalysts

Boronic Acid Catalysis: an Atom-Economical Platform for Direct
Activation and Functionalization of Carboxylic Acids and Alcohols

SIGMA-ALDRICH

Programmed heterocycle synthesis through aryne precursors

Add **Aldrich**

A variety of substituted heterocycles can be readily synthesized by the trapping of heterocyclic arynes. Aldrich offers indolyne and pyridyne precursors based on chemistry developed in Professor Neil Garg's laboratory.


Advantages of aryne precursors include:

- Arynes are generated in situ under mild conditions
- Reactions are run at ambient temperature
- Reactions are compatible with halides and other cross-coupling partners




Brought to you in collaboration with Professor Neil Garg


For more information, visit
Aldrich.com/garg




L511625





L511668



L511633



L511641



References:

- (1) Im, G.-Y. J.; Bronner, S. M.; Goetz, A. E.; Paton, R. S.; Cheong, P. H.-Y.; Houk, K. N.; Garg, N. K. *J. Am. Chem. Soc.* **2010**, *132*, 17933
- (2) Goetz, A. E.; Garg, N. K. *Nature Chem.* **2013**, *5*, 54.
- (3) Saito, N.; Nakamura, K.; Shibano, S.; Ide, S.; Minami, M.; Sato, Y. *Org. Lett.* **2013**, *15*, 386.

Aldrich Chemical Co. LLC
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FRX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-325-5812
SAFC* 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 314-771-5765
24-Hour Emergency 314-776-6555
Website sigma-aldrich.com
Email aldrich@siatal.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 2060, Milwaukee, WI 53201, USA
sharbil.firsan@siatal.com

Subscriptions

Request your FREE subscription to the
Aldrichimica Acta, through:

IPad® App: Aldrich.com/acta-app
Web: Aldrich.com/acta
Email: sams-esa@siatal.com
Phone: 800-325-3010 (USA)
Mail: Attn: Marketing Services
Sigma-Aldrich Co. LLC
1050 Spruce Street
St. Louis, MO 63103

The entire *Aldrichimica Acta* archive is available at Aldrich.com/acta.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich, a member of the Sigma-Aldrich Group.

©2014 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA, SAFC, SIGMA-ALDRICH, ALDRICH, and SUPELCO are trademarks of Sigma-Aldrich Co. LLC, registered in the US and other countries. Aldrich is a trademark of Sigma-Aldrich Co. LLC. Aldrich brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the products for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.

IPad is a registered trademark of Apple Inc.






"PLEASE BOTHER US."

Dear Fellow Chemists,

Professor E. J. Corey of Harvard University kindly suggested that we consider offering *N*-trifluoroethylanilines, which should prove increasingly useful for the synthesis of highly valuable fluorinated drugs and agrochemicals. We are pleased to offer a number of *N*-trifluoroethylanilines, three of which are featured below.

(1) Muller, K.; Faeh, C.; Diederich, F. *Science* **2007**, 317, 1881. (2) Jin, Z.; Hammond, G. B.; Xu, B. *Aldrichimica Acta* **2012**, 45, 67. (3) Al-Maharik, N.; O'Hagan, D. *Aldrichimica Acta* **2011**, 44, 65. (4) O'Hagan, D. *Chem. Soc. Rev.* **2008**, 37, 308.



L512699	2-Methoxy- <i>N</i> -(2,2,2-trifluoroethyl)aniline	1 g
L512702	4-[<i>N</i> -(2,2,2-Trifluoroethyl)amino]benzonitrile	250 mg
L512710	4-Ethoxycarbonyl- <i>N</i> -(2,2,2-trifluoroethyl)aniline	1 g

We welcome your product ideas. Do you need an *N*-trifluoroethylaniline that is not featured on our website or in our catalog? Ask Aldrich! For more than 60 years, your research needs and suggestions have shaped the Aldrich product offering. Email your suggestion to techserv@siatal.com.


Scott Batcheller
Manager
Chemistry R&D and Product Management

TABLE OF CONTENTS

Highly Enantioselective Hydroformylation of Alkenes by Rhodium-Diazaphospholane Catalysts	29
<i>Gene W. Wong and Clark R. Landis,* University of Wisconsin-Madison</i>	
Boronic Acid Catalysis: an Atom-Economical Platform for Direct Activation and Functionalization of Carboxylic Acids and Alcohols	41
<i>Hongchao Zheng and Dennis G. Hall,* University of Alberta</i>	

ABOUT OUR COVER

Lake Albano (oil on canvas, 121.9 x 170.4 cm) was completed in 1762 by Richard Wilson (1712/1714–1782). Born and raised in Wales, Wilson left for London at about age 15 to study portraiture under Thomas Wright. After about six years of training, he struck out on his own as a portraitist for almost two decades, but without attaining the recognition and financial success he had hoped for. In 1752, he made a trip to Italy that proved to be a turning point in his career as an artist. In Rome, he came to admire and was heavily influenced by the works of celebrated landscape painters, such as Claude Lorrain and Gaspard Dughet, who had lived and worked in the city a century earlier. In 1757 or 1758, he returned to Britain, where he devoted himself to painting primarily Welsh and Italian (often idealized) sceneries and to training a new generation of artists. The fame and success that eluded him as a portrait artist he attained as a landscape painter. Considered by many to be a pioneer of British landscape painting, Wilson has been an acknowledged strong influence on later British landscape artists, in particular John Constable and Joseph Turner.



Detail from **Lake Albano**. Photo courtesy National Gallery of Art, Washington, DC.

Lake Albano* embodies many of the elements that characterize Wilson's landscapes: small human figures; a body of water; a building or two, often in the middle ground of the painting; bright and generally clear skies; and delicate trees and foliage, often more prominent in the foreground. His landscapes, whether actual or idealized, depict charming and serene scenes, with human and animal figures inserted to give depth and scale but not distract from the beauty or majesty of the surroundings. While his landscapes tend to capture the general appearances of nature, his rendering of light and distance validate his keen and delicate observations of the natural world.

This painting is part of the Paul Mellon Collection at the National Gallery of Art, Washington, DC.

* Actual **Lake Albano** was a favorite setting for a number of landscape artists of the 18th and 19th centuries. To find out who else has painted it, visit Aldrich.com/acta472

Highly enantioselective allylic allylation enabled by Landis ligands


Add **Aldrich**

Chiral Ligands for Catalysis


Highly active and highly selective ligands developed by Professor Clark Landis and co-workers are ideally suited for asymmetric allylic allylation reactions, which are essential for accessing a variety of valuable chiral intermediates.



Brought to you in collaboration
with Professor Clark Landis



685859
(*R,R,S*)-DiazaPhos-PPE



685240
(*S,S,S*)-DiazaPhos-PPE

For more information about
our available Landis ligands, visit
Aldrich.com/landisligands

References:

- (1) Clark, T. P.; Landis, C. R. *J. Am. Chem. Soc.* **2003**, *125*, 11792.
- (2) Thomas, P. J.; Axtell, A. T.; Klosin, J.; Peng, W.; Rand, C. L.; Clark, T. P.; Landis, C. R.; Abboud, K. A. *Org. Lett.* **2007**, *9*, 2665.

- (20) Arnold, K.; Batsanov, A. S.; Davies, B.; Whiting, A. *Green Chem.* **2008**, *10*, 124.
- (21) (a) Tang, P. *Org. Synth.* **2005**, *81*, 262. (b) Mylavarapu, R. K.; Kondaiah, G. C. M.; Kolla, N.; Veeramalla, R.; Koilkonda, P.; Bhattacharya, A.; Bandichhor, R. *Org. Process Res. Dev.* **2007**, *11*, 1065. (c) Levonis, S. M.; Bornaghi, L. F.; Houston, T. A. *Aust. J. Chem.* **2007**, *60*, 821.
- (22) Maki, Y.; Ishihara, K.; Yamamoto, H. *Org. Lett.* **2005**, *7*, 5043.
- (23) Latta, R.; Springsteen, G.; Wang, B. *Synthesis* **2001**, 1611.
- (24) Arnold, K.; Davies, B.; Heralut, D.; Whiting, A. *Angew. Chem., Int. Ed.* **2008**, *47*, 2673.
- (25) Wipf, P.; Wang, X. *J. Comb. Chem.* **2002**, *4*, 656.
- (26) Al-Zoubi, R. M.; Marion, O.; Hall, D. G. *Angew. Chem., Int. Ed.* **2008**, *47*, 2876.
- (27) Gemigon, N.; Al-Zoubi, R. M.; Hall, D. G. *J. Org. Chem.* **2012**, *77*, 8386.
- (28) (a) Kalgutkar, A. S.; Marnett, A. B.; Crews, B. C.; Rimmel, R. P.; Marnett, L. J. *J. Med. Chem.* **2000**, *43*, 2860. (b) Moth, C. W.; Prusakiewicz, J. J.; Marnett, L. J.; Lybrand, T. P. *J. Med. Chem.* **2005**, *48*, 3613.
- (29) Yuchun, X.; Huizhou, L.; Jiayong, C. *Int. J. Pharm.* **2000**, *196*, 21.
- (30) Gemigon, N.; Zheng, H.; Hall, D. G. *Tetrahedron Lett.* **2013**, *54*, 4475.
- (31) Marcelli, T. *Angew. Chem., Int. Ed.* **2010**, *49*, 6840.
- (32) Maki, T.; Ishihara, K.; Yamamoto, H. *Org. Lett.* **2005**, *7*, 5047.
- (33) Yamashita, R.; Sakakura, A.; Ishihara, K. *Org. Lett.* **2013**, *15*, 3654.
- (34) (a) Sakakura, A.; Ohkubo, T.; Yamashita, R.; Akakura, M.; Ishihara, K. *Org. Lett.* **2011**, *13*, 892. (b) Sakakura, A.; Yamashita, R.; Ohkubo, T.; Akakura, M.; Ishihara, K. *Aust. J. Chem.* **2011**, *64*, 1458.
- (35) Fotsop, D. F.; Roussi, F.; Leverrier, A.; Bret  ch  , A.; Gu  ritte, F. *J. Org. Chem.* **2010**, *75*, 7412.
- (36) Zheng, H.; Hall, D. G. *Tetrahedron Lett.* **2010**, *51*, 3561.
- (37) Zheng, H.; McDonald, R.; Hall, D. G. *Chem.—Eur. J.* **2010**, *16*, 5454.
- (38) (a) Alvarez, R.; Velazquez, S.; San-Felix, A.; Aquaro, S.; de Clercq, E.; Perno, C.-F.; Karlsson, A.; Balzarini, J.; Camarasa, M. *J. Med. Chem.* **1994**, *37*, 4185. (b) Salameh, B. A.; Leffler, H.; Nilsson, U. J. *Bioorg. Med. Chem. Lett.* **2005**, *15*, 3344.
- (39) Hamada, M.; Roy, V.; McBrayer, T. R.; Whitaker, T.; Urbina-Blanco, C.; Nolan, S. P.; Balzarini, J.; Snoeck, R.; Andrei, G.; Schinazi, R. F.; Agrofoglio, L. A. *Eur. J. Med. Chem.* **2013**, *67*, 398.
- (40) Childs, R. F.; Mulholland, D. L.; Nixon, A. *Can. J. Chem.* **1982**, *60*, 801.
- (41) McCubbin, J. A.; Hosseini, H.; Krokhn, O. V. *J. Org. Chem.* **2010**, *75*, 959.
- (42) McCubbin, J. A.; Krokhn, O. V. *Tetrahedron Lett.* **2010**, *51*, 2447.
- (43) McCubbin, J. A.; Nassar, C.; Krokhn, O. V. *Synthesis* **2011**, 3152.
- (44) Zheng, H.; Hall, D. G. University of Alberta, Edmonton, AB, Canada. Unpublished work, 2012.
- (45) Zheng, H.; Ghanbari, S.; Nakamura, S.; Hall, D. G. *Angew. Chem., Int. Ed.* **2012**, *51*, 6187.
- (46) (a) Tius, M. A. *Eur. J. Org. Chem.* **2005**, 2193. (b) Pellissier, H. *Tetrahedron* **2005**, *61*, 6479. (c) Frontier, A. J.; Collison, C. *Tetrahedron* **2005**, *61*, 7577. (d) Nakanishi, W.; West, F. G. *Curr. Opin. Drug Discovery Dev.* **2009**, *12*, 732.
- (47) For recent examples, see: (a) Hastings, C. J.; Pluth, M. D.; Bergman, R. G.; Raymond, K. N. *J. Am. Chem. Soc.* **2010**, *132*, 6938. (b) Narayan, R.; Fr  hlich, R.; W  rthwein, E.-U. *J. Org. Chem.* **2012**, *77*, 1868.
- (48) Zheng, H.; Lejkowski, M.; Hall, D. G. *Tetrahedron Lett.* **2013**, *54*, 91.
- (49) (a) Cadierno, V.; Crochet, P.; Gimeno, J. *Synlett* **2008**, 1105. (b) Ahlsten, N.; Bartoszewicz, A.; Martin-Matute, B. *Dalton Trans.* **2012**, *41*, 1660. (c) Garc  a-  lvarez, J.; Garc  a-Garrido, S. E.; Crochet, P.; Cadierno, V. *Curr. Top. Catal.* **2012**, *10*, 35.
- (50) (a) Engel, D. A.; Dudley, G. B. *Org. Biomol. Chem.* **2009**, *7*, 4149. (b) Cadierno, V.; Crochet, P.; Garc  a-Garrido, S. E.; Gimeno, J. *Dalton Trans.* **2010**, *39*, 4015.
- (51) Zheng, H.; Lejkowski, M.; Hall, D. G. *Chem. Sci.* **2011**, *2*, 1305.
- (52) (a) Morrill, C.; Grubbs, R. H. *J. Am. Chem. Soc.* **2005**, *127*, 2842. (b) Morrill, C.; Beutner, G. L.; Grubbs, R. H. *J. Org. Chem.* **2006**, *71*, 7813.
- (53) (a) Pennell, M. N.; Unthank, M. G.; Turner, P.; Sheppard, T. D. *J. Org. Chem.* **2011**, *76*, 1479. (b) Pennell, M. N.; Turner, P. G.; Sheppard, T. D. *Chem.—Eur. J.* **2012**, *18*, 4748.
- (54) Debache, A.; Boumound, B.; Amimour, M.; Belfaitah, A.; Rhouati, S.; Carboni, B. *Tetrahedron Lett.* **2006**, *47*, 5697.
- (55) (a) Debache, A.; Boulcina, R.; Belfaitah, A.; Rhouati, S.; Carboni, B. *Synlett* **2008**, 509. (b) Nemouchi, S.; Boulcina, R.; Carboni, B.; Debache, A. *Comptes Rendus Chimie* **2012**, *15*, 394.
- (56) (a) Goswami, S. V.; Thorat, P. B.; Bhusare, S. R. *Heterocycl. Commun.* **2012**, *18*, 245. (b) Goswami, S. V.; Thorat, P. B.; Bhusare, S. R. *Tetrahedron Lett.* **2012**, *53*, 6771.
- (57) Goswami, S. V.; Thorat, P. B.; Chakrawar, A. V.; Bhusare, S. R. *Mol. Divers.* **2013**, *17*, 33.
- (58) Aelvoet, K.; Batsanov, A. S.; Blatch, A. J.; Grosjean, C.; Patrick, L. G. F.; Smethurst, C. A.; Whiting, A. *Angew. Chem., Int. Ed.* **2008**, *47*, 768.
- (59) (a) Arnold, K.; Batsanov, A. S.; Davies, B.; Grosjean, C.; Sch  tz, T.; Whiting, A.; Zawatzky, K. *Chem. Commun.* **2008**, 3879. (b) Georgiou, I.; Whiting, A. *Org. Biomol. Chem.* **2012**, *10*, 2422.
- (60) Li, M.; Yang, T.; Dixon, D. J. *Chem. Commun.* **2010**, *46*, 2191.
- (61) Rao, G.; Philipp, M. *J. Org. Chem.* **1991**, *56*, 1505.
- (62) Goswami, S. V.; Thorat, P. B.; Kadam, V. N.; Khiste, S. A.; Bhusare, S. R. *Chin. Chem. Lett.* **2013**, *24*, 422.
- (63) Goswami, S. V.; Thorat, P. B.; Bhusare, S. R. *J. Chem. Sci.* **2013**, *125*, 745.
- (64) Li, Y.; Molina de La Torre, J. A.; Grabow, K.; Bentrup, U.; Junge, K.; Zhou, S.; Br  ckner, A.; Beller, M. *Angew. Chem., Int. Ed.* **2013**, *52*, 11577.
- (65) Tale, R. H.; Patil, K. M.; Dapurkar, S. E. *Tetrahedron Lett.* **2003**, *44*, 3427.
- (66) Tale, R. H.; Patil, K. M. *Tetrahedron Lett.* **2002**, *43*, 9715.
- (67) Tale, R. H.; Sagar, A. D.; Santan, H. D.; Adude, R. N. *Synlett* **2006**, 415.

Trademarks. Taxol[®] (Bristol-Myers Squibb Company).

About the Authors

Hongchao Zheng was born in Heilongjiang, People's Republic of China. He received his B.S. degree in chemistry in 2001 and his M.Sc. degree in 2005 from Tsinghua University. After immigrating to Canada, Zheng earned a second M.Sc. degree in 2007 from McMaster University (Canada), where he worked under the guidance of Prof. Graham A. McGibbon. In the same year, he joined Professor Dennis Hall's group at the University of Alberta (Canada) to pursue his Ph.D. research in the field of boronic acid catalysis. He obtained his Ph.D. degree in 2012, and is currently conducting research in the field of gold catalysis as a postdoctoral associate in the laboratory of Professor Michel R. Gagn   at the University of North Carolina at Chapel Hill (USA).

Dennis Hall obtained his Ph.D. degree in 1995 from Universit   de Sherbrooke, Canada, under the direction of Prof. Pierre Deslongchamps. He was then an NSERC Postdoctoral Fellow in the laboratory of Prof. Peter G. Schultz in the Department of Chemistry at UC Berkeley. He initiated his independent career in 1997 at the University of Alberta (Edmonton, Canada), where he is presently University Professor of Chemistry. The unifying theme of his research program is the development of new synthetic and biological applications of organoboronic acid derivatives. Hall served as Editor for the successful monograph *Boronic Acids* (Wiley-VCH), now in its 2nd edition. For his contribution of reference 27, he was awarded the *Journal of Organic Chemistry* 2013 Outstanding Author of the Year Award (ACS Organic Division).  

ALDRICH
Chemistry

P.O. Box 14508
St. Louis, MO 63178
USA

SAFETY AND COMPLIANCE, SIMPLIFIED


Subscribe today for instant access to up-to-date
Chemical Safety and Regulatory Data (CSRD)
from Sigma-Aldrich.

A simple upload of Sigma-Aldrich CSRD content to your system will:

- Ensure compliant handling and transfer of chemicals
- Help you easily generate safety and compliance reports
- Enable you to make informed decisions regarding your chemical inventory

For more information, contact us at
ehsservices@sial.com, or visit
sigma-aldrich.com/ehs-services

iPad is a registered trademark of Apple Inc.



SIGMA-ALDRICH

VOL. 48, NO. 1 • 2015

ALDRICH
Chemistry

Aldrichimica ACTA

GREEN CHEMISTRY ISSUE • B. H. LIPSHUTZ, GUEST EDITOR



CONTRIBUTORS TO THIS ISSUE

Anastas • Constable • Gallou and Hamann • Gladysz • Jessop •
Koenig • Krische • Lipshutz • Ritter • Stahl • Tucker • Warner

ADDITIONAL ONLINE CONTRIBUTORS

Rogers • Sheldon

SIGMA-ALDRICH

TABLE OF CONTENTS

Paul T. Anastas: Green Chemistry Next: Moving from Evolutionary to Revolutionary	3
Stephen K. Ritter: Taking a Measure of Green Success	5
David J. C. Constable: Moving Toward a Green Chemistry and Engineering Design Ethic	7
Kelsey C. Miles and Shannon S. Stahl*: Practical Aerobic Alcohol Oxidation with Cu/Nitroxyl and Nitroxyl/NO ₂ Catalyst Systems	8
Fabrice Gallou and Lawrence G. Hamann: Toward a Sustainability Mindset in the Practice of Pharmaceutical Chemistry—from Early Discovery to Manufacturing	12
Inji Shin, T. Patrick Montgomery, and Michael J. Krische*: Catalytic C–C Bond Formation and the Hendricksonian Ideal: Atom- and Redox-Economy, Stereo- and Site-Selectivity	15
John L. Tucker: Pharmaceutical Green Chemistry at Amgen: Seeing with New Eyes	16
Philip G. Jessop: Switchable Solvents as Media for Synthesis and Separations	18
Stefan G. Koenig*, Marian C. Bryan, and Kristi L. Budzinski: The Leading Edge of Green Chemistry at Genentech	22
Bruce H. Lipshutz* and Subir Ghorai: Green Chemistry in the Introductory Organic Laboratory	23
Tathagata Mukherjee and John A. Gladysz*: Fluorous Chemistry Meets Green Chemistry: A Concise Primer	25
John C. Warner: Where We Should Focus Green Chemistry Efforts	29

Additional Online Contributors

Steven P. Kelley and Robin D. Rogers*: A Practical Overview of Organic Synthesis in Ionic Liquids *Aldrich.com/acta*

Roger A. Sheldon: Biocatalysis and Biomass Conversion in Ionic Liquids *Aldrich.com/acta*

ALDRICH
Chemistry

Choose the “greener” alternative

Add **Aldrich**

Aldrichimica Acta—Three “Green” Options

Receive timely, topical reviews by leaders in chemical synthesis and reduce paper waste. All options are complimentary.

- Download to your iPad®
- Receive your subscription via email
- Review issues online

You can also view the previous 47 years of the *Acta* via our online archive.

Switch to a “greener” *Acta*; visit
Aldrich.com/acta

iPad is a registered trademark of Apple Inc.



SIGMA-ALDRICH

Aldrich Chemical Co. LLC
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-325-5832
SAFC 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 314-771-5765
24-Hour Emergency 314-776-6555
Website sigma-aldrich.com
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 2060, Milwaukee, WI 53201, USA
sharbil.firsan@sial.com

Subscriptions

Request your FREE subscription to the
Aldrichimica Acta, through:

iPad® App: Aldrich.com/acta-app
Web: Aldrich.com/acta
Email: sams-usa@sial.com
Phone: 800-325-3010 (USA)
Mail: Attn: Marketing Services
Sigma-Aldrich Co. LLC
3050 Spruce Street
St. Louis, MO 63103

The entire *Aldrichimica Acta* archive is available at Aldrich.com/acta.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich,
a member of the Sigma-Aldrich Group.

©2015 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA,
SAFC, SIGMA-ALDRICH, ALDRICH, and SUPELCO are trademarks
of Sigma-Aldrich Co. LLC, registered in the US and other
countries. Add Aldrich is a trademark of Sigma-Aldrich Co.
LLC. Aldrich brand products are sold through Sigma-Aldrich,
Inc. Purchaser must determine the suitability of the products
for their particular use. Additional terms and conditions may
apply. Please see product information on the Sigma-Aldrich
website at www.sigmaaldrich.com and/or on the reverse side
of the invoice or packing slip.

iPad is a registered trademark of Apple Inc.



"PLEASE BOTHER US."

Dear Fellow Chemists,

Professor Bruce H. Lipshutz and co-workers at UC Santa Barbara have developed three generations of surfactants (PTS, TPGS, and Nok) that enable transition-metal-catalyzed cross-coupling reactions to be run in the absence of organic solvents. To further adoption of this technology, Professor Lipshutz has suggested that we offer individual kits containing all materials needed to carry out such transformations: premeasured catalyst, surfactant, and additives—all in a reaction vessel with a stir bar; just add your substrates! Sigma-Aldrich now offers ten such kits for rapid and convenient olefin metathesis, as well as C–C, C–B, C–CF_n, and C–N bond constructions, in water at room temperature.



For more information about the technology and our full listing in this area, visit Aldrich.com/lipshutzkits

802271	Suzuki–Miyaura coupling in a box, complete with Lipshutz surfactant in water	1 Kit
802298	Buchwald–Hartwig amination in a box, complete with Lipshutz surfactant in water	1 Kit

We welcome your product ideas. Do you need a kit that is not featured on our website or in our catalog? Ask Aldrich! For more than 60 years, your research needs and suggestions have shaped the Aldrich product offering. Email your suggestion to techserv@sial.com

Scott Batcheller
Manager
Chemistry R&D and Product Management

ABOUT OUR COVER

A self-sustaining country estate*—what a fitting subject for the cover of our *Acta* issue that is dedicated to the topic of Green Chemistry and sustainability in chemical synthesis! **Wivenhoe Park, Essex** (oil on canvas, 56.1 × 101.2 cm) was painted in 1816 by John Constable (1776–1837), one of the great British landscapists of the 19th century, who is credited with elevating the status of this genre. Constable's interest in painting started early, but faced initial disapproval from his family; nevertheless, he was able to eventually enroll at London's Royal Academy of Art, where he received his formal training. While he never traveled abroad, he had the opportunity to study the works of earlier Dutch and French landscapists. He struggled in his early years to get recognition from the British art establishment, which he did not receive until he was into his forties. His artistic influence was greatest after his death and outside of England, particularly in France.



Detail from **Wivenhoe Park, Essex**. Photo courtesy National Gallery of Art, Washington, DC.

Constable's delight in, and profound connection to, nature is obvious in his many paintings of the countryside of his native Suffolk and other counties in southeast and southwest England. Nowhere is this more apparent than in this painting with its almost photographic quality. Here, Constable aims to capture on canvas a fleeting moment in the life of this idyllic setting and to convey to the viewer the same sensation of balance and harmony between man and nature that it evoked in him. The precise and crisp brushwork, attention to detail, and the use of billowing, brightly lit clouds to communicate movement, are characteristic of Constable's style.

This painting is part of the Widener Collection at the National Gallery of Art, Washington, DC.

* *Wivenhoe Park* was a working, self-sustaining country estate when Constable painted it. Can you identify some of the elements in the painting that indicate this sustainability? To find out, visit Aldrich.com/acta481

Try your reaction in water—just add your substrates!

Add **Aldrich**

Aldrich Chemistry offers a single-use kit for the facile setup of commonly used cross-couplings and other reactions to be done in water via the assistance of Lipshutz Designer Surfactant, TPGS-750-M, on a 0.5 mmol scale. Written protocols are available on our website.

Kit Components

- Reaction vial with stir bar
- Premeasured additives and precatalysts
- TPGS-750-M in water
- Syringe and septum cap
- Fume hood magnet with standard instructions

Coming Soon! Reaction Kits

- Suzuki–Miyaura Coupling (802271)
- Buchwald–Hartwig Amination (802298)
- Sonogashira Coupling (802301)
- Heck Coupling (802328)
- Stille Coupling (802336)
- Negishi Coupling (802344)
- Miyaura Borylation (802352)
- Olefin Metathesis (802360)
- Trifluoromethylation (802379)
- Amide Coupling (802387)



For more information or to order:
Aldrich.com/lipshutzkits



Where We Should Focus Green Chemistry Efforts

Insight from the Co-Author of the 12 Principles of Green Chemistry

John C. Warner

Warner Babcock Institute for Green Chemistry, 100 Research Drive, Wilmington, MA 01887, USA • Email: John.Warner@warnerbabcock.com

I am often asked to identify what I feel are the most important technologies that need to be improved from a Green Chemistry perspective. I am sure that I am not alone in receiving this type of question. One can understand why this question is so prevalent. It is certainly no secret that today's funding resources for scientific research are scarce. It makes sense to want to make progress in a fewer number of focused areas, rather than be spread too thin making imperceptible headway across too many pathways. It is logical to want to identify the technologies that have the most negative impact on human health and the environment, and cross-reference them with the most promising early-stage ideas. In so doing, the assumption follows that these early successes will provide confidence and momentum to expand and take on more difficult challenges.

My reaction to this question is to become quite uncomfortable: who am I to be so presumptuous as to pick what is most important? In choosing the focus areas, should we identify classes of materials: phthalates, brominated flame retardants, parabens, or plastics? Perhaps we should identify negative impacts: endocrine disruption, global climate change, carcinogenesis, or fossil fuel consumption? Maybe we should focus on areas of research: catalysis, solvents, nanotechnology, or feedstock? It seems like we need help in every area. But there isn't really any "one size fits all" winner. I would imagine different industry sectors have different top priority needs. What the petroleum industry needs is probably different from what the pharmaceutical industry needs; and people working on the front lines of patient care at a cancer hospital will likely see different needs than farmers working in the fields of a vegetable farm.

And then there is the question of how innovation should be "caused to happen" in any given area once it is selected. I often say that Green Chemistry has three long-term requirements: (i) It must have less impact on human health and the environment than some incumbent technology. But this is not enough. (ii) It also must perform better than the incumbent technology in order to be successful in replacing it. No one is going to use a cleaner that doesn't clean, just because it is greener. (iii) It also must have appropriate costs. Society has not demonstrated an overwhelming willingness to pay too much of a premium for a green technology. Reflecting upon this, Green Chemistry really provides a "holy grail" of sorts: better performance, better cost, and "Oh, by the way... it is better for human health and the environment." Given these criteria, what market barrier exists for adopting a greener product save its invention in the first place?


So the question then becomes, how does identifying a focus area that "needs improvement" catalyze innovation? How do scientists learn how to "make better products"? The United States Environmental Protection Agency has given out 5 Presidential Green Chemistry Challenge awards every year since 1996. Other international organizations have various award programs as well. Exemplifying and disseminating case studies of successful discoveries and implementations of Green Chemistry is critical. But is it enough? Does listening to a symphony teach someone to play a musical instrument? Does watching the Boston Marathon train someone how to run 26.2 miles? While illustrative examples are important, something deeper is necessary.

We really need to get at the roots of how we train chemists in the first place. We need to reflect upon the words and semantics we use in our scientific "disciplines". We have a body of knowledge called "organic chemistry". We have another body of knowledge called "physical

chemistry". There is "analytical chemistry", "inorganic chemistry", and several other subdisciplines. But, it is important to realize that these subdisciplines are merely intellectual constructs. They allow us to apply the reductionist approach while learning and growing the body of knowledge. We place boundaries on a given subdiscipline for bureaucratic reasons: in order to define the scope of a journal, the syllabus of a class, or the table of contents of a textbook. But, at the end of the day, there really is no such thing as an "organic chemist". One cannot possibly function in the world of chemistry doing "only" organic chemistry. Organic chemists function simultaneously as analytical, physical, theoretical, and other types of chemist. They choose to identify organic chemistry as their "specialty"; but, in no way could they restrict themselves to doing "only" organic chemistry. All of the other subdisciplines work in the same way. It is an illusion to believe that any of these subdisciplines have any independent reality of their own. In my opinion, Green Chemistry should work the same way; and the field is growing as such. There now are many journals, classes, and textbooks in Green Chemistry. The only problem is that Green Chemistry is still seen as something of an "elective" in chemistry. Can you imagine having a student graduate with a degree in chemistry after only being shown "illustrative examples" of physical chemistry? It is inconceivable that a chemistry degree program at a university suggest that, instead of having several semester-long classes in organic chemistry, the initiated should look at examples of papers in the *Journal of Organic Chemistry*. While these hypothetical examples are laughable, we are still at a point in our evolution where Green Chemistry is treated as such. This must change.

Of course, there are moral and ethical issues surrounding this. It is an inescapable truth that the field of chemistry has a certain obligation to society to make sure that future practitioners learn some fundamental principles regarding the making of materials and products that have reduced impact on human health and the environment. But, this is also about innovation and economic competitiveness. Chemists need to treat Green Chemistry as simply a part of the fundamentals of chemistry—nothing more, and nothing less.

This is where the Green Chemistry Commitment program of the nonprofit organization Beyond Benign (www.beyondbenign.org) comes into play. Obviously, it is difficult to introduce Green Chemistry into the basic chemistry curriculum. If no one has had these classes before, who is going to teach them? Beyond Benign seeks to create a community of chemistry departments to share best practices for what uniquely works for them. If a chemistry department wants to develop a standalone class in Green Chemistry, it should do so. If a department wants to integrate it across various existing classes, it should do so. If a department wants to address Green Chemistry in the labs, it should do so. By sharing best practices, the community can grow and move from illustrative examples to basic pedagogy.

So when I am asked to identify what I feel are the most important technologies that need to be improved from a Green Chemistry perspective, my answer is instantaneous: **education**. Imagine a world where all chemists had fundamental training in how to design technologies that had reduced impact in human health and the environment. Instead of merely picking the technologies, materials, and endpoints that we need to improve upon, let's focus on how we train the scientists in the first place; and we can consequently improve *all*. 

ALDRICH[®]
Chemistry


P.O. Box 14508
St. Louis, MO 63178
USA

GREEN CHEMISTRY ISSUE • B. H. LIPSHUTZ, GUEST EDITOR

Introducing Our New **Greener Alternatives Platform**

Using Greener Alternatives Doesn't Mean Making Compromises

Our new platform gives you more choices to reduce the environmental impact of your research, while still delivering the quality and efficacy you need.

- Search for "Greener Alternatives" on our website to view all compliant products
- Look for the Greener Alternatives product identifier  on our website
- Use DOZN™, our industry-first quantitative green chemistry analysis tool, to evaluate products
- View documentation supporting product greenness



Discover how you can minimize your environmental footprint at
sigma-aldrich.com/green



SIGMA-ALDRICH

VOL. 48, NO. 2 • 2015

ALDRICH
Chemistry

Aldrichimica **ACTA**



Catalytic Asymmetric Hydrogenation of α -Substituted Ketones and Aldehydes
via Dynamic Kinetic Resolution: Efficient Approach to Chiral Alcohols

SnAP Reagents for a Cross-Coupling Approach to the One-Step Synthesis of
Saturated N-Heterocycles

Development of Solvent Selection Guides

SIGMA-ALDRICH

Diversify your lead

Add Aldrich™

Late Stage Functionalization (LSF) allows synthetic and medicinal chemists to structurally modify one Lead Compound into a small library of analogs for further screening and eventual optimization. Diversification achieved directly from a compound *already in hand* obviates the need for ground-up syntheses of each analog; a task often hampered by time, cost and the availability of appropriate starting materials.

The Aldrich Late Stage Functionalization Toolkit (ALD00444) provides ample quantities of our most popular and effective LSF reagents, presenting an opportunity to diversify one heterocycle-containing lead compound into a minimum of twelve value-added analogs. In collaboration with Professor Phil Baran and his research group, innovative LSF is accomplished with the Baran Diversinates™ and Palau'Chlor® found in this toolkit.

ALDRICH LATE STAGE FUNCTIONALIZATION TOOLKIT (ALD00444)

Late Stage Functionalization

Late Stage Functionalization (LSF) allows synthetic and medicinal chemists to structurally modify one lead compound into a small library of analogs for further screening and eventual optimization. Diversification achieved directly from a compound *already in hand* obviates the need for ground-up syntheses of each analog; a task often hampered by time, cost and the availability of appropriate starting materials.

The Aldrich Late Stage Functionalization Toolkit provides ample quantities of our most popular and effective LSF reagents to date. Armed with this Toolkit, you can diversify one lead compound into a minimum of twelve value-added analogs.

Intelligent Diversification

The components of the LSF Toolkit



For a full listing of available Diversinates, visit:

Aldrich.com/diversinates

Aldrich Chemical Co. LLC.
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)
or 414-438-2199
Mail P.O. Box 2060
Milwaukee, WI 53201, USA

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-325-5832
SAFC® 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 314-771-5765
24-Hour Emergency 314-776-6555
Website sigma-aldrich.com
Email aldrich@sial.com

General Correspondence

Editor: Sharbil J. Firsan, Ph.D.
P.O. Box 2060, Milwaukee, WI 53201, USA
sharbil.firsan@sial.com

Subscriptions

Request your FREE subscription to the
Aldrichimica Acta, through:

iPad® App: Aldrich.com/acta-app
Web: Aldrich.com/acta
Email: sams-usa@sial.com
Phone: 800-325-3010 (USA)
Mail: Attn: Marketing Services
Sigma-Aldrich Corporation
3050 Spruce Street
St. Louis, MO 63103

The entire *Aldrichimica Acta* archive is available at Aldrich.com/acta.

Aldrichimica Acta (ISSN 0002-5100) is a publication of Aldrich, a member of the Sigma-Aldrich Group.

©2015 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA, SAFC, SIGMA-ALDRICH, ALDRICH, and SUPELCO are trademarks of Sigma-Aldrich Co. LLC, registered in the US and other countries. Add Aldrich is a trademark of Sigma-Aldrich Co. LLC. Aldrich brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the products for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.

iPad is a registered trademark of Apple Inc.



"PLEASE BOTHER US."

Dear Fellow Chemists,

Amos B. Smith, III, Rhodes-Thompson Professor of Chemistry at the University of Pennsylvania, kindly suggested that we consider offering 3,6-dichloro-1,2,4,5-tetrazine, which has been employed to staple cysteine residues of peptide chains, and can either be released upon photolysis or used in inverse-electron-demand Diels–Alder reactions to incorporate various chemical probes.

(1) Brown, S. P.; Smith, A. B., III. *J. Am. Chem. Soc.* **2015**, *137*, 4034. (2) Courter, J. R.; Abdo, M.; Brown, S. P.; Tucker, M. J.; Hochstrasser, R. M.; Smith, A. B., III. *J. Org. Chem.* **2014**, *79*, 759. (3) Tucker, M. J.; Abdo, M.; Courter, J. R.; Chen, J.; Brown, S. P.; Smith, A. B., III; Hochstrasser, R. M. *Proc. Natl. Acad. Sci. U.S.A.* **2013**, *110*, 17314.



792594 3,6-Dichloro-1,2,4,5-tetrazine 1 mg
10 g

We welcome your product ideas. Do you need a product that is not featured on our website or in our catalog? Ask Aldrich! For more than 60 years, your research needs and suggestions have shaped the Aldrich product offering. Email your suggestion to techserv@sial.com.

Scott Batcheller
Manager
Chemistry R&D and Product Management

TABLE OF CONTENTS

Catalytic Asymmetric Hydrogenation of α -Substituted Ketones and Aldehydes via Dynamic Kinetic Resolution: Efficient Approach to Chiral Alcohols..... 33
Jian-Hua Xie and Qi-Lin Zhou, Nankai University

SnAP Reagents for a Cross-Coupling Approach to the One-Step Synthesis of Saturated N-Heterocycles 43
Michael U. Luescher, Kimberly Geoghegan, Paula L. Nichols, and Jeffrey W. Bode, ETH-Zürich*

Development of Solvent Selection Guides 51
Harry E. Eastman, Craig Jamieson, and Allan J. B. Watson, University of Strathclyde*

ABOUT OUR COVER

Johann Georg von Dillis (1759–1841), a German master painter of the late 18th and early 19th centuries, painted and signed *A Royal Party Admiring the Sunset atop the Hesselberg Mountain* in 1801. He received his first drawing lessons while attending the Gymnasium in Munich. He then studied art at the Munich Zeichnungsschule under the guidance of F. I. Oefele and J. J. Dörner the Elder. His art career began ca. 1786 and, until his retirement in the late 1830s, consisted of official appointments by the courts of Maximilian I and Ludwig I, commissioned sketches and drawings, and a professorship of landscape painting at the Munich Royal Academy of Fine Arts. He travelled frequently and widely throughout Europe on Bavarian State business, and it was during these travels that he became acquainted with, and influenced by, the work of P.-H. de Valenciennes, S. Denis, J.-X. Bidauld, W. Allston, and, especially, J. M. W. Turner. His artwork influenced such artists as C. E. F. Blechen, and his gallery work had a profound influence on the world of art in Germany, especially in Bavaria.



Detail from *A Royal Party Admiring the Sunset atop the Hesselberg Mountain*. Photo courtesy National Gallery of Art, Washington, DC.

Von Dillis excelled at landscapes, and was instrumental in moving the genre from the classical tradition of idealized pastorals to a new, realistic form, whereby the artist draws from life, not only what he perceives with his eyes upon close examination of nature, but also what he feels within himself. This romantic approach to landscape painting values emotions, which are aroused nowhere better than in the countryside with its powerful pull on the artist, as this work demonstrates. Not only is the average person drawn to nature, but even sophisticated, urban dwellers such as this elegantly attired group of royals* and their aides can appreciate the wonders of nature. This composition (watercolor, gouache, and pen and gray ink over graphite on laid paper) measures 37 x 42.7 cm and is von Dillis's best-known work. With its crisp colors, decidedly finished look, and the artist's meticulous attention to detail, it clearly highlights von Dillis's complete immersion in the work.

This painting is part of the Wolfgang Ratjen Collection at the National Gallery of Art, Washington, DC.


*Who could the "royals" depicted in this painting be? To help solve the mystery, visit Aldrich.com/acta482

Rediscover Aldrich catalysis with the repriced Buchwald portfolio


Add **Aldrich**TM

Easily activated Palladium Precatalysts for facile C-N and C-C bond formations.


Buchwald palladium precatalysts are highly reactive, air-stable, and crystalline white solids that generate the active monoligated Pd(0) species in situ, providing a more efficient and easily accessible catalytic system.




G1 = 1st Generation Precatalysts



G2 = 2nd Generation Precatalysts




G3 = 3rd Generation Precatalysts




To see the complete list of available Buchwald Precatalysts and Ligands, visit Aldrich.com/Buchwald

Now the most competitively priced Buchwald offering!


Experiment with our Buchwald portfolio of products, which we are now offering at up to 69% off of our former list prices without compromising the high quality that you have come to expect from Aldrich. Examples of the new competitive price reductions are shown below.




730998-500MG
tBuBrettPhos
New Price—reduced 69%




638072-5G
SPhos
New Price—reduced 38%




718742-100MG
BrettPhos
New Price—reduced 37%




753009-5G
SPhos Pd G2
New Price—reduced 34%




761605-5G
BrettPhos Pd G3
New Price—reduced 35%



763403-250MG
RuPhos Pd G3
New Price—reduced 21%



638080-5G
tBuXPhos
New Price—reduced 41%



708739-1G
tBuXPhos Pd G1
New Price—reduced 35%

GCI Pharmaceutical Roundtable: Collaboration to Deliver a Solvent Selection Guide for the Pharmaceutical Industry, 2008. ACS Green Chemistry Institute® Web site. <http://www.acs.org/content/dam/acsorg/greenchemistry/industriainnovation/roundtable/solvent-selection-guide.pdf> (accessed April 13, 2015).

- (9) MacMillan, D. S.; Murray, J.; Sneddon, H. F.; Jamieson, C.; Watson, A. J. B. *Green Chem.* **2012**, *14*, 3016.
- (10) (a) Taygerly, J. P.; Miller, L. M.; Yee, A.; Peterson, E. A. *Green Chem.* **2012**, *14*, 3020. (b) Peterson, E. A.; Dillon, B.; Raheem, I.; Richardson, P.; Richter, D.; Schmidt, R.; Sneddon, H. F. *Green Chem.* **2014**, *16*, 4060.
- (11) MacMillan, D. S.; Murray, J.; Sneddon, H. F.; Jamieson, C.; Watson, A. J. B. *Green Chem.* **2013**, *15*, 596.
- (12) McGonagle, F. I.; MacMillan, D. S.; Murray, J.; Sneddon, H. F.; Jamieson, C.; Watson, A. J. B. *Green Chem.* **2013**, *15*, 1159.
- (13) Skowerski, K.; Biatacki, J.; Tracz, A.; Olszewski, T. K. *Green Chem.* **2014**, *16*, 1125.
- (14) For a recent review of LCA, see Jiménez-González, C.; Overcash, M. R. *Green Chem.* **2014**, *16*, 3392.
- (15) Neher, R. In *Thin Layer Chromatography*; Marini-Bettolo, G. B., Ed.; Elsevier: Amsterdam, The Netherlands, 1964; pp 75–86.
- (16) Roughley, S. D.; Jordan, A. M. *J. Med. Chem.* **2011**, *54*, 3451.
- (17) Cooper, T. W. J.; Campbell, I. B.; Macdonald, S. J. F. *Angew. Chem., Int. Ed.* **2010**, *49*, 8082.
- (18) Vidal, C.; García-Álvarez, J.; Hernán-Gómez, A.; Kennedy, A. R.; Hevia, E. *Angew. Chem., Int. Ed.* **2014**, *53*, 5969.
- (19) (a) Adams, J. P.; Alder, C. M.; Andrews, I.; Bullion, A. M.; Campbell-Crawford, M.; Darcy, M. G.; Hayler, J. D.; Henderson, R. K.; Oare, C. A.; Pendrak, I.; Redman, A. M.; Shuster, L. E.; Sneddon, H. F.; Walker, M. D. *Green Chem.* **2013**, *15*, 1542. (b) Henderson, R. K.; Hill, A. P.; Redman, A. M.; Sneddon, H. F. *Green Chem.* **2015**, *17*, 945.

Trademarks. ACS Green Chemistry Institute® (American Chemical Society); Isopar™ (Exxon Mobil Corp.).

About the Authors

Harry Eastman was born in 1992 in Portsmouth, England, U.K. He studied chemistry at the University of Surrey, working for a year with Professor A. M. P. Koskinen at Aalto University, and completed an M.Sci. degree in Chemistry program in 2014. He began his Ph.D. studies with Dr. Allan Watson in 2014, and is currently focused on the development of new methods for organic synthesis.

Craig Jamieson earned his B.Sc. (Hons) in Chemistry degree at the University of Glasgow in 1996. His Ph.D. studies were carried out under the direction of Professor R. Ramage at the University of Edinburgh (1999). Following postdoctoral research under the supervision of Professor S. V. Ley at the University of Cambridge, he was appointed in 2001 as Principal Scientist in GlaxoSmithKline's Discovery Medicinal Chemistry group, where he worked on a range of exploratory medicinal chemistry programs. In 2004, Craig joined Organon Laboratories (later Merck Research Labs) as a Group Leader in the Medicinal Chemistry Department, with responsibility for hit-to-clinical candidate optimization. In August 2010, he accepted an appointment as John Anderson Research Lecturer in Chemical Biology in the Department of Pure and Applied Chemistry at the University of Strathclyde.

Allan J. B. Watson obtained his M.Sci. degree in 2004 from the University of Strathclyde and his Ph.D. degree (with Professor W. J. Kerr) in 2008. This was followed by postdoctoral work with Professor D. W. C. MacMillan at Princeton University, and an industrial postdoctoral research position at GlaxoSmithKline. He returned to a lecturer position at the University of Strathclyde in 2011. His research interests are synthetic methodology, boron chemistry, medicinal chemistry, and sustainable synthesis. ☉


PRODUCT HIGHLIGHT

Aldrich® Micro Photochemical Reactor

This newly designed apparatus enables researchers to process up to 16 samples at a time. Safe LED light is an energy-efficient source for performing photoreactions.

Features

- Aluminum sample base with riser feet; holds up to 16 samples and sits atop a hot plate stirrer (not included)
- LED light ring is IP68 double-density 12 vdc water-proof blue
- Uses 4 mL sample vials with either PTFE-lined solid-top caps or holed caps with PTFE face rubber septa (purchase separately)
- Power supply is wall plug power supply; 500 mA 5–6 watts, CE compliant




Aldrich Prod. No. ALDKIT001

For more information, visit sigma-aldrich.com

ALDRICH
Chemistry

P.O. Box 14508
St. Louis, MO 63178
USA

SEARCH FOR ANSWERS, NOT JUST STRUCTURES



INTRODUCING BUILDING BLOCK EXPLORER

A comprehensive collection of building blocks
backed by scientific expertise, accurate inventory
and reliable shipping — all in one place.

START EXPLORING | sigma-aldrich.com/explore

©2015 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA, SIGMA-ALDRICH, SUPELCO, and FLUKA are trademarks of Sigma-Aldrich Co. LLC or its affiliates, registered in the US and other countries.

SIGMA | **ALDRICH** | **Fluka** | **SUPELCO** | **SAFC**



83599-h
RPU
83599-512540
1085

SIGMA-ALDRICH